

In-Gel Trypsin Digestion

Reagents: (prepare freshly right before the digestion)

(Use HPLC grade solvents, highest possible grade reagents and MilliQ water for all preparations)

Wash solution: 50% methanol and 5% acetic acid.

25 mM NH_4HCO_3 in 50% acetonitrile (ACN)

25 mM NH_4HCO_3

5% formic acid in water

5% formic acid in 50% ACN

12.5 ng/ μl trypsin in 25 mM NH_4HCO_3 (Sequencing Grade Modified Trypsin, Promega Cat# V5111)

1.5-ml microcentrifuge tubes: Low protein binding, siliconized tubes are recommended. Rinse tubes in USP-grade methanol and dry before the use.

Procedure:

1. Excise protein spot/band (from SYPRO Ruby, Coomassie, or MS compatible silver stained gels), cut gel slice into small pieces ($\sim 1 \text{ mm}^3$) and place them in a 1.5-ml microcentrifuge tube (siliconized, methanol washed).
2. Add 200 μl Wash solution to the gel pieces for 4 hr to overnight. Use gel loading pipet tips to remove the supernatant and discard. Repeat once for 2-3 hr.
3. Add 200 μl of 25 mM NH_4HCO_3 in 50% ACN and vortex 10 min. Repeat this step once or twice.
4. Dry the gel pieces in a Speed Vac. (~ 30 min)
5. Reduction and Alkylation (Optional, although this has been done before 2nd dimension PAGE, for better digestion and recovery this step is recommended).
 - a.) Add 40 μl (or enough to cover) of 10 mM DTT in 25 mM NH_4HCO_3 and incubate at 56°C for 1 h. Remove and discard supernatant. Let the tube cool to room temp.
 - b.) Add the same volume of 55 mM iodoacetamide in 25 mM NH_4HCO_3 . Incubate at room temp for 45 min in the dark with occasional vortex. Remove and discard supernatant.
 - c.) Wash with ~ 100 μl of 25 mM NH_4HCO_3 and vortex 10 min. Remove and discard supernatant.
 - d.) Dehydrate gel pieces with ~ 100 μl of 25 mM NH_4HCO_3 in 50% ACN, vortex 5 min. Remove and discard supernatant.
 - e.) Repeat steps c. and d. for one more time.
 - f.) Dry the gel in a Speed Vac.

6. Add 25 μl (or enough for gel to absorb) of 12.5 ng/ μl Trypsin in 25 mM NH_4HCO_3 to the dry gel pieces, wait for ~ 15 min until the trypsin solution has been absorbed by the gel pieces, remove excess solution and discard. Add a minimum amount of 25 mM NH_4HCO_3 without trypsin to keep gel pieces immersed throughout digestion. Incubate at 37°C for 12-16 h.
7. To recover the peptides from the gel.
 - a.) Add 30 μl of 5% formic acid in 50% ACN, agitate at rt for 30-60 min. Remove the solution into a clean siliconized tube.
 - b.) Add 15 μl of 5% formic acid, incubate at rt for 10 min, and then add 15 μl of 100% ACN, agitate for 30-60 min.
 - c.) Combine two extract solutions and dry in a Speed Vac.
8. Dissolve the dry pellet in ~ 10 -20 μl of 0.1% formic acid, and store the peptide solution at -20°C .

Reference:

1. Shevchenko, A., Wilm, M., Vorm, O., Mann, M. 1996. Mass spectrometric sequencing of proteins in silver-stained polyacrylamide gels. *Anal. Chem.* **68**, 850-858
2. Jiménez, C. R., Huang, L., Qiu, Y., and Burlingame, A. L. 2003. In-gel digestion of proteins for MALDI-MS fingerprint mapping. in *Current Protocols in Protein Science*, John E. Coligan, Ben M. Dunn, David W. Speicher, and Paul T. Wingfield (eds.), John Wiley & Sons, Inc. Unit 16.4