



# OLYMPUS FV3000 基本操作手冊




- 本機
- 資源回收筒
- FV  
FV31S-DT
- FV  
FV31S-SW
- CS  
cellSens  
Dimensi...
- Temp
- FV3000  
Manual



點選畫面上的  
"FV31S-SW"  
啟動軟體

**OLYMPUS**



請記得勾選核取方塊  
以操作電動載物臺

按下OK,進行下一步


FV31S-SW  
FLUOVIEW

Initializing...




Copyright 2013-2021 OLYMPUS CORPORATION All rights reserved.


OLYMPUS



如搭配IX3-SSU電動載物臺時，請點選”No”，載物臺清潔請由工程師處理



按下"Layout"，叫出介面配置選單，選取"OLYMPUS standard" 按下"Apply"套用



掃描器控制

偵測器控制

Z stack & Time lapse 設定

存檔路徑  
檔案名稱  
執行掃描

即時影像/已掃描影像

OLYMPUS FV31S-SW

File Tool Window Tools Window Help

Layout Acquisition Viewer cellSens

LSM Imaging x Ocular

Scan Settings

Scanner:

Type:  Galvano  Resonant

Mode:  OneWay  Roundtrip

Interlace:  OFF  2x  4x

Speed: 2.0 us/pixel

Image Size:

Aspect Ratio:  1:1  4:3

Scan Size: 512x512

High voltage: 2000 V

Pixel: 2.0 us/line 2000 ms/frame 2.17 sec

Area Settings

Rotation: 0.0 deg

Pan X: 0.00 um

Pan Y: 0.00 um

Zoom: 1.00

PMT Setting

PMT

Mode:  VBF  Lambda

Average:  None  Line  Frame 2 Times

Sequential Scan:  None  Line  Frame

Dye & Detector Select

Phase1

Confocal Aperture: Auto 117 um

Laser ND Filter:  None  10%

P1 CH1 Alexa Fluor 405 HSD1 430 -500 nm G-1

405 0.9 %

HV 700 V

Gain 1.000 x

Live Graph Table Analysis

Analysis ROI:

Live Livex2 Livex4 1:1 100 % LUT: Hi-Lo Auto Contrast Reset

Title Single

Series

LSM

Time:  ON  OFF

Z:  ON  OFF

Series: 2 min 50 sec

Time Lapse

LSM Total: 0:02:50.84

Interval: FreeRun

Scan: 0:02:50.84

Cycle: 5

Z Section

Motor:  Start/End  Range

Current: 3746.61

Origin: 200.00

Register

Move

Slices: 76

Step Size: 3.99

Start: 3746.61

End: 4045.86

Acquisition

Normal Sync Sequence MATL

Imaging

LSM Start

Stimulation Start

SERIES DONE Append

APPEND NEXT

LSM D:\CFC

Biofilm

current scan condition


Total scanning time: 0:02:50.84


remaining time: -

Next scan start at: -

Windows taskbar: 上午 10:02 2021/7/20

# XY單張取像 Step by step





先選取Ocular視窗，點選”Ocular”  
使用螢光(EPI)或是穿透光(DIA)找尋  
目標樣本，找到目標後點選”LSM”  
將光路切換回掃描器



Mode :  VBF

Average:  None

Sequential Scan:  None

**Dye & Detector Select**

Phase 1

Confocal Aperture  Auto 117

Airy Disk x1.00

HSD1 HSD2

Group 1 CH1

Group 2 CH2

Group 3

Laser ND Filter:  None

P1 CH1 Alexa Fluor

405 0.9 %

HV 700 V

Gain 1.000 x

Scan: 0:02:50.84 Rest: 0:00:00.00

Cycle: 5

Z Section

Motor:  Start/End  Range

Current: 2499.38 Origin: 200

Slices: Step Size:

Start: 3746.61 End:

FV Dye and Detector Select

Dye Channel

Selected Dye Spectrum Assigned Dye Spectrum

400 500 600 700 800

Dye List

Recent Dyes :

Dye	Excitation	Emission
Alexa Fluor 568	579	603
DAPI	359	461
Alexa Fluor 405	401	422
EGFP	488	510
Yellow	561	586
Cy2	484	500

Registered Dyes :

Dye	Excitation	Emission
Violet	405	430
Cyan	445	470
Blue	488	513
Emerald	514	539
Yellow	561	586
Orange	604	610

Assigned Dye Add Phase

Phase 1

Channel	Dye	Detector
CH1	Alexa Fluor 405	HSD1
CH2	EGFP	HSD2
CH3	Yellow	SD3

Add Remove All Clear



選取偵測器控制視窗中的  
 "Dye & Detector Select"，選取樣本  
 中所使用的螢光種類，按下"OK"軟體  
 會自動套用相對應的PMT與雷射設定

**OK** Cancel

Acquisition Viewer cellSens

上午 10:10 2021/7/20

# 掃描器設定視窗



掃描模式 **Oneway: 一般 Roundtrip: 往復**

往復掃描速度選擇

掃描速度 數字愈大速度愈慢

影像比例 **1:1 或是 4:3**

影像尺寸 下拉選單中可選擇不同尺寸  
一般掃描器建議使用 **1024x1024**  
高速掃描器尺寸固定為 **512x512**

掃描器可自由  
選轉角度

偵測器電壓校正,建議選**"ON"**,避免調整  
掃描速度時產生過曝現象

Zoom : 3.00 **光學變焦**



# 掃描器設定視窗上半部

Sequential Scan:  None  Line  Frame

Dye & Detector Select

Phase 1

Confocal Aperture: Auto 120  $\mu$ m

Airy Disk x1.00

Mode:  VBF  Lambda

偵測模式，請選VBF

Average:  None  Line  Frame

2 Times

影像平均，可提高S/N ratio  
次數愈多效果愈好，掃描時間相對較久

Sequential Scan:  None  Line  Frame

序列掃描，為避免漏光發生建議一定要開啟，選擇Line模式可同時更新所有螢光頻道，Frame模式一次僅更新一個頻道

Confocal Aperture: Auto 120  $\mu$ m

針孔光圈大小，請啟用Auto功能

Airy Disk x1.00

針孔光圈大小調整，建議使用1.00

Group 1 CH1





# 掃描器設定視窗下半部

序列掃描群組設定

雷射減光片，選取10%時，雷射最小強度可設定為0.01%

螢光種類


所使用PMT

雷射波長與強度


HV: PMT電壓, 代表靈敏度, 建議 < 750

Gain: 放大畫面上所有訊號

Offset: 過濾畫面上背景及雜訊




開始預覽即時影像時，建議影像尺寸  
使用512 x 512 ①，雷射使用0.5% ②偵測  
器HV設定為550，按下”Live x 2” ③預掃  
同時按下”Hi-Lo” ④將即時影像  
改為黑白模式，更容易找到樣本影像




1

3

2




找到樣本後，調整雷射強度及HV 1，讓即時影像上面出現些許紅色訊號(過曝部分) 2 即可，一般而言，雷射以及HV設定數值愈低，對樣本傷害愈小，同時畫質也愈好，調整好後按下"Stop" 3 停止預掃




將影像尺寸調成1024x1024 **1**，選擇  
存檔路徑以及設定檔案名稱 **2**，確認  
Time以及Z **3**均選取"OFF"後，按下  
LSM Start **4**即可開始掃描XY影像，  
掃描完畢後軟體會自動儲存檔案


# XYZ(Z stack)取像 Step by step







開始預覽即時影像時，建議影像尺寸  
使用512 x 512 ①，雷射使用0.5% ②偵測  
器HV設定為550，按下”Live x 2” ③預掃  
同時按下”Hi-Lo” ④將即時影像  
改為黑白模式，更容易找到樣本影像



先確認Z選項選取"ON" 1


調整調節輪移動焦距至樣本的頂端(或底端) 2

此時即時影像應該幾乎無訊號 2

按下"Register" 3 紀錄此時的高度

3

2



調整調節輪移動焦距至樣本的另一端  
此時即時影像應該幾乎無訊號 1

按下"Register" 2 紀錄此時的高度


按下"Optimize" 3 讓軟體自動決定最佳間隔

2  
Register


Optimize 3

1

CH1 Violet, CH2 Alexa Fluor 568



按下"LSM Start" 1 軟體便會開始掃描Z stack影像  
此處會顯示掃描進度與Z軸位置 2



1

2

3

掃描進行同時 1，可以點選已掃檔案 2

選擇ThreeSides 3 可觀察XZ及YZ切面

**Layout**

LSM Imaging x Ocular

▼ Scan Settings

Scanner:

Type:  Galvano  Resonant

Mode:  OneWay  Roundtrip

Interlace:  OFF  2x  4x

Speed: 2.0 us/pixel

Image Size:

Aspect Ratio:  1:1  4:3

Scan Size: 800x800

High voltage correction:  ON  OFF

Pixel: 2.0 usec Line: 2.690 msec Frame: 8.611 sec

▼ Area Settings

Rotation: 0.0 deg

Pan X: 0.00 um

Pan Y: 0.00 um

**PMT Setting**

▼ PMT

Mode:  VBF  Lambda

Average:  None  Line  Frame 2 Times

Sequential Scan:  None  Line  Frame

Dye & Detector Select

Confocal Aperture Auto 120 um

Airy Disk x1.00

HSD1 HSD2

Group 1 CH1

Group 2 CH2

Laser ND Filter:  None  10%

P1  CH1 Violet HSD1 430-470 nm G-1

405 0.3 %

HV 521 V

Gain 1.000 x

Offset 4 %

**Series**

LSM

Time:  ON  OFF

Z:  ON  OFF

Series: 1 min 44 sec

LSM total: 0:01:44.24

Interval: FreeRun

Scan: 0:01:44.24 Rest: 0:00:00.00

Cycle: 5

▼ Z Section

Motor:  Start/End  Range

Current: 3695.99 Origin: 200.00

Slices: 12

Step Size: 0.83

Start: 3703.46 End: 3694.33

**\*Acquire**

Normal Sync Sequence MATL

Imaging

SERIES DONE Append

APPEND NEXT 1

Bleach

LSM DACFC Neuron

Current scan condition

Z position : 9 / 12

Total scanning time : 0:01:44.24

Remaining time : 0:00:22

Next scan start at : -

Live Graph Table Analysis Neuron Neuron\_0001 x

Analysis ROI:   Loop Start: 1 - End: 12

Tile Single ThreeSides Volume

CH1 Violet CH2 Alexa Fluor...

CH1 Violet, CH2 Alexa Fluor 568


Size:800x800 0.265[um/pixel] xy=(8,14) Int:266, 922


**Image List**

- Neuron.oir 07/20/2021 10:22:57.448 AM XY
- Neuron\_0001.oir

也可點選Volume 1 即時觀察3D立體影像

# XYT(Time Lapse)取像 Step by step





Time :  ON  OFF 1

LSM Total: 0:00:24.33 4

Interval: 0:00:05.0 2


Cycle: 5 3

先確認Time選項選取"ON" 1


Interval 2 設定拍攝間隔，以秒為單位  
例如想要五分鐘拍一張，請輸入"300"

Cycle 3 設定拍攝循環數，輸入數值後  
這邊會顯示總拍攝時間 4





掃描完畢後，點選檔案名稱 1  
從下拉選單中選取播放速度後 2  
按下 3 按鈕即可開始播放影片



XYT檔案可以由Export選單中直接輸出成WMV影片檔

Export dialog box details:

- File name: Neuron\_0002.wmv
- Save as type: WMV (\*.wmv)
- CH / Range: CH1 Violet
- Range selection: Current Frame
- Table:

	Start	End	Step	Total
Lambda	0	0	0	0/0
Z Series	0	0	0	0/0
T Series	1	10	1	10/10
- OutputSelection: LSM
- Information: Bit depth 24 Bit/Pixel
- ROI overlay: All ROI, Overlay Color scale checked
- OutputFormat: RGB Color with Merge, 24bit Full Color, Amount
- TIFF: Compression None, LZW
- JPEG: CompressionRatio 70 % \* 1% (lowest) - 100% (highest), 70% (default)
- Movie: Frame rate 3.0 Frame/sec \* 30-0.1
- File: Save properties as ASCII text checked
- Comment: (empty)

File menu and Acquisition settings:

- File menu: Open, Save, Save As, Export (Ctrl+E), Exit
- Acquisition settings: Type (Galvano, Resonant), Mode (OneWay, Roundtrip), Interlace (OFF, 2x, 4x), Speed (2.0 us/pixel), Image Size (512x512), Aspect Ratio (1:1, 4:3), Scan Size (512x512), High voltage correction (ON, OFF), Pixel (2.0 usec Line: 2.114 msec Frame: 4.332 sec)

PMT Setting and Acquire dialog:

- PMT Setting: Mode (VBF), Average (None), Sequential Scan (None), Dye & Detector Select (Phase1), Confocal Aperture (Auto, 120), Airy Disk (x1.00), HSD1, HSD2, Group 1 (CH1), Laser ND Filter (None), P1 CH1 Violet, HV (521 V), Gain (1.000 x), Offset (4 %)
- Acquire dialog: Normal, Sync, Sequence, Imaging (LSM Start, SERIES DONE, APPEND NEXT), LSM (D\CFC, Neuron)

Series and Time Lapse settings:

- Series: LSM Time (ON, OFF), Z (ON, OFF), Series: 4.332 sec
- Time Lapse: Interval (FreeRun), Scan (0:00:04.33), Rest (0:00:00.00), Cycle (10)
- Z Section: Motor (Start/End, Range), Current, Origin (200.00), Slices (12), Step Size (0.83), Optimize



Current scan condition and Image List:


- Current scan condition: Total scanning time (0:00:43.31), Remaining time (-), Next scan start at (-)
- Image List: Neuron.oir (07/20/2021 10:22:57.448 AM), XY, Neuron\_0001.oir

Viewer and cellSens interface:

- Viewer: Analysis ROI, Fluor 568
- cellSens: Image display area showing a fluorescence microscopy image of neurons with red and blue channels.

# 共振式掃描器(Resonant scanner) XYT(Time Lapse)取像 Step by step





點選“Resonant”<sup>1</sup> 可將掃描器切換至高速模式  
選擇“OneWay”<sup>2</sup> 時，速度可達每秒鐘15張  
此時影像尺寸最大值會固定在512x512<sup>3</sup>

LSM Imaging x Ocular

Scan Settings

Scanner:

Type:  Galvano  Resonant

Mode:  OneWay  Roundtrip **1**

Interlace:  OFF  2x  4x

Speed: 0.067 us/pixel

Image Size:

Aspect Ratio:  1:1  4:3

Scan Size: 512x512

High voltage correction:  ON  OFF

Pixel: 0.067 usec Line: 0.063 msec Frame: 33.333 msec

Area Settings

Rotation: 0.0 deg Reset

Pan X: 0.00 um Reset

Pan Y: 0.00 um Reset

PMT Setting x

PMT

Mode:  VBF  Lambda

Average:  None  Line  Frame 2 Times

Sequential Scan:  None  Line  Frame

Dye & Detector Select

Phase1

Confocal Aperture Auto 120 um

Airy Disk x1.00

HSD1 HSD2

Group 1 CH1 CH2

Laser ND Filter:  None  10%

P1 CH1 Violet HSD1 430 -470 nm G-1

405 0.3 %

HV 521 V

Gain 1.000 x

Offset 4 %

P1 CH2 Alexa Fluor 568 HSD2 570 -670 nm G-1

Series x

LSM

Time:  ON  OFF

Z:  ON  OFF

Series: 33.333 msec

Time Lapse

LSM Total: 0:00:00.03

Interval: FreeRun

Scan: 0:00:00.03 Rest: 0:00:00.00

Cycle: 10

Z Section

Motor:  Start/End  Range

Current: Origin: 200.00

Register Move

Slices: 12

Step Size: 0.83

Optimize

Start: 3703.46 End: 3694.33

Register

Acquire x

Normal Sync Sequence MATL

Imaging

LSM Start

SERIES DONE Append

APPEND NEXT 1

Bleach

Simulation Start

LSM DACFC Neuron

Current scan condition

Total scanning time: 0:00:00.03

Remaining time: -

Next scan start at: -

Live Graph Table Analysis Neuron x Neuron\_0001 Neuron\_0002

Analysis ROI:

LUT: Hi-Lo Auto Contrast Reset

Tile Single

CH1 Violet CH2 Alexa Fluor...

CH1 Violet, CH2 Alexa Fluor 568

Size: 1024x1024 0.207[um/pixel] xy=(0,142) Int:221, 238

Image List x


Neuron\_0001.oir

07/20/2021 10:22:57.448 AM


XY

Neuron\_0001.oir


**選擇“RoundTrip” **1** 時，速度可達每秒鐘30張**





若想進一步提高掃描速度，可以點選此處 **1**，在畫面上畫出掃描區域 **2**，即可獲得更高的fps，將尺寸縮小至512x32時，可達極速438 fps



如圖中所示，掃描此影像範圍可達100 fps  
拍攝300張連續圖像只需要3秒鐘時間





如果想要重複上次拍攝條件，請先開啟舊檔案  
然後按下此按鍵，系統就會喚回拍攝參數