

Quality Control of Poly (A)⁺ RNA

The quality of the poly (A)⁺ RNA is done by evaluation of the quality of the 1st strand cDNA. That is done by gel electrophoresis analysis of labeled 1st strand cDNAs.

©*Materials and Reagents*

10X Klenow reaction buffer – Invitrogen Cat.18012-021

Klenow enzyme – Invitrogen Cat.18012-021

Cy3 dUTP or Cy5-dUTP (25 nmole) – Amersham Cat. PA53022 / PA55022

Oligo dT-V (A/C/G) (2 µg/µl)

5X Superscript II reaction buffer - Invitrogen Cat. 10864014

Superscript II reverse transcriptase (200 U/µl) - Invitrogen Cat. 10864014

50X dNTPs (25mM each except for dTTP at 10 mM) - Amersham Cat. 27-2035-01

10X dNTPs (0.25 mM each except dTTP at 0.09 mM – Amersham Cat. 27-2035-01

0.1M DTT - Invitrogen Cat. 10864014

PCR clean up kit – Qiagen Cat. 28106

TE pH 8.0

2X loading buffer (0.08% bromophenol blue, 14% (w/v) sucrose in ultra pure water)

1X TAE buffer

Lambda *Hind* III DNA marker

©*Labeling 1st strand cDNA with Cy-dye*

1. Mix 2 µl RNA (80 ng mRNA) with 0.5 µl oligo dT-V (2 µg/µl)
2. Heat 10 min at 70°C and quickly chill on ice.

3. Add

5X Superscript II buffer	1.0 μ l
0.1 M DTT	0.5 μ l
50X dNTPs (25 mM dNTPs except dTTP at 10 mM)	0.1 μ l
Cy3- or Cy5-dUTP	0.5 μ l
Superscript II reverse transcriptase (200 U/ μ l)	0.4 μ l

4. Incubate at 42°C for 2 hrs.

5. Clean up reaction and elute with 50 μ l EB buffer from Qiagen PCR clean up kit.

6. Dry the eluted DNA down and resuspend in 3 μ l TE pH 8.0.

7. Mix 3 μ l of labeled 1st strand cDNA (labeled 1st strand cDNA from 80 ng mRNA) with 3 μ l of 2X loading buffer.

8. Load 6 μ l into each well.

9. Run a 1.2% agarose gel in 1X TAE buffer at 100mA, 30 min. A labeled DNA molecular weight marker is used for comparison. DNA molecular weight marker is labeled as follows:

©***Labeling Lambda Hind III with Cy-dye***

1. Mix the following components and incubate at 37°C for 1hr

Klenow 10X buffer	4 μ l
10X dNTP (0.25 mM each except dTTP at 0.09 mM)	4 μ l
λ Hind III (0.5 μ g / μ l)	4 μ l
Cy3-dUTP/or Cy5-dUTP	1 μ l
Klenow enzyme	2 μ l

Ultra pure water

25 μ l

2. Clean up reaction with Qiagen PCR clean up kit and elute with 50 μ l EB buffer.
3. Dry the eluted DNA down and resuspend in 10 μ l TE.
4. Mix 3 μ l of the labeled DNA with 3 μ l of 2X loading buffer.
5. Load 6 μ l into each well. Run together with your labeled sample on a 1.2% agarose gel in 1X TAE buffer at 100mA, 30 min.

©QC gel analysis

Scan gel using a Fluorescence Imager (e.g., Typhoon, Variable Mode Imager, Molecular Dynamics, NJ). Generally we look for a fairly bright smear between 500-2000 bp in size.