

Thermo

Proteome Discoverer

User Guide

Software Version 2.1

XCALI-97751 Revision A July 2015

Thermo
SCIENTIFIC

© 2015 Thermo Fisher Scientific Inc. All rights reserved.

Proteome Discoverer, Orbitrap Elite, Orbitrap Velos, Orbitrap Fusion Lumos, and Q Exactive are trademarks; and LTQ, Orbitrap, Thermo Scientific, and Xcalibur are registered trademarks of Thermo Fisher Scientific Inc. in the United States.

The following are registered trademarks in the United States and possibly other countries:

iTRAQ is a registered trademark of AB SCIEX. UniProt is a registered trademark of European Molecular Biology Laboratory Incorporated Association. NIST is a registered trademark of the National Institute of Standards and Technology. Applied Biosystems is a registered trademark of Applied Biosystems, LLC. SEQUEST is a registered trademark of the University of Washington. Intel and Xeon are registered trademarks of Intel Corporation.

Mascot is a registered service mark of Matrix Science Ltd.

TMT is a registered trademark of Proteome Sciences plc in the United Kingdom.

Excel, Microsoft, and Windows are registered trademarks of Microsoft Corporation in the United States and other countries.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history: Release A, July 2015

Software version: Thermo Proteome Discoverer version 2.1, Microsoft Windows 7 with Service Pack 1, Mascot Server 2.1

For Research Use Only. Not for use in diagnostic procedures.

Contents

	Preface	xi
	Related Documentation	xi
	System Requirements	xii
	Cautions and Special Notices	xiv
	Contacting Us	xiv
Chapter 1	Introduction	1
	Features	1
	Proteome Discoverer Application	2
	Study Management	3
	Search Engines	3
	Workflow Editor	4
	Quantification	4
	ProteinCenter Access	5
	Graphical Views	5
	FASTA Databases and Indexes	6
	Data Filtering and Validation	6
	Export Formats	6
	Qual Browser Application	6
	Fragmentation Methods	6
	Peptides and Fragment Ions	7
	MudPIT Experiments	8
	Workflow	8
	Inputs and Outputs	11
	FASTA Databases	11
	Inputs	11
	Outputs	12
	New Features in This Release	12
	Changes in Quantification	13
	Protein FDR Validator Node	15
	New PSM Rank Filter	15
	Parallel Job Execution	15
	Displaying Species Names and Species Maps	15
	Displaying Master Proteins by Default	16
	Exporting Data	16

Chapter 2	Getting Started	17
	Opening the Proteome Discoverer Application	17
	Closing the Proteome Discoverer Application	18
	Using the Start Page	18
	Opening the Start Page	19
	Opening Result Files on the Start Page	19
	Retaining Study Names or Result File Names on the Start Page	19
	Deleting a Study Name or a Result File Name from the Start Page	20
	Clearing the Start Page	20
	Closing the Start Page	21
	Configuring Search Engine Parameters	21
	Configuring the Sequest HT Search Engine	22
	Configuring the Mascot Search Engine	25
	Configuring Parallel Job Execution Parameters	32
	Studies in the Proteome Discoverer Application	33
	Analyses in the Proteome Discoverer Application	33
	Workflow Editor	34
	Performing a Search	36
	Before Performing a Search	36
	Using Studies	36
	Using Analyses	70
	Specifying Quantification Ratios from Selected Sample Groups	82
	Performing the Search	91
	Working with the Search Results	95
	Using the Workflow Editor to Create Workflows	103
	Using Workflow Templates	119
	Creating Specific Types of Workflows	127
	Searching Multiple Sequence Databases with Mascot	134
	Creating a Multiconsensus Report	135
Chapter 3	Using the Proteome Discoverer Daemon Utility	143
	Starting the Proteome Discoverer Daemon Utility in a Window	143
	Selecting the Server	144
	Running the Proteome Discoverer Daemon Utility from the Window	145
	Monitoring Job Execution in the Proteome Discoverer Daemon Utility	148
	Logging in to a Remote Server	149

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System	150
Before You Start	151
Running the Proteome Discoverer Daemon Utility from a Parameter File.	153
Creating a Processing Method That Calls the Proteome Discoverer Daemon Utility	154
Batch Processing with a Processing Method That Calls the Proteome Discoverer Daemon Utility	157
Batch Processing with Multiple Processing Methods	160
Processing MudPIT Samples by Using a Processing Method	163
Running the Proteome Discoverer Daemon Utility from the Command Line.	165
Syntax	166
Examples	168
Chapter 4 Searching for Data	171
Using FASTA Databases	171
Displaying FASTA Files	172
Adding FASTA Files to the Proteome Discoverer Application	174
Exporting FASTA Files	180
Deleting FASTA Files	180
Adding Protein Sequences and References to a FASTA Database File	181
Finding Protein Sequences and References	182
Compiling a FASTA Database.	187
Excluding Individual Protein References and Sequences from a FASTA Database.	191
FASTA Database Utilities Dialog Box Parameters	192
Managing FASTA Indexes.	196
Adding or Modifying FASTA Parsing Rules	209
Identifying Contaminants During Searches	217
Displaying Species Names for Proteins and Peptide Groups	219
Searching Spectrum Libraries	220
Displaying Spectrum Libraries.	220
Adding a Spectrum Library	221
Deleting a Spectrum Library	224
Searching Spectrum Libraries with the MSPepSearch Node	224
Visually Verifying Spectrum Library Matches	225

Defining Chemical Modifications	227
Dynamic Modifications.	227
Static Modifications.	227
Opening the Chemical Modifications View.	228
Adding Chemical Modifications	231
Adding Amino Acids	232
Deleting Amino Acids	233
Deleting Chemical Modifications	234
Importing Chemical Modifications	234
Using the Qual Browser Application	236
Customizing Cleavage Reagents	237
Opening the Cleavage Reagents View	238
Adding a Cleavage Reagent	243
Deleting a Cleavage Reagent	243
Modifying a Cleavage Reagent.	244
Filtering Cleavage Reagent Data	244
Chapter 5 Filtering Data	245
Filtering Results.	245
Filtering PSMs	246
Filtering PSMs with the Percolator Node	247
Filtering PSMs with the Fixed Value PSM Validator Node	248
Filtering PSMs with the Target Decoy PSM Validator Node.	248
Filtering PSMs with the MSF Files Node	248
Filtering PSMs with the Peptide and Protein Filter Node	254
Filtering Proteins.	254
Filtering Phosphorylation Site Probabilities.	254
Filtering with Display Filters	254
Filtering the Data on a Results Report Page.	255
Example 1	261
Example 2	262
Example 3	264
Example 4	266
Filtering with Filter Sets	267
Using Display Filters to Filter Numerical Values by a Specified Precision Value.	270
Finding Common or Unique Proteins in Multiple Searches	272
Applying Filters Specific to Different Searches in Multiconsensus Reports.	276
Chapter 6 Validating Results	277
Target FDRs	278
Peptide Confidence Indicators.	278
Setting Up FDRs.	279

	Calculating FDRs for PSMs, Peptide Groups, Proteins, and Protein Groups	281
	Calculating False Discovery Rates for PSMs	281
	Calculating False Discovery Rates for Peptide Groups	281
	Calculating False Discovery Rates for Proteins	281
	Calculating False Discovery Rates for Protein Groups	282
Chapter 7	Grouping Peptides and Proteins	283
	Grouping Proteins	283
	Master Proteins	283
	Protein Grouping Algorithm	283
	Protein Grouping Workflow	286
	Number of Unique Peptides Column on the Proteins Page	287
	PSMs Identified by Multiple Workflow Nodes	287
	Grouping PSMs	288
Chapter 8	Obtaining Protein Annotation Information	289
	ProteinCenter	290
	Gene Ontology (GO) Database Annotation	291
	Pfam Database Annotation	291
	Entrez Gene Database Annotation	292
	Ensembl Genome Database Annotation	292
	UniProt Database Annotation	293
	Configuring the Proteome Discoverer Application for Protein Annotation	293
	Downloading Files from ProteinCenter	296
	Exporting Files Downloaded from ProteinCenter	296
	Creating a Protein Annotation Workflow	296
	Displaying the Annotated Protein Results	299
	Displaying Predefined GO Protein Annotation Results	299
	Displaying Annotation Aspects	304
	Annotation Aspects View Parameters	317
	Displaying GO Accessions	317
	Displaying Protein Family (Pfam) Annotations	320
	Displaying Entrez Gene Database Identifications	321
	Displaying Ensembl Genome Database Annotations	322
	Displaying UniProt Annotations	323
	Uploading Results to ProteinCenter	323
	Accessing the ProteinCard Page	326

ProteinCard Page Parameters.	328
General Page	328
Keys Page	330
Features Page.	331
Molecular Functions Page	333
Cellular Components Page	334
Biological Processes Page	335
Diseases Page	337
External Links Page	337
GO Slim Categories	338
GO Slim Categories for Molecular Functions	339
GO Slim Categories for Cellular Components	340
GO Slim Categories for Biological Processes	343
Chapter 9 Searching for Post-Translational Modifications	347
Using the ptmRS Node	347
How the ptmRS Node Calculates Sequence and Site Probabilities.	348
Creating a PTM Analysis Workflow with the ptmRS Node.	351
Columns in the Results Report	355
Methylation.	358
Using the Peptide in Protein Annotation and the PSM Grouper Nodes	359
Peptide Positions in Proteins	359
Flanking Residues	360
Modifications of Proteins on the PSMs Page	362
Modifications of Proteins on the Proteins Page	363
Modifications of Peptides on the Peptide Groups Page	364
Viewing PTM Information on the Protein Identification Details View	369
Filtering Phosphorylation Site Probabilities.	374
Example 1	376
Example 2	377
Example 3	378
Chapter 10 Performing Quantification	381
Precursor Ion Quantification.	381
Chemical Elements Supported in Precursor Ion Quantification.	383
SILAC 2plex Methods.	383
SILAC 3plex Methods.	385
Dimethylation 3plex Method	385
18O Labeling Method.	385
Reporter Ion Quantification	385
TMT Quantification.	386
iTRAQ Quantification	388
Precursor Ion Area Detection	389

Performing Quantification	389
Performing Precursor Ion Quantification	389
Performing Reporter Ion Quantification	396
Performing Precursor Ion Area Detection	402
Using a Quantification Method.	408
Setting Up the Quantification Method	408
Specifying the Quantification Channels.	410
Correcting Reporter Ion Quantification Results for Isotopic Impurities . . .	419
Checking the Quantification Method	424
Changing a Quantification Method.	424
Removing a Quantification Method	426
Importing a Quantification Method	426
Exporting a Quantification Method	427
Restoring Quantification Method Template Defaults	427
Searching for Quantification Modifications with Mascot	428
Generating Quantification Ratios	431
Calculating Quantification Results	431
Quantification Methodology	431
Calculating PSM Abundances	432
Classifying PSMs.	436
Calculating Peptide Group Abundances	442
Classifying Peptide Groups	442
Calculating Protein Abundances	443
Normalizing Peptide Groups and Protein Abundances	443
Sample Information Used to Calculate and Display Quantification Results	444
Calculating Group Abundances and Ratios	446
Calculating Standard Errors.	447
Calculating Protein Group Ratios	447
Calculating Areas.	447
Calculating emPAI Values.	448
Filtering Quantification Data with Signal-to-Noise Values	449
Using Signal-to-Noise Values as Quantification Channel Values.	449
Filtering Quantification Data with Average Reporter Ion Signal-to-Noise Values	450
Identifying Isotope Patterns in Precursor Ion Quantification.	451
Viewing Quantification Results.	453
Using the Quantification Page of the Result Summaries	453
Displaying Data Distribution Maps for Quantification	453
Displaying the Quantification Channel Values Chart	459
Displaying the Quantification Spectrum Chart	466
Displaying the Quan Spectra Page.	474
Displaying Quantification Ratio Charts.	475
Treatment of Missing Reporter Ions for Quantification	479
Troubleshooting Quantification	486

Appendix A	FASTA Reference	491
	FASTA Databases	491
	NCBI	491
	UniRef100	492
	SwissProt and TrEMBL	492
	Custom Database Support	493
	Custom Parsing Rule A	493
	Custom Parsing Rule B	493
	Custom Parsing Rule C	493
Appendix B	Chemistry References	495
	Amino Acid Mass Values	495
	Enzyme Cleavage Properties	496
	Fragment Ions	497
Appendix C	Technical and Biological Replicates	499
	Definitions	499
	Specifying Replicates	499
	Index	503

Preface

This guide describes how to use the Thermo Proteome Discoverer™ 2.1 application for peptide and protein mass spectrometry analyses.

Contents

- [Related Documentation](#)
- [System Requirements](#)
- [Cautions and Special Notices](#)
- [Contacting Us](#)

Related Documentation

The Proteome Discoverer application includes complete documentation. In addition to this guide, you can also access the following document as a PDF file from the data system computer:

- *Proteome Discoverer Quick Start*

❖ To view the product manual

From the Microsoft™ Windows™ taskbar, choose **Start > All Programs > Thermo Proteome Discoverer 2.1 > Proteome Discoverer 2.1 User Guide**.

❖ To download user documentation from the Thermo Scientific website

1. Go to www.thermoscientific.com.
2. In the Search box, type the product name and press Enter.
3. In the left pane, select **Documents & Videos**, and then under Refine By Category, click **Operations and Maintenance**.

4. (Optional) Narrow the search results or modify the display as applicable:
 - For all related user manuals and quick references, click **Operator Manuals**.
 - For installation and preinstallation requirements guides, click **Installation Instructions**.
 - For documents translated into a specific language, use the Refine By Language feature.
 - Use the Sort By options or the Refine Your Search box (above the search results display).
5. Download the document as follows:
 - a. Click the document title or click **Download** to open the file.
 - b. Save the file.

For access to the application Help, follow this procedure.

❖ **To view the Proteome Discoverer Help**

From the application window, choose **Help > Proteome Discoverer Help**.

- If information about setting parameters is available for a specific view, page, or dialog box, click **Help** or press the F1 key for information about setting parameters.

System Requirements

The Proteome Discoverer application requires a license. In addition, ensure that the system meets these minimum requirements.

Table 1. System requirements

System	Minimum requirements
Computer	<ul style="list-style-type: none"> • 2 GHz processor • (Recommended) Two Intel™ Xeon™ 6-core processors, 2.4 GHz • 2 GB RAM • 64-bit operating system for the Proteome Discoverer application; 32- or 64-bit operating system for the Proteome Discoverer Daemon utility • Video card and monitor capable of 1280 × 1024 resolution (XGA) • Screen resolution of 96 dpi • 1 TB available on drive C • NTFS format
Software	<ul style="list-style-type: none"> • Adobe™ Reader™ 10 • Microsoft Windows 7 Professional with Service Pack 1 • Microsoft .NET Framework 4.5.1
Mascot SM Server	<ul style="list-style-type: none"> • Mascot Server 2.1 <ul style="list-style-type: none"> – Mascot servers running version 2.1 are usable, but retrieving the result files (protein sequences) from the servers can be a lengthy process because you can only retrieve the protein sequences one at a time. – Mascot servers running version 2.1 should have all available updates, patches, or both from Matrix Science already installed. Ensure that you have installed a patch that enables MIME format for the result files; otherwise, the Proteome Discoverer application cannot receive the search results from the Mascot server. • Mascot Server 2.2: Proteome Discoverer 2.1 does not support error-tolerant searches. • Mascot Server 2.3: Proteome Discoverer 2.1 does not support error-tolerant searches and searches against multiple-sequence databases. • Mascot Server 2.4: Proteome Discoverer 2.1 does not support error-tolerant searches.

Note Ensure that port 28199 is not blocked by firewalls.

Note Before installing the Proteome Discoverer application, ensure that the Windows operating system has the latest Microsoft .NET Framework and Windows updates installed.

Cautions and Special Notices

Make sure you follow the cautions and special notices presented in this guide. Cautions and special notices appear in boxes; those concerning safety or possible system damage also have corresponding caution symbols.

This guide uses the following types of cautions and special notices.



CAUTION Highlights hazards to humans, property, or the environment. Each CAUTION notice is accompanied by an appropriate CAUTION symbol.






IMPORTANT Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the system.

Note Highlights information of general interest.

Tip Highlights helpful information that can make a task easier.

Contacting Us

There are several ways to contact Thermo Fisher Scientific for the information you need. You can use your smartphone to scan a QR code, which opens your email application or browser.

Contact us	Customer Service and Sales	Technical Support
	(U.S.) 1 (800) 532-4752	(U.S.) 1 (800) 532-4752
	(U.S.) 1 (561) 688-8731	(U.S.) 1 (561) 688-8736
	us.customer-support.analyze@thermofisher.com 	us.techsupport.analyze@thermofisher.com 

Contact us

Customer Service and Sales

Technical Support



❖ **To find global contact information or customize your request**

1. Go to www.thermoscientific.com.
2. Click **Contact Us**, select the **Using/Serviceing a Product** option, and then type the product name.
3. Use the phone number, email address, or online form.



❖ **To find product support, knowledge bases, and resources**

Go to www.thermoscientific.com/support.

❖ **To find product information**

Go to www.thermoscientific.com/lc-ms.

Note To provide feedback for this document:

- Send an email message to Technical Publications (techpubs-lcms@thermofisher.com).
 - Complete a survey at www.surveymonkey.com/s/PQM6P62.
-

Introduction

This chapter introduces you to the Proteome Discoverer application and describes its features and functionality.

Contents

- [Features](#)
- [Workflow](#)
- [Inputs and Outputs](#)
- [New Features in This Release](#)

Features

The Proteome Discoverer application is a client-server application that uses workflows to process and report mass spectrometry data. It compares the raw data taken from mass spectrometry or spectral libraries to the information from a selected FASTA database and identifies proteins from the mass spectra of digested fragments. You can use this application to analyze spectral data from all Thermo Scientific and other mass spectrometers. Specifically, the application does the following:

- Works with peak-finding search engines such as Sequest™ HT and Mascot to process all data types collected from low- and high-mass-accuracy mass spectrometry (MS) instruments. The peak-finding algorithm searches the raw mass spectrometry data and generates a peak list and relative abundances. The peaks represent the fragments of peptides for a given mass and charge.
- Produces complementary data from a variety of dissociation methods and data-dependent stages of tandem mass spectrometry.
- Combines, filters, and annotates results from several database search engines and from multiple analysis iterations. The search engines correlate the uninterrupted tandem mass spectra of peptides with databases, such as FASTA. (See [“Using FASTA Databases”](#) on [page 171](#).)

Proteome Discoverer Application

The Proteome Discoverer application consists of the main client, called Proteome Discoverer; the Magellan server; the Proteome Discoverer Daemon utility; and the different client applications that you can connect to the Proteome Discoverer application.

Proteome Discoverer Client

In general, you interact with the main window of the Proteome Discoverer client. You use it for creating and submitting the workflows to be processed; monitoring job progress; configuring the server; administering repository data, such as FASTA files, spectral libraries, and modifications; and displaying the results.

You create workflows in the application's Workflow Editor, create processing and consensus analysis sequences within a study, and then submit them for execution to the Magellan server running a local or remote process.

Magellan Server

The Magellan server is a command-line application that performs all workflow-based data processing. It can run on the same computer or on a remote machine. It is responsible for registering and managing processing nodes that you use within the data processing workflows, processing workflow jobs, storing data, and general data maintenance. The server is extendable with new processing nodes. The Magellan server publishes client gateways by .NET remoting, and multiple clients can connect to these gateways to communicate with the server. The Magellan server maintains a repository of the registered processing nodes, their configurations, and metadata about the available nodes. It also maintains repositories of registered data such as FASTA files, spectral libraries, chemical modifications, and so forth.

The Magellan server interprets user-specified workflows and converts them to jobs that can be processed. It controls the execution of data processing jobs and provides real-time status updates to connected clients.

Proteome Discoverer Daemon Utility

You use the Proteome Discoverer Daemon utility to automatically start workflow jobs from a remote system. You can use this utility to start batches of files to process. It can perform multiple searches on multiple raw data files at any given time. You can use it to perform searches on multiple raw data files taken from multiple samples or replicates from the same sample. Through its command line, you can automate the submission of workflow jobs to the Magellan server. You can also use the Proteome Discoverer Daemon utility to submit jobs manually. For more information on the Proteome Discoverer Daemon utility, see [“Using the Proteome Discoverer Daemon Utility”](#) on page 143.

The Proteome Discoverer application includes the following additional features.

Study Management

The Proteome Discoverer application uses a study management system to organize your experiment data and to automate processing. The study management system comprises studies and analyses. A study contains your input files, the information about the samples contained in these files, information about the treatment of the samples, the workflows used to process the input files, and the results of the processing. An analysis contains the processing and consensus workflows that you assemble in the Workflow Editor to process your data. Each analysis is associated with a study.

Search Engines

The Proteome Discoverer application supports the Sequest HT and Mascot search engines; each produces complementary data. These search engines are available as nodes in the Workflow Editor.

Sequest HT Search Engine

The Sequest HT search engine is distributed by Thermo Fisher Scientific. It can analyze different data types:

- Electron-transfer dissociation (ETD)
- Electron transfer dissociation with HCD or CID activation (EThcD)
- Electron-capture dissociation (ECD)
- Collision-induced dissociation (CID)
- High-energy collision-induced dissociation (HCD)
- Pulsed Q collision-induced dissociation (PQD)

ETD and ECD generate primarily c and z fragment ions with preferences for precursor ion charge states of +3 or higher. CID and HCD generate primarily b and y fragment ions with preferences for precursor ion charge states of +3 or lower. PQD and HCD do not exhibit a low-mass cutoff and are good for reporter-ion experiments.

Frequently, peptides identified by CID, PQD, or HCD are not observed with ETD or ECD, and vice versa, so that combining results from, for example, CID and ETD can enhance sequence coverage. Many times CID and ETD identify the same peptides, often with different precursor ion charge states. Combining ETD and CID results improves confidence in identifications.

The Sequest HT search engine calculates XCorr scores for peptide matches and provides the peptide matches having the best XCorr score for each spectrum. It calculates a preliminary SpScore score and uses it to filter peptide candidates. It calculates XCorr values for peptide spectrum matches (PSMs) only if they pass the SpScore filter. The Sequest HT node calculates the XCorr value for every peptide candidate.

Mascot Search Engine

The Mascot search engine, created by Matrix Science, uses mass spectrometry data to identify proteins from primary sequence databases. For more details on Mascot, visit www.matrixscience.com.

Workflow Editor

For ease of use and greater flexibility, the Proteome Discoverer application offers a dual-workflow search capability to search for matching proteins and peptides.

You can use the application's Workflow Editor to create customized searches and customized results reports. You create two workflows: a customized processing workflow that generates the primary search results from the input data files and a customized consensus workflow to display the results. The Workflow Editor can accept single or multiple input raw data files for its processing workflow and single or multiple MSF files for its consensus workflow. It can search with multiple algorithms and merges results from multiple fragmentation methods. For detailed information about the Workflow Editor, see “Workflow Editor” on page 34.

Quantification

You can perform precursor ion quantification (for example, SILAC), reporter ion quantification (for example, iTRAQ™ and Tandem Mass Tag™ [TMT]), and precursor ion area detection with the Proteome Discoverer application. For details, see “Precursor Ion Quantification” on page 381, “Reporter Ion Quantification” on page 385, and “Precursor Ion Area Detection” on page 389.

SILAC is an isotopically labeled quantification method that uses in-vivo metabolic labeling to detect differences in the abundance of proteins in multiple samples. SILAC uses the Precursor Ions Quantifier node in the Workflow Editor.

iTRAQ and TMT are very similar isobarically labeled quantification methods that use external reagents, or tags, to chemically label proteins and peptides to detect differences in abundances. TMT quantification offers default 2plex, 6plex, and 10plex quantification methods, and iTRAQ offers 4plex and 8plex quantification methods. You can use these methods to create your own quantification templates. iTRAQ and TMT use the Reporter Ions Quantifier node in the Workflow Editor.

The application also offers precursor ion area detection, which you can use to determine the area for any quantified peptide. This type of quantification uses the Precursor Ions Area Detector node. For more information about precursor ion area detection, see “Precursor Ion Area Detection” on page 389.

ProteinCenter Access

The Proteome Discoverer application gives you access to annotation information from ProteinCenter, which is a web-based application that you can use to download biologically enriched annotation information for a single protein. You can download information such as the following:

- Molecular functions, cellular components, and biological processes from the Gene Ontology (GO) database
- Classification information for protein families from the Protein Family (Pfam) database from the Wellcome Trust Sanger Institute (WTSI)
- Gene identifications from the Entrez Gene database maintained by the National Center for Biotechnology Information (NCBI)
- Genomic annotations of genetically sequenced organisms from the Ensembl genome database, which is a joint project of the European Bioinformatics Institute (EBI) and the WTSI
- Post-translational modification information from the UniProt™ database

The data in ProteinCenter is updated biweekly. You can use this information to annotate the proteins in your results report (.pdResult) file if you have an account with ProteinCenter. For information, see [“Obtaining Protein Annotation Information”](#) on page 289. You can also upload search results directly from the Proteome Discoverer application to ProteinCenter.

Graphical Views

The application includes a number of graphical tools to help you analyze your data.

- Data tables in the results report are organized by proteins, protein groups, PSMs, peptide groups, and MS/MS spectrum information. Depending on the nodes that you use in the workflows, you can optionally display a result statistics table and tables on decoy proteins, decoy protein groups, decoy PSMs, and decoy peptide groups. In addition, you can display tables on annotated modifications, found modifications, and unknown modifications.
- Several graphical views contain detailed information about the selected peptides and proteins. You can display more than one view to perform a comparative analysis of your selected peptide or proteins. Refer to the Help.
- Data distribution maps display the areas and ratios calculated by the Peptide and Protein Quantifier node as data distribution maps that show the distribution of values across the available files, samples, and sample groups, or across ratios and ratio groups.

FASTA Databases and Indexes

You can download FASTA databases that you can use when searching for and analyzing data in the Proteome Discoverer application and its indexes. See [“Using FASTA Databases”](#) on [page 171](#).

Data Filtering and Validation

A number of protein and peptide filtering, validation, and grouping options are available to help you sort and filter your data. For information on Proteome Discoverer’s filtering capabilities, see [“Filtering Data”](#) on [page 245](#). For information on validating your data, see [“Validating Results”](#) on [page 277](#). For information on grouping, see [“Grouping Peptides and Proteins”](#) on [page 283](#).

Export Formats

You can export the spectra from protein and peptide reports in standard spectrum data formats, such as MZDATA, DTA, MZML, and MGF. You can also export search results to XML and tab-delimited TXT files. In addition, you can export annotated spectra for selected peptides into a ZIP file that includes an HTML page with peptide information and links to spectrum images. The Help describes how to export your data to these and other formats.

Qual Browser Application

The Proteome Discoverer application can start the Qual Browser application, if Xcalibur is installed on the same computer. You can use the Qual Browser application to view the entire ion chromatogram and browse individual precursor and MSⁿ data. You can filter the results in a variety of ways, for example, to produce a selected ion chromatogram. The Qual Browser application automatically displays the elemental composition, theoretical mass, delta values, and ring and double-bond (RDB) equivalents for your high-resolution data. For more information about the Qual Browser application, see [“Using the Qual Browser Application”](#) on [page 236](#).

Fragmentation Methods

The Proteome Discoverer application supports the following fragmentation types:

- CID – Uses the collision-induced dissociation method of fragmentation, where molecular ions are accelerated to high kinetic energy in the vacuum of a mass spectrometer and then allowed to collide with neutral gas molecules such as helium, nitrogen, or argon. The collision breaks the bonds and fragments the molecular ions into smaller pieces.

- ECD – Uses the electron capture dissociation method of fragmentation, where multiply protonated molecules are introduced to low-energy free electrons. Capture of the electrons releases electric potential energy and reduces the charge state of the ions by producing odd-electron ions, which easily fragment.
- HCD – Uses the high-energy collision-induced dissociation method of fragmentation, where the projectile ion has laboratory-frame translational energy higher than 1 keV. HCD produces a highly abundant series of reporter ions for TMT and iTRAQ quantification.
- ETD – Uses the electron transfer dissociation method of fragmentation, where singly charged reagent anions transfer an electron to multiply protonated peptides within an ion trap mass analyzer to induce fragmentation. ETD cleaves along the peptide backbone while side chains and modifications such as phosphorylation are left intact. This method is used to fragment peptides and proteins. ETD activation includes the ETciD activation type.
- EThcD – Uses the electron transfer dissociation method of fragmentation but produces additional b and y ions.
- IRMPD – With the infrared multi-photon dissociation method of fragmentation, an infrared laser is directed at the ions in the vacuum of the mass spectrometer. The target ions absorb multiple infrared photons until they reach more energetic states and begin to break bonds, resulting in fragmentation.
- PQD – Uses the pulsed Q collision-induced dissociation method of fragmentation, where precursor ions are activated at a high value, a time delay occurs to allow the precursor to fragment, and then a rapid pulse is applied at a low value where all fragment ions are trapped. The product ions can then be scanned out of the ion trap and detected. PQD fragmentation produces precise, reproducible fragmentation and has been used for iTRAQ peptide quantification on the LTQ™ mass spectrometer using both electrospray and MALDI source ionization.


Peptides and Fragment Ions

The types of fragment ions observed in an MS/MS spectrum depend on several factors, such as the primary sequence, the energy source, and the charge state.

Fragment ions of peptides are produced by a CID process in which a peptide ion is fragmented in a collision cell. Low-energy CID spectra are generated by MS/MS and ESI, and are sequence-specific. The fragment ion spectra contain peaks of the fragment ions formed by cleavage of the peptide bond and are used to determine the amino acid sequence. A fragment must have at least one charge for it to be detected. If this charge is retained on the N terminal fragment, the ion is classed as a, b, or c. If the charge is retained on the C terminal fragment, the ion type is x, y, or z. A subscript indicates the number of residues in the fragment.

In addition to the proton carrying the charge, c ions and y ions abstract an additional proton from the precursor peptide, as shown in [Figure 1](#).

Figure 1. Structures of six singly charged sequence ions



MudPIT Experiments

Multidimensional Protein Identification Technology (MudPIT) experiments investigate complex proteomes by applying multidimensional chromatography to the samples before acquisition in the mass spectrometer. Typically, this process results in several dozen or even a few hundred fractions that are separately analyzed by LC/MS, resulting in one raw data file per sample fraction. Analyzing gel slices or performing in-depth follow-up acquisitions also results in multiple fractions. Because all these fractions belong to the same sample, the Proteome Discoverer application can process all raw data files from these fractions as one contiguous input file and generate a single result file. For detailed information about processing MudPIT samples, see [“Using the Proteome Discoverer Daemon Utility”](#) on [page 143](#).

Workflow

Using the standard Proteome Discoverer workflows involves the following steps when you process, analyze, and interpret mass spectrometry data. [Figure 2](#) on [page 10](#) shows these steps graphically.


1. If you are using Sequest HT, download a FASTA database, if necessary.
2. Create a study and an analysis. Identify the raw input file or files containing the data in the sample.
3. In the Workflow Editor, create a processing workflow and a consensus workflow, and set the parameters of these nodes.

4. Begin a search of the raw data. The Proteome Discoverer application initiates a search of the FASTA database.

Note You can filter data according to false discovery rates that you define through the use of decoy databases that you specify in the workflow.

5. Sort and filter the search report, generate graphs and views, and interpret the search results.
6. (Optional) Review the quantification results and change parameters.
7. (Optional) Reanalyze the quantification results.
8. (Optional) Export results and data to other applications.

Figure 2. The Proteome Discoverer workflow



Inputs and Outputs

The Proteome Discoverer application can accept several different file formats as input and can export data in several formats.

FASTA Databases

The Proteome Discoverer application includes FASTA databases, including several example FASTA databases and example RAW data files. Use these files when exploring and learning how to use the application. For a detailed description of the different types of FASTA databases and their purpose, see [“Using FASTA Databases”](#) on page 171.

Inputs

You can use the application with data from the following mass spectrometers:

- Ion traps, such as the LTQ Velos
- Orbitraps, such as the Q Exactive™
- Hybrid mass spectrometers, which combine ion traps and Fourier transform ion cyclotron resonance (FT-ICR) or Orbitrap™ analyzers, such as the LTQ Orbitrap Velos™ or the Orbitrap Fusion Lumos™.
- Quadrupole time-of-flight (QTOF) mass spectrometers
- Triple quadrupole mass spectrometers

The application accepts the following file types as input to the processing workflow:

- Xcalibur raw data files contain raw data collected from a mass spectrometer.
- Mascot Generic Format (MGF) files are mass spectral files produced during Mascot analysis. They contain a list of precursor ions, their fragments, and the masses of the fragments.
- MZDATA files are common data format files developed by the Human Proteome Organization (HUPO) for proteomics mass spectrometry data. These files are in version 1.05 format. They are exported with XML indentation enabled so that the different XML tags are broken into multiple lines instead of merged into one line.
- MZXML files are standard 2.x mass spectrometer data format files developed at the Seattle Proteome Center at the Institute for Systems Biology (ISB) that contain a list of precursor ions, their fragments, and the masses of each fragment.
- MZML files are a combination of .mzData and .mzXML formats developed by the Human Proteome Organization Standard Initiative (HUPO-PSI) and the ISB. The Proteome Discoverer application supports version 1.1.0 of the MZML format.

The application accepts MSF files as input to the consensus workflow.

Outputs

The Proteome Discoverer application creates MSF files as output to its processing workflow.

The application creates .pdResult files as standard output to its consensus workflow.

You can also export spectra in the following formats:

- DTA Archive (DTA) files are files containing MSⁿ data for single or grouped scans.
- Mascot Generic Format (MGF) files are mass spectral files produced during Mascot analysis. They contain a list of precursor ions, their fragments, and the masses of the fragments.
- MZDATA files are common data format files developed by the HUPO-PSI for proteomics mass spectrometry data. These files are in version 1.05 format. They are exported with XML indentation enabled so that the different XML tags are broken into multiple lines instead of merged into one line.
- MZML files are a combination of MZDATA and MZXML formats developed by the HUPO-PSI and the ISB. The Proteome Discoverer application supports version 1.1.0 of the MZML format.

You can export peptide and protein search results in the following formats:

- ProtXML files contain protein identifications from MS/MS-derived peptide sequence data.
- PepXML files contain peptides that are included in the results of searches performed by the Sequest HT and Mascot search engines. They are in PepXML format version 1.14, which is an open data format developed by SPC/Institute for Systems Biology for storing, exchanging, and processing peptide sequence assignments from MS/MS scans.
- Tab-delimited TXT files are in a simple text format that stores tabular data and is widely used to exchange data between different computer programs.

New Features in This Release

The 2.1 release of the Proteome Discoverer application focuses on improvements to the quantification process but includes other new features as well.

Changes in Quantification

The 2.1 release makes a number of changes to the way the application quantifies data.

New Method of Calculating Quantification Results

The Proteome Discoverer application uses a new methodology to calculate the quantification results. When you include the Peptide and Protein Quantifier node in the consensus workflow, the node performs these basic steps to calculate the quantification results:

1. Calculates peak area values for precursor area quantification and signal-to-noise (S/N) values (for the Orbitrap analyzer) or abundance values for reporter ion quantification.
2. Sums the quantification channel values of the PSMs to the peptide groups and proteins.
3. Optionally normalizes the abundance values so that all channels have the same total abundance.
4. Optionally scales the normalized abundance values for each protein and peptide so that their sum is 100 times the number of quantification channels.
5. Groups the optionally normalized and scaled abundance values according to the sample grouping that you specified when you set up the analysis.
6. Calculates the quantification ratios from the grouped (averaged) abundance values according to the ratios selected when you set up the analysis.

The node also calculates areas from the data provided by the Precursor Ions Quantifier node and the Precursor Ions Area Detector node. In addition, it calculates the emPAI values.

For detailed information about the method that the Peptide and Protein Quantifier node uses to calculate quantification results, see [“Quantification Methodology”](#) on [page 431](#).

New Quantification Columns in the .pdResult File

The new quantification methodology adds the following new columns to the Proteins and Peptide Groups pages of the .pdResult file:

- Abundances
- Abundance Counts
- Abundance Ratios
- Abundance Ratios (log2)
- Abundance Variability [%]
- Normalized Abundances
- Grouped Abundances
- Grouped Abundance Counts
- Grouped Abundance Standard Errors [%]
- # Razor Peptides
- Protein Quan Usage
- Peptide Quan Usage

For information on these columns on the Proteins page, see the Help. For information on these columns on the Peptide Groups page, see the Help.

New Parameters for the Peptide and Protein Quantifier Node

To support the new quantification methodology, the Peptide and Protein Quantifier node in the consensus workflow has added the following new parameters:

- Peptides to Use
- Reporter Abundance Based On
- Average Reporter S/N Threshold
- Normalization Mode
- Proteins for Normalization
- Scaling Mode
- Show Quan Value Counts
- Show Quan Ratios As
- Replace Missing Values with Minimum Value
- Reject Quan Results with Missing Channels

These parameters are explained in the Help.

Adding Impurity Correction Factors to TMT 10plex Kits

You can now add isotopic impurity correction factors and variant relations as they are stated in the product data sheets (the Certificates of Analysis) to the TMT 10plex kits. You must create a new quantification method to apply these factors. For instructions, see [“Excluding PSMs with High Levels of Coisolation”](#) on page 435.

Filtering Quantification Data with Signal-to-Noise Values

For Orbitrap data, you can filter reporter ion quantification spectra that have too much variability in the intensities of the reporter ions by using signal-to-noise values instead of intensities. The Peptide and Protein Quantifier node offers two new parameters to offset the effects of the variability that results in Orbitrap data when the mass peak arises from a smaller population of ions.

When the number of ions is too small in Orbitrap data, you can optionally use the Reporter Abundance Based On parameter to select signal-to-noise (S/N) values of the reporter peaks instead of their intensities as quantification channel values. For more information on this type of filtering, see [“Using Signal-to-Noise Values as Quantification Channel Values”](#) on page 449.

You can filter quantification channel values by the sum of the S/N values of the reporter ion peaks. Use the Average Reporter S/N Threshold parameter to specify an average reporter S/N threshold value that determines which PSMs the application excludes from quantification. The application filters quantification values by the sum of the S/N values of the reporter ion peaks and excludes all PSMs that have an average S/N value smaller than this threshold value.

This parameter adds a new column, Average Reporter S/N Threshold, to the PSMs and Quan Spectra pages of the .pdResult report. For more information on this type of filtering, see [“Filtering Quantification Data with Average Reporter Ion Signal-to-Noise Values”](#) on page 450.

New TMT Quantification Methods

This release adds the low-resolution TMTe 6plex quantification method, which you must use when you acquire MS/MS data in an ion trap or in an instrument with low resolution. In this method, different channel impurities can affect each other (mass differences between ^{15}N and ^{13}C are not resolved).

The release also includes the low-resolution iodo TMTe 6plex quantification method.

Protein FDR Validator Node

The Protein FDR Validator node now bases its validation on protein scores calculated from the posterior error probability (PEP) values of the peptides, if these values are available. If no PEP values are available, the node calculates the protein score from the search engine scores as it did in the Proteome Discoverer 2.0 application. For detailed information on this node and the methodology that it uses, see the Help.

New PSM Rank Filter

The Percolator, Fixed Value PSM Validator, and Target Decoy PSM Validator nodes in the processing workflow have a new rank filter called Maximum Rank, which specifies a maximum search engine rank threshold and filters out all PSMs with a search engine rank higher than this value. [“Filtering PSMs with the MSF Files Node”](#) on page 248 describes this filter.

Parallel Job Execution

To avoid overloading your data system computer, you can specify how many processing and consensus workflows to execute at the same time by setting two new configuration parameters, Max. Number of Processing Workflows in Parallel Execution and Max. Number of Consensus Workflows in Parallel Execution. You must restart the application after you set these parameters. For details, see [“Configuring Parallel Job Execution Parameters”](#) on page 32.

Displaying Species Names and Species Maps

To assist you in identifying the species present in a sample, the application extracts from the FASTA database the species names of the proteins. It displays and annotates them as colored entries in a distribution map column or as semicolon-separated strings in a text column on the Proteins and Peptide Groups pages of the .pdResult file. See [“Displaying Species Names for Proteins and Peptide Groups”](#) on page 219.


Displaying Master Proteins by Default

You can display only master proteins on the Proteins page by default. Use the Display Filter node in the consensus workflow, process the workflow, and open the Proteins page of the .pdResult file. If you want the Proteins page to display other proteins by default, you can use a filter set to display the appropriate information. For more information, see the Help.

Exporting Data

This release offers three new ways to export data.

Exporting a FASTA File Downloaded from ProteinCenter

The FASTA Files view now includes a new icon,  **Export**, that you can use to export a FASTA file that you downloaded from ProteinCenter. For instructions on exporting a FASTA file, see “Exporting FASTA Files” on [page 180](#).

Exporting Search Results as Excel Workbooks

You can export data in the .pdResult file as Microsoft Excel™ workbooks, including hierarchical data in three layers as proteins (layer 1), their connected peptides (layer 2), and their connected PSMs (layer 3). You can select which data to export on each level from all the visible results on a .pdResult file page. To export search results from the .pdResult file to an Excel workbook, see the Help.

You can copy the following items on the .pdResult file to the Clipboard so that you can paste them into a third-party application file such as Excel:

- Rows
- Rows with headers
- Ranges of cells

Exporting Search Results to a Tab-Delimited Text File

When you export search results to a tab-delimited text file, the Export to Text (tab-delimited) dialog box includes a new option, Generate R-Friendly Headers. This option replaces the special characters in the headers of the columns on the pages of the .pdResult file with the appropriate characters instead of dots when you export data to R, the open source environment for statistical computing and graphics. For details, see the Help.

Getting Started

This chapter describes how to set up a search in the Proteome Discoverer application.

Contents


- [Opening the Proteome Discoverer Application](#)
- [Closing the Proteome Discoverer Application](#)
- [Using the Start Page](#)
- [Configuring Search Engine Parameters](#)
- [Configuring Parallel Job Execution Parameters](#)
- [Studies in the Proteome Discoverer Application](#)
- [Analyses in the Proteome Discoverer Application](#)
- [Workflow Editor](#)
- [Performing a Search](#)
- [Creating a Multiconsensus Report](#)

To start performing a search immediately, see [“Performing a Search”](#) on [page 36](#).

Opening the Proteome Discoverer Application


Open the Proteome Discoverer application by using the Start menu or by clicking the desktop icon.

❖ To open the Proteome Discoverer application

From the Start menu, choose **Programs > Thermo Proteome Discoverer 2.1** or click the **Proteome Discoverer** icon, , on your desktop.

The Proteome Discoverer Start Page opens, as shown in [Figure 3](#).

Figure 3. Proteome Discoverer Start Page



For information on the features of the Proteome Discoverer interface and how to customize them, refer to the Help.

To open a .pdResult file, refer to the Help.

To use the Start Page, see [“Using the Start Page”](#) on [page 18](#).

Closing the Proteome Discoverer Application

Save your changes before you exit the Proteome Discoverer application, because it does not prompt you.

To close a results report file before you close the application, refer to the Help.

❖ To close the Proteome Discoverer application

1. (Optional) Save any changes that you made to the report, study, or workflow.
2. Choose **File > Exit**.

—or—

Click the **X** in the upper right corner of the main Proteome Discoverer window.

The Proteome Discoverer application closes.

Using the Start Page

The Start Page provides links to the most important commands used in working with the Proteome Discoverer application. It features three areas:

- Start, which contains the commands required to create a new study or analysis, open studies, and open results
- Recent Studies, which lists the 20 most recently opened studies

- Recent Results, which lists the 20 most recently opened .pdResult files

The links on the Start Page correspond to File menu commands.

- [Opening the Start Page](#)
- [Opening Result Files on the Start Page](#)
- [Retaining Study Names or Result File Names on the Start Page](#)
- [Deleting a Study Name or a Result File Name from the Start Page](#)
- [Clearing the Start Page](#)
- [Closing the Start Page](#)

Opening the Start Page

❖ To open the Start Page

Start the Proteome Discoverer application, as described in “[Opening the Proteome Discoverer Application](#)” on page 17.

–or–

Choose **View > Start Page** or click the **Start Page** icon, .

Opening Result Files on the Start Page

For instructions on opening .pdResult files on the Start Page, refer to the Help.

Retaining Study Names or Result File Names on the Start Page


When the number of study or file names exceeds 20 on the Start Page, the Proteome Discoverer application begins removing individual files, oldest first, to keep the number at 20. However, you might want to continually work with some studies or result files over a period of time. You can prevent their names from disappearing from the lists of recent items by “pinning” them.

❖ To retain study names or .pdResult file names on the Start Page

1. On the Start Page, hold the mouse pointer over the name of the study or .pdResult file that you want to retain.
2. Click the pin symbol that appears to the left of the entry.

[Figure 4](#) shows the Start Page with pinned entries. A pinned entry still changes its position in the list, depending on its age, but it does not disappear from the list. You can remove an item from its “pinned” status by clicking the pin.

Figure 4. Pinned items on the Start Page



Deleting a Study Name or a Result File Name from the Start Page

You can delete the name of an individual study name or a .pdResult file name from the Start Page.

❖ To delete an individual study or results file name from the Start Page

1. Right-click the name of the .pdResult file that you want to delete.
2. Choose **Remove from List** on the shortcut menu.

—or—

1. In Windows Explorer, delete the study or .pdResult file.
2. On the Start Page, click the name of the study or .pdResult file that you deleted.

An error message box appears.

3. Click **OK**.

Clearing the Start Page


You can delete all the study names or file names from the Start Page at once.

❖ To delete all study names or results file names from the Start Page

1. Choose **File > Recent Studies > Clear Recent Study List** or **File > Recent Results > Clear Recent Result List**.

The Clear Recent Study List dialog box appears for studies, as shown in [Figure 5](#), or the Clear Recent Result List dialog box appears for result files.

Figure 5. Clear Recent Study List



2. Click **OK**.

All the study or results file names disappear from the Start Page and the File menu.

Closing the Start Page


❖ To close the Start Page

Click the **X** on the Start Page tab.

Configuring Search Engine Parameters


Before you begin your search, you can configure certain search parameters for the Sequest HT and Mascot search engines.

❖ To configure search parameters

1. Choose **Administration > Configuration**, or click the **Edit Configuration** icon, .

The Administration page changes to the Configuration view. By default, this view is set to the Mascot configuration view shown in [Figure 6](#).

Figure 6. Default (Mascot) configuration view of the Administration page



2. Follow these procedures:

- [Configuring the Sequest HT Search Engine](#)
- [Configuring the Mascot Search Engine](#)

Configuring the Sequest HT Search Engine

Before using the Sequest HT search engine, set the parameters that configure the SequestHT node and the Sequest protein-scoring calculation.

❖ To configure the Sequest HT search engine

1. On the Administration page, expand the **Sequest** folder under Processing Settings in the Configuration area.
2. Click **Sequest HT Node**.

All of the SequestHT node configuration parameters shown in [Figure 7](#) appear.

Figure 7. Sequest HT node configuration parameters

1. Workload Level	
Automatic	True
Number of Spectra Processed At Once	3000
Number of Parallel Tasks	0
2. XCorr Confidence Thresholds (low-resolution data)	
z=1: High Confidence XCorr	1.5
z=1: Medium Confidence XCorr	0.7
z=2: High Confidence XCorr	2
z=2: Medium Confidence XCorr	0.9
z=3: High Confidence XCorr	2.5
z=3: Medium Confidence XCorr	1.2
z>=4: High Confidence XCorr	3
z>=4: Medium Confidence XCorr	1.5
3. XCorr Confidence Thresholds (high-resolution data)	
z=1: High Confidence XCorr	1.2
z=1: Medium Confidence XCorr	0.7
z=2: High Confidence XCorr	1.9
z=2: Medium Confidence XCorr	0.8
z=3: High Confidence XCorr	2.3
z=3: Medium Confidence XCorr	1
z>=4: High Confidence XCorr	2.6
z>=4: Medium Confidence XCorr	1.2

For information on these parameters, refer to the Help

3. In the Automatic box, specify whether you want the Proteome Discoverer application to automatically estimate the workload level.

The default is True, which means that the application automatically estimates the workload level.

4. (Optional) If you set the Automatic parameter to False, do the following:
 - a. In the Number of Spectra Processed at Once box, specify the maximum number of spectra that the Sequest HT search engine can process at once.

The minimum value is 1000, and there is no maximum. The default is 3000.

The larger the value, the more memory is required.

- b. In the Number of Parallel Tasks box, specify the number of search tasks that Sequest HT can perform at the same time.

The minimum value is 0, and there is no maximum. The default is 0.

If you set this parameter to 0, this search engine performs as many parallel tasks as the number of available CPUs can handle.

5. If you are using the Sequest HT search engine to search low-resolution data, set the XCorr confidence thresholds under the XCorr Confidence Thresholds (low-resolution data) parameter.


The application uses these values only when you include the Fixed Value PSM Validator node in a workflow. The Target Decoy PSM Validator node and the Percolator node assign false discovery rates differently.

The default values appear in [Figure 7](#).

6. If you are using the Sequest HT search engine to search high-resolution data, set the XCorr confidence thresholds under the XCorr Confidence Thresholds (high-resolution data) parameter.


The application considers these values only when the workflow includes the Fixed Value PSM Validator node.

The default values appear in [Figure 7](#) on [page 23](#).


7. To apply any changes to the settings of the node configuration parameters, click the **Apply** icon, .

The message box shown in [Figure 8](#) appears:

Figure 8. Administration message box



8. Click **OK**.


Note Click the **Reset** icon, , to return to the default values.

❖ To configure the Sequest HT protein-scoring calculation

1. On the Administration page, expand the **Sequest** folder under Processing Settings.
2. Click **Sequest Protein Score**.

The Protein Scoring Options configuration parameter shown in [Figure 9](#) appears.

Figure 9. Sequest Protein Scoring Options configuration parameter



For information on this parameter, refer to the Help.

3. In the Peptide Relevance Factor box, specify a factor to apply to the protein score.


The Proteome Discoverer application calculates the protein score for a protein as follows:


$$\text{protein_score} = (\text{sum_of_all_cross-correlation_factors_of_0.8_or_above}) + (\text{peptide_charge} \times \text{peptide_relevance_factor})$$

The minimum value is 0.0, and the maximum value is 0.8. The default is 0.4.

You can review the setting of the Sequest HT protein-scoring configuration parameters by choosing View > Result Summaries, clicking the Configuration tab, clicking the Consensus WF tab, and viewing the Configuration for Protein Scorer section (see Figure 10) on the Consensus WF subpage. For information on the Result Summaries, refer to the Help.


Figure 10. Protein-scoring configuration parameter settings in the Protein Scorer section of the Search Summary



- To apply any changes to the settings of the Sequest protein-scoring configuration parameters, click the **Apply** icon,  **Apply**.

The message box shown in Figure 8 on page 24 appears.

- Click **OK**.

Note Click the **Reset** icon,  **Reset**, to return to the default values.

Configuring the Mascot Search Engine

Before using the Mascot search engine, you must direct the Proteome Discoverer application to the location of the Mascot server and configure the parameters that control access to the server. If your Mascot search fails, the troubleshooting guidelines can help you check for server problems.

- [Directing the Proteome Discoverer Application to the Mascot Server Location](#)
- [Configuring Mascot Parameters](#)
- [Troubleshooting Failed Mascot Searches](#)

Directing the Proteome Discoverer Application to the Mascot Server Location

Open a web browser and try to access the Mascot server through its URL. If you cannot, the Mascot server might not be running, or the URL might not be correct. In this case, contact your system administrator to assist you.

To test the connection between the Proteome Discoverer application and the Mascot server, follow this procedure.

❖ To test the connection to the Mascot server

1. Start the Mascot Communication Tester by choosing **Tools > Mascot Communication Tester**.

The Mascot Communication Tester window appears, as shown at the top of [Figure 11](#). The Proteome Discoverer application automatically populates the fields of the window, including the URL in the Mascot Server Url box.

2. If the URL does not automatically appear or you want to enter it manually, enter the URL of the Mascot server in the Mascot Server Url box.


This is the same URL that you entered in the Mascot Server URL box in the Configuration view of the Proteome Discoverer application.

If the Web server that hosts the Mascot server enables any user authentication, provide a user name and a password for this authentication in the corresponding boxes in the Mascot Communication Tester window.

3. Click **Test Communication**.

An automated communication test starts, as shown in [Figure 11](#).

Figure 11. Test running in the Mascot Communication Tester window



The upper pane displays the progress of the communication test, and the lower pane provides detailed information. When errors occur, the lower pane displays the particular test that failed, which helps you determine why the connection failed. You can click the blue highlighted links shown in the lower pane to test whether the communication works from within a Web browser.

If this information does not resolve the problem, use the Copy Info button to copy the log information into an e-mail and send it to Thermo Fisher Scientific Customer Service at the e-mail address given in “[Contacting Us](#)” on [page xiv](#).

Configuring Mascot Parameters

Before using the Mascot search engine, set the parameters that configure the Mascot server, the Mascot node, and the Mascot protein-scoring calculation.

❖ To configure the Mascot server

1. Choose **Administration > Configuration**.
2. On the Administration page, expand the **Mascot** folder under Processing Settings in the Configuration section.


3. Click **Mascot Server**.

The parameters shown in [Figure 12](#) appear.

The Proteome Discoverer application generates an MGF file that contains the search settings and all mass spectral information. It submits this file to the Mascot server through a web server, which might have a file size limitation. A search that generates large amounts of data—for example, a search with multiple raw data files—could create an MGF file that exceeds this limitation. The Max. MGF File Size [MB] parameter avoids this limitation by performing several separate Mascot searches and merging the results.

4. To split the MGF file and avoid any potential file-size limitations on the web server, enter the maximum size that the MGF file can be in the Max. MGF File Size [MB] box, in megabytes, as shown in [Figure 12](#).

Figure 12. Maximum MGF file size for Mascot



1. Mascot Server	
Max. MGF File Size [MB]	500
Mascot Server URL	
Number of attempts to submit the search	20
Time interval between attempts to submit a search [sec]	90
2. Mascot Server Authentication	
Mascot User Name	
Mascot Server Password	
3. Web Server Authentication	
User Name	
Password	

Make sure the size is less than the file size permitted by the web server.

The minimum file size is 20, and there is no maximum. The default file size is 500 megabytes.

For information on these parameters, refer to the Help.

5. In the Number of Attempts to Submit the Search box, specify the number of times that the application tries to submit the search when the Mascot server is busy.

The minimum value is 0, and there is no maximum value. The default is 20.


6. In the Time Interval Between Attempts to Submit a Search [sec] box, specify the interval of time, in seconds, that elapses between attempts to submit a search when the Mascot server is busy.

The minimum value is 20, and there is no maximum value. The default is 90 seconds.

7. If you are accessing a Mascot server through your own network and security for that server is enforced, enter your user name and password in the boxes under Mascot Server Authentication.


8. If you are accessing a Mascot server through the web and security for that server is enforced, enter your user name and password in the boxes under Web Server Authentication.

- If you changed any settings of the Mascot server configuration parameters, click




The message box shown in [Figure 13](#) appears:

Figure 13. Administration message box



- Click **OK**.

Note Click the **Reset** icon,  **Reset**, to return to the default values.

❖ **To configure the Mascot node**

- On the Administration page, expand the **Mascot** folder under Processing Settings in the Configuration section.
- Click **Mascot Node**.


The configuration parameters for the Mascot node now appear, as shown in [Figure 14](#).

Figure 14. Configuration parameters for the Mascot node

1. Default Confidence Thresholds	
Significance High	0.01
Significance Middle	0.05


For information on these parameters, refer to the Help.

- Set the Default Confidence Thresholds parameters as follows:
 - Significance High: Calculates the thresholds for high-confidence peptides. The application automatically sets this value to the calculated relaxed significance when it performs a decoy search. The minimum value is 0.0, and the maximum value is 1.0. The default is 0.01.
 - Significance Middle: Calculates the thresholds for medium-confidence peptides. The application automatically sets this value to the calculated relaxed significance when it performs a decoy search. The minimum value is 0.0, and the maximum value is 1.0. The default is 0.05.
- To apply any changes to the settings of the Mascot node configuration parameters, click




The message box shown in [Figure 13](#) appears:

Figure 15. Administration message box



5. Click **OK**.


Note Click the **Reset** icon,  **Reset**, to return to the default values.

❖ **To configure the Mascot protein-scoring calculation**

1. On the Administration page, expand the **Mascot** folder under Processing Settings.
2. Click **Mascot Protein Score**.

The Protein Scoring Options configuration parameters shown in [Figure 16](#) appear.

Figure 16. Protein Scoring Options configuration parameters



For information on these parameters, refer to the Help.

You can review these parameter settings in the Protein Scorer section of the Consensus Workflow page of the Result Summaries. ([Figure 10](#) on [page 25](#) shows this section for the equivalent Sequest HT parameters.) For information on the Result Summaries, refer to the Help.


3. Specify whether to use MudPIT scoring or normal scoring.

MudPIT scoring, an option for Mascot searches, is used for generating the protein score from data sets with large numbers of peptide matches, as in a complex proteome digest. Mascot stores 10 possible matches for each spectrum. The protein score is normally a sum of all possible matches. However, this might result in a high “protein” score that comprises numerous low-scoring erroneous matches. The Use MudPIT Scoring parameter prevents low-scoring peptide matches from contributing to a falsely reported high-scoring protein.

- (Default) Automatic: Automatically determines whether to use MudPIT scoring.
- True: Always uses MudPIT scoring.
- False: Never uses MudPIT scoring.


For information about MudPIT scoring, consult the Mascot Help.

4. If you changed any settings of the Mascot protein-scoring configuration parameters, click




The message box shown in [Figure 17](#) appears:

Figure 17. Administration message box



5. Click **OK**.

Note Click the **Reset** icon,  **Reset**, to return to the default values.

Troubleshooting Failed Mascot Searches

If all your searches using the Mascot search engine fail, follow these instructions to locate the problem.

❖ To troubleshoot failed Mascot searches

1. Verify that the Mascot server is running and accessible from the computer that is running the Proteome Discoverer application. For details, see [“Directing the Proteome Discoverer Application to the Mascot Server Location”](#) on page 26.
2. If your Mascot server is running, verify that it is operating properly by submitting a simple search from the Mascot web interface. Do one of the following:
 - If the search from the Mascot web interface is successful, go to [step 3](#).
 - If the search fails, contact your system administrator. There might be a problem with the Mascot server itself.
3. If your Mascot server is operating properly and you can access it from the Proteome Discoverer application, try to perform a very simple search using the Mascot node. Do one of the following as applicable:
 - If simple searching fails, there might be a general problem in the interaction between the Proteome Discoverer application and the Mascot server. In this case, file an error report.
 - If you can perform simple Mascot searches, investigate your failing searches more closely:

Does the search finish successfully on the Mascot server according to the Mascot search log?

Do the process messages sent to the job queue during the search indicate the problem?
4. If these suggestions still do not resolve the search problem, file an error report.

Configuring Parallel Job Execution Parameters

To avoid overloading your machine, you can specify how many processing and consensus workflows to execute at the same time. (For definitions of processing and consensus workflows, see “[Workflow Editor](#)” on [page 34](#).) You must restart the application after you set the appropriate configuration parameters.

IMPORTANT Use caution when setting these parameters. Workflows share resources other than processors, such as hard disks or network connections, so running many jobs in parallel could cause Proteome Discoverer, other applications, or the operating system to become unresponsive.

❖ To specify the number of workflows to execute in parallel

1. Choose **Administration > Configuration**.
2. In the Configuration area in the lower left corner of the Proteome Discoverer window, click the + symbol to the left of Server Settings.
3. Click **Parallel Job Execution**.
4. In the Max. Number of Processing Workflows in Parallel Execution box, specify the number of processing workflows to execute at the same time.

Default: 1

Range: 1 to half the number of dual-processor cores on your machine.


For example, if your machine has a dual-core processor with hyperthreading, it has four logical cores. The maximum value for the parallel execution of processing workflows would therefore be half of four, or two.

5. In the Max. Number of Consensus Workflows in Parallel Execution box, specify the number of consensus workflows to execute at the same time.

Default: 1


Range: 1 to half the number of dual-processor cores on your machine

For example, if your machine has a dual-core processor with hyperthreading, it has four logical cores. The maximum value for the parallel execution of consensus workflows would therefore be half of four, or two.

6. Click the **Apply** icon,  **Apply** .

The Settings Applied confirmation box opens, as shown in [Figure 18](#).

Figure 18. Settings Applied confirmation box



7. Save any appropriate files, and restart the Proteome Discoverer application.

Studies in the Proteome Discoverer Application

The first step in conducting a search in the Proteome Discoverer application is to create a study. A study is primarily a container that you use to structure and annotate your data with meta information, such as how you treated the samples, how you performed the chromatography, or how you set up the data acquisition. A study also contains the input files, the information about the samples contained in those files, and their treatment information. It also keeps track of your processed results.

You can use a study to organize all the quality control files that you run regularly; to track all input and result files when you test different sample preparations, chromatography, or acquisition settings; or to organize true biological experiments with technical and biological replicates, different treatments, time points, and so forth.

To create a study, see [“Using Studies”](#) on [page 36](#).

Analyses in the Proteome Discoverer Application

An analysis in the Proteome Discoverer application is the actual data processing in the context of a study. It consists of the different data processing and reporting steps to perform, their associated workflows, their parameter settings, their input files, and other information related to how your data is being processed. A study can contain multiple analyses.

An analysis consists of multiple data-processing steps in the form of a workflow. The processing step and the consensus step are described in the next topic, [“Workflow Editor.”](#)

To create an analysis, see [“Using Analyses”](#) on [page 70](#).

Once you have set up an analysis, you can save it as a .pdAnalysis file, which you can then use as a template to set up new analyses.

Workflow Editor

The Workflow Editor is both a stand-alone window in the Proteome Discoverer application and a window included in a study. It is a flexible, complex tool that you can use to create customized searches and a customized results report. Use it to create one workflow to perform the search and another workflow to produce the results. Each workflow consists of a string of nodes that you can choose from a group of nodes unique to that workflow. You can further customize each workflow by setting parameters for each node.

In the Workflow Editor, you create the following two workflows:


- A processing workflow, which processes data in one or more RAW data files, generates the primary search results (for example, the results of a sequence database search by an engine such as Sequest HT) or extracting the raw quantification values (for example, the extracted reporter peak intensities where samples use isobaric tags). The processing workflow accepts spectrum container files such as RAW and MZML files as input and produces MSF files as output.
- A consensus workflow, which collects and assembles the data from one or more processing workflow results. The consensus workflow accepts an MSF file or files as input and produces a .pdResult file as output. The application can create different .pdResult files from the same MSF file and can combine different MSF files into one .pdResult file called a multiconsensus report.

You can create a reusable processing and consensus workflow template by saving your design to load and use at another time.

As a stand-alone window, the Workflow Editor is intended for use when you just want to edit or create a workflow template for later reuse.

[Figure 19](#) shows the stand-alone Workflow Editor.

Figure 19. Workflow Editor



The Workflow Editor consists of the following features:

- Workflow Nodes pane—Contains the nodes available for a processing workflow on one page, Processing Workflow, and the nodes available for a consensus workflow on another page, Consensus Workflow.
- Parameters pane—Displays the parameters available for the selected node.
- Workflow Tree pane—Where you connect nodes together to create a processing or consensus workflow.
- Post-Processing Nodes pane—Where you place the Post-Processing nodes, which process the results generated by the entire workflow. This pane is only visible when you create a consensus workflow.

For detailed instructions on using the Workflow Editor to create a workflow, see [“Using the Workflow Editor to Create Workflows”](#) on page 103.

Performing a Search

Follow the steps outlined in this topic to perform a search.

- [Before Performing a Search](#)
- [Using Studies](#)
- [Using Analyses](#)
- [Specifying Quantification Ratios from Selected Sample Groups](#)
- [Performing the Search](#)
- [Working with the Search Results](#)
- [Using the Workflow Editor to Create Workflows](#)
- [Using Workflow Templates](#)
- [Creating Specific Types of Workflows](#)
- [Searching Multiple Sequence Databases with Mascot](#)
- [Creating a Multiconsensus Report](#)

Before Performing a Search

Follow these steps before using the Workflow Editor either as a stand-alone tool or in the context of a study to create a workflow:

- Configure the search engine that you will use. For information on this procedure, see [“Configuring the Sequest HT Search Engine”](#) on [page 22](#) or [“Configuring the Mascot Search Engine”](#) on [page 25](#).
- If you are using the Sequest HT search engine, download a FASTA file or files, if necessary. See [“Adding FASTA Files to the Proteome Discoverer Application”](#) on [page 174](#).

If you are using the Mascot search engine, you must ask your administrator to download a FASTA file or files to the Mascot server.

- Make spectrum source files available as RAW, MGF, MZDATA, MZXML, or MZML files.

Using Studies

This topic explains how to create a study and how to use the Proteome Discoverer application to perform a database search. For illustrative purposes, it uses an example study called Bailey_2014, which uses a publicly available data set from the Chorus Project (<https://chorusproject.org>)¹. This project resides under the Elution Order Algorithm project and includes data about the following.

- Two sets of biological replicates


Four mice (replicates) were sacrificed and dissected. Individual organs of interest from them were homogenized, and the proteins were extracted from them and labeled. Then the differentially labeled organ-specific proteomes were mixed together if they came from the same mouse. Each pooled mouse sample was then run twice, using different acquisition method parameters.

For information on replicates, see “[Technical and Biological Replicates](#)” on [page 499](#).

- The TMT 8plex quantification method
- Label switching
- One biological factor: different tissues
- One technical factor: different acquisition methods

[Figure 20](#) shows how different tissue samples are distributed over four biological mouse replicates.

Figure 20. Distribution of different tissue samples over four biological mouse replicates



See the following topics:

- [Creating a Study](#)
- [Opening an Existing Study](#)
- [Adding a Quantification Method to the Study](#)
- [Adding the Study Factors](#)
- [Adding Input Files](#)
- [Importing MSF and .pdResult Files](#)

¹ Bailey, D.J.; McDevitt, M.T.; Westphall, M.S.; Pagliarini, D. J; Coon, J. J. Intelligent data acquisition blends targeted and discovery methods. *Journal of Proteome Research*, **2014**, 13 (4): 2152–2161.

- [Specifying the Quantification Method for Multiple Input Files](#)
- [Setting the Factor Values for the Samples](#)
- [Setting Values for Multiple Samples at the Same Time on the Input Files Page](#)
- [Setting Values for Multiple Samples at the Same Time on the Samples Page](#)
- [Saving a Study](#)
- [Copying a Study to Another Computer](#)
- [Study Window Parameters](#)

Creating a Study

Note If you want to use a custom quantification method, create it before you create or open a study. For instructions, see [“Adding a Quantification Method to the Study”](#) on page 41.

The first step is to create a study or open an existing study. To open an existing study, see [“Opening an Existing Study”](#) on page 40.

❖ To create a study

1. On the Start Page, click **New Study/Analysis**.

–or–


Choose **File > New Study/Analysis**.

–or–

Click the **Create New Study/Analysis** icon, .

The New Study and Analysis dialog box opens, as shown in [Figure 21](#).

Figure 21. New Study and Analysis dialog box



In this example, you only specify the name of the study and a root directory to save the study in.


2. In the Study Name box, specify the mandatory study name.

The application generates a default study name by searching for the common part of the file names when you add multiple files at once and using this common part as the suggestion for the name of the new study.

3. In the Study Root Directory box, specify the folder where you will store the study folder. Click the **Browse** button (...), and in the Select Folder dialog box, specify the folder and click **Select Folder**.
4. Click **OK**.

The application creates a new study folder in the folder that you specified as the root directory and opens a new page with the study name (Study: Bailey_2014 in the example), as shown in [Figure 22](#). It appends a .pdStudy extension to the study file name.

Figure 22. New Study: *Study_name* page



On the Study Definition page, you add a description of your study, select the quantification method or methods to use in the study, and set up the new factors to use to describe and distinguish your samples.

Opening an Existing Study


❖ To open an existing study

In the Recent Studies area of the Start Page, click the study name.

—or—

In the Recent Studies area of the Start Page, right-click the study name and choose **Open**.

—or—

Choose **File > Open Study** or click the **Open Existing Study** icon, , and browse to the study folder and then select and open the desired .pdStudy file.

—or—

1. In the Start area of the Start Page, click **Open Study**.
2. In the Open Study dialog box, browse to and select the name of the study, and click **Open**.

Adding a Description

You can optionally add a description of the study by typing it in the Study Description area of the Study Definition page.

Adding a Quantification Method to the Study


A quantification method contains the specification of the available quantification channels. The Proteome Discoverer application currently supports precursor ion-based quantification methods and MS/MS reporter ion-based quantification methods. It also supports precursor ion area detection. You can specify a quantification method for each of the input files.

The example used in this topic uses a custom TMT 8plex method, which you must create.

❖ To create a quantification method


1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon, .

The Quantification Methods view opens. It lists all of the available methods for both precursor ion and reporter ion quantification.

2. Click the **Add** icon,  **Add**.

The Create New Quantification Method dialog box appears, as shown in [Figure 23](#).


Figure 23. Create New Quantification Method dialog box



3. Select the **From Existing Method** option, and select **TMT 10plex** from the adjacent list.
4. Click **Create**.

The Quantification Method Editor opens, as shown in [Figure 24](#).

Figure 24. Quantification Method Editor



- In the Active column of the Quantification Method Editor, clear the check boxes for the following two channels:
 - 128N
 - 130N
- (Optional) Enter the isotopic impurities from the certificate of analysis for TMT reagent labels used for experiment. For instructions, see [“Correcting Reporter Ion Quantification Results for Isotopic Impurities”](#) on page 419.
- Click **OK**.
- In the Save Quantification Method dialog box, type the name of the quantification method that you want to create: **TMT 8plex**.
- Click **Save**.

The application adds the TMT 8plex method to the Quantification Methods view, as shown in [Figure 25](#).

Figure 25. New custom TMT 8plex method in the Quantification Methods view

Status	Method Name	Description	Is Active
▶	Dimethylation 3plex (C2H6, C2H2D4, 13C2D6)	Dimethylation 3plex (C2H4, C2D4, 13C2D4) Method	<input checked="" type="checkbox"/>
▶	Full 18O Labeling (O2 18O2)	18O labeling method for fully labeled samples	<input checked="" type="checkbox"/>
▶	Incomplete 18O Labeling (O2 O18O + 18O2)	18O labeling method for incompletely labeled samples	<input checked="" type="checkbox"/>
▶	iodo TMT6plex	Method for cysteine-reactive 6-plex Tandem Mass Tag® of ProteomeS...	<input checked="" type="checkbox"/>
▶	ITRAQ 4plex	Method for ITRAQ™ 4-plex mass tags by Applied Biosystems	<input checked="" type="checkbox"/>
▶	ITRAQ 8plex	Method for ITRAQ™ 8-plex mass tags by Applied Biosystems	<input checked="" type="checkbox"/>
▶	Low Resolution iodo TMT 6plex	Method for low resolution cysteine-reactive 6-plex Tandem Mass Tag...	<input checked="" type="checkbox"/>
▶	Low Resolution TMTe 6plex	Method for low resolution 6-plex Tandem Mass Tag® of Proteome Scie...	<input checked="" type="checkbox"/>
▶	SILAC 2plex (Arg10, Lys6)	SILAC 2plex (Arg10, Lys6) Method	<input checked="" type="checkbox"/>
▶	SILAC 2plex (Arg10, Lys8)	SILAC 2plex (Arg10, Lys8) Method	<input checked="" type="checkbox"/>
▶	SILAC 2plex (Ile6)	SILAC 2plex (Ile6) Method	<input checked="" type="checkbox"/>
▶	SILAC 3plex (Arg6, Lys4 Arg10, Lys8)	SILAC 3plex (Arg6, Lys4 Arg10, Lys8) Method	<input checked="" type="checkbox"/>
▶	SILAC 3plex (Arg6, Lys6 Arg10, Lys8)	SILAC 3plex (Arg6, Lys6 Arg10, Lys8) Method	<input checked="" type="checkbox"/>
▶	TMT 10plex	Method for 10-plex Tandem Mass Tag® of Proteome Sciences plc	<input checked="" type="checkbox"/>
▶	TMT 10plex Corr		<input checked="" type="checkbox"/>
▶	TMT 2plex	Method for 2-plex Tandem Mass Tag® of Proteome Sciences plc	<input checked="" type="checkbox"/>
▶	TMT 8plex	Method for 10-plex Tandem Mass Tag® of Proteome Sciences plc	<input checked="" type="checkbox"/>
▶	TMT 6plex	Method for 6-plex Tandem Mass Tag® of Proteome Sciences plc	<input checked="" type="checkbox"/>

New custom TMT
8plex quantification
method

For more information on adding a quantification method, see “” on page 419.

❖ **To select the quantification method to use in the study**

1. Click the **Study Definition** tab in the study, if it is not already selected.


The Quantification Methods area of the Study Definition page lists all the quantification methods that are currently available.

2. Select the check box corresponding to the quantification method or methods that you want to use.

In this example, the samples are labeled with TMT 8plex, so you would select the TMT 8plex check box, as shown in Figure 26.

If the Quantification Methods pane does not include the TMT 8plex method, choose **File > Save All**, and close and reopen the study.

Figure 26. Quantification method selected



3. After you add the input files, specify the quantification method for each input file (see “Adding Input Files” on page 48). For instructions, see “Specifying the Quantification Method for Multiple Input Files” on page 53.

Adding the Study Factors

In this step, you add the study factors that you want to use for your samples.

A factor is a single biological or technical parameter that you control, for example, genotype, diet, environmental stimulus, age, column length, spray voltage, or collision energy.

In theory, you can track every parameter as a factor in your study, but normally you only track the parameters that actually differentiate samples from each other. For example, the gradient that you use for the chromatography is a parameter in your experiment. However, if you use the same gradient for all your samples, you would not add this as an explicit factor to your study.


The experiment in the example dataset was performed by using eight different tissues and two different acquisition methods, so you will add two categorical factors, “Acquisition” and “Tissue,” to the study.

❖ To add categorical study factors

1. Add the first categorical (non-numeric) factor (Acquisition in the example):
 - a. In the Study Factors area of the Study Definition page, choose **Add > Categorical Factor**.

The categorical factor dialog box shown in Figure 27 appears.

Figure 27. Categorical factor dialog box



[new factor] is highlighted.

- b. Type a name over [new factor] for the new categorical factor, for example, **Acquisition**. See [Figure 28](#).

Note If the full categorical box becomes compressed, click **Edit** to restore it to its original size.

- c. In the box to the left of the Add and Delete buttons, type the name of the first value and click **Add**. See [Figure 28](#).


In this example, the value is the acquisition method, and the first acquisition method is intelligent data acquisition (IDA).

- d. In the same box, type the name of the second acquisition method and click **Add**.

In this example, the second acquisition method is data-dependent acquisition (DDA). See [Figure 28](#).

[Figure 28](#) shows the completed categorical factor dialog box.


Figure 28. Completed categorical factor dialog box



- e. Click **Apply** in the categorical factor dialog box.

Confirm that the Study Definition page resembles [Figure 29](#).

Figure 29. First categorical factor added to the study



2. Add the second categorical factor (Tissue in the example):
 - a. In the Study Factors area of the Study Definition page, choose **Add > Categorical Factor**.
The dialog box shown in [Figure 27](#) on [page 45](#) appears.
[new factor] is highlighted.
 - b. Type a name over [new factor] for the new factor (see [Figure 28](#) on [page 45](#)), for example, Tissue.
 - c. In the box to the left of the Add and Delete buttons (see [Figure 28](#) on [page 45](#)), type the name of the first type of tissue and click **Add**.
In this example, the first type of tissue is Cerebellum.
 - d. In the same box, type the name of any additional types of tissue and click **Add** after each one.


The example adds the following types of tissue to the study:

- Cerebellum
- Cerebrum
- Heart
- Kidney
- Liver

- Lung
 - Muscle
 - Spleen
- e. Click **Apply** in the Tissue dialog box.

Confirm that the Study Definition page resembles [Figure 30](#).

Figure 30. Second categorical factor added to the study




❖ **To add numeric study factors**

Note The example used throughout this topic does not use numeric study factors.

1. In the Study Factors area of the Study Definition page, choose **Add > Numeric Factor**.
The numeric factor dialog box shown in [Figure 31](#) appears.


Figure 31. Numeric factor dialog box



2. Backspace over [new factor] and type a name for the new numeric factor, for example, **time point**.
3. In the Factor Unit box, type the factor unit, for example, **hour**.
4. In the box to the left of the Add and Delete buttons, type the name of the first value, for example, **1**, and click **Add**.
5. In the same box, type the name of other values, if any, and click **Add**.
6. Click **Apply** in the numeric factor dialog box.

Figure 32 shows the completed numeric factor dialog box.

Figure 32. Completed numeric factor box




Adding Input Files

Add the input files from the example data set to your study. (You can add input files to the study at any point.) For the types of input files supported, see “Inputs” on page 11. You can add individual input files, multiple unrelated input files, or fractions.

Adding a Single Input File or Multiple Unrelated Input Files

❖ To add a single input file or multiple unrelated input files

1. Click the **Add Files** icon,  .
2. In the Add Files dialog box, browse to the location of the input files, and select them.

Note You can also add existing MSF files to a study. The application adds the MSF file as a new result to the Analysis Results page and adds the raw data files that were processed to create the MSF files to the Input Files page.


3. Click **Open**.

–or–

Drag the input file or files from Windows Explorer and drop them onto the Input Files page.

The input files appear on the Input Files page, as shown in [Figure 33](#). Each file on the page receives a unique identifier: F1, F2, ..., F*n*. The Proteome Discoverer application adds each file as a single study file.

Figure 33. Input files on the Input Files page




ID	Name	File Type	Quan Method	Sample Information
F1	29May3013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F2	31May3013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F3	31May3013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F4	31May3013_DJB_mouse_tmt8_BR3_unfrac_165min_mae15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F5	31May3013_DJB_mouse_tmt8_BR4_unfrac_165min_mae15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F6	29May3013_DJB_mouse_tmt8_BR1_unfrac_165min_dda15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F7	29May3013_DJB_mouse_tmt8_BR2_unfrac_165min_dda15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]
F8	29May3013_DJB_mouse_tmt8_BR3_unfrac_165min_dda15_1	.raw		Sample Type: [Sample], Acquisition: [n/a], Tissue: [n/a]

Adding Fractions


Fractions are raw data files derived from a sample that was separated into multiple vials, with each vial subjected to different treatment. For example, a sample could be divided into eight vials, with a different pH level applied to each vial.

One type of fraction is a MudPIT file. MudPIT experiments investigate complex proteomes by applying multidimensional chromatography to the samples before acquisition in the mass spectrometer. Typically, this process results in several dozen or even a few hundred fractions that are separately analyzed by LC/MS, resulting in one raw data file per sample fraction.

Analyzing gel slices or performing in-depth follow-up acquisitions also results in multiple fractions. Because all these fractions belong to the same sample, the Proteome Discoverer application can process all raw data files from these fractions as one contiguous input file, and it generates only one result file.

You can search MudPIT fractions in the Proteome Discoverer application by using the Add Fractions icon,  **Add Fractions**, or in the Proteome Discoverer Daemon utility by using the MudPIT button. However, Thermo Fisher Scientific recommends that you use the Proteome Discoverer application to process local batch or MudPIT samples, and use the Proteome Discoverer Daemon utility only when you need remote access and automatic processing from the command line—for example, when you use the Xcalibur data system to send data to a remote Proteome Discoverer computer (host) for acquisition. For information on processing MudPIT and batch samples with the Proteome Discoverer Daemon utility, see [“Using the Proteome Discoverer Daemon Utility”](#) on page 143.

❖ **To add fractions**


1. Click the **Add Fractions** icon,  **Add Fractions**.
2. In the Add Fractions dialog box, browse to the location of the input files, and select them.
3. Click **Open**.

The Proteome Discoverer application adds all selected files as a single study file.

Creating Subsets of Fractions


You can optionally use only parts of a file set in a fractionated sample by re-adding them as a subset on the Input Files page.

❖ **To create subsets of fractions**


1. Add the original input files as follows:
 - a. On the Input Files page, click the **Add Fractions** icon,  **Add Fractions**.
 - b. In the Add Fractions dialog box, browse to the directory where the fraction files are located, select the fraction files, and click **Open**.

The fraction files appear as a file set called F29 on the Input Files page, as shown in [Figure 34](#).

Figure 34. Fraction files added to the Input Files page as file set 29




2. Re-add some of the fraction files as follows:

- a. Click the **Add Fractions** icon,  **Add Fractions**.
- b. In the Add Fractions dialog box, browse to the same directory where the original fraction files are located, select the fraction files that you want to re-add, and click **Open**.

The Proteome Discoverer application adds the selected fraction files to the original file set as a subset. [Figure 35](#) shows three of the original fraction files added to the file set as a subset of F29 called F30.

Figure 35. Three fraction files re-added as a subset of F29



You can create any number of subsets.

When adding fractions to a study, the application first checks to see whether the study already contains any new files whose location changed within the study. If it finds these kinds of files, it updates the file locations but *does not* create new sample files. For example, suppose that the study contains a file called *drive:\xxx\file1*. You add file1 through file5. The application updates the path of the file1 file but does not create a new sample file.

If all selected files are new—that is, if the study does not already contain the files—the application creates a new sample file set that contains all the selected files. For example, suppose that you want to add file1 through file5, and all these files are new. The application creates a new sample file containing files file1 through file5.

If the study contains at least one new file, the application checks to see how many existing study files contain new files that are not subsets. For example, suppose that the existing sample file consists of all new files. You add file1 through file5. The study contains a sample file set called F1 consisting of file1 through file5 and no other files. In this case, the application does nothing.

As another example, suppose that an existing sample file set called F1 contains all new files but also some others, for instance, file1 through file8. In this case, the application creates a new subset containing file1 through file5 and names it [Subset of F1] F2... .

As a another example, suppose that an existing file set called F1 contains only new files, such as file1 through file3. The application creates a new file set, F2, which consists of file1 through file5 and is actually a superset of F1. But it adds F2 as a subset to the original F1 file set. The F2 subset added to F1 is called [Subset of F2] F1... .

As a final example, suppose that an existing file set called F1 contains only new files, such as file1 through file3. You want the application to create a new file set, F2, which consists of file1 through file3 and file6, as a superset of F1. But because file4 and file5 are missing, the application issues an error message and takes no action.

The application also does not support multiple file sets containing new files. When it encounters this situation, it issues an error message and takes no action.

Importing MSF and .pdResult Files

You can add existing MSF and .pdResult files to the study.

When you add an MSF file, the application adds the input files of the MSF file as fractions. If it cannot, you cannot add the MSF file to the study.

When you add a .pdResult file, the application adds the input files from each MSF file of the .pdResult file as one fraction. For example, suppose that a .pdResult file was created from two MSF files, msf1 and msf2. Msf1 was created from file1 through file4, and msf2 was created from file5 through file8. The application adds file1 through file4 as one file set and file5 through file8 as a second file set. If the application cannot do this, you cannot add the .pdResult file to the study.


Specifying the Quantification Method for Multiple Input Files

In this step, you specify the quantification method that was used for each of the files. In the example data set, all samples were labeled with the TMT 8plex custom quantification method.

❖ To set the quantification method for each of the input files

1. Click the **Input Files** tab, if it is not already selected.
2. In each sample row, click the Quan Method column and select the quantification method (for this example, **TMT 8plex**) from the list. See [Figure 36](#).

Figure 36. Setting the quantification method for each raw data file



Setting the Factor Values for the Samples

When you select a quantification method for a file, the Proteome Discoverer application generates a sample placeholder for each quantification channel.

The term “sample” refers to a distinct measurement of an analyte. At a minimum, a RAW data file contains at least a single sample, or in the case of sample multiplexing, it uses isobaric or metabolic labeling for multiple samples. If you measure the analyte from the same vial in two acquisitions—that is, if you measure two technical replicates—the application considers them two different samples. If you split the analyte, label each with a different isobaric tag, mix them, and acquire them into a single RAW data file, the application considers the two differently labeled parts as two different samples.

Each sample is associated with a sample type. Currently only quantification uses sample types. The sample type can be one of the following:

- Sample: A specimen from a larger biological entity
- Control: A sample typically used as a reference sample in a quantification experiment
- Blank: A sample consisting only of solvent and no sample mixture
- Standard: A sample consisting of a standard quality-control peptide mixture

You can create ad hoc relative quantitative ratios between any study factors or sample types, for example, blank/standard, control/sample, blank/control, and standard/blank. You can leave all of the channels as samples or make all of the channels any combination of sample, control, blank, or standard.

The Sample Type column on the Samples page of the study displays the sample type of each sample. The default sample type is Sample. If a file has samples for different quantification channels, mark one of the samples Control. In this example, the channel that was used to label the mouse kidney tissues is used as the control for generating the quantification ratios.

You can also use the control channels in scaling. For more information, see “[Normalizing Peptide Groups and Protein Abundances](#)” on page 443.

Each sample is associated with a quantification channel shown in the Quan Channel column and with values for each of the factors that you specified for your study. Previously, you specified a factor for the acquisition method used and a factor for the tissue that was extracted and labeled. You now set the correct factor values for each of the samples in the study.


Each sample has an automatically generated sample name composed of the raw data file name and the appended name of the quantification channel. You can change this name, but the name must be unique among all samples in the study.

❖ To view the samples

1. On the Input Files page, click the gray arrow to the left of a sample to display its constituent file entries, as shown in [Figure 37](#).

A hierarchical view opens, showing the samples contained in a raw data file. For each of the raw data files in the example, there are eight samples for the eight quantification channels of the TMT 8plex method. See [Figure 37](#).

Figure 37. Samples in a raw data file



❖ To set the factor values for the samples


1. Click the **Input Files** tab, if it is not already selected.
2. Click the gray arrow next to the first sample to expand the information about the sample.
3. For the first factor (in the example, Acquisition), set the value for each sample in each raw data file by selecting the down arrow in the factor column and then selecting the value from the list.

In the example, select **IDA** in the Acquisition column, as shown in [Figure 38](#).

- For the second factor (in the example, Tissue), set the value for each sample in each raw data file by selecting the down arrow in the factor column and then selecting the value from the list.

Figure 38 shows this process for the Tissue factor.

Figure 38. Setting the factor values for the Tissue factor



- Set the same values in the Acquisition and the Tissue columns for the rest of the samples.

After you finish setting the factor values for each sample, the Input Files page resembles Figure 39. (In the example data set, you must set 128 factor values for eight files with eight samples each with two factors each.)

Note For instructions on changing the values for factors or other study variables for multiple samples at once, see “Setting Values for Multiple Samples at the Same Time on the Input Files Page” on page 58.

Figure 39. Completed sample factors

ID	Name	File Type	Quan Method	Sample Information		
F4	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [IDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]		
Sample	Sample Identifier	Sample Type	Quan Channel	Control Chann	Acquisition	Tissue
4	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Control	126	-	IDA	Kidney
30	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	127_N	126	IDA	Cerebellum
31	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	127_C	126	IDA	Muscle
32	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	128_C	126	IDA	Cerebrum
33	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	129_N	126	IDA	Lung
34	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	129_C	126	IDA	Liver
35	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	130_C	126	IDA	Heart
36	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	Sample	131	126	IDA	Spleen
File:						
ID	Name	Date Modified	Size			
F4.1	WUSSJO-BGIBS-PC\c\Program Files\Proteome Discoverer source files\Studies\mouse_tmt8_data\29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1.raw	5/29/2014 9:26:36 PM	328754564 [Byte]			
F5	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [DDA] Tissue: [n/a, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]		
Sample	Sample Identifier	Sample Type	Quan Channel	Control Chann	Acquisition	Tissue
5	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Control	126	-	DDA	n/a
37	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	127_N	126	DDA	Cerebellum
38	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	127_C	126	DDA	Muscle
39	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	128_C	126	DDA	Cerebrum
40	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	129_N	126	DDA	Lung
41	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	129_C	126	DDA	Liver
42	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	130_C	126	DDA	Heart
43	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	Sample	131	126	DDA	Spleen
File:						
ID	Name	Date Modified	Size			
F5.1	WUSSJO-BGIBS-PC\c\Program Files\Proteome Discoverer source files\Studies\mouse_tmt8_data\31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1.raw	5/29/2014 9:31:46 PM	299616516 [Byte]			
F6	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [DDA, IDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]		
Sample	Sample Identifier	Sample Type	Quan Channel	Control Chann	Acquisition	Tissue
6	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	Control	126	-	DDA	Kidney
44	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	Sample	127_N	126	DDA	Cerebellum
45	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	Sample	127_C	126	DDA	Muscle
46	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	Sample	128_C	126	DDA	Cerebrum
47	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	Sample	129_N	126	DDA	Lung

6. Click the red down arrow next to each sample to compress the information displayed.

The Sample Information column summarizes the information about the samples contained in a file, as shown in Figure 40.

Figure 40. Sample Information Column

ID	Name	File Type	Quan Method	Sample Information
F1	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [IDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F2	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [IDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F3	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [IDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F4	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [IDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F5	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [DDA] Tissue: [n/a, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F6	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [DDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F7	31May2013_DJB_mouse_tmt8_BR3_unfrac_165min_mae15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [DDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]
F8	31May2013_DJB_mouse_tmt8_BR4_unfrac_165min_mae15_1	.raw	TMT 8plex	Sample Type: [Control, Sample] Acquisition: [DDA] Tissue: [Kidney, Cerebellum, Muscle, Cerebrum, Lung, Liver, Heart, Spleen]

7. (Optional) Click the **Samples** tab.

Figure 41 shows the Samples page. It displays the same sample information as the Input files page.

Figure 41. Samples page

Sample	Sample Identifier	Sample Type	Acquisition	Tissue
S1	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Control	IDA	Kidney
S2	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Control	IDA	Kidney
S3	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Control	IDA	Kidney
S4	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Control	IDA	Kidney
S5	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae1	Control	DDA	n/a
S6	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae1	Control	DDA	Kidney
S7	31May2013_DJB_mouse_tmt8_BR3_unfrac_165min_mae1	Control	DDA	Kidney
S8	31May2013_DJB_mouse_tmt8_BR4_unfrac_165min_mae1	Control	DDA	Kidney
S9	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Cerebellum
S10	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Muscle
S11	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Cerebrum
S12	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Lung
S13	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Liver
S14	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Heart
S15	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda1!	Sample	IDA	Spleen
S16	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Cerebellum
S17	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Muscle
S18	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Cerebrum
S19	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Lung
S20	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Liver
S21	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Heart
S22	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda1!	Sample	IDA	Spleen
S23	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Cerebellum
S24	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Muscle
S25	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Cerebrum
S26	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Lung
S27	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Liver
S28	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Heart
S29	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda1!	Sample	IDA	Spleen
S30	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Cerebellum
S31	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Muscle
S32	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Cerebrum
S33	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Lung
S34	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Liver
S35	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Heart
S36	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda1!	Sample	IDA	Spleen

Setting Values for Multiple Samples at the Same Time on the Input Files Page

Highly multiplexed data are results obtained from processing several samples from one raw data file mixed and analyzed together in one LC/MS analysis, where isotopic and isobaric labels were used in quantification to distinguish the contribution of the individual samples. If you have highly multiplexed data, several files and samples with many study variables to set (sample type, quantification channel, study factors), or both, you can set values for study variables for multiple samples at once. You can use either the mouse or the keyboard to set these values on the Input Files page.


You can also use multicell editing on the Samples page to enter factor values for several samples. For example, if you have acquired data from several different tissues and the tissue is encoded in the file and sample names, you can use the filter option on the Samples page to filter all samples that contain “liver” and then set the Tissue factor to “liver” for all these samples.

❖ To set values for multiple sample cells at the same time by using the mouse

1. Click the **Input Files** tab.
2. To select multiple samples for one study variable, select the first cell and drag the cursor to select the remaining cells.


Figure 42 shows multiple cells being selected in the Acquisition column.

Figure 42. Selecting multiple cells at one time




3. Press the F2 key to enter multicell editing mode.
4. Select the new value from the list in the last cell that you selected, as shown in [Figure 43](#).
5. Press the RETURN or ENTER key, or click elsewhere in the application.

Figure 43. Selecting a new value to apply to all cells in multicell editing mode



In the example, the values in the cells of the Acquisition column change from IDA to DDA, as shown in Figure 44.

Figure 44. Multiple values changed at one time




❖ **To set values for multiple sample cells at the same time by using the keyboard**

1. Click the **Input Files** tab.
2. Click the first cell to edit.
3. Press the ESC key twice until you are no longer in cell editing mode.
4. Holding down the SHIFT key, use the up and down arrow keys on your keyboard to select the cells that you want to edit.
5. Press the F2 key to enter multicell editing mode.
6. In the last cell selected, use the up and down arrow keys to select a new cell value from the list.
7. Press the RETURN or ENTER key to apply your change to all selected cells.

Setting Values for Multiple Samples at the Same Time on the Samples Page

Use the following procedure to change values for multiple samples on the Samples page.


❖ **To set values for multiple samples at the same time on the Samples page**

1. Click the **Samples** tab.
2. Click the **Contains** icon, , if necessary, in the Acquisition column.


3. In the box next to the icon, type the name of the filter.

For example, suppose you want to change the acquisition method for a number of samples from IDA to DDA. If you want to display only samples that contain “DDA,” you would type **dda** in the filter box. Figure 45 shows the results.

Figure 45. Multicell editing of samples using filters



4. Select the first cell to change, and then drag the cursor to select the remaining cells that you want to change.
5. Press the F2 key to enter multicell editing mode.
6. Select the new value from the list in the last cell that you selected.
7. Press the RETURN or ENTER key, or click elsewhere in the application.


To return the samples to their unfiltered state, select the **Clear All Filters** icon, .


Saving a Study

You can save a study manually or automatically.


You can save a study manually at any time. A change in a study that requires you to save it is indicated with an asterisk (*) in the tab after the study name.

Note Studies and analyses in the Proteome Discoverer application are separate and must be saved separately. Saving a study does not save an analysis, and saving an analysis does not save a study.

You can set an option to have the application save studies automatically when you click the Run icon,  Run . This option also saves the results generated in the study. It saves the analysis containing the workflow, but when you close the study and reopen it, you must access the analysis by clicking the Analysis Results tab, and then doing one of the following:

Click the **Reprocess** icon,  Reprocess ▾ , and choose **All Analysis Steps** (to open both the processing and consensus workflow) or **Last Consensus Step** (to open just the consensus workflow).

–or–


Click the **Show Details** icon,  Show Details .

❖ To save a study manually

Choose **File > Save**.

The Proteome Discoverer application saves the study in the *study_name*.pdStudy file in the study directory.

❖ To save a study automatically

1. At some point before you click the Run icon,  Run , select **Tools > Options**.
2. In the Options dialog box, select **Study Options** in the left pane.
3. Select the **Auto Save When Starting Analysis** check box.
4. Click **OK**.

The Auto Save When Starting Analysis option does not take effect for any open studies. The application automatically saves only newly opened studies.

❖ To save all open studies

Choose **File > Save All**.

The Proteome Discoverer application saves all studies open in the application.

Copying a Study to Another Computer

If your computer has run out of storage space or no longer works, you might want to move your study data to a new computer to continue working with it. You can copy a study folder with all of its data to another computer and continue to use it on another installation of the Proteome Discoverer application.

Your computer system must meet the following prerequisites to effect this transfer:

- All file references in a study must not point to fixed locations on a specific computer. You must be able to update the storage paths to a new location.


The study management system automatically tries to locate files in a study folder or any of its subfolders. You can resolve all missing file issues in one of the following ways:

- Copy all missing input files or missing result files into the study folder or a subfolder. For example, within the study folder you can create a folder called Input Files and copy all raw data files into the Input Files folder.
 - Update the file location by re-adding the files to the study. If the study can no longer find a *drive:\data\abc.raw* file because you have moved the file to *drive:\data\project xy\abc.raw*, you can add the *abc.raw* file again from the new location. The study management system stores the new location for the file. It also retains all quantification methods, factors, and other study variables set for the file.
- The quantification methods stored in a study do not keep their connection to the initially connected processing server from which they were received.

When you go to the Study Definition page of an open study, you find all quantification methods currently available on the connected processing server listed in the Quantification Methods box in the upper right corner. By clicking the check box next to a listed quantification method, you can add this method to the study. When you save a study, the study management system stores the added quantification methods in the study file.

The quantification methods stored in the study are kept independent of the processing server. However, once you add a quantification method to a study, the method is independent of the connected processing server, and any changes you make to the quantification methods do not appear in the quantification methods already added to your study. You must be able to display and edit the quantification methods within the study itself. If you hold the mouse pointer over one of the check boxes representing a quantification method, two small icons for editing a method and exporting a method appear, as shown in [Figure 46](#).

Figure 46. Icons for editing and exporting a quantification method



Clicking the Edit icon opens the Quantification Method Editor.

Note If a quantification method is already set for one or more input files, you cannot change the number and names of the quantification channels. In this case, you can only change the modifications set for the method or the channels. Click the **Export** icon to save the quantification method to a METHOD file that you can then import into another Proteome Discoverer instance.

The quantification methods are now automatically set for the quantification nodes contained in the workflow to avoid any accidental mismatches between the quantification methods specified for an input file and the quantification methods set in a workflow. The analysis validation detects any inconsistencies.

You can also use the Proteome Discoverer Daemon utility to select a quantification method. For information on this procedure, see [“Running the Proteome Discoverer Daemon Utility from the Window”](#) on page 145.

Study Window Parameters

Table 2 describes the parameters available in pages of a study.

Table 2. Study page parameters (Sheet 1 of 6)





Parameter	Description
 Add Files	<p>Opens the Add Files dialog box so that you can select the files to submit to the processing workflow. You can add RAW, MSF, MGF, MZXML, MZDATA, MZML, and .pdResult files.</p> <p>This icon is available only when the Input Files page is selected.</p>
 Add Fractions	<p>Opens the Add Fraction dialog box so that you can select a series of input spectrum files in any format that are fractions. The formats can be RAW, MGF, MZXML, MZDATA, or MZML.</p> <p>This icon is available only when input files are selected on the Input Files page.</p>
 Remove Files	<p>Removes the selected input files from the Input Files page.</p> <p>This icon is available only when the Input Files page is selected.</p>
 Open Containing Folder	<p>Opens the folder or folders containing the .pdResult file and the MSF file corresponding to a selected analysis result on the Analysis Result page. This command opens the folder or folders in Windows Explorer.</p>

Table 2. Study page parameters (Sheet 2 of 6)



Parameter	Description
 New Analysis	Opens a new Analysis window.
 Open Analysis Template	Opens the Open Analysis Template dialog box so that you can select a custom .pdAnalysis file to use as a template.
Study Definition page	Displays basic study information, the quantification method, and study factors.
Study Summary	Displays information about the study.
Study Name	Displays the name of the study that you entered in the Study Name box of the Create New Study/Analysis dialog box.
Study Directory	Displays the directory where the study is located. You entered this location in the Study Root Directory box of the Create New Study/Analysis dialog box.
Last Changed	Displays the date and time that you saved the most recent change.
Creation Date	Displays the time and date that you clicked OK in the New Study and Analysis dialog box to create the study.
Study Description	Displays your description of the study.
Quantification Methods	Displays the quantification methods available for the selected input files in the study.
Study Factors	Displays information about the factors that distinguish two or more samples. Examples of factors are drug treatment; time of drug application; differences in tissue, organism, patient, and so forth; differences in sample preparation, chromatography settings, acquisition parameters; or differences in the isobaric or metabolic labels used.
Paste (Study Factors area)	Pastes the copied study factor box into the Study Factors area and assigns each copy the study factor name and an incrementing number: <i>study_factor1</i> , <i>study_factor2</i> , . . . <i>study_factor_n</i> . This command is only available after you select the Copy command.
Copy (Study Factors area)	Copies the selected study factor box.

Table 2. Study page parameters (Sheet 3 of 6)

Parameter	Description
Add (Study Factors area)	<p>Specifies the type of study factor to add to the study. It can be one of the following:</p> <ul style="list-style-type: none"> • Categorical Factor: Adds a new study factor box for a factor that belongs to a non-numerical class. See Figure 28 on page 45. • Numeric Factor: Adds a new study factor box for a factor that is numerical in nature. See Figure 32 on page 48.
Input Files page	Lists the input files submitted to the Proteome Discoverer application for the analysis.
ID	Displays the identification number that the Proteome Discoverer application assigns to the raw data file when you add the file to the study. File names begin with an F, and sample names begin with an S.
Name	Displays the name of the raw data file.
File Type	Displays the type of the input file, for example, a raw data file.
Quan Method	Displays the type of quantification performed on the samples. A quantification method specifies the available quantification channels. The Proteome Discoverer application currently supports precursor ion quantification methods using stable isotopes and MS/MS reporter ion quantification methods using isobaric tags. You can specify a quantification method for each of the input files.
Sample Information	Displays information about the samples in the input files, such as the sample type, the acquisition method, and the tissue type.
Samples page	Displays information about the samples contained in the input files.
Sample	Specifies the number of the sample.
Sample Identifier	Specifies the name of the sample.

Table 2. Study page parameters (Sheet 4 of 6)




Parameter	Description
Sample Type	<p>Displays the sample type associated with a sample for quantification purposes. The sample type can be one of the following:</p> <ul style="list-style-type: none"> • Sample: A specimen from a larger biological entity • Control: A sample typically used as a reference sample in a quantification experiment • Blank: A sample consisting only of solvent and no sample mixture • Standard: A sample consisting of a standard quality-control peptide mixture
<i>Study_variables</i>	<p>Display the factors that distinguish two or more samples, for example, drug treatment; time of drug application; differences in tissue, organism, or patient; differences in sample preparation, chromatography settings, or acquisition settings; or differences in isobaric or metabolic labels.</p>
Analysis Results page	<p>Displays the results of the search.</p>
 Open Result	<p>Opens the .pdResult file containing the search results for the selected file.</p>
 Show Details	<p>Opens an Analysis Sequence Details window, which displays the Workflow Editor and the Analysis window containing the processing and consensus workflows used to generate the selected results file.</p>
 Reprocess ▾	<p>Repeats the analysis if you change the input files or the workflow used to process the RAW data files. You can change the workflows before rerunning the processing or consensus workflows.</p> <ul style="list-style-type: none"> • All Analysis Steps: Re-performs all the analysis steps. • Last Consensus Step: Re-performs just the last step in the consensus workflow. • Use Results to Create (Multi) Consensus: Uses a new consensus workflow to reprocess MSF results associated with one or more .pdResult files.

Table 2. Study page parameters (Sheet 5 of 6)

Parameter	Description
Search	Specifies the string to search for in the search results listed on the Analysis Results page. Specify the string in the Search box. By default, the application searches only in the File Name and Description columns. If you select the column name from the Search For menu to the right of the Search box, the application searches for the string only in that column.
Search for	Specifies the name of the column to search for the string specified in the Search box.
Execution State	Displays the status of the search producing the results. These status states are the same as those displayed by the job queue. For information on these states, refer to the Help.
Creation Date	Displays the time and date that the Proteome Discoverer application began the analysis.
File Name	Displays the name of the input file being processed or completed. The application uses the name of the first sample as a default if you do not specify an analysis name.
Description	Displays your brief description of the analysis. This description has no relation to any other descriptions already introduced on any other level of the study.
Workflow page	For information on the parameters of the Workflow page, see Table 4 on page 118 .

Table 2. Study page parameters (Sheet 6 of 6)

Parameter	Description
Grouping & Quantification page	Displays information about how the application groups quantification ratios and samples.
Study Variables	Specifies the categorical study factors to group samples and quantification ratios by. This area contains check boxes for all the study factors that you selected in the Study Factors Area of the Study Definition page. <ul style="list-style-type: none">Files: Groups samples and quantification ratios by files.Quan Channels: Groups samples and quantification ratios by quantification channels.<i>Study_factors</i>: Groups samples and quantification ratios by user-defined study factors.Sample Types: Groups factors that include Sample, Control, Standard, and Blank. Variables displayed in italics contain only a single value.
Manual Ratio Generation	Contains menus where you choose the numerator and denominator for each quantification ratio. <ul style="list-style-type: none">Numerator: Specifies the numerator of the quantification ratio.Denominator: Specifies the denominator of the quantification ratio.
Bulk Ratio Generation	Displays the study factor values to use as the denominators for semiautomatically generated quantification ratios.
Generated Sample Groups	Displays samples grouped by the values set for the selected study factors. The application ignores any other differences among the samples that might be present.
Generated Ratios	Displays the quantification ratios generated.

Using Analyses


When you have finished setting up a study, you can create an analysis. See these topics:

- [Creating an Analysis](#)
- [Opening an Existing Analysis](#)
- [Adding Input Files to an Analysis](#)

- Using Multiple Processing Steps in an Analysis
- Adding or Deleting a Processing Step
- Creating the Workflows
- Saving an Analysis
- Analysis Window Parameters


Creating an Analysis

❖ To create an analysis

On the Study: *Study_name* page, click the **New Analysis** icon,  **New Analysis**.

An Analysis window opens on the right side of the Study: *Study_name* page, as shown in Figure 47.

Figure 47. Analysis window



Analysis window containing a Consensus Step box and a Processing Step box

An Analysis window contains the following items:

- A Consensus Step box, which represents the consensus workflow step of the data processing
- A Processing Step box, which represents the processing workflow step of the data processing
- A Child Steps bar, which contains an Add button that you can use to add another Processing Step box. Multiple Processing Step boxes are useful when you want to process the same data in different ways—for example, by using different nodes or different node settings.


In addition, two new tabs appear on the Study: *Study_name* page: the Workflows tab and the Grouping & Quantification tab, as shown in [Figure 47](#).

If you switch to an existing study that includes an Analysis window containing an existing consensus workflow, processing workflow, or both, and you want to open a new analysis, you must close the existing Analysis window and open a new Analysis window. To close the window, click **X** in the upper right corner.


Opening an Existing Analysis

You can open an existing analysis through a template or through a results file.

❖ To open an existing analysis through a template

1. In an open study, click the **Open Analysis Template** icon,  **Open Analysis Template**.
2. In the Open Analysis Template dialog box, browse to and select the .pdAnalysis file that you want to open.
3. Click **Open**.

❖ To open an existing analysis through a results file

1. In an open study, click the **Analysis Results** tab in the study.
2. Select the appropriate .pdResults file.
3. Click the **Reprocess** icon,  **Reprocess**.
4. To open an Analysis window containing both the processing and consensus workflows, choose **All Analysis Steps** in the Reprocess list.

–or–


To open an Analysis window in order to execute just the consensus workflow, choose **Last Consensus Step** in the Reprocess list.

Adding Input Files to an Analysis

❖ To add the input files to an analysis

Select and drag the files from the Input Files page to the Input Files box of the Processing Step box in the analysis window.

The input files are listed in the Input Files area of the Processing Step box, as shown in [Figure 48](#).

Figure 48. Input files in the Input Files area of the Processing Step box of the Analysis window

Note You can remove a file from the Input Files area of the Processing Step box by clicking the X to the left of the file name.

Using Multiple Processing Steps in an Analysis

An analysis can contain more than one processing step. You can use additional processing steps when you want to process one set of input files differently from another set of input files. For example, you might want to run a Percolator validation separately on a set of files, such as the replicates of the control group and the replicates of the treatment group.

❖ To add multiple processing steps

In an open analysis window, click **Add** on the Child Steps bar.

A new processing step without a workflow appears.

❖ To duplicate an existing processing step, including its workflow

In an open analysis window, click **Clone** on the processing step bar.

Another window with a child step appears. You can now define a different processing workflow for the same file or a different file.

Adding or Deleting a Processing Step

You might want to perform only the processing step in a workflow. For example, you only need the processing step if you want to export peak lists as an MGF file with the Spectrum Exporter node. You can remove the Consensus step and re-add it later.

❖ **To perform only the processing step in an analysis**

In the title bar of the Consensus Step box, click the **X** in the upper right corner to remove the consensus step from the analysis.

The Consensus Step box disappears, and the Add icon,  **Add**, appears in the Analysis title bar.

❖ **To re-add a deleted consensus step in an analysis**

In the Analysis title bar, click the **Add** icon, .

The Consensus Step box reappears, and the Add icon disappears.

Creating the Workflows


The next step in performing a search is to create the workflows to use for the processing and consensus steps. This topic gives a brief overview of the steps required to create processing and consensus workflows. For detailed information on creating workflows, see [“Using the Workflow Editor to Create Workflows”](#) on page 103. For information on the Workflow Editor, see [“Workflow Editor”](#) on page 34.

❖ **To create the processing workflow**

1. Click the **Workflows** tab.

Note The Workflows tab does not appear until you add or open an analysis.

The Workflow Editor opens.

2. Click the **Show Workflow** icon, , on the title bar of the Processing Step box to indicate that you want to create a processing workflow.

The Workflow Nodes pane lists the nodes available for use in the processing workflow. You might need to click the Workflow Nodes tab to see this pane.

3. Create the appropriate processing workflow in the Workflow Tree pane of the Workflow Editor. For detailed information on this procedure, see [“Creating a Processing Workflow”](#) on page 104.

Default processing workflows are available for several instruments. You can use these default workflows as is or modify them to suit your needs. For information on accessing them, see [“Using Common Workflow Templates”](#) on page 119. You must check the node parameters in the workflow to make sure they are appropriate for your data and add modifications to the search for the quantification to be successful.

The example shown in [Figure 50](#) on page 77 includes the following nodes:

- Spectrum Files node
- Event Detector node
- Spectrum Selector node

- Sequest HT node
- Percolator node
- Precursor Ions Area Detector node
- Reporter Ions Quantifier node

For information about these nodes, refer to the Help.

4. Connect the nodes together, as needed (see the next figure).

In this example, the only connections that you must make are from the Spectrum Selector node to the Sequest HT node to the Percolator node.

5. Set the appropriate parameters for each node as follows:


- a. Click the node.
- b. (Optional) Click **Show Advanced Parameters** in the Parameters pane to the left to display all parameters.
- c. Set the appropriate parameters.

6. For this example, set the parameters of the Sequest HT node as follows:

- Set the Protein Database parameter to an appropriate database, for example, **SwissProt**.
- Set the N-Terminal Modification parameter under Dynamic Modifications (Peptide Terminus), for example, **TMT 6plex/+229.163 Da**.
- Set a Dynamic Modification parameter, for example, **TMT6plex /+229.163 Da (K)**.
- Set a Static Modification parameter, for example, **Carbamidomethyl/+57.021 Da (C)**.

Figure 49 shows these parameter settings.

Figure 49. Example of parameters set for the Sequest HT node





Parameters	
Show Advanced Parameters	
1. Input Data	
Protein Database	Swissprot2.fasta
Enzyme Name	Trypsin (Full)
Max. Missed Cleavage Sites	2
Min. Peptide Length	6
Max. Peptide Length	144
2. Tolerances	
Precursor Mass Tolerance	10 ppm
Fragment Mass Tolerance	0.6 Da
Use Average Precursor Mass	False
Use Average Fragment Mass	False
3. Spectrum Matching	
Use Neutral Loss a Ions	True
Use Neutral Loss b Ions	True
Use Neutral Loss y Ions	True
Use Flanking Ions	True
Weight of a Ions	0
Weight of b Ions	1
Weight of c Ions	0
Weight of x Ions	0
Weight of y Ions	1
Weight of z Ions	0
4. Dynamic Modifications	
Max. Equal Modifications Per Peptide	3
1. Dynamic Modification	TMT6plex / +229.163 Da (K)
2. Dynamic Modification	None
3. Dynamic Modification	None
4. Dynamic Modification	None
5. Dynamic Modification	None
6. Dynamic Modification	None
5. Dynamic Modifications (peptide terminus)	
1. N-Terminal Modification	TMT6plex / +229.163 Da (N-Terminus)
2. N-Terminal Modification	None
3. N-Terminal Modification	None
1. C-Terminal Modification	None
2. C-Terminal Modification	None
3. C-Terminal Modification	None

Figure 50 shows a complete example of the processing workflow.


Figure 50. Processing workflow example

The screenshot shows the Proteome Discoverer software interface. The main window is titled 'Study: Bailey_2014 *'. The 'Workflow Nodes' pane on the left lists various processing steps, including 'Event Detector', 'Reporter Ions Quantifier', 'Spectrum Selector', 'Precursor Ions Area Detector', 'Sequest HT', and 'Percolator'. The 'Workflow Tree' pane in the center displays a flowchart of the workflow. The 'Analysis' pane on the right shows the details of the selected 'Processing Step', including the workflow name, result file, and a list of input files.


7. (Optional) Save the processing workflow:

- a. In the Workflow box above the Workflow Tree pane, type a name for the processing workflow.
- b. (Optional) In the Description box, type a brief description of the processing workflow.
- c. In the Workflow Editor, click the **Save** icon,  **Save**.
- d. In the Save Workflow dialog box, do the following:
 - i. Select the file to save the workflow in, or type a file name in the File Name box. You can save the workflow in the study folder or in the Common Templates folder (select the **Save Common** icon, , in this case), or in a separate folder of workflows.
 - ii. Click **Save**.

The application saves the workflow in the *file_name*.pdProcessingWF file.

Note A yellow triangle containing an exclamation mark in the upper right corner of the Processing Step box () usually indicates that a node in the workflow is obsolete, a parameter is missing, a required node is missing, or the output file name is invalid or missing. Hold the mouse pointer over the triangle to display further details about what is missing.

❖ **To create the consensus workflow**

1. Click the **Show Workflow** icon, , on the title bar of the Consensus Step box to indicate that you want to create a consensus workflow.

The Workflow Nodes pane lists the nodes available for use in the consensus workflow. You might need to click the Workflow Nodes tab to see this pane.

Create the appropriate consensus workflow in the Workflow Tree pane. For detailed information on creating a consensus workflow, see [“Creating a Consensus Workflow”](#) on [page 112](#).

Default consensus workflows are available for several instruments. You can use these default workflows as is or modify them to suit your needs. For information on accessing them, see [“Using Common Workflow Templates”](#) on [page 119](#). You must check the node parameters in the workflow to make sure they are appropriate for your data and add modifications to the search for the quantification to be successful.

The example workflow in [Figure 51](#) includes the following nodes in the Workflow Tree pane:

- MSF Files node
- PSM Grouper node
- Peptide Validator node
- Peptide and Protein Filter node
- Protein Scorer node
- Protein Grouping node
- Peptide and Protein Quantifier node

The example workflow in [Figure 51](#) includes the following nodes in the Post-Processing Nodes area:

- Result Statistics node
- Data Distributions node

For information about these nodes, refer to the Help.

2. Set the appropriate parameters for each node.


[Figure 51](#) shows an example of the consensus workflow.

Figure 51. Consensus workflow example


The screenshot displays the Proteome Discoverer software interface. The main window is titled 'Study: Bailey_2014'. The 'Workflow Editor' is active, showing a 'Workflow Tree' with the following nodes: MSF Files (0), PSM Grouper (1), Peptide Validator (2), Peptide and Protein Filter (3), Protein Scorer (4), Result Statistics (7), and Data Distributions (8). The 'Analysis' panel on the right shows the 'Consensus Step' configuration. The workflow is named 'Bailey_2014' and the result file is '29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda15_1.pdResult'. The 'Input Files' list includes:

File ID	File Name	Sample Type
F1	29May2013_DJB_mouse_tmt8_BR1_unfrac_165min_dda15_1	TMT 8plex Sample Type: [Control. S
F2	29May2013_DJB_mouse_tmt8_BR2_unfrac_165min_dda15_1	TMT 8plex Sample Type: [Control. S
F3	29May2013_DJB_mouse_tmt8_BR3_unfrac_165min_dda15_1	TMT 8plex Sample Type: [Control. S
F4	29May2013_DJB_mouse_tmt8_BR4_unfrac_165min_dda15_1	TMT 8plex Sample Type: [Control. S
F5	31May2013_DJB_mouse_tmt8_BR1_unfrac_165min_mae15_1	TMT 8plex Sample Type: [Control. S
F6	31May2013_DJB_mouse_tmt8_BR2_unfrac_165min_mae15_1	TMT 8plex Sample Type: [Control. S
F7	31May2013_DJB_mouse_tmt8_BR3_unfrac_165min_mae15_1	TMT 8plex Sample Type: [Control. S
F8	31May2013_DJB_mouse_tmt8_BR4_unfrac_165min_mae15_1	TMT 8plex Sample Type: [Control. S

3. (Optional) Save the consensus workflow:

- In the Workflow box, type a name for the consensus workflow.
- (Optional) In the Description box, type a brief description of the consensus workflow.
- In the Workflow Editor, click the **Save** icon,  **Save**.
- In the Save Workflow dialog box, do the following:
 - Browse to the study folder, and select the file to save the workflow in, or type a file name in the File Name box.
 - In the Save As Type box, select **Consensus Workflow File (*.pdConsensusWF)**.
 - Click **Save**.

The application saves the workflow in a *file_name.pdConsensusWF* file.


Note A yellow triangle containing an exclamation mark in the upper right corner of the Consensus Step box () usually indicates that a node in the workflow is obsolete, a parameter is missing, a required node is missing, or the output file name is invalid or missing. Hold the mouse pointer over the triangle to display details about what is missing.

Saving an Analysis

To use an analysis as a template for later reuse, you can save it as a .pdAnalysis template file.

Note Studies and analyses in the Proteome Discoverer application are separate, so you must save them separately. Saving a study does not save an analysis, and saving an analysis does not save a study.

❖ To save an analysis as a template for later reuse

1. In the upper right corner of the Analysis window, click the **Save** icon,  **Save** .
2. In the Save Analysis Template dialog box, browse to the location where you want to store the template.
3. In the File Name box, browse to the study folder, and type or select the template file name.
4. In the Save as Type box, select **Analysis Templates (*.pdAnalysis)**.
5. Click **Save**.

The application saves the analysis in a file with a .pdAnalysis extension.

The .pdAnalysis template file saves the processing and consensus workflows. It saves neither the input files nor the study variables that were selected to group the samples and quantification ratios.

Viewing the Workflow and the Analysis After the Search

To view the workflow and the analysis that produced a results file after you perform a search, see [“Viewing the Workflow and the Analysis from the Results”](#) on page 97.

Analysis Window Parameters

Table 3 describes the parameters available in the Analysis window.

Table 3. Analysis window parameters (Sheet 1 of 2)









Parameter	Description
As Batch	Determines whether the application performs an analysis in batch mode. For information, see “Performing a Search in Batch Mode” on page 94.
 Run	Performs an analysis.
 Save	Opens the Save Analysis Template dialog box so that you can save the analysis in a .pdAnalysis template file.

Table 3. Analysis window parameters (Sheet 2 of 2)

Parameter	Description
Consensus Step box	Displays the name of the consensus workflow, the name of the consensus workflow results file, and the processing workflow that generated the MSF file that was submitted to the consensus workflow as an input file.
	Indicates that you cannot change the consensus workflow.
	Opens the Workflow Tree window containing the consensus workflow.
	Indicates that the consensus workflow contains an error.
Workflow	Displays the name of the consensus workflow.
Result File	Displays the name of the file containing the results of the consensus workflow.
Child Steps	Displays the processing workflow.
Add	Adds an empty Processing Step box.
Processing Step box	Displays the name of the processing workflow, the name of the processing workflow results file, and the input files that were submitted to the processing workflow as input files.
	Indicates that you cannot change the processing workflow.
	Opens the Workflow Tree window containing the processing workflow.
	Indicates that the processing workflow contains an error.
Reprocess	Opens a locked processing workflow so that you can make changes to it.
Clone	Opens a child step window so that you can create a different processing workflow for the same file or a different file.
Workflow	Displays the name of the processing workflow.
Result File	Displays the name of the file containing the processing workflow results.
Input Files	Lists the input raw data files that were submitted to the processing workflow.

Specifying Quantification Ratios from Selected Sample Groups

After you set up the workflow to use for the analysis, you can specify the ratios to report for the quantification and how to group your samples in relation to the specified factor values.

You can generate arbitrary ratios from selected sample group abundances without being restricted to predefined sample and control specifications. Previous versions of the Proteome Discoverer application automatically generated quantification ratios that relied on the specification of sample or control sample types. This feature is appropriate for some experiments but too restrictive for experiments in which the same sample should occur in the denominator for some ratios and in the numerator for others.

For example, in a typical mutant-versus-wild-type experiment, you usually would monitor ratios like this to study the effect of the treatment on a mutant:

$$\frac{wt_treat}{wt}, \frac{m_treat}{m}, \frac{m_treat}{wt_treat}$$

where:

- *wt_treat* is a treated wild type.
- *wt* is a wild type.
- *m_treat* is a treated mutant.
- *m* is a mutant.

Also, in some experiments there is no clear control available, and all pair-wise comparisons of available samples are of interest.

On the Grouping & Quantification page of a study, you can generate custom quantification ratios that do not restrict you to predefined sample and control specifications.

In the Grouping & Quantification view, you first select the study variables that were used to group your samples and select the numerators and denominators of your ratios. Then, you manually specify the quantification ratios to generate or semiautomatically create all possible quantification ratios against an ad-hoc selection of denominator values.


❖ To open the Grouping & Quantification page

On the Study: *Study_name* page, click the **Grouping & Quantification** tab.

Note The Grouping & Quantification tab does not appear until you add input files that are associated with a single quantification method.

Figure 52 shows the Grouping & Quantification page of the study.

Figure 52. Grouping and Quantification page



The page contains the following areas:

- **Study Variables area:** At the top left is a list of the selectable study factors, or study variables, that you specified in the Study Factors area of the Study Definition page. A study factor is anything that captures the difference between two samples, for example, drug treatment, the time of drug application, or differences in tissue, organisms, or patients. It can be differences in sample preparation, chromatography settings, or acquisition parameters; or differences in the isobaric or metabolic labels used. Study factors can vary from sample to sample. In the example used throughout this topic, the study factors are the quantification channels associated with each sample and the factors for the acquisition method and tissue. Selecting the study factors to use for grouping provides an effective means of sorting the samples and quantification ratios into replicate and treatment groups.
 - **Files:** Groups samples and quantification ratios by files.
 - **Quan Channels:** Groups samples and quantification ratios by quantification channels.
 - **Study factors:** Groups samples and quantification ratios by user-defined study factors.
 - **Sample Types:** Groups factors that include Sample, Control, Standard, and Blank. Variables displayed in italics contain only a single value.

Study factors are user-defined, but the File, Quan Channel, and Sample Type variables always appear on the Grouping & Quantification page for all studies.

- **Manual Ratio Generation area:** Contains menus where you can select the numerator and denominator for each quantification ratio.
- **Bulk Ratio Generation area:** Displays the study factor values to use as the denominators for semiautomatically generated quantification ratios.
- **Generated Sample Groups area:** Displays samples grouped by the values set for the selected study factors. The application ignores any other differences among the samples that might be present.

Selecting these variables indicates that you want to group your samples and quantification ratios according to the specified study variables and ignore any other difference between the samples that might be present. For example, if you group by tissue, you put into the same sample group all samples that share the same tissue factor. You also group all quantification ratios together that have the same tissue in the numerator or denominator. However, you ignore any differences in the acquisition method used.


The grouping of samples affects how identifications are displayed in the distribution maps (see “[Sample Information Used to Display Identifications and Quantifications Among Files and Samples](#)” on page 455) and the calculation of areas from precursor ion quantification. For the latter, the application reports separate area values for each sample or sample group.

When you group ratios, the application calculates an averaged ratio for a group of ratios that are similar with respect to the selected study variable in the numerator and denominator.

Note The application first groups samples to average abundance values from replicate samples and then calculates ratios from the grouped abundances. The ratios are freely defined in the Grouping and Quantification view, without being restricted to predefined control samples.

If the application does not use all available sample group values in at least one ratio, it displays a warning in the Generated Sample Groups area, like that shown in [Figure 53](#), and marks the unused sample groups.

Figure 53. Warning in the Generated Sample Groups area



- **Generated Ratios area:** Displays the ratios generated.

❖ **To generate custom quantification ratios semiautomatically**

1. In the Study Variables area, select the check box of the study factors, or variables, that you want to use to group your samples and from which you want to draw the numerators and denominators of the ratios. For the example used throughout this topic, select the **Tissue** check box to indicate that the samples and quantification ratios are grouped by the values set for the tissue factor.

After you select the study factors, the Generated Sample Groups area displays the generated sample groups. When performing the quantification, the application calculates abundance values for each sample and averages the abundance values of all samples in a sample group.

This step creates sample groups using the selectable study factors available for the files to be analyzed. The application places all samples that have the same combination of study factor values in the same group. In the example shown in [Figure 58](#) on [page 88](#), the application groups the samples by Acquisition and Tissue. It also groups all samples together that have the same combination of values for acquisition and tissue.


The order of study factors is relevant for the semiautomatic generation of ratios.

2. (Optional) To change the placement of a study factor in the list of study factors, do the following:

- a. In the Study Variables area, select the check box for a study factor.

A placement handle in the form of a green rectangle appears to the left of the selected check box, as shown in [Figure 54](#).

Figure 54. Placement handle in the Study Variables area






- b. Hold the cursor over the placement handle.

White up and down arrows now appear on the handle.

- c. Drag the cursor up or down to move the variable to its new place in the list of variables, or click the up or down arrows to move the study factor.


3. (Optional) To sort study factors in the Bulk Ratio Generation area and sample groups in the Generated Sample Groups area, click one of the following to the right of each study factor in the Study Variable area.

- To sort these items in descending order, click the **Sort Descending** icon, .
- To sort these items in ascending order, click the **Sort Ascending** icon, .
- To leave these items unsorted, click the **No Sorting** icon, .

4. In the Bulk Ratio Generation area, select the check box for the type of tissue to use in the denominator of the ratio. For the example, select **Tissue: Kidney**.

The Bulk Ratio Generation area displays a list of the denominator values for this type of study factor. If you select only one study factor, it displays a list of the available denominator values for this factor, as shown in [Figure 55](#).

Figure 55. Denominator values available for a single study factor




If you select multiple study factors, the Bulk Ratio Generation area displays the denominator values available for each factor.

5. (Optional) To select the same study factor for all the denominators, do the following:
 - a. Hold the cursor over a denominator value.

An icon containing four check boxes in a square appears on the left side of that item, as shown in [Figure 56](#).
 - b. Click the icon.

The application selects the same study factor for all denominators (see [Figure 56](#)).


Figure 56. Icon for multiple selection of denominators in the Bulk Ratio Generation area




6. Click **Add Ratios**.

The application generates all possible ratios against the selected denominator values and adds them to the Generated Ratios area. [Figure 57](#) shows the generated quantification ratios and ratio groups in the Generate Ratios area after the selection of Tissue as the study variable to group by and Kidney as the denominator to use.


Figure 57. Quantification ratios and ratio groups generated by grouping the samples and ratios by tissue



7. (Optional) Select another variable or variables. In the example, select the **Acquisition** check box in the Study Variables area.
When you add a second variable, the information in the Generated Ratios area becomes invalid and appears in a gold color.
8. In the Generate Ratios area, click the **Clear All** icon,  **Clear All**, to delete the previous ratios.
9. In the Bulk Ratio Generation area, select the check box for the type of tissue to use in the denominator. For the example, select **Tissue: Kidney**.
10. Click **Add Ratios**.

[Figure 58](#) shows the generated quantification ratios and ratio groups on the Generated Ratio Groups page after selecting Acquisition and Tissue as the study variable to group by.

Figure 58. Quantification ratios and ratio groups generated by grouping the samples and ratios by tissue and acquisition



❖ **To generate custom quantification ratios manually**

1. On the Study: *Study_name* page, click the **Grouping & Quantification** tab.


Note The Grouping & Quantification tab does not appear until you add or open an analysis and specify Input Files for the analysis.

2. In the Study Variables area, select the check box of the study factor or factors that you want to use to group your samples and from which you want to draw the numerators and denominators of the ratios. For the example used throughout this chapter, select the **Tissue** check box to indicate that the samples and quantification ratios are grouped by the values set for the tissue factor.
3. In the Manual Ratio Generation area, select the numerator value from the Numerator list.
4. Select the denominator value from the Denominator list.

5. Click **Add Ratio**.

The generated ratio appears in the Generated Ratios area, as shown in [Figure 59](#).

Figure 59. Ratio manually generated for one study factor



6. (Optional) Select another study factor or factors in the Study Variables area, for example, **Acquisition**.

When you add a second factor, the information in the Generated Ratios area becomes invalid and appears in a gold color.

7. In the Generated Ratios area, click the **Clear All** icon, , to delete the previous ratios.


8. In the Manual Ratio Generation area, select the numerator value from the Numerator list, for example, **(Kidney, DDA)**.

9. Select the denominator value from the Denominator list, for example, **(Kidney, IDA)**.

10. Click **Add Ratio**.

[Figure 60](#) shows the ratios generated in the Generated Ratios area.

Figure 60. Ratios manually generated for two study factors



❖ To generate custom quantification ratios based on channels


1. In the Study Variables area, select the **Quan Channel** check box.
2. To generate quantification ratios semiautomatically, follow the procedure “To generate custom quantification ratios semiautomatically” on page 84.

—or—

To generate quantification ratios manually, follow the procedure “To generate custom quantification ratios manually” on page 88.

Figure 61 shows the ratios generated in the Generated Ratios area.

Figure 61. Ratios generated by quantification channel



❖ To save the settings on the Grouping & Quantification page

Note The application does not save the settings on the Grouping & Quantification page with a study or with an analysis. Instead, it associates the settings with search results, so you must load them from data sets that have already been processed within the study or recreate them from scratch.

1. Click the **Analysis Results** tab of the study.
2. Select the result on the Analysis Results page, and choose either **Reprocess > All Analysis Steps** or **Reprocess > Last Consensus Step**.


If you choose Use Results to Make New (Multi) Consensus, the Grouping and Quantification page no longer displays ratios.

Performing the Search

You can perform a search in group mode or batch mode.

Performing a Search in Group Mode


❖ To perform a search in group mode

In the upper right corner of the Analysis window, click the **Run** icon, .

The Proteome Discoverer application validates the analysis setup before it starts processing and issues error or warning messages in the Analysis Validation Issues dialog box, shown in [Figure 62](#), if it finds errors. For example, it might issue an error message if not all the input files have the same quantification method. Or, it might issue a warning message if you added several input files with quantification but did not set any of the study variables to group your samples and quantification ratios.


You can ignore warnings, which are marked by an exclamation mark inside a yellow triangle. Because warnings are only hints that the analysis might not be set up correctly, you can click Ignore in the Analysis Validation Issues dialog box. [Figure 62](#) shows an example of the warnings in this dialog box.

Figure 62. Warnings in the Analysis Validation Issues dialog box



You cannot ignore validation errors, which are marked with a red exclamation mark, as shown in [Figure 63](#). You must resolve them.

Figure 63. Error in the Analysis Validation Issues dialog box




The validation cannot detect every potential problem but can check for these specific problems or inconsistencies:

- The set workflows have major errors because of missing mandatory node types such as validator nodes, missing connections, or missing parameters. In most cases, you catch these errors earlier because you cannot click the Run button.
- When you use quantification nodes or have assigned quantification methods to your input files, the validation checks for the following:
 - All quantification nodes used in the workflow of a particular processing step have the same kind of quantification method; for example, you cannot have mixed precursor and reporter ion quantification nodes in a workflow.
 - Not all input files of a particular processing step have the same quantification method assigned; for example, you have mixed precursor and reporter ion quantification files in a processing step.
 - The kind of quantification node used in the workflow of a particular processing step is not compatible with the quantification method set for the input files of this step; for example, you try to process reporter quantification files with a precursor quantification node or vice versa.
 - All input files have a quantification method set, but no quantification node is used in the workflow of the processing step.
 - A quantification node is used in the workflow of a processing step, but none of the input files has a quantification method set.
 - You have not selected any of the study variables for grouping on the Grouping & Quantification page, which is only shown if you process more than one file. If you use only one file, you usually do not want to group. However, you could create biological replicates using different channels of a quantification method such as iTRAQ in an experiment on a single raw data file in which you want to group by tissue, matrix or experimental condition.
 - Some of the processing steps provide quantification values, but the workflow of the consensus step does not contain any node that processes quantification values.
 - The quantification method used for the input files of a particular processing step is not valid.

Once the application validates the analysis, it begins processing it. The job queue opens so that you can monitor the progress of the job, as shown in [Figure 64](#).

Figure 64. Search progress in the job queue



The job is done when “Completed” appears in the Execution State column for the processing and consensus workflows.

The application uses the name of the first raw data file or sample as the default name of the results file as a whole.

Performing a Search in Batch Mode


You might want to process each file in a set of files with the same processing workflow and the same consensus workflow. Processing a set of files in this way is called batch mode. Batch mode is only available if there is more than one input file and if the analysis has just one processing step—that is, if there is just one Processing box in the Analysis window.


❖ To perform a search in batch mode

1. Set up or open a study, as described in “Using Studies” on page 36, and add the input files, as described in “Adding Input Files” on page 48.
2. Set up or open an analysis, as described in “Using Analyses” on page 70.
3. Drag the input files to the Processing Step box.

4. In the title bar of the Analysis window, select the **As Batch** check box, as shown in Figure 65.

Figure 65. As Batch check box selected



5. Click the **Run** icon, .

The application processes all of the input files separately and generates one .pdResult file for each of the input files. The result files have the same name as the single input file.

Working with the Search Results

The following topics explain how to open, close, view, convert, and reprocess search results.

- [Opening the Search Results](#)
- [Deleting Search Results](#)
- [Viewing the Workflow and the Analysis from the Results](#)

- [Converting Results](#)
- [Reprocessing an Existing Analysis](#)

Opening the Search Results

For information on opening the search results, refer to the Help.


Figure 66 shows the search results for the example search used in this chapter.

Figure 66. Search results

Checked	Protein Group ID	# Proteins	# Unique Peptides	# Peptides	# PSMs	Group Description
<input type="checkbox"/>	440	1	42	42	152	[Master Protein] Titin (EC 2.7.11.1) (Connectin) (Rhabdomyosarcoma antigen MU-RMS-40.14) - Homo sapiens (Human)
<input type="checkbox"/>	463	4	18	18	91	[Master Protein] Spectrin alpha chain, brain (Spectrin, non-erythroid alpha chain) (Alpha-II spectrin) (Fodrin alpha chain) - Mus musculus (Mouse)
<input type="checkbox"/>	72	4	17	17	77	[Master Protein] Plectin-1 (Plectin-6) (PLTN) (PCN) - Mus musculus (Mouse)
<input type="checkbox"/>	774	7	13	13	123	[Master Protein] Clathrin heavy chain - Rattus norvegicus (Rat)
<input type="checkbox"/>	513	8	12	18	238	[Master Protein] Myosin-6 (Myosin heavy chain 6) (Myosin heavy chain, cardiac muscle alpha isoform) (MyHC-alpha) - Mus musculus (Mouse)
<input type="checkbox"/>	909	38	12	14	111	[Master Protein] Vimentin - Mus musculus (Mouse)
<input type="checkbox"/>	82	8	10	11	58	[Master Protein] Neurofilament triplet L protein (68 kDa neurofilament protein) (Neurofilament light polypeptide) (NF-L) - Mus musculus (Mouse)
<input type="checkbox"/>	227	9	10	10	133	[Master Protein] Aconitate hydratase, mitochondrial precursor (EC 4.2.1.3) (Citrate hydro-lyase) (Aconitase) - Bos taurus (Bovine)
<input type="checkbox"/>	320	5	10	10	69	[Master Protein] Spectrin beta chain, brain 1 (Spectrin, non-erythroid beta chain 1) (Beta-II spectrin) (Fodrin beta chain) - Mus musculus (Mouse)
<input type="checkbox"/>	899	4	10	10	68	[Master Protein] Pyruvate carboxylase, mitochondrial precursor (EC 6.4.1.1) (Pyruvic carboxylase) (PCB) - Rattus norvegicus (Rat)
<input type="checkbox"/>	294	5	9	12	120	[Master Protein] Alpha-actinin-3 (Alpha actinin skeletal muscle isoform 3) (F-actin cross-linking protein) - Mus musculus (Mouse)
<input type="checkbox"/>	365	26	8	8	105	[Master Protein] ATP synthase subunit beta, mitochondrial precursor (EC 3.6.3.14) - Mus musculus (Mouse)
<input type="checkbox"/>	532	13	8	8	43	[Master Protein] Hexokinase-1 (EC 2.7.1.1) (Hexokinase type I) (HK I) (Hexokinase, tumor isozyme) - Mus musculus (Mouse)
<input type="checkbox"/>	482	14	7	13	123	[Master Protein] Glycogen phosphorylase, muscle form (EC 2.4.1.1) (Myophosphorylase) - Mus musculus (Mouse)
<input type="checkbox"/>	709	23	7	9	81	[Master Protein] Desmin - Mus musculus (Mouse)
<input type="checkbox"/>	58	7	7	7	62	[Master Protein] Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A2 / hnRNP B1) - Mus musculus (Mouse)
<input type="checkbox"/>	923	11	7	7	54	[Master Protein] Catalase (EC 1.11.1.6) - Mus musculus (Mouse)
<input type="checkbox"/>	395	6	7	7	41	[Master Protein] Synaptosomal-associated protein 25 (SNAP-25) (Synaptosomal-associated 25 kDa protein) (Super protein) (SUP) - Homo sapiens (Human)
<input type="checkbox"/>	593	8	7	7	26	[Master Protein] 2-oxoglutarate dehydrogenase E1 component, mitochondrial precursor (EC 1.2.4.2) (Alpha-ketoglutarate dehydrogenase) - Rattus norvegicus (Rat)
<input type="checkbox"/>	328	6	6	13	170	[Master Protein] Myosin-1 (Myosin heavy chain 1) (Myosin heavy chain 2x) (MyHC-2x) (Myosin heavy chain, skeletal muscle, adult 1) - Mus musculus (Mouse)
<input type="checkbox"/>	315	40	6	9	86	[Master Protein] Heat shock cognate protein HSP 90-beta - Gallus gallus (Chicken)
<input type="checkbox"/>	661	22	6	8	83	[Master Protein] Sodium/potassium-transporting ATPase alpha-1 chain precursor (EC 3.6.3.9) (Sodium pump 1) (Na ⁺ /K ⁺ ATPase 1) - Rattus norvegicus (Rat)
<input type="checkbox"/>	270	12	6	7	90	[Master Protein] Alpha-internexin (Alpha-Inx) - Rattus norvegicus (Rat)
<input type="checkbox"/>	410	4	6	6	65	[Master Protein] Ubiquinol-cytochrome-c reductase complex core protein 1, mitochondrial precursor (EC 1.10.2.2) (Core I protein) - Mus musculus (Mouse)
<input type="checkbox"/>	380	4	6	6	63	[Master Protein] Carbamoyl-phosphate synthase [ammonia], mitochondrial precursor (EC 6.3.4.16) (Carbamoyl-phosphate synthetase I) (CPSase I) - Mus musculus (Mouse)
<input type="checkbox"/>	509	3	6	6	57	[Master Protein] Selenium-binding protein 2 (56 kDa acetaminophen-binding protein) (AP56) - Mus musculus (Mouse)
<input type="checkbox"/>	555	3	6	6	47	[Master Protein] Lamin-B1 - Mus musculus (Mouse)
<input type="checkbox"/>	876	11	6	6	39	[Master Protein] Glutamate dehydrogenase 1, mitochondrial precursor (EC 1.4.1.3) (GDH) (Memory-related protein 2) (MRG-2) - Rattus norvegicus (Rat)
<input type="checkbox"/>	129	12	6	6	31	[Master Protein] Moesin (Membrane-organizing extension spike protein) - Mus musculus (Mouse)
<input type="checkbox"/>	86	6	6	6	30	[Master Protein] Transitional endoplasmic reticulum ATPase (TER ATPase) (15S Mg ²⁺ -ATPase p97 subunit) (Valosin-containing protein) (VCF) - Mus musculus (Mouse)
<input type="checkbox"/>	772	2	6	6	29	[Master Protein] Spectrin beta chain, brain 2 (Spectrin, non-erythroid beta chain 2) (Beta-III spectrin) (SPNB-3) (Beta Spill sigma 1) (Spectrin-like) - Mus musculus (Mouse)
<input type="checkbox"/>	346	7	6	6	26	[Master Protein] Plastin-2 (L-plastin) (Lymphocyte cytosolic protein 1) (LCP-1) (65 kDa macrophage protein) (pp65) - Mus musculus (Mouse)
<input type="checkbox"/>	174	3	6	6	23	[Master Protein] Peroxisomal multifunctional enzyme type 2 (MFE-2) (D-bifunctional protein) (DBP) (17-beta-hydroxysteroid dehydrogenase 4) (17-OHSD4) - Mus musculus (Mouse)
<input type="checkbox"/>	664	3	6	6	21	[Master Protein] Prohibitin-2 (B-cell receptor-associated protein BAP37) (BAP-37) - Rattus norvegicus (Rat)
<input type="checkbox"/>	686	3	6	6	13	[Master Protein] Laminin subunit beta-2 precursor (S-laminin) (Laminin chain B3) - Rattus norvegicus (Rat)
<input type="checkbox"/>	878	11	5	7	103	[Master Protein] Sodium/potassium-transporting ATPase alpha-3 chain (EC 3.6.3.9) (Sodium pump 3) (Na ⁺ /K ⁺ ATPase 3) (Alpha(III)) - Homo sapiens (Human)
<input type="checkbox"/>	123	6	5	5	54	[Master Protein] Aspartate aminotransferase, cytoplasmic (EC 2.6.1.1) (Transaminase A) (Glutamate oxaloacetate transaminase 1) - Mus musculus (Mouse)
<input type="checkbox"/>	929	7	5	5	54	[Master Protein] Electron transfer flavoprotein subunit beta (Beta-ETF) - Mus musculus (Mouse)
<input type="checkbox"/>	275	1	5	5	46	[Master Protein] Liver carboxylesterase 31 precursor (EC 3.1.1.1) (ES-Male) (Esterase-31) - Mus musculus (Mouse)
<input type="checkbox"/>	222	6	5	5	42	[Master Protein] Aldehyde dehydrogenase, mitochondrial precursor (EC 1.2.1.3) (ALDH class 3) (ALDH M3) (ALDH E3) - Mus musculus (Mouse)

Deleting Search Results

You can delete search results from a study.


1. In an open study, click the **Analysis Results** tab if it is not already selected.
2. Select the row of results that you want to delete, and select the **Remove Files** icon,  **Remove Files**, or press the DELETE key.

3. In the Remove Analysis Result dialog box, click **Yes**.

Viewing the Workflow and the Analysis from the Results


After you perform a search, you can see the workflow and the analysis that produced the search results.

❖ To access the workflow and the analysis from the results

1. In an open study, click the **Analysis Results** tab.
2. Select the name of the results file.
3. Click the **Show Details** icon,  Show Details .

The Analysis Sequence Details window opens, as shown in [Figure 67](#). It displays the Workflow Editor and the Analysis window that contains the processing and consensus workflows used to generate the selected results file.

Figure 67. Analysis Sequence Details window



Searching Analysis Results

You can search for a specified string in the search results on the Analysis Results page of a study. Filtering the analysis results in studies enables you to reprocess previous workflows more quickly.


❖ To search for a specified string

1. In an open study, click the **Analysis Results** tab.
2. In the Search box, shown in [Figure 68](#), type the string that you want to search for.

By default, the application searches for the string in the File Name and Description columns.

3. (Optional) To search in a single column, select the name of the column from the Search For menu to the right of the Search box.

Figure 68. Search box and Search For menu on the Analysis Results page



Converting Results

You can convert MSF files to .pdResult files or use .pdResult files to obtain a different results file. You can perform this conversion in two ways:


- Add result files to a new study.
- Add result files to an existing study.

You can process these files as a batch and generate a .pdResult file for each input file or one .pdResult file for multiple input files.

Adding Result Files to a New Study

You can add existing MSF or .pdResult files to a new study.


❖ To add MSF or .pdResult files to a new study

1. Create a study. For instructions, see “Creating a Study” on page 38.
2. In the New Study and Analysis dialog box, click the **Add Files** icon,  **Add Files** .
3. In the Add Input File dialog box, select the MSF files, .pdResult files, or both that you want to import.
4. Click **Open**.
5. In the New Study and Analysis dialog box, click **OK**.

An Analysis window opens that contains a processing step for each file added to the study.

6. Create a consensus workflow. For instructions, see “Creating a Consensus Workflow” on page 112.
7. (Optional) To generate a .pdResults file for each input file, click the **As Batch** check box in the Analysis window title bar.


When you select the As Batch check box, the application creates a .pdResult file for each input file. When you do not select it, it creates one .pdResult file for all input files.

8. Click the **Run** icon,  **Run** .

Adding Result Files to an Existing Study


You can add MSF or .pdResult files to an existing study.


❖ To add result files to an open study

1. In an open study, click the **Add Files** icon,  **Add Files** .
2. In the Add Files dialog box, browse to the location of the MSF or .pdResult files that you want to add to the existing results in the study, select them, and click **Open**.

In [Figure 69](#), an MSF file is added to the Analysis Results page, which contains a .pdResult file.

Figure 69. Existing results file and the added MSF file on the Analysis Results page




1. Click the **Analysis Results** tab.
2. On the Analysis Results page, select the existing result files and the added result files.
3. Click the **Reuse** icon,  **Reuse** .

An Analysis window opens that contains a processing step for each result file in the study.

4. Create a consensus workflow. For instructions, see “[Creating a Consensus Workflow](#)” on [page 112](#).
5. (Optional) To generate a .pdResults file for each input file, click the **As Batch** check box in the Analysis window title bar.

When you select the As Batch check box, the application creates a .pdResult file for each input file. When you do not select it, it creates one .pdResult file for all input files.


6. Click the **Run** icon,  **Run** .

Reprocessing an Existing Analysis

You can change the input files or the workflow used to process the RAW data files and then reprocess the existing results. You can decide to reprocess all analysis steps or only the last consensus step. You can also use a new consensus workflow to reprocess MSF results associated with one or more .pdResult files.


Using the Reprocess command with the Reuse command, you can select multiple analysis results at once to reuse for a new consensus from existing analysis results.

❖ To reprocess an existing search by using different parameters

1. In an open study, click the **Analysis Results** tab if it is not already selected.
2. Select the search that you want to reprocess.
3. Click the **Reprocess** icon,  **Reprocess** ▾ .
4. To rerun all the analysis steps, choose **All Analysis Steps**.


–or–

To rerun just the last step in the consensus workflow, choose **Last Consensus Step**.

5. Change the appropriate workflows, settings, parameters, input files, and so forth.
6. In the Analysis window, click the **Run** icon,  **Run** .


The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

❖ **To reprocess existing MSF results by using a new consensus workflow**


1. Click the **Analysis Results** tab.
2. Select the first .pdResult file associated with the MSF file that you want to reprocess with a new consensus workflow.
3. Click the **Reprocess** icon,  **Reprocess** ▾ .
4. Choose **Use Results to Create New (Multi) Consensus**.

The application adds the analysis associated with the selected .pdResult file, as shown in [Figure 70](#).

Figure 70. Consensus and processing steps associated with a .pdResult file




The lock symbol, , indicates that you cannot change the processing workflow. (If you want to change the processing workflow click **Reprocess** on the Processing Step title bar.)


5. (Optional) Select any other .pdResult files associated with additional MSF files that you want to reprocess with a new consensus workflow, click the **Reprocess** icon,  **Reprocess** ▾ , and choose **Use Results to Create New (Multi) Consensus**.
6. Click the **Workflows** tab, and then click the **Consensus** tab, if it is not already open.
7. Create a new consensus workflow. For instructions, see [“Creating a Consensus Workflow”](#) on [page 112](#).

—or—


Open an existing workflow by doing one of the following:

- In the Workflow Editor within the analysis, click the **Open** icon,  **Open** , browse to the location where you stored the consensus workflow, select the .pdConsensusWF file, and click **Open**.

—or—

- If you stored the consensus workflow in the Common Templates directory (see “Using Custom Workflow Templates” on page 123), click the **Open Common** icon,  **Open Common**, select the .pdConsensusWF file, and click **Open**.

The template appears in the Workflow Tree pane of the Workflow Editor.

8. Click the **Run** icon,  **Run**.

Using the Workflow Editor to Create Workflows

This topic gives detailed information about creating processing and consensus workflows.

- [Opening the Stand-alone Workflow Editor](#)
- [Closing the Workflow Editor](#)
- [Creating a Processing Workflow](#)
- [Creating a Consensus Workflow](#)
- [Incorporating an Existing Workflow into a Study and an Analysis](#)
- [Workflow Editor Parameters](#)


Opening the Stand-alone Workflow Editor

You can open the Workflow Editor in the context of an analysis, or you can open the stand-alone Workflow Editor outside of an analysis.

❖ To open the Workflow Editor from inside an analysis


1. In a study, create an analysis. See “Creating an Analysis” on page 71.
2. Click the **Workflows** tab.

–or–

To create a consensus workflow, click the **Show Workflow** icon, , on the title bar of the Consensus Step box. To create a processing workflow, click the same icon on the title bar of the Processing Step box.

The Workflow Editor opens in the middle of the Study: *Study_name* page.

❖ To open the stand-alone Workflow Editor

Choose **View > Workflow Editor** or click the **Workflow Editor** icon, .

The Workflow Editor opens on a separate page called Workflow Editor.

Closing the Workflow Editor

To close the Workflow Editor in the context of an analysis, you must close the analysis itself.

❖ **To close the Workflow Editor inside an analysis**

Click the **X** on the Analysis title bar.

❖ **To close the stand-alone Workflow Editor**


Click the **X** on the Workflow Editor tab.

Creating a Processing Workflow

This topic expands on the basic procedure outlined in “[To create the processing workflow](#)” on [page 74](#). This procedure is the same whether you use the Workflow Editor as a stand-alone tool or in the context of an analysis.

❖ **To create a new processing workflow**

1. To use the Sequest HT search engine or the Mascot search engine in the workflow, configure it first by following the instructions in “[Configuring the Sequest HT Search Engine](#)” on [page 22](#) or “[Configuring the Mascot Search Engine](#)” on [page 25](#).


2. Choose **View > Workflow Editor** or click the **Workflow Editor** icon,  .

The Workflow Editor opens, as shown in [Figure 19](#) on [page 35](#).

3. Click the **Workflow Nodes** tab in the Workflow Nodes pane.
4. Click the **Processing Workflow** tab on the Workflow Nodes page.
5. From the Data Input area of the Workflow Nodes pane, drag the **Spectrum Files** node to the Workflow Tree pane.

The Spectrum Files node is required. For its location, see [Figure 71](#), and for more information, refer to the Help.

Figure 71. Selected data input node in the Processing Workflow window



6. Drag the **Spectrum Selector** node to the workspace and place it beneath the Spectrum Files node.

Selecting the Spectrum Selector node in the Workflow Tree pane displays its parameters in the right pane.

Note You can set the Spectrum Selector node to select which precursor mass to use for a given MSⁿ scan, such as choosing the precursor from the parent scan.

7. Drag a search engine node (**Mascot** or **Sequest HT**) to the workspace and place it below the Spectrum Selector node.
8. Drag a PSM Validation node (**Fixed Value PSM Validator**, **Percolator**, or **Target Decoy PSM Validator**) to the workspace and place it below the search engine node.
9. Depending on your data needs, drag any other appropriate nodes from the Workflow Nodes pane to the Workflow Tree pane.

For a description of the nodes that you can select, refer to the Help.

You can also add third-party nodes that are in your installation but are not documented in this manual. For further information on those nodes, consult the relevant third-party documentation.

The numbers that appear on each workflow node indicate the order in which the application processes the nodes.

You cannot drag workflow nodes into the Workflow Tree pane that cannot logically be added. For example, you cannot add the Reporter Ions Quantifier node and the Precursor Ions Quantifier node to the workspace together because they are mutually exclusive.

10. Organize the nodes to reflect a procedural order from top to bottom so that the Spectrum Files node remains on top as the root node.


Delete a node by selecting the node in the Workflow Tree pane and pressing DELETE or by right-clicking the node and choosing **Cut** (or **CTRL+X**) from the shortcut menu.

You can use the Cut command and the Paste (or CTRL+P) command on the shortcut menu to move a node to another place in the workspace or use the Copy (or CTRL+C) and Paste commands to duplicate a node in the workspace.

You can paste copied or cut nodes into other workflows.

11. Connect the nodes if they do not automatically connect:
 - a. Click the top node to enable a blue handle at the bottom center of the node, as shown in [Figure 72](#).


Figure 72. Activated node



Joining the nodes together creates a sequence of steps for the application to follow.


- b. Drag the blue handle down to the node below it, as shown in [Figure 73](#).

Figure 73. Joining two nodes




IMPORTANT If the next node appears with a red edge at this point, you cannot connect it to the previous node.

If the Workflow Editor prevents you from connecting two nodes, the workflow is erroneous.

- c. Link all the nodes to develop a workflow.
12. After you join all your chosen nodes, align them by clicking the **Auto Layout** icon,  **Auto Layout**, or right-clicking a node and choosing **Auto Layout** from the shortcut menu.
 13. (Optional) To renumber the nodes, right-click and choose **Auto Number** from the shortcut menu.

[Figure 74](#) shows a basic processing workflow.

Figure 74. Basic processing workflow



14. Set the parameters for each node in the Parameters pane as follows:
 - a. Click the node.
 - b. (Optional) Click **Show Advanced Parameters** so that you can view all the node's parameters.

The available parameters for the node appear in the Parameters pane, as shown in the example for the Spectrum Selector node in [Figure 75](#).

Figure 75. Parameters for Spectrum Selector node in the Parameters pane


Parameters	
<input type="button" value="Hide Advanced Parameters"/>	
1. General Settings	
Precursor Selection	Use MS1 Precursor
Use New Precursor Reevaluation	True
Use Isotope Pattern in Precursor	True
2. Spectrum Properties Filter	
Lower RT Limit	0
Upper RT Limit	0
First Scan	0
Last Scan	0
Ignore Specified Scans	
Lowest Charge State	0
Highest Charge State	0
Min. Precursor Mass	350 Da
Max. Precursor Mass	5000 Da
Total Intensity Threshold	0
Minimum Peak Count	1
3. Scan Event Filters	
Mass Analyzer	(Not specified)
MS Order	Is Not MS1
Activation Type	(Not specified)
Min. Collision Energy	0
Max. Collision Energy	1000
Scan Type	Is Full
Polarity Mode	(Not specified)
4. Peak Filters	
S/N Threshold (FT-only)	1.5
5. Replacements for Unrecogn	
Unrecognized Charge Replacement	Automatic
Unrecognized Mass Analyzer Replacement	ITMS
Unrecognized MS Order Replacement	MS2
Unrecognized Activation Type Replacement	CID
Unrecognized Polarity Replacement	+
Unrecognized MS Resolution@21	60000
Unrecognized MSn Resolution@:	30000
6. Precursor Pattern Extraction	
Precursor Clipping Range Before	2.5 Da
Precursor Clipping Range After	5.5 Da

- c. Set each node's parameters. Complete this step for each node that you select.

With the Mascot and Sequest HT nodes, you can specify multiple FASTA protein databases to search at one time.

From the dropdown list of the Sequest HT node's Protein Database parameter, select the check boxes next to the names of the appropriate FASTA files, as shown in [Figure 76](#). There is no limit to the number of FASTA files that you can select.


Figure 76. Selecting multiple FASTA databases for the Sequest node



To search multiple databases with Mascot server 2.3 or later, see [“Searching Multiple Sequence Databases with Mascot”](#) on page 134.

Figure 77 shows the parameters set for the Sequest HT node.

Figure 77. Setting parameters for the Sequest HT node



When you click some parameters, two lists appear, as shown in [Figure 78](#).

Figure 78. Settings and filters

3. Scan Event Filters	
Mass Analyzer	(Not specified)
MS Order	Is Not MS1
Activation Type	Is Not MS1
Min. Collision Energy	Any
Max. Collision Energy	Is
Scan Type	Is Not
Polarity Mode	
4. Peak Filters	
S/N Threshold (FT-only)	1.5
5. Replacements for Unrecognized Properties	
Unrecognized Charge Replacements	Automatic
Unrecognized Mass Analyzer Replacements	ITMS
Unrecognized MS Order Replacements	MS2
Unrecognized Activation Type Replacements	CID

You can use the list on the left to apply a filter to the setting that you select in the list on the right. The list on the left consists of three options:



- **Is:** Applies the setting selected in the list on the right. In the example in [Figure 78](#), “Is” means that the workflow processes data from the CID activation type.
- **Is Not:** Applies all settings in the list on the right except the selected setting. In the example in [Figure 78](#), “Is Not” means that the workflow processes data from all activation types except CID.
- **(Default) Any:** Applies all settings available for the parameter in the list on the right. In the example in [Figure 78](#), “Any” means that the workflow processes data from any activation type available in the list on the right.

You can filter input data before searching the database to remove lower-quality spectral peak lists from your analysis. This step might help to decrease search times and false positive identifications. The Spectrum Filters area of the Workflow Nodes pane for the processing workflow provides three types of spectrum filters to use for your search. Use these pre-analysis filters to streamline your search results. For information about these nodes, refer to the Help.

Use the Scan Event Filter node for high-mass-accuracy data, such as Mascot analysis and Sequest HT analysis of mixed fragmentation-mode-type data (CID and ETD). It can filter information according to fragmentation type, mass analyzer identity, and other parameters. For information about the Scan Event Filter node, refer to the Help.

15. (Optional) In the Name box, type a name for the workflow.
16. (Optional) In the Description box, type a brief description of the workflow.


Providing a workflow description is highly recommended.

17. (Optional) Save the processing workflow as a template in a .pdProcessingWF file:
 - a. Click the **Save** icon,  **Save** .
 - b. In the Save Workflow dialog box, do the following:
 - i. Select the file to save the workflow in, or type a file name in the File Name box. You can save the workflow in the study folder or in the Common Templates folder (select the **Save Common** icon,  **Save Common** , in this case), or in a separate folder of workflows.
 - ii. Click **Save**.

The application saves the workflow in a file with a .pdProcessingWF extension.


For more information on using templates, see “[Using Workflow Templates](#)” on [page 119](#).

The application returns you to the Processing Workflow window.

18. In the Analysis window, click the **Run** icon,  **Run** .

The job queue appears, as shown in [Figure 79](#), displaying the status of your search.

Figure 79. Status of the processing workflow in the job queue



Use the job queue to check the status of your search.

For information about the job queue, refer to the Help.


The application creates an MSF file containing the results. The Processing Files pane lists the name of the output MSF file, which is the same as the name of the input file but with an .msf extension.

Creating a Consensus Workflow

You can create a new consensus workflow that generates a results report for a new processing workflow. Or, you can create a stand-alone consensus workflow if you have a preexisting MSF file for which you want to generate a results report.

For a description of all the nodes available in the consensus workflow, refer to the Help.

❖ To create a consensus workflow


1. Choose **View > Workflow Editor** or click the **Workflow Editor** icon, .

The Workflow Editor opens, as shown in [Figure 19](#) on [page 35](#).

2. Create a processing workflow by following the instructions in “[Creating a Processing Workflow](#)” on page 104.
3. Click the **Consensus Workflow** tab in the Workflow Nodes pane.
4. From the Data Input area of the Workflow Nodes pane, drag the **MSF Files** node to the Workflow Tree pane.

The MSF Files node is required. [Figure 80](#) shows the location of this node.

Figure 80. Selected data input node in the Consensus Workflow window



5. Drag the following nodes to the workspace and place them beneath the MSF Files node.

IMPORTANT Although these nodes are not required, they generate meaningful results. Reports generated from a workflow consisting of only the MSF Files node contain only the Proteins, PSMs, and MS/MS Spectrum Info pages. The PSMs and proteins have no confidences, and the proteins have no scores.

- PSM Grouper node: For information, refer to the Help.
 - Peptide Validator node: For information, refer to the Help.
 - Peptide and Protein Filter node: For information, refer to the Help.
 - Protein Scorer node: For information, refer to the Help.
 - Protein Grouping node: For information refer to the Help.
6. If you performed quantification in the Processing workflow, drag the **Peptide and Protein Quantifier** node to the Workflow Tree pane.

This node is required if you performed quantification. For more information on this node, refer to the Help.

7. Depending on your data needs, drag any other appropriate nodes from the Workflow Nodes pane to the Workflow Tree pane.

For a description of the nodes that you can select, refer to the Help.

The numbers that appear on each workflow node indicate the order in which the application processes the nodes.

You can also add third-party nodes that are in your installation but are not documented in this manual. For further information on those nodes, consult the relevant third-party documentation.

You cannot drag workflow nodes into the Workflow Tree pane that cannot logically be added.


8. Organize the nodes to reflect a procedural order from top to bottom so that the MSF Files node remains on top as the root node.

Delete a node by selecting the node in the Workflow Tree pane and pressing DELETE, or by right-clicking the node and choosing **Cut** (or **CTRL+X**) from the shortcut menu.

You can use the Cut command and the Paste (or CTRL+P) command on the shortcut menu to move a node to another place in the workspace, or use the Copy (or CTRL+C) and Paste commands to duplicate a node in the workspace.


You can paste copied or cut nodes into other workflows.

The application automatically connects the nodes in the consensus flow. If two adjacent nodes do not connect, there is an error in the flow, and you must rearrange the nodes until they connect.

9. Align the nodes by clicking the **Auto Layout** icon,  **Auto Layout**, or right-clicking a node and choosing **Auto Layout** from the shortcut menu.
10. (Optional) To renumber the nodes, right-click and choose **Auto Number** from the shortcut menu.
11. Set the parameters for each node in the Workflow Tree pane, following the instructions in [step 14 on page 107](#) of “[Creating a Processing Workflow](#).”

[Figure 81](#) shows a basic consensus workflow that does not involve quantification.

Figure 81. Basic consensus workflow





12. Drag any appropriate nodes from the Post-Processing category in the Workflow Nodes pane to the Post-Processing Nodes pane.

You cannot place nodes from the Post-Processing category in the Workflow Tree pane. Similarly, you cannot place nodes from any other category in the Post-Processing Nodes pane.

13. (Optional) In the Name box, type a name for the workflow.
14. (Optional) In the Description box, type a brief description of the workflow.

Providing a workflow description is highly recommended.


15. (Optional) Save the consensus workflow as a template in a .pdConsensusWF file as follows:

- a. Click the **Save** icon,  **Save**.
- b. In the Save Workflow dialog box, do the following:
 - i. Select the file to save the workflow in, or type a file name in the File Name box. You can save the workflow in the study folder or in the Common Templates folder (select the **Save Common** icon,  **Save Common**, in this case), or in a separate folder of workflows.
 - ii. Click **Save**.

The application saves the workflow with a .pdConsensusWF extension.


For more information on using templates, see “Using Workflow Templates” on page 119.

The application returns you to Consensus Workflow window.

16. In the Analysis window, click the **Run** icon,  **Run** .

The job queue appears, as shown in Figure 82, displaying the status of the consensus flow.

Figure 82. Status of the consensus workflow



Use the job queue to check the status of the consensus workflow.





For information about the job queue, refer to the Help.

The Proteome Discoverer application generates a .pdResult file. For information on opening and using the .pdResult file, refer to the Help.


Incorporating an Existing Workflow into a Study and an Analysis

You might want to create a workflow outside of the context of a study or an analysis and incorporate it into a study and an analysis later.


❖ **To incorporate an existing workflow into a study and an analysis**

1. Click **View > Workflow Editor** or click the **Workflow Editor** icon, .
The stand-alone Workflow Editor opens on a separate page.
2. Create the processing workflow. For instructions, see “[Creating a Processing Workflow](#)” on [page 104](#).
3. Click the **Save** icon,  **Save**, browse to the location where you want to save the .pdProcessingWF file, and click **Save**.
4. Create the consensus workflow. For instructions, see “[Creating a Consensus Workflow](#)” on [page 112](#).
5. Click the **Save** icon,  **Save**, browse to the location where you want to save the .pdConsensusWF file, and click **Save**.
6. Create or open a study and an analysis:
 - To create a study, see “[Creating a Study](#)” on [page 38](#).
 - To open an existing study, see “[Opening an Existing Study](#)” on [page 40](#).
 - To create an analysis, see “[Creating an Analysis](#)” on [page 71](#).
 - To open an existing analysis, see “[Opening an Existing Analysis](#)” on [page 72](#).
7. Click the **Workflows** tab to open the Workflow Editor within the analysis.
8. In the Workflow Editor within the analysis, click the **Open** icon,  **Open**, browse to the location where you stored the processing workflow, select the .pdProcessingWF file, and click **Open**.


–or–

If you stored the processing workflow in the Common Templates directory (see “[Using Custom Workflow Templates](#)” on [page 123](#)), click the **Open Common** icon,  **Open Common**, select the .pdProcessingWF file, and click **Open**.

The template appears in the Workflow Tree pane of the Workflow Editor.

9. In the Workflow Editor within the analysis, click the **Open** icon,  **Open**, browse to the location where you stored the consensus workflow, select the .pdConsensusWF file, and click **Open**.

–or–

If you stored the consensus workflow in the Common Templates directory (see “[Using Custom Workflow Templates](#)” on [page 123](#)), click the **Open Common** icon,  **Open Common**, select the .pdConsensusWF file, and click **Open**.

The template appears in the Workflow Tree pane of the Workflow Editor.

Workflow Editor Parameters

Table 4 describes the parameters available in the Workflow Editor.

Table 4. Workflow Editor parameters (Sheet 1 of 2)







Parameter	Description
Workflow Nodes pane	Lists all the nodes available for constructing a processing and a consensus workflow and the parameters for each node. This pane includes the Processing Workflow page and the Consensus Workflow page. For a detailed description of each node shown in these pages, refer to the Help.
Processing Workflow page	Lists all the nodes available for constructing a processing workflow.
Consensus Workflow page	Lists all the nodes available for constructing a consensus workflow.
Parameters page	Lists all the parameters that you can set for a selected node.
Show/Hide Advanced Parameters	Determines whether more advanced parameters are displayed or hidden in the Parameters pane.
 Open	Opens the Open Workflow dialog box so that you can select an existing workflow.
 Open Common	Opens the Open Workflow dialog box in the Common Templates area so that you can select a processing or consensus workflow template.
 Save	Opens the Save Workflow dialog box so that you can save a processing or consensus workflow. By default, the dialog box opens in the study directory. The type of workflow that it saves depends on whether you selected the Processing Workflow tab or the Consensus Workflow tab in the Workflow Nodes pane.
 Save Common	Opens the Save Workflow dialog box in the Common Templates area so that you can save a processing or consensus workflow template to the Common Templates directory.
 Auto Layout	Automatically adjusts and aligns the connecting arrows and nodes in the Workflow Tree pane in the Workflow Editor. For information about using this command, see “Creating a Processing Workflow” on page 104 and “Creating a Consensus Workflow” on page 112 .

Table 4. Workflow Editor parameters (Sheet 2 of 2)

Parameter	Description
 Clear	Removes the contents of the Workflow Tree pane.
Workflow	Displays the name of the current workflow in the Workflow Tree pane.
Description	Displays a brief description of the current workflow in the Workflow Tree pane.
Workflow Tree pane	Provides a workspace for placing all nodes except the post-processing nodes into a workflow layout.
Post-Processing Nodes pane	Provides a workspace for placing the post-processing nodes in a consensus workflow.

Using Workflow Templates

A workflow template is an XML file that contains information about a processing workflow or a consensus workflow. For processing workflows, this file has a .pdProcessingWF suffix. For consensus workflows, it has a .pdConsensusWF suffix.

The Proteome Discoverer application offers standard workflow templates, or you can create your own.

See these topics:

- [Using Common Workflow Templates](#)
- [Using Custom Workflow Templates](#)
- [Correcting Template Errors](#)
- [Saving a Workflow as a Common Template](#)
- [Saving a Workflow as a Custom Template](#)
- [Deleting a Workflow Template](#)
- [Opening a Workflow from a .pdResult File](#)

Using Common Workflow Templates

When you install the Proteome Discoverer application, it creates a Common Templates folder containing subfolders of factory default templates for processing workflows, consensus workflows, analyses, and filter sets. The processing workflows are divided into subfolders of templates appropriate for Q Exactive, LTQ Orbitrap, and Orbitrap Fusion output. You and other users of the system can also use this folder to store common workflow templates that you can then share.

The application installs the Common Templates folder in the following location:


C:\Users\Public\Public Documents\Thermo\Proteome Discoverer 2.1\Common Templates\

You can use either of the following two methods to access the Common Templates folder and copy your preferred workflow templates to and from it.

❖ **To open the Common Templates folder through the Documents library in Windows Explorer**


1. Open Windows Explorer.
2. Choose **Libraries > Documents > Public Documents > Thermo > Proteome Discoverer 2.1 > Common Templates**, as shown in [Figure 83](#).

Figure 83. Location of the Common Templates folder through the Documents library



❖ **To open the Common Templates folder directly in Windows Explorer**

1. Open Windows Explorer.
2. Choose **Drive C: > Users > Public > Public Documents > Thermo > Proteome Discoverer 2.1 > Common Templates**, as shown in [Figure 84](#).

Figure 84. Location of the Common Templates folder through the local disk directory

Note The C:\Users\Public\Documents\ folder is a special folder in Windows. The absolute path to this folder, as well as the folder names shown in Windows Explorer, might differ on different systems, depending on your Windows version. Programmatically, always access these special folders by their symbolic names. To find out their exact location on your system, you can refer to the PUBLIC environment variable that is set to the following:

```
PUBLIC = C:\Users\Public
```

or to the following registry entry:


```
HKEY_LOCAL_MACHINE\SOFTWARE\Microsoft\Windows\CurrentVersion\Explorer\Shell Folders\Common Documents
```

The Documents library is often a concatenation of two different folders. To see them, do the following:

1. Open Windows Explorer.
2. Choose **Libraries > Documents**.
3. Click **2 locations**.

The Documents Library Locations dialog box, shown in [Figure 85](#), opens to show the actual folders that are joined together.

Figure 85. Document Library Locations dialog box




Follow these instructions to open a template in the Common Templates folder.

❖ **To open a common workflow template**

1. Create an analysis (for instructions, see “[Creating an Analysis](#)” on page 71), or open an existing analysis (for instructions, see “[Opening an Existing Analysis](#)” on page 72).

The Workflows tab appears.

2. Click the **Workflows** tab.
3. Click the Consensus Step bar or the Processing Step bar in the Analysis window.
4. In the Workflow Editor window, click the **Open Common** icon,  **Open Common**.

If you clicked the Processing Step bar, the Open Workflow dialog box displays .pdProcessingWF files in the Common Templates folder. If you clicked the Consensus Step bar, the dialog box displays .pdConsensusWF files in the Common Templates folder.


5. In the Open Workflow dialog box, select a .pdProcessingWF file for a processing workflow or a .pdConsensusWF file for a consensus workflow from the Common Templates folder.
6. Click **Open**.

The workflow appears in the Workflow Tree area. For information on correcting any template errors displayed, see “[Correcting Template Errors](#)” on page 123.


Using Custom Workflow Templates

Custom workflow templates are templates that you created yourself and stored. If you saved a workflow as part of an analysis, you can access the workflow by opening the analysis template.

❖ To open a custom workflow template for a new analysis

1. Create an analysis. For instructions, see “Using Analyses” on page 70.
2. Click the **Open Analysis Template** icon,  **Open Analysis Template** .
3. In the Open Analysis Template dialog box, select the .pdAnalysis file that you want to use as a template.
4. Click **Open**.

❖ To open a custom workflow template in an existing analysis

1. Open an existing analysis.
The Workflows tab appears.
2. Click the **Workflows** tab.
3. Click the Consensus Step or the Processing Step title bar in the Analysis window.
4. In the **Workflow Editor** window, click the **Open** icon,  **Open** .
5. In the Open Workflow dialog box, browse to the .pdProcessingWF or .pdConsensusWF file that the custom template is stored in, and click **Open**.

The workflow appears in the Workflow Tree area. For information on correcting any template errors displayed, see “Correcting Template Errors.”

Correcting Template Errors

When you open a template file, either standard or custom, the Proteome Discoverer application validates parameter settings and displays the selected workflow in the Workflow Tree area of the Workflow Editor. The name of the workflow appears in the Workflow box, and a description of the workflow appears in the Description box. The application uses warning symbols to indicate outdated nodes and to display error information in the Workflow Failures pane.

If the application cannot load the selected file, it displays a message box with information about the issue. It cannot load files that cannot be updated, are read-only or invalid, or were created with a newer version of the Proteome Discoverer application.

When you open an existing workflow template, some of the nodes in the Workflow Tree pane might show a triangular yellow warning symbol, shown in [Figure 86](#). This symbol indicates that the version of the node that was used when the template was created has been superseded by a later version in the current Proteome Discoverer application. When you save or resave the workflow, the application automatically updates the node to the latest version.

Figure 86. Warning symbol indicating an outdated node version




Figure 87 shows a round blue warning symbol containing an exclamation point, which indicates that one or more of the parameter settings for the node are incorrect or outdated. Click on the node and reset the parameters in the Parameters pane.

Figure 87. Warning symbol indicating incorrect or outdated parameter settings





Figure 88 shows a blocked gear (or settings) symbol, which indicates that the flow contains a node that is no longer available for the given workflow because it is outdated, not installed, installed but not licensed, or has incorrect parameter settings.


Figure 88. Warning symbol indicating a node no longer available for the given workflow



When you open a workflow that contains an unavailable node, a Load Failures pane opens beneath the Workflow Nodes pane, as shown in Figure 89.

If this warning symbol is attached to a node, the application cannot update it because it is not licensed or no longer licensed, it has been replaced by a different node, or the node was moved to the consensus workflow. In these cases, it does not permit you to save the workflow.

Figure 89. Load Failures pane



Load Failures pane

The Load Failures pane has three columns:

- **Error Information:** Displays information about the problem that the application encountered in the workflow.
- **Parameter:** Displays the name of the node parameter that has an erroneous setting.
- **Value:** Displays the erroneous setting of the node parameter.


When a warning symbol is attached to a node, the application automatically updates the node with the correct version and saves the previous parameter values in the updated node. It does not include node parameters that are no longer available, and it adds any new parameters set to their default values.

If the Parameter and Value columns indicate a problem with the parameter settings, enter the correct parameter settings in the Parameters pane.

Saving a Workflow as a Common Template

If you want to share a changed standard template or a new custom workflow with others, you can save it as a standard template in the Common Templates folder.

❖ To save a workflow as a common template

1. To save a common workflow under a new file name or to save a new custom workflow, type a new name in the Workflow box of the Workflow Editor and a new description in the Description box.
2. Click the **Save Common** icon,  Save Common .

By default, the Save Workflow dialog box opens in the Common Templates directory. For the location of this folder, see “[Using Common Workflow Templates](#)” on [page 119](#).

3. In the Save Workflow dialog box, do the following:

- a. In the File Name box, type the file name.


If you altered a standard workflow template, the name of this file must be different from the name of the original template.

- b. Click **Save**.

Saving a Workflow as a Custom Template

You can save a new workflow as a custom template for your own use.

❖ To save a workflow as a custom template

1. Type a new name in the Workflow box of the Workflow Editor and a new description in the Description box.
2. Click the **Save** icon,  Save .

By default, the Save Workflow dialog box opens in the last directory used.

3. In the Save Workflow dialog box, do the following:

- a. If you do not want to save the workflow template in the study directory, browse to the directory where you want to save the custom template.
- b. In the File Name box, type the file name.
- c. Click **Save**.

Deleting a Workflow Template

You can delete a workflow template.

❖ To delete an existing workflow template


1. In Windows Explorer, navigate to the folder containing the template.
2. Select the template to delete, and either press the DELETE key or right-click and choose **Delete**.
3. In the Delete File dialog box, click **Yes**.

Opening a Workflow from a .pdResult File


From a .pdResult file, you can open the workflow that generated the file.

❖ To open the workflow that generated a .pdResult file

Do one of the following:

- a. If you recently generated a .pdResult file, select **Administration > Show Job Queue** to open the job queue.
- b. Select the appropriate search.
- c. Click the **Open Study** icon,  **Open Study** .

–or–

- a. In the Workflow Editor window of an open analysis, click the **Open** icon,  **Open** .
The Open Workflow dialog box opens.
- b. In the File Type box of the Open Workflow dialog box, select **Proteome Discoverer Result File (*.pdResult)**.
- c. Browse to the location of the appropriate .pdResult file, select the file, and click **Open**.

You can also obtain information about the workflow that generated a .pdResult file by using the file's Search Summary. For information on the Search Summary, refer to the Help.

Creating Specific Types of Workflows

Follow these procedures to create specific types of workflows.

- [Creating Quantification Workflows](#)
- [Creating Protein Annotation Workflows](#)
- [Creating PTM Analysis Workflows](#)

- [Creating Parallel Workflows](#)
- [Creating HCD/EThcD Workflows](#)
- [Creating Spectra Exportation Workflows](#)
- [Creating a Marked Contaminants Workflow](#)

Creating Quantification Workflows

To perform quantification, you must run quantification processing and consensus workflows. A quantification processing workflow includes one of three quantification nodes found in the quantification section of the Workflow Nodes pane of the Processing Workflow window. A quantification consensus workflow includes the Peptide and Protein Quantifier node, found in the Workflow Nodes pane of the Consensus Workflow window.

[Table 5](#) lists these nodes and where to find information about creating a quantification workflow for each.

Table 5. Quantification nodes in the processing and consensus workflows

Quantification node	Use	For more information
Processing workflow		
Precursor Ions Quantifier node	For precursor ion quantification (for example, SILAC)	See “ Creating a Processing Workflow for Precursor Ion Quantification ” on page 390.
Reporter Ions Quantifier node	For reporter ion quantification (for example, iTRAQ and TMT)	See “ Creating a Processing Workflow for Reporter Ion Quantification ” on page 397.
Precursor Ions Area Detector node	For precursor ion area detection	See “ Creating a Processing Workflow for Precursor Ion Area Detection ” on page 403.
Consensus workflow		
Peptide and Protein Quantifier node	For all types of quantification	See “ Creating a Processing Workflow for Precursor Ion Quantification ” on page 390, “ Creating a Processing Workflow for Reporter Ion Quantification ” on page 397, and “ Creating a Processing Workflow for Precursor Ion Area Detection ” on page 403.

In the processing workflow, you must attach the Precursor Ions Area Detector node or the Precursor Ions Quantifier node directly to the Event Detector node.

Creating Protein Annotation Workflows

To create a workflow that uses the ProteinCenter Annotation node to retrieve annotation information, see [“Creating a Protein Annotation Workflow”](#) on page 296. This node loads information from the GO Slim, Pfam, Entrez Gene, Ensembl Genome, and UniProt databases from ProteinCenter and installs it in the Proteome Discoverer results files.


Creating PTM Analysis Workflows

To create a workflow that enables you to focus on studying the biologically relevant post-translational modifications of proteins, you can create a PTM analysis workflow to display phosphorylation modifications in the .pdResult results. For instructions, see [“Creating a PTM Analysis Workflow with the ptmRS Node”](#) on page 351.

Creating Parallel Workflows

Parallel workflows are workflows that search the same raw data file and the same part of the spectrum but specify different criteria, different search nodes for the search, or both. The processing part of these workflows resembles the example in [Figure 90](#). You can use parallel workflows to conduct two or more searches that use two or more search engines on the same raw data and, at the same time, compare the results of these two searches. For example, you might want to search both CID and ETD data from the same raw data file to increase the chances of finding a match. CID data contains b and y ions, and ETD data contains b, c, and z ions, so the two types of data are complementary. You can also use a parallel workflow for quantification.


Figure 90. Parallel processing workflow



The following instructions show you how to create the basic parallel workflow shown in [Figure 90](#).

❖ **To create a parallel processing workflow**

1. Create or open a study and an analysis:
 - To create a study, see [“Creating a Study”](#) on [page 38](#).
 - To open an existing study, see [“Opening an Existing Study”](#) on [page 40](#).
 - To create an analysis, see [“Creating an Analysis”](#) on [page 71](#).
 - To open an existing analysis, see [“Opening an Existing Analysis”](#) on [page 72](#).
2. Follow the general instructions for creating a processing workflow in [“Creating a Processing Workflow”](#) on [page 104](#).
3. Drag the **Spectrum Files** node to the Processing Workflow Tree pane, and specify the name and path of the raw data file in the Parameters pane.
4. Drag the **Spectrum Selector** node to the Processing Workflow Tree pane and place it directly under the Spectrum Files node. Set the parameters.

5. Drag two **Scan Event Filter** nodes to the Processing Workflow Tree pane and place them side by side beneath the Spectrum Selector node. In the Parameters pane, set the Activation Type parameter to **CID** for one node and to **ETD** for the other node.
6. Drag the a search engine node such as **Sequest HT** to the Processing Workflow Tree pane and place it beneath the Scan Event Filter node set to the **CID** activation type.
7. Drag another search engine node such as **Mascot** to the Processing Workflow Tree pane and place it beneath the Scan Event Filter node set to the **ETD** activation type.
8. Drag a PSM Validator node such as **Fixed Value PSM Validator** beneath one of the search engine nodes.
9. Connect the nodes as shown in [Figure 90](#) on [page 130](#).
10. Create a consensus workflow, following the instructions in “[Creating a Consensus Workflow](#)” on [page 112](#).
11. Save the analysis. See “[Saving an Analysis](#)” on [page 80](#).
12. Save the study. See “[Saving a Study](#)” on [page 62](#).
13. Click the **Run** icon,  , in the upper right corner of the Analysis window.

The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

Creating HCD/EThcD Workflows

To improve the quality of ETD fragmentation, newer mass spectrometers such as the Orbitrap Fusion can supplement the electron transfer dissociation with HCD or CID activation. These supplemental activations are called EThcD and ETciD activations in recent instrument firmware.

- ETciD activation, which was formerly known as “supplemental activation,” produces the same fragment ions as ETD-only activation, so you do not need to change the search parameters. ETciD scans or ETD scans with “sa” in their names are simply recognized as ETD scans.
- EThcD activation generates additional b and y ions, so you must change the search parameters. These parameters are different from those for ETD- or HCD-only activation.

If the acquisition method used both HCD and EThcD activations simultaneously, the raw data file contains spectra from both activations. In these cases, the search workflow splits the scan by activation type and searches HCD and EThcD spectra separately with the appropriate settings.

A typical processing workflow for raw data files containing both HCD scans and EThcD scans that are acquired in parallel uses Scan Event Filter nodes to split the scan by activation type. The Proteome Discoverer application then searches the HCD scans with b and y ions and the EThcD scans with c and z ions and a certain number of b and y ions.

Figure 91 shows a typical processing workflow for searching a raw data file that contains both HCD and EThcD spectra. In this workflow, the Activation Type parameter of one Scan Event Filter node is set to HCD, and the same parameter for the other Scan Filter node is set to EThcD.

Figure 91. Typical search workflow for the simultaneous search of EThcD and HCD spectra





Figure 92 shows the results of the workflow presented in Figure 91.

Figure 92. Search results from the simultaneous search of HCD and EThcD spectra in the same raw data file



Creating Spectra Exportation Workflows

You can create a workflow that exports data to DTA, MGF, MZDATA, and MZML files. For information on creating this workflow, refer to the Help.

Creating a Marked Contaminants Workflow

You can mark as contaminants all proteins in the results file that are listed in a designated FASTA file or files. To create this workflow, see [“Identifying Contaminants During Searches”](#) on page 217.


Searching Multiple Sequence Databases with Mascot

If you connect the Proteome Discoverer application to Mascot server 2.3 or later, you can search more than one sequence database in the same search.

❖ To search multiple sequence databases with Mascot

1. Connect to Mascot server 2.3 or later. For instructions, see “[Configuring the Mascot Search Engine](#)” on page 25.
2. Create or open a study and an analysis:
 - To create a study, see “[Creating a Study](#)” on page 38.
 - To open an existing study, see “[Opening an Existing Study](#)” on page 40.
 - To create an analysis, see “[Creating an Analysis](#)” on page 71.
 - To open an existing analysis, see “[Opening an Existing Analysis](#)” on page 72.
3. Create a processing workflow. Follow the instructions in “[Creating a Processing Workflow](#)” on page 104.
4. For the search engine node, drag the **Mascot** node to the Workflow Tree pane in the Processing Workflow window.
5. To set the parameters for the Mascot node, select two or more databases in the Protein Database box, as shown in [Figure 93](#).

Figure 93. Selecting multiple sequence databases for the Mascot node



6. Continue with the workflow given in “Creating a Processing Workflow” on page 104.

When you run a search against multiple sequence databases, the results include proteins from both or all databases, as shown in Figure 94.

Figure 94. Result of performing a Mascot search against multiple sequence databases

	Checked	Accession	Description	ΣCoverage	Σ# Peptides	Σ# PSMs	#
1	<input type="checkbox"/>	IPI01012004.1	cDNA PSEC0175 fis, clone OVARC1000169, highly similar t	0.00 %	0	1	1
2	<input type="checkbox"/>	IPI01009528.1	cDNA PSEC0148 fis, clone PLACE1007202, highly similar t	0.00 %	0	1	1
3	<input type="checkbox"/>	P86235	Protein disulfide-isomerase A3 (Fragments) OS=Mesocricet	0.00 %	0	1	1
4	<input type="checkbox"/>	IPI00893541.1	14 kDa protein	0.00 %	0	1	1
5	<input type="checkbox"/>	IPI00847663.2	cDNA FLJ53558, highly similar to Protein disulfide-isomeras	0.00 %	0	1	1
6	<input type="checkbox"/>	P27773	Protein disulfide-isomerase A3 OS=Mus musculus GN=Pdis	0.00 %	0	1	1
7	<input type="checkbox"/>	P30101	Protein disulfide-isomerase A3 OS=Homo sapiens GN=PDIs	0.00 %	0	2	2
8	<input type="checkbox"/>	P11598	Protein disulfide-isomerase A3 OS=Rattus norvegicus GN=I	0.00 %	0	1	1
9	<input type="checkbox"/>	Q4V1	Protein disulfide-isomerase A3 OS=Chlorocebus aethiops C	0.00 %	0	1	1

Creating a Multiconsensus Report

A multiconsensus report is a .pdResult file generated from a workflow that includes multiple search engines or multiple input files. For the input files, you can use any of the input files listed in “Inputs and Outputs” on page 11 or MSF files.

The .pdResult report displays the combined results in columns that are labeled by search engine type, for example, Score Mascot and Score Sequest HT.

- On the Proteins page, the #PSMs, # Peptides, Coverage, and Score columns display information side by side for each protein. When the proteins referenced in the individual searches have different accession numbers, the application displays the accession number of the first search in the result set.
- On the PSMs page, the Search ID column displays the order in which the search was submitted to the job queue. A multiconsensus report resulting from a workflow that includes both Sequest HT and Mascot contains XCorr (Sequest HT) or Ions Score (Mascot) columns. These columns score the number of fragment ions that are common to two different peptides with the same precursor mass and calculate the cross-correlation score for all candidate peptides queried from the database. The default setting shows the top matches per peptide and search engine.

The identification information from the individual search nodes appear side by side on the PSMs and Peptide Groups pages. The grouped peptides do not represent actual matches found during the search but rather the unified information from all matches found for this particular peptide sequence.

You can access the peptide information by examining both the Proteins page for all peptides associated with a protein and the PSMs page for all peptides, including those not associated with any protein. You can view the associated peptides on the Proteins page. Select a protein row and click **Show Associated Tables** at the bottom of the page.

To filter multiconsensus reports for specific results, see [“Finding Common or Unique Proteins in Multiple Searches”](#) on page 272 and [“Applying Filters Specific to Different Searches in Multiconsensus Reports”](#) on page 276.

❖ **To create a multiconsensus report from a search using multiple input files**


See [“Adding Input Files”](#) on page 48.

❖ **To create a multiconsensus report by using multiple search engines**

1. Create or open a study and an analysis:
 - To create a study, see [“Creating a Study”](#) on page 38.
 - To open an existing study, see [“Opening an Existing Study”](#) on page 40.
 - To create an analysis, see [“Creating an Analysis”](#) on page 71.
 - To open an existing analysis, see [“Opening an Existing Analysis”](#) on page 72.
2. Add the appropriate input file or files.
3. Create a processing workflow by following the instructions in [“Creating a Processing Workflow”](#) on page 104.

In the workflow, include at least two of the search engine nodes listed under Sequence Database Search in the Workflow Nodes pane.


Figure 95. Workflow with multiple search engines




4. Create a consensus workflow by following the instructions in [“Creating a Consensus Workflow”](#) on page 112.
5. For multiconsensus reports that include multiple occurrences of the same type of search nodes, set the Merge Mode parameter of the MSF Files node, as shown in [Figure 96](#), to determine how the application merges the data. Use one of the following settings:
 - Globally by Search Engine Type: See [“Globally by Search Engine Type Setting”](#) on page 139.
 - Per File and Search Engine Type: See [“Per File and Search Engine Type Setting”](#) on page 140.
 - Do Not Merge: See [“Do Not Merge Setting”](#) on page 141.

For more information on the MSF Files node, refer to the Help.

Figure 96. Setting the Merge Mode parameter of the MSF Files node



6. Click the **Run** icon, , in the upper right corner of the Analysis window.

The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

Merging Results by Search Engine Node

When you generate a .pdResult file from a workflow containing multiple search engines, you can merge the search results of each search engine type (such as Sequest HT or Mascot) to prevent too many result columns from appearing in the report. In addition, you might want to divide the results by search engine—for example, if someone wants to review the differences between instrument methods used in the different files or compare search engine settings.

You might want to merge search results in the following cases:

- The application creates a report from many MSF files (for example, 50), and each contains at least one search node.

- The application creates a report from only a few MSF files (for example, three), but each workflow contains 10 search nodes.

In both cases, you might prefer to reduce the number of data columns to obtain a clear overview of the results.

When you create a results report generated by a workflow that contains MSF files generated by several search nodes, you can use the following settings for the Merge Mode parameter of the MSF Files node to display the search-engine-specific data on the Protein Groups and the Proteins pages:

- Globally by Search Engine Type
- Per File and Search Engine Type setting
- Do Not Merge

The following topics illustrate the Merge Node parameter settings. They are based on the processing of two files with a processing workflow that contains one Mascot node and two Sequest HT nodes that use different sequence databases.

For more information on the MSF Files node, refer to the Help.

Globally by Search Engine Type Setting

The default setting of the Merge Mode parameter, Globally by Search Engine Type, merges the results of all search nodes of the same type. It displays the best PSM values on the Peptide Groups page and shows the combined protein score that was calculated from all PSMs identified by the search engine type.

[Figure 97](#) displays the merged results of all Sequest HT and Mascot nodes and one protein score by node type on the Peptide Groups page. The results from the two different search engine types appear together in two different column groups. The top illustration shows them grouped on the Peptide Groups page, and the bottom illustration shows them grouped on the Proteins page.

2 Getting Started

Performing a Search

Figure 97. Results grouped by search engine type on the Peptide Groups page and the Proteins page

Protein Groups										Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info
#	Checked	Confidence	Sequence	# Protein Groups	# Proteins	XCorr Sequest HT	Confidence Sequest HT	Ions Score Mascot	Confidence Mascot				
1	<input type="checkbox"/>		ILIANTGMDTDK	1	7	2.89		54					
2	<input type="checkbox"/>		ESEDKPEIEDVGSDEEEK	1	6	3.46		50					
3	<input type="checkbox"/>		HIYYITGETK	1	15	2.62		55					
4	<input type="checkbox"/>		TAVVVGITDDVR	1	14	2.93		58					
5	<input type="checkbox"/>		IDNDGDGFVTTEELK	1	4	2.64		83					
6	<input type="checkbox"/>		AGTGVDNVLEAATR	1	10	3.62		57					
7	<input type="checkbox"/>		GAEAAVVTGGVVPVQGSK	1	8	4.32		96					
8	<input type="checkbox"/>		QDLPNAMNAEITDK	1	22	3.13		64					
9	<input type="checkbox"/>		SQIFSTASDNQPTVIK	1	9	3.09		54					

Protein Groups										Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info
#	Checked	Master	Accession	# Peptides	# PSMs	# AAs	Score Sequest HT	# Peptides Sequest HT	Score Mascot	# Peptides Mascot			
1	<input type="checkbox"/>		P08238	18	155	724	139.81	16	660	15			
2	<input type="checkbox"/>	<input type="checkbox"/>	B4DMA2	16	84	686	136.22	16					
3	<input type="checkbox"/>	<input type="checkbox"/>	B4DGL0	15	76	714	122.43	15					
4	<input type="checkbox"/>		P07900	15	128	732	112.43	13	529	14			
5	<input type="checkbox"/>	<input type="checkbox"/>	Q2VPJ6	12	63	585	103.56	12					
6	<input type="checkbox"/>		Q71U34	11	94	646	100.21	10	345	9			
7	<input type="checkbox"/>	<input type="checkbox"/>	Q53G26	10	56	646	100.21	10					
8	<input type="checkbox"/>	<input type="checkbox"/>	Ensembl	10	74	627	100.21	10					
9	<input type="checkbox"/>		P31327	14	92	1500	94.46	14	292	8			
10	<input type="checkbox"/>	<input type="checkbox"/>	Ensembl	14	76	1506	94.46	14					
11	<input type="checkbox"/>	<input type="checkbox"/>	A8K134	14	58	1500	94.46	14					

Per File and Search Engine Type Setting

The Per File and Search Engine Type setting of the Merge Mode parameter merges the results of all search engine nodes of the same type for each input file.

In [Figure 98](#), the results report displays one group of search-engine-specific values per file and the search engine type. There are two groups of Sequest HT columns and two groups of Mascot columns. The top illustration shows the data grouped on the Peptide Groups page, and the bottom illustration shows it grouped on the Proteins page.

Figure 98. Results grouped by search engine type for each input file on the Peptide Groups page and the Proteins page

Protein Groups										Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info
#	Checked	Confidence	Sequence	# Protein Groups	# Proteins	XCorr A Sequest HT	Confidence A Sequest HT	Ions Score A Mascot	Confidence A Mascot	XCorr B Sequest HT	Confidence B Sequest HT	Ions Score B Mascot	Confidence B Mascot
1	<input type="checkbox"/>		ILIANTGMDTDK	1	7	2.89		54		2.89		54	
2	<input type="checkbox"/>		VGLQVAVK	1	12	3.35		70		3.35		70	
3	<input type="checkbox"/>		VDNDEHQLSLR	1	8	3.91		83		3.91		83	
4	<input type="checkbox"/>		VDATAETDLAK	1	1	2.45		69		2.45		69	
5	<input type="checkbox"/>		EDMAALEK	2	32	2.60		54		2.60		54	
6	<input type="checkbox"/>		TAFQALDAAGDK	1	4	3.39		85		3.39		85	
7	<input type="checkbox"/>		ISSDLGHPVPK	1	12	4.10		60		4.10		60	
8	<input type="checkbox"/>		SVTEGGAELSNEER	1	11	3.85		94		3.85		94	
9	<input type="checkbox"/>		GATQQLDEAER	1	10	3.53		67		3.53		67	
10	<input type="checkbox"/>		GYDVIQAQSGTGK	1	29	2.94		71		2.94		71	

Protein Groups										Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info	
#	Checked	Master	Accession	# Peptides	# PSMs	# AAs	Score A Sequest HT	# Peptides A Sequest HT	Score A Mascot	# Peptides A Mascot	Score B Sequest HT	# Peptides B Sequest HT	Score B Mascot	# Peptides B Mascot
1	<input type="checkbox"/>		P08238	18	155	724	127.69	16	572	14	139.81	16	641	14
2	<input type="checkbox"/>	<input type="checkbox"/>	B4DMA2	16	84	686	124.10	16			136.22	16		
3	<input type="checkbox"/>	<input type="checkbox"/>	B4DGL0	15	76	714	110.32	15			122.43	15		
4	<input type="checkbox"/>		P07900	15	128	732	103.04	13	501	13	112.43	13	496	12
5	<input type="checkbox"/>	<input type="checkbox"/>	Q2VPJ6	12	63	585	94.17	12			103.56	12		
6	<input type="checkbox"/>		P63261	7	83	375	85.99	7	232	4	84.04	6	219	4
7	<input type="checkbox"/>	<input type="checkbox"/>	Q4R561	7	83	375	85.99	7	232	4	84.04	6	219	4
8	<input type="checkbox"/>	<input type="checkbox"/>	B4DVQ0	7	53	333	85.99	7			84.04	6		

Do Not Merge Setting

The Do Not Merge setting of the Merge Mode parameter does not merge nodes. It displays the search-engine-specific values for each search engine. The File Limit for Automatic Merge parameter of the MSF Files Node sets an upper limit on the number of files that the report can display for unmerged data.

Figure 99 shows four groups of Sequest HT columns and two groups of Mascot columns. The top illustration shows the data grouped on the Peptide Groups page, and the bottom illustration shows it grouped on the Proteins page.

Figure 99. Results in one group of search-engine-specific values per search node on the Peptide Groups page and the Proteins page

Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info														
#	Checked	Confidence	Sequence	# Protein Groups	# Proteins	XCorr A3 Seq	Confidence A3 Seq	XCorr A4 Seq	Confidence A4 Seq	Ions Score A0	Confidence A0 Seq	XCorr B3 Seq	Confidence B3 Seq	XCorr B4 Seq	Confidence B4 Seq	Ions Score B0 Ma	Confidence B0 Ma		
1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LIANTGMDTDK	1	7	2.89	<input checked="" type="checkbox"/>	2.89	<input checked="" type="checkbox"/>	54	<input checked="" type="checkbox"/>	2.89	<input checked="" type="checkbox"/>	2.89	<input checked="" type="checkbox"/>	54	<input checked="" type="checkbox"/>		
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SVTEGGAELSNEER	1	11	3.85	<input checked="" type="checkbox"/>	3.85	<input checked="" type="checkbox"/>	94	<input checked="" type="checkbox"/>	3.85	<input checked="" type="checkbox"/>	3.85	<input checked="" type="checkbox"/>	94	<input checked="" type="checkbox"/>		
3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	GATGQLDEAER	1	10	3.53	<input checked="" type="checkbox"/>	3.53	<input checked="" type="checkbox"/>	67	<input checked="" type="checkbox"/>	3.53	<input checked="" type="checkbox"/>	3.53	<input checked="" type="checkbox"/>	67	<input checked="" type="checkbox"/>		
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	GYDVIAGQSGTGK	1	29	2.94	<input checked="" type="checkbox"/>	2.94	<input checked="" type="checkbox"/>	71	<input checked="" type="checkbox"/>	2.94	<input checked="" type="checkbox"/>	2.94	<input checked="" type="checkbox"/>	71	<input checked="" type="checkbox"/>		
5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	VDAQFGQSDQGR	1	2	2.68	<input checked="" type="checkbox"/>	2.68	<input checked="" type="checkbox"/>	53	<input checked="" type="checkbox"/>	2.68	<input checked="" type="checkbox"/>	2.68	<input checked="" type="checkbox"/>	53	<input checked="" type="checkbox"/>		
6	<input type="checkbox"/>	<input checked="" type="checkbox"/>	VEIINDQGQNR	7	151	2.63	<input checked="" type="checkbox"/>	2.63	<input checked="" type="checkbox"/>	58	<input checked="" type="checkbox"/>	2.63	<input checked="" type="checkbox"/>	2.63	<input checked="" type="checkbox"/>	58	<input checked="" type="checkbox"/>		
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	HFVALSTNTTK	1	9	3.13	<input checked="" type="checkbox"/>	3.13	<input checked="" type="checkbox"/>	62	<input checked="" type="checkbox"/>	3.13	<input checked="" type="checkbox"/>	3.13	<input checked="" type="checkbox"/>	62	<input checked="" type="checkbox"/>		
8	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LATQLTGPVMPVIR	1	5	2.97	<input checked="" type="checkbox"/>	2.97	<input checked="" type="checkbox"/>	68	<input checked="" type="checkbox"/>	2.97	<input checked="" type="checkbox"/>	2.97	<input checked="" type="checkbox"/>	68	<input checked="" type="checkbox"/>		
9	<input type="checkbox"/>	<input checked="" type="checkbox"/>	GVEEEEEDGEMRE	1	1	3.09	<input checked="" type="checkbox"/>	3.09	<input checked="" type="checkbox"/>	51	<input checked="" type="checkbox"/>	3.09	<input checked="" type="checkbox"/>	3.09	<input checked="" type="checkbox"/>	51	<input checked="" type="checkbox"/>		
10	<input type="checkbox"/>	<input checked="" type="checkbox"/>	ELGSSVALYSR	1	12	2.56	<input checked="" type="checkbox"/>	2.56	<input checked="" type="checkbox"/>	53	<input checked="" type="checkbox"/>	2.56	<input checked="" type="checkbox"/>	2.56	<input checked="" type="checkbox"/>	53	<input checked="" type="checkbox"/>		

Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info														
#	Checked	Master	Accession	# Peptides	# PSMs	# AAs	Score A3 Seq	# Peptides A3 Seqs	Score A4 Seqs	# Peptides A4 Seqs	Score A0 Masc	# Peptides A0 Ma	Score B3 Sequet	# Peptides B3 Seq	Score B4 Sequet H	# Peptides B4 Seq	Score B0 Masc	# Peptides B0 Ma	
1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B4DMA2	16	84	686	124.10	16					136.22	16					
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	P08238	18	155	724	124.10	16	49.15	11	572	14	136.22	16	51.96	12	641	14	
3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B4DGL0	15	76	714	110.32	15					122.43	15					
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	P07900	15	128	732	99.45	13	42.69	10	501	13	108.85	13	45.50	11	496	12	
5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q2VPJ6	12	63	585	94.17	12					103.56	12					
6	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q63G76	7	53	375	85.99	7					84.04	6					
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q63GK6	7	53	375	85.99	7					84.04	6					
8	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B4DVQ0	7	53	333	85.99	7					84.04	6					
9	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q8WVW5	7	53	363	85.99	7					84.04	6					
10	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q63G99	7	53	375	85.99	7					84.04	6					
11	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q4R61	7	83	375	85.99	7	29.73	4	232	4	84.04	6	29.73	4	219	4	

2 Getting Started

Performing a Search

Using the Proteome Discoverer Daemon Utility

This chapter describes the Proteome Discoverer Daemon utility, which you can use to perform a search either locally or remotely by using its own application window, the command line, or a parameter file that contains command-line commands. You can use it to start workflows, monitor job execution on the configured server, perform batch processing, and process Multidimensional Protein Identification Technology (MudPIT) samples. The Proteome Discoverer Daemon utility can perform multiple searches on multiple raw data files taken from multiple samples or from one sample.

Contents

- [Starting the Proteome Discoverer Daemon Utility in a Window](#)
- [Selecting the Server](#)
- [Running the Proteome Discoverer Daemon Utility from the Window](#)
- [Running the Proteome Discoverer Daemon Utility from a Parameter File](#)
- [Monitoring Job Execution in the Proteome Discoverer Daemon Utility](#)
- [Logging in to a Remote Server](#)
- [Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System](#)
- [Running the Proteome Discoverer Daemon Utility from the Command Line](#)

For information about MudPIT and creating a MudPIT workflow, see [“Adding Input Files”](#) on [page 48](#).


Starting the Proteome Discoverer Daemon Utility in a Window

You can start the Proteome Discoverer Daemon utility on the command line or in a window. To run it on the command line, see [“Running the Proteome Discoverer Daemon Utility from the Command Line”](#) on [page 165](#).

3 Using the Proteome Discoverer Daemon Utility

Selecting the Server

❖ To start the utility in a window

1. Start the utility in Windows by choosing **Start > All Programs > Thermo Proteome Discoverer *release_number* > Proteome Discoverer Daemon *release_number*** or by clicking the **Daemon** icon, , on your desktop.
2. After the utility window appears, connect to a computer that is running the Proteome Discoverer application.

Selecting the Server

The Proteome Discoverer Daemon utility can connect to a remote computer running the Magellan server so that you can perform searches on multiple raw data files from multiple samples or a single sample on a remote computer. It can also connect to a local server.

❖ To specify the server to connect to


1. Click the **Configuration** tab in the Proteome Discoverer Daemon utility window.
2. From the Host list, select the name of the server that you want to use, or type the server name.

You must connect the utility to a computer running the Magellan server. Your local host is the default server, that is, the computer that you are working on. To connect remotely to the Magellan server, see [“Logging in to a Remote Server”](#) on [page 149](#).

3. In the User box, type the user’s login name on the server.

The Configuration page now resembles [Figure 100](#).

Figure 100. Configuration page of the Proteome Discoverer Daemon utility



4. Click **Apply** to activate the new settings.
5. To return to the previous settings, click **Reset**.

Running the Proteome Discoverer Daemon Utility from the Window

You can use the Proteome Discoverer Daemon utility to run processing and consensus workflows on a remote server. These workflows are saved in templates that you created in the Workflow Editor. You must set the values of the required parameters for the workflow templates before you run the utility. The specified values—for example, FASTA databases—must be available on the server where they will be processed. You cannot change the workflows by using the Proteome Discoverer Daemon interface.


You can choose quantification methods from predefined quantification methods on the server or from a method file that you load with the Browse button.

❖ To start a workflow for batch processing or MudPIT processing

1. Click the **Start Jobs** tab.

The Start Jobs page appears, as shown in [Figure 101](#).

Figure 101. Start Jobs page of the Proteome Discoverer Daemon utility



2. Click the **Load Files** tab if it is not already selected.
3. Click **Add**.
4. In the Open dialog box, locate the folder that contains the raw data files, select the spectrum (RAW) files that you want to load, and click **Open**.

The selected spectrum files appear on the Load Files page.

To remove a file from the Load Files page, select the file and click **Remove**.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Window

5. To specify the type of processing, select the **Batch Processing** or **MudPIT** option.
 - Batch processing (the default): Executes the workflow once for each spectrum file.
 - MudPIT: Feeds all spectrum files into one workflow.

When you select the MudPIT option, the Treat as Replicates check box and the Output Filename box become available.

6. (Optional) For the MudPIT option, select the **Treat as Replicates** check box to group the samples by quantification channel and to average the abundance values for each channel across files.

This option is only available when you select the MudPIT option.

If you do not select the Treat as Replicates option, the application groups the samples by quantification channel and by file. It creates the sample groups and ratios only within files, not across files.

7. Select a processing template from template files located on the Daemon server.
 - In the Workflows and Quan Method area of the Start Jobs page, do one of the following:

Click the **Browse** button (...), select the processing template from the Open dialog box, and click **Open**.

–or–

Select a processing workflow template from the Processing menu.

The menu box lists only the templates that you selected in the current run. When you restart the application, the menu box does not list any templates. When you connect the application to another server, the application verifies that the templates are still valid and can be run by the connected server.

The template files for processing workflows have a .pdProcessingWF extension.

8. Select a consensus template.

- In the Workflows and Quan Method area of the Start Jobs page, do one of the following:

Click the **Browse** button (...), select a consensus workflow template from the Open dialog box, and click **Open**.

–or–

Select a consensus workflow template from the Consensus menu.

The menu box lists only the templates that you selected in the current run. When you restart the application, the menu box does not list any templates. When you connect the application to another server, the application verifies that the templates are still valid and can be run by the connected server.


The template files for consensus workflows have a .pdConsensusWF extension.

9. To obtain quantification results, select a quantification method from the Quan Method menu in the Workflows and Quan Method area of the Start Jobs page.
10. (Optional) To specify the control channel to use in creating ratios when you process quantification samples, select the control channel from the Control Channel list.
11. If you selected the MudPIT option, in the Output Filename box, type the name of the output file where you want to store the results of the search.

The Start Jobs page should now resemble [Figure 102](#) for batch processing or [Figure 103](#) for MudPIT processing.

12. Click **Start** to execute the job.


Figure 102. Start Jobs page of the Proteome Discoverer Daemon utility for batch processing



3 Using the Proteome Discoverer Daemon Utility

Monitoring Job Execution in the Proteome Discoverer Daemon Utility

Figure 103. Start Jobs page of the Proteome Discoverer Daemon utility for MudPIT processing



Monitoring Job Execution in the Proteome Discoverer Daemon Utility


You can use the Job Queue page in the Proteome Discoverer Daemon utility window to monitor the execution of the jobs that you submit. It performs the same function as the job queue on the Proteome Discoverer interface. For information about the features of the job queue on the Proteome Discoverer interface, refer to the Help.

A progress bar displays the progress of the overall batch processing. This progress bar is visible only if you have started batch jobs.

❖ To monitor the job execution

Click the **Job Queue** tab of the Proteome Discoverer Daemon utility window.


Figure 104 shows the completed job for batch processing.

Figure 104. Job Queue page of the Proteome Discoverer Daemon utility for batch processing


Logging in to a Remote Server

The searches started by the Proteome Discoverer application consume memory and can potentially cause the data-acquiring computer to crash and lose the sample in the mass spectrometer. To avoid this outcome, Thermo Fisher Scientific recommends that you connect the Proteome Discoverer Daemon utility to a remote computer running the Magellan server before data acquisition.

❖ To log in to a remote server

1. Start the Proteome Discoverer application on the remote machine.
2. If you want to store the output files in a location other than the default folder, do the following:
 - a. Choose **Administration > Configuration > Server Settings > Discoverer Daemon**.
The PublicFiles folder is the default folder displayed in the Current File Directory box, as shown in [Figure 105](#).
 - b. In the New Directory box, browse to the location where you want to store the output files.
 - c. Click  **Apply** .


The Proteome Discoverer Daemon utility issues a message informing you that the Proteome Discoverer application will apply the change the next time that you start it.

To return to the default directory, click  **Reset** .

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

Figure 105. Proteome Discoverer Daemon area of the Configuration view



3. Start the Proteome Discoverer Daemon utility on the local machine.

A message box informs you that the utility cannot connect to the server.

4. Click **OK** in the message box.

The Proteome Discoverer Daemon utility opens with the Configuration page selected.

5. In the Host box, type the name of the remote computer.

6. In the User box, type the login name of the remote server.

7. Click **Apply**.

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

You can use the parameter file created in the Proteome Discoverer Daemon utility to call the application from the Xcalibur data system. You can start the Proteome Discoverer Daemon utility by adding a parameter file to the processing method specified in the Xcalibur injection sequence.

These topics describe how to run the utility from the Xcalibur data system:

- [Before You Start](#)
- [Running the Proteome Discoverer Daemon Utility from a Parameter File](#)
- [Creating a Processing Method That Calls the Proteome Discoverer Daemon Utility](#)
- [Batch Processing with a Processing Method That Calls the Proteome Discoverer Daemon Utility](#)
- [Batch Processing with Multiple Processing Methods](#)
- [Processing MudPIT Samples by Using a Processing Method](#)

Before You Start

Before you run the Proteome Discoverer Daemon utility from the Xcalibur data system, perform the following steps to set up the interface between the Proteome Discoverer Daemon utility and the Xcalibur data system.

❖ **To prepare to run the Proteome Discoverer Daemon utility from the Xcalibur data system**

1. Before you start the Proteome Discoverer Daemon utility, install the Proteome Discoverer application on a remote computer to decouple data processing from data acquisition.


Thermo Fisher Scientific strongly recommends that you perform data analysis and data recording on two different computers to avoid disturbing the data acquisition by resource-consuming data processing.

2. Install the Proteome Discoverer Daemon utility on the same computer where the Xcalibur data system is running.
3. Start the Proteome Discoverer application.
4. In the Proteome Discoverer application, prepare the workflow to be used by the Proteome Discoverer Daemon utility, as shown in [Figure 106](#). Save this workflow.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

Figure 106. Basic workflow used for the samples




After you install the Proteome Discoverer Daemon utility, the Proteome Discoverer application places the directory where it saves the raw data files and stores the results in the following folder:

Drive:\ProgramData\Thermo\Proteome Discoverer <release_number>\PublicFiles

This directory might be invisible to you because the *drive:\Documents and ProgramData* directory is hidden. To display hidden directories, choose **Tools > Folder Options > View > Hidden files and folders > Show hidden files and folders** in Windows Explorer.

5. (Optional) To change this directory for easier data access, open the Proteome Discoverer application, choose **Administration > Configuration**, click **Discoverer Daemon** beneath Server Settings in the Configuration area on the left side of the Administration view, and change the directory in the New Directory box, shown in [Figure 107](#).

The settings are applied after you restart the Proteome Discoverer application.

Figure 107. Change in the destination directory where results from the Proteome Discoverer Daemon utility are stored

Running the Proteome Discoverer Daemon Utility from a Parameter File

In the Proteome Discoverer Daemon utility, you can create a parameter file that you can use to call the application from the Xcalibur data system. The application automatically translates the options that you set on the utility interface and in the workflow used for the search into text commands in the parameter file.


❖ To create a parameter file that calls the Proteome Discoverer Daemon utility

1. Set up the search according to the instructions in [“Running the Proteome Discoverer Daemon Utility from the Window”](#) on [page 145](#). You do not have to have files loaded to create a parameter file.
2. Click the **Export Parameter File** tab, shown in [Figure 108](#), on the Start Jobs page.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

Figure 108. Export Parameter File page



3. In the Number of Rawfiles box for a MudPIT search, select the number of files that will appear in the Xcalibur Sequence Setup dialog box.

The Number of Rawfiles option is not available when you select the Batch Processing option.

4. Click **Export**.

The Save a Parameter File dialog box appears.

5. Specify the path and name of the parameter file, and click **Save**.

The Proteome Discoverer application writes the parameter file in XML format to the specified directory.

To call the Proteome Discoverer Daemon utility through the parameter file, see [“Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System”](#) on page 150.

Creating a Processing Method That Calls the Proteome Discoverer Daemon Utility

The following procedure describes how to create a processing method that calls the Proteome Discoverer Daemon utility. It assumes that you have already created an appropriate processing method for your raw data files. Processing methods have a .pmd file extension.

❖ **To add a processing method that calls the Proteome Discoverer Daemon utility to a processing method**

1. Choose **Start > All Programs > Thermo Xcalibur > Xcalibur** to start the Xcalibur data system.

The Roadmap view of the Xcalibur Home Page window opens.

2. In the Roadmap view, do one of the following:

- Choose **GoTo > Processing Setup**.

The Processing Setup window opens.

3. Open the processing method that you want to modify as follows:


- a. Choose **File > Open**.
- b. Browse to the location of the processing method file and select the file.
- c. Click **Open**.

The selected processing method opens in the Processing Setup window.

4. Open the Programs view of the Processing Setup window as follows:

- a. Choose **View > View Bar**.

The view bar appears on the left side of the dialog box.

- b. On the view bar, click the **Programs** icon, .

The Programs view of the Processing Setup window opens, as shown in [Figure 109](#).

Figure 109. Programs view with an empty table



Programs:									
----- Sample type -----									
	Enable	Std	QC	Unk	Other	Action	Program or Macro Name	Sync	Parameters
*	<input type="checkbox"/>	Yes	Yes	Yes	Yes	Run Program		Yes	

5. If the Programs view contains an empty table, right-click the table and choose **Insert Row** from the shortcut menu.

A new row appears above the placeholder row, as shown in [Figure 110](#). An asterisk to the left side of a table row defines the row as a placeholder row.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

Figure 110. Programs view with an unedited table row

	Enable	Std	QC	Unk	Other	Action	Program or Macro Name	Sync	Parameters
1	<input checked="" type="checkbox"/>	Yes	Yes	Yes	Yes	Run Program		Yes	
*		Yes	Yes	Yes	Yes	Run Program		Yes	

6. In the added table row, specify the name and location of the parameter file as follows:
 - a. In the Enable column, select the check box.
 - b. In the Action list column, select **Run Program**.
 - c. Right-click the **Program or Macro Name** column and choose **Browse** from the shortcut menu, as shown in Figure 111.

Figure 111. Programs view with the shortcut menu displayed

	Enable	Std	QC	Unk	Other	Action	Program or Macro Name	Sync	Parameters
1	Yes	Yes	Yes	Yes	Yes	Run Program			
*		Yes	Yes	Yes	Yes	Run Program			

The Browse for Program dialog box opens.

- d. Browse to the following executable, and click **Open**:

C:\Program Files\Thermo\Proteome Discoverer Daemon 2.1\System\Release\System\Release\DiscovererDaemon.exe

Note If the following warning appears, click **OK**:

The file 'DiscovererDaemon' does not exist on this computer.

- e. In the Parameters column, type the location of the parameter file containing the commands that will execute the Proteome Discoverer Daemon utility:

-p *path_to_parameter_file\parameter_filename* %R

IMPORTANT If the name of the parameter file contains a space, you must enclose the name in quotation marks, as in this example:

-p "C:\Xcalibur\methods\batch processing.param" %R

7. In the Std, QC, Unk, Other, and Sync columns, accept the default settings or modify them according to your requirements. For information about setting the sample types to

be sent to the Proteome Discoverer Daemon utility, see [“To specify the sample types to be sent to the Proteome Discoverer Daemon utility.”](#)

To send all sample types to the Proteome Discoverer Daemon utility, make sure that all of the sample type columns are set to **Yes**, as shown in [Figure 112](#).

Figure 112. Program table with a call to the Proteome Discoverer Daemon utility

Programs:									
----- Sample type -----									
	Enable	Std	QC	Unk	Other	Action	Program or Macro Name	Sync	Parameters
1	Yes	Yes	Yes	Yes	Yes	Run Program	C:\Program Files\Thermo\Discoverer	Yes	-p "C:\Daemon\
*		Yes	Yes	Yes	Yes	Run Program		Yes	

-p "C:\Daemon\data\daemon.param" %R

8. Click **OK** to save the changes to the processing method.

9. Choose **File > Save**.

❖ **To specify the sample types to be sent to the Proteome Discoverer Daemon utility**

1. If the processing method that you want to modify is not open, open it and make sure that the parameter file and its location are specified as described in [“To add a processing method that calls the Proteome Discoverer Daemon utility to a processing method,”](#) on page 155.
2. In the Std, QC, Unk, and Other columns, do the following:
 - To send a sample to the Daemon application, make sure that “Yes” appears in the column for its sample type.
 - To avoid processing a sample with the Proteome Discoverer Daemon utility, clear the column for its sample type.

Tip Use the Other column for the Blank sample type. For example, if you do not want to send blank samples to the Proteome Discoverer Daemon utility for further processing, clear the Other column.

3. Save the processing method.

Batch Processing with a Processing Method That Calls the Proteome Discoverer Daemon Utility

To inject samples and to acquire and process data files with the Xcalibur data system, you must create one or more instrument methods, one or more processing methods, and a sequence that defines the sample injection set.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

For information about creating an instrument method for your LC/MS system, refer to the Help for the LC devices and the Help for the mass spectrometer. For information about creating processing methods and sequences, refer to the Xcalibur Help.

Tip For a typical LC/MS experiment, an autosampler automates the sample injection process, and the position nomenclature depends on the autosampler tray type.

For information about specifying the autosampler tray type and the position nomenclature for the specified tray type, refer to the Help for the autosampler.

For some autosamplers, you can change the tray type from the Sequence Setup view by choosing Change > Tray Name and then selecting a different tray type.


To start the Proteome Discoverer Daemon utility from the Xcalibur data system version 2.10 or later, you must add a processing method that calls the application to the sequence.

❖ To set up and run an injection sequence with a processing method that starts the Proteome Discoverer Daemon utility

1. From the Home Page window of the Xcalibur data system, do one of the following:

- Click the **Sequence View** icon, , on the Home Page window toolbar.

–or–

- Click the **Sequence Setup** icon, , on the Roadmap view.

The Sequence Setup view opens with an empty sequence table. Refer to the Xcalibur – Sequence Setup view Help for information about filling out the sequence table.

2. In the Proc Meth column, select a processing method with a parameter file that calls the Proteome Discoverer Daemon utility as follows:


- Type the file location and name of the processing method.

–or–

- Double-click the column to open the Select Processing Method dialog box, where you can browse to and select the processing method.

You can now start the sequence without first saving it, or you can save the sequence for later use.


3. In the sequence table, select the row or rows that you want to run.

4. Choose **Actions > Run Sequence** or click the **Run Sequence** icon, .

If you have changed the instrument configuration in Foundation platform after the previous sequence run, the Change Instruments In Use dialog box opens. Otherwise, the Run Sequence dialog box opens, as shown in [Figure 113](#).

For an LC/MS system, the autosampler (or device with an autosampler) is specified as the start instrument. When the autosampler makes an injection, it triggers the mass spectrometer to begin data acquisition.

Figure 113. Run Sequence dialog box



5. Select the Programs check box, if it is not already selected.

You must select this check box to start the Proteome Discoverer Daemon utility.

6. Click **OK**.

If you have not already saved the sequence, the File Summary Information dialog box opens.

7. Save the sequence as follows:
 - a. In the File Summary Information box, click **OK**.
 - b. In the File Name box, type a unique name for the sequence.
 - c. In the Save In list, select the appropriate folder location for the sequence.
 - d. Click **Save**.

The Xcalibur data system adds the sequence to the acquisition queue.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

For each sequence row, after the data system acquires a raw data file, it sends the processing method and the raw data file to the Proteome Discoverer application, which stores the raw data file and the .pdResult file in the server output directory specified in the Server Output Directory box of the Export Parameter File page of the Start Jobs page. All the search results of the batch processing are stored in the same directory. If the same directory name is used for the results of another batch process, the application appends the date and an incremental index number to the folder name.

Batch Processing with Multiple Processing Methods

In some cases, you might need to use more than one processing method in the sequence. For example, the `sequest.pmd` method runs the Proteome Discoverer Daemon utility with a parameter file containing a simple Sequest HT workflow, and the `export.pmd` method runs the utility with an export workflow.

❖ To use more than one processing method in a sequence

1. In the Sequence Setup view, choose **File > New**.
The New Sequence Template dialog box opens.
2. Enter the appropriate values in each of the boxes.
3. In the Bracket Type area, select the **None** option, as shown in [Figure 114](#).

With this bracket type, you can change the processing methods individually for each sample.

Figure 114. New Sequence Template with the selection of None for the bracket type




Figure 115 shows a sequence using two different processing methods.

Figure 115. Sequence with two different processing methods

	Sample Type	File Name	Sample ID	Path	Inst Meth	Proc Meth	Position	Inj Vol
1	Unknown	BSA1	1	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\batch	A1	10.00
2	Unknown	BSA2	2	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\batch_export	A2	10.00
3	Unknown	BSA3	3	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\batch	A3	10.00
*								0.00

4. Click **OK**.

In this example, the Xcalibur data system starts two different workflows (performing a Sequest HT search and exporting a raw data file) for the recorded raw data files in the Proteome Discoverer application, as shown in Figure 116.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

Figure 116. Two workflows in the job queue started by two different processing methods

The screenshot shows the 'Discoverer Daemon' application window with the 'Job Queue' tab selected. The interface includes a toolbar with 'Pause', 'Resume', 'Abort', 'Remove', 'Refresh', and 'Display Verbose Messages' options. The main area displays a table of jobs in the queue, with two jobs highlighted in green, indicating they are completed.

Execution State	Progress	Type	Name	Submitted at	Study
Completed	100 %	Consensus	BSA3	1/16/2015 2:40 PM	BSA
Completed	100 %	Processing	BSA3	1/16/2015 2:40 PM	BSA

Below the job list, a detailed log for the BSA3 job is shown, including the following messages:

- 2:40 PM ProcessingJob ----- Total Job execution took: 4.5 s. -----
- 2:40 PM ProcessingJob Finished D:\Work\BSA\BSA3.msf
- 2:40 PM (4):Fixed ValuePSM Valida... -- Total execution of Fixed Value PSM Validator (4) took 0.3 s --
- 2:40 PM (4):Fixed ValuePSM Valida... Evaluating peptides of Sequest HT (2) started
- 2:40 PM (2):Sequest HT -- Total search time was 1.9 s --
- 2:40 PM (2):Sequest HT Stored 150 PSMs for 521 spectra
- 2:40 PM (2):Sequest HT Performing target search containing 521 spectra took 1.5 s.
- 2:40 PM (2):Sequest HT Start Sequest HT target search for 521 spectra.
- 2:40 PM (2):Sequest HT There is already an adequate target FASTA index.
- 2:40 PM (3):Spectrum Selector -- Total execution of Spectrum Selector (3) took 0.9 s --
- 2:40 PM (3):Spectrum Selector Sent 521 spectra from 1 files.
- 2:40 PM (3):Spectrum Selector Sent 521 spectra from file 1.
- 2:40 PM (2):Sequest HT Sequence Database: bovine.fasta
- 2:40 PM (3):Spectrum Selector Reading from file 1 of 1:C:\DevelopmentTestData\work\BSA3.raw (1330 spectra total)
- 2:40 PM ProcessingJob Processing D:\Work\BSA\BSA3.msf

A second job list is shown below, with two completed jobs for BSA2:

Execution State	Progress	Type	Name	Submitted at	Study
Completed	100 %	Consensus	BSA2	1/16/2015 2:39 PM	BSA
Completed	100 %	Processing	BSA2	1/16/2015 2:39 PM	BSA

The log for the BSA2 job includes the following messages:

- 2:39 PM ProcessingJob ----- Total Job execution took: 1.3 s. -----
- 2:39 PM ProcessingJob Finished D:\Work\BSA\BSA2.msf
- 2:39 PM (5):Spectrum Exporter -- Total execution of Spectrum Exporter (5) took 0.2 s --
- 2:39 PM (5):Spectrum Exporter Exported 521 spectra to "D:\Work\BSA\trypMyo.mgf"
- 2:39 PM (3):Spectrum Selector -- Total execution of Spectrum Selector (3) took 0.9 s --
- 2:39 PM (3):Spectrum Selector Sent 521 spectra from 1 files.
- 2:39 PM (3):Spectrum Selector Sent 521 spectra from file 1.
- 2:39 PM (5):Spectrum Exporter Start exporting of spectra to "D:\Work\BSA\trypMyo.mgf"
- 2:39 PM (3):Spectrum Selector Reading from file 1 of 1:C:\DevelopmentTestData\work\BSA2.raw (1330 spectra total)
- 2:39 PM ProcessingJob Processing D:\Work\BSA\BSA2.msf

Processing MudPIT Samples by Using a Processing Method


You can process MudPIT samples by using a processing method.

❖ To process MudPIT samples

1. Start the Proteome Discoverer Daemon utility and export a parameter file for MudPIT processing. For information about exporting a parameter file, see [“Running the Proteome Discoverer Daemon Utility from a Parameter File”](#) on page 153.

Figure 117 shows how to configure the Export Parameter File page in the utility to export a parameter file. In the following example, the parameter file is saved in C:\Xcalibur\methods.

Figure 117. Selecting MudPIT processing in the Start Jobs page



This example features two MudPIT samples. Each one consists of two raw data files (for a total of four raw data files).

2. Define a processing method (see [“Creating a Processing Method That Calls the Proteome Discoverer Daemon Utility”](#) on page 154) using the parameter file exported in step 1, and select the method as the processing method in the Proc Meth column, as shown in Figure 118.

Figure 118. Sequence used for MudPIT processing


	Sample Type	File Name	Sample ID	Path	Inst Meth	Proc Meth	Position	Inj Vol	Level
1	Unknown	BSA1	1	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\mudPIT	1	10.00	
2	Unknown	BSA2	2	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\mudPIT	2	10.00	
3	Unknown	BSA3	3	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\mudPIT	3	10.00	
4	Unknown	BSA4	3	C:\Projects\Daemon	C:\Xcalibur\methods\Daemon	C:\Xcalibur\methods\mudPIT	3	10.00	
*								0.00	

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Xcalibur Data System

3. Start processing the MudPIT samples in the Run Sequence dialog box, as shown in Figure 119.

Figure 119. Starting the processing of the MudPIT samples





The Proteome Discoverer Daemon utility processes the two samples as MudPIT, as shown in Figure 120.

Figure 120. Two MudPIT samples processed in the Proteome Discoverer Daemon utility

Execution State	Progress	Type	Name	Submitted at	Study	Data Source	Description	U
Completed	100 %	Consensus	mudpitResult	1/16/2015 4:14 PM	Daemon Temporary Study	C:\ProgramData\Thermo\Proteome Discoverer 2.0\PublicFiles\DiscovererDaemon...	Result filtered fo...	Then
Completed	100 %	Processing	mudpitResult	1/16/2015 4:14 PM	Daemon Temporary Study	Files:\DiscovererDaemon\SpectrumFiles\CWF_Basic_20150116_1\BSA3.raw\Dis...	Basic processing...	Then
Completed	100 %	Consensus	mudpitResult	1/16/2015 4:13 PM	Daemon Temporary Study	C:\ProgramData\Thermo\Proteome Discoverer 2.0\PublicFiles\DiscovererDaemon...	Result filtered fo...	Then
Completed	100 %	Processing	mudpitResult	1/16/2015 4:13 PM	Daemon Temporary Study	Files:\DiscovererDaemon\SpectrumFiles\CWF_Basic_20150116\BSA1.raw\Disco...	Basic processing...	Then

The Proteome Discoverer application saves the data in the two MudPIT samples in two directories, each one containing the raw data files of one MudPIT sample (in this example, two raw data files), as shown in Figure 121.

Figure 121. Two directories where the raw data files of each MudPIT group are saved



Running the Proteome Discoverer Daemon Utility from the Command Line

You can run the Proteome Discoverer Daemon utility from the command line or from the interface window.

❖ **To run the Proteome Discoverer Daemon utility on the command line**

1. Open a command shell and use the `cd` command to move to **Program Files > Thermo > Proteome Discoverer Daemon 2.1 > System > Release**.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Command Line

When you run the Proteome Discoverer Daemon utility from the command line, you must start it in the Proteome Discoverer Daemon 2.1 folder, not in the Proteome Discoverer folder.

2. Type **DiscovererDaemon** and any of the following options on the command line:

DiscovererDaemon

*[-e folderName FileCount workflows quanMethodFile
ParameterAssignment]*

[-c folderName]

[-a folderName spectrumFile]

[-l serverName userName]

[-r outputFilename]

[-p parameterFile rawFile]

[-f folderName]

For more information, see the following topics:

- [Syntax](#)
- [Examples](#)

Syntax

The Proteome Discoverer Daemon utility command-line syntax includes the following parameters:

- *[-e folderName filecount workflows quanMethodFile parameterAssignment]*

Executes the workflow on the server using these specified parameters:

- *folderName*: Specifies the location where the raw data files are stored. You can give it any name, for example, RawFiles or Fractions.
- *FileCount*: Specifies the number of spectrum files that must be included before the workflow is executed. This parameter is intended to be used with MudPIT experiments and acquisition on several machines. If the workflow should be executed regardless of the number of files contained in the file collection, use ANY instead of a number.
- *workflows*: Specifies a list of the template workflow files to run. You must have created these files in the Proteome Discoverer application by choosing Workflow Editor > *processing_workflow|consensus_workflow* > Save > *filename*.

You must list the workflows in the order that they are to be executed. Separate the workflow names by semicolons. For example, to run a processing workflow, C:\test.pdProcessingWF, and a consensus workflow, test.pdConsensusWF, use this syntax:

```
test.pdProcessingWF;test.pdConsensusWF
```

- *quanMethodFile*: Specifies the name of the file that contains the quantification method.
- *ParameterAssignment*: Specifies the name and value of a parameter in the format of *parameter=value*. Some examples follow.

This example sets the FASTA database for any node to equine.fasta:

```
FastaDatabase=equine.fasta
```

The next example sets the FASTA database for all Mascot nodes to equine.fasta:

```
Mascot.FastaDatabase=equine.fasta
```

The last example sets the FASTA database for Mascot nodes having 4 as the processing node number to equine.fasta. It is equivalent to [4].FastaDatabase=equine.fasta because the processing node numbers are unique.

```
Mascot[4].FastaDatabase=equine.fasta
```

- [-c *folderName*]
 - Remote server: Creates a user-named folder in the PublicFiles folder on the server where you store output files. The PublicFiles folder is the default folder in the Current File Directory box in the view displayed in the Proteome Discoverer application when you select Administration > Configuration > Server Settings > Discoverer Daemon. The -c option automatically appends the name to the date and, if the directory already exists, an incremental index number.

You can create a folder only in the directory configured in the view opened by the Administration > Configuration > Server Settings > Discoverer Daemon command on the remote server. If you attempt to create a folder in a location other than the PublicFiles folder in the Current File Directory box, the Proteome Discoverer Daemon utility issues a message informing you that it will apply the change the next time that you start the application.

This option performs the same function as the -f *foldername* option, except that you can use the name of the folder more than once. When you use the name more than once, the Proteome Discoverer Daemon utility appends the date and an incremental index number to the name.

 - Local server: Does nothing.
- [-a *folderName SpectrumFile*]
 - Remote server: Uploads the spectrum file to the location that you specified on the configured server. *SpectrumFile* is the name of the spectrum file.
 - Local server: Does nothing.
- [-l *serverName userName*]: Connects the Proteome Discover Daemon application to the specified local or remote host machine.
 - *serverName*: Specifies the name of the local or remote host.
 - *userName*: Specifies the name to log in.

3 Using the Proteome Discoverer Daemon Utility

Running the Proteome Discoverer Daemon Utility from the Command Line

- [-r *outputFilename*]: Specifies the name of the output file. You must use this option with the -e option, as in this example:

```
DiscovererDaemon.exe -e sfcid any mascot3.xml -r silac1noMT_AS4DE.msfc
```

- [-p *parameterFile rawFile*]: Processes the specified raw data file with all the parameters given in the parameter file, including the connection to the server.
 - *parameterFile*: Specifies the name of the parameter file. In addition to the parameters, the parameter file also contains the workflows to run.
 - *rawFile*: Specifies the name of the raw data file.

In the following example of the -p syntax, the Proteome Discoverer Daemon utility processes the 9mix_LysC_monolith.raw file with the parameters given in the parameter file called C:\Xcalibur\methods\batchprocessing.param.

```
DiscovererDaemon.exe -p C:\Xcalibur\methods\batchprocessing.param  
9mix_LysC_monolith.raw
```

- [-f *folderName*]: On a remote server, this option creates a user-named folder in the PublicFiles folder of the server where the local version of the raw data file and the result files are stored. If the directory already exists, the Proteome Discoverer Daemon utility issues an error message, and the process returns with exit code -1 (standard exit code 0).

If you attempt to create a file other than in the PublicFiles folder in the Current File Directory box, the Proteome Discoverer Daemon utility issues a message informing you that the Proteome Discoverer application will apply the change the next time that you start it.

This option performs the same function as the -c *foldername* option, except that you cannot use the name of the folder more than once.

On a local server, this option does nothing.

Examples

The following are examples of the Proteome Discoverer Daemon utility command-line syntax.

Example 1

This example constructs the spectrum file collection called Rawfiles, adds the TrypMyo.raw file to the collection, and executes the SequestEquine workflow using the raw data file in the Rawfiles directory:

```
DiscovererDaemon.exe -c Rawfiles -a Rawfiles  
C:\Rawfiles\TrypMyo.raw -e Rawfiles ANY C:\Workflows\SequestEquine.xml
```


Example 2

In the following example, the Proteome Discoverer Daemon utility evaluates several fractions in a single workflow:

```
DiscovererDaemon.exe -f Fractions
DiscovererDaemon.exe -a Fractions C:\rawfiles\fraction1.raw
DiscovererDaemon.exe -a Fractions C:\rawfiles\fraction2.raw
DiscovererDaemon.exe -a Fractions C:\rawfiles\fractionN.raw
DiscovererDaemon.exe -e Fractions ANY
pathToFile.pdProcessingWF;pathToOtherFile.pdConsensusWF
```

Example 3

The next example demonstrates that you can start several workflows with one invocation of the Proteome Discoverer Daemon utility.

```
DiscovererDaemon.exe
-f RawFile
-a RawFile C:\Rawfiles\TrypMyo.raw
-e RawFile ANY pathToFile.pdProcessingWF;pathToOtherFile.pdConsensusWF
-a RawFile C:\Rawfiles\BSADigest.raw
-e RawFile ANY pathToFile.pdProcessingWF;pathToOtherFile.pdConsensusWF
```

Example 4

The following example runs the Proteome Discoverer Daemon utility on a remote host called protlab2, uploads the iTRA_BSA_3ITMS2_3HCD.raw spectrum file to the server, and executes the workflows in C:\Workflows:

```
DiscovererDaemon.exe -l protlab2 leo_davinci -c sfcid -a sfcid
iTRA_BSA_3ITMS2_3HCD.raw -e sfcid any C:\workflows\MascotEcoli\
pathToFile.pdProcessingWF;pathToOtherFile.pdConsensusWF
```

Example 5

In this example, the following sequence of commands submits multiple raw data files for processing on a remote server. The 020110303 notation indicates the date that the Daemon command was issued.

```
DiscovererDaemon.exe -c AllTrypMyo
DiscovererDaemon.exe -a AllTrypMyo_020110303
C:\DaemonTest\mudpit4\Tryp_Myo.raw
DiscovererDaemon.exe -a AllTrypMyo_020110303
C:\DaemonTest\mudpit4\Tryp_Myo_1.raw
DiscovererDaemon.exe -a AllTrypMyo_020110303
C:\DaemonTest\mudpit4\Tryp_Myo_2.raw
```

```
DiscovererDaemon.exe -e AllTrypMyo_020110303 3  
C:\DaemonTest\mudpit4\pathToFile.pdProcessingWF;  
pathToOtherFile.pdConsensusWF
```

Example 6

In this example, the next sequence of commands submits multiple raw data files for processing on a local server:

```
DiscovererDaemon.exe -a AllTrypMyo C:\DaemonTest\mudpit4\Tryp_Myo.raw  
DiscovererDaemon.exe -a AllTrypMyo  
C:\DaemonTest\mudpit4\Tryp_Myo_1.raw  
DiscovererDaemon.exe -a AllTrypMyo  
C:\DaemonTest\mudpit4\Tryp_Myo_2.raw  
DiscovererDaemon.exe -e AllTrypMyo 3  
C:\DaemonTest\mudpit4\pathToFile.pdProcessingWF;  
pathToOtherFile.pdConsensusWF
```

Example 7

If you want to perform quantification, the commands resemble this example:

```
DiscovererDaemon.exe -c iTRAQInput  
DiscovererDaemon.exe -a iTRAQInput C:\testData\iTraq\  
myo_8plex_hcd_pqd.raw  
DiscovererDaemon.exe -e iTRAQInput 1  
C:\Users\Public\Documents\Thermo\Proteome_Discoverer_2.1\  
Common_Templates\iTRAQ_8plex_Sequest_HT_equine.pdProcessingWF;  
C:\Users\Public\Documents\Thermo\Proteome_Discoverer_2.1\  
Common_Templates\WF_standard_automatic_validation_with_Quan  
.pdConsensusWF;  
C:\Users\Public\Documents\Thermo\Proteome_Discoverer_2.1\Com  
mon_Templates\iTRAQ_8plex.method
```

The Proteome Discoverer Daemon utility appends a time stamp to each file when it processes the file on a remote server.

Searching for Data

This chapter describes the features that you can use when searching for and analyzing data in the Proteome Discoverer application.

Contents

- [Using FASTA Databases](#)
- [Identifying Contaminants During Searches](#)
- [Searching Spectrum Libraries](#)
- [Defining Chemical Modifications](#)
- [Using the Qual Browser Application](#)
- [Customizing Cleavage Reagents](#)

Using FASTA Databases

You can use the FASTA database utilities to add, delete, and find protein references and sequences. You can also extract information from an existing FASTA file, place it into a new FASTA file, and compile it for availability in the Proteome Discoverer application.

For more information about FASTA databases, see [“FASTA Reference”](#) on [page 491](#).

- [Displaying FASTA Files](#)
- [Adding FASTA Files to the Proteome Discoverer Application](#)
- [Deleting FASTA Files](#)
- [Adding Protein Sequences and References to a FASTA Database File](#)
- [Finding Protein Sequences and References](#)
- [Compiling a FASTA Database](#)
- [Excluding Individual Protein References and Sequences from a FASTA Database](#)
- [FASTA Database Utilities Dialog Box Parameters](#)

- [Managing FASTA Indexes](#)
- [Adding or Modifying FASTA Parsing Rules](#)
- [Identifying Contaminants During Searches](#)

Displaying FASTA Files

You can list all the FASTA files that you have downloaded from other sources onto your hard drive and registered.

❖ To list the available FASTA files

Choose **Administration > Maintain FASTA Files** or click the **Maintain FASTA Files** icon, , either in the toolbar or on the Administration page.

The FASTA Files view shown in [Figure 122](#) appears. It lists all the FASTA files that you have downloaded from other sources and registered. It displays the processed FASTA file properties, such as the file name, file size, and the number of proteins stored. The Proteome Discoverer application analyzes each protein entry to determine if the FASTA file meets the application requirements for use in a spectra search. It processes the FASTA file and makes it available for use.

Figure 122. FASTA Files view

The screenshot shows the Proteome Discoverer interface with the FASTA Files view. The top toolbar contains several icons, each labeled with a function:

- Add icon
- Download icon
- Update icon
- Remove icon
- Cancel icon
- Export icon
- Refresh icon
- Check for Updates icon

The main window displays a table of FASTA files with the following columns: Name, Protein Database, Taxonomy ID, Version, File Size [kB], # Sequences, # Residues, Status, Last Modified, and Update Available. The table lists various species and their corresponding FASTA files, such as Homo_sapiens_Ref, equine.fasta, and ecoli.fasta.

Name	Protein Database	Taxonomy ID	Version	File Size [kB]	# Sequences	# Residues	Status	Last Modified	Update Available
Homo_sapiens_Ref	Custom			22657	31237	19719930	Availa...	11/16/2012	<input type="checkbox"/>
equine.fasta	Custom			149	346	75453	Availa...	08/02/2004	<input type="checkbox"/>
ecoli.fasta	Custom			6916	4939	1563388	Availa...	06/23/2006	<input type="checkbox"/>
cRAP.fasta	Custom			40	115	38188	Availa...	10/02/2013	<input type="checkbox"/>
PD_Contaminants...	Custom			155	298	130376	Availa...	11/11/2014	<input type="checkbox"/>
Mus_musculus_Uni...	Custom			38540	73353	3182568	Availa...	11/04/2014	<input type="checkbox"/>
Chlamydomonas_r...	Custom			8488	15066	6958843	Availa...	11/04/2014	<input type="checkbox"/>
Arabidopsis_thalia...	Custom			25510	51667	20424187	Availa...	11/04/2014	<input type="checkbox"/>
yeast.fasta	Custom			14516	15178	6101720	Availa...	06/23/2006	<input type="checkbox"/>
bovine.fasta	Custom			670	1057	242627	Availa...	06/23/2006	<input type="checkbox"/>
Lasiorhinus krefftii...	TrEMBL	30667	2014-10-01	0	1	30	Availa...	11/28/2014	<input type="checkbox"/>
Squatina squatina...	GI	263718	2014-10-19	3	6	1312	Availa...	11/28/2014	<input type="checkbox"/>
Ardops nichollsi (Tr...	TrEMBL	148034	2014-10-01	4	10	3429	Availa...	11/28/2014	<input type="checkbox"/>
Delphinus delphis (...)	SwissProt	9728	2014-10-01	2	5	1298	Availa...	11/28/2014	<input type="checkbox"/>
Xiphias gladius (Tr...	TrEMBL	8245	2014-10-01	120	211	61708	Availa...	11/28/2014	<input type="checkbox"/>
Hippopotamus am...	SwissProt	9833	2014-10-01	8	21	4841	Availa...	11/28/2014	<input type="checkbox"/>
Corvus hawaiiensis...	TrEMBL	134902	2014-10-01	0	2	448	Availa...	11/28/2014	<input type="checkbox"/>
Papio (SwissProt T...	SwissProt	9554	2014-10-01	112	195	66278	Availa...	11/28/2014	<input type="checkbox"/>
Cynomys parviden...	SwissProt	99827	2014-10-01	0	1	379	Availa...	11/28/2014	<input type="checkbox"/>
Sciurus carolinensi...	SwissProt	30640	2014-10-01	1	5	1054	Availa...	11/28/2014	<input type="checkbox"/>
Quercus robur (Swi...	SwissProt	38942	2014-10-01	0	1	504	Availa...	11/28/2014	<input type="checkbox"/>
Fagus sylvatica (S...	SwissProt	28930	2014-10-01	1	3	525	Availa...	11/28/2014	<input type="checkbox"/>
Meleagris gallopav...	SwissProt	9103	2014-10-01	24	59	13257	Availa...	11/28/2014	<input type="checkbox"/>
Mammuthus primig...	SwissProt	37349	2014-10-01	12	14	4097	Availa...	11/28/2014	<input type="checkbox"/>
Elephas maximus (...)	TrEMBL	9783	2014-10-01	122	183	43600	Availa...	11/28/2014	<input type="checkbox"/>
Canis lupus familia...	SwissProt	9615	2014-10-01	923	857	381701	Availa...	11/28/2014	<input type="checkbox"/>

Adding FASTA Files to the Proteome Discoverer Application

You must add a FASTA file to the Proteome Discoverer application before you can conduct a search with the Sequest HT node. You must add a FASTA file to the server that is running Mascot before you can conduct a search with the Mascot search engine. You can download FASTA files from two sources:

- You can download a controlled protein database directly from ProteinCenter and add it to a FASTA file. These controlled databases offer you access to the latest protein information. The ProteinCenter database service provides extensive information about proteins, peptides, and genes. In addition, it gives you the ability to select proteins of a specified taxonomy to download and use as a FASTA file.
- You can add a FASTA file that you have downloaded from other sources onto your hard drive.

❖ To download a FASTA file from ProteinCenter

1. Configure the Proteome Discoverer application for protein annotation. For instructions, see “[Configuring the Proteome Discoverer Application for Protein Annotation](#)” on [page 293](#).

2. Choose **Administration > Maintain FASTA Files** or click the **Maintain FASTA Files** icon, .


The Administration page appears with the FASTA files view, shown in [Figure 122](#) on [page 173](#).

3. Click the **Check for Updates** icon,  **Check for Updates**.

The Proteome Discoverer application updates the available databases in the job queue.


4. Click **FASTA Files** under Content Management on the Configuration view to return to the FASTA Files view.

The Download icon on the FASTA Files view becomes available.

5. Click the **Download** icon,  **Download**.


The Download from ProteinCenter dialog box appears, as shown in [Figure 123](#).

Figure 123. Download from ProteinCenter dialog box




6. In the Taxonomy ID box, do the following:
 - a. Type the taxonomy identification number of the appropriate organism-specific sequence database in ProteinCenter.

The taxonomy identification number is a unique number identifying a biological species, a special subspecies, or a bacterial strain. You can find the organism's taxonomy identification number on the UniProt™ website (<http://www.uniprot.org/taxonomy/>) or at the NCBI.

- i. In the box to the right of the **Taxonomy** menu on the UniProt website, type the name of the organism that you are interested in, for example, **Baker's yeast**.
- ii. Click the **Search** icon, .
- iii. (Optional) Under the Taxon heading, click the name of the subspecies that you are interested in, for example, **Saccharomyces cerevisiae**.

The taxonomy identification number appears on the Taxon Identifier line, as shown in [Figure 124](#).

Figure 124. Taxonomy identification on the UniProt website




- b. (Optional) To include data for a subspecies or subcategory of the selected species in the downloaded database, select the **Include All Subcategories** check box in the Download from ProteinCenter dialog box. For more information on this option, see the Help.
- c. From the list in the Database box of the Download from ProteinCenter dialog box, select the name of the original source database to download the proteins from.

The default is SwissProt.

Figure 125 shows the completed Download from ProteinCenter dialog box.

Figure 125. Completed Download from ProteinCenter dialog box



- d. Click **Import**.

The application now displays the download as a job running in the job queue.

- When the job queue displays “Completed” in the Execution State column, click **FASTA Files** under Content Management on the Configuration view to return to the FASTA Files view.

The downloaded database appears in the FASTA Files view. It might take several minutes to appear.

- If you do not see the downloaded database after a few minutes, click the **Refresh** icon,




Figure 126 shows the *Saccharomyces cerevisiae* (4932) species database downloaded from the SwissProt database.

Figure 126. Database downloaded from ProteinCenter

Name	Protein Database	Taxonomy ID	Version	File Size [kB]	# Sequences	# Residues	Status	Last Modified	Update Available
Mus musculus (SwissProt TaxID=10090)	SwissProt	10090	2015-04-29	16661	24742	14437410	Available	07/08/2015	<input type="checkbox"/>
Saccharomyces cerevisiae (SwissProt TaxID=4932)	SwissProt	4932	2015-04-29	23	41	20162	Available	07/24/2015	<input checked="" type="checkbox"/>

Database downloaded from ProteinCenter


If you select the SwissProt database, all proteins of the specified organism that are present in the SwissProt database are downloaded. The FASTA file is automatically given the scientific name of the organism and the database source, with the taxonomy identification in parentheses, for example, *Saccharomyces cerevisiae* (SwissProt TaxID =4932). It also gives the date the database was updated on the server, for example, (SwissProt - 2015-04-29).

If no databases appear in the list, follow the procedure in “[To update a FASTA file from ProteinCenter.](#)”


❖ To download a FASTA file from sources other than ProteinCenter

- Choose **Administration > Maintain FASTA Files** or click the **Maintain FASTA Files** icon, .

The Administration page appears with the FASTA files view, shown in [Figure 122](#) on [page 173](#).

- Click the **Add** icon, .

- In the Open dialog box that appears, browse to and select the FASTA file that you want to register, and then click **Open**.

The FASTA file that you selected appears as a job in the job queue. To cancel the addition of this file, click the **Abort** icon, .

When you see “Completed” in the Execution State column, the database has finished downloading.

4. To add another FASTA file, wait until the Execution State column indicates that the addition of the FASTA file is completed, click **FASTA Files** in the left pane of the Administration page under Content Management, and then click **Add** to add the next file.



The amount of time that it takes to import a FASTA file depends on the file size. When the application finishes importing a FASTA file, it displays “Available” in the Status column. The FASTA file is now available to use for a protein or peptide search with the Proteome Discoverer application.

❖ To update a FASTA file from ProteinCenter


Note You can only use this procedure to update FASTA databases that you originally downloaded directly from ProteinCenter.

1. Choose **Administration > Maintain FASTA Files** or click the **Maintain FASTA Files** icon,  .

The Administration page appears with the FASTA files view, shown in [Figure 122](#) on [page 173](#).

2. Select the database that you want to check for the availability of updated information. Click  at the beginning of a database row to select the row.
3. Click the **Check for Updates** icon,  **Check for Updates** .

The application now displays the search for updates as a job running in the job queue.


4. When the job queue displays Completed in the Execution State column, close the Proteome Discoverer application and reopen it.
5. Choose **Administration > Maintain FASTA Files** or click the **Maintain FASTA Files** icon,  .

If an update is available for the selected database, a check mark appears in the Update Available column of the FASTA Files view, as shown in [Figure 127](#).


Figure 127. Check marks indicating update availability

The screenshot shows the Proteome Discoverer Administration interface. The main window displays a table of FASTA files with columns for Name, Protein Database, Taxonomy ID, Version, File Size [kB], # Sequences, # Residues, Status, Last Modified, and Update Available. Checkmarks in the Update Available column indicate that updates are available for several files, including those from SwissProt and TrEMBL databases.

Name	Protein Database	Taxonomy ID	Version	File Size [kB]	# Sequences	# Residues	Status	Last Modified	Update Available
Homo_sapiens...	Custom			22657	31237	19719930	Availa...	11/16/2012	<input type="checkbox"/>
equine.fasta	Custom			149	346	75453	Availa...	08/02/2004	<input type="checkbox"/>
ecoli.fasta	Custom			6916	4939	1563388	Availa...	06/23/2006	<input type="checkbox"/>
cRAP.fasta	Custom			40	115	38188	Availa...	10/02/2013	<input type="checkbox"/>
PD_Contamina...	Custom			155	298	130376	Availa...	11/11/2014	<input type="checkbox"/>
Mus_musculus...	Custom			38540	73353	31825568	Availa...	11/04/2014	<input type="checkbox"/>
Chlamydomon...	Custom			8488	15066	6958843	Availa...	11/04/2014	<input type="checkbox"/>
Arabidopsis_th...	Custom			25510	51667	20424187	Availa...	11/04/2014	<input type="checkbox"/>
yeast.fasta	Custom			14516	15178	6101720	Availa...	06/23/2006	<input type="checkbox"/>
bovine.fasta	Custom			670	1057	242627	Availa...	06/23/2006	<input type="checkbox"/>
Lasiorhinus kre...	TrEMBL	30667	2014-10-01	0	1	30	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Squatina squati...	GI	263718	2014-10-19	3	6	1312	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Ardops nichols...	TrEMBL	148034	2014-10-01	4	10	3429	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Delphinus dep...	SwissProt	9728	2014-10-01	2	5	1298	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Xiphias gladius...	TrEMBL	8245	2014-10-01	120	211	61708	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Hippopotamus...	SwissProt	9833	2014-10-01	8	21	4841	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Corvus hawaii...	TrEMBL	134902	2014-10-01	0	2	448	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Papio (SwissPr...	SwissProt	9554	2014-10-01	112	195	66278	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Cynomys parvi...	SwissProt	99827	2014-10-01	0	1	379	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Sciurus carolin...	SwissProt	30640	2014-10-01	1	5	1054	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Quercus robur...	SwissProt	38942	2014-10-01	0	1	504	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Fagus sylvatica...	SwissProt	28930	2014-10-01	1	3	525	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Meleagris gallo...	SwissProt	9103	2014-10-01	24	59	13257	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Mammuthus pr...	SwissProt	37349	2014-10-01	12	14	4097	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Elephas maxim...	TrEMBL	9783	2014-10-01	122	183	43600	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Canis lupus fa...	SwissProt	9615	2014-10-01	923	857	381701	Availa...	11/28/2014	<input checked="" type="checkbox"/>
Escherichia coli...	SwissProt	83333	2014-10-01	2323	4387	1381889	Availa...	12/08/2014	<input type="checkbox"/>

6. Select the database with the check mark in the Update Available column. Click  at the beginning of a row to select the row.

The Update icon becomes available.

7. Click the **Update** icon,  **Update**.

The application now displays the update as a job running in the job queue.

8. (Optional) When the job queue displays Completed in the Execution State column, click **FASTA Files** under Content Management to return to the FASTA Files view.

The data on the ProteinCenter server is updated every two weeks. When the ProteinCenter server contains new data, the downloaded FASTA files might change.

An update on the ProteinCenter server does not mean that new proteins for the FASTA files are automatically available. When you can access the new proteins depends on the timing of submissions to a specific database and that database's release cycle. For example, the UniProt databases are updated every month, so a new ProteinCenter version after two weeks does not necessarily contain new proteins.

Exporting FASTA Files

You can export FASTA files from the view activated by Administration > Maintain FASTA Files, including FASTA files uploaded to the Proteome Discoverer application and files downloaded from ProteinCenter.

The exported FASTA file might be different from the original FASTA file. The application removes invalid and duplicate protein sequences during importation and cannot export them, although it retains protein title lines. The exported FASTA file lists the duplicate proteins with their title lines in consecutive order.

❖ To export a FASTA file

1. Choose **Administration > Maintain FASTA Files**.
2. Select the FASTA file that you want to export.
3. Click the **Export** icon,  **Export** .
4. In the Save As dialog box, browse to the directory where you want to save the FASTA file, and click **Save**.

The job queue opens and displays a status of “Completed” in the Execution State column when the application finishes exporting the file.



Deleting FASTA Files


You can delete a FASTA file from the application.

❖ To delete a FASTA file

1. Choose **Administration > Maintain FASTA Files** or click the **Maintain FASTA Files** icon,  .

The Administration page appears with the FASTA files view, shown in [Figure 122](#) on [page 173](#).

2. Click  at the beginning of a row to select the row to delete.
3. Click the **Remove** icon,  **Remove** .
4. In the Remove FASTA Databases dialog box, click **OK**.

The FASTA file that you selected appears as a job in the job queue. After you start the deletion of the file, you cannot cancel the deletion. You can remove the completed job from the job queue by clicking the Remove icon,  **Remove** , and then clicking **OK** in the Delete Jobs dialog box.

Adding Protein Sequences and References to a FASTA Database File


You can add protein sequences and protein references to a registered FASTA database file. The protein sequence refers to the sequence of amino acids that constitute the protein, and the protein reference refers to the name or reference of the protein.

❖ To add a protein sequence and reference to a FASTA database file

1. Choose **Tools > FASTA Database Utilities**.
2. In the FASTA Database Utilities dialog box, click the **Add Protein References** tab.
The Add Protein References page of the dialog box appears.
3. Click the **Browse** button (...) next to the FASTA File box.
4. In the Save/Add to FASTA File dialog box, select the FASTA database that you want to add the protein sequence and reference to, and click **Save**.
5. In the Enter Description box of the FASTA Database Utilities dialog box, type a description of the protein sequence that you are adding.
6. In the Enter Protein Sequence box, type the protein sequence that you want to add to the FASTA database.

The Add Protein References page should resemble the illustration in [Figure 128](#).

Figure 128. Add Protein References page of the FASTA Database Utilities dialog box



7. Click **Add Entry** to add the protein sequence.

For information on the Add Protein References page parameters, refer to the Help.

Finding Protein Sequences and References

You can find a protein sequence or reference in an existing FASTA database file.


- [To find a protein sequence or reference](#)
- [To filter a protein reference search](#)
- [To refine a filtered protein reference search](#)
- [To delete conditions in filtered protein reference searches](#)

❖ To find a protein sequence or reference

1. Choose **Tools > FASTA Database Utilities**.
2. In the FASTA Database Utilities dialog box, click the **Find Protein References** tab.

The Find Protein References page appears, as shown in [Figure 129](#).

Figure 129. Find Protein References page of the FASTA Database Utilities dialog box




3. Click the **Browse** button (...) next to the FASTA Database box to locate the FASTA file of interest.
4. In the Please Select a FASTA Database dialog box, select the FASTA file, and click **Open**.
5. In the Search For box of the Find Protein References page, type an amino acid sequence or a protein reference search string.
6. In the Search In area, select the **References** option (the default) to search for a reference or the **Sequences** option to search for a sequence.
 - **References:** Searches for the search string in the protein references.
 - **Sequences:** Searches for the specified amino acid sequence within the protein sequences.

You can further refine the results by using filters either before or after you run the search. For instructions on filtering, see [“To filter a protein reference search”](#) on [page 184](#).

7. In the Maximum Number of Matches Reported box, select the maximum number of references or sequences to report.
8. Click **Start Search**.

Results appear if the search parameters match the data, as shown in [Figure 130](#). Click a protein row to see the amino acid sequences that constitute that protein.
9. To suspend the search, click **Stop Search**.

Figure 130. Find Protein References page of the FASTA Database Utilities dialog box



- (Optional) To save a protein result row in another FASTA database, select the protein row, click **Save/Add Selected to Database**, select the database in the Save/Add to FASTA File dialog box, and click **Save**.


For information on the Find Protein References page parameters, refer to the Help.

❖ To filter a protein reference search

- On the Find Protein References page of the FASTA Database Utilities dialog box, click the line below “Reference” in the middle of the page to access a list of operators that you can use to filter the references. (The default operator is “Starts with.”) For a list of all operators, refer to the Help.
- In the line below the operator that you selected, type the search string or condition that you want the operator to apply to.

The example in [Figure 131](#) filters out those protein references that contain “fragment.”

Figure 131. Filtering out protein references containing “fragment”




❖ **To refine a filtered protein reference search**

1. Select the **Custom** option from the list in the line below the operator.

To make the Custom option available, click the down arrow in the line below the operator, as shown in [Figure 132](#).


Figure 132. Selecting the Custom option



Click this down arrow.

The Custom option opens the Custom Filter dialog box, shown in [Figure 133](#), so you can add multiple conditions.

Figure 133. Custom Filter dialog box



2. Click **Add**.

A new line appears in the Operator (left) and Operand (right) lists.

3. Select an operator from the Operator list.

4. Type an operand on the line in the Operand column.

5. In the Filter Based On list, do one of the following:


Select the **All** option to indicate whether the search algorithm should search for protein references that meet both conditions.

–or–

Select the **Any** option to indicate whether the search algorithm should search for protein references that meet only one of the conditions.


Figure 134 gives an example of a search for protein references that meet both of the conditions.

Figure 134. Custom filter with a second condition



6. Click **OK**.

❖ **To delete conditions in filtered protein reference searches**

- To delete a condition in the Custom Filter dialog box, select the check box to the left of the appropriate condition in the Operator column, and click **Delete**.
- To delete the condition in the Reference area on the Find Proteins References page, click the **Clear Reference Filter Criteria** icon, , in the line below the operator.

- To delete all conditions in both the Custom Filter dialog box and the Reference area on the Find Proteins References page, click the **Clear All Filter Criteria** icon, , to the left of the filters.

Compiling a FASTA Database

You can extract information from an existing FASTA file and place it into a new FASTA file, replace an existing FASTA file, or append it to an existing FASTA file. Then you must compile the new or changed FASTA file to make it available in the Proteome Discoverer application.


❖ To compile a FASTA database

1. Choose **Tools > FASTA Database Utilities**.
2. In the FASTA Database Utilities dialog box, click the **Compile FASTA Database** tab.
The Compile FASTA Database page appears.
3. In the Original box, browse to the FASTA file that you are taking the information from, or type its path and name.
4. In the Please Select a FASTA Database dialog box, click **Open**.
5. In the Target box, browse to the FASTA file that you are placing the extracted information into, or type its path and name.
6. In the Save/Add to FASTA File dialog box, select the file, verify that the file extension is .fasta, and click **Save**.
7. In the Target Database Options area, select one of the following options to indicate what you want to do with the extracted information:
 - **Create/Replace**: Creates a new FASTA file for storing the information or overwriting an existing FASTA file. This option is the default.
 - **Append**: Adds the extracted information to an existing FASTA file.
8. In the Search In area, specify whether the Proteome Discoverer application should search for the search string in the protein references or sequences.
 - **References**: Searches for the search string in the protein references.
 - **Sequences**: Searches for the specified amino acid sequence within the protein sequences.
9. To disregard the case of the information to be extracted, select the **Ignore case of reference strings** check box.

10. Specify the information to be extracted:


- a. Click  above the Step 1: String(s) to Include box.

A line enabling you to specify the first set of conditions appears in the box.

- b. Click the first line in the Select Operator column, and select the operator to apply to the information to be extracted. You can select from the following:
- Starts With: Extracts information that begins with this string.
 - Does Not Start With: Extracts information that does not begin with this string.
 - Ends With: Extracts information that ends with this string.
 - Does Not End With: Extracts information that does not end with this string.
 - Contains: Extracts information that includes this string.
 - Does Not Contain: Extracts information that does not includes this string.
- c. Click the first line in the Condition column, and type the condition that the information must meet in order to be extracted.
- d. Repeat [step a](#) through [step c](#) to add more sets of conditions for the information to be extracted.
- e. To delete a set of conditions, in the Active column select the line that you want to delete and click .

The Compile FASTA Database page should now resemble the example in [Figure 135](#).

Figure 135. Compile FASTA Database page of the FASTA Database Utilities dialog box




11. Click **Compile Database**.

Click **Stop** to halt the compilation.


12. After the compilation, click **Start Search** on the Find Protein References page to view the results of the extraction, as shown in the example in [Figure 136](#).

You do not have to enter information into the Search For box.

Figure 136. Results of search




13. (Optional) To specify any information that you want to exclude from the extracted results, follow these steps:

- a. Click  above the Step 2: String(s) to Exclude From the Results of Step 1 box on the Compile FASTA Database page.

A line enabling you to specify the first set of conditions now appears in the box.

- b. Click the first line in the Select Operator column, and select the operator to apply to the information from the list. You can choose from the following:
 - Starts With: Excludes information that begins with this string.
 - Does Not Start With: Excludes information that does not begin with this string.
 - Ends With: Excludes information that ends with this string.
 - Does Not End With: Excludes information that does not end with this string.

- Contains: Excludes information that includes this string.
 - Does Not Contain: Excludes information that does not include this string.
- c. Click the first line in the Condition column, and type the condition that the information must meet in order to be excluded.
 - d. Repeat [step a](#) through [step c](#) to add more sets of conditions for the information that you want to exclude.
 - e. To delete a set of conditions, in the Active column select the line that you want to delete and click .
14. Click **Compile Database**.
 15. Click **Start Search** on the Find Protein References page to view the results of the extraction, as shown in the example in [Figure 136](#).



You do not have to enter information into the Search For box.

For information on the Compile FASTA Database page parameters, refer to the Help.

Excluding Individual Protein References and Sequences from a FASTA Database

You can exclude individual entries from a FASTA file.

❖ To exclude individual protein references and sequences from a FASTA file

1. Choose **Tools > FASTA Database Utilities**.
2. In the FASTA Database Utilities dialog box, click the **Compile FASTA Database** tab.
3. In the Original box, browse to the FASTA database that contains the protein that you want to remove, or type its path and name. In the Please Select a FASTA Database dialog box, click **Open**.
4. In the Target box, browse to the output FASTA file or type its path and name. In the Save/Add to FASTA File dialog box, select the file, verify that the file extension is .fasta, and click **Save**.
5. Select the **Ignore Case of References Strings** check box.
6. Click  above the Step 1: String(s) to Include box.
A line enabling you to specify the first set of conditions now appears in the box.
7. Click the first line in the Select Operator column, and select **Contains**, if it is not already selected. Leave the first line in the Condition column blank.
8. Click  above the Step 2: String(s) to Exclude From the Results of Step 1 box.
A line enabling you to specify the first set of conditions now appears in the box.
9. Click the first line in the Select Operator column, and select **Contains**.

10. In the first line of the Condition column, type the protein reference or sequence that you want to remove.

11. Click **Compile Database**.

The compiling process creates the target FASTA file that excludes protein entries that match the condition.

FASTA Database Utilities Dialog Box Parameters

The following topics list the parameters on the pages of the FASTA Database Utilities dialog box.

- [Add Protein References Page](#)
- [Compile FASTA Database Page](#)
- [Find Protein References Page](#)

Add Protein References Page

Table 6 describes the parameters on the Add Protein References page of the FASTA Database Utilities dialog box.

Table 6. Add Protein References page parameters

Parameters	Description
FASTA File	Specifies the name of the FASTA file to modify.
Enter Description	Describes the protein sequence that you are appending to the selected FASTA file. The description should include the name of the protein and the source tissue.
Enter Protein Sequence	Specifies the amino acid sequence of the protein that you are appending to the selected FASTA file.
Add Entry	Appends and saves the protein sequence and reference to the currently selected FASTA file.

Compile FASTA Database Page

Table 7 describes the parameters on the Compile FASTA Database page of the FASTA Database Utilities dialog box.

Table 7. Compile FASTA Database page parameters (Sheet 1 of 3)

Parameters	Description
Original	Specifies the path and name of the source FASTA database that you are taking the information from. It must be an existing FASTA file.
Target	Specifies the path and name of the FASTA file that you are placing or adding the extracted information to.
Target Database Options	<p>Determines whether the extracted information is placed in a new FASTA file or is appended to an existing FASTA file.</p> <ul style="list-style-type: none"> • Create/Replace (Default): Creates a new FASTA file where you will store the extracted information. • Append: Adds the extracted information to an existing FASTA file.
Search In	<p>Determines whether to search for the search string in the protein references or sequences.</p> <ul style="list-style-type: none"> • References: Searches for the search string in the protein references. • Sequences: Searches for the specified amino acid sequence within the protein sequences.
Ignore Case of Reference Strings	Disregards the case of the text in the Step 1: String(s) to Include box and the Step 2: String(s) to Exclude From the Results of Step 1 box.
Step 1: String(s) to Include:	Specifies the information to be extracted from the original file.
Active	Specifies the number of the condition set and allows you to select it.

Table 7. Compile FASTA Database page parameters (Sheet 2 of 3)





Parameters	Description
Select Operator	<p>Specifies the action to apply to the condition selected in the Condition column:</p> <ul style="list-style-type: none"> • Starts With: Extracts information that begins with this string. • Does Not Start With: Extracts references that do not begin with this string. • Ends With: Extracts information that ends with this string. • Does Not End With: Extracts references that do not end with this string. • Contains: Extracts information that includes this string. • Does Not Contain: Extracts references that do not contain this string.
Condition	Specifies the condition that the information must meet in order to be extracted.
	Adds a line to the Step 1: String(s) to Include box.
	Deletes a line from the Step 1: String(s) to Include box.
Step 2: String(s) to Exclude From the Results of Step 1	Specifies the information to be excluded from the extracted results.
Active	Specifies the number of the condition set and allows you to select it.
Select Operator	<p>Specifies the action to apply to the condition selected in the Condition column:</p> <ul style="list-style-type: none"> • Starts With: Excludes information that begins with this string. • Does Not Start With: Excludes information that does not begin with this string. • Ends With: Excludes information that ends with this string. • Does Not End With: Excludes information that does not end with this string. • Contains: Excludes information that includes this string. • Does Not Contain: Excludes information that does not include this string.

Table 7. Compile FASTA Database page parameters (Sheet 3 of 3)

Parameters	Description
Condition	Specifies the condition that the information must meet in order to be excluded.
	Adds a line to the Step 2: String(s) to Exclude From the Results of Step 1 box.
	Deletes a line from the Step 2: String(s) to Exclude From the Results of Step 1 box.
Compile Database	Compiles and saves the changes to the FASTA file specified in the Target box.
Stop	Suspends the compilation.



Find Protein References Page

[Table 8](#) describes the parameters on the Find Protein References page of the FASTA Database Utilities dialog box.

Table 8. Find Protein References page parameters (Sheet 1 of 2)

Parameters	Description
FASTA Database	Specifies the name of the FASTA file that you want to search for sequences or references in.
Search For	Specifies the amino acid sequence to search for.
Search In	Determines whether to search for the search string in the protein references or sequences. <ul style="list-style-type: none"> References: Searches for the search string in the protein references. Sequences: Searches for the specified amino acid sequence within the protein sequences.
Maximum Number of Matches Reported	Specifies the maximum number of matching references or sequences to report.
Start Search	Begins the search for the reference or sequence.
Save/Add Selected to Database	Adds a selected protein result row to a specified database.
Stop Search	Suspends the active search.

Table 8. Find Protein References page parameters (Sheet 2 of 2)

Parameters	Description
Reference: List of Operators	Specifies the operator for the condition: <ul style="list-style-type: none"> • Equals • Does Not Equal • Less Than • Less Than or Equal To • Greater Than • Greater Than or Equal To • Like • Matches Regular Expression • Starts With • Contains • Ends With • Does Not Start With • Does Not Contain • Does Not End With • Does Not Match • Not Like
Reference: List of Condition Settings	Custom: Opens the Enter Filter Criteria for Reference dialog box so that you can define additional conditions.
Clear Reference Filter icon ()	Deletes the condition under “Reference” on the Find Protein References page.
Clear All Filters icon ()	Deletes all conditions in the Enter Filter Criteria for Reference dialog box and the condition under “Reference” on the Find Protein References page.

Managing FASTA Indexes

A FASTA index is a type of lookup table containing masses, theoretical peptide sequences, and associated proteins, which minimizes search time. The index lists all possible amino acid sequences that can be produced when an enzyme digests a protein or peptide. The peptide fragments are listed by molecular weight. The index stores information about every nominal mass, every peptide that has that mass, every protein that contains this peptide, and the location of its protein description in the FASTA file. Rather than read all protein sequences from the FASTA file, digest them in silico with the specified enzyme, calculate the mass of each peptide, and compare it to the given precursor mass, the Proteome Discoverer application looks for the specific mass in the FASTA index and uses it to find the peptides that have this mass and the associated proteins that contain the peptides.

For full enzymatic searches, the Proteome Discoverer application automatically creates FASTA indexes as they are needed. It does not automatically create FASTA indexes during semi-enzymatic or no-enzyme searches because these searches usually consume a large amount of space on a computer's hard disk. However, you can manually create FASTA indexes for these types of searches.

- [Specifying the Location and Number of FASTA Indexes Stored](#)
- [Displaying the FASTA Indexes View](#)
- [Specifying the Columns to Display](#)
- [Manually Creating a FASTA Index](#)
- [Controlling Automatic FASTA Index Removal](#)
- [Deleting a FASTA Index](#)
- [Changing Number and Location of Stored FASTA Indexes](#)
- [Removing FASTA Indexes When a FASTA File Is Deleted](#)

Specifying the Location and Number of FASTA Indexes Stored


If you do not want to store the FASTA indexes in the default directory shown in [Figure 137](#), you can specify an alternate directory in the FASTA Indexes configuration view. You can also change the maximum number of FASTA indexes stored.

❖ To specify the location and number of the FASTA indexes stored

1. Choose **Administration > Server Settings > FASTA Indexes**.


The configuration view shown in [Figure 137](#) appears.

Figure 137. FASTA Indexes configuration view




2. In the New Directory box, browse to the location of the folder to store the FASTA indexes in.
3. In the New Maximum Number of FASTA Indexes box, select the maximum number of FASTA indexes to store.

If you generate more FASTA indexes than the number to store in the New Maximum Number of FASTA Indexes box, the Proteome Discoverer application discards the difference from the oldest FASTA indexes the next time that you restart the application.


4. If you changed any settings, click the **Apply** icon,  .

The message box shown in [Figure 138](#) appears.

Figure 138. Administration message box



5. Click **OK**.

Note Click  **Reset** to return to the default values.

6. Restart your machine.

Displaying the FASTA Indexes View


You can access FASTA indexes through the FASTA Indexes view.

❖ To display the FASTA Indexes view

1. Choose **Administration > Maintain FASTA Indexes**.


The FASTA Indexes view appears, as shown in [Figure 139](#).

Figure 139. FASTA Indexes view



- Click the plus (+) sign to the left of a database name to vertically display the settings for that database, as shown for the uniprot.fasta database in Figure 140.

Figure 140. Database settings in the FASTA Indexes view



Specifying the Columns to Display

Use the Field Chooser to specify the columns that you want to display.

❖ To set the columns that you want to display



1. Click the **Field Chooser** icon, .
2. In the Field Chooser dialog box, shown in [Figure 141](#), select the check boxes corresponding to the columns that you want to display in the FASTA Indexes view.

Figure 141. Field Chooser dialog box in the FASTA Indexes view




The Proteome Discoverer application instantly makes the selected columns visible and the cleared columns invisible. For a description of these columns, see the Help.

Manually Creating a FASTA Index

The Proteome Discoverer application automatically creates FASTA indexes for a full enzymatic digestion during a Sequest HT search, if an adequate FASTA index does not already exist. You can manually create a FASTA index for a semi-enzymatic or non-specific digestion (see “[Controlling Automatic FASTA Index Removal](#)” on [page 205](#)).


You can only create a specific FASTA index once.

❖ To manually create a FASTA index

1. Choose **Administration > Maintain FASTA Indexes**.
2. Click the **Add** icon,  **Add**.

The FASTA Index Creator dialog box appears, as shown in [Figure 142](#).

Figure 142. FASTA Index Creator dialog box



3. In the General section, specify whether the available FASTA indexes will be removed from memory after the number of indexes reaches the specified maximum.

- True (Default): Automatically removes the FASTA indexes from memory.
- False: Keeps the FASTA indexes in memory.

For information about how the Proteome Discoverer application removes FASTA indexes after the maximum has been reached, see [“Controlling Automatic FASTA Index Removal”](#) on page 205. For instructions on specifying the maximum number of indexes, see [“Changing Number and Location of Stored FASTA Indexes”](#) on page 207.

4. In the Input Data section, specify the basic information that the Proteome Discoverer application needs to create the index:

- FASTA File: Select the FASTA database to be indexed from the list.
- Enzyme Name: Select the enzyme used in the digestion from the list on the left (the enzymes on this list are set in the Cleavage Reagents window) and the type of digestion from the list on the right:
 - Full: Specifies a full enzymatic digestion.
 - Semi: Specifies semi-enzymatic digestion.
 - Unspecific: Specifies a non-specific digestion.
 - No Cleavages: Specifies that no cleavages occur.

- **Maximum Missed Cleavage Sites:** Specifies the maximum number of internal cleavage sites per peptide fragment that is acceptable for an enzyme to miss when cleaving peptides during digestion. Normally the digestion time is too short to enable the enzyme to cleave the peptide at all positions, so you must specify the number of missed positions in one resulting peptide fragment where the enzyme could cleave but did not.

The minimum value is 0, and the maximum value is 12. The default is 2.

5. In the Mass Range Settings section, set the limits of the mass range of the singly charged precursor ion to be processed:

- **Minimum Precursor Mass:** Specifies the minimum mass of the precursor ion. The minimum value is 0.0 Da, and the maximum value is 10000.0 Da. The default is 350 Da.
- **Maximum Precursor Mass:** Specifies the maximum mass of the precursor ion. The minimum value is 0.0 Da, and the maximum value is 10000.0 Da. The default is 5000 Da.
- **Use Average Precursor Mass:** Determines whether the average mass is used to match the precursor ion.
 - True: Uses the average mass to match the precursor ion.
 - False (Default): Uses the monoisotopic mass to match the precursor ion, which is the mass of the most abundant isotope of the protein, peptide, or fragment ion.


6. In the Static Modifications area, specify the static modifications that occur on the amino acid:


- **Peptide N-Terminus:** Select the static modification that occurs on the N terminus of the peptide.
- **Peptide C-Terminus:** Select the static modification that occurs on the C terminus of the peptide.
- **Static Modification:** Select the static modification that occurs on the amino acid side chain.

7. Click **OK**.

The Proteome Discoverer application starts creating the FASTA index, and the job queue appears, as shown in [Figure 143](#).


Figure 143. Creation of the FASTA index displayed in the job queue



- When the job finishes, choose **Administration > Maintain FASTA Indexes** to display the FASTA Indexes view.
- In the FASTA Indexes view, click the **Refresh** icon,  Refresh .

The new FASTA index appears in the FASTA Indexes view on the Administration page, as shown in Figure 144. The Proteome Discoverer application creates an index for the specified FASTA file and the decoy version of the FASTA file.

Figure 144. FASTA Indexes view




Controlling Automatic FASTA Index Removal

After the number of FASTA indexes reaches the specified maximum, the Proteome Discoverer application automatically removes from memory the number of FASTA indexes over the maximum. It first removes the oldest indexes (that is, the ones with the earliest access time). However, you can mark specific FASTA indexes so that they will not be removed from memory, even after the maximum is reached.


❖ To deactivate automatic FASTA index removal

1. In the FASTA Indexes view on the Administration page, clear the **Auto Remove** check box.

The Apply icon now becomes available.

2. Click the **Apply** icon,  **Apply**.
3. In the Remove FASTA indexes confirmation box, click **OK**.

❖ **To activate automatic FASTA index removal**

1. Select the **Auto Remove** check box.
2. Click the **Apply** icon,  Apply.
3. In the Remove FASTA indexes confirmation box, click **OK**.



Deleting a FASTA Index

You can only delete FASTA indexes that have the Auto Remove check box selected.

❖ **To delete a FASTA index**


1. Be sure that the Auto Remove check box is selected for the index that you want to delete.
2. Select the index that you want to delete by clicking the box next to the plus (+) sign.

The box now changes to the  arrow.



3. Click the  arrow.
4. Click the **Remove** icon,  Remove .
5. Click **OK** in the Remove FASTA indexes confirmation box.

The name of the deleted index disappears from the FASTA Indexes table and reappears in a separate table called Deleted FASTA Indexes, as shown in [Figure 145](#). It no longer appears in the FASTA Indexes table. However, because the FASTA index might be used in some calculations, its removal from the application only takes place the next time that the server starts.

Figure 145. Deleted FASTA Indexes table



❖ To restore a deleted FASTA index

1. In the Deleted FASTA Indexes table, select the deleted index by clicking the  arrow.
2. Click the **Restore** icon, .
3. In the Restore FASTA indexes confirmation box, click **OK**.

The restored index appears in the FASTA Indexes table and disappears from the Deleted FASTA Indexes table.

Changing Number and Location of Stored FASTA Indexes


You can specify a new directory for storing the FASTA indexes and change the maximum number of FASTA indexes stored. The Proteome Discoverer application counts all FASTA indexes, even the indexes that cannot be automatically removed with the Auto Remove option.

❖ **To change the number and location of stored FASTA indexes**

1. Click the **Options** icon,  **Options** .

The FASTA Indexes Options dialog box appears, as shown in [Figure 146](#).

Figure 146. FASTA Indexes Options dialog box



Note Another way to access these options is to choose Administration > Configuration and click FASTA Indexes in the Server Settings area.

The FASTA Indexes Options dialog box contains two read-only parameters:

- The FASTA Index Directory box displays the name of the current directory where the FASTA indexes are saved.
 - The Maximum Number of FASTA Indexes box displays the current maximum number of FASTA indexes allowed.
2. In the New Directory box, browse to the directory where you want to store the FASTA indexes.

You can change the directory only if the server runs on the local machine.
 3. In the New Maximum Number of FASTA Indexes box, type the new maximum number of FASTA indexes allowed.
 4. Click **OK**.
 5. In the FASTA index settings confirmation box, click **OK**.

After you confirm the changes, the Proteome Discoverer application saves them, but the changes are only executed the next time that the server starts. You can undo the changes made since the last time that the server started and before the next time that the server starts, even though you clicked OK in the FASTA Indexes Options dialog box and closed it. For example, when you change the location of the directory in the FASTA Indexes Options dialog box, click OK, and close the dialog box, the server moves all FASTA indexes to the new target directory when the server restarts. But if you reinvoke the dialog box and click Reset before restarting the server, the changes that you made previously are deleted, and the directory reverts to its previous location.

❖ **To reset the changes made in a previous FASTA index session**

1. Click the **Options** icon,  Options .

The FASTA Indexes Options dialog box appears, as shown in [Figure 146](#).

2. Click **Reset**.

Removing FASTA Indexes When a FASTA File Is Deleted

When you or the Proteome Discoverer application deletes a FASTA file, the application removes the FASTA indexes belonging to the deleted FASTA file the next time that the server starts.

Adding or Modifying FASTA Parsing Rules

The protein descriptions and accessions shown on the Proteins page of the .pdResult file are taken from the title lines of the uploaded FASTA file. The rules for extracting these values are defined as regular expressions (see http://en.wikipedia.org/wiki/Regular_expression). If you want to use unsupported FASTA files that might contain descriptions or accessions that are difficult to read—for example, FASTA files from the Saccharomyces Genome Database (SGD) or the Arabidopsis Information Resource (TAIR) web pages—you might need to add new parsing rules to the system or modify existing parsing rules.


If you want to add or modify a FASTA parsing rule, download a FASTA file from an appropriate website, add or change the parsing rule, test the parsing rule, correct the parsing rule until it meets your needs, and then apply it.

❖ **To add or modify FASTA parsing rules**

Choose **Administration > Maintain FASTA Parsing Rules**.

The FASTA Parsing Rules view opens, as shown in [Figure 147](#).


Figure 147. FASTA Parsing Rules view




The FASTA Parsing Rules view includes the following features:

- Parsing rule category selection area: Displays the four categories (Title Line Rules, Accession Rules, Taxonomy Rules, and Avoid Expression Rules) into which the Proteome Discoverer application groups parsing rules.
- Parsing rule text area: Displays the regular expression of the parsing rule. Each line corresponds to a single regular expression. These expressions are tested as alternatives (“or” connected). Figure 148 shows regular expressions for SwissProt accessions in the parsing rules text area.
- List of parsing rules in the selected category: Displays all available parsing rules in the selected category. This list corresponds to the available values of the appropriate parameter of the MSF Files node. If you select a single entry, the FASTA Parsing Rules view displays the parsing rule in the parsing rule area on the right.
- Test area: Loads the title lines of a sample FASTA file to test the matching of the expression.

Figure 148. Regular expressions for SwissProt accessions in the parsing rule text area




❖ **To add a regular expression**

1. Click the **Add** icon,  .
2. Enter the new parsing rule as a regular expression in the parsing rule text area, and change the preset rule name, if needed.

You can find a guide explaining the syntax of regular expressions on the Microsoft Developer Network (MSDN) home page (Regular Expression Language - Quick Reference page) at <http://msdn.microsoft.com/en-us/library/az24scfc.aspx>. The title line rules and the accession rules use named capture groups to capture the accession and description of the protein.

- Title line rules (Title Line Rules area) must specify a named capture group, AC1, for the protein accession and a group, Desc1, for the protein description.
- Accession rules (Accession Rules area) require a capture group, AC, for the accession.
- Taxonomy rules (Taxonomy Rules area) require a matched taxonomy string.
- Avoid expression rules (Avoid Expression Rules area) require a matched string.

❖ **To change or rename a regular expression**

Select the parsing rule in the list, and click on the **Edit** icon,  .


You can now edit the rule in the parsing rule text area.

❖ **To test a parsing rule**

1. Click **Load FASTA File** to load the FASTA file for which the rule is designed.

The first and last five title lines of the file are displayed in the text box of the test area, as shown in [Figure 149](#).

Figure 149. Loaded FASTA file from UniProt with first and last five title lines in the test area



2. Click **Test Rule** to test the selected parsing rule.

The Proteome Discoverer application applies the newly defined parsing rule to the small sample set of FASTA title lines read from the file. It highlights the matched values with a colored background, as shown in the following examples. If the application cannot match anything in the FASTA title line, it displays [No Match].

Test of the standard title line rule:

```
>sp|P31946|1433B_HUMAN 14-3-3 protein beta/alpha OS=Homo sapiens GN=YWHAB PE=1 SV=3
>sp|P62258|1433E_HUMAN 14-3-3 protein epsilon OS=Homo sapiens GN=YWHAE PE=1 SV=1
```

Test of the SwissProt accession rule:

```
>sp|P31946|1433B_HUMAN 14-3-3 protein beta/alpha OS=Homo sapiens GN=YWHAB PE=1 SV=3
>sp|P62258|1433E_HUMAN 14-3-3 protein epsilon OS=Homo sapiens GN=YWHAE PE=1 SV=1
```

Test of the human taxonomy rule:


```
>sp|P31946|1433B_HUMAN 14-3-3 protein beta/alpha OS=Homo sapiens GN=YWHAB PE=1 SV=3
>sp|P62258|1433E_HUMAN 14-3-3 protein epsilon OS=Homo sapiens GN=YWHAE PE=1 SV=1
```

Test of the refseq-id accession rule with no match:

```
[No Match] >sp|P31946|1433B_HUMAN 14-3-3 protein beta/alpha OS=Homo sapiens GN=YWHAB PE=1 SV=3
[No Match] >sp|P62258|1433E_HUMAN 14-3-3 protein epsilon OS=Homo sapiens GN=YWHAE PE=1 SV=1
```

4. Correct the parsing rule until it meets your needs.

❖ **To apply the changes**

To apply the changes, click the **Apply** icon,  **Apply**.

After you apply the changes, the changed or newly defined parsing rules are available in the MSF Files node. For more information on the MSF Files node, refer to the Help.

Using the MSF Files Node to Change Title Line Parsing

You can use the MSF Files node to determine how the Proteome Discoverer application parses the FASTA title lines of found proteins. When the application parses these title lines, it applies a set of predefined parsing rules to extract the accession and description of the protein. If it finds a protein in more than one results file from which it generates a report, and if the FASTA files of the two reports differ, the application displays the first available description and accession on the Proteins page. You can view all other accessions and descriptions when you move the mouse over the cells in the corresponding columns, as shown in [Figure 150](#).

Figure 150. Available accessions and descriptions

Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info
	Accession	Description			
1	17543850	Ribosomal protein, Large subunit, Acidic (P1) family membe			
2	17531551	C06A8.3 [Caenorhabditis elegans]			
3	25150244	Ribosomal protein, Large subunit, Acidic (P1) family membe			
4	25141400	C37A2.7 [Caenorhabditis elegans]			
5	17556226	ECTin family member (lec-6) [Caenorhabditis elegans]			
6	17542014	ribosomal Protein, Small subunit family member (rps-8) [Caenorhabditis elegans]			
7	1754310	Malate DeHydrogenase family member (mdh-1) [Caenorhabditis elegans]			
8	17536383	ENOLase family member (enol-1) [Caenorhabditis elegans]			

Accession

Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info					
	Accession	Description				Σ Coverage	Σ # Peptides	Σ # PSMs	# AAs	MW [kDa]
1	17543850	Ribosomal protein, Large subunit, Acidic (P1) family membe				70.00 %	4	6	110	10.9
2	17531551	C06A8.3 [Caenorhabditis elegans]				46.36 %	3	7	151	15.9
3	25150244	Ribosomal protein, Large subunit, Acidic (P1) family membe				45.95 %	2	4	111	11.3
4	25141400	C37A2.7 [Caenorhabditis elegans]				42.99 %	2	5	107	10.8
5	17556226	galactose 4-epimerase [Caenorhabditis elegans]								16.0
6	17542014	Ribosomal protein, Large subunit, Acidic (P1) family member (rla-1) [Caenorhabditis elegans]								23.7
7	17554310	Malate DeHydrogenase family member (mdh-1) [Caenorhabditis elegans]								35.1
8	17536383	ENOLase family member (enol-1) [Caenorhabditis elegans]				32.72 %	8	14	434	46.6

Description

The MSF Files node has advanced parameters that you can set to change the title line parsing. Changing these settings affects how the Proteome Discoverer application displays protein accession and description information in the report. You can use the node's Title Line Rule parameter to define an alternative parsing rule to use to parse the accession and description. If the specified rule does not match anything in the FASTA title line, the Proteome Discoverer application tries the standard parsing rules to match the protein accession and description.

In cases where multiple accessions and descriptions are available for a given protein (see [Figure 150 on page 213](#)), you can determine the order of the displayed accessions and descriptions by using three other node parameters. The application applies them in this order:

Note The Proteome Discoverer application applies the following rules only if multiple equally scored proteins are returned for the same set of PSMs, which rarely happens. Otherwise, it marks the highest-scoring protein reference as the master protein and displays it as such on the Protein Groups page. Even then, the protein marked as master is typically the longest of the group of equally scored protein references (the others are alternate master protein candidates).

1. Preferred Accession, which allows you to select a parsing rule to extract the preferred protein accession from the FASTA entry. If the Proteome Discoverer application finds a preferred accession, it displays it instead of the primary accession.

If you select a rule for the Preferred Accession parameter that matches one of the accessions, the Proteome Discoverer application moves this accession and description to the first position in the list of available accessions and descriptions.

2. Preferred Taxonomy, which allows you to select a parsing rule to extract the preferred taxonomy from the FASTA entry. If the Proteome Discoverer application finds a preferred taxonomy, it displays the accession and description of this entry, except when an entry containing a preferred accession is better than an entry containing preferred taxonomy and no preferred accession.

If you performed the search without a preferred taxonomy and the application identifies proteins with the same sequence from different species, you can select a rule with the Preferred Taxonomy parameter to display accession and description from the right species.

Here is an example showing the precedence of the Preferred Taxonomy parameter. Suppose that you have one protein with more than one accession and description. Both descriptions contain some common keywords, for example:

Description 1: xxxxxxxxxx human abc

Description 2: yyyyyyyyyy human abc

You set the Preferred Taxonomy parameter to a rule that includes “human.” Then you set the Avoid Expressions parameter to a rule that includes “human.”

In this case, the Preferred Taxonomy parameter has a higher precedence than the Avoid Expression parameter.

3. Avoid Expressions, which allows you to select the terms that the Proteome Discoverer application should avoid when parsing the protein description. If more than one description is available, the application prefers the description containing none of the specified terms.

Some of the publicly available protein databases like UniProt or NCBI collect experimental verified and curated proteins as well as unverified proteins, for example, the output from bioinformatic algorithms predicting potential proteins in a sequenced genome. In many cases, these unverified proteins contain words like “predicted” or “hypothetical” in the description. The Avoid Expressions parameter matches these words in the description. Where there are two different descriptions from two different databases used in a complex search, the Proteome Discoverer application displays the description that does not contain a word that you selected with the Avoid Expressions parameter.

The Proteome Discoverer application only applies the Avoid Expressions parameter if there are different FASTA title lines for the same protein that conform to the other title line rules described in this topic.

Here is an example showing how the Avoid Expressions parameter works. Suppose that you have one protein with more than one accession and description. Both descriptions contain some common keywords, for example:

Description 1: xxxxxxxxxxx human abc

Description 2: yyyyyyyyyy human abc

You set the Avoid Expressions parameter to a rule that includes “abc.”

In this case, the Proteome Discoverer application takes the first accession in the list.

The application saves the parsing rules for the FASTA title lines in the FastaTitleParsingRules.xml file, which is stored in the C:\ProgramData\Thermo\Proteome Discoverer 2.1\MagellanDBs folder or equivalent. This file contains the parsing rules, along with name and meta information. A parsing rule is a list of regular expressions. If the application uses the parsing rule, it applies all regular expressions in the list in the order of the list, starting with the first. It uses the first rule that matches the title line to read out the accession and, if declared, the description.

The file contains a section for each of the four parsing rule parameters of the MSF Files node. The basic *.xml document with the four sections without rules looks like this:

```
<?xml version="1.0" encoding="utf-16"?>
<FastaTitlelineRules
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xmlns:xsd="http://www.w3.org/2001/XMLSchema">
  <TitlelineRules>
    ... enter rules here ...
  </TitlelineRules>
  <AccessionRules>
    ... enter rules here ...
  </AccessionRules>
  <TaxonomyRules>
```

```
... enter rules here ...  
</TaxonomyRules>  
<AvoidExpressionRules>  
... enter rules here ...  
</AvoidExpressionRules>  
</FastaTitleLineRules>
```

The parsing rules with their lists of regular expressions are defined as follows as in this example showing an accession rule to match SwissProt accessions:

```
<ParsingRule name="swissprot" isVisible="true" changeable="true">  
<RuleParts>  
<RulePart>\|SWISS-PROT:(?&lt;AC&gt;[A-N,R-Z][\d][A-Z][A-Z,\d][A-Z,\d][\d])</RulePart>  
<RulePart>\|SWISS-PROT:(?&lt;AC&gt;[O,P,Q][\d][A-Z,\d][A-Z,\d][A-Z,\d][\d])</RulePart>  
<RulePart>sp\|(?&lt;AC&gt;[A-N,R-Z][\d][A-Z][A-Z,\d][A-Z,\d][\d])\|  
</RulePart>  
  
<RulePart>sp\|(?&lt;AC&gt;[O,P,Q][\d][A-Z,\d][A-Z,\d][A-Z,\d][\d])\|  
</RulePart>  
</RuleParts>  
</ParsingRule>
```

The first line of the rule defines the name that is displayed in the node parameter. It specifies whether the rule is visible in the list for the parameter. For new rules, always set `isVisible` and `changeable` to `true`; otherwise, it is impossible to apply the rule.

The body of the rule is a list of statements containing one or more elements named `RulePart`. All regular expressions are “or” connected in the final parsing rule. For accession rules, a named capture group, `(?<AC>)`, must match the accession. The Proteome Discoverer application evaluates it to extract the accession for display. You must use `<` and `>` instead of `<` and `>` because the `<` and `>` are not allowed in the XML entry. The title line rules contain two named capture groups, `?<AC1>` and `?<DESC1>`, for accession and description, respectively.

For information on adding or modifying FASTA parsing rules, see [“Adding or Modifying FASTA Parsing Rules”](#) on page 209.

Identifying Contaminants During Searches

In many cases, a proteomics sample contains contaminant proteins from sample preparations that are identified in the search. These contaminants have no impact on the scientific interpretation of the data, so you want to remove them from the displayed results file. You might also want to remove proteins from the sample preparation (for example, antibodies used in immunoprecipitation or proteins from affinity purification) and proteins that are known carryovers from the previous experiment (for example, peptides that still adhere to the HPLC column after the washing procedure).

Suppose that a core facility receives a sample from a microbiological laboratory and a request to search the data against the proteins of all known bacteria. If the preceding experiment involved a proteome of *E. coli*, it is likely that the HPLC column retains some *E. coli* peptides and that proteins are detected on the basis of these peptides. It is helpful if you can specify the *E. coli* FASTA file as an additional contaminant database and have the proteins highlighted in an additional column.

You can use a consensus workflow that includes the Protein Marker node to mark as contaminants all proteins in the results file that are listed in a designated FASTA file or files as contaminants. You can specify an original FASTA file and up to three additional FASTA files containing contaminants. You can also specify the name of the column in the results report that displays an X symbol marking the proteins.

Before you begin creating a workflow, upload the FASTA file or files containing the proteins marked as contaminants to the Proteome Discoverer application. For instructions, see [“Adding FASTA Files to the Proteome Discoverer Application”](#) on page 174.

❖ To create a marked contaminants processing workflow

1. Create or open a study and an analysis:
 - To create a study, see [“Creating a Study”](#) on page 38.
 - To open an existing study, see [“Opening an Existing Study”](#) on page 40.
 - To create an analysis, see [“Creating an Analysis”](#) on page 71.
 - To open an existing analysis, see [“Opening an Existing Analysis”](#) on page 72.
2. Follow the general instructions for creating a processing workflow with the Workflow Editor. See [“Creating a Processing Workflow”](#) on page 104.

The basic processing workflow for identifying contaminants is the same as the basic processing workflow shown in [Figure 74](#) on page 107.


❖ To create a marked contaminants consensus workflow

1. Create a basic consensus workflow as described in [“Creating a Consensus Workflow”](#) on page 112, but attach the Protein Grouping node to the ProteinScorer node, and attach the Protein Marker node directly to the Peptide and Protein Filter node, as shown in [Figure 151](#).

4 Searching for Data

Identifying Contaminants During Searches

Figure 151. Basic consensus workflow for marking proteins as contaminants



2. Set the parameters for the nodes, including the Protein Marker node. For information on the Protein Marker node parameters, see the Help.
3. Run the workflow.

The results report contains a Contaminant column on the Proteins, Peptide Groups, and PSMs pages that displays an X symbol next to the proteins marked as contaminants in the searched FASTA file or files. [Figure 152](#) shows the Contaminant column on the Proteins page.

Figure 152. Contaminant column in results file

Protein Groups															
Proteins															
Peptide Groups															
PSMs															
MS/MS Spectrum Info															
	Checked	Accession	Description	Σ Coverage	Σ # Peptides	Σ # PSMs	# AAs	MW [kDa]	calc. pI	Is Master Protein	Coverage A2	# PSMs A2	# Peptides A2	Score A2	Contaminant
1	<input type="checkbox"/>	16131810	protein chain elongation factor EF-Tu (duplicate of tufA) [Escherichia coli K12]	5%	1	2	394	43.3	5.45	✓	5%	2	1	4.64	
2	<input type="checkbox"/>	90111711	ketoacid-binding protein [Escherichia coli K12]	10%	1	1	128	13.6	5.50	✓	10%	1	1	3.42	X
3	<input type="checkbox"/>	16131632	ketol-acid reductoisomerase [Escherichia coli K12]	2%	1	2	491	54.0	5.31	✓	2%	2	1	3.22	
4	<input type="checkbox"/>	16131678	5-methyltetrahydropteroyl-triglutamate-homocysteine methyltransferase [Escherichia coli K12]	1%	1	1	753	84.6	5.92	✓	1%	1	1	2.56	X
5	<input type="checkbox"/>	16130546	CP4-57 prophage; predicted protein [Escherichia coli K12]	3%	1	1	538	62.0	5.91	✓	3%	1	1	2.44	
6	<input type="checkbox"/>	24054939	acyltransferase for 30S ribosomal subunit protein S18 [Shigella flexneri]	8%	1	2	148	16.6	4.97	✓	8%	2	1	2.34	

Contaminant column

Displaying Species Names for Proteins and Peptide Groups

In some cases, a proteomics sample might contain more than one species, either intentionally or accidentally, as a result of sample preparation. These different species can affect the scientific interpretation of the data, so you might want to see them or filter them in the result report. You can display the species names of the proteins and the peptide groups that they belong to as distribution maps or as text strings separated by semicolons in columns on the Proteins and Peptide Groups pages of the .pdResult file. The application extracts these names from the FASTA title lines of all proteins found in the sequence database search, removes duplicates, and assigns a unique color to each species name if you choose to display the names as a distribution map. [Figure 153](#) shows the protein and peptide group names displayed in the Species Map and the Species columns on the Proteins page of the .pdResult report.

Figure 153. Species names as distribution maps and as text columns

The figure shows two screenshots of the Proteome Discoverer interface. The top screenshot displays the 'Proteins' tab with a table of protein groups. The columns include 'Checked', 'Description', '# Peptides', '# PSMs', 'Species Map', and 'Species'. The 'Species Map' column shows a grid of colored squares representing different species. The 'Species' column lists the species names. The bottom screenshot displays the 'Peptide Groups' tab with a table of peptide groups. The columns include 'Checked', 'Confidence', 'Sequence', '# Protein Groups', '# Proteins', '# PSMs', 'Species Map', and 'Species'. The 'Species Map' column shows a grid of colored squares representing different species. The 'Species' column lists the species names.

Protein Groups	Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info
1525	Heat shock cognate 71 kDa protein (Fragment) OS=Homo	3	4	Homo sapiens
1526	plastin-3 isoform 4 [Homo sapiens]	1	1	Homo sapiens
1527	Activated RNA polymerase II transcription cofactor 4 (Frag	1	1	Homo sapiens
1528	PREDICTED: Y-box-binding protein 2 isoform X3 [Mus mus	2	2	Mus musculus
1529	adenylate kinase 2, mitochondrial isoform b [Mus musculus	1	1	Mus musculus
1530	Uncharacterized protein OS=Bos taurus GN=KRT2 PE=3 S	1	1	Bos taurus
1531	Malate synthase A [Escherichia coli K-12]	1	1	Escherichia coli K-12
1532	Protease 1 OS=Achromobacter lyticus PE=1 SV=1	1	2	Achromobacter lyticus
1533	Tax_id=9606 Gene_Symbol= Similar to Keratin, type II cyt	2	2	Homo sapiens; Rattus norvegicus; Gallus gallus; Mus musculus; Bo
1534	actin, cytoplasmic 1 [Homo sapiens]	1	1	Homo sapiens; Rattus norvegicus; Gallus gallus; Mus musculus; Bo
1535	Glyceraldehyde-3-phosphate dehydrogenase OS=Escheric	1	1	Escherichia coli; Escherichia coli K-12

Protein Groups	Peptide Groups	PSMs	MS/MS Spectrum Info	Annotated Modifications	
1	AVLLGPPGAGK	1	18	1	Rattus norvegicus; Homo sapiens; Gallus gallus; Mus musculus

Protein Groups	Proteins	PSMs	MS/MS Spectrum Info
1	adenylate kinase 2, mitochondrial isoform b [Rattus norvegi	1	1
2	adenylate kinase 2, mitochondrial isoform b [Homo sapiens]	1	1
3	cDNA FLJ55097, highly similar to Adenylate kinase isoerzy	1	1
4	adenylate kinase 2, mitochondrial isoform a [Homo sapiens]	1	1
5	adenylate kinase 2, mitochondrial isoform c [Homo sapiens]	1	1
6	PREDICTED: adenylate kinase 2, mitochondrial isoform X7	1	1
7	adenylate kinase 2, mitochondrial isoform b [Mus musculus]	1	1
8	Adenylate kinase 2, mitochondrial OS=Homo sapiens GN=K	1	1
9	Adenylate kinase 2, isoform CRA_a OS=Homo sapiens GN=	1	1

The application can display up to 30 species names in a distribution map. If it extracts more names, it issues a warning and displays the additional names as text strings separated by semicolons.

Use the As Species Map and As Species Names parameters of the ProteinMarker node to display species names in the result file. For a description of these parameters, see the Help. Use the node in the same workflow as that given in [“Identifying Contaminants During Searches”](#) on [page 217](#). It connects to the Peptide and Protein Filter node.

The annotation is restricted to FASTA database formats of UniProt (Swiss-Prot, TrEMBL), NCBI (RefSeq), and ProteinCenter downloads.

Searching Spectrum Libraries

Spectrum library search is a different search approach from the sequence database search commonly used in shotgun proteomics. The main difference between a database search and a spectrum library search is in the origin of the spectra that the measured spectra from your experiments are compared to. Sequence database searches use theoretical spectra generated from peptide sequences, but spectrum libraries are libraries of measured (consensus) spectra from actual previous experiments. Using a library of already well-identified peptides avoids identifying already known peptides over and over again by a time-consuming database search. Restricting the library to previously identified peptides also drastically reduces the search space and therefore the search time. In addition, comparisons that use consensus spectra consider the measured peak intensities, increasing the selectivity and making the identification more accurate.

You can use the MSPepSearch node to search large spectrum libraries downloaded from the NIST or the PeptideAtlas home page.

- [Displaying Spectrum Libraries](#)
- [Adding a Spectrum Library](#)
- [Deleting a Spectrum Library](#)
- [Searching Spectrum Libraries with the MSPepSearch Node](#)
- [Visually Verifying Spectrum Library Matches](#)

Displaying Spectrum Libraries


You can display a list of all the spectrum libraries that you registered in the Proteome Discoverer application.

❖ To list the available spectrum libraries

Choose **Administration > Maintain Spectrum Libraries**.

The Spectrum Libraries view shown in [Figure 154](#) appears. It lists all the spectrum libraries that you downloaded from NIST or the Peptide Atlas home page and registered. It displays the processed spectrum library properties, such as the file name, file size, the number of proteins stored, and the library type, which determines the search node to use. The Proteome Discoverer application processes the spectrum library and makes it available for use.

Figure 154. Spectrum Libraries view




Adding a Spectrum Library

You must add a spectrum library to the Proteome Discoverer application before you can conduct a search with the MSPepSearch node. In the registration process, the Proteome Discoverer application automatically recognizes the type of the spectral library. The type determines the search node that you can use the library with. Adding the spectrum libraries is similar to the procedure for adding FASTA files.

❖ **To add a spectrum library for searching with the MSPepSearch node**

1. Download the appropriate spectrum libraries from the National Institute of Standards and Technology (NIST) at <http://peptide.nist.gov> or from Peptide Atlas at <http://www.peptideatlas.org/speclib>.

The Proteome Discoverer application recognizes the following file formats for searching spectrum libraries with the MSPepSearch node:


- *.zip/*.gz files from NIST or PeptideAtlas. You can find these files in the *_nist.tar.gz file on the library download site at NIST or the *_nist.zip file on the PeptideAtlas home page. The file must contain a complete spectrum library in MSPepSearch. If files are missing, the Proteome Discoverer application does not add the library.
2. In the Proteome Discoverer application, choose **Administration > Maintain Spectrum Libraries**.
 3. Click the **Add** icon, .
 4. In the Select a Spectrum Library dialog box, do the following:
 - a. In the list box in the lower right corner of the Select a Spectrum Library dialog box, select **All Spectrum Library Files** (*.gz, *.msp, *.zip) or Zip archives (*.gz, *.zip).
 - b. Browse to the location of the spectrum library where you downloaded and unpacked the *_nist.tar.gz file.
 - c. Select the *filename.gz* file.
 - d. Click **Open**.

When you add a spectrum library file, the Proteome Discoverer application takes the following steps:

- Extracts the archive file.
- Extracts spectra for visualization.

During library creation, the job queue in the Administration view displays each step, as shown in [Figure 155](#).

Figure 155. Adding a spectrum library for searching with the MSPepSearch node



When the Proteome Discoverer application finishes adding the spectrum library, the spectrum library file appears in the Spectrum Libraries view, as shown in [Figure 154](#) on [page 221](#).

Now you are ready to search the spectrum library. To search with the MSPepSearch node, see [“Searching Spectrum Libraries with the MSPepSearch Node.”](#) For more information on the MSPepSearch node, refer to the Help.

Deleting a Spectrum Library


You can delete a spectrum library from the application.

❖ To delete a spectrum library

1. Choose **Administration > Maintain Spectrum Libraries**.

The Administration page appears with the Spectrum Libraries view.

2. Click  at the beginning of a row to select the row.
3. Click the **Remove** icon,  **Remove**.
4. In the Remove Spectrum Libraries Databases dialog box, click **OK**.

The Spectrum Libraries file that you selected appears as a job in the job queue. After you start the deletion of the file, you cannot cancel the deletion. You can remove the completed job from the job queue by clicking the Remove icon, , and then clicking OK in the Delete Jobs dialog box.

Searching Spectrum Libraries with the MSPepSearch Node

You can use the MSPepSearch node as an alternative to a search node such as Sequest HT.

The spectrum library search reports the three scores shown in [Table 9](#). Dot score and reversed dot score are secondary scores, and their values are not shown by default.

Table 9. Scores generated by the MSPepSearch node


Score	Description
MSPepSearch	Is the main score of MSPepSearch.
Dot score	Is the score from a cross-correlation computed between two spectra.
Reverse dot score	Is the reversed spectral dot product.

❖ To create a processing workflow for searching spectrum libraries with the MSPepSearch node

Follow the instructions in “[Creating a Processing Workflow](#)” on [page 104](#), but substitute the MSPepSearch node for the search engine node.

[Figure 156](#) shows the basic processing workflow for searching spectrum libraries with the MSPepSearch node.

Figure 156. Basic processing workflow using the MSPepSearch node to search spectrum libraries



The Fixed Value PSM Validator node is the only possible peptide validator for the MSPepSearch node. It is impossible to perform a decoy search because there is no proper decoy spectrum library.

For a description of the parameters available in the MSPepSearch node, refer to the Help.


❖ **To create a consensus workflow for searching a spectrum libraries with the MSPepSearch node**

Follow the instructions in [“Creating a Consensus Workflow”](#) on page 112.

Visually Verifying Spectrum Library Matches

You can visually verify matches between measured spectra from your experiment and the reference spectra in the spectrum library for peptides identified with the MSPepSearch node. In the Peptide Spectrum Match Identification Details view, you can display a mirror plot of the matching peptides, as shown in [Figure 157](#). You can use the reference spectrum with the fragment match settings (refer to the Help).

Figure 157. Mirror plot in the Peptide Spectrum Match Identification Details view



The Proteome Discoverer application displays the reference spectrum using intensities multiplied by -1 in the same plot as the measured spectrum. In the reference spectrum, it also labels peaks of the a, b, c, ion series and the x, y, and z ion series, as well as the peaks from the precursor peptide. It does not display labels for all fragments with a mass difference, isotope peaks, and “?” peaks in the spectrum library.

❖ **To generate a mirror plot**

1. Open the .pdResult file for the results of the spectrum library search performed with the MSPepSearch node.
2. Follow the instructions for generating a Peptide Spectrum Match Identification Details view given in the Help.

Defining Chemical Modifications

You can update the chemical modifications that you use to conduct a peptide identification search. The available modifications are defined in the Chemical Modifications view on the Administration page that opens when you choose Administration > Maintain Chemical Modifications. Use this view to customize the chemical modifications that you use to do your search. You can import a new list or the latest UNIMOD list. You can also modify the chemical modification list by adding amino acids to the modifications, creating new modifications, or activating or deactivating existing modifications.

Note A modification must be active to be usable during a search.

- [Dynamic Modifications](#)
- [Static Modifications](#)
- [Opening the Chemical Modifications View](#)
- [Adding Chemical Modifications](#)
- [Adding Amino Acids](#)
- [Deleting Amino Acids](#)
- [Importing Chemical Modifications](#)

The Proteome Discoverer application offers two types of modifications, dynamic and static.

Dynamic Modifications

Dynamic modifications, also known as variable amino acid modifications, are modifications that might or might not be present. They are mainly used for determining post-translational modifications (PTMs). For example, some phosphorylated peptide serines are modified, and some are not modified.

Static Modifications

Static modifications apply the same specific mass to all occurrences of that named amino acid, as in an exhaustive chemical modification.

A static modification might result from derivatization or isotopic labeling of an amino acid. For example, carbamidomethylated cysteine has a delta mass of 57.021464, which is added to each cysteine residue appearing in a protein.


In static searches, the Proteome Discoverer application assumes that every amino acid residue will be modified in that way. Constant mass is changed. The search engines perform static modification searches by adding the specified constant value to the mass of the specified amino acid.

Opening the Chemical Modifications View

The Chemical Modifications view is an advanced feature of the Proteome Discoverer application. You use it to build and maintain the static and dynamic modifications data that is available when you define your search settings.

In the Chemical Modifications view, you can explore the default types of modifications and their corresponding amino acids. It contains the modification's delta mass, amino acids, and substitutions. By using the Chemical Modifications view, you can add amino acids to existing modifications and create new modifications.

❖ To open the Chemical Modifications view

1. Choose **Administration > Maintain Chemical Modifications**, or click the **Maintain Chemical Modifications** icon, , either on the toolbar or on the Administration page.

The Chemical Modifications view appears on the Administration page, as shown in [Figure 158](#). The amino acids listed are those where the modifications can appear.

Figure 158. Chemical Modifications view

Is Active	Modification	Abbreviation	Delta Mass	Delta Average Mass	Substitution	Leaving Group	Pos
<input checked="" type="checkbox"/>	Acetyl	Acetyl	42.010565	42.0367	H(2) C(2) O		Any
<input checked="" type="checkbox"/>	Acetyl	Acetyl	42.010565	42.0367	H(2) C(2) O		Protei
<input checked="" type="checkbox"/>	Acetyl	Acetyl	42.010565	42.0367	H(2) C(2) O		Any_↑
<input checked="" type="checkbox"/>	Amidated	Amidated	-0.984016	-0.9848	H N O(-1)		Protei
<input checked="" type="checkbox"/>	Amidated	Amidated	-0.984016	-0.9848	H N O(-1)		Any_↓
<input checked="" type="checkbox"/>	Carbamidomethyl	Carbamidomethyl	57.021464	57.0513	H(3) C(2) N O		Any
<input checked="" type="checkbox"/>	Carbamidomethyl	Carbamidomethyl	57.021464	57.0513	H(3) C(2) N O		Any_↑
<input checked="" type="checkbox"/>	Carbamyl	Carbamyl	43.005814	43.0247	H C N O		Any
<input checked="" type="checkbox"/>	Carbamyl	Carbamyl	43.005814	43.0247	H C N O		Any_↑
<input checked="" type="checkbox"/>	Carbamyl	Carbamyl	43.005814	43.0247	H C N O		Protei
<input checked="" type="checkbox"/>	Carboxymethyl	Carboxymethyl	58.005479	58.0361	H(2) C(2) O(2)		Any_↑
<input checked="" type="checkbox"/>	Carboxymethyl	Carboxymethyl	58.005479	58.0361	H(2) C(2) O(2)		Any
<input checked="" type="checkbox"/>	Deamidated	Deamidated	0.984016	0.9848	H(-1) N(-1) O		Any
<input checked="" type="checkbox"/>	Deamidated	Deamidated	0.984016	0.9848	H(-1) N(-1) O		Protei
<input checked="" type="checkbox"/>	dHex(1)Hex(5)HexN...	dHex(1)Hex(5)HexN...	2059.734933	2060.8689	dHexHex(5) He...		Any
<input checked="" type="checkbox"/>	dHex(1)Hex(5)HexN...	dHex(1)Hex(5)HexN...	2350.83035	2352.1234	dHexHex(5) He...		Any
<input checked="" type="checkbox"/>	Dimethyl	Dimethyl	28.0313	28.0532	H(4) C(2)		Any
<input checked="" type="checkbox"/>	Dimethyl	Dimethyl	28.0313	28.0532	H(4) C(2)		Any_↑
<input checked="" type="checkbox"/>	Dimethyl	Dimethyl	28.0313	28.0532	H(4) C(2)		Protei
<input checked="" type="checkbox"/>	Dimethyl:2H(4)	Dimethyl:2H(4)	32.056407	32.0778	2H(4) C(2)		Protei
<input checked="" type="checkbox"/>	Dimethyl:2H(4)	Dimethyl:2H(4)	32.056407	32.0778	2H(4) C(2)		Any_↑
<input checked="" type="checkbox"/>	Dimethyl:2H(4)	Dimethyl:2H(4)	32.056407	32.0778	2H(4) C(2)		Any
<input checked="" type="checkbox"/>	Dimethyl:2H(6)13C(2)	Dimethyl:2H(6)13C(2)	36.07567	36.0754	H(-2) 2H(6) 13C...		Any
<input checked="" type="checkbox"/>	Dimethyl:2H(6)13C(2)	Dimethyl:2H(6)13C(2)	36.07567	36.0754	H(-2) 2H(6) 13C...		Any_↑
<input checked="" type="checkbox"/>	dNIC	dNIC	109.048119	109.1205	H 2H(3) C(6) N O		Any
<input checked="" type="checkbox"/>	EDEDTIDVFQQQTGG	EDEDTIDVFQQQTGG	1662.700924	1663.6508	H(102) C(69) N(...)		Any
<input checked="" type="checkbox"/>	EEEDVIEVYQEQTGG	EEEDVIEVYQEQTGG	1705.73189	1706.7153	H(107) C(72) N(...)		Any
<input checked="" type="checkbox"/>	Formyl	Formyl	27.994915	28.0101	C O		Protei
<input checked="" type="checkbox"/>	Formyl	Formyl	27.994915	28.0101	C O		Any
<input checked="" type="checkbox"/>	Formyl	Formyl	27.994915	28.0101	C O		Any_↑
<input checked="" type="checkbox"/>	Gln->pyro-Glu	Gln->pyro-Glu	-17.026549	-17.0305	H(-3) N(-1)		Any_↑
<input checked="" type="checkbox"/>	Hex(5)HexNac(4)Ne...	Hex(5)HexNac(4)Ne...	1913.677025	1914.7277	Hex(5) HexNac(...)		Any
<input checked="" type="checkbox"/>	Hex(5)HexNac(4)Ne...	Hex(5)HexNac(4)Ne...	2204.772441	2205.9822	Hex(5) HexNac(...)		Any

- Click + to the left of each modification row to see the amino acids that the modification is found on, the letter abbreviation of this amino acid, and the modification type or category. Figure 159 shows an example of the information given for the Acetyl modification. Table 10 lists the available modification categories.

Figure 159. Modification information for Acetyl displayed in columns




Table 10. Available modification categories (Sheet 1 of 2)

Classification	Description
Post-translational	Protein modification after translation (<i>in vivo</i>)
Co-translational	Amino acid modified in translation (for example, myristyl glycine)
Pre-translational	Amino acid modified before integration into a protein (for example, formyl methionine)
Chemical derivative	Chemically induced modification (for example, during sample preparation)
Artifact	Modification made during sample preparation
N-linked glycosylation	Glycosylation (<i>in vivo</i>)

Table 10. Available modification categories (Sheet 2 of 2)


Classification	Description
O-linked glycosylation	Glycosylation (in vivo)
Other glycosylation	Glycosylation (in vivo)
Synthetic peptide protection group	Protection group used in chemical peptide synthesis (for example, trityl (triphenylmethyl))
Isotopic label	Label for quantification
Non-standard residue	Amino acid derivative like selenomethionine
Multiple	More than one classification possible
AA substitution	Amino acid replaced by another amino acid (mutation)
Other	Modification not fitting into another category

The Proteome Discoverer application automatically imports the classifications from unimod.org, the protein modifications online database for mass spectrometry applications. You can also manually define your own classifications.

Adding Chemical Modifications

You can create new chemical modifications and add them to the Chemical Modifications view. For example, you might have a new or experimental label that you want to add to the list of chemical modifications.

❖ To add a new chemical modification

1. Choose **Administration > Maintain Chemical Modifications**, or click the **Maintain Chemical Modifications** icon, , either on the toolbar or on the Administration page.


The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

2. Click the **Add a Modification** heading.

An empty row appears, as shown in [Figure 160](#).

Figure 160. Adding a row in the Chemical Modifications view

Is Active	Modification	Abbreviation	Delta Mass	Delta Average Mass	Substitution	Leaving Group	Position	Unimod Accession No.
<input checked="" type="checkbox"/>			1	1			Any	0
<input checked="" type="checkbox"/>	Acetyl	Acetyl	42.010565	42.0367	H(2) C(2) O		Any	1
<input checked="" type="checkbox"/>	Acetyl	Acetyl	42.010565	42.0367	H(2) C(2) O		Protein_N...	1
<input checked="" type="checkbox"/>	Acetyl	Acetyl	42.010565	42.0367	H(2) C(2) O		Any_N_Te...	1
<input checked="" type="checkbox"/>	Amidated	Amidated	-0.984016	-0.9848	H N O(-1)		Protein_C...	2
<input checked="" type="checkbox"/>	Amidated	Amidated	-0.984016	-0.9848	H N O(-1)		Any_C_Te...	2
<input checked="" type="checkbox"/>	Carbamidomethyl	Carbamidomethyl	57.021464	57.0513	H(3) C(2) N O		Any	4
<input checked="" type="checkbox"/>	Carbamidomethyl	Carbamidomethyl	57.021464	57.0513	H(3) C(2) N O		Any_N_Te...	4
<input checked="" type="checkbox"/>	Carbamyl	Carbamyl	43.005814	43.0247	H C N O		Any	5
<input checked="" type="checkbox"/>	Carbamyl	Carbamyl	43.005814	43.0247	H C N O		Any_N_Te...	5

3. In the empty row, enter the name of the modification, the delta masses, the chemical substitution, the chemical group that is leaving, the position, and the abbreviations of the modifications.
4. If you select Any in the Position column, a message box opens to inform you that you must specify which amino acids (target amino acids) will possibly have the modification. For instructions on this procedure, see [“Adding Amino Acids.”](#)
5. To accept the new modifications, click the **Apply** icon, .
6. Add an amino acid to the modifications. See [“Adding Amino Acids.”](#)

❖ **To update an existing chemical modification**

1. Choose **Administration > Maintain Chemical Modifications**.

The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

2. In the Modification column, click the cell that you want to update.
3. Type your changes for the delta masses, the substitution, the group that it is leaving, the position, or the abbreviations of the modifications.

For chemical modifications that you add yourself, you can edit any column except the Unimod Accession No. column. The Unimod Accession No. column identifies these modifications by a zero. For chemical modifications that you import from UNIMOD, you can edit only the Modification and Abbreviation columns. UNIMOD chemical modifications are identified by a number greater than zero in the Unimod Accession No. column.

Columns that you can edit activate an edit button when you click them. Columns that you cannot edit display a gray background.

4. To accept the changes, click the **Apply** icon, .

Adding Amino Acids

You can add amino acids to a modification that has been set up for any position.

❖ **To add an amino acid to a modification**

1. Choose **Administration > Maintain Chemical Modifications**.

The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

2. Click + to the left of the modification row that you want to update.

The row must display Any in the Position column.

The list of classifications now appears, as shown in [Figure 159](#) on [page 230](#).

3. Click the **Add a Modification** line below the list of amino acids.


Figure 161 shows this line.

Figure 161. Adding an amino acid to a modification

Amino Acid Name	One Letter Code	Classification
Cysteine	C	Chemical derivative
Aspartic Acid	D	Artefact
Glutamic Acid	E	Artefact
Histidine	H	Artefact
Lysine	K	Artefact

* Add a Modification...

An empty row appears.

- In the empty row, select the amino acid from the list in the Amino Acid Name column.
The amino acid and the one-letter abbreviation appear.
- From the list in the Classification column, select the type of modification.
- To save the modifications, click the **Apply** icon,  **Apply**.

When you reimport data from unimod.org, the Proteome Discoverer application retains the modification that you added. However, if you want to change the classification of an amino acid, you must do so before reimporting the Unimod data. After you import the Unimod data, the only way to change the classification is to delete the amino acid and re-add it with another classification.

Deleting Amino Acids

You can also delete amino acids from chemical modifications.

❖ To delete an amino acid from a chemical modification

- Choose **Administration > Maintain Chemical Modifications**.

The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

- Click **+** to the left of the modification row that you want to delete.

The row expands and the associated amino acids appear.

- Select the amino acid row that you want to delete.

- Click the **Remove** icon,  **Remove**.

- In the Delete Row dialog box, click **Yes**.

The row is removed from the chemical modifications table.


Deleting Chemical Modifications

You can remove chemical modifications from the Chemical Modifications view.

❖ To delete a modification

1. Choose **Administration > Maintain Chemical Modifications**.

The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

2. Select the row of the modification that you want to delete.
3. Click the **Remove** icon,  **Remove**.
4. In the Delete Row dialog box, click **Yes**.

The row is removed from the chemical modifications table.

Importing Chemical Modifications


You can import chemical modifications from a local file or obtain an updated version from unimod.org, a public domain database.

When you install the Proteome Discoverer application, it automatically imports accessions from unimod.org as chemical modifications.

❖ To import chemical modifications from a local file

1. Choose **Administration > Maintain Chemical Modifications**.

The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

2. Click the **Import** icon, .
3. In the Import From list of the Import Modifications dialog box, select **Local File**.
4. In the adjacent box, click the **Browse** button (...) to browse to your file, or type the name and path of the file in the box.
5. To overwrite an existing upload, select the **Overwrite Existing** check box.
6. Click **Import**.

A status message appears.

7. When the upload is complete, click **Close**.

❖ To import chemical modifications from unimod.org


1. Choose **Administration > Maintain Chemical Modifications**.

The Chemical Modifications view appears, as shown in [Figure 158](#) on [page 229](#).

2. Click the **Import** icon, .

The Import Modifications dialog box appears, as shown in [Figure 162](#).

Figure 162. Import Modifications dialog box



3. In the Import From list, select **Unimod**.

The UNIMOD URL appears in the adjacent box.

4. To overwrite an existing upload, select the **Overwrite Existing** check box.
5. Click **Import**.

A status message appears.

6. When the upload is complete, click **Close**.

For chemical modifications imported from unimod.org, you can only edit the Is Active, Modification, and Abbreviation columns. You do not have access to the Delta Mass, Delta Average Mass, Substitution, Leaving Group, Position, and UNIMOD Accession No. columns. Chemical modifications imported from unimod.org have a number greater than zero in the Unimod Access No. column.

If you select the Overwrite Existing check box, the Proteome Discoverer application does the following when it imports chemical modifications from unimod.org:

- Updates the columns that are inaccessible to you.
- Updates the names and the abbreviations of the modifications.
- Adds any new amino acids found in unimod.org.
- Adds any amino acids that you removed if they are defined in unimod.org.
- Removes any amino acids that you added if they are defined in unimod.org.

If you do not select the Overwrite Existing check box, the Proteome Discoverer application performs the same tasks as it does during installation:

- Updates the columns that are inaccessible to you.
- Leaves the modification name and abbreviation unchanged.
- Adds any new amino acids found in unimod.org.
- Adds any amino acids that you removed if they are defined in unimod.org.
- Leaves unchanged any amino acids that you added.

Using the Qual Browser Application

The Proteome Discoverer application provides a link to the Qual Browser application from selected PSMs if the application was already installed on the computer. With the Qual Browser application, you can examine spectra and chromatograms in detail, view the entire ion chromatogram, and browse individual precursor and MSⁿ data. You can filter the results in a variety of ways, for example, to produce a selected ion chromatogram. The Qual Browser application automatically displays the elemental composition, theoretical mass, delta values, and ring and double-bond (RDB) equivalents for your high-resolution data. (RDB equivalents measure the number of unsaturated bonds in a compound and limit the calculated formulas to only those that make sense chemically.)

You must have the Xcalibur data system installed to use the Qual Browser application. For information about using the Qual Browser application, refer to the *Thermo Xcalibur Qual Browser User Guide*.

You must also have a search results file open and a specific peptide or search input row selected before the Qual Browser application becomes available. If you are viewing the Administration page, the Qual Browser application does not open a raw data file.


❖ To open the Qual Browser application

1. Open a .pdResults file.
2. Choose **Tools > Open QualBrowser**, choose **Tools > Open QualBrowser** or press CTRL+SHIFT+B to open the Spectrum window.

Note You must have a search results (.pdResult) file open and selected before the Open QualBrowser command becomes available on the Tools menu. In addition, the Open QualBrowser command is available only when you select at least a single peptide or a search input item first. You cannot use QualBrowser if the original raw data file or files are missing. The .pdResult file and the raw data file must reside in the same directory.

The Qual Browser application opens, as shown in [Figure 163](#).

Figure 163. Qual Browser application window



3. Right-click and choose **Display Options** from the shortcut menu.
4. To automatically annotate your peaks with the elemental composition, theoretical mass, RDB equivalent, or mass delta, click the **Composition** tab and select the labels for display.

Customizing Cleavage Reagents

In the Cleavage Reagents view, you can explore the default types of reagents and their corresponding settings. You can also add, remove, and modify the reagents and their corresponding settings. The Cleavage Reagents view contains the cleavage sites, cleavage inhibitors, abbreviations, and cleavage specificities.

For detailed information about enzyme cleavage, see “[Enzyme Cleavage Properties](#)” on [page 496](#).

- [Opening the Cleavage Reagents View](#)
- [Adding a Cleavage Reagent](#)
- [Deleting a Cleavage Reagent](#)
- [Modifying a Cleavage Reagent](#)
- [Filtering Cleavage Reagent Data](#)

Opening the Cleavage Reagents View

❖ **To display the Cleavage Reagents view**

Choose **Administration > Maintain Cleavage Reagents**.

The Cleavage Reagents view appears, as shown in [Figure 164](#).

Figure 164. Cleavage Reagents view

Name	Cleavage Sites	Cleavage Inhibitors	Offset	Abbreviation	Cleavage Specificities
* Click here to add a new record...					
Trypsin_R	R	P	1	-	Full; Semi
Trypsin_K	K	P	1	-	Full; Semi
Trypsin(KRLNH)	HKLNR	-	1	Try_a	Full; Semi
Trypsin	KR	P	1	Try	Full; Semi
Staph_Protease	E	-	1	-	Full
Proline_Endopept	P	-	1	-	Full
No-Enzyme	-	-	0	-	Unspecific
No-Cleavage	-	-	0	-	Semi; Semi (N-Term); Semi (...)
LysC	K	-	1	-	Full
IodosoBenzoate	W	-	1	-	Full
GluC	DE	-	1	-	Full
Elastase/Tryp/Ch...	AFIKLRVWY	P	1	-	Full
Elastase	AILV	P	1	-	Full
Cyanogen_Bromi...	M	-	1	-	Full
Clostripain	R	-	1	-	Full
Chymotrypsin(F...	FWY	P	1	ChTr	Full
Chymotrypsin	FLWY	-	1	ChyTr	Full
AspN	D	-	0	-	Full

Cleavage Enzyme Specificities

The Proteome Discoverer application supports the following cleavage specificities for enzymes:

- Full: Specifies a full enzymatic digestion.
- Semi: Specifies a semi-enzymatic digestion. The cleavage is specific at one terminal end.
- Semi (N-Term): Specifies that the cleavage be specific at the amino terminal end and unspecific at the carboxyl terminal end.
- Semi (C-Term): Specifies that the cleavage be specific at the carboxyl terminal end and unspecific at the amino terminal end.
- Enzymatic Unspecific: Specifies an unspecific digestion with an extra weight for enzymatically digested peptides.

The list of cleavage reagents includes the No-Enzyme and No-Cleavage reagents:

- The No-Enzyme reagent digests a sequence at every amino acid. For example, if the FASTA file contains the ACDEFGHIK sequence, No-Enzyme cleaves it into peptides A, AC, ACD, ..., C, CD, CDE, ..., D, DE, DEF.

The application supports the following specificity for the No-Enzyme reagent:

- Unspecific: Specifies a non-specific digestion.
- The No-Cleavage reagent does not cleave a sequence. For example, if the FASTA file contains the ACDEFGHIK sequence, the full sequence is the only possible sequence.

The application supports the following cleavage specificities for the No-Cleavage reagent:

- No-Cleavage: Specifies that no cleavages occur, so intact proteins result.
- Semi, Semi-N, Semi-C: Specifies terminal signal peptides that are cleaved off in vivo.

Cleavage Reagents View Parameters

[Table 11](#) describes the parameters on the Cleavage Reagents view.

Table 11. Cleavage Reagents view parameters (Sheet 1 of 4)

Parameter	Description
Remove	Deletes the highlighted cleavage reagent.
Apply	Saves any changes made to the cleavage reagents.
Name	Specifies the name of the reagent used for the protein digestion. <ul style="list-style-type: none"> • All: Returns the filtered search results to the results that were first loaded. • Custom: Opens the Custom Filter dialog box, shown in Figure 165 on page 244. • Blanks: Filters out rows that have data-filled cells in the column whose funnel icon you clicked. • NonBlanks: Filters out rows that have empty cells in the column whose funnel icon you clicked. • List of enzymes: Lists the specific names of all available enzymes that act as reagents.

Table 11. Cleavage Reagents view parameters (Sheet 2 of 4)

Parameter	Description
Cleavage Sites	<p>Specifies the position (amino acid) at which to cleave the sequence.</p> <ul style="list-style-type: none"> • All: Returns the filtered search results to the results that were first loaded. • Custom: Opens the Custom Filter dialog box, shown in Figure 165 on page 244. • Blanks: Filters out rows that have data-filled cells in the column whose funnel icon you clicked. • NonBlanks: Filters out rows that have empty cells in the column whose funnel icon you clicked. • List of amino acids: Lists the names of the amino acids at which the enzyme cleaves the peptide.
Cleavage Inhibitors	<p>Specifies the amino acids that block cleavage when adjacent to the cleavage site.</p> <ul style="list-style-type: none"> • All: Returns the filtered search results to the results that were first loaded. • Custom: Opens the Custom Filter dialog box, shown in Figure 165 on page 244. • Blanks: Filters out rows that have data-filled cells in the column whose Filter icon you clicked. • NonBlanks: Filters out rows that have empty cells in the column whose Filter icon you clicked. • -: No amino acids inhibit cleavage. • P: Adjacent proline amino acids inhibit cleavage.

Table 11. Cleavage Reagents view parameters (Sheet 3 of 4)


Parameter	Description
Offset	<p>Specifies at which side of the given amino acids the protease should cleave:</p> <ul style="list-style-type: none">• All: Returns the filtered search results to the results that were first loaded.• Custom: Opens the Custom Filter dialog box, shown in Figure 165 on page 244.• Blanks: Filters out rows that have data-filled cells in the column whose Filter icon you clicked.• NonBlanks: Filters out rows that have empty cells in the column whose Filter icon you clicked.• 0: The protease cleaves at the N-terminal side of the specified cleavage sites.• 1: The protease cleaves at the C-terminal side of the specified cleavage sites.
Abbreviation	<p>Specifies the user-defined abbreviation.</p> <ul style="list-style-type: none">• All: Returns the filtered search results to the results that were first loaded.• Custom: Opens the Custom Filter dialog box, shown in Figure 165 on page 244.• Blanks: Filters out rows that have data-filled cells in the column whose Filter icon you clicked.• NonBlanks: Filters out rows that have empty cells in the column whose Filter icon you clicked.• List of abbreviated enzymes: Lists the existing user-defined abbreviations.

Table 11. Cleavage Reagents view parameters (Sheet 4 of 4)

Parameter	Description
Cleavage Specificities	<p>Selects the cleavage specificity for enzymes. You can select more than one cleavage specificity. The available cleavage specificities are the following:</p> <ul style="list-style-type: none"> • Full: Specifies a full enzymatic digestion. Sequences are not cleaved. • Semi: Specifies a semi-enzymatic digestion. The cleavage is specific at one terminal end. • Semi (N-Term): Specifies that the cleavage be specific at the amino terminal end and unspecific at the carboxyl terminal end. • Semi (C-Term): Specifies that the cleavage be specific at the carboxyl terminal end and unspecific at the amino terminal end. • Enzymatic Unspecific: Specifies an unspecific digestion with an extra weight for enzymatically digested peptides.
Click here to add a new record	Opens a row in the reagents table so you can enter the name and information about a new cleavage reagent.

Adding a Cleavage Reagent

❖ To add a new cleavage reagent

1. Click the Name column cell and click **Click Here To Add a New Record**.
2. Modify the default values in the row of that new reagent.
3. Click the **Apply** icon,  **Apply** .


Deleting a Cleavage Reagent

❖ To delete a cleavage reagent

1. Click the box in the * column next to the row that you want to delete.
2. Click **Delete**.
3. Click **Yes** in the confirmation box that appears.

Modifying a Cleavage Reagent

❖ To modify a cleavage reagent

1. Click in the column for the reagent you want to modify, select the current contents, and enter the new information.
2. Click the **Apply** icon,  **Apply** .

Filtering Cleavage Reagent Data

❖ To filter cleavage reagent data



1. Click the **Filter** icon,  , next to the header of the column.
2. Select one of the following:
 - All: Returns the filtered search results to the results that were first loaded.
 - Custom: Opens the Custom Filter dialog box, shown in [Figure 165](#).

Figure 165. Custom Filter dialog box



For information about using this type of dialog box, see “[Filtering with Display Filters](#)” on [page 254](#).

- Blanks: Filters out rows that have data-filled cells in the column whose Filter icon you clicked.
- NonBlanks: Filters out rows that have empty cells in the column whose Filter icon you clicked.

Filtering Data

A single or multiconsensus .pdResult file displays a list of PSMs and proteins identified by the search engine that you specify. This chapter explains how to sort and filter the data in your Proteome Discoverer results report.

Contents

- [Filtering Results](#)
- [Filtering PSMs](#)
- [Filtering Proteins](#)
- [Filtering Phosphorylation Site Probabilities](#)
- [Filtering with Display Filters](#)

Filtering Results

You can filter the results of an analysis during processing or afterward in the opened results report.

- In the processing workflow, you can use the Maximum Delta Cn parameter of the PSM validation nodes (Percolator, Target Decoy PSM Validator, Fixed Value PSM Validator) to reduce the number of PSMs stored and to affect the PSM FDR estimation. For information on these nodes, refer to the Help.
- In the consensus workflow, you can use the filtering parameters of the MSF Files node, the Peptide and Protein Filter node, and the Protein Grouping node to determine the set of protein groups, proteins, peptide groups, and PSMs displayed in a .pdresult file. For information on these nodes, refer to the Help.
- In the .pdResult file, you can apply an extensive set of editable filters, called display filters, to any of the open pages to simplify and streamline the data to review. You can use display filters to show or hide data with particular properties. Applying display filters does not affect calculations performed by workflows, such as confidence levels, coverage, or rank.

Filtering PSMs

Search engines often provide multiple possible matching peptides as explanations for the same spectrum. Most of the time you can clearly distinguish the top-scoring match from the other PSMs, but sometimes, especially in the presence of dynamic modifications, the best-scoring matches of the same spectrum have very similar scores. In this case, you can filter the results to select the best-scoring PSMs and the matches that have very similar scores by using the Maximum Delta Cn parameter or the Maximum Rank parameter of the following nodes:

- Processing workflow:
 - [Filtering PSMs with the Percolator Node](#)
 - [Filtering PSMs with the Fixed Value PSM Validator Node](#)
 - [Filtering PSMs with the Target Decoy PSM Validator Node](#)
- Consensus workflow:
 - [Filtering PSMs with the MSF Files Node](#)

The ΔCn value displays the normalized score difference between the currently selected PSM and the highest-scoring PSM for that spectrum:

$$\Delta Cn(\text{rank } i) = \frac{\text{score}_{\text{rank } 1} - \text{score}_{\text{rank } i}}{\text{score}_{\text{rank } 1}}$$

The Maximum Delta Cn parameter filters out all PSMs with a ΔCn score larger than the specified value.

On the PSMs page of the .pdResult file, the ΔCn column displays the ΔCn values. For example, [Figure 166](#) shows how the score of a peptide ranked 2 compares to other high-confidence peptides from the same spectrum.

Figure 166. Δ Cn scores for multiple high-confidence peptides from the same spectrum

The screenshot displays the Proteome Discoverer interface. The top panel shows a list of PSMs with columns for Checked, Best PSM Ambiguity, Mass Analyzer, Activation Type, MS Order, # PS, Isolation Interference [%], Ion Inject Time [ms], Precursor m/z [Da], Precursor MH+ [Da], Precursor Charge, RT [min], First Scan, and Spectrum File. The bottom panel shows a detailed view of the top 10 PSMs, including identifying nodes, peptide sequences, modifications, protein group counts, master protein accessions, protein accessions, missed cleavages, charge, Δ Score, Δ Cn, Rank, and Search Engine Rank.


#	Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions	# Missed Cleavages	Charge	Δ Score	Δ Cn	Rank	Search Engine Rank
1	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	S18(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0032	0.0000	1	1
2	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	T17(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0032	0.0032	2	2
3	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	Y23(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0127	0.0127	3	3
4	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	S22(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0127	0.0127	3	3
5	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	S10(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0223	0.0223	5	5
6	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	Y13(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0254	0.0254	6	6
7	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	S11(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0350	0.0350	7	7
8	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTsGYAGLLSSAYGGLTSP	S3(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0413	0.0413	8	8
9	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	TTSGYAGLLSSAYGGLTSP	T2(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0429	0.0429	9	9
10	<input checked="" type="checkbox"/>	●	Sequest HT (A19)	Unambiguous	ITSGYAGLLSSAYGGLTSP	T1(Phospho)	1	2	IP101021414.1	IP101021414.1; IP100554648.	0	3	0.0429	0.0429	9	9

Filtering PSMs with the Perculator Node

Figure 167 shows the parameters available for the Perculator node (for information about these parameters, refer to the Help). Setting its Maximum Delta Cn parameter to 0 ensures that the application consider only rank 1 PSM matches for further processing with Perculator. When you set this parameter to a slightly higher value, the application might consider PSMs with a rank higher than 1.

A slightly higher value can be a good option in case of phosphopeptide analysis, where the first two or more matches often refer to the same peptide sequence and differ solely in the potential location of the phosphorylation site. Their scores are usually very close, if not identical, so they generate very low Delta Cn values. The default setting of the Maximum Delta Cn parameter, 0.05, is intended to accommodate such cases. You can specify a maximum value of 0.1 for the Maximum Delta parameter. Allowing PSMs with even higher Delta Cn values negatively affects Perculator's performance. In addition, it hugely increases the volume of data to be processed.

Figure 167. Percolator node parameters



Filtering PSMs with the Fixed Value PSM Validator Node

For information on filtering with the Fixed Value PSM Validator node in the processing workflow, refer to the Help.

Filtering PSMs with the Target Decoy PSM Validator Node

For information on filtering with the Target Decoy PSM Validator node in the processing workflow, refer to the Help.

Filtering PSMs with the MSF Files Node

You can also use four parameters of the MSF Files node in the consensus workflow to filter PSMs: the Maximum Delta Cn, the Maximum Rank, the Maximum Delta Mass, and the Score/Threshold parameters of the individual search node.

❖ **To filter PSMs with the Maximum Delta Cn parameter**

1. Create a consensus workflow that contains the MSF Files node.

For information about creating a basic consensus workflow, see “[Creating a Consensus Workflow](#)” on [page 112](#).

2. Set any parameters for the MSF Files node.

For information on these parameters, refer to the Help.


3. Click **Show Advanced Parameters** in the Parameters pane of the Workflow Editor, and specify a value for the Maximum Delta Cn parameter to set a threshold that determines which PSMs the application transfers to the report file.

[Figure 168](#) shows this parameter.

The threshold can be within a range of 0 to 0.1. The default is 0.05.

The application transfers only those PSMs with a ΔCn value greater than the specified threshold to the report file.

Figure 168. MFS Files node parameters



❖ **To filter PSMs with the Maximum Rank parameter**

1. Create a consensus workflow that contains the MSF Files node.
For information, see [“Creating a Consensus Workflow”](#) on page 112.
2. Set any parameters for the MSF Files node.
For information on these parameters, refer to the Help.
3. Click **Show Advanced Parameters** in the Parameters pane of the Workflow Editor.
4. Specify a value for the Maximum Rank parameter to set a threshold that determines which PSMs the application transfers to the report file.

Figure 168 shows this parameter.

You can set the threshold to any positive value. The default is 0, which means that the Maximum Rank filter is not applied.

The application transfers to the report file only those PSMs with a rank better than or equal to the specified threshold.

❖ **To filter PSMs with the Maximum Delta Mass parameter**

1. Create a consensus workflow that contains the MSF Files node.

For information about creating a basic consensus workflow, see [“Creating a Consensus Workflow”](#) on page 112.

2. Set any parameters for the MSF Files node.

For information on these parameters, refer to the Help.

3. Specify a value for the Maximum Delta Mass parameter to indicate whether the application should apply a delta mass filter and exclude from the final result peptide matches with a larger mass difference between the theoretical and the found peptide m/z .

Use a value within the range of 0.0–5.0 Da or 0.0–5000 ppm. The default is 0 ppm. If you set this parameter to zero, the application does not apply a delta mass filter.

❖ **To filter PSMs with the Score and Threshold parameters**

1. Create a consensus workflow that contains the MSF Files node.

For information about creating a basic consensus workflow, see [“Creating a Consensus Workflow”](#) on page 112.


2. Set any parameters for the MSF Files node.

For information about these parameters, refer to the Help.

3. Set the filter parameters for the node:

- a. In the 1. Score box under PSM Filters, select the name of the score to filter by from the list, as shown in [Figure 169](#).

Figure 169. Selecting the score to filter by



If you do not see the name of the filter that you want to filter by, see [“To specify the list of score names to filter by”](#) on [page 253](#) for instructions on adding score names to the list.

- b. In the 1. Threshold box, type the score to filter by.
- c. Repeat the previous two steps for any additional filters.

You can specify up to 10 PSM scores and thresholds to filter by.


For each score type, the application excludes scores lower than the specified threshold.

❖ **To specify the list of score names to filter by**

1. Choose **Administration > Configuration > MSF Files**.
2. Click the **Browse** button (...) next to the PSM Scores parameter.

The Edit Parameter Text for PSM Scores dialog box appears, as shown in [Figure 170](#).

Figure 170. Edit Parameter Text for PSM Scores dialog box



3. Beneath the last line, type the name of the score type that you want to add in the format *search_node_name:score_name*.

You can find the search node name in the Workflow Editor or in the Result Summaries when you open the .pdResult file. The column on the PSMs page that shows the score also displays the search node name.

The specified string for the score type name must be in the format *search_node_name:score_name*. If a string does not match this format, it is not available in the score parameters of the filter node. If the search node name and score names are not correct, the filter node cannot filter the PSMs on the basis of that score. If the node cannot find a specified score, it issues a warning message.

4. Click **OK**.
5. Close the study and reopen it.

The added score type is now available in the list in the Edit Parameter Text dialog box.

Filtering PSMs with the Peptide and Protein Filter Node

You can also filter PSMs by using the Peptide and Protein Filter node in the consensus workflow. This node filters by score for the type of score specified. For information on this node, refer to the Help.

Filtering Proteins


You can filter proteins with the Peptide and Protein Filter node. For information on this node, refer to the Help.

When you place the Display Filter node in the consensus workflow, the Proteins page displays only master proteins by default. If you want to filter the proteins on the Proteins page with a filter set other than the default, follow the procedure in [“To apply stored filters to a .pdResult file”](#) on [page 268](#).

Filtering Phosphorylation Site Probabilities

Filtering phosphorylation site probabilities involves both display filters and the ptmRS node in the processing workflow. For information on filtering phosphorylation modification data, refer to [“Filtering Phosphorylation Site Probabilities”](#) on [page 374](#).

Filtering with Display Filters

You can use the View > Display Filter command or the Display Filter icon, , to display only those items that meet a particular condition that you define. Refining your search results in this way reduces the complexity of the displayed results, allowing for a quicker and more productive analysis. You can sort and filter your results by a number of criteria, such as charge state, modifications, or even peptide probability. You can also create and apply multiple filters to your search results. You can save such filter sets to disk and load them in into the same or other report files.


- [Filtering the Data on a Results Report Page](#)
- [Example 1](#)
- [Example 2](#)
- [Example 3](#)
- [Example 4](#)
- [Filtering with Filter Sets](#)
- [Using Display Filters to Filter Numerical Values by a Specified Precision Value](#)
- [Finding Common or Unique Proteins in Multiple Searches](#)

- Applying Filters Specific to Different Searches in Multiconsensus Reports
-

Filtering the Data on a Results Report Page


You can use the application's display filters to refine the results in the .pdResult file.

❖ To filter the data on a .pdResult data page

1. Open a .pdResult file, and click on the tab of the page that you want to view.
2. Choose **View > Display Filter**, or click the **Display Filter** icon, .

The dockable Display Filter pane opens at the top of the results file, as shown in Figure 171.

Figure 171. Display Filters pane



The screenshot shows the Proteome Discoverer software interface. At the top, a dockable 'Display Filter' pane is open, containing a toolbar with 'Load', 'Save', 'Clear', and 'Apply' buttons. Below the toolbar are radio buttons for 'Protein Groups', 'Proteins', 'Peptide Groups', 'PSMs', 'MS/MS Spectrum Info', and 'Quan Spectra'. The 'PSMs' section is active, showing an 'AND' filter configuration with 'Add group' and 'Add property' buttons.

Below the filter pane is a table with the following columns: Checked, Confidence, Identifying Node, PSM Ambiguity, Annotated Sequence, Modifications, Contaminant, Standard, Decoy, # Protein Groups, # Proteins, Master Protein Accessions, and Protein Accessions. The table contains 28 rows of data, each representing a search result. The status bar at the bottom indicates: 'Ready 792 Protein Groups; 2170 Proteins; 2911 Peptide Groups; 9775 PSMs; 72125 MS/MS Spectrum Info; 72155 Quan Spectra'.

Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	Contaminant	Standard	Decoy	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	IVRPEVDVMcTAFHDNEETF	N-Term(TMT6plex); C10(Car	X	X		1	1	P02768	P02768
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSSYGSGG	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSSYGSGG	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	IVRPEVDVMcTAFHDNEETF	N-Term(TMT6plex); C10(Car	X	X		1	1	P02768	P02768
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSSYGSGG	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSSYGSGG	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	nQEAGHEAAGEEAAEASGF	N-Term(TMT6plex); K22(TM				1	1	ENSMUSP00000109741	ENSMUSP000001097
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	nQEAGHEAAGEEAAEASGF	N-Term(TMT6plex); K22(TM				1	1	ENSMUSP00000109741	ENSMUSP000001097
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gLADGVDkDAVNAAIQAI	N-Term(TMT6plex); K9(TMT				1	1	ENSMUSP00000059678	ENSMUSP000000596
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVLPGMVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP000000424
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Selected	TKEQANAVSEAVSSVNTV	K22(TMT6plex)		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Selected	TKEQANAVSEAVSSVNTV	K22(TMT6plex)		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Rejected	eQANAVSEAVSSVNTVATI	N-Term(TMT6plex); K20(TM		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVLPGMVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP000000424
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVLPGMVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP000000424
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Rejected	eQANAVSEAVSSVNTVATI	N-Term(TMT6plex); K20(TM		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gLADGVDkDAVNAAIQAI	N-Term(TMT6plex); K9(TMT				1	1	ENSMUSP00000059678	ENSMUSP000000596
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	sGGGGGGGSGGGGVS	N-Term(TMT6plex); C10(Car	X	X		1	1	P13645	P13645
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	eGTGSTATSSGAGGAVG	N-Term(TMT6plex); K19(TM				1	1	ENSMUSP00000022925	ENSMUSP000000229
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Selected	TKEQANAVSEAVSSVNTV	K22(TMT6plex)		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	tKEQANAVSEAVSSVNTVA	N-Term(TMT6plex); K2(TMT		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVLPGMVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP000000424
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vRPeVYYGGAEIGQQIR	N-Term(TMT6plex); C4(Car				1	1	ENSMUSP0000000804	ENSMUSP000000080
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	sTSSSSSSSNQQTKE	N-Term(TMT6plex); K17(TM	X	X		2	2	P51965	P51965; ENSMUSP00
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	ISScGGGGGSGGGGFGS	N-Term(TMT6plex); C4(Car	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	tKEQANAVSEAVSSVNTVA	N-Term(TMT6plex); K2(TMT		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVLPGMVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP000000424
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	sGGGGGGGSGGGGVS	N-Term(TMT6plex); C10(Car	X	X		1	1	P13645	P13645


The left pane of the filter area lists the names of the filter groups. These names correspond to the names of the pages in the results report, such as Protein Groups, Proteins, PSMs, Decoy Peptide Groups, Quan Spectra, and so forth.

For convenience, the application automatically highlights a filter group when you navigate through different pages in the results report.

3. Select the filter category, for example, **PSMs**.
4. Set the first condition:
 - a. Click the Add Property box in the right pane, which has a blue background, and select the property that you want from the menu.

An operator box and a values box appear to the right of the element, as shown in [Figure 172](#). The Apply button is now highlighted.

Figure 172. Elements of the Display Filter pane



- b. Click the operator box, and choose an operator for the condition from the menu.
For numerical values, the default is Is Equal To. For textual values, the default is Contains.
 - c. In the values box, type a value.

You can use simple wildcards to filter for specific patterns in strings. Choose the Contains Pattern operator and type a question mark (?) in the values box for any single character or an asterisk (*) for a sequence of arbitrary characters of any length.


5. (Optional) Set the additional conditions, if applicable, by clicking in the Add Property box and repeating the previous steps.

Note Multiple conditions are arranged as siblings of an And box, which means that both conditions must apply to display a particular item. Change this box to Or if only one of the conditions applies.

- (Optional) To add another filter, click **Add Property**.

Another level of filtering appears above the original level, and a new Add Property box appears under the original, as shown in [Figure 173](#).


Figure 173. Display Filter pane with two filters




- From the menu to the right of the selected filter group, choose the level that you want the filter to apply to, as shown in [Figure 174](#):

- **Affects First Level Only:** Applies the filter to the items on a page, but not to the items on any associated pages.
- **Affects All Levels:** Applies the filter to both the items on a page and to the items on any associated pages at all levels.


Figure 174. Selecting the level of items to apply the filter to




- To apply the filters in the filter set to the results report, click the **Apply** icon, .

For examples showing how to build increasingly complex filters, see “[Example 1](#)” on page 261, “[Example 2](#)” on page 262, and “[Example 3](#)” on page 264.

❖ **To temporarily disable filters**

- In the left pane of the Display Filter pane, click the **On/Off** icon, , next to the filter category that you want to disable temporarily.

The icon changes to an Off icon.

- Click the **Apply** icon, .


5 Filtering Data

Filtering with Display Filters



Using the example in [Figure 174](#), the PSM filter is temporarily disabled, and [Figure 175](#) shows the results. The filters for Protein Groups, Proteins, Peptide Groups, MS/MS Spectrum Info, and Quan Spectra are still enabled.

3. To re-enable the filter categories, click the **On/Off** icon and click the **Apply** icon.



Figure 175. Temporarily disabled Protein Groups and PSMs filters




❖ To clear conditions and filters

- To remove a condition, click **Remove** next to it.
- To clear the filter on the currently selected page, click the **Clear** icon,  **Clear**, and then click the **Apply** icon,  **Apply**.

For example, if the PSMs page is currently selected, the Clear icon clears only the filter on the PSMs page.

- To clear the filters on all pages of the open results report, click the **Clear All** icon,  **Clear All**, and then click the **Apply** icon,  **Apply**.

❖ To save a filter set

1. Click the **Save** icon,  **Save**, to save the filter set in a file on the disk.
2. In the Export Filter dialog box, browse to the file to save the filter set in, or type the name in the File Name box.
3. Click **Save**.

The file has a .filterset suffix.

When you reopen the .pdResult file containing the saved filter set, the “filtered” notation appears in parentheses on the tabs of the pages that have filters set so that you know immediately that the results on that page are filtered. In [Figure 176](#), the results on the PSMs page are filtered.

Figure 176. .pdResult file containing filtered results for the PSMs page

PSMs tab indicating that the page contains filtered results

	Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	Contaminant	Standard	Decoy	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions
1	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gEGQLSAAER	N-Term(TMT6plex)				1	1	ENSMUSP00000046101	ENSMUSP0000000461
2	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	aQASAPAQAPk	N-Term(TMT6plex); K11(TMT				1	4	ENSMUSP00000080203	ENSMUSP000001297
3	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	ISSPATLNSR	N-Term(TMT6plex)	X	X		1	1	P00761	P00761
4	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gEGQLSAAER	N-Term(TMT6plex)				1	1	ENSMUSP00000046101	ENSMUSP0000000461
5	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	dkEAcVhk	N-Term(TMT6plex); K2(TMT	X	X		1	1	P02787	P02787
6	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	tVSLGAGAk	N-Term(TMT6plex); K9(TMT				1	4	ENSMUSP00000075067	ENSMUSP000000750
7	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	sPVVSGDTSPr	N-Term(TMT6plex)	X	X		1	6	P10636	ENSMUSP000001026
8	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	aQASAPAQAPk	N-Term(TMT6plex); K11(TMT				1	4	ENSMUSP00000080203	ENSMUSP000001297
9	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	siINNAEK	N-Term(TMT6plex); K7(TMT				1	4	ENSMUSP00000119975	ENSMUSP000001160
10	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	eIEVGGGR	N-Term(TMT6plex)				1	2	ENSMUSP00000073880	ENSMUSP000001005
11	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vTGGAAASK	N-Term(TMT6plex); K8(TMT				1	2	ENSMUSP00000099843	ENSMUSP000001312
12	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vATVSLPR	N-Term(TMT6plex)	X	X		1	1	P00761	P00761
13	<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vATVSLPR	N-Term(TMT6plex)	X	X		1	1	P00761	P00761

❖ **To load a filter set**

1. Choose **View > Display Filter**.

The Display Filter pane appears, as shown in [Figure 171](#) on [page 255](#).

2. Click the **Load** icon, **Load...**

3. In the Load Filter Set dialog box, browse to the filter set (*name.filterset*) file that you want to load, or type the name of the file in the File Name box.

4. Click **Open**.

The selected filter set now appears in the Display Filter pane.

When you load filter sets that do not apply to the type of pages that the report contains, the application still preserves the filters but makes them unavailable. When you load filter sets with conditions that refer to unavailable properties (columns) in the report, application shows the condition as unavailable with a line through it.

For examples of inapplicable but preserved filters, see [“Example 4”](#) on [page 266](#).


5. To remove conditions that are temporarily disabled, click **Remove** next to the condition. To remove all such unused conditions from all filters, click **Remove Unused Items**.

❖ **To view the filter settings in the Result Summaries**

1. Open a .pdResult file.
2. Choose **View > Summaries**.
3. In the Result Summaries page, click the **Validation** tab.
4. Click the **Filters** tab in the upper right corner of the Validation page.

The filters that you set in the Display Filter pane are reflected in the Applied Display Filters area of the Filters and Counts page of the Validation page, as shown in [Figure 177](#).


Figure 177. Filters shown in the Display Filter pane and in the Result Summaries pane



Following are some examples showing how to use display filters.


Example 1

Suppose that you would like to view only PSMs that are associated with a protein that contains the number “076070” in their corresponding protein accessions.

1. Open a .pdResult file.
2. Click the **Display Filter** icon, , on the main toolbar.
3. In the left pane of the Display Filter pane, select **PSMs**.
4. In the right pane of the window, click in the Add Property box and select **Protein Accessions**.
5. Click in the operator box, and choose **Contains** if it is not already selected.
6. In the values box to the right of the operator box, type **O76070**.

The Display Filter pane now resembles [Figure 178](#).

Figure 178. Filtering protein accessions



7. Click the **Apply** icon,  **Apply**.

[Figure 179](#) shows only entries that contain O76070 in their protein accessions.

Figure 179. Entries on the PSMs page that contain 076070 in their protein accessions

Protein accessions containing 076070

The screenshot shows the Proteome Discoverer software interface. The 'Display Filter' pane is active, showing a filter for 'Protein Accessions contains 076070'. The main table displays the results of this filter, showing various PSMs with columns for Contaminant, Standard, Decay, # Protein Groups, # Proteins, Master Protein Accessions, Protein Accessions, # Missed Cleavages, Charge, ΔScore, ΔCn, Rank, Search Engine Rank, m/z [Da], MH+ [Da], and ΔM [ppm].

Contaminant	Standard	Decay	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions	# Missed Cleavages	Charge	ΔScore	ΔCn	Rank	Search Engine Rank	m/z [Da]	MH+ [Da]	ΔM [ppm]
	X		1	1	076070	076070	1	2	0.0290	0.0000	1	1	1231.68005	2462.35283	9.29
	X		1	1	076070	076070	1	2	0.0292	0.0000	1	1	1231.68005	2462.35283	9.29
ex): K20(TM	X		1	1	076070	076070	0	2	0.0290	0.0000	2	2	1231.68005	2462.35283	1.05
ex): K20(TM	X		1	1	076070	076070	0	2	0.0292	0.0000	2	2	1231.68005	2462.35283	1.05
	X		1	1	076070	076070	1	2	0.0323	0.0000	1	1	1231.68005	2462.35283	9.29
ex): K2(TMT	X		1	1	076070	076070	1	3	0.6867	0.0000	1	1	974.22461	2920.65927	1.18
ex): K2(TMT	X		1	1	076070	076070	1	3	0.6408	0.0000	1	1	974.22461	2920.65927	1.18
ex): K2(TMT	X		1	1	076070	076070	1	3	0.6252	0.0000	1	1	974.22369	2920.65653	0.24
ex): K2(TM	X		1	1	076070	076070	0	2	0.0323	0.0000	2	2	1231.68005	2462.35283	1.05
ex): K2(TMT	X		1	1	076070	076070	1	3	0.6658	0.0000	1	1	974.22461	2920.65927	1.18
ex): K2(TMT	X		1	1	076070	076070	1	3	0.5997	0.0000	1	1	974.22363	2920.65635	0.18
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	2	3	0.0029	0.0000	1	1	769.45703	2306.36654	10.14
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	1	3	0.0029	0.0000	2	2	769.45703	2306.36654	1.35
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	2	3	0.0088	0.0000	1	1	769.45551	2306.35196	8.16
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	1	3	0.0088	0.0000	2	2	769.45551	2306.35196	-0.64
	X		1	1	076070	076070	1	3	0.0403	0.0000	1	1	821.45789	2462.35910	11.84
	X		1	1	076070	076070	1	3	0.0393	0.0000	1	1	821.45789	2462.35910	11.84
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	2	3	0.0045	0.0000	1	1	769.45648	2306.35489	9.43
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	1	3	0.0045	0.0000	2	2	769.45648	2306.35489	0.63
ex): K2(TMT	X		1	1	076070	076070	0	2	0.3379	0.0000	1	1	952.02136	1903.03545	1.08
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	2	3	0.0093	0.0000	1	1	769.45703	2306.36654	10.14
	X		1	1	076070	076070	1	3	0.0418	0.0000	1	1	821.45789	2462.35910	11.84
ex): K20(TM	X		1	1	076070	076070	0	3	0.0403	0.0000	2	2	821.45789	2462.35910	3.59
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	1	3	0.0093	0.0000	2	2	769.45703	2306.36654	1.35
ex): K2(TMT	X		1	1	076070	076070	0	2	0.3375	0.0000	1	1	952.02136	1903.03545	1.08
ex): K20(TM	X		1	1	076070	076070	0	3	0.0393	0.0000	2	2	821.45789	2462.35910	3.59
ex): K20(TM	X		1	1	076070	076070	1	3	0.0523	0.0000	1	1	821.45782	2462.35892	11.76
ex): K2(TMT	X		1	2	076070	076070; ENSMUSP0000002	2	3	0.0063	0.0000	1	1	769.45551	2306.35196	8.16

Example 2

The following example extends the filter defined in [Example 1](#) by adding a second condition that requires a top search rank.

1. In the Display Filter pane shown in [Figure 178](#) on [page 261](#), click in the Add Property box, and select **Search Engine Rank**.
2. Leave the default **Is Equal To** setting in the operator box.
3. In the values box to the right of the operator box, type **1**.

The Display Filter pane now resembles [Figure 180](#).

Figure 180. Search Engine Rank added as a second filter




4. Click the **Apply** icon,  **Apply**.

Figure 181 shows the results of filtering PSMs by protein accessions that contain the word “076070” and a search engine rank of “1.”

Figure 181. Entries on the PSMs page that contain 076070 in their protein accessions and 1 in their search engine rank

Search engine rank of 1

	Contaminant	Standard	Decoy	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions	# Missed Cleavages	Charge	ΔScore	ΔCn	Rank	Search Engine Rank	m/z [Da]	MH+ [Da]	ΔM [ppm]
		X		1	1	O76070	O76070	1	2	0.0290	0.0000	1	1	1231.68005	2462.35283	9.29
		X		1	1	O76070	O76070	1	2	0.0292	0.0000	1	1	1231.68005	2462.35283	9.29
		X		1	1	O76070	O76070	1	2	0.0323	0.0000	1	1	1231.68005	2462.35283	9.29
ex): K2(TMT)		X		1	1	O76070	O76070	1	3	0.6867	0.0000	1	1	974.22461	2920.65927	1.18
ex): K2(TMT)		X		1	1	O76070	O76070	1	3	0.6408	0.0000	1	1	974.22461	2920.65927	1.18
ex): K2(TMT)		X		1	1	O76070	O76070	1	3	0.6252	0.0000	1	1	974.22369	2920.65653	0.24
ex): K2(TMT)		X		1	1	O76070	O76070	1	3	0.6658	0.0000	1	1	974.22461	2920.65927	1.18
ex): K2(TMT)		X		1	1	O76070	O76070	1	3	0.5997	0.0000	1	1	974.22363	2920.65635	0.18
ex): K2(TMT)		X		1	2	O76070	O76070; ENSMUSP0000002	2	3	0.0029	0.0000	1	1	769.45703	2306.35654	10.14
ex): K2(TMT)		X		1	2	O76070	O76070; ENSMUSP0000002	2	3	0.0088	0.0000	1	1	769.45551	2306.35196	8.16
		X		1	1	O76070	O76070	1	3	0.0403	0.0000	1	1	821.45789	2462.35910	11.84
		X		1	1	O76070	O76070	1	3	0.0393	0.0000	1	1	821.45789	2462.35910	11.84
ex): K2(TMT)		X		1	2	O76070	O76070; ENSMUSP0000002	2	3	0.0045	0.0000	1	1	769.45648	2306.35489	9.43
ex)		X		1	1	O76070	O76070	0	2	0.3379	0.0000	1	1	952.02136	1903.03545	1.08
ex): K2(TMT)		X		1	2	O76070	O76070; ENSMUSP0000002	2	3	0.0093	0.0000	1	1	769.45703	2306.35654	10.14
		X		1	1	O76070	O76070	1	3	0.0418	0.0000	1	1	821.45789	2462.35910	11.84
ex)		X		1	1	O76070	O76070	0	2	0.3375	0.0000	1	1	952.02136	1903.03545	1.08
		X		1	1	O76070	O76070	1	3	0.0523	0.0000	1	1	821.45782	2462.35892	11.76
ex): K2(TMT)		X		1	2	O76070	O76070; ENSMUSP0000002	2	3	0.0063	0.0000	1	1	769.45551	2306.35196	8.16
		X		1	2	O76070	O76070; ENSMUSP0000002	1	2	0.0349	0.0000	1	1	924.52850	1848.04973	11.99
		X		1	1	O76070	O76070	1	3	0.0462	0.0000	1	1	821.45551	2462.35196	8.93
		X		1	1	O76070	O76070	1	3	0.0546	0.0000	1	1	821.45520	2462.35105	8.56
ex)		X		1	1	O76070	O76070	0	2	0.2782	0.0000	1	1	952.01990	1903.03252	-0.46
ex)		X		1	1	O76070	O76070	0	2	0.2534	0.0000	1	1	952.02228	1903.03728	2.05
ex)		X		1	1	O76070	O76070	0	2	0.2589	0.0000	1	1	952.02246	1903.03765	2.24
		X		1	2	O76070	O76070; ENSMUSP0000002	1	2	0.0528	0.0000	1	1	924.52881	1848.05034	12.32
		X		1	2	O76070	O76070; ENSMUSP0000002	2	3	0.0262	0.0000	1	1	693.06763	2077.18833	8.72
		X		1	1	O76070	O76070	1	3	0.0602	0.0000	1	1	821.45587	2462.35306	9.38

Show Associated Tables

Ready 792 Protein Groups; 2170 Proteins; 2911 Peptide Groups; 63/9775 PSMs; 72125 MS/MS Spectrum Info; 72155 Quan Spectra

Both conditions must be met for a data row to pass the filter because they are connected by the AND Boolean operator. If you want only one condition to apply, click in the Boolean operator box, which displays AND in [Figure 181](#), and select **OR**.


Example 3

Suppose that you want to build a more complex filter from the examples shown in “[Example 1](#)” on [page 261](#) and “[Example 2](#)” on [page 262](#). The new filter corresponds to “protein accession contains O76070” and “search rank is equal to 1” or “protein accession does not contain O76070.”

1. In the Display Filter pane shown in [Figure 180](#) on [page 263](#), click **Add Group**.


The Display Filter pane displays a new group, as shown in [Figure 182](#).

Figure 182. Adding a group



2. In the Boolean operator box of the new condition, select **OR**, as shown in [Figure 183](#).

Figure 183. Selecting the OR Boolean operator




3. Click in the Add Property box at the lowest level, as shown in [Figure 184](#), and select **Protein Accessions**.

You must add the new condition to the outer level corresponding to OR rather than to the inner level corresponding to AND.

4. In the operator box, select **Does Not Contain**.
5. In the values box to the right of the operator box, type **O76070**.

The Display Filter pane resembles [Figure 184](#).

Figure 184. OR condition added



6. Click the **Apply** icon, **Apply**.

Figure 185 shows the results of the filtering. In this case, all proteins were filtered out.

Figure 185. Entries on the PSMs page that contain O76070 in their protein accessions and 1 in their search engine rank or do not contain O76070 in their protein accessions



Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	Contaminant	Standard	Decoy	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	IVRPEVDVMcTAFHDNEETF	N-Term(TMT6plex); C10(Car	X	X		1	1	P02768	P02768
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSYSYSGC	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSYSYSGC	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	IVRPEVDVMcTAFHDNEETF	N-Term(TMT6plex); C10(Car	X	X		1	1	P02768	P02768
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSYSYSGC	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gGGGGYGSGSSYSYSGC	N-Term(TMT6plex)	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	nQEAGHEAAGEEAEASGF	N-Term(TMT6plex); K22(TM				1	1	ENSMUSP00000109741	ENSMUSP0000010
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	nQEAGHEAAGEEAEASGF	N-Term(TMT6plex); K22(TM				1	1	ENSMUSP00000109741	ENSMUSP0000010
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gLADGVdKDAVNAAIQQA	N-Term(TMT6plex); K9(TMT				1	1	ENSMUSP00000059678	ENSMUSP0000005
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVlKPGMVVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP0000004
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Selected	TKEQANAVSEAVVSSNTV	K22(TMT6plex)		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Selected	TKEQANAVSEAVVSSNTV	K22(TMT6plex)		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVlKPGMVVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP0000004
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVlKPGMVVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP0000004
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	gLADGVdKDAVNAAIQQA	N-Term(TMT6plex); K9(TMT				1	1	ENSMUSP00000059678	ENSMUSP0000005
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	sGGGGGGGGcGGGGVS	N-Term(TMT6plex); C10(Car	X	X		1	1	P13645	P13645
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	eGTGSTSSGSAGGAVGk	N-Term(TMT6plex); K19(TM				1	1	ENSMUSP00000022925	ENSMUSP0000002
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Selected	TKEQANAVSEAVVSSNTV	K22(TMT6plex)		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	hKEQANAVSEAVVSSNTVA	N-Term(TMT6plex); K2(TMT		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVlKPGMVVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP0000004
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vRPeVYGGAEIGQQIR	N-Term(TMT6plex); C4(Car				1	1	ENSMUSP0000000804	ENSMUSP0000000
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	aSTSSSSSSSNQqTEk	N-Term(TMT6plex); K17(TM	X	X		1	2	P51965	P51965; ENSMUSP
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	fSScGGGGGSGFAGGGFg	N-Term(TMT6plex); C4(Car	X	X		1	1	P04264	P04264
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	hKEQANAVSEAVVSSNTVA	N-Term(TMT6plex); K2(TMT		X		1	1	O76070	O76070
<input checked="" type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	vETGVlKPGMVVTFAPVNV	N-Term(TMT6plex); K7(TMT				1	1	ENSMUSP00000042457	ENSMUSP0000004

Ready 792 Protein Groups; 2170 Proteins; 2911 Peptide Groups; 9750/9775 PSMs; 72125 MS/MS Spectrum Info; 72155 Quan Spectra

Example 4


You can save and load filter sets to and from the disk, even though they might not entirely apply to another report. As an example, assume that you create a filter set that only displays protein groups with more than four associated proteins and proteins that are connected to at least ten identified peptides. [Figure 186](#) shows these filters.

Figure 186. Specifying protein group filters and protein filters




When you load this filter set into a different results report that does not contain the Grouped Protein Groups page, the protein groups filter is still preserved but is unavailable, as shown in [Figure 187](#).


Figure 187. Preserved but unavailable protein groups filter



When you select the unavailable filter set, the message shown in [Figure 188](#) appears.


Figure 188. Selecting the preserved but unavailable protein groups filter



To remove all filters that refer to unavailable pages, click the **Clear** icon,  **Clear**.


Now suppose that you define the filter for the heavy-light ratio of a quantification report that only shows ratios above 2.5. [Figure 189](#) shows this filter.


Figure 189. Specifying a quantification filter



When you load this filter into a report that does not contain quantification information and therefore has no Heavy/Light column, this condition is unavailable, so it appears with a line drawn through it, as shown in [Figure 190](#).

Figure 190. Loading a quantification filter into a report containing no quantification data



The application still preserves such conditions in the internal representation of the filter but ignores them for the current report. You can remove a temporarily disabled condition by clicking the Remove button next to the condition. To remove all such unused conditions from all filters, click the **Clear** icon,  **Clear**.

Filtering with Filter Sets

You can save a set of your favorite display filters as the default filters of a .pdResult file. You set the display filters and save them to a .filterSet file, and then apply them in a newly generated .pdResult file by using the Display Filter node in the consensus workflow.

Saving a Filter Set in a .pdResult File

❖ **To save a filter set in a .pdResult file**

1. Open the .pdResult file and set the filters in the Display Filter pane, following the instructions in [“Filtering with Display Filters”](#) on [page 254](#).
2. Save the filter set as described in [“To save a filter set”](#) on [page 258](#).
3. (Optional) Choose **File > Save**.

This command saves the current layout and filter set in the .pdResult file and in the .filterSet file. It also saves the filter set in the following directory in the Common Templates area:

```
C:\Users\Public\Public Documents\Thermo\Proteome Discoverer 2.1\  
Common Templates\FilterSets
```

You can apply this filter set as a template when you generate other .pdResult files by using the Display Filter node.

Applying Stored Filters to a .pdResult File


Use the following procedure to apply stored filters to a .pdResult file or to temporarily override the display of only master proteins on the Proteins page of the .pdResult file.

❖ To apply stored filters to a .pdResult file

1. In an open study, create a consensus workflow. For instructions, see [“Creating a Consensus Workflow”](#) on page 112.
2. Drag the Display Filter node to the Post-Processing Nodes pane.
3. Click the **Display Filter** node.
4. Click the **Browse** button (...) next to the Filter Set parameter.
5. In the Parameter Text for Filter Set (Read Only) dialog box, shown in [Figure 191](#), click **Load File**.

Note The default filter set in this dialog box displays only master proteins by default on the Proteins page.

Figure 191. Parameter Text for Filter Set (Read Only) dialog box



6. In the Open dialog box, browse to the .filterSet file that you want to load, and click **Open**.
7. In the Parameter Text for Filter Set dialog box, click **OK**.
The new filter set replaces the default filter set.
8. Run the consensus workflow.
9. (Optional) Open the .pdResult file that contains the filter set.

Changing the Default Filter for the Proteins Page

By default, the Proteins page displays only master proteins. If you want the Proteins page to display a different set of proteins by default, follow these steps.

❖ To change the default filter set for the Proteins page

1. Choose **Administration > Configuration**.
2. Click the **Browse (...)** button next to the Default Filter Set parameter.
3. In the Parameter Text for Filter Set (Read Only) dialog box, shown in [Figure 191](#), click **Load File**.

4. In the Open dialog box, browse to the new default .filterSet file that you want to load, and click **Open**.
5. In the Parameter Text for Filter Set dialog box, click **OK**.

The new default filter set replaces the factory default filter set.


Using Display Filters to Filter Numerical Values by a Specified Precision Value

In reporting, the application displays most numerical values with a significantly lower precision than it uses to store the values internally. The display filters can filter with high precision, but you might not want to choose such precision during filtering. In addition, the precision might be data-dependent; for example, filtering the delta mass requires greater precision than is required in filtering probabilities.

If an operator tests for the equality of two values, you can control the numerical precision with which the application filters data that are represented as floating-point values. The filter automatically extends any value that you enter to at least two digits for the fractional part of the value.

To use a higher precision, enter a value with more digits. The example in [Figure 192](#) uses a precision of five digits (3.22528) in the filter. The filter for such values automatically rounds the internal data before filtering according to the precision of the given value.

Figure 192. Automatic rounding when filtering high-precision data



Filtering numerical data inherently involves a precision problem. Numbers categorized as “double” or “float” follow the IEEE 754 format, which specifies that they should always be compared to a tolerance interval, which is effectively equivalent to rounding. This comparison is especially important in the Proteome Discoverer application, which stores data internally with full precision but rounds them for user convenience in the displayed pages. For example, although the data shown in the XCorr column in Figure 193 is only displayed with two digits of precision (3.23), the internal precision of the value is higher (3.2253262996673584), as shown in the ToolTip.

5 Filtering Data

Filtering with Display Filters

Figure 193. Example data with floating-point values

ProteinCenter Distribution Marker Quan report of Sequest HT TMT 10-plex search of HF Q Exactive data

Display Filter: AND XCorr is greater than 3.22528


Ion Inference [%]	Ion Inject Time [ms]	RT [min]	First Scan	Spectrum File	Ions Matched	XCorr	Reporter Quan Result ID	Quan Info	Quan Usage	(127_N) / (127_C) / (128_N) / (128_C)
9	8.167	126.9514	50882	141007KSAM00279TMT20-f3.raw	0/0	3.59	103349	Unique	Used	3.074 4.024 32.352
0	5.553	135.3280	54479	141007KSAM00279TMT20-f3.raw	0/0	3.59	106774	Unique	Used	1.881 3.752 3.174
0	50.000	96.9789	37992	141007KSAM00279TMT20-f3.raw	0/0	3.59	91106	Unique	Used	1.465 1.773 5.671
0	50.000	198.3035	81306	141007KSAM00279TMT20-f3.raw	0/0	3.59	131957	Unique	Used	1.004 0.765 8.713
9	33.001	122.0221	48762	141007KSAM00279TMT20-f3.raw	0/0	3.59	101330	Unique	Used	1.100 0.732 7.830
16	9.251	91.3153	35552	141007KSAM00279TMT20-f3.raw	0/0	3.59	88782	Unique	Used	1.717 4.877 17.760
0	22.220	97.8687	38374	141007KSAM00279TMT20-f3.raw	0/0	3.56	91469	Unique	Used	0.775 0.736 0.999
0	50.000	77.2144	29482	141007KSAM00279TMT20-f3.raw	0/0	3.55	83001	Unique	Used	1.413 0.957 5.202
0	8.744	113.8851	45262	141007KSAM00279TMT20-f3.raw	0/0	3.52	97994	Unique	Used	0.965 0.820 1.008
0	50.000	103.2807	40703	141007KSAM00279TMT20-f3.raw	0/0	3.47	93688	Unique	Used	0.960 0.925 2.622
9	9.202	110.9069	43986	141007KSAM00279TMT20-f3.raw	0/0	3.46	96799	Unique	Used	0.841 0.809 1.016
0	49.928	120.8771	48270	141007KSAM00279TMT20-f3.raw	0/0	3.44	100861	Unique	Used	1.169 0.910 3.534
0	50.000	70.1461	26463	141007KSAM00279TMT20-f3.raw	0/0	3.44	80126	Unique	Used	1.256 1.327 6.123
0	50.000	80.6061	30936	141007KSAM00279TMT20-f3.raw	0/0	3.42	84386	No Quan Labels	Not Used	1.606 1.340 4.484
0	50.000	100.8772	39668	141007KSAM00279TMT20-f3.raw	0/0	3.41	92702	Unique	Used	1.374 1.108 6.571
0	27.502	100.5271	39517	141007KSAM00279TMT20-f3.raw	0/0	3.41	92558	Unique	Used	1.180 0.876 3.574
0	50.000	177.6078	72540	141007KSAM00279TMT20-f3.raw	0/0	3.39	123859	Unique	Used	1.508 2.531 7.543
0	21.782	68.7110	25853	141007KSAM00279TMT20-f3.raw	0/0	3.38	79545	Unique	Used	0.855 1.053 2.878
0	50.000	72.0326	27266	141007KSAM00279TMT20-f3.raw	0/0	3.36	80890	Unique	Used	3.695 5.945 21.255
0	50.000	69.1354	26033	141007KSAM00279TMT20-f3.raw	0/0	3.33	79711	Unique	Used	0.995 0.812 3.498
0	50.000	81.2021	31192	141007KSAM00279TMT20-f3.raw	0/0	3.31	84629	Unique	Used	1.137 0.821 4.157
0	50.000	60.4115	22265	141007KSAM00279TMT20-f3.raw	0/0	3.29	76128	Unique	Used	1.567 1.088 6.308
0	50.000	70.5791	26647	141007KSAM00279TMT20-f3.raw	0/0	3.25	80301	Unique	Used	2.165 0.806 6.908
0	45.226	112.4828	44661	141007KSAM00279TMT20-f3.raw	0/0	3.23	97443	Unique	Used	0.896 1.121 1.311
0	50.000	115.3900	45908	141007KSAM00279TMT20-f3.raw	0/0	3.23	98612	Unique	Used	1.305 0.876 4.606

Ready | 792 Protein Groups; 2170 Proteins; 2911 Peptide Groups; 9773/9775 PSMs; 72125 MS/MS Spectrum Info; 72155 Quan Spectra

Finding Common or Unique Proteins in Multiple Searches

A multiconsensus report can help you find common or unique proteins between two or more searches. You can use the display filters of the Proteins, PSMs, or MS/MS Spectrum Info page to automate the task of creating a list of all proteins unique to one search or search node or a list of all proteins common to all searches or search nodes.


❖ To find proteins common to two or more searches


1. Choose **View > Display Filter**, or click the **Display Filter** icon, .
2. Select **Proteins** in the left pane of the Display Filters pane.
3. In the right pane, click in the Add Property box and select **Coverage *first_search_engine***, for example, **Coverage Mascot**.

4. In the operator box, select **Is Greater Than**.
5. In the values box, select **0**.
6. Click in the Add Property box again, and select **Coverage *second_search_engine***, for example, **Coverage Sequest HT**.
7. In the operator box, select **Is Greater Than**.
8. In the values box, select **0**.

The filter in the Display Filters pane now resembles [Figure 194](#).


Figure 194. Filtering to select common proteins




9. For each additional search engine, repeat [step 3](#) through [step 5](#) earlier in this procedure.
10. Click the **Apply** icon,  **Apply**.

The Coverage *Search_engine* columns in the results report display the proteins found by both searches. The Coverage Mascot and Coverage Sequest HT columns display these proteins in [Figure 195](#).

Figure 195. Multiconsensus report showing proteins common to Mascot and Sequest HT search in results report





❖ To find proteins unique to each search

1. Choose **View > Display Filter**, or click the **Display Filter** icon, .
2. Select **Proteins** in the left pane of the Display Filter pane.
3. In the right pane, click in the Add Property box, and select **# PSMs first_search_engine**, for example, **# PSMs Mascot**.
4. In the operator box, select **Is Not Set**.

The filter in the Display Filter pane now resembles [Figure 196](#).

Figure 196. Filtering to select unique proteins



5. Click the **Apply** icon,  **Apply** .

Zeros in the Coverage *Search_engine* column in the results report indicate the proteins found only by the selected search. Figure 197 displays the proteins unique to the search with Mascot.

Figure 197. Multiconsensus report showing proteins unique to the Mascot search

Checked	Master	Accession	Description	Coverage	Score Mascot	Coverage Mascot	Peptides Mascot	# PSMs Mascot	Score Sequest	Coverage Sequest HT	#
<input type="checkbox"/>		IPI00795257.1	32 kDa protein	20%		0%			78.97	20%	
<input type="checkbox"/>		IPI00789134.1	Glyceraldehyde 3-phosphate dehydrogenase	22%		0%			78.97	22%	
<input type="checkbox"/>		IPI00382470.3	heat shock protein 90kDa alpha (cytosolic), class A member	15%		0%			69.91	15%	
<input type="checkbox"/>		IPI00037070.3	Isoform 2 of Heat shock cognate 71 kDa protein	19%		0%			63.04	19%	
<input type="checkbox"/>		IPI00889534.1	carbamoyl-phosphate synthetase 1 isoform a precursor	6%		0%			58.55	6%	
<input type="checkbox"/>		IPI00894498.1	Beta actin variant (Fragment)	11%		0%			57.75	11%	
<input type="checkbox"/>		IPI00894365.2	cDNA FLJ52842, highly similar to Actin, cytoplasmic 1	12%		0%			57.75	12%	
<input type="checkbox"/>		IPI00910870.1	cDNA FLJ59163, highly similar to Heat shock cognate 71 kD	25%		0%			43.17	25%	
<input type="checkbox"/>		IPI00908770.1	cDNA FLJ53063, highly similar to Tubulin beta-7 chain	16%		0%			41.85	16%	
<input type="checkbox"/>		IPI00788737.1	39 kDa protein	13%		0%			39.81	13%	
<input type="checkbox"/>		IPI00797221.5	Glyceraldehyde-3-phosphate dehydrogenase	17%		0%			39.81	17%	
<input type="checkbox"/>		IPI00794605.1	18 kDa protein	7%		0%			39.16	7%	
<input type="checkbox"/>		IPI00909140.1	cDNA FLJ53012, highly similar to Tubulin beta-7 chain	17%		0%			38.80	17%	
<input type="checkbox"/>		IPI00604607.2	Hsp99-alpha-delta-N	15%		0%			38.70	15%	
<input type="checkbox"/>		IPI00759806.1	Isoform MBP-1 of Alpha-enolase	14%		0%			37.36	14%	
<input type="checkbox"/>		IPI00790339.1	22 kDa protein	15%		0%			36.56	15%	
<input type="checkbox"/>		IPI00796881.1	16 kDa protein	20%		0%			36.56	20%	
<input type="checkbox"/>		IPI00908605.1	cDNA FLJ59940, highly similar to Tubulin beta-2C chain	15%		0%			35.69	15%	
<input type="checkbox"/>		IPI00908398.1	cDNA FLJ52029, highly similar to Tubulin beta-7 chain	15%		0%			35.60	15%	
<input type="checkbox"/>		IPI00375531.2	non-metastatic cells 1, protein (NM23A) expressed in isofon	24%		0%			34.55	24%	
<input type="checkbox"/>		IPI00911039.1	cDNA FLJ54408, highly similar to Heat shock 70 kDa protei	11%		0%			33.40	11%	
<input type="checkbox"/>		IPI00845339.1	cDNA FLJ54392, highly similar to Heat shock 70 kDa protei	10%		0%			33.40	10%	
<input type="checkbox"/>		IPI00910742.1	cDNA FLJ53341, highly similar to Tubulin beta-4 chain	8%		0%			32.52	8%	
<input type="checkbox"/>		IPI00797270.2	cDNA FLJ52069, highly similar to Triosephosphate isomera	30%		0%			31.73	30%	
<input type="checkbox"/>		IPI00027107.5	Tu translation elongation factor, mitochondrial precursor	11%		0%			31.27	11%	
<input type="checkbox"/>		IPI00910482.1	cDNA FLJ54407, highly similar to Heat shock 70 kDa protei	9%		0%			30.25	9%	
<input type="checkbox"/>		IPI00909560.1	cDNA FLJ53645, highly similar to Pyruvate kinase isozyme	14%		0%			30.05	14%	
<input type="checkbox"/>		IPI00220644.8	Isoform M1 of Pyruvate kinase isozymes M1/M2	13%		0%			30.05	13%	

You can repeat the procedure using # PSMs *Search_engine* to find the proteins unique to the *Search_engine* search.


In general, if one search or search node uniquely identifies a protein, the number of identified peptides for this protein is non-zero for this search or search node but zero in all other searches or search nodes. On the other hand, if a protein is common to two or more searches or search nodes, the number of identified peptides is non-zero for all considered searches or search nodes.

Applying Filters Specific to Different Searches in Multiconsensus Reports

You cannot apply filters specific to different searches in multiconsensus reports, but you can apply filters that are specific to the different search engines used (Sequest HT and Mascot). The available score-based filters automatically change the list of scores that you can select for filtering, depending on the search engines used for the loaded reports.


For a Sequest HT search report, you can apply the XCorr and Charge display filters, as shown in [Figure 198](#).

Figure 198. XCorr and Charge filters applied to Sequest HT multiconsensus reports



For a Mascot search report, you can apply the Ions Score display filter, as shown in [Figure 199](#).

Figure 199. Ions Score filter applied to Mascot multiconsensus reports



Validating Results

This chapter explains how to validate your results by calculating false discovery rates.

Contents

- [Target FDRs](#)
- [Peptide Confidence Indicators](#)
- [Setting Up FDRs](#)
- [Calculating FDRs for PSMs, Peptide Groups, Proteins, and Protein Groups](#)

The false discovery rate (FDR), or the false positive rate, is a statistical value that estimates the number of false positive identifications among all identifications found by a peptide identification search. It is a measure of the certainty of the identification. You can use the Proteome Discoverer decoy database search feature to determine FDRs.

You can use FDRs to validate MS/MS searches of large data sets, but they are not effective on searches of a small number of spectra or searches against a small number of protein sequences, because the number of matches will likely be too small to give a statistically meaningful estimate.

A decoy database strategy gives a probability value to identifiers and the percentage of false discoveries that you can expect. A one percent false discovery rate is a typical target for searches.

A good decoy database should contain entries that resemble real proteins but do not contain genuine peptide sequences. The simplest approach to obtaining such a decoy database is to reverse all protein sequences, which is the scheme that the Proteome Discoverer application uses by default. However, you can request that the search be performed against the concatenated database, if required. It is a suitable approach for enzymatic MS/MS searches.

CAUTION Reversing the database is not suitable for peptide mass fingerprinting or no-enzyme MS/MS searches, especially for dynamic modifications. You might see mass shifts at each end of a peptide sequence that transform a genuine y-series match into a false b-series match or vice versa.

You can perform the decoy database search in two ways:

- Perform two separate searches, one against the non-decoy database and one against the decoy database. Then count the number of matches from both searches to determine the false discovery rates. This approach is the more conservative approach.
- Create a concatenated database from the non-decoy and the decoy database and then perform the search against this concatenated database.

The difference between the two approaches becomes clear when you find two significant matches for a given spectrum. The first match is from the non-decoy database, and the second one is from the decoy database. Because the Proteome Discoverer application considers only the top matches when calculating the FDRs, finding two significant matches for a given spectrum is not considered a false positive in the concatenated database approach, but it counts in the separate databases approach. The latter case is more conservative and is the default approach that the application currently uses. You can select the validation methodology by setting the Concatenated FDR Calculation parameter of the Peptide Validator node.

To calculate the FDR, the application counts the matches that pass a given set of filter thresholds from the decoy database and from the non-decoy database. It counts only the top match per spectrum, assuming that for any given spectrum only one peptide can be the correct match.

For detailed information on how the application calculates the FDRs for PSMs, peptide groups, proteins, and protein groups, refer to the Help.

Target FDRs

If you set an FDR target value for a decoy database search, the application determines and applies filter thresholds to identified matches so that the resulting FDR is not higher than the set target value. The application distributes the confidence indicators applied to each peptide match according to these calculated filter thresholds.


You must specify two target values for a decoy database search: a strict target FDR and a more relaxed FDR. [Figure 201](#) on [page 280](#) shows the decoy search setting with target FDRs of one percent and five percent, respectively. After completing the search, the system automatically determines two sets of filter settings so that the resulting separate FDRs do not exceed their corresponding target value.

Peptide Confidence Indicators

The application uses the filter settings that determine FDRs to distribute the confidence indicators for the peptide matches; these are the green, yellow, and red circles attached to each peptide match. Whenever you perform a decoy database search during the database search and apply filter settings to achieve the specified target FDRs, the same filters are used to

distribute the confidence indicators. Peptide matches that pass the filter associated with the strict FDR are assigned a green high-confidence indicator, peptide matches that pass the filter associated with the relaxed FDR but do not pass the filter associated with the strict FDR are assigned a yellow medium-confidence indicator, and all other peptide matches receive a red low-confidence indicator. [Figure 200](#) gives an example of these confidence indicators.

Figure 200. Decoy database search results



	Checked	Confidence	Sequence	PSM Ambiguity	# Proteins	# Protein Groups	Protein Accessions	Modifications	Activation Type	ΔScore	ΔCn	Rank	Search Engine Rank	Concatenate
1	<input type="checkbox"/>	●	NSLDDPNLk	Unambiguous	2	1	Q10265; O59855	K9(TMT6plex)	HCD	0.4549	0.0000	1	1	1
2	<input type="checkbox"/>	●	QTVAVGVik	Unambiguous	2	1	Q10119; P50522	K9(TMT6plex)	HCD	0.5438	0.0000	1	1	1
3	<input type="checkbox"/>	●	aDASAVEAFR	Unambiguous	1	1	O60101	N-Term(TMT6plex)	HCD	0.6227	0.0000	1	1	1
4	<input type="checkbox"/>	●	QLIVAVNk	Unambiguous	2	1	Q10119; P50522	K8(TMT6plex)	HCD	0.3397	0.0000	1	1	1
5	<input type="checkbox"/>	●	qTVAVGVik	Unambiguous	2	1	Q10119; P50522	N-Term(TMT6plex); K9(TMT)	HCD	0.3799	0.0000	1	1	1
6	<input type="checkbox"/>	●	eLDYFSk	Unambiguous	1	1	O60101	N-Term(TMT6plex); K7(TMT)	HCD	0.2148	0.0000	1	1	1
7	<input type="checkbox"/>	●	aADEFLLk	Unambiguous	1	1	P40370	N-Term(TMT6plex); K8(TMT)	HCD	0.5028	0.0000	1	1	1
8	<input type="checkbox"/>	●	IPLQDVYk	Unambiguous	2	1	Q10119; P50522	N-Term(TMT6plex); K8(TMT)	HCD	0.6168	0.0000	1	1	1
9	<input type="checkbox"/>	●	aLEELLk	Unambiguous	1	1	O94489	N-Term(TMT6plex); K7(TMT)	HCD	0.1429	0.0000	1	1	1
10	<input type="checkbox"/>	●	vVDLVAYTAAk	Unambiguous	1	1	P78958	N-Term(TMT6plex); K11(TMT)	HCD	0.6107	0.0000	1	1	1
11	<input type="checkbox"/>	●	QTVAVGVik	Unambiguous	2	1	Q10119; P50522	K9(TMT6plex)	HCD	0.5105	0.0000	1	1	1
12	<input type="checkbox"/>	●	aDASAVEAFR	Unambiguous	1	1	O60101	N-Term(TMT6plex)	HCD	0.5280	0.0000	1	1	1
13	<input type="checkbox"/>	●	aADEFLLk	Unambiguous	1	1	P40370	N-Term(TMT6plex); K8(TMT)	HCD	0.3565	0.0000	1	1	1
14	<input type="checkbox"/>	●	QTVAVGVik	Selected	2	1	Q10119; P50522	K9(TMT6plex)	HCD	0.0331	0.0000	2	2	2
15	<input type="checkbox"/>	●	IPLQDVYk	Unambiguous	2	1	Q10119; P50522	N-Term(TMT6plex); K8(TMT)	HCD	0.1609	0.0000	1	1	1

Setting Up FDRs

You can set up FDRs with any of the following nodes in the workflow.

To validate the identifications from different search engines, the Proteome Discoverer application offers the Percolator node (refer to the Help) and the Target Decoy PSM Validator node (refer to the Help) in the processing workflow. Percolator is a superior validation algorithm that uses a machine learning approach, but it requires a sufficient number of target and decoy matches that are not always available. In these cases, you can use the Target Decoy PSM Validator node. This node triggers a target and decoy search and calculates score thresholds to achieve the specified target false discovery rate (FDR). The derived score thresholds for the strict and relaxed FDR separate the identified PSMs into high-, medium-, and low-confidence identifications.

[Figure 201](#) illustrates the basic processing workflow to use for setting up a decoy database search. For more information on creating a processing workflow, see [“Creating a Processing Workflow”](#) on page 104.

Note In this workflow, you can use the Percolator node in place of the Target Decoy PSM Validator node.

To work properly, Percolator needs a sufficient number of PSMs from the target and the decoy search. If the search identified fewer than 200 target or decoy PSMs, or if fewer than 20 percent decoy PSMs are available compared to the number of target matches, Percolator rejects them for processing and displays an appropriate message in the Proteome Discoverer job queue or in the Search Summary of an open report.

Figure 201. Basic processing workflow for setting up a decoy database search





Figure 202 illustrates the basic consensus workflow to use for setting up a decoy database search. For more information on creating a consensus workflow, see “[Creating a Consensus Workflow](#)” on page 112.

Figure 202. Basic consensus workflow for setting up a decoy database search

Calculating FDRs for PSMs, Peptide Groups, Proteins, and Protein Groups

This topic directs you to sources of detailed information on calculating FDRs for PSMs, peptide groups, proteins, and protein groups.

Calculating False Discovery Rates for PSMs

Use the Percolator node to calculate the FDRs for PSMs. For information on this node, refer to the Help and [“Filtering PSMs with the Percolator Node”](#) on page 247.

Calculating False Discovery Rates for Peptide Groups

Use the quality algorithm in the Percolator node to calculate the FDRs for peptide groups. For information on this node, refer to the Help.

Calculating False Discovery Rates for Proteins

Use the Protein FDR Validator node or the FidoCT Protein Validator node to calculate the FDRs for proteins. For detailed information on these nodes, refer to the Help.

6 Validating Results

Calculating FDRs for PSMs, Peptide Groups, Proteins, and Protein Groups

Calculating False Discovery Rates for Protein Groups

Use the FidoCT Protein Validator node to calculate the FDRs for protein groups. For information on this node, refer to the Help.

Grouping Peptides and Proteins

This chapter describes how the Proteome Discoverer application groups proteins and peptides.

Contents

- [Grouping Proteins](#)
- [Grouping PSMs](#)

Grouping Proteins

The application groups proteins through the Protein Grouping node in the consensus workflow. This topic explains the algorithm that it uses and the workflow that you should use to group proteins.

Master Proteins


For information on master proteins, see the Help.

Protein Grouping Algorithm

The application uses a protein grouping inference process to group proteins when you use the Protein Grouping node in the consensus workflow. [Figure 203](#) shows the steps involved in this process.

The application initially collects all peptide spectrum matches (PSMs) that were not filtered out by the Peptide and Protein filter node.

Figure 203. Protein grouping inference process in the Proteome Discoverer application



1. The application creates preliminary protein groups from the PSMs collected. It combines all proteins into one protein group that contains the same subset of peptides.

The application takes the next steps in the protein grouping process if you set the Apply Strict Parsimony Principle parameter of the Protein Grouping node in the consensus workflow to True.

2. The application removes all protein groups that have no unique peptides among the peptides that it considers for the protein grouping process. If a protein group does not contain at least one unique peptide, the application also includes all of the peptides included by other protein groups, so there is no supporting evidence for the existence of

this protein group. At this point, the application explicitly retains all protein groups that form circular rings of overlapping shared peptides. For example, suppose a circular ring comprises these protein groups:

- ABCD (identified by peptides a, b, c, and d)
- CDEF (identified by peptides c, d, e, and f)
- EFAB (identified by peptides e, f, a, and b)

To explain all identified peptides, only two of the three protein groups are needed, but at this point it is not clear which to take and which to reject. The application postpones the resolution of this issue until step 5.

3. The application first collects all spectra with more than one peptide match to consider for the protein grouping process. It then resolves these ambiguous cases and selects one of the PSMs to use for the protein grouping process while rejecting the remaining peptide matches of a spectrum. In cases where the application considers more than one PSM for a spectrum, it resolves this ambiguity by selecting the PSM that is connected to the “best” protein group and rejecting the other PSMs. The “best” protein group is the group with the highest number of unambiguous and unique peptides.
4. The application resolves the cases where protein groups form circular rings of overlapping identified peptides. This step is the last step of the protein group inference process, resulting in the final list of protein groups that the application reports in the Proteins page of the .pdResult file.

The PSM Ambiguity column on the PSMs and MS/MS Spectrum Info pages can help you understand the process of selecting PSMs for the protein group. This column is available for every PSM, every search input entry (representing the searched spectra), and every peptide group. For the search input entries and the peptide groups, this column displays the best PSM ambiguity from all connected PSMs. Refer to the Help for a description of the categories of ambiguity in this column.

Consider the example shown in [Figure 204](#), where the application identifies eight different PSMs for search input 20. Even though only seven PSMs ranked 1 through 7 are of high confidence, all eight PSMs meet the specified protein grouping criteria because, on the basis of user-specified criteria in the Peptide and Protein Filter node, all PSMs for the top-scored proteins were retained. However, because the search input cannot be assigned unambiguously to a single protein, the PSM ambiguity is set to Ambiguous.

7 Grouping Peptides and Proteins

Grouping Proteins

Figure 204. PSMs shown for search input


Protein Group ID	# Proteins	# Unique Peptides	# Peptides	# PSMs	Group Description
1	1612	29	0	15	2164 [Master Protein] Myoglobin (Horse Heart) Mutant With Ser 92 Replaced By Asp (S92d) [MASS=16979]
2	3721	27	1	15	1768 [Master Protein] Myoglobin (Horse Heart) Mutant With Leu 104 Replaced By Asn (L104n) [MASS=16952]
3	1838	27	1	15	1723 [Master Protein] Myoglobin - horse
4	2486	27	1	14	1679 [Master Protein] Myoglobin (Horse Heart) Mutant With His 64 Replaced By Tyr (H64y) [MASS=16978]
5	1187	27	1	14	1638 [Master Protein] Myoglobin Mutant With His 64 Replaced By Thr (H64t)
6	3815	1	3	8	562 [Master Protein] MYOGLOBIN
7	4214	1	4	8	558 [Master Protein] Myoglobin OS=Homo sapiens GN=MB PE=1 SV=2

Checked	Best PSM Ambiguity	Mass Analyzer	Activation Type	MS Order	# PS	Isolation Interference [%]	Ion Inject Time [ms]	Precursor m/z [Da]	Precursor MH+ [Da]	Precursor Charge	RT [min]	First Scan	Spectrum File
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	78	300.000	539.87659	1078.74590	2	35.6968	1874	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	48	300.000	539.87555	1078.74382	2	37.2210	1958	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	42	300.000	539.87555	1078.74382	2	42.6514	2261	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	73	300.000	412.53067	1647.10085	4	31.5065	1632	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	8	71	150.000	477.76028	954.51329	2	44.2932	2356	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	16	300.000	526.82050	1052.63371	2	80.1368	4510	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	79	300.000	413.92282	1239.75391	3	26.5586	1363	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	8	71	150.000	477.76028	954.51329	2	44.2932	2356	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	FTMS	HCD	MS2	8	15	300.000	627.36310	1880.07474	3	56.3918	3092	myo_Splex_hcd_pod_raw
<input checked="" type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	8	2	114.176	526.82245	1052.63762	2	68.4475	3805	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	8	7	150.000	526.82233	1052.63738	2	74.9539	4194	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	8	25	150.000	526.82294	1052.63860	2	83.3743	4692	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	7	2	114.176	526.82245	1052.63762	2	68.4475	3805	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	7	21	150.000	470.77271	1800.06899	4	46.0509	2454	myo_Splex_hcd_pod_raw
<input type="checkbox"/>	Ambiguous	ITMS	PGD	MS2	7	54	150.000	478.51993	1911.05788	4	55.3775	3026	myo_Splex_hcd_pod_raw

Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions	# Missed Cleavages	Charge	ΔScore	ΔCn	Rank	Search Engine R
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	DTDMLAASK		1	1	1171681	1171681	0	2	0.0000	0.0000	1	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	aLELFR	N-Term(TRAQ@plex)	8	11	999881; 1942750; 70561; 494711; 2554649; 127638; P02144	999881; 1942750; P68082; 7	0	2	0.0000	0.0000	2	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	aLELFR	N-Term(TRAQ@plex)	1	1	E1X7N1	E1X7N1	0	2	0.0000	0.0000	2	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	mSSLQGKMR	M1(Oxidation)	1	1	114701	114701	1	2	0.0145	0.0145	4	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	TTTASSKAAK		1	1	124110	124110	1	2	0.0290	0.0290	5	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	aLGFPAK	N-Term(TRAQ@plex)	1	4	E1X8C2	E1X8C2; G4QLT9; A5U1B1;	0	2	0.0435	0.0435	6	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	YGGRRPLDK		1	1	7450552	7450552	0	2	0.0483	0.0483	7	
<input type="checkbox"/>	●	Sequest HT (A2)	Ambiguous	YITEQLSAK		1	2	1055211	994791; 1055211	0	2	0.0966	0.0966	8	

Protein Grouping Workflow

To group proteins, use the Protein Grouping node in the consensus workflow. Figure 205 shows the basic workflow to use to group proteins with the Protein Grouping node.

Figure 205. Basic workflow for grouping proteins

For information on the parameters of the Protein Grouping node, refer to the Help.

Number of Unique Peptides Column on the Proteins Page

The application counts the number of peptides that are only contained in a protein group and displays the value in the # Unique Peptides column on the Proteins page. The Proteome Discoverer application counts only the peptides that display a status of Selected or Unambiguous in the PSM Ambiguity column, because assessing the uniqueness of peptides that were not used to form protein groups has no relevance.

PSMs Identified by Multiple Workflow Nodes

In search results where the application identifies PSMs by multiple search nodes within a single workflow, the protein grouping algorithm selects one of the PSMs identified for the same spectrum for building the protein groups.

In search results where the application identifies PSMs by multiple search nodes from multiple workflows (multiconsensus report), it treats PSMs and spectra from the different workflows as being separate, even if it searched the same raw data files and therefore the same spectra. In this case, it is difficult to determine whether the application searched the exact same spectra, because they might have changed in the different workflows.

Grouping PSMs

You can group PSMs into peptide groups by using the PSM Grouper node. For information on this node, refer to the Help.

Obtaining Protein Annotation Information

This chapter explains how the Proteome Discoverer application retrieves annotation information from ProteinCenter. Annotation information is metadata for investigated proteins from available public databases and includes widely used annotations like the GO (Gene Ontology) terms for each protein, Pfam (Protein Families) classification of proteins, Entrez Gene identifiers, Ensembl gene identifications from the Ensembl reference genomes of genetically sequenced organisms, and UniProt information about post-translational modifications (PTMs).

Contents

- [ProteinCenter](#)
- [Gene Ontology \(GO\) Database Annotation](#)
- [Pfam Database Annotation](#)
- [Entrez Gene Database Annotation](#)
- [Ensembl Genome Database Annotation](#)
- [UniProt Database Annotation](#)
- [Configuring the Proteome Discoverer Application for Protein Annotation](#)
- [Downloading Files from ProteinCenter](#)
- [Exporting Files Downloaded from ProteinCenter](#)
- [Creating a Protein Annotation Workflow](#)
- [Displaying the Annotated Protein Results](#)
- [Uploading Results to ProteinCenter](#)
- [Accessing the ProteinCard Page](#)
- [ProteinCard Page Parameters](#)
- [GO Slim Categories](#)

ProteinCenter

ProteinCenter is a web-based application that you can use to download biologically enriched annotation information for a single protein. You can download the following types of information:

- Molecular functions, cellular components, and biological processes from the GO database
- Classification information for protein families from the Pfam database
- Gene identifications from the Entrez Gene database
- Genomic annotations of genetically sequenced organisms from the Ensembl genome database
- Post-translational modification information from the UniProt database.

The data in ProteinCenter is updated biweekly.

You can also download FASTA databases from ProteinCenter (see [“Adding FASTA Files to the Proteome Discoverer Application”](#) on page 174). You can also upload search results directly from the Proteome Discoverer application to ProteinCenter.

The application gives you access to ProteinCenter in two ways:

- You can create a consensus workflow that includes the ProteinCenter Annotation node, which retrieves GO, Pfam, Entrez, Ensembl, and UniProt database information from ProteinCenter and stores it in the Proteome Discoverer results files. The application displays this information in columns on the Proteins page of the .pdResult file. For information on setting up a protein annotation workflow to achieve these results, see [“Configuring the Proteome Discoverer Application for Protein Annotation”](#) on page 293 and [“Creating a Protein Annotation Workflow”](#) on page 296.
- The ProteinCard page available for each protein displays the annotation data available in ProteinCenter and displays it on a page of the Protein Identification Details view (see [“Accessing the ProteinCard Page”](#) on page 326). You can display this information for the following proteins:
 - Proteins on the Proteins page of the .pdResult file
 - Proteins associated with identified peptides

You can access the ProteinCard page for each protein by double-clicking its row in the .pdResult report; or clicking its row, choosing View > Protein Details, and then clicking the ProteinCard tab of the Protein Identification Details view. The ProteinCard page itself is divided into separate pages representing different aspects of that protein: General, Keys, Features, Molecular Functions, Cellular Components, Biological Processes, Diseases, and External Links. You can display a ProteinCard page for every identified protein whose accession is tracked in ProteinCenter. For information on the ProteinCard page, see [“Accessing the ProteinCard Page”](#) on page 326 and [“ProteinCard Page Parameters”](#) on page 328.

You can download organism-specific sequence databases from ProteinCenter. For instructions, see “Adding FASTA Files to the Proteome Discoverer Application” on page 174.

You can also upload protein results directly from the application to ProteinCenter. For information, see “Uploading Results to ProteinCenter” on page 323.

Gene Ontology (GO) Database Annotation

The Gene Ontology (GO) database is a collaborative effort, incorporating community input from database and genome annotation groups to address the need for consistent descriptions of gene products in different databases. The GO project has developed three structured, controlled vocabularies (ontologies) that describe gene products in a species-independent way.

biological processes

cellular components

molecular functions

Each gene ontology is divided into categories and subcategories called GO terms, which define the protein in more specific terms. For example, chloroplast, a term in the cellular component ontology, is subdivided as follows:

chloroplast

[p] *chloroplast* envelope

[p] *chloroplast* membrane

[i] *chloroplast* inner membrane

[i] *chloroplast* outer membrane

You can find more information on the GO Ontology website at www.geneontology.org/.

GO Slim annotations are a collection of annotations assembled by ProteinCenter from various sources such as the European Bioinformatics Institute (EBI) and the GO consortium. These annotations are displayed in the Molecular Function, Cellular Component, and Biological Process columns of the results report (.pdResult) file. In addition, you can define your own categories of GO Slim annotations.

Pfam Database Annotation

In addition to GO and GO Slim annotations, you can also use ProteinCenter annotations from the Pfam database at the Wellcome Trust Sanger Institute (WTSI) (<http://pfam.xfam.org/>). These are annotations of protein families, which are proteins with similar sequences and similar biological functions. A special sequence comparison algorithm

called the Hidden Markov Model groups proteins into the families by comparing the sequences. Each family has its own identification number that starts with Pf The Proteins page of the .pdResult file displays this number in the Pfam IDs column. You can use the Pfam identification number to search in the Pfam database to obtain more details about the protein family.

The Pfam annotation system is an alternative to GO annotations. You might want to use the Pfam system to filter your proteins when you want the results to be traceable, scored, and uniformly grouped. You might also consider its computationally based data more reliable. However, it might be easier to use the hierarchy and grouping of the GO system to help you interpret results.

Table 12 compares the features of the GO and Pfam databases.

Table 12. Comparison of GO and Pfam features

GO features	Pfam features
Proteins grouped in biologically meaningful categories	Proteins grouped by similarity
Deep hierarchical order of terms	Few hierarchies
Data input by experts with different confidence levels and differing opinions	Computational data input with no human influence or expert knowledge

Entrez Gene Database Annotation

The database maintained by the National Center for Biotechnology Information (NCBI) includes the Entrez Gene database, which contains gene-specific information. Each gene stored in the Entrez Gene database has a unique identification. All proteins derived from the same gene have the same gene identification. When you use the Proteome Discoverer application to retrieve the Entrez Gene identifications from ProteinCenter, the Proteins page of the results report displays these identifications in the Entrez Gene ID column. You can use this information to group or cluster together the proteins that are biologically meaningful. Because not all genes are stored in the Entrez Gene database, some proteins do not have a valid gene identification. In this case, the application displays a blank or empty cell in the Entrez Gene ID column on the Proteins page of the results file.

Ensembl Genome Database Annotation

The Ensembl genome database is a joint project of the European Bioinformatics Institute (EBI) and the Wellcome Trust Sanger Institute (WTSI). It stores and automatically annotates the reference genomes of the genetically sequenced organisms. Because this database focuses on completed gene-sequencing projects, it contains fewer species than the large protein databases like the UniProt and NCBI databases, which also collect data from partially sequenced organisms.

UniProt Database Annotation


From ProteinCenter, you can retrieve information on known PTMs from the UniProt database and compare it with information on found PTMs. For details on this feature, see “[Viewing PTM Information on the Protein Identification Details View](#)” on [page 369](#).

Configuring the Proteome Discoverer Application for Protein Annotation

To retrieve protein annotations during report processing, downloading FASTA files, and using the ProteinCard page, you do not need to configure the Proteome Discoverer application. The application automatically provides the URL during installation.

However, to upload results to your local installation of ProteinCenter, you must configure the application in the Administration view and in the Options dialog box.

❖ To configure the application for uploading results to ProteinCenter in the Administration view


1. Choose **Administration > Configuration**, or click the **Edit Configuration** icon, .
2. Under Processing Settings in the Configuration area of the left pane, click **ProteinCenter** if it is not already selected.

The ProteinCenter view appears, as shown in [Figure 206](#).

8 Obtaining Protein Annotation Information

Configuring the Proteome Discoverer Application for Protein Annotation

Figure 206. ProteinCenter view




1. In the Upload URL box, specify the URL of the ProteinCenter location to upload a data set to.
2. In the User Name box, type the user name that allows you to access the specified upload server.

You were given this user name when you obtained a ProteinCenter license.


3. In the Password box, type the password that allows you to access the specified upload server.

You were given this password when you obtained a ProteinCenter license.

4. If you changed any settings, click the **Apply** icon,  **Apply** .

The message box shown in [Figure 207](#) appears:

Figure 207. Administration message box




5. Click **OK**.

❖ **To configure the application for uploading results to ProteinCenter in the Options dialog box**

1. Choose **Tools > Options**.
2. In the left pane of the Options dialog box, click **ProteinCenter**.

The ProteinCenter page of the dialog box opens, as shown in [Figure 208](#).

Figure 208. ProteinCenter page of the Options dialog box



- a. In the ProteinCenter URL box, type the URL of the ProteinCenter web server to use.
- b. In the Upload URL box, type the URL of the location to upload the files to.
- c. In the User Name box, type the user name of your ProteinCenter user account.
- d. In the Password box, type the password of your ProteinCenter user account.
- e. Click **OK**.

A message box appears with the following message:

Settings of Protein Center changed. Do you want to save your changes?

3. Click **Yes**.

For instructions on completing the upload to ProteinCenter, see [“Uploading Results to ProteinCenter”](#) on [page 323](#).

Downloading Files from ProteinCenter

To download files from ProteinCenter, see [“Adding FASTA Files to the Proteome Discoverer Application”](#) on page 174.

Exporting Files Downloaded from ProteinCenter

To export FASTA files that you downloaded from ProteinCenter, see [“Exporting FASTA Files”](#) on page 180.

Creating a Protein Annotation Workflow

You can retrieve annotations of all identified proteins from ProteinCenter by using the ProteinCenter Annotation node in a workflow. This node can retrieve the following information:

- Gene Ontology (GO) annotations, which are displayed in the GO Accessions column of the results report (.pdResult) file. These also include GO Slim annotations, which are displayed in the Molecular Function, Cellular Component, and Biological Process columns of the results report (.pdResult) file.

In addition, you can define your own categories of GO Slim annotations. In the Proteome Discoverer application, the categories containing user-defined GO Slim terms are called annotation aspects. For information on defining new annotation aspects, see [“Defining Annotation Aspects”](#) on page 304.

- Protein family (Pfam) annotations, which are displayed in the Pfam IDs column of the results report (.pdResult) file.
- Gene identifications from the Entrez Gene database, which are displayed in the Entrez Gene ID column of the results report (.pdResult) file.
- Gene identifications from the reference genomes of the genetically sequenced organisms in the Ensembl genome database, which are displayed in the Ensemble Gene ID and Chromosome columns of the Proteins page of the .pdResult file.
- UniProt PTM modifications documented in the UniProt database, which are displayed on the Protein Identification Details view in the Proteins page of the .pdResult file.

The data in ProteinCenter is updated biweekly.


❖ To create a processing protein annotation workflow

1. Follow the general instructions for creating a processing workflow with the Workflow Editor. See [“Creating a Processing Workflow”](#) on page 104.

The basic processing workflow for protein annotation is the same as the basic processing workflow shown in [Figure 74](#) on page 107.

2. Include the following nodes in the processing workflow:
 - The Spectrum Files node
 - The Spectrum Selector node
 - A search engine node (Mascot or Sequest HT)
 - The Fixed Value PSM Validator node, Percolator node, or Target Decoy PSM Validator node
3. Add any other appropriate nodes.
4. Connect the nodes together.
5. Set the parameters for each node.

For information about all the parameters that you can set for each node, refer to the Help.

6. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the processing workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the processing workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save**.
 - d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Processing Workflow File (*.pdProcessingWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdProcessingWF file.

❖ To create a consensus protein annotation workflow

1. If you are going to use the ProteinCenter Annotation node in the consensus workflow, configure this node by following the instructions in [“Configuring the Proteome Discoverer Application for Protein Annotation”](#) on page 293.
2. Create a consensus workflow by following the general instructions in [“Creating a Consensus Workflow”](#) on page 112.
3. Include the following nodes in the consensus workflow:
 - MSF Files node
 - PSM Grouper node
 - Peptide Validator node

8 Obtaining Protein Annotation Information


Creating a Protein Annotation Workflow


- Peptide and Protein Filter node
- ProteinScorer node
- ProteinCenter Annotation node
- Protein Grouping node

Connect the ProteinCenter Annotation node to the ProteinScorer node if the application does not automatically connect it.

Figure 209 shows the basic protein annotation consensus workflow.


Figure 209. Basic protein annotation consensus workflow



4. (Optional) Add any other appropriate nodes.
5. Set the parameters for each node.
6. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the consensus workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the consensus workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save** .

- d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Consensus Workflow File (*.pdConsensusWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdConsensusWF file.

7. Save the analysis. See “Saving an Analysis” on page 80.
8. Save the study. See “Saving a Study” on page 62.
9. Select the input files.
10. Click the **Run** icon, , in the upper right corner of the Analysis window.

The job queue appears, as shown in Figure 79 on page 112, displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

Displaying the Annotated Protein Results

The Proteome Discoverer application retrieves GO, Pfam, Entrez Gene, Ensembl, and UniProt PTM annotation data from ProteinCenter when it finishes processing all search nodes. You can display both predefined and user-defined annotated protein results in the results report (.pdResult file).

Note The application cannot retrieve annotations from searches conducted in the UniRef FASTA database because of the prefix appended to the accession number.

- [Displaying Predefined GO Protein Annotation Results](#)
- [Displaying Annotation Aspects](#)
- [Displaying GO Accessions](#)
- [Displaying Protein Family \(Pfam\) Annotations](#)
- [Displaying Entrez Gene Database Identifications](#)
- [Displaying Ensembl Genome Database Annotations](#)
- [Displaying UniProt Annotations](#)

Displaying Predefined GO Protein Annotation Results

Follow these procedures to display predefined GO protein categories, called annotation aspects, in the results report.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

To display user-defined annotation aspects, see “Displaying Annotation Aspects” on page 304.

❖ To display predefined GO protein annotation aspects in the results report


1. Generate the .pdResult file by following the instructions in “Creating a Protein Annotation Workflow” on page 296.
2. If the Molecular Function, Cellular Component, and Biological Processes columns are not visible on the Proteins page (they are visible by default), select them in the Field Chooser dialog box.

For information on the Field Chooser dialog box, refer to the Help.

The application displays the results on the Proteins page of the .pdResult file as colored boxes similar to those shown in ProteinCenter. Figure 210 gives an example.



















































Figure 210. GO Slim category boxes for the protein groups shown in the results of an annotation search

GO Slim category boxes




The application represents each aspect of the annotation (biological processes, cellular components, and molecular functions) in a separate column. Each box represents a GO Slim category, which is a selected subset of the Gene Ontology annotations. If the protein annotation is included in one of these subsets, the corresponding box is highlighted by a color specific to this GO Slim category. Figure 211 provides the column names and shows the meaning of the GO Slim category colors.

Figure 211. GO Slim category colors

Molecular Function		Cellular Component		Biological Process	
	antioxidant activity		cell surface		cell communication
	catalytic activity		chromosome		cell death
	DNA binding		cytoplasm		cell differentiation
	enzyme regulator activity		cytoskeleton		cell division
	metal ion binding		cytosol		cell growth
	motor activity		endoplasmic reticulum		cell organization and biogenesis
	nucleotide binding		endosome		cell proliferation
	protein binding		extracellular		cellular component movement
	receptor activity		golgi		cellular homeostasis
	RNA binding		membrane		coagulation
	signal transducer activity		mitochondrion		conjugation
	structural molecule activity		nucleus		defense response
	transcription regulator activity		organelle lumen		development
	translation regulator activity		proteasome		metabolic process
	transporter activity		ribosome		regulation of biological process
			spliceosome		reproduction
			vacuole		response to stimulus
					transport


If you move the mouse over the GO Slim category boxes, a ToolTip showing the category name appears, as shown in the Molecular Function column in [Figure 212](#).

Figure 212. ToolTip identifying the annotation category



Ready 40 Protein Groups; 50 Proteins; 127 Peptide Groups; 287 PSMs; 11805 MS/MS Spectrum Info

❖ **To filter the identified proteins by GO Slim categories**

1. In the .pdResult file, choose **View > Display Filter** or click the **Display Filter** icon, .
2. In the Display Filters pane, select **Proteins** in the left column.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

3. Click **Add Property**, and choose a property that represents an annotation aspect, such as **Molecular Function**, **Cellular Component**, or **Biological Process**.

The example in [Figure 213](#) shows the selection of Biological Processes.


Figure 213. Filtering annotation aspects

Checked	Master	Accession	Description	Contaminant	Standard	Decoy	cRAP	Coverage
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P10599	Thioredoxin					
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P00167	Cytochrome P450 2C9		X		X	46%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P00709	Alpha-lactalbumin		X		X	27%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P08758	Annexin A2		X		X	41%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P01112	GTPase		X		X	22%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P55957	BH3-interacting domain protein		X		X	22%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P02741	C-reactive protein		X		X	7%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P09211	Glutathione S-transferase		X		X	55%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P01133	Pro-epitope		X		X	1%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P12081	Histidine		X		X	20%
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P02753	Retinol-binding protein 4		X		X	24%

The application creates an expression with an “is true” operator.

4. In the blank target value box in the dialog box, select a predefined Go Slim value for the annotation, as shown in [Figure 214](#).


Figure 214. Selecting a predefined value for the annotation



5. Click the **Apply** icon, **Apply** .

Only those proteins that feature the selected aspect are now visible in the .pdResult file, as shown in Figure 215.

Figure 215. Proteins filtered by Go Slim category



You can combine several aspect filters by using AND or OR as with any other numerical or textual data type.

Displaying Annotation Aspects

Although 50 GO Slim terms are available in the Molecular Function, Cellular Component, and Biological Processes categories, these annotation categories, called annotation aspects, might not be specific enough for your needs. You might want to define your own aspects to display in the results report. You can define these aspects in the Administration view.


To display predefined aspects such as Molecular Function, see “[Displaying Predefined GO Protein Annotation Results](#)” on page 299.

- [Defining Annotation Aspects](#)
- [Using User-Defined Annotation Aspects](#)
- [Editing User-Defined Annotation Aspects](#)
- [Removing User-Defined Annotation Aspects](#)
- [Exporting User-Defined Annotation Aspects](#)
- [Importing User-Defined Annotation Aspects](#)
- [Defining New Annotation Aspects for the First Time After Installation](#)
- [Annotation Aspect Editor Parameters](#)
- [Annotation Aspects View Parameters](#)


Defining Annotation Aspects

The Annotation Aspects view, shown in [Figure 216](#), displays the currently available annotation aspects. When you click + to the left of an aspect, the application displays the GO Slim terms, or annotation groups, associated with that aspect. It also displays the color code used in the results report for each group and the group name used in the ToolTip that appears when you move the mouse over the GO Slim category boxes, as shown in [Figure 212](#) on page 301. [Figure 216](#) shows these groups for the Biological Processes category.

Figure 216. Annotation Aspects view



❖ **To define an annotation aspect**


1. Choose **Administration > Maintain Annotation Aspects**.
2. In the Annotation Aspect view on the Administration page, click the **Add** icon,  **Add**, in the toolbar.

The Annotation Aspect Editor opens and loads the hierarchical definitions of the GO terms and the Pfam terms in the Available Annotations area in the lower right part of the Editor, as shown in [Figure 217](#).

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

Figure 217. Annotation Aspect Editor



The Annotation Aspects Editor consists of three parts:

- On the top is the aspect definition area, which includes the name and description of the annotation aspect and colored boxes representing the annotation groups of the aspect.
 - On the lower left is the Annotation Group area, which contains the definitions of a single annotation group.
 - On the lower right is the Available Annotations area, which lists the available annotations.
3. In the Aspect Name box in the aspect definition area, enter the aspect name, for example, **Metal_Ion_Binding**.

You can use letters, numbers, and underscores (_) in the aspect name. Spaces are not allowed.


4. In the Description box, type a brief definition of the new aspect, for example, “Chemical element bound to the protein.”

5. Click **New Group** in the Annotation Group area of the Editor.

A red box appears in the aspect definition area of the Editor beneath the Description box.

6. In the Name box in the Annotation Group area, enter the name of the group, for example, “Magnesium” (see [Figure 218](#)).

7. Specify the fill color of the group rectangle used if a protein belongs to the new group (see [Figure 218](#)):

- a. Click the button  next to Color to open the Color dialog box.
- b. Click the color that you want, or define a custom color. For information on defining a custom color, refer to the Help.
- c. Click **OK**.

8. (Optional) If you specify more than one group, change the position of a rectangle in the aspect definition area.

Starting from the left, the rectangles in the aspect definition area are assigned a sequential position number in the Position box. To change the position of a rectangle in the aspect definition area, select the rectangle and enter the new position number into the Position box. You can also use the up and down arrows in the Position box to set the new position number. For example, if the rectangle corresponding to potassium is the fourth from the left, and you want it to be the seventh rectangle from the left, enter 7 in the Position box.

The Available Annotations area lists all GO annotations with the accession number and a description. When you click + to the left of a GO annotation, the Editor displays all child terms, if any, for the annotation.

9. In the Annotation Database box, select the database of annotation to use, either **Gene Ontology** or **Protein Families (Pfam)**.

For information on the Gene Ontology (GO) annotation, see “[Gene Ontology \(GO\) Database Annotation](#)” on [page 291](#). For information on Pfam annotation, see “[Pfam Database Annotation](#)” on [page 291](#).


10. Scroll down the list of GO or Pfam terms to find the appropriate term, in this case, “metal ion binding.” For faster access, use the column filters in the gray row below the column headers. Type the name you are looking for in the filter cell. The entries in the table are automatically restricted to values containing the given name.


[Figure 218](#) shows a partial list of the available annotations and the child terms for metal ion binding.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

Figure 218. List of child terms for metal ion binding



11. To associate GO terms with the group, select the terms to include in the list of available annotations, and click the left arrow button, , to move them to the Included Annotations pane.

All GO terms are listed with an accession number and a description. The proteins in the results report (.pdResult file) reflect all GO or Pfam terms, including the appended child terms, that you move to the group definition. For example, when you define a group with 17 different GO or Pfam terms, the results report displays a colored rectangle if the corresponding protein is annotated with one of these 17 terms. Otherwise, the rectangle remains empty.

Figure 219 shows the complete definition of an aspect containing one group, which is highlighted for all proteins that are annotated with “GO:0000287- magnesium ion binding.” With this annotation aspect, you can filter the results report (.pdResult file) for all proteins containing magnesium.

Figure 219. Annotation aspect with completed definition for one group





Figure 220 shows the metal ion binding aspect with more group definitions. You can define a maximum of 25 groups per aspect.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

Figure 220. Annotation aspect with eleven groups




To remove a group, select the colored box corresponding to the group and click **Delete Group** in the Annotation Group area.

❖ To save user-defined aspect annotations

Click **Save** in the lower right corner of the Annotation Aspect Editor.

Before the application saves an aspect, it verifies that the aspect name and the group names are defined and valid. It also verifies that you set the color and the list of group definitions. If it finds missing mandatory values, it opens a message box listing the missing values. You must set all mandatory values before the application saves the annotation aspect. After the application saves the annotation aspect, it lists the aspect in the Administration view, as shown in [Figure 221](#). Click the plus sign, +, to the left of the row to display the single groups of the aspect, as shown in [Figure 221](#).



Figure 221. User-defined aspect in the Annotation Aspects view



Using User-Defined Annotation Aspects

By default, a user-defined aspect is active in the Annotation Aspects view and therefore usable in the ProteinCenter Annotation node.

❖ To deactivate a user-defined annotation aspect


1. In the table in the Annotation Aspects view, click the check box in the Active column for the appropriate annotation aspect.
2. Click the **Apply** icon,  **Apply**.
3. To make an annotation aspect active, select its check box in the Active column, and click the **Apply** icon,  **Apply**.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

You can select user-defined aspects in the ProteinCenter Annotation node in a consensus workflow. As shown in [Figure 222](#), the ProteinCenter Annotation node has up to six annotation aspects available for the Annotation Aspects parameter that you can apply to the annotation. When you click on each numbered aspect, a dropdown list appears from which you can choose any predefined and user-defined aspect. The list contains all aspect marked “Active” in the Administration view.

Figure 222. User-defined annotation aspect in the Parameters pane of the Workflow Editor



The application uses the new annotation aspects in its consensus workflow and displays the results of the annotation in the .pdResult file, as shown in [Figure 223](#).

Figure 223. User-defined annotation aspects in the results report

#	Checked	Accession	Description	I Coverage	F PSMs	# AAs	MW [kDa]	calc: pI	Is Master Protein	Biological Process	Cellular Component	Molecular Function	Metal_Ion_Binding	Pfam IDs	GO Accessions	Entrez Gene ID
1	<input type="checkbox"/>	P12081	Histidine-tRNA ligase, cytoplasmic OS=Homo sapiens GN=	15.72 %	30	509	57.4	5.88	✓					PF00458, PF00587, PF0	(GO:000166); GO:00: 3035	
2	<input type="checkbox"/>	P01008	Antithrombin-III OS=Homo sapiens GN=SERPINC1 PE=1 S	12.07 %	25	464	52.6	6.71	✓					PF00079	(GO:0002020); GO:00: 462	
3	<input type="checkbox"/>	P02768	Serum albumin OS=Homo sapiens GN=ALB PE=1 SV=2	12.81 %	17	609	69.3	6.28	✓					PF00273	(GO:0002576); GO:00: 213	
4	<input type="checkbox"/>	P04040	Catalase OS=Homo sapiens GN=CAT PE=1 SV=3	12.52 %	13	527	59.7	7.39	✓					PF01199, PF06628	(GO:0003002); GO:00: 847	
5	<input type="checkbox"/>	P02787	Serotransferrin OS=Homo sapiens GN=TF PE=1 SV=3	10.46 %	16	698	77.0	7.12	✓					PF00405, PF12974	(GO:0002576); GO:00: 7018	
6	<input type="checkbox"/>	P01127	Platelet-derived growth factor subunit B OS=Homo sapiens	17.01 %	12	241	27.3	9.16	✓					PF00341, PF04692	(GO:0001139); GO:00: 5155	
7	<input type="checkbox"/>	P55957	BH3-interacting domain death agonist OS=Homo sapiens G	13.33 %	15	195	22.0	5.44	✓					PF06393	(GO:0001836); GO:00: 637	
8	<input type="checkbox"/>	P08758	Annexin A5 OS=Homo sapiens GN=ANXA5 PE=1 SV=2	15.94 %	12	320	35.9	5.05	✓					PF00191	(GO:0004859); GO:00: 308	
9	<input type="checkbox"/>	Q06830	Peroxisomal oxidase OS=Homo sapiens GN=PRDX1 PE=1 SV+	7.64 %	10	199	22.1	8.13	✓					PF00578, PF08534, PF1	(GO:0003002); GO:00: 5052	
10	<input type="checkbox"/>	P06396	Gelsolin OS=Homo sapiens GN=GSN PE=1 SV=1	3.84 %	11	782	85.6	6.28	✓					PF00626	(GO:0001726); GO:00: 2934	
11	<input type="checkbox"/>	P05413	Fatty acid-binding protein, heart OS=Homo sapiens GN=FA	3.08 %	9	133	14.8	6.80	✓					PF00061	(GO:0005251); GO:00: 2170	
12	<input type="checkbox"/>	P00167	Cytochrome b5 OS=Homo sapiens GN=CYB5A PE=1 SV=2	8.1 %	7	134	15.3	4.96	✓					PF00173	(GO:0004033); GO:00: 1528	
13	<input type="checkbox"/>	P16083	Ribosylhydroxymethyltransferase [quinone] OS=H	19.05 %	9	231	25.9	6.29	✓					PF02525, PF03358	(GO:0001512); GO:00: 4835	
14	<input type="checkbox"/>	P69905	Hemoglobin subunit alpha OS=Homo sapiens GN=HBAT P1	10.56 %	9	142	15.2	8.68	✓					PF00042	(GO:0004601); GO:00: 3040	
15	<input type="checkbox"/>	P00709	Alpha-lactalbumin OS=Homo sapiens GN=LALBA PE=1 SV	15.49 %	6	142	16.2	5.00	✓					PF00062	(GO:0004461); GO:00: 3906	
16	<input type="checkbox"/>	P10543	Nucleoside diphosphate kinase 2 OS=Homo sapiens GN=	7.28 %	4	81	9.1	8.43	✓					PF02240, PF11976	(GO:0005515); GO:00: 4738	
17	<input type="checkbox"/>	P02768	Lactotransferrin OS=Homo sapiens GN=LTF PE=1 SV=6	4.93 %	4	710	78.1	8.12	✓					PF00405	(GO:0004252); GO:00: 4057	
18	<input type="checkbox"/>	P06732	Creatine kinase M-type OS=Homo sapiens GN=CKM PE=1	6.56 %	5	381	43.1	7.25	✓					PF00217, PF02807	(GO:0001166); GO:00: 1158	
19	<input type="checkbox"/>	P02144	Myoglobin OS=Homo sapiens GN=MB PE=1 SV=2	2.73 %	5	154	17.2	7.68	✓					PF00042	(GO:0001666); GO:00: 4151	
20	<input type="checkbox"/>	P68871	Hemoglobin subunit beta OS=Homo sapiens GN=HBB PE-	15.65 %	4	147	16.0	7.28	✓					PF00042	(GO:0004601); GO:00: 3043	
21	<input type="checkbox"/>	P00915	Carbonic anhydrase 1 OS=Homo sapiens GN=CA1 PE=1 S	13.79 %	4	261	28.9	7.12	✓					PF00194	(GO:0004089); GO:00: 759	
22	<input type="checkbox"/>	P02769	Serum albumin OS=Bos taurus GN=ALB PE=1 SV=4	3.62 %	4	607	69.2	6.18	✓					PF00273	(GO:0003677); GO:00: 280717	
23	<input type="checkbox"/>	O00762	Ubiquitin-conjugating enzyme E2 C OS=Homo sapiens GN-	13.97 %	4	179	19.6	7.37	✓					PF00179, PF05743	(GO:0000075); GO:00: 11065	
24	<input type="checkbox"/>	P01112	GTase HRas OS=Homo sapiens GN=HRAS PE=1 SV=1	11.11 %	4	189	21.3	5.31	✓					PF00009, PF00071, PF0	(GO:0001139); GO:00: 3265	
25	<input type="checkbox"/>	P02741	C-reactive protein OS=Homo sapiens GN=CRP PE=1 SV=1	7.14 %	4	224	25.0	5.63	✓					PF00354	(GO:0001666); GO:00: 1401	
26	<input type="checkbox"/>	P01344	Insulin-like growth factor II OS=Homo sapiens GN=IGF2 PE	10.56 %	4	180	20.1	9.32	✓					PF00049, PF08365	(GO:0001501); GO:00: 3481	
27	<input type="checkbox"/>	P08263	Glutathione S-transferase A1 OS=Homo sapiens GN=GSTA	6.76 %	4	222	25.6	8.88	✓					PF00043, PF02798	(GO:0003002); GO:00: 2938	
28	<input type="checkbox"/>	P09211	Glutathione S-transferase P OS=Homo sapiens GN=GSTP1	5.24 %	3	210	23.3	5.64	✓					PF00043, PF02798	(GO:0003002); GO:00: 2950	
29	<input type="checkbox"/>	P10543	Microtubule-associated protein tau OS=Homo sapiens GN-	3.30 %	3	758	78.9	6.71	✓					PF02418	(GO:0002276); GO:00: 4137	
30	<input type="checkbox"/>	R00432	Catalase OS=Bos taurus GN=CAT PE=1 SV=3	3.81 %	3	527	59.9	7.28	✓					PF01199, PF06628	(GO:0003002); GO:00: 531682	
31	<input type="checkbox"/>	P01031	Complement C5 OS=Homo sapiens GN=C5 PE=1 SV=4	0.78 %	3	1676	188.2	6.52	✓					PF00207, PF01759, PF0	(GO:000187); GO:00: 727	
32	<input type="checkbox"/>	P05019	Insulin-like growth factor I OS=Homo sapiens GN=IGF1 PE	7.69 %	2	195	21.8	9.72	✓					PF00049	(GO:0001501); GO:00: 3479	
33	<input type="checkbox"/>	P62987	Ubiquitin-60S ribosomal protein L40 OS=Homo sapiens GN	14.06 %	3	128	14.7	8.83	✓					PF00240, PF01020, PF1	(GO:0000075); GO:00: 7311	
34	<input type="checkbox"/>	P0CG47	Polyubiquitin-B OS=Homo sapiens GN=UBB PE=1 SV=1	3.58 %	3	229	25.7	7.43	✓					PF00240, PF11976	(GO:0000075); GO:00: 7314	
35	<input type="checkbox"/>	P62979	Ubiquitin-40S ribosomal protein S27a OS=Homo sapiens G	11.54 %	3	156	18.0	9.64	✓					PF00240, PF01599, PF1	(GO:0000075); GO:00: 6233	
36	<input type="checkbox"/>	P0CG48	Polyubiquitin-C OS=Homo sapiens GN=UBC PE=1 SV=3	3.65 %	3	685	77.0	7.66	✓					PF00240, PF11976	(GO:0000075); GO:00: 7316	
37	<input type="checkbox"/>	P62937	Peptidyl-prolyl cis-trans isomerase A OS=Homo sapiens GA	7.27 %	2	165	18.0	7.81	✓					PF00160	(GO:0000413); GO:00: 5478	
38	<input type="checkbox"/>	P63165	Small ubiquitin-related modifier 1 OS=Homo sapiens GN-S	6.93 %	2	101	11.5	5.52	✓					PF00240, PF11976	(GO:0005515); GO:00: 7341	
39	<input type="checkbox"/>	P00441	Superoxide dismutase [Cu-Zn] OS=Homo sapiens GN=SOC	8.44 %	3	154	15.9	6.13	✓					PF00080	(GO:000187); GO:00: 6647	
40	<input type="checkbox"/>	P02062	Hemoglobin subunit beta Equus caballus GN=HBB PE-	6.85 %	2	146	16.0	7.02	✓					PF00042	(GO:0005344); GO:00: 100054109	
41	<input type="checkbox"/>	P02042	Hemoglobin subunit delta OS=Homo sapiens GN=HBD PE-	6.80 %	2	147	16.0	8.05	✓					PF00042	(GO:0005344); GO:00: 3045	
42	<input type="checkbox"/>	P00706	Lyszyme C-3 OS=Anas platyrhynchos PE=1 SV=1	4.65 %	2	129	14.5	8.82	✓					PF00062	(GO:0003796); GO:00: 0	
43	<input type="checkbox"/>	P61626	Lyszyme C OS=Homo sapiens GN=LYZ PE=1 SV=1	4.05 %	2	146	16.5	9.16	✓					PF00062	(GO:0003796); GO:00: 4069	
44	<input type="checkbox"/>	P00268	Lyszyme C OS=Callinectes sapidus GN=Lyz PE=1 SV=1	4.08 %	2	147	16.3	9.07	✓					PF00062	(GO:0003796); GO:00: 982318	

Editing User-Defined Annotation Aspects


You can change an existing user-defined annotation aspect.

❖ To edit user-defined annotation aspects

1. In the table in the Annotation Aspects view on the Administration page, double-click the annotation aspect that you want to edit.

—or—

Do the following:

- a. Select the row containing the user-defined annotation aspect that you want to edit. Click the row in the column to the left of the ID column to select the aspect.
- b. Click the **Edit** icon,  **Edit**, in the toolbar.

The Annotation Aspects Editor opens.

If you double-click a predefined aspect like Biological Process, the Annotation Aspect Editor opens, but you cannot change any of the settings. You can only edit the settings for user-defined annotation aspects.

2. Make the changes that you want, and click **Save** in the lower right-hand corner of the Annotation Aspects Editor.

8 Obtaining Protein Annotation Information


Displaying the Annotated Protein Results

The Annotation Aspect Editor closes, leaving the Annotation Aspects view displayed on the Administration page.

Removing User-Defined Annotation Aspects

You can remove a user-defined annotation aspect.

❖ To remove user-defined annotation aspects


1. In the table on the Annotation Aspects view on the Administration page, select the row containing the user-defined annotation aspect that you want to remove. Click the row in the column to the left of the ID column to select the aspect.
2. Click the **Remove** icon, , on the toolbar.

The annotation aspect no longer appears in the table in the Annotation Aspects view.

Exporting User-Defined Annotation Aspects

You can export user-defined annotation aspects to an XML file. You might want to use this XML file to share user-defined annotation aspects between different installations of the Proteome Discoverer application, or use it as a backup if you install the application on a new computer.


❖ To export a user-defined annotation aspect

1. In the table on the Annotation Aspects view on the Administration page, select the row containing the user-defined annotation aspect that you want to export. Click the row in the column to the left of the ID column to select the aspect.
2. Click the **Export** icon, , in the toolbar.
3. In the Export Annotation Aspect dialog box, enter the name of the XML file where you want to save the annotation aspects.
4. Click **Save**.

Importing User-Defined Annotation Aspects

You can import user-defined annotation aspects from an XML file.

❖ To import a user-defined annotation aspect

1. In the table in the Annotation Aspects view of the Administration page, select the row containing the user-defined annotation aspect that you want to import. Click the row in the column to the left of the ID column to select the aspect.
2. Click the **Import** icon, , in the toolbar.


3. In the Import Annotation Aspect dialog box, select the XML file containing the annotation aspects that you want to import.
4. Click **Open**.

The imported annotation aspect now appears in the table in the Annotation Aspects view.

Defining New Annotation Aspects for the First Time After Installation


After you install the Proteome Discoverer application, only an initial version of the annotation database is available, and it does not contain any annotation terms. You must have an updated annotation database available before you can define new aspects. The application updates the annotation database during the first use of the ProteinCenter Annotation node in a consensus workflow. If you define new aspects before using the ProteinCenter Annotation node for the first time or after a long period of no usage in a consensus workflow (in this case, the local annotation database would be outdated), manually update the annotation database by following this procedure.

❖ To update the annotation aspect database for the first time after installation

Click the **UpdateDB** icon,  **UpdateDB**, in the Annotation Aspects view.

The view changes to the job queue view, as shown in [Figure 224](#).

Figure 224. Job queue showing update of annotation aspect database



Annotation Aspect Editor Parameters

[Table 13](#) describes the parameters available in the Annotation Aspect Editor.

Table 13. Annotation Aspect Editor parameters (Sheet 1 of 2)

Parameter	Description
Aspect Name	Specifies the name of the annotation aspect.
Description	Briefly describes the annotation aspect.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results






Table 13. Annotation Aspect Editor parameters (Sheet 2 of 2)

Parameter	Description
Aspect definition area	Displays colored rectangles corresponding to the GO or Pfam terms (called groups) assigned to the annotation aspect.
Annotation Group	Contains parameters that define a group.
Name	Specifies the name of the group.
New Group	Places a rectangle in the aspect definition area to represent the new group being defined.
Position	Specifies the position of a rectangle representing a group in the aspect definition area.
Color	Specifies the color that the group displays in the results report (.pdResult file).
Delete Group	Deletes a group.
Included Annotations	Lists the annotations assigned to a group.
Accession	Displays the accession number of the annotation.
Description	Briefly describes the annotation.
# Child Terms	Displays the number of child terms that the annotation has.
<	Assigns a GO or Pfam term to an annotation group.
>	Removes a GO or Pfam term from an annotation group.
Available Annotations	Lists the annotations in the selected database.
Annotation Database	Selects the database of annotations to display.
Accession	Displays the accession number of the annotation.
Description	Briefly describes the annotation.
# Child Terms	Displays the number of child terms that the annotation has.
Cancel	Discards any changes to the annotation aspect and closes the Annotation Aspect Editor.
Saves	Saves the annotation aspect.

Annotation Aspects View Parameters

Table 14 describes the parameters available in the Annotation Aspects view on the Administration page.

Table 14. Annotation Aspects view parameters

Parameter	Description
 Apply	Implements any changes that you make to the ID, Read Only, or Active columns.
 Add	Opens the Annotation Aspects Editor so that you can define a new annotation aspect.
 Remove	Removes an annotation aspect from the table on the Annotation Aspects view.
 Edit	Opens the Annotation Aspects Editor so that you can change the settings of an existing annotation aspect.
Import	Imports an XML file containing annotation aspects and displays it in the table on the Annotation Aspects view.
Export	Places the selected annotation aspects in an XML file.
 UpdateDB	Updates an annotation database before you define a new annotation aspect for the first time after installation of the Proteome Discoverer application.
Active	Determines whether the annotation aspect is displayed in the ProteinCenter Annotation node.
Aspect Name	Displays the name of the annotation aspects.
Position	Specifies the position of the colored rectangle representing a group assigned to the annotation aspect.
Color	Specifies the color that a group of an annotation aspect displays in the results report (.pdResult file).
Description	Briefly describes the group of an annotation aspect.

Displaying GO Accessions

Gene ontology terms are related in hierarchical graphs. The GO term annotated to a special protein is always part of a complex directed graph. All ancestor elements—that is, the elements between the annotated GO term and one of the three top-level terms (Molecular Function, Cellular Component, and Biological Process)—are additional less-specific

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

descriptions of the annotated value. For example, the “iron ion binding (GO:0005506)” term contains in its graph the “metal ion binding (GO:0046872)” value, which is less specific. All GO terms contained in the graph of the annotated GO term of the protein are represented in the GO Accessions column on the Proteins page.

❖ To display GO accessions

1. Open the generated .pdResult file by following the instructions in the Help.
2. Click the **Proteins** tab.
3. In the Field Chooser dialog box of the proteins page, select the **GO Accessions** column.

For information on the Field Chooser dialog box, refer to the Help.

The application displays the protein’s GO terms contained in the graph of the annotated GO term on the Proteins page of the .pdResult report in the GO Accessions column, as shown in [Figure 225](#).

Figure 225. Go Accessions column on the Proteins page of the .pdResult file

Checked	Accession	Description	Is Master Protein	Biological Process	Cellular Component	Molecular Function	Pfam IDs	GO Accessions	Entrez Gene ID
<input type="checkbox"/>	P00167	Cytochrome b5 OS=Homo sapiens GN=CYB5A F	✓				PF00173	[GO:0004033]; GO:001	1528
<input type="checkbox"/>	P05413	Fatty acid-binding protein, heart OS=Homo sapie	✓				PF00061	[GO:0005215]; GO:001	2170
<input type="checkbox"/>	Q06830	Peroxioredoxin-1 OS=Homo sapiens GN=PRDX1	✓				PF00578, PF08534, PF	[GO:0000302]; GO:001	5052
<input type="checkbox"/>	P0CG48	Polyubiquitin-C OS=Homo sapiens GN=UBC PE=	✓				PF00240, PF11976	[GO:0000075]; GO:001	7316
<input type="checkbox"/>	P0CG47	Polyubiquitin-B OS=Homo sapiens GN=UBB PE=	✓				PF00240, PF11976	[GO:0000075]; GO:001	7314
<input type="checkbox"/>	P02144	Myoglobin OS=Homo sapiens GN=MB PE=1 SV=	✓				PF00042	[GO:0001666]; GO:001	4151
<input type="checkbox"/>	P16083	Ribosylidihydroxynicotinamide dehydrogenase [quin	✓				PF02525, PF03358	[GO:0001512]; GO:001	4835
<input type="checkbox"/>	Q15843	NEDD8 OS=Homo sapiens GN=NEDD8 PE=1 S'	✓				PF00240, PF11976	[GO:0005515]; GO:001	4738
<input type="checkbox"/>	P01127	Platelet-derived growth factor subunit B OS=Homo	✓				PF00341, PF04692	[GO:0000139]; GO:001	5155
<input type="checkbox"/>	P10145	Interleukin-8 OS=Homo sapiens GN=IL8 PE=1 S	✓				PF00048	[GO:0001525]; GO:001	3576
<input type="checkbox"/>	P08758	Annexin A5 OS=Homo sapiens GN=ANXA5 PE=	✓				PF00191	[GO:0004859]; GO:001	308
<input type="checkbox"/>	P12081	Histidine-tRNA ligase, cytoplasmic OS=Homo sa	✓				PF00458, PF00587, PF	[GO:0000166]; GO:001	3035
<input type="checkbox"/>	P68871	Hemoglobin subunit beta OS=Homo sapiens GN=	✓				PF00042	[GO:0004601]; GO:001	3043
<input type="checkbox"/>	P00709	Alpha-lactalbumin OS=Homo sapiens GN=LALB	✓				PF00062	[GO:0004461]; GO:001	3906
<input type="checkbox"/>	P62987	Ubiquitin-60S ribosomal protein L40 OS=Homo s	✓				PF00240, PF01020, PF	[GO:0000075]; GO:001	7311
<input type="checkbox"/>	O00762	Ubiquitin-conjugating enzyme E2 C OS=Homo sa	✓				PF00179, PF05743	[GO:0000075]; GO:001	11065
<input type="checkbox"/>	P00915	Carbonic anhydrase 1 OS=Homo sapiens GN=C.	✓				PF00194	[GO:0004089]; GO:001	759
<input type="checkbox"/>	P55957	BH3-interacting domain death agonist OS=Homo	✓				PF06393	[GO:0001836]; GO:001	637
<input type="checkbox"/>	P02768	Serum albumin OS=Homo sapiens GN=ALB PE=	✓				PF00273	[GO:0002576]; GO:001	213
<input type="checkbox"/>	P04040	Catalase OS=Homo sapiens GN=CAT PE=1 SV=	✓				PF00199, PF06628	[GO:0000302]; GO:001	847
<input type="checkbox"/>	P01008	Antithrombin-III OS=Homo sapiens GN=SERPIN	✓				PF00079	[GO:0002020]; GO:001	462
<input type="checkbox"/>	P62979	Ubiquitin-40S ribosomal protein S27a OS=Homo	✓				PF00240, PF01599, PF	[GO:0000075]; GO:001	6233
<input type="checkbox"/>	P01112	GTPase HRas OS=Homo sapiens GN=HRAS PE	✓				PF00009, PF00071, PF	[GO:0000139]; GO:001	3265
<input type="checkbox"/>	P61769	Beta-2-microglobulin OS=Homo sapiens GN=B2I	✓				PF07654	[GO:0000139]; GO:001	567
<input type="checkbox"/>	P69905	Hemoglobin subunit alpha OS=Homo sapiens GN=	✓				PF00042	[GO:0004601]; GO:001	3040
<input type="checkbox"/>	P01344	Insulin-like growth factor II OS=Homo sapiens GN=	✓				PF00049, PF08365	[GO:0001501]; GO:001	3481
<input type="checkbox"/>	P10599	Thioredoxin OS=Homo sapiens GN=TXN PE=1 S	✓				PF00085	[GO:0000122]; GO:001	7295
<input type="checkbox"/>	P02787	Serotransferrin OS=Homo sapiens GN=TF PE=1	✓				PF00405, PF12974	[GO:0002576]; GO:001	7018
<input type="checkbox"/>	P00441	Superoxide dismutase [Cu-Zn] OS=Homo sapien	✓				PF00080	[GO:0000187]; GO:001	6647
<input type="checkbox"/>	P05019	Insulin-like growth factor I OS=Homo sapiens GN=	✓				PF00049	[GO:0001501]; GO:001	3479
<input type="checkbox"/>	P62937	Pantothelinol ribitolane isomerase A OS=Homo s	✓				PF00160	[GO:0000113]; GO:001	5478

The square parentheses around a GO term indicate that the term is a child term of the GO terms without parentheses on the same line. For example, on the seventh line in [Figure 226](#), [GO:0001512] is the child term of the higher-level terms GO:0016491, GO:0003824, and GO:0003674. The latter three terms represent a higher level in the biological context of the protein.

Figure 226. Child and higher-level GO terms

Accession	Description	Biological Process	Cellular Component	Molecular Function	GO Accessions
P00167	Cytochrome b5 OS=Homo sapiens GN=CYP11B	[GO:0004033]	[GO:0004033]	[GO:0004033]	[GO:0004033]; [GO:0016616]; [GO:0016614]; [GO:0016491]; [GO:0003824]; [GO:0003674]
P05413	Fatty acid-binding protein, heart OS=Homo sapiens GN=ABHD11	[GO:0005215]	[GO:0005215]	[GO:0005215]	[GO:0005215]; [GO:0003674]
Q06830	Peroxisomal protein 1 OS=Homo sapiens GN=PRX1	[GO:0000302]	[GO:0000302]	[GO:0000302]	[GO:0000302]; [GO:0006979]; [GO:0010035]; [GO:0006950]; [GO:0042221]; [GO:0050896]; [GO:0008150]
P0CG48	Polyubiquitin-C OS=Homo sapiens GN=UBI1	[GO:0000075]	[GO:0000075]	[GO:0000075]	[GO:0000075]; [GO:0007050]; [GO:0071156]; [GO:0022402]; [GO:0045786]; [GO:0010564]; [GO:0007049]; [GO:000099]
P0CG47	Polyubiquitin-B OS=Homo sapiens GN=UBI1	[GO:0000075]	[GO:0000075]	[GO:0000075]	[GO:0000075]; [GO:0007050]; [GO:0071156]; [GO:0022402]; [GO:0045786]; [GO:0010564]; [GO:0007049]; [GO:000099]
P02144	Myoglobin OS=Homo sapiens GN=MB PE=	[GO:0001666]	[GO:0001666]	[GO:0001666]	[GO:0001666]; [GO:0006950]; [GO:0036293]; [GO:0050896]; [GO:0070482]; [GO:0008150]; [GO:0009628]
P16083	Ribosylidihydroxynicotinamide dehydrogenase	[GO:0001512]	[GO:0001512]	[GO:0001512]	[GO:0001512]; [GO:0016491]; [GO:0003824]; [GO:0003674]
Q15843	NEDD8 OS=Homo sapiens GN=NEDD8 PE=	[GO:0005515]	[GO:0005515]	[GO:0005515]	[GO:0005515]; [GO:0005488]; [GO:0003674]
P01127	Platelet-derived growth factor subunit B OS=	[GO:0000139]	[GO:0000139]	[GO:0000139]	[GO:0000139]; [GO:0012505]; [GO:0031090]; [GO:0044431]; [GO:0044464]; [GO:0016020]; [GO:0043227]; [GO:00444]
P10145	Interleukin-8 OS=Homo sapiens GN=IL8 PE=	[GO:0001525]	[GO:0001525]	[GO:0001525]	[GO:0001525]; [GO:0048514]; [GO:0048646]; [GO:0001568]; [GO:0009653]; [GO:0032502]; [GO:0001944]; [GO:00488]
P08758	Annexin A5 OS=Homo sapiens GN=ANXA5	[GO:0004859]	[GO:0004859]	[GO:0004859]	[GO:0004859]; [GO:0055102]; [GO:0004857]; [GO:0030234]; [GO:0003674]
P12081	Histidine-tRNA ligase, cytoplasmic OS=Homo sapiens GN=HIST1	[GO:0000166]	[GO:0000166]	[GO:0000166]	[GO:0000166]; [GO:0036094]; [GO:1901265]; [GO:0005488]; [GO:0097159]; [GO:0003674]
P68871	Hemoglobin subunit beta OS=Homo sapiens GN=HBB	[GO:0004601]	[GO:0004601]	[GO:0004601]	[GO:0004601]; [GO:0016209]; [GO:0016684]; [GO:0003674]; [GO:0016491]; [GO:0003824]
P00709	Alpha-lactalbumin OS=Homo sapiens GN=LALBA	[GO:0004461]	[GO:0004461]	[GO:0004461]	[GO:0004461]; [GO:0035250]; [GO:0008194]; [GO:0008378]; [GO:0016757]; [GO:0016758]; [GO:0016740]; [GO:00038]
P62987	Ubiquitin-60S ribosomal protein L40 OS=Homo sapiens GN=UBI1	[GO:0000075]	[GO:0000075]	[GO:0000075]	[GO:0000075]; [GO:0007050]; [GO:0071156]; [GO:0022402]; [GO:0045786]; [GO:0010564]; [GO:0007049]; [GO:000099]
O00762	Ubiquitin-conjugating enzyme E2 C OS=Homo sapiens GN=UBE2C	[GO:0000075]	[GO:0000075]	[GO:0000075]	[GO:0000075]; [GO:0007050]; [GO:0071156]; [GO:0022402]; [GO:0045786]; [GO:0010564]; [GO:0007049]; [GO:000099]
P00915	Carbonic anhydrase 1 OS=Homo sapiens GN=CA1	[GO:0004089]	[GO:0004089]	[GO:0004089]	[GO:0004089]; [GO:0016836]; [GO:0016835]; [GO:0016829]; [GO:0003824]; [GO:0003674]
P55957	BH3-interacting domain death agonist OS=Homo sapiens GN=BIH1	[GO:0001836]	[GO:0001836]	[GO:0001836]	[GO:0001836]; [GO:0008637]; [GO:0097190]; [GO:0006915]; [GO:0007005]; [GO:0007165]; [GO:0097194]; [GO:00125]
P02768	Serum albumin OS=Homo sapiens GN=ALB	[GO:0002576]	[GO:0002576]	[GO:0002576]	[GO:0002576]; [GO:0006887]; [GO:0016192]; [GO:0032940]; [GO:0006810]; [GO:0009987]; [GO:0046903]; [GO:00516]
P04040	Catalase OS=Homo sapiens GN=CAT PE=	[GO:0000302]	[GO:0000302]	[GO:0000302]	[GO:0000302]; [GO:0006979]; [GO:0010035]; [GO:0006950]; [GO:0042221]; [GO:0050896]; [GO:0008150]
P01008	Antithrombin-III OS=Homo sapiens GN=SERPINC1	[GO:0002020]	[GO:0002020]	[GO:0002020]	[GO:0002020]; [GO:0019899]; [GO:0005515]; [GO:0005488]; [GO:0003674]
P62979	Ubiquitin-40S ribosomal protein S27a OS=Homo sapiens GN=UBI1	[GO:0000075]	[GO:0000075]	[GO:0000075]	[GO:0000075]; [GO:0007050]; [GO:0071156]; [GO:0022402]; [GO:0045786]; [GO:0010564]; [GO:0007049]; [GO:000099]
P01112	GTPase HRas OS=Homo sapiens GN=HRAS	[GO:0000139]	[GO:0000139]	[GO:0000139]	[GO:0000139]; [GO:0012505]; [GO:0031090]; [GO:0044431]; [GO:0044464]; [GO:0016020]; [GO:0043227]; [GO:00444]
P61769	Beta-2-microglobulin OS=Homo sapiens GN=MB21D2	[GO:0000139]	[GO:0000139]	[GO:0000139]	[GO:0000139]; [GO:0012505]; [GO:0031090]; [GO:0044431]; [GO:0044464]; [GO:0016020]; [GO:0043227]; [GO:00444]
P69905	Hemoglobin subunit alpha OS=Homo sapiens GN=HBA1	[GO:0004601]	[GO:0004601]	[GO:0004601]	[GO:0004601]; [GO:0016209]; [GO:0016684]; [GO:0003674]; [GO:0016491]; [GO:0003824]
P01344	Insulin-like growth factor II OS=Homo sapiens GN=IGF2	[GO:0001501]	[GO:0001501]	[GO:0001501]	[GO:0001501]; [GO:0048731]; [GO:0007275]; [GO:0048856]; [GO:0032501]; [GO:0032502]; [GO:0008150]
P10699	Thioredoxin OS=Homo sapiens GN=TXNIP	[GO:0000323]	[GO:0000323]	[GO:0000323]	[GO:0000323]; [GO:0000357]; [GO:0000365]; [GO:0045892]; [GO:0000365]; [GO:0000361]; [GO:0010639]; [GO:000513]

4. Move the cursor over the GO Accessions column.

The application displays the annotated GO term and all ancestor terms associated with a protein, as shown in Figure 227. It shows the term annotated to the protein in brackets, followed by their ancestor terms. Each annotated GO term starts on a new line. If you want all proteins to have a higher-level annotation that is not provided by the Molecular Function, Cellular Component, and Biological Process annotation columns, you can filter for the GO term in this column.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

Figure 227. The complete list of GO terms associated with a protein

Checked	Accession	Description	Is Master Protein	Biological Process	Cellular Component	Molecular Function	Pfam IDs	GO Accessions	Entrez Gene
<input type="checkbox"/>	P00167	Cytochrome b5 OS=Homo sapiens GN=CYB5A PE=1 SV=2	✓				Pf00173	[GO:0004033]; GO:000	1528
<input type="checkbox"/>	P05413	Fatty acid-binding protein, heart OS=Homo sapiens GN=FA	✓				Pf00061	[GO:0005215]; GO:000	2170
<input type="checkbox"/>	Q06830	Peroxiredoxin-1 OS=Homo sapiens GN=PRDX1 PE=1 SV=	✓				Pf00578, Pf08534, Pf1	[GO:0000302]; GO:000	5052
<input type="checkbox"/>	P0CG48	Polyubiquitin-1 (GO:0005215); GO:0003674							
<input type="checkbox"/>	P0CG47	Polyubiquitin-2 (GO:0005324); GO:0015245; GO:0005319; GO:0022892; GO:0005215; GO:0003674							
<input type="checkbox"/>	P02144	Myoglobin OS (GO:0005504); GO:0008289; GO:0033293; GO:0005488; GO:0031406; GO:0003674							
<input type="checkbox"/>	P16083	Ribosylidihydroxy (GO:0005515); GO:0005488; GO:0003674							
<input type="checkbox"/>	Q15843	NEDD8 OS=H (GO:0005737); GO:0044424; GO:0005622; GO:0044464; GO:0005575; GO:0005623							
<input type="checkbox"/>	P01127	Platelet-derived (GO:0005829); GO:0044444; GO:0005737; GO:0044424; GO:0005622; GO:0044464; GO:0005575; GO:0005623							
<input type="checkbox"/>	P10145	Interleukin-8 (GO:0006631); GO:0032787; GO:0044255; GO:0019752; GO:0006629; GO:0044237; GO:0043436; GO:0044238; GO:0008152; GO:0009987; GO:0006082; GO:0042180; GO:0008150; GO:0044281							
<input type="checkbox"/>	P08758	Annexin A5 (GO:0006810); GO:0051234; GO:0008150; GO:0051179							
<input type="checkbox"/>	P12081	Histidine-tRNA (GO:0008092); GO:0005515; GO:0005488; GO:0003674							
<input type="checkbox"/>	P68871	Hemoglobin beta (GO:0008285); GO:0008283; GO:0042127; GO:0048523; GO:0008150; GO:0050794; GO:0009987; GO:0048519; GO:0050789; GO:0065007							
<input type="checkbox"/>	P00709	Alpha-lactalbumin (GO:0008289); GO:0005488; GO:0003674							
<input type="checkbox"/>	P62987	Ubiquitin-60S (GO:0015909); GO:0015908; GO:0006869; GO:0015718; GO:0010876; GO:0071702; GO:0046942; GO:0033036; GO:0006810; GO:0015849; GO:0051179; GO:0051234; GO:0008150							
<input type="checkbox"/>	O00762	Ubiquitin-conjugase (GO:0016528); GO:0005737; GO:0044424; GO:0005622; GO:0044464; GO:0005575; GO:0005623							
<input type="checkbox"/>	P00915	Carbonic anhydrase (GO:0032868); GO:0043434; GO:0009725; GO:0009719; GO:0010033; GO:0050896; GO:0042221; GO:0008150							
<input type="checkbox"/>	P55957	BH3-interacting domain (GO:0042493); GO:0042221; GO:0050896; GO:0008150							
<input type="checkbox"/>	P02768	Serum albumin OS=Homo sapiens GN=ALB PE=1 SV=2	✓				Pf00273	[GO:0002576]; GO:000	213
<input type="checkbox"/>	P04040	Catalase OS=Homo sapiens GN=CAT PE=1 SV=3	✓				Pf00199, Pf06628	[GO:0000302]; GO:000	847
<input type="checkbox"/>	P01008	Antithrombin-III OS=Homo sapiens GN=SERPINC1 PE=1 S	✓				Pf00079	[GO:0002020]; GO:000	462
<input type="checkbox"/>	P62979	Ubiquitin-40S ribosomal protein S27a OS=Homo sapiens G	✓				Pf00240, Pf01599, Pf1	[GO:0000075]; GO:000	6233
<input type="checkbox"/>	P01112	GTPase HRas OS=Homo sapiens GN=HRAS PE=1 SV=1	✓				Pf00009, Pf00071, Pf0	[GO:0000139]; GO:000	3265
<input type="checkbox"/>	P61769	Beta-2-microglobulin OS=Homo sapiens GN=B2M PE=1 SV	✓				Pf07654	[GO:0000139]; GO:000	567
<input type="checkbox"/>	P69905	Hemoglobin subunit alpha OS=Homo sapiens GN=HBA1 PE	✓				Pf00042	[GO:0004601]; GO:000	3040
<input type="checkbox"/>	P01344	Insulin-like growth factor II OS=Homo sapiens GN=IGF2 PE	✓				Pf00049, Pf08365	[GO:0001501]; GO:000	3481
<input type="checkbox"/>	P10599	Thioredoxin OS=Homo sapiens GN=TXN PE=1 SV=3	✓				Pf00085	[GO:0000122]; GO:000	7295
<input type="checkbox"/>	P02787	Serotransferrin OS=Homo sapiens GN=TF PE=1 SV=3	✓				Pf00405, Pf12974	[GO:0002576]; GO:000	7018
<input type="checkbox"/>	P00441	Superoxide dismutase [Cu-Zn] OS=Homo sapiens GN=SOD	✓				Pf00080	[GO:0000187]; GO:000	6647
<input type="checkbox"/>	P05019	Insulin-like growth factor I OS=Homo sapiens GN=IGF1 PE	✓				Pf00049	[GO:0001501]; GO:000	3479

Displaying Protein Family (Pfam) Annotations

As noted in “[Pfam Database Annotation](#)” on [page 291](#), you can retrieve Pfam annotations from the Pfam database as an alternative to GO annotations.

❖ To display Protein Family (Pfam) annotations

1. Open the results report (.pdResult) file by following the instructions in the Help.
2. Click the **Proteins** tab.
3. If the Pfam IDs column is not visible on the Proteins page, select it in the Field Chooser dialog box.

For information on the Field Chooser dialog box, refer to the Help.

[Figure 228](#) shows the Pfam IDs column on the Proteins page.

Figure 228. Pfam ID annotations for the protein groups shown in the results of an annotation search

Protein Groups	Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info	Is Master Protein	Biological Process	Cellular Component	Molecular Function	Pfam IDs	GO Accessions	Entrez Gene ID
1	P00167	Cytochrome b5 OS=Homo sapiens GN=CYB5A f	✓						Pf00173	[GO:0004033]; GO:000	1528
2	P05413	Fatty acid-binding protein, heart OS=Homo sapie	✓						Pf00061	[GO:0005215]; GO:000	2170
3	Q06830	Peroxiredoxin-1 OS=Homo sapiens GN=PRDX1	✓						Pf00578, Pf08534, Pf1	[GO:0000302]; GO:000	5052
4	P0CG48	Polyubiquitin-C OS=Homo sapiens GN=UBC PE=	✓						Pf00240, Pf11976	[GO:0000075]; GO:000	7316
5	P0CG47	Polyubiquitin-B OS=Homo sapiens GN=UBB PE=	✓						Pf00240, Pf11976	[GO:0000075]; GO:000	7314
6	P02144	Myoglobin OS=Homo sapiens GN=MB PE=1 SV=	✓						Pf00042	[GO:0001666]; GO:000	4151
7	P16083	Ribosyl-dihydroxynicotinamide dehydrogenase [quin	✓						Pf02525, Pf03358	[GO:0001512]; GO:000	4835
8	Q15843	NEDD8 OS=Homo sapiens GN=NEDD8 PE=1 S'	✓						Pf00240, Pf11976	[GO:0005515]; GO:000	4738
9	P01127	Platelet-derived growth factor subunit B OS=Horr	✓						Pf00341, Pf04692	[GO:0000139]; GO:000	5155
10	P10145	Interleukin-8 OS=Homo sapiens GN=IL8 PE=1 S	✓						Pf00048	[GO:0001525]; GO:000	3576
11	P08758	Annexin A5 OS=Homo sapiens GN=ANXA5 PE=	✓						Pf00191	[GO:0004859]; GO:000	308
12	P12081	Histidine-tRNA ligase, cytoplasmic OS=Homo sz	✓						Pf00458, Pf00587, Pf0	[GO:0000166]; GO:000	3035
13	P68871	Hemoglobin subunit beta OS=Homo sapiens GN=	✓						Pf00042	[GO:0004601]; GO:000	3043
14	P00709	Alpha-lactalbumin OS=Homo sapiens GN=LALB	✓						Pf00062	[GO:0004461]; GO:000	3906
15	P62987	Ubiquitin-60S ribosomal protein L40 OS=Homo s	✓						Pf00240, Pf01020, Pf1	[GO:0000075]; GO:000	7311
16	O00762	Ubiquitin-conjugating enzyme E2 C OS=Homo sz	✓						Pf00179, Pf05743	[GO:0000075]; GO:000	11065
17	P00915	Carbonic anhydrase 1 OS=Homo sapiens GN=C-	✓						Pf00194	[GO:0004089]; GO:000	759
18	P55957	BH3-interacting domain death agonist OS=Homo	✓						Pf06393	[GO:0001836]; GO:000	637
19	P02768	Serum albumin OS=Homo sapiens GN=ALB PE=	✓						Pf00273	[GO:0002576]; GO:000	213
20	P04040	Catalase OS=Homo sapiens GN=CAT PE=1 SV=	✓						Pf00199, Pf06628	[GO:0000302]; GO:000	847
21	P01008	Antithrombin-III OS=Homo sapiens GN=SERPIN	✓						Pf00079	[GO:0002020]; GO:000	462
22	P62979	Ubiquitin-40S ribosomal protein S27a OS=Homo	✓						Pf00240, Pf01599, Pf1	[GO:0000075]; GO:000	6233
23	P01112	GTPase HRas OS=Homo sapiens GN=HRAS PE	✓						Pf00009, Pf00071, Pf0	[GO:0000139]; GO:000	3265
24	P61769	Beta-2-microglobulin OS=Homo sapiens GN=B2I	✓						Pf07654	[GO:0000139]; GO:000	567
25	P69905	Hemoglobin subunit alpha OS=Homo sapiens GN	✓						Pf00042	[GO:0004601]; GO:000	3040
26	P01344	Insulin-like growth factor II OS=Homo sapiens GN	✓						Pf00049, Pf08365	[GO:0001501]; GO:000	3481
27	P10599	Thioredoxin OS=Homo sapiens GN=TXN PE=1 S	✓						Pf00085	[GO:0000122]; GO:000	7295
28	P02787	Serotransferrin OS=Homo sapiens GN=TF PE=1	✓						Pf00405, Pf12974	[GO:0002576]; GO:000	7018
29	P00441	Superoxide dismutase [Cu-Zn] OS=Homo sapien	✓						Pf00080	[GO:0000187]; GO:000	6647
30	P05019	Insulin-like growth factor I OS=Homo sapiens GN	✓						Pf00049	[GO:0001501]; GO:000	3479
31	P62937	Peptidyl-prolyl cis-trans isomerase A OS=Homo s	✓						Pf00160	[GO:0000413]; GO:000	5478

Displaying Entrez Gene Database Identifications

Entrez Gene database identifications are unique identifications assigned to all genes stored in the Entrez Gene database, NCBI's database of gene-specific information. The Proteome Discoverer application displays these identifications in the Entrez Gene ID column on the Proteins page, as shown in [Figure 229](#). All proteins derived from the same gene have the same gene identification. You can use this information to group or cluster biologically meaningful proteins together. Because not all genes are stored in the Entrez Gene database, some proteins do not have a valid gene identification. In this case, the column is empty. For more information on the Entrez Gene identifications, see [“Entrez Gene Database Annotation”](#) on [page 292](#).

❖ To display Entrez Gene database identifications

1. Open the results report (.pdResult) file by following the instructions in the Help.
2. Click the **Proteins** tab.
3. In the Field Chooser dialog box of the Proteins page, select the **Entrez Gene ID** column.

For information on the Field Chooser dialog box, refer to the Help.

8 Obtaining Protein Annotation Information

Displaying the Annotated Protein Results

The application displays the gene identifications in the Entrez Gene ID column of the Proteins page of the .pdResult) file, as shown in [Figure 229](#). Because not all genes are stored in the Entrez Gene database, some proteins do not have a valid gene identification. In this case, the column displays a value of 0.

Figure 229. Entrez Gene ID column on the Proteins page

Protein Groups	Checked	Accession	Description	Is Master Protein	Biological Process	Cellular Component	Molecular Function	Pfam IDs	GO Accessions	Entrez Gene ID
1	<input type="checkbox"/>	P00167	Cytochrome b5 OS=Homo sapiens GN=CYE	✓				PF00173	[GO:0004033]; GO:000	1528
2	<input type="checkbox"/>	P05413	Fatty acid-binding protein, heart OS=Homo s	✓				PF00061	[GO:0005215]; GO:000	2170
3	<input type="checkbox"/>	Q06830	Peroxisomal protein 1 OS=Homo sapiens GN=PR	✓				PF00578, PF08534, PF1	[GO:0000302]; GO:000	5052
4	<input type="checkbox"/>	P0CG48	Polyubiquitin-C OS=Homo sapiens GN=UBC	✓				PF00240, PF11976	[GO:0000075]; GO:000	7316
5	<input type="checkbox"/>	P0CG47	Polyubiquitin-B OS=Homo sapiens GN=UBE	✓				PF00240, PF11976	[GO:0000075]; GO:000	7314
6	<input type="checkbox"/>	P02144	Myoglobin OS=Homo sapiens GN=MB PE=1	✓				PF00042	[GO:0001666]; GO:000	4151
7	<input type="checkbox"/>	P16083	Ribosylhydronicotinamide dehydrogenase	✓				PF02525, PF03358	[GO:0001512]; GO:000	4835
8	<input type="checkbox"/>	Q15843	NEDD8 OS=Homo sapiens GN=NEDD8 PE=	✓				PF00240, PF11976	[GO:0005515]; GO:000	4738
9	<input type="checkbox"/>	P01127	Platelet-derived growth factor subunit B OS=	✓				PF00341, PF04692	[GO:0000139]; GO:000	5155
10	<input type="checkbox"/>	P10145	Interleukin-8 OS=Homo sapiens GN=IL8 PE=	✓				PF00048	[GO:0001525]; GO:000	3576
11	<input type="checkbox"/>	P08758	Annexin A5 OS=Homo sapiens GN=ANXA5	✓				PF00191	[GO:0004859]; GO:000	308
12	<input type="checkbox"/>	P12081	Histidine-tRNA ligase, cytoplasmic OS=Hon	✓				PF00458, PF00587, PF0	[GO:0000166]; GO:000	3035
13	<input type="checkbox"/>	P68871	Hemoglobin subunit beta OS=Homo sapiens	✓				PF00042	[GO:0004601]; GO:000	3043
14	<input type="checkbox"/>	P00709	Alpha-lactalbumin OS=Homo sapiens GN=L	✓				PF00062	[GO:0004461]; GO:000	3906
15	<input type="checkbox"/>	P62987	Ubiquitin-60S ribosomal protein L40 OS=Ho	✓				PF00240, PF01020, PF1	[GO:0000075]; GO:000	7311
16	<input type="checkbox"/>	O00762	Ubiquitin-conjugating enzyme E2 C OS=Hor	✓				PF00179, PF05743	[GO:0000075]; GO:000	11065
17	<input type="checkbox"/>	P00915	Carbonic anhydrase 1 OS=Homo sapiens G	✓				PF00194	[GO:0004089]; GO:000	759
18	<input type="checkbox"/>	P55957	BH3-interacting domain death agonist OS=H	✓				PF06393	[GO:0001836]; GO:000	637
19	<input type="checkbox"/>	P02768	Serum albumin OS=Homo sapiens GN=ALB	✓				PF00273	[GO:0002576]; GO:000	213
20	<input type="checkbox"/>	P04040	Catalase OS=Homo sapiens GN=CAT PE=1	✓				PF00199, PF06628	[GO:0000302]; GO:000	847
21	<input type="checkbox"/>	P01008	Antithrombin-III OS=Homo sapiens GN=SEF	✓				PF00079	[GO:0002020]; GO:000	462
22	<input type="checkbox"/>	P62979	Ubiquitin-40S ribosomal protein S27a OS=H	✓				PF00240, PF01599, PF1	[GO:0000075]; GO:000	6233
23	<input type="checkbox"/>	P01112	GTPase HRas OS=Homo sapiens GN=HRA	✓				PF00009, PF00071, PF0	[GO:0000139]; GO:000	3265
24	<input type="checkbox"/>	P61769	Beta-2-microglobulin OS=Homo sapiens GN	✓				PF07654	[GO:0000139]; GO:000	567
25	<input type="checkbox"/>	P69905	Hemoglobin subunit alpha OS=Homo sapier	✓				PF00042	[GO:0004601]; GO:000	3040
26	<input type="checkbox"/>	P01344	Insulin-like growth factor II OS=Homo sapier	✓				PF00049, PF08365	[GO:0001501]; GO:000	3481
27	<input type="checkbox"/>	P10599	Thioredoxin OS=Homo sapiens GN=TXN PE=	✓				PF00085	[GO:0000122]; GO:000	7295
28	<input type="checkbox"/>	P02787	Serotransferrin OS=Homo sapiens GN=TF F	✓				PF00405, PF12974	[GO:0002576]; GO:000	7018
29	<input type="checkbox"/>	P00441	Superoxide dismutase [Cu-Zn] OS=Homo sa	✓				PF00080	[GO:0000187]; GO:000	6647
30	<input type="checkbox"/>	P05019	Insulin-like growth factor I OS=Homo sapien	✓				PF00049	[GO:0001501]; GO:000	3479
31	<input type="checkbox"/>	P62937	Peptidyl-prolyl cis-trans isomerase A OS=Hc	✓				PF00160	[GO:0000413]; GO:000	5478

Displaying Ensembl Genome Database Annotations

The Ensembl genome database stores and automatically annotates the reference genomes of the genetically sequenced organisms. For fully sequenced species, the Ensembl gene identification and chromosome information is available if you retrieve protein annotations from ProteinCenter. If this information is available, it appears in the Ensembl Gene ID and Chromosome columns on the Proteins page of the .pdResult file, as shown in [Figure 230](#). As with the Entrez Gene identifications, you can use this information to group or cluster proteins by these values.

❖ To display Ensembl genome database annotations

1. Open the results report (.pdResult) file by following the instructions in the Help.
2. Click the **Proteins** tab.
3. In the Field Chooser dialog box of the proteins page, select the **Ensembl Gene ID** and the **Chromosome** column.

For information on the Field Chooser dialog box, refer to the Help.

The application displays the genome and chromosome identifications in the Ensembl Gene ID and Chromosome columns, respectively, of the Proteins page of the .pdResult file, as shown in Figure 230.

Figure 230. Columns for Ensembl Genome database identifications and chromosome number

Checked	Accession	Description	Coverage	# Peptides	# PSMs	# AAs	Ensembl Gene ID	Chromosome	Biological Process	Cellular Component	Molecular Function
<input type="checkbox"/>	P07900	Heat shock protein HSP 90-alpha OS=H	13%	8	21	732	ENSG00000080824	14			
<input type="checkbox"/>	B4DMA2	cDNA FLJ54023, highly similar to Heat	7%	5	20	686		6			
<input type="checkbox"/>	P08238	Heat shock protein HSP 90-beta OS=Hc	7%	5	20	724	ENSG000000096384	6			
<input type="checkbox"/>	Q2VPJ6	HSP90AA1 protein (Fragment) OS=Hon	14%	7	19	585		14			
<input type="checkbox"/>	B4DGL0	cDNA FLJ53619, highly similar to Heat	6%	4	19	714		6			
<input type="checkbox"/>	Q8TBA7	HSP90AA1 protein (Fragment) OS=Hon	11%	6	19	638		14			
<input type="checkbox"/>	Q6UYC3	Lamin A/C OS=Homo sapiens GN=LMN	13%	7	15	614	ENSG00000160789	1			
<input type="checkbox"/>	Q5I6Y4	Lamin A/C transcript variant 1 OS=Hom	12%	7	15	664		1			
<input type="checkbox"/>	P02545	Prelamin-A/C OS=Homo sapiens GN=L	12%	7	15	664	ENSG00000160789	1			
<input type="checkbox"/>	Q5I6Y6	Lamin A/C transcript variant 1 OS=Hom	12%	7	15	664		1			
<input type="checkbox"/>	Q58FF7	Putative heat shock protein HSP 90-bet	8%	4	18	597		4			
<input type="checkbox"/>	Q6PJ43	ACTG1 protein (Fragment) OS=Homo s	9%	2	19	263		17			
<input type="checkbox"/>	B4E335	cDNA FLJ52842, highly similar to Actin.	7%	2	19	351		7			
<input type="checkbox"/>	Q53G99	Beta actin variant (Fragment) OS=Homc	6%	2	19	375		7			
<input type="checkbox"/>	Q8WVW5	Putative uncharacterized protein (Fragr	6%	2	19	363		17			
<input type="checkbox"/>	P63261	Actin, cytoplasmic 2 OS=Homo sapiens	6%	2	19	375	ENSG00000267807; E	HG21_PATCH: 17			
<input type="checkbox"/>	B4DVQ0	cDNA FLJ58286, highly similar to Actin.	7%	2	19	333		17			
<input type="checkbox"/>	Q53GK6	Beta actin variant (Fragment) OS=Homc	6%	2	19	375		7			
<input type="checkbox"/>	B3KVVQ3	cDNA FLJ43573 fis, clone RECTM2001	9%	2	19	253		7			
<input type="checkbox"/>	Q1KLZ0	HCG15971, isoform CRA_a OS=Homo	6%	2	19	375	ENSG00000075624	7			
<input type="checkbox"/>	B7ZAP6	cDNA, FLJ79260, highly similar to Actin	8%	2	19	294		7			
<input type="checkbox"/>	B4E3A4	cDNA FLJ57283, highly similar to Actin.	6%	2	19	356		17			
<input type="checkbox"/>	Q53G76	Beta actin variant (Fragment) OS=Homc	6%	2	19	375		7			
<input type="checkbox"/>	B4DW52	cDNA FLJ55253, highly similar to Actin.	7%	2	19	347		7			
<input type="checkbox"/>	Q5HY54	Filamin-A OS=Homo sapiens GN=FLNA	4%	8	14	2607	ENSG00000196924; E	X; HG1497_PATCH			
<input type="checkbox"/>	P21333	Filamin-A OS=Homo sapiens GN=FLNA	4%	8	14	2647	ENSG00000269329; E	HG1497_PATCH: X			
<input type="checkbox"/>	Q60FE5	Filamin A OS=Homo sapiens GN=FLNA	4%	8	14	2620		X			
<input type="checkbox"/>	Q60FE6	Filamin A OS=Homo sapiens GN=FLNA	4%	8	14	2612		X			
<input type="checkbox"/>	Q86U12	Full-length cDNA clone CS0CAP007YF:	14%	5	17	413	ENSG00000080824	14			
<input type="checkbox"/>	Q3BDU5	Rhabdomyosarcoma antigen MU-RMS-	13%	6	13	487		1			
<input type="checkbox"/>	Q8N519	LMNA protein OS=Homo sapiens PE=2	14%	6	13	465		1			
<input type="checkbox"/>	Q5I6Y5	Lamin A/C transcript variant 1 OS=Hom	14%	6	13	480		1			
<input type="checkbox"/>	Q5TCJ4	Lamin A/C OS=Homo sapiens GN=LMN	14%	6	13	465		1			
<input type="checkbox"/>	Q5TCI8	Lamin A/C OS=Homo sapiens GN=LMN	12%	5	12	491	ENSG00000160789	1			
<input type="checkbox"/>	B4DI32	cDNA FLJ56081, highly similar to Lamin	10%	5	12	574		1			
<input type="checkbox"/>	Q5TCI9	Lamin A/C (Fragment) OS=Homo sapien	11%	5	12	513		1			
<input type="checkbox"/>	D6RAQ3	Prelamin-A/C OS=Homo sapiens GN=L	10%	5	12	565	ENSG00000160789	1			
<input type="checkbox"/>	E7EUI9	Prelamin-A/C OS=Homo sapiens GN=L	10%	5	12	574	ENSG00000160789	1			
<input type="checkbox"/>	E7EUY0	DNA-dependent protein kinase catalytic	3%	9	14	4096	ENSG00000253729	8			
<input type="checkbox"/>	P78527	DNA-dependent protein kinase catalytic	3%	9	14	4128		8			

Displaying UniProt Annotations

For information on displaying UniProt PTM annotations, refer to the Help.

Uploading Results to ProteinCenter

If you have a user account on a ProteinCenter server, you can upload search results directly from the Proteome Discoverer application to ProteinCenter. You can upload your results to ProteinCenter by exporting the results through a node in the consensus workflow or by exporting the results in an open .pdResult file.

❖ To upload search results to ProteinCenter through the consensus workflow

1. Configure the Proteome Discoverer application for uploading search results to ProteinCenter. For instructions, see [“To configure the application for uploading results to ProteinCenter in the Administration view”](#) on page 293.
2. Create or open a study and an analysis:
 - To create a study, see [“Creating a Study”](#) on page 38.
 - To open an existing study, see [“Opening an Existing Study”](#) on page 40.
 - To create an analysis, see [“Creating an Analysis”](#) on page 71.
 - To open an existing analysis, see [“Opening an Existing Analysis”](#) on page 72.


3. Open the Workflow Editor and create a processing and a consensus workflow.

For information on creating a processing workflow, see [“Creating a Processing Workflow”](#) on page 104.

For information on creating a consensus workflow, see [“Creating a Consensus Workflow”](#) on page 112.

4. In the consensus workflow, place the To ProteinCenter node in the Post-Processing Nodes pane.

For information on creating a quantification processing workflow, see [“Creating a Processing Workflow for Precursor Ion Quantification”](#) on page 390 or [“Creating a Processing Workflow for Reporter Ion Quantification”](#) on page 397.

5. In the Target Name parameter of the To ProteinCenter node, specify the name of the data set to upload.
6. Save the analysis. See [“Saving an Analysis”](#) on page 80.
7. Save the study. See [“Saving a Study”](#) on page 62.
8. Click the **Run** icon,  **Run**, in the upper right corner of the Analysis window.

The job queue appears, as shown in [Figure 79](#) on page 112, displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

The To ProteinCenter node exports the search results from a .pdResult file to ProtXML and then uploads them to the ProteinCenter server.

If the node encounters errors, it tries several times to upload the file. The number of times that it tries is specified by the Number of Attempts to Submit the Annotation Request and the Time Interval Between Attempts to Submit the Annotation Request [sec] parameters in the ProteinCenter view shown in [Figure 206](#) on page 294. If you are not authorized to submit the data, the node tries only once.

If an error occurs during upload to ProteinCenter, the workflow does not fail, but the job queue displays an error message.


❖ **To upload search results to ProteinCenter from an open .pdResult file**

1. Configure the Proteome Discoverer application for uploading search results to ProteinCenter. For instructions, see “[To configure the application for uploading results to ProteinCenter in the Administration view](#)” on [page 293](#).
2. Open a .pdResult file in the Proteome Discoverer application. Refer to the Help.
3. To export only the result data from selected protein groups, in the .pdResult file select the Checked box in the row of the protein groups that you want to include.
4. Choose **Tools > To ProteinCenter**.

The Export to ProteinCenter dialog box opens.

5. In the Destination box, specify the name of the data set to upload to ProteinCenter, as shown in [Figure 231](#).

Figure 231. Export to ProteinCenter dialog box



6. Select the **Checked Protein Groups** check box in the Export to ProteinCenter dialog box.

If you do not select Checked Protein Groups, the application exports the result data of all protein groups.


7. Click **Export**.

After the application exports the data set to ProteinCenter, you can log in to your ProteinCenter account. The uploaded data set appears under the Incoming node in the ProteinCenter window, as shown in [Figure 232](#).

8 Obtaining Protein Annotation Information

Accessing the ProteinCard Page

Figure 232. Uploaded data set in ProteinCenter window



Accessing the ProteinCard Page

You can access the data in ProteinCenter through the ProteinCard page of the Protein Identification Details view for each protein. In the ProteinCard page, a protein is considered a specific amino acid sequence in a given species.

❖ To open the ProteinCard page

1. Double-click a grid cell on the Proteins page of the .pdResult file.

—OR—

Select a cell and choose **View > Protein Details**.

—OR—

Click the **Protein Details** icon, .

2. In the Protein Identification Details view, click the **ProteinCard** tab.

You might experience a short delay as the Proteome Discoverer application accesses the URL.

After loading data from the ProteinCenter server, the Proteome Discoverer application displays the data in the ProteinCard page. By default, the ProteinCard page shows the General page, shown in [Figure 234](#) on [page 329](#).


3. Click on the tab of the page containing the information that you are seeking:

- [General Page](#)
- [Keys Page](#)
- [Features Page](#)
- [Molecular Functions Page](#)
- [Cellular Components Page](#)
- [Biological Processes Page](#)
- [Diseases Page](#)
- [External Links Page](#)

4. Click **OK** to close the Protein Identification Details dialog box.

If the entire protein is not found in ProteinCenter but a protein with the same sequence exists, the ProteinCard page displays a warning that the displayed information is from a protein with different accession, as shown in [Figure 233](#). If there is more than one protein with the same sequence but from different organisms, an additional list box appears so that you can select the correct species.

Figure 233. Warning displayed for protein with different accession



ProteinCard Page Parameters

The ProteinCard page of the Protein Identification Details dialog box contains the following pages.

- [General Page](#)
- [Keys Page](#)
- [Features Page](#)
- [Molecular Functions Page](#)
- [Cellular Components Page](#)
- [Biological Processes Page](#)
- [Diseases Page](#)
- [External Links Page](#)

General Page

The General page of the ProteinCard page, shown in [Figure 234](#), displays information about the protein: its name, its description, its function, the keywords that produce it in a database search, and the gene that ultimately directs the protein's synthesis through RNA.

Figure 234. General page of the ProteinCard page




Table 15 lists the parameters on the General page of the ProteinCard page.

Table 15. Parameters on the General page of the ProteinCard page (Sheet 1 of 2)

Parameter	Description
Top area	Displays the protein name in bold font on the first line. The second line in bold font is the official symbol of the gene that ultimately directs the synthesis of the protein through RNA, and the text following it is the alternative name or names of the gene.
Top right area	Displays the name of the species that contains the gene that ultimately directs the synthesis of this protein through RNA, the number of the chromosome that the gene resides on, and the location of the chromosome that the gene resides on. The name of the species links to the National Center for Biotechnology Information (NCBI) taxonomy browser.

8 Obtaining Protein Annotation Information

ProteinCard Page Parameters

Table 15. Parameters on the General page of the ProteinCard page (Sheet 2 of 2)

Parameter	Description
Gene Details area	Displays information about the gene that directs the synthesis of the protein. If no information about the gene is available, the application provides a link to the Entrez Gene website.
Protein Details area	Lists the keywords that produce this protein in a database search, the functions of the protein, and a description of the protein.

Keys Page

The Keys page of the ProteinCard page, shown in [Figure 235](#), lists all the accession keys for a given protein.

Figure 235. Keys page of the ProteinCard page

Primary Key	Src	Secondary Key	Src	Description
ENSP00000434657	ESBL			-
ENSP00000434412	ESBL			-
ENSP00000350894	ESBL			-
333360851	GI	NP_001193943.1	REF	serpin H1 precursor
		hsa:871	KEGG	
193783823	GI	BAG53805.1	DBJ	unnamed protein product
119595385	GI	EAW74979.1	GB	serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1), isoform CRA_a
119595383	GI	EAW74977.1	GB	serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1), isoform CRA_a
119595381	GI	EAW74975.1	GB	serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1), isoform CRA_a
119595380	GI	EAW74974.1	GB	serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1), isoform CRA_a
53831040	GI	AAU95378.1	GB	arsenic-transactivated protein 3
37589973	GI	AAH36298.2	GB	Serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1)
33090237	GI	AAP93914.1	GB	proliferation-inducing gene 14
32454741	GI	NP_001226.2	REF	serpin H1 precursor
30583027	GI	AAP35758.1	GB	serine (or cysteine) proteinase inhibitor, clade H (heat shock protein 47), member 2
20141241	GI	P50454.2	SP	Serpin H1
15779117	GI	AAH14623.1	GB	Serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1)
8574447	GI	BAA96789.1	DBJ	rheumatoid arthritis related antigen RA-A47
8574445	GI	BAA96788.1	DBJ	rheumatoid arthritis related antigen RA-A47
IPI00032140.4	IPI			Serpin H1
P50454	SP	Q8IY96	UNI	Serpin H1
		SERPH_HUMAN	UNI	
		SPH2_HUMAN	UNI	

Table 16 lists the parameters on the Keys page of the ProteinCard page.

Table 16. Parameters on the Keys page of the ProteinCard page

Parameter	Description
Primary Key	Lists the accession key of the database that the sequence was imported from. It is linked to the original database records in the source database, such as Ensembl, SGD, NRDB, IPI or UniProt. This parameter emphasizes the preferred type of accession.
Src	Specifies the abbreviation of the primary source database.
Secondary Key	Lists the secondary accession key, which is either an alternative key used in the source database or the key of the original database.
Src	Specifies the abbreviation of the secondary source database.
Description	Displays the original description for the original database entry.

An exclamation mark flags outdated protein keys, and the keys are linked to the outdated history in their respective source database.

Features Page

The Features page of the ProteinCard page, shown in Figure 236, includes a selection of sequence features from UniProt, from various conserved domain predictions, and from the computational enrichment undertaken by ProteinCenter. (Computational enrichment refers to information that has no experimental evidence but was found by using a computer prediction program.) The features are sorted according to their start positions in the protein sequence.

8 Obtaining Protein Annotation Information

ProteinCard Page Parameters

Figure 236. Features page of the ProteinCard page




Table 17 lists the parameters on the Features page of the ProteinCard page.


Table 17. Parameters on the Features page of the ProteinCard page

Parameter	Description
Source	Specifies the name of the database that provided information about the feature: <ul style="list-style-type: none"> • InterPro • Tmap (computational enrichment) • PrediSi (computational enrichment) • Pfam (computational enrichment) • UniProt
Category	Displays the type of information that UniProt, InterPro, and Tmap include for each row. For example, UniProt might include “CARBOHYD” as one of its types of information, and InterPro might include “SSF57184” as one of its types of information.
From	Specifies the start position of the amino acid.
To	Specifies the end position of the amino acid.
Acc	Specifies the accession identifier for the domain linked to InterPro or Pfam.
Description	Describes the feature.

Molecular Functions Page

The Molecular Functions page of the ProteinCard page, shown in [Figure 237](#), summarizes information about the function of the protein. It consolidates GO data and Enzyme Category (EC) information. The EC designation indicates whether a protein has been categorized with a certain enzyme function.

Figure 237. Molecular Functions page of the ProteinCard page



[Table 18](#) lists the parameters on the Molecular Functions page of the ProteinCard page.

Table 18. Parameters on the Molecular Functions page of the ProteinCard page (Sheet 1 of 2)

Parameter	Description
GO Id	Lists the GO code for each of the protein's molecular functions. Each code is linked to the QuickGO browser of the European Bioinformatics Institute (EBI), which hosts several databases and services.
Evidence Codes	Lists the evidence codes for each of the protein's molecular functions for GO annotation. Evidence codes describe how the GO information was proven—for example, by computer prediction or by experiment.
PMIDs	Lists the molecular function codes in the PubMed database, which is maintained by the U.S. National Library of Medicine (NLM) and the National Institutes of Health (NIH). Each code is linked to the PubMed browser.

Table 18. Parameters on the Molecular Functions page of the ProteinCard page (Sheet 2 of 2)


Parameter	Description
Go Slim	Specifies the basic GO Slim category for the GO term. GO Slim categories are reduced versions of the GO ontologies containing a subset of the terms in the entire GO database. They give a broad overview of the ontology content without the detail of the specific fine-grained terms. Table 21 on page 339 provides the Go Slim categories for molecular functions.
Name	Describes the molecular function for a GO term. This description is created by the GO consortium.

Enzymes with an EC number for IUBMB Enzyme Nomenclature are displayed with links to detailed information at the International Union of Biochemistry and Molecular Biology.

Cellular Components Page

The Cellular Components page of the ProteinCard page, shown in [Figure 238](#), summarizes information about where the protein carries out its function in the cell.

Figure 238. Cellular Components page of the ProteinCard page



[Table 19](#) lists the parameters on the Cellular Components page of the ProteinCard page.

Table 19. Parameters on the Cellular Components page of the ProteinCard page

Parameter	Description
GO Id	Lists the GO code for each of the protein's molecular functions. Each code is linked to the QuickGO browser of the EBI, which hosts a number of databases and services.
Evidence Codes	Lists the evidence codes for each of the protein's cellular components for GO annotation. Evidence codes describe how the GO information was proven—for example, by computer prediction or by experiment.
PMIDs	Lists the cellular component codes in the PubMed database, which is maintained by the NLM and the NIH. Each code is linked to the PubMed browser.
Go Slim	Specifies the basic GO Slim category for the GO term. GO Slim categories are reduced versions of the GO ontologies containing a subset of the terms in the entire GO database. They give a broad overview of the ontology content without the detail of the specific fine-grained terms. Table 22 on page 340 provides the Go Slim categories for cellular components.
Name	Describes the cellular component for a GO term. This description is created by the GO consortium.

Enzymes with an EC number for IUBMB Enzyme Nomenclature are displayed with links to detailed information at the International Union of Biochemistry and Molecular Biology.

Biological Processes Page

The Biological Processes page of the ProteinCard page, shown in [Figure 239](#), summarizes information about the biological processes that the protein is a part of.

Figure 239. Biological Processes page of the ProteinCard page




Table 20 lists the parameters on the Biological Processes page of the ProteinCard page.

Table 20. Parameters on the Biological Processes page of the ProteinCard page


Parameter	Description
GO Id	Lists the GO code for each of the protein's molecular functions. Each code is linked to the QuickGO browser of the EBI, which hosts a number of databases and services.
Evidence Codes	Lists the evidence codes for each of the protein's biological processes for GO annotation. Evidence codes describe how the GO information was proven—for example, by computer prediction or by experiment.
PMIDs	Lists the biological process codes in the PubMed database, which is maintained by the NLM and the NIH. Each code is linked to the PubMed browser.
Go Slim	Specifies the basic GO Slim category for the GO term. GO Slim categories are reduced versions of the GO ontologies containing a subset of the terms in the entire GO database. They give a broad overview of the ontology content without the detail of the specific fine-grained terms. Table 23 provides the Go Slim categories for biological components.
Name	Describes the biological process for a GO term. This description is created by the GO consortium.

Enzymes with an EC number for IUBMB Enzyme Nomenclature are displayed with links to detailed information at the International Union of Biochemistry and Molecular Biology.

Diseases Page

The Diseases page of the ProteinCard page, shown in [Figure 240](#), lists the diseases that the selected protein is associated with.

Figure 240. Diseases page of the ProteinCard page




External Links Page

The External Links page of the ProteinCard page, shown in [Figure 241](#), lists the web links to resources containing information about the protein.

8 Obtaining Protein Annotation Information

GO Slim Categories

Figure 241. External Links page of the ProteinCard page



Click the appropriate link to open the browser for the database. The external links contain links to resources containing information about the respective protein.

GO Slim Categories

This topic defines the GO Slim terms for molecular functions, cellular components, and biological processes.

GO Slim Categories for Molecular Functions

Table 21 describes the GO Slim categories for molecular functions.

Table 21. GO Slim categories for molecular functions (Sheet 1 of 2)

GO Slim molecular function	Description
Antioxidant activity	Inhibition of the reactions brought about by dioxygen (O ₂) or peroxides. Usually the antioxidant is effective because it can be more easily oxidized than the substance protected. The term is often applied to components that can trap free radicals, breaking the chain reaction that normally leads to extensive biological damage.
Catalytic activity	Catalysis of a biochemical reaction at physiological temperatures. In biologically catalyzed reactions, the reactants are known as substrates, and the catalysts are naturally occurring macromolecular substances known as enzymes. Enzymes possess specific binding sites for substrates and are usually composed wholly or largely of protein.
DNA binding	Selective interaction with DNA (deoxyribonucleic acid).
Enzyme regulator activity	Modulation of an enzyme.
Metal ion binding	Selective interaction with any metal ion.
Motor activity	Catalysis of movement along a polymeric molecule such as a microfilament or microtubule, coupled to the hydrolysis of a nucleoside triphosphate.
Nucleotide binding	Selective interaction with a nucleotide, which is any compound consisting of a nucleoside that is esterified with (ortho)phosphate or an oligophosphate at any hydroxyl group on the ribose or deoxyribose moiety.
Protein binding	Selective interaction with any protein or protein complex (a complex of two or more proteins that may include other nonprotein molecules).
Receptor activity	The mediation by protein or gene products of a signal from the extracellular environment to an intracellular messenger.
RNA binding	Selective interaction with an RNA molecule or a portion of it.
Signal transducer activity	Mediation of the transfer of a signal from the outside to the inside of a cell by means other than the introduction of the signal molecule itself into the cell.

Table 21. GO Slim categories for molecular functions (Sheet 2 of 2)

GO Slim molecular function	Description
Structural molecule activity	The action of a molecule that contributes to the structural integrity of a complex or assembly within or outside a cell.
Transcription regulator activity	Activity that plays a role in regulating transcription; it might bind a promoter or enhancer DNA sequence or interact with a DNA-binding transcription factor.
Translation regulator activity	The initiation, activation, perpetuation, repression, or termination of polypeptide synthesis at the ribosome.
Transporter activity	Activity that enables the directed movement of substances (such as macromolecules, small molecules, ions) into, out of, within, or between cells.

GO Slim Categories for Cellular Components

Table 22 describes the GO Slim categories for cellular components.

Table 22. GO Slim categories for cellular components (Sheet 1 of 4)

GO Slim cellular component	Description
Cell surface	Proteins that are attached to the external part of the cell wall, cell membrane, or both.
Chromosome	A structure composed of a very long molecule of DNA and associated proteins (for example, histones) that carry hereditary information.
Cytoplasm	All of the contents of a cell excluding the plasma membrane and nucleus but including other subcellular structures.
Cytoskeleton	Any of the various filamentous elements that form the internal framework of cells and that typically remain after treatment of the cells with mild detergent to remove membrane constituents and soluble components of the cytoplasm. The term embraces intermediate filaments, microfilaments, microtubules, the microtrabecular lattice, and other structures characterized by a polymeric filamentous nature and long-range order within the cell. The various elements of the cytoskeleton not only serve in the maintenance of cellular shape but also have roles in other cellular functions, including cellular movement, cell division, endocytosis, and movement of organelles.
Cytosol	That part of the cytoplasm that does not contain membranous or particulate subcellular components.

Table 22. GO Slim categories for cellular components (Sheet 2 of 4)

GO Slim cellular component	Description
Endosome	A membrane-bound organelle that carries materials newly ingested by endocytosis. It passes many of the materials to lysosomes for degradation.
Endoplasmatic reticulum	The irregular network of unit membranes, visible only by electron microscopy, that occurs in the cytoplasm of many eukaryotic cells. The membranes form a complex meshwork of tubular channels, which are often expanded into slit-like cavities called cisternae. The endoplasmatic reticulum takes two forms, rough (or granular), with ribosomes adhering to the outer surface, and smooth, with no ribosomes attached.
Extracellular	The space external to the outermost structure of a cell. For cells without external protective or external encapsulating structures, this term refers to the space outside of the plasma membrane. It only applies to proteins that are not attached to the cell surface. It covers the host cell environment outside an intracellular parasite.
Golgi	A compound membranous cytoplasmic organelle of eukaryotic cells consisting of flattened, ribosome-free vesicles arranged in a more or less regular stack. The Golgi apparatus differs from the endoplasmic reticulum in often having slightly thicker membranes, appearing in sections as a characteristic shallow semicircle so that the convex side (cis or entry face) abuts the endoplasmic reticulum, secretory vesicles emerging from the concave side (trans or exit face). In vertebrate cells, there is usually one such organelle, but in invertebrates and plants, where they are known usually as dictyosomes, there may be several scattered in the cytoplasm. The Golgi apparatus processes proteins produced on the ribosomes of the rough endoplasmic reticulum. Such processing includes modification of the core oligosaccharides of glycoproteins and the sorting and packaging of proteins for transport to a variety of cellular locations.

Table 22. GO Slim categories for cellular components (Sheet 3 of 4)

GO Slim cellular component	Description
Membrane	<p>Double layer of lipid molecules that encloses all cells, and, in eukaryotic cells, many organelles. The membrane can be a single or double lipid bilayer. It also includes associated proteins.</p> <p>Note This term is not restricted to the plasma membrane but applies to all types of membranes present in the cell, that is, nuclear membranes and mitochondrial membranes.</p>
Mitochondrion	A semiautonomous, self-replicating organelle that occurs in varying numbers, shapes, and sizes in the cytoplasm of virtually all eukaryotic cells. It is notably the site of tissue respiration.
Nucleus	A membrane-bounded organelle of eukaryotic cells where chromosomes are housed and replicated. In most cells, the nucleus contains all of the cell's chromosomes except the organelle chromosomes and is the site of RNA synthesis and processing. In some species or in specialized cell types, RNA metabolism or DNA replication might be absent.
Spliceosome	A ribonucleoprotein complex containing RNA and small nuclear ribonucleoproteins (snRNPs), which is assembled during the splicing of messenger RNA primary transcript to excise an intron.
Protein complex	Any protein group composed of two or more subunits, which may or may not be identical. Protein complexes might have other associated non-protein prosthetic groups, such as nucleic acids, metal ions, or carbohydrate groups.
Ribosome	An intracellular organelle, about 200 angstroms in diameter, consisting of RNA and protein. It is the site of protein biosynthesis resulting from translation of messenger RNA (mRNA).

Table 22. GO Slim categories for cellular components (Sheet 4 of 4)

GO Slim cellular component	Description
Vacuole	A closed structure found only in eukaryotic cells, completely surrounded by unit membrane and containing liquid material. Cells contain one or several vacuoles that might have different functions from each other. Vacuoles have a diverse array of functions. They can act as a storage organelle for nutrients or waste products, as a degradative compartment, as a cost-effective way of increasing cell size, and as a homeostatic regulator controlling both the turgor pressure and the pH of the cytosol.
Organelle lumen	The volume enclosed by the membranes of a particular organelle, for example, endoplasmic reticulum lumen or the space between the two lipid bilayers of a double membrane surrounding an organelle (for example, nuclear membrane lumen).

GO Slim Categories for Biological Processes

[Table 23](#) describes the GO Slim categories for biological processes.

Table 23. GO Slim categories for biological processes (Sheet 1 of 3)

Go Slim biological process	Description
Cell communication	Any process that mediates interactions between a cell and its surroundings. Cell communication encompasses interactions such as signaling or attachment between one cell and another cell, between a cell and an extracellular matrix, or between a cell and any other aspect of its environment.
Cell death	The specific activation or halting of processes within a cell so that its vital functions markedly cease, rather than simply deteriorate gradually over time, which culminates in cell death.
Cell differentiation	The process where relatively unspecialized cells—for example, embryonic or regenerative cells—acquire specialized structural features, functional features, or both that characterize the cells, tissues, or organs of the mature organism or some other relatively stable phase of the organism's life history. Differentiation includes the processes involved in commitment of a cell to a specific fate.

Table 23. GO Slim categories for biological processes (Sheet 2 of 3)

Go Slim biological process	Description
Cell division	The processes resulting in the physical partitioning and separation of a cell into daughter cells.
Cell growth	The process by which a cell irreversibly increases in size over time by accretion and biosynthetic production of matter similar to that already present.
Cell homeostasis	The processes involved in the maintenance of an internal equilibrium at the level of the cell.
Cell motility	Any process involved in the controlled movement of a cell.
Cell organization and biogenesis	A process that is carried out at the cellular level and that results in the formation, arrangement of constituent parts, or disassembly of a cellular component. The process includes the plasma membrane and any external encapsulating structures, such as the cell wall and cell envelope.
Cell proliferation	The multiplication or reproduction of cells, resulting in the rapid expansion of a cell population.
Coagulation	The process by which a fluid solution, or part of it, changes into a solid or semisolid mass.
Conjugation	The union or introduction of genetic information from compatible mating types that results in a genetically different individual. Conjugation requires direct cellular contact between the organisms.
Defense response	Reactions triggered in response to the presence of a foreign body or the occurrence of an injury, which result in restriction of damage to the organism attacked or prevention and recovery from the infection caused by the attack.
Development	The biological process whose specific outcome is the progression of an organism over time from an initial condition (for example, a zygote or a young adult) to a later condition (for example, a multicellular animal or an aged adult).
Metabolic process	Processes that cause many of the chemical changes in living organisms, including anabolism and catabolism. Metabolic processes typically transform small molecules but also include macromolecular processes such as DNA repair and replication, and protein synthesis and degradation.

Table 23. GO Slim categories for biological processes (Sheet 3 of 3)

Go Slim biological process	Description
Regulation of biological process	Any process that modulates the frequency, rate, or extent of a biological process. Biological processes are regulated by many means, for example, control of gene expression, protein modification, or interaction with a protein or substrate molecule.
Reproduction	The production by an organism of new individuals that contain some portion of their genetic material inherited from that organism.
Response to stimulus	A change in state or activity of a cell or an organism (in terms of movement, secretion, enzyme production, gene expression, and so forth) as a result of a stimulus.
Transport	The directed movement of substances (such as macromolecules, small molecules, ions) into, out of, within, or between cells.

8 Obtaining Protein Annotation Information

GO Slim Categories

Searching for Post-Translational Modifications

This chapter describes how to search for peptides containing amino acids with post-translation modifications, such as phosphorylation or glycosylation, and how to interpret the results.

Contents

- [Using the ptmRS Node](#)
- [Using the Peptide in Protein Annotation and the PSM Grouper Nodes](#)
- [Viewing PTM Information on the Protein Identification Details View](#)
- [Filtering Phosphorylation Site Probabilities](#)

Using the ptmRS Node

The ptmRS node in the processing workflow provides a confidence measure for the localization of phosphorylation and any other modifications in peptide sequences. You can use it with all common fragmentation techniques, such as CID, ETD, and HCD, and peptide identifications from all available database search engines. This node uses parallelized calculations with multiple processor cores, if available, to improve performance. It also determines the optimal number of peaks to consider for localization of PTM sites for each m/z window individually, which increases the sensitivity of site localization for CID data. Depending on the applied fragmentation technique, the ptmRS node uses different fragment ion types for scoring to provide the highest sensitivity. For CID data, it scores only singly and doubly charged b and y ions. For analysis of HCD spectra, it also considers neutral loss ions. In contrast, when localizing PTM sites in ETD spectra, the node exclusively considers singly charged c, z, and y+H ions.


The input to ptmRS is a list of filtered MS/MS fragment spectra with corresponding peptide identifications.

How the ptmRS Node Calculates Sequence and Site Probabilities

The calculation of both ptmRS sequence probabilities and ptmRS site probabilities is based on the assumption that if the sequence of amino acids assigned to the respective MS/MS spectrum is correct and the MS/MS spectrum acquired did not result from fragmentation of two co-eluting peptides with similar m/z ratios, one of the putative isoforms must be the true assignment. Therefore, the sum of all sequence probabilities must equal 100 percent. If one or both of these assumptions is not correct, the calculated probabilities might not correlate with the true probabilities for correct modification localization.

[Figure 242](#) shows the basic internal workflow of the ptmRS node for phosphorylation.

Figure 242. Internal ptmRS workflow example for phosphorylation



The ptmRS node performs the following steps:

1. On the basis of all potential modification sites in the peptide sequence, ptmRS calculates all putative phosphorylation isoforms.
2. ptmRS divides the fragment spectrum into 100-Th windows and extracts the *i* most intense peaks for each window, where *i* is the peak depth used for the filtered fragment spectrum.

- For each phosphorylation isoform, S , and each peak depth, i , ptmRS calculates the list of theoretical fragment ions and matches them to the extracted peaks of the experimental spectrum.

The probability, p , of matching a fragment peak purely by chance is calculated as follows:

$$p = \frac{N_{Peaks} \times d}{w}$$

where:

- N_{Peaks} is the number of extracted peaks.
 - d is the mass tolerance for matching peaks to the theoretical fragment ions.
 - w is the extracted mass range.
- From the number of theoretical fragment ions, n , the number of matched fragment peaks, k_S , and the probability, p , ptmRS calculates the probability, $P_{S,i}$ of matching k_S or more peaks purely by chance as the cumulative binomial probability of matching k_S or more peaks in n attempts^{1,2}:

$$P_{S,i}(X \geq k_S) = \sum_{k=k_S}^n \binom{n}{k} p^k \times (1-p)^{(n-k)}$$

- ptmRS reports the binomial peptide score of a specific isoform at peak depth, i , as follows:

$$\text{pRS}_{S,i} = 10 \times \log_{10}(P_{S,i}(X \geq k_S))$$

- ptmRS determines the optimal peak depth $i_{optimal}$ as the peak depth with the largest difference between the best-scoring phosphorylation isoform and the second-best-scoring isoform¹.
- From the binomial peptide score at the optimal peak depth, ptmRS calculates an isoform confidence sequence probability and modification site probabilities. It calculates the ptmRS sequence probability as follows²:

$$\text{ptmRS sequence probability} = \frac{P_{S,i_{optimal}}^{-1}}{\sum_S P_{S,i_{optimal}}^{-1}}$$

¹ Beausoleil S.A., et al. *Nat. Biotech.* **2006**, Volume 24: 1285-1292.

² Olsen J.V., et al. *Cell.* **2006**, 127: 635-648.

Table 24 gives an example showing how ptmRS calculates sequence probabilities.

Table 24. Calculating ptmRS sequence probabilities

Putative isoform	pRS score	1/P value	pRS sequence probability
AM(pS)PALGVM(pS)FSGVQ	121	1.26×10^{12}	$1.26 \times 10^{12} / 1.89 \times 10^{12=0.67}$
AM(pS)PALGVMSF(pS)GVQ	118	6.31×10^{12}	0.33
AMSPALGVM(pS)F(pS)GVQ	59	7.94×10^{12}	0.00
		$\Sigma = 1.89 \times 10^{12}$	$\Sigma = 1.00$

8. The ptmRS node calculates the ptmRS site probability, for example, for a particular phosphorylation site by summing up the ptmRS sequence probabilities of those isoforms where the respective site is phosphorylated³. Table 25 gives an example showing how ptmRS calculates site probabilities.

Table 25. Calculating ptmRS site probabilities

pRS sequence probability	Putative isoform							
0.67	AM	(pS)	PALGVM	(pS)	F	S	GVQ	
0.33	AM	(pS)	PALGVM	S	F	(pS)	GVQ	
0.00	AM	S	PALGVM	(pS)	F	(pS)	GVQ	
ptmRS site probabilities			1.00			0.67	0.33	

Creating a PTM Analysis Workflow with the ptmRS Node

If you want to focus on studying the biologically relevant post-translational modifications of proteins, you can create a PTM analysis workflow. The processing workflow must include the ptmRS node (refer to the Help), which generates all PSM-specific modification site data. The Peptide in Protein Annotation node (refer to the Help) in the consensus workflow provides additional information about the modifications on the protein level. For in-depth PTM analysis, use both nodes together.


³ Olsen, J., et al. *Cell*. **2006**, *127*: 635–648.

The *ptmRS* node calculates modification site probabilities and makes them available in the Protein Identification Details view when you choose View > Protein Details on the Proteins page of the .pdResult file. This view color-codes the found modification above the amino acid sequences to indicate the probability of the modification being found on those portions of the amino acid. The PTM Site Probabilities area to the left of the sequence table displays a legend explaining the color-coding. For more information on this view, see [“Viewing PTM Information on the Protein Identification Details View”](#) on page 369.



For information on the parameters of the *ptmRS* node, refer to the Help.

❖ To create a basic PTM analysis processing workflow

1. Create or open a study and an analysis:
 - To create a study, see [“Creating a Study”](#) on page 38.
 - To open an existing study, see [“Opening an Existing Study”](#) on page 40.
 - To create an analysis, see [“Creating an Analysis”](#) on page 71.
 - To open an existing analysis, see [“Opening an Existing Analysis”](#) on page 72.
2. Create a basic processing workflow, as described in [“Creating a Processing Workflow”](#) on page 104, by using the following nodes:
 - Spectrum Files
 - Spectrum Selector
 - Search engine node, such as Sequest HT or Mascot
 - PSM Validation node, such as Fixed Value PSM Validator, Percolator, or Target Decoy PSM Validator
3. Attach a *ptmRS* node to the PSM validator node (in this case the Target Decoy PSM Validator node), as shown in [Figure 243](#).

Figure 243. Basic PTM analysis processing workflow

4. (Optional) Save the workflow:

- a. In the Name box above the Workflow Tree pane, type a name for the processing workflow.
- b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the processing workflow.
- c. In the Workflow Editor, click the **Save** icon,  **Save** , or the **Save Common** icon,  **Save Common** .
- d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Processing Workflow File (*.pdProcessingWF)**.
 - iii. Click **Save**.


The application saves the workflow in the *file_name*.pdProcessingWF file.



❖ To create a basic PTM analysis consensus workflow

1. Follow the general instructions for creating a consensus workflow in the Workflow Editor. See [“Creating a Consensus Workflow”](#) on page 112.
2. Include the following nodes in the consensus workflow:
 - MSF Files node
 - PSM Grouper node
 - Peptide Validator node
 - Peptide and Protein Filter node
 - Protein Scorer node
 - Protein Grouping node
 - Peptide in Protein Annotation node


Figure 244 shows the nodes required to create a basic PTM analysis consensus workflow.

Figure 244. Basic PTM analysis consensus workflow



3. (Optional) Add any other appropriate nodes.
4. Connect the nodes together. Connect the Peptide in Protein Annotation node to the Protein Grouping node, if it not already connected.
5. Set the parameters for each node.
6. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the processing workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the processing workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save** , or the **Save Common** icon,  **Save Common** .
 - d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Consensus Workflow File (*.pdConsensusWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdConsensusWF file.

7. Save the study. See “Saving a Study” on page 62.
8. Click the **Run** icon,  **Run** , in the upper right corner of the Analysis window.

The job queue appears, as shown in Figure 79 on page 112, displaying the status of your search as the search progresses. For information about the job queue, refer to the Help.

Columns in the Results Report

The PSMs page of the .pdResults report displays the primary results of the ptmRS processing. The node adds the following new columns to the report. These columns report the probabilities for single modifications, score, and isoform confidences.

- **ptmRS Best Site Probabilities** column: Displays the most likely positions of the modifications and their site score for each PSM. For each modification site, this value is an estimate of the probability (0–100%) that the site is truly modified. Any ptmRS site probabilities above 75% indicate that a site is truly modified. For peptide groups, this column shows the best site probabilities of the first PSM. This column is visible by default.
- **ptmRS *Modification* Site Probabilities**: Displays the modification site probabilities for all possible modification positions for a single modification. In the column title, *Modification* is replaced by the name of the actual modification, for example, Oxidation Site Probabilities.

9 Searching for Post-Translational Modifications

Using the ptmRS Node

- **ptmRS Binomial Peptide Score:** Reflects the quality of the match between a modification position isoform and the respective tandem spectrum. This peptide score is based on the cumulative binomial probability that the observed match is a random event. The value of the binomial peptide score heavily depends on the data scored, but usually scores above 50 indicate a good match. You can find details on the calculation of this score elsewhere.⁴
- **ptmRS Isoform Confidence Probability:** Displays the estimate of the probability (0–100%) that an individual isoform is correct. This estimate assumes that the search engine correctly identified the corresponding peptide in terms of amino acid sequence, as well as the number and identity of the modifications.⁴
- **ptmRS Isoform Group Report:** Displays the confidence estimate of an individual isoform. This estimate assumes that the search engine correctly identified the corresponding peptide in terms of amino acid sequence, as well as the number and identity of the modifications.

Figure 245 shows the columns that the ptmRS node adds to the PSMs page of the .pdResult report.

Figure 245. Columns added to the PSMs page by the ptmRS node

The screenshot displays the Proteome Discoverer software interface, specifically the PSMs page. The table shows various protein entries with their corresponding PSMs and the additional columns provided by the ptmRS node. The columns are: ptmRS: Binomial Peptide Score, ptmRS: Isoform Confidence Probability, ptmRS: Isoform Group Report, ptmRS: Best Site Probabilities, ptmRS: Methyl Site Probabilities, ptmRS: Oxidation Site Probabilities, and ptmRS: Phospho Site Probabilities. The table lists protein descriptions such as Histone H2A-Ill - Gallus gallus (Chicken), Histone H2A type 3 - Rattus norvegicus (Rat), and Histone H4.3 (HM4) - Zea mays (Maize). The status bar at the bottom indicates 45/136 Protein Groups, 1052 Proteins, 340 Peptide Groups, 674 PSMs, and 2893 MS/MS Spectrum Info.

⁴ Taus, T. et al., *J. Proteome Res*, 2011, 5354–5362

Figure 246 shows another example of the results of ptmRS processing on the PSMs page. The top of Figure 246 shows the ptmRS Best Site Probabilities column. The bottom of the figure shows all the ptmRS Modification Site Probabilities columns for oxidation, phosphorylation, and acetylation, along with the ptmRS Binomial Peptide Score column, the ptmRS Isoform Confidence Probability column, the ptmRS Isoform Group report column, and the ptmRS Best Site Probabilities column.

Figure 246. ptmRS modification site probability score columns on the PSMs page

Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	Master Protein Accessions	Protein Accessions	Activation Type	ptmRS [4]: Best Site Probabilities
<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	EPTSSEGGGLEGSSAAGI		Q04323	Q04323	HCD	
<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	NmGGPYGGGNYGPGGSG	M2(Oxidation)	O88569	O88569; Q2HJ60; Q8TTV2; I	HCD	M2(Oxidation): 100
<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	ESEDKPEIEDVGSDEEEEK*	S13(Phospho)	P07900	P07900	HCD	S13(Phospho): 100
<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	NmGGPYGGGNYGPGGSG	M2(Oxidation)	O88569	O88569; Q2HJ60; Q8TTV2; I	HCD	M2(Oxidation): 100
<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	SLAGSSGPGASSGTSGDHK		P29692	P29692	HCD	
<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	GPPQEEEEEEEEEEATK		P51858	P51858	HCD	

ptmRS Modification Site Probabilities columns

ptmRS [4]: Binomial Peptide Score	ptmRS [4]: Isoform Confidence Probability	ptmRS [4]: Isoform Group Report	ptmRS [4]: Best Site Probabilities	ptmRS [4]: Phospho Site Probabilities	ptmRS [4]: Oxidation Site Probabilities	ptmRS [4]: Acet
257	100.0 %		M2(Oxidation): 100		M(2): 100.0	
205	100.0 %	y1+--H(3) O(4) P(15): 1714.77: y	S13(Phospho): 100	S(2): 0.0; S(13): 100.0		
254	100.0 %		M2(Oxidation): 100		M(2): 100.0	

On the Peptide Groups page, the PSM Grouper node places the modification site scores in the Modifications column, as shown in Figure 247. The application uses the modification site scores to select the best modification sites.

Figure 247. ptmRS modification site probability scores in the Modifications column on the Peptide Groups page

Checked	Confidence	PSM Ambiguity	Annotated Sequence	Modifications	Modifications (all positions)	Modifications in Master Proteins
<input type="checkbox"/>	●	Selected	[K] LATLAASRSNCMK [R]	1xPhospho [T3(99.8)]	1xPhospho [T3(99.8); S7(0.2); S9(0)]	P06773 1xPhospho [T172(99.8); S176(0.2); S178(0)]
<input type="checkbox"/>	●	Rejected	[K] AFFESHFAPSAER [T]			
<input type="checkbox"/>	●	Unambiguous	[K] AAEAAAAPAESAAPAAAGEEPSK			
<input type="checkbox"/>	●	Unambiguous	[R] NTTQNTGYSSGTQANYPVR [I]			
<input type="checkbox"/>	●	Unambiguous	[K] KLGEMWSEQSAK [D]			
<input type="checkbox"/>	●	Unambiguous	[K] ESEDKPEIEDVGSDEEEEEK [D]	1xPhospho [S13(100)]	1xPhospho [S2(0); S13(100)]	P07900 1xPhospho [S252(0); S263(100)]
<input type="checkbox"/>	●	Unambiguous	[K] DVTNNVHYENYR [CS]			
<input type="checkbox"/>	●	Unambiguous	[K] TPAQYDASELK [A]			
<input type="checkbox"/>	●	Unambiguous	[K] YGLQDSDEEEEHPSK [T]	1xPhospho [S6(100)]	1xPhospho [Y1(0); S6(100); S15(0)]	P52948 1xPhospho [Y883(0); S888(100); S897(0)]
<input type="checkbox"/>	●	Unambiguous	[R] LESGMQNSIHTK [T]			
<input type="checkbox"/>	●	Unambiguous	[R] DGMNDQGGYGSVGR [M]			
<input type="checkbox"/>	●	Unambiguous	[K] GVVEVTHDLQK [H]			
<input type="checkbox"/>	●	Unambiguous	[K] QEYDESGPSIVHR [K]			
<input type="checkbox"/>	●	Unambiguous	[R] ILEQEEEEQAGKPGPSK [K]			
<input type="checkbox"/>	●	Unambiguous	[R] KPLLESGLTGK [G]			
<input type="checkbox"/>	●	Unambiguous	[K] VLQHYQESDKGEELGPGNVQK			
<input type="checkbox"/>	●	Unambiguous	[K] ITDSAGHILYSK [E]			
<input type="checkbox"/>	●	Unambiguous	[R] SSGSPYGGGYSGGGGGYGYG			
<input type="checkbox"/>	●	Unambiguous	[R] SSGSPYGGGYSGGGGGYGYG	1xPhospho [S4(100)]	1xPhospho [S1(0); S2(0); S4(100); Y6(0); Y10(0); S12(0); S16(0)]	Q8BG05 1xPhospho [S356(0); S357(0); S359(100); Y36
<input type="checkbox"/>	●	Unambiguous	[R] STAGDTHLGGEDFNR [LM]			
<input type="checkbox"/>	●	Selected	[KR] KEEELQAALAR [VLA]			
<input type="checkbox"/>	●	Rejected	[R] KSSELKSPK [Q]	3xAcetyl [S2(99.9); S3(100); S7	3xAcetyl [K1(0.1); S2(99.9); S3(100); K6(0.1); S7(99.9); K10(0)]	Q6R2W3 3xAcetyl [K798(0.1); S799(99.9); S800(100); K
<input type="checkbox"/>	●	Unambiguous	[K] HTVDDGLDIR [K]			

On the Proteins page, the Peptide in Protein Annotation node creates modification site scores in the Modifications column, as shown in Figure 248.

9 Searching for Post-Translational Modifications

Using the ptmRS Node

Figure 248. ptmRS modification site probability scores in the Modifications column on the Proteins page

Checked	Master	Accession	Description	Coverage	# Peptides	# PSMs	# Unique Peptides	# Protein Groups	# AAs	MW [kDa]	calc. pI	Modifications
<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q9TTV2	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A	7%	1	3	1	1	341	35.9	8.98	Oxidation [M315(100)]
<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q5RBU8	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A	7%	1	3	1	1	353	37.4	8.97	Oxidation [M327(100)]
<input type="checkbox"/>	<input checked="" type="checkbox"/>	O88569	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A	7%	1	3	1	1	353	37.4	8.95	Oxidation [M327(100)]
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P22626	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A	7%	1	3	1	1	353	37.4	8.95	Oxidation [M327(100)]
<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q2HJ60	Heterogeneous nuclear ribonucleoproteins A2/B1 (hnRNP A	7%	1	3	1	1	341	36.0	8.65	Oxidation [M315(100)]
<input type="checkbox"/>	<input checked="" type="checkbox"/>	P07900	Heat shock protein HSP 90-alpha (HSP 86) (Renal carcinom	6%	3	4	3	1	732	84.6	5.02	Phospho [S252(0);S263(100)]

In the Protein Identification Details view, the application displays the modification site probabilities calculated by the ptmRS node as color-coded modification abbreviations above the sequence, as shown in Figure 249.

Figure 249. Modification site probabilities shown as color-coded modification abbreviations above the sequence in the Protein Identification Details view

Protein Identification Details

Coverage: 47.06%

Found Modifications:

- D Dimethyl (K,R)
- M Methyl (K,R)
- O Oxidation (M)

PTM Site Probabilities:

Threshold: 25%

25 - 45%
45% - 75%
75% - 99%
99% - 100%

Sequence: M A R T K Q T A R K S T G G K A P R K Q L A T K A A R K S A P T T G G V G G P H R Y R P C T V A L R E I R R Y Q K S T E L L I R K L P F Q R L V R E I A Q D F K T D L R F Q S A A I G A L Q E A S E A Y

Modifications: Q10453

101

Modifications: Q10453

L V L G L F E D T N L C A I H A K R V T I M P K D I Q L A R R I R G E R A

Methylation

Methylation is a PTM that can occur more than once on certain amino acids. For example, lysine might be modified by one, two, or three methyl groups. Search engines like Sequest HT or Mascot use three different modifications (methyl, dimethyl, and trimethyl) to denote these three modifications.

When there is more than one possible modification target, the methyl groups might be distributed differently over the sequence. For example, a peptide with a dimethyl modification on lysine might actually carry two single methyl groups on two lysine residues next to each other (both peptides have the same precursor mass in the MS1 scan). The ptmRS node

calculates site probabilities for these different modification isoforms. After the *ptmRS* modification site scoring, a peptide with a dimethyl modification identified by the search engine might have two single methyl modifications. The *ptmRS* scores are the preferred criteria for the modification display for peptide groups and proteins.

Using the Peptide in Protein Annotation and the PSM Grouper Nodes

Proteomics results contain data that is determined by the position of an identified peptide in a protein sequence, such as the position of the identified peptide in a protein itself, derived data like the flanking residues, and the positions of modifications in the protein sequence.

If you search for proteins with a given modification, you must have the modification information (that is, type and position) directly available on the Proteins page of the .pdResult file. If your research focuses on PTMs, you might work mainly with the Peptide Groups page. In this case, you might want to know the position of the modification in the proteins containing the peptide. You might also need to identify the flanking residues when the digestion of the protein in the sample preparation breaks potential modification sites into two pieces that you can no longer identify.

Consider an example involving flanking residues. Enzymes add biologically relevant PTMs to the protein inside the cell. These enzymes bind to the target protein only on a defined short sequence called the motif. Suppose that you are looking for potential glycosylation sites with the NxS/T sequence motif. The glycosylation occurs only in the local sequence *Asn-Xaa-Ser/Thr/Cys*, where *Xaa* is a random amino acid except proline. *Xaa* might be arginine or lysine and therefore a cleavage site of trypsin. If the glycosylated peptide contains two asparagine residues, one in the middle of the sequence and one at the position before the last amino acid, you must use the flanking residues of the peptide to decide whether the asparagine residue at the end of the sequence is a valid glycosylation site. If so, the flanking residue must be either serine, threonine, or cysteine.

You can use the Peptide in Protein Annotation node in the consensus workflow to add columns to the Proteins and Peptide Groups pages of the results report that display the positions of modifications in the protein sequence. For detailed information about the parameters of these nodes, refer to the Help.

Peptide Positions in Proteins

The Protein Positions for Peptides parameter of the Peptide in Protein Annotation node in the consensus workflow determines whether the Proteome Discoverer application displays the peptide or PSM position for all proteins or only for the master proteins. When you set this parameter to For All Proteins, the application adds a Position in Proteins column to the Peptide Groups page that displays the position of the peptides or PSMs for all proteins that it finds the peptide or PSM in. It displays these positions in the following form:


protein accession [*start - end*]

9 Searching for Post-Translational Modifications

Using the Peptide in Protein Annotation and the PSM Grouper Nodes

When you set the parameter to Only for Master Proteins, the application adds a Positions in Master Proteins column to the Peptide Groups page that displays the position of the peptides or PSMs only for the master protein of each protein group. [Figure 250](#) shows this column.

Figure 250. Positions in Master Proteins column on the Peptide Groups page



	Checked	Confidence	Annotated Sequence	Modifications (all PTM positions)	Modifications in Master Proteins	Master Protein Accessions	Positions in Master Proteins	XCorr Sequest HT	Confidence Sequest HT
1	<input type="checkbox"/>	●	[K]LDATVHGEVSSK [F]			Q17770	Q17770 [81-92]	2.38	●
2	<input type="checkbox"/>	●	[K]GHYTEGAELVDNVLDIR [RK]			P52275	P52275 [104-121]	6.09	●
3	<input type="checkbox"/>	●	[R]VFEVSLGDLNNSEADFR [K]			P48154	P48154 [64-80]	2.54	●
4	<input type="checkbox"/>	●	[K]FTDVPPVFLDGK [K]			P34255	P34255 [233-244]	2.92	●
5	<input type="checkbox"/>	●	[K]VLGTDELYEAR [Y]			P18334	P18334 [247-259]	2.95	●
6	<input type="checkbox"/>	●	[K]FDGLVLSGLPSASSELSQSR [I]			P05690	P05690 [34-53]	4.55	●
7	<input type="checkbox"/>	●	[R]AVAGLSLDAILAK [R]			O01868	O01868 [84-96]	2.88	●
8	<input type="checkbox"/>	●	[K]NAQEELAENVVGLRR [E]			P02566	P02566 [1498-1512]	3.12	●
9	<input type="checkbox"/>	●	[K]AYLDFVVR [H]			Q18115	Q18115 [314-321]	2.08	●
10	<input type="checkbox"/>	●	[R]VLAQGIHFSASNAFLYETVR [R]			P18948	P18948 [841-861]	3.54	●
11	<input type="checkbox"/>	●	[K]ALEILPTLVVK [E]			Q9N358	Q9N358 [138-148]	2.24	●
12	<input type="checkbox"/>	●	[R]IHFPLAAYTPLISAEK [A]			P34690	P34690 [263-278]	3.29	●
13	<input type="checkbox"/>	●	[K]ITGDFEDHVQSDIVAFNKI [-]			Q9U2H9	Q9U2H9 [244-263]	4.78	●

On the Peptide Groups page, there is often more than one protein for each peptide. You can click Show Associated Tables to display tables that show distinct peptide-protein relations, so the position of the peptide in the protein is unambiguous and displayed without protein accession. If a peptide occurs more than once in a protein, the tables list both ranges.

Flanking Residues

The flanking residues of a peptide are the amino acids that border a cleavage site. The Annotate Flanking Residues of the Peptide parameter of the Peptide in Protein Annotation node determines whether the Proteome Discoverer application annotates the amino acids before and after the peptide in the protein sequence. If you set this parameter to True, the application replaces the Sequence column on the Peptide Groups page by an Annotated Sequence column that displays the one-letter code representing the flanking residues of a peptide in a protein. [Figure 251](#) shows the Annotated Sequence column on the Peptide Groups page.

Figure 251. Annotated Sequence column on the Peptide Groups page

#	Checked	Confidence	Annotated Sequence	Modifications	Modifications in Proteins	# Protein Groups	# Proteins	# PSMs	Master Protein Accessions	Positions in Master Protei
1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] LGEMWNTAADDKQPYEK [K]			1	5	1	P10103	P10103 [129-146]
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] HPDSSVNFAEFSK [K]			1	4	1	P17741	P17741 [31-43]
3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] KTVTAMDVYALKR [Q]			1	19	4	Q6LAF1; P84043	Q6LAF1 [80-93]; P84043
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] KTVTAMDVYALKR [Q]		P02310 1×Oxidation [M85(100)]; Q41811 1×Oxidati	1	19	6	Q6LAF1; P84043	Q6LAF1 [80-93]; P84043
5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] KNPEVPVNFAEFSK [C]			1	4	2	P40618	P40618 [30-44]
6	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] ADREDESSPYAALAAQDVAQR [C]		Q9CAX6 1×Oxidation [M74(100)]; P19950 1×Oxidati	1	6	1	P13471	P13471 [64-84]
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] SKEEAEALYHSKYEEELQVTVGR [H]			1	1	1	P35908	P35908 [369-390]
8	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] KTVTAMDVYALKR [Q]			1	1	2	P27996	P27996 [80-93]
9	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] LLELDPPPKDVEDR [G]			1	1	1	P19661	P19661 [57-70]
10	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] KEEDPKTASQLLVNLR [A]			1	2	1	Q60611	Q60611 [411-427]
11	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] SAGKYDYYLNPQGG [A]			1	4	1	Q00566	Q00566 [116-130]
12	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] VQASLAANTFTITGHAEIK [Q]			1	2	1	P20290	P20290 [123-141]
13	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] SKAEAEALYHSKYEEELQVTVGR [H]			1	1	1	P04264	P04264 [365-386]
14	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] NQPGKYSQLVETIR [R]			1	1	2	Q92522	Q92522 [43-57]
15	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] FADEHVFGSPFTVK [I]			1	2	1	O75369	O75369 [2075-2088]
16	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] IKDDPEDDGYFAPPKEDIKPLK [R]	1×Methyl [K2(100)]	Q04750 1×Methyl [K119(100)]	0	1	1	Q04750	Q04750 [118-139]
17	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] LILPGELAKHAVSEGTK [A]			2	88	6	P07794; P10854; Q96A08; POC1H3; Q8CGP0	P07794 [99-115]; P10854
18	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] AQTLPSTVVTITSESSPGKR [E]			1	2	1	Q01082	Q01082 [2326-2345]
19	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] TPEVDEALEKFDK [A]			1	2	1	P02755	P02755 [143-156]
20	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] KESYSIYYKLVK [Q]			1	14	1	POC1H3; Q8CGP0	POC1H3 [35-47]; Q8CGP
21	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] LVREIAQDFKTLR [F]			2	32	1	P68429; Q10453; P59169; P68432	P68429 [71-84]; Q10453
22	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] LAREIAQDFKTLR [F]	1×Dimethyl [R3(100)]	Q5RCC9 1×Dimethyl [R73(100)]	1	1	1	Q5RCC9	Q5RCC9 [71-84]
23	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] TYFPHFDLSHGSGQIK [A]			0	1	1	P07414	P07414 [41-56]
24	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] TYFPHFDLSHGSGQVK [GA]	1×Methyl [K16(100)]	P01954 1×Methyl [K56(100)]; P01936 1×Methyl [K5	0	2	1	P01936	P01936 [41-56]
25	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] TYFPHFDLSHGSAQVK [GAD]			1	70	1	P01966	P01966 [42-57]
26	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] SLETENSALQLQVTEREEVVR [G]			1	3	1	P14733	P14733 [53-72]
27	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] R] TVTAMDVYALKRQGR [T]		P02310 1×Oxidation [M85(100)]; Q41811 1×Oxidati	1	20	1	Q6LAF1; P84043	Q6LAF1 [81-96]; P84043
28	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] AQNDLIMWDELKIK [A]			1	1	1	P18493	P18493 [243-257]
29	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] ITRFRPDSADGMLLYNGQK [R]		P98160 1×Oxidation [M3700(100)]; Q05793 1×Oxide	1	2	1	P98160	P98160 [3690-3707]
30	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] LAAAILGGVDQHHIKPGAK [V]			1	3	1	P22509	P22509 [150-168]
31	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] VDVLVNLPGHEHGQR [L]			1	3	2	P04937; P07589	P04937 [938-952]; P0758
32	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] LASVPAGGAVAVSAAPGSAAPAAGSAI			1	4	1	P05387	P05387 [62-94]
33	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] IKDEPEDDGYFAPPKEDIKPLK [R]			1	1	1	Q9WUL0	Q9WUL0 [118-139]
34	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[R] GENLVSMTVEGPPPKDTGIAR [V]		P27048 1×Oxidation [M80(100)]; P17136 1×Oxidati	1	10	1	P14678	P14678 [74-94]
35	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] VKEPSVGEATSTSDILK [V]			1	2	1	Q5R5Q2	Q5R5Q2 [230-246]
36	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] KLGEWNTAADDKQPYEK [K]		Q9UGV6 1×Oxidation [M132(100)]; P63158 1×Oxide	1	5	2	P10103	P10103 [128-146]
37	<input type="checkbox"/>	<input checked="" type="checkbox"/>	[K] KHPDASVNFAEFSK [C]			1	4	4	P10103	P10103 [130-141]

You can click Show Associated Tables to display pages that show the positions of the peptide in the proteins and the flanking residues. You can use these pages to obtain detailed information about the flanking residues of the peptide in different proteins. To specify the number of flanking residues to display in the associated pages, use the Number Flanking Residues in Connection Tables parameter of the Peptide in Protein Annotation node. You can display up to five flanking amino acids in the associated pages.

The application displays the flanking residues in brackets before and after the peptide sequence and separates them by a dot in the sequence, as shown in the first sequence in Figure 252. If the application finds a peptide more than once in a protein or if it finds a peptide in more than one protein, it might display more than one flanking residue. In these cases, it displays all possible residues within a pair of brackets, which are commonly used to display a sequence consensus.


Figure 252 shows different flanking residues for two proteins in brackets. The sequences compared in the bottom of the view show the residues, which are identical, for the two proteins, in green. The flanking residues are different on both sides of the peptide, so there are gaps (since different residues are not aligned) or alternative shading to indicate the difference between the two sequences.

9 Searching for Post-Translational Modifications

Using the Peptide in Protein Annotation and the PSM Grouper Nodes

Figure 252. Peptide group consensus sequence with different flanking residues from two proteins

Consensus sequence showing flanking residues in brackets



Full protein sequences of the different flanking residues shown at top in brackets

Modifications of Proteins on the PSMs Page

The Proteome Discoverer application displays the positions of the identified modifications in a protein the same way that it displays the positions of a peptide in a protein. The way it displays peptide modifications is described in the Help.

You can display modification positions for all proteins or only for master proteins.

If you select the Peptide in Protein Annotation node in the consensus workflow and set its Protein Modifications Reported parameter to For All Proteins, the application adds a Modification in Proteins column to the Peptides Group page that displays the positions of the modifications in all proteins.

If you set the node's Protein Modifications Reported parameter to For Master Proteins, the application adds a Modification in Master Proteins column to the Peptides Group page that displays the modifications in only the master proteins.

The application always creates two different modification columns. One modification column contains only the best possibilities for the modifications. The other contains all possible positions. The Modification Sites Shown parameter of the PSM Grouper node determines whether the application shows one or both of these columns. The column not shown is still available, and you can display it by using the Field Chooser (refer to the Help). For detailed information on the parameters of the PSM Grouper node, refer to the Help and the decision tree for the modification display (Figure 255 on page 366) and the accompanying examples.

You can also use the PSM Grouper node's Site Probability Threshold parameter to determine whether the application should display only sites with a minimal site probability score.

Modifications of Proteins on the Proteins Page

The Modification Sites Reported parameter of the Peptide in Protein Annotation node determines which modifications the application annotates on the Proteins page of the .pdResult file. These are the options:

All Combined—The application annotates all modifications and displays the annotations in a single column called Modifications.

Only Specific—The application annotates only the modifications specified by the modification parameters that reside below the Modification Sites Reported parameter in the same section and displays two columns for each specified modification as follows:

- The first column, named by the modification name, shows the positions of all identified modifications of this kind in the protein.
- The second column, also named by the modification name plus an additional count, contains the total count of this kind of modification in each protein.

All and Specific—The application annotates all identified modifications of the protein and displays the annotations in a single column called Modifications. It also annotates the specified modifications and displays the positions and count in single columns.

You can restrict the concatenated list of modifications in the Modification column on the Proteins page. Use the Site Probability Threshold parameter of the PSM Grouper node and the Report Only PTMs parameter of the Peptide in Protein Annotation node to show only modification sites above a given probability threshold or modifications that are PTMs.

You can also use the Peptide in Protein Annotation node's Minimum PSM Confidence parameter to restrict the modifications shown to those with a minimal confidence level of the PSMs that identified the modifications.

Depending on the protein, the list of modifications can become long. If you are interested in only one type of modification, use the Modification parameter (under Modifications in Protein) to display a set of five (maximum) selected modifications in extra columns. For each of these specified modifications, the application creates two new columns. One column contains the modifications and all sites where the modification is found, and the second column shows how many times this modification occurs in the protein.

9 Searching for Post-Translational Modifications

Using the Peptide in Protein Annotation and the PSM Grouper Nodes

Using phosphorylation as an example, [Figure 253](#) shows in the Modifications column all modifications and the sites where the modifications are found. It also shows the number of modifications in the Phospho Count column and the separate modifications in the Phospho Positions column.

Figure 253. Phosphorylation modifications on the Proteins page

Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info							
	Checked	Master	Accession	Description	Coverage	# Peptides	Modifications	Phospho(S.T.Y) Count	Phospho(S.T.Y) Positions	Score A2	# Peptides A2	
1	<input type="checkbox"/>	<input type="checkbox"/>	B4DMA2	cDNA FLJ54023, highly similar to Heat shock protein HSP 90-beta OS=Homo sapiens GN=HS	32%	23	Phospho [S188.S217]	2	Phospho [S188.S217]	169.65	23	
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	P08238	Heat shock protein HSP 90-beta OS=Homo sapiens GN=HS	30%	23	Phospho [S226.S255]	2	Phospho [S226.S255]	169.65	23	
3	<input type="checkbox"/>	<input type="checkbox"/>	B4DGLD	cDNA FLJ53619, highly similar to Heat shock protein HSP 90-alpha OS=Homo sapiens GN=H	29%	21	Phospho [S216.S245]	2	Phospho [S216.S245]	164.71	21	
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	P07900	Heat shock protein HSP 90-alpha OS=Homo sapiens GN=H	28%	19	Phospho [S231]	1	Phospho [S231]	119.38	19	
5	<input type="checkbox"/>	<input type="checkbox"/>	Q6PK50	HSP90AB1 protein (Fragment) OS=Homo sapiens GN=HSF	39%	13	Phospho [S226.S255]	2	Phospho [S226.S255]	118.00	13	
6	<input type="checkbox"/>	<input type="checkbox"/>	Q2VVPJ6	HSP90AA1 protein (Fragment) OS=Homo sapiens GN=HSF	33%	18	Phospho [S231]	1	Phospho [S231]	114.10	18	
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Q5R208	Carbamoylphosphate synthetase I OS=Homo sapiens GN=H	14%	19	Phospho [S913]	1	Phospho [S913]	101.93	19	
8	<input type="checkbox"/>	<input type="checkbox"/>	Q5R210	Carbamoylphosphate synthetase I OS=Homo sapiens GN=H	13%	18		0		99.79	18	
9	<input type="checkbox"/>	<input type="checkbox"/>	P31327	Carbamoyl-phosphate synthase [ammonia], mitochondrial C	13%	18		0		99.79	18	
10	<input type="checkbox"/>	<input type="checkbox"/>	Q5R206	Carbamoylphosphate synthetase I OS=Homo sapiens GN=H	13%	18		0		99.79	18	
11	<input type="checkbox"/>	<input type="checkbox"/>	J3KQL0	Carbamoyl-phosphate synthase [ammonia], mitochondrial C	13%	18		0		99.79	18	
12	<input type="checkbox"/>	<input type="checkbox"/>	A8K134	cDNA FLJ75245, highly similar to Homo sapiens carbamoyl	13%	18		0		99.79	18	
13	<input type="checkbox"/>	<input type="checkbox"/>	Q59HF8	Carbamoyl-phosphate synthetase 1 variant (Fragment) OS=	13%	18		0		99.79	18	
14	<input type="checkbox"/>	<input type="checkbox"/>	Q8TBA7	HSP90AA1 protein (Fragment) OS=Homo sapiens GN=HSF	25%	14	Phospho [S137]	1	Phospho [S137]	94.45	14	
15	<input type="checkbox"/>	<input type="checkbox"/>	Q1KLZ0	HCG15971, isoform CRA_a OS=Homo sapiens GN=PS1TF	31%	11		0		85.30	11	

Modifications of Peptides on the Peptide Groups Page

When you use the PSM Grouper node in the consensus workflow, the Proteome Discoverer application adds columns of information about modifications to the Peptide Groups page of the .pdResult file. For information on the Peptide Groups page, refer to the Help.

The content of the Modifications column shown in the Peptide Groups page is a merged list of all the modifications found in the PSMs belonging to the peptide group. [Figure 254](#) shows the format of each kind of modification.

Figure 254. Modification format

2xPhospho[S4(98);T7(97.6)]

Number of modifications of this type in the peptide

Type of modification as shown in the Chemical Modifications view

Target amino acid and position in the peptide

Site probability, if available

The PSM Grouper node creates two different modification columns. One column contains all possible modification sites, and the other one contains only the best possible site for each modification. If a modification has no unique “best” position, the application displays the modification with the target amino acid but without specifying the best position. In this case, you must review the PSMs of the peptide group by clicking Show Associated Tables on the results pages. Depending on the setting of the Modification Sites Shown parameter of the


node, the application displays one or both of these columns. If one of the columns is invisible, you can display it by using the Field Chooser in the upper left corner of the result page. You can restrict the display of modifications in the Peptide Groups page through the parameter settings of the PSM Grouper node. For information on these parameters, refer to the Help.

[Figure 255](#) shows the decisions that determine how the application displays a modification on the Peptide Groups page.

9 Searching for Post-Translational Modifications

Using the Peptide in Protein Annotation and the PSM Grouping Nodes




Figure 255. Decision tree determining the display of modifications on the Peptide Groups page



The following examples illustrate how the application displays different modifications.

Example 1

In this example, there is one modification with different positions and different confidences.

PSM number	PSM confidence level	Modified PSM sequence
PSM 1		APEPTIDEWITHM(O _x)ANYMETHIQNIN
PSM 2		APEPTIDEWITHM(O _x)ANYMETHIQNIN
PSM 3		APEPTIDEWITHMANYM(O _x)ETHIQNIN

The application displays the modification as 1x Oxidation[M13; M17] in the Modifications (All Possible Sites) column of the Peptide Groups page and 1x Oxidation[M13] in the Modifications column, which displays the best modification sites.

Example 2

In this example, the site probabilities are available, and enough scores are above the threshold.

The peptide has two phosphorylations (from the precursor mass).

PSM number	Modified PSM sequence
PSM 1	APEPTIDES(Phospho 80.0)WITHS(Phospho 80.0)ERANDT(Phospho 40.0)HR Binomial Score: 112
PSM 2	APEPTIDES(Phospho 98.0)WITHS(Phospho 76.0)ERANDT(Phospho 36.0)HR Binomial Score: 42

The application displays the modifications as 2x Phospho[S9(80.0) and S14(80.0)] in the Modifications column and 2xPhospho [S9(80.0), S14(80.0), T20(40)] in the Modifications (All Possible Sites) column.

Example 3

In this example, the site probabilities are available, but not enough scores are above the threshold.

A peptide has two phosphorylations (from the precursor mass). The Site Probability Threshold parameter of the PSM Grouper node is set to 75.0. See the next table.

9 Searching for Post-Translational Modifications

Using the Peptide in Protein Annotation and the PSM Grouping Nodes




PSM number	Modified PSM sequence
PSM 1	APEPTIDES(Phospho 80.0)WITHS(Phospho 60.0)ERANDT(Phospho 60.0)HR Binomial Score: 112
PSM 2	APEPTIDES(Phospho 98.0)WITHS(Phospho 66.0)ERANDT(Phospho 35.0)HR Binomial Score: 42

The application displays the modifications as 2x Phospho[S9(80.0); S/T] in the Modifications column.

One phosphorylation has an 80% probability of being located on serine 9. The other phosphorylation might be located on any serine or threonine in the sequence. If no threonine is present, the string becomes phospho[S9(98.0); S]. The Modifications (All Possible Sites) column contains 2xPhospho[S9(80);S14(60);T20(60)].

Example 4

In this example, the position of the modifications with the highest PSM confidence is not unambiguous.

PSM number	PSM confidence level	Modified PSM sequence
PSM 1		APEPTIDEWITHM(O _x)ANYMETHIQNIN
PSM 2		APEPTIDEWITHMANYM(O _x)ETHIQNIN
PSM 3		APEPTIDEWITHMANYM(O _x)ETHIQNIN




The sites with highest PSM probability are M13 and M17.

The mass of the precursor ion indicates that the peptide contains one modification. In the list of PSMs are two high-confidence (green) PSMs where the modification is in different positions. The application cannot determine which position of the modification is the correct one, so it reports no position.

The application displays the modifications as Oxidation[M] in the Modifications column and Oxidation[M13; M17] in the Modifications (All Possible Sites) column.

Example 5

In this example, the position of modifications with the highest PSM confidence is unambiguous.

PSM number	PSM confidence level	Modified PSM sequence
PSM 1		APEPTIDEWITHM(O _x)ANYMETHIQNIN
PSM 2		APEPTIDEWITHM(O _x)ANYMETHIQNIN
PSM 3		APEPTIDEWITHMANYM(O _x)ETHIQNIN

The site with the highest PSM probability is M13.

The mass of the precursor ion indicates that the peptide contains one modification. All high-confidence PSMs report the modification in the same position. In this case, there is a reasonable certainty that the modification is on the position identified in the high-confidence matches, so the application reports it.

The application displays the modifications as 1x Oxidation[M13] in the Modifications column and 1x Oxidation [M13; M17] in the Modifications (All Possible Sites) column.

Viewing PTM Information on the Protein Identification Details View

To view information on the Proteins page about the modifications found in the search, use the Protein Identification Details view. The following topics explain how to use this view to investigate PTMs. For basic information about the Protein Identification Details view, refer to the Help.

❖ To display a protein sequence

1. Open the .pdResult file.
2. Double-click a grid cell on the Proteins page.

–or–

Select a grid cell and choose **View > Show Protein Details**.

–or–


Click the **Show Protein Details** icon, .

The Protein Identification Details view appears, as shown in [Figure 256](#).

9 Searching for Post-Translational Modifications

Viewing PTM Information on the Protein Identification Details View

Figure 256. Protein Identification Details dialog box



Protein Identification Details

Coverage ProteinCard

Heat shock protein HSP 90beta

Annotate PTMs reported in Uniprot
 Show only PTMs
 Include PSMs that are filtered Out

Coverage: 64.50%

Found Modifications:

C Carbamidomethyl (C)
O Oxidation (M)
P Phospho (S)
T TMT6plex (N-term.K)

Sequence Modification List

Modifications	1	11	21	31	41	51	61	71	81	91
IPI00414676.6	MPEEVHGGEE	EVETFAFQAE	IAQLMSLIIN	TFYSNKEIFL	RELISNASDA	LDKIRVESLT	DPSKLDSCKE	LKIDIIPNPO	ERTLTUVDTG	IGMTKADLIN
IPI00414676.6	NLGTIARSGT	KAFMEALQAG	ADISMIGQFG	VGFSAYLVA	EKVVVITIKHN	DDEQVAMESS	AGGSFVTRAD	HGEPIGRGTV	VILHLKEDQT	EYLEERRVKE
IPI00414676.6	VVKKHSQFIG	YPITLYLEKE	REKESIDDEA	EEENGEKEEE	DRDDEEKPKI	EDVGSDEEDD	SGKDRKKKTK	KIKERKIDQE	ELNKTQPIWT	RNPDDITQEE
IPI00414676.6	YGEFYKSLTN	DMEDHLAVKH	FSVEQQLQFR	ALLFIPRRAP	FDLFFENKWK	NNIKLYVRV	FIMDSCELI	PEYLNFRIGV	VDSDELPLNI	SREMLQQSKI
IPI00414676.6	LKVIRKNIUK	KCLELFSELA	EDKENYKIFY	EAFSKNLKLG	IHEDSTNRRR	LSELLRYHTS	QSGDEMTSLS	EYVSRMKEIQ	KSIYYITGES	KEQVANSASFV
IPI00414676.6	ERVVRKRGFEV	VVMTEPIDEY	CVQQLKEFDG	KSLVSVTKEG	LLEPEDEEEK	KHMEESKAKF	ENLCKLMKEI	LDKKVEKVII	SNRLVSSPCC	IVTSTYGVTA
IPI00414676.6	NMERIMKAQA	LRDNSTMGYM	MAKKHLEINP	DHPIVETLRQ	KAEADKNDKA	VKDLVLLFE	TALLSSGFSL	EDPQTHSNRI	YRMIKLGLGI	DEDEVAAREP
IPI00414676.6	NAAVPDEIPP	LEGEDASRM	EEVD							

OK Help

The Protein Identification Details view features two pages, Coverage and ProteinCard.

- For detailed information on the Coverage page, refer to the Help.
- For detailed information on the ProteinCard page, see [“Accessing the ProteinCard Page”](#) on page 326.

❖ To display a protein's PTMs

1. In the Protein Identification Details view, click the **Coverage** tab if it is not already selected.
2. Click the **Modification List** tab.

The Modification List page displays all identified modifications in the sequence view. An example of the Modification List page is shown in the Help.

For information on the columns in this page, refer to the Help.

3. In the Number of Amino Acids Before and After Modification box, specify the length of the sequence motif before and after the site of the PTM.

❖ **To show only PTMs**

1. In the Protein Identification Details view, click the **Coverage** tab if it is not already selected.
2. Select the **Show Only PTMs** check box.

This option displays only the PTMs on the Sequence and Modification List pages, so it helps you focus on biologically relevant modifications. It excludes all modifications classified as artifacts, chemical derivatives, isotopic labels, and synthetic peptide protection groups from the display.

❖ **To display filtered-out PSMs**

1. In the Protein Identification Details view, click the **Coverage** tab if it is not already selected.
2. Select the **Include PSMs That Are Filtered Out** check box.

This option filters out PSMs that have been found but filtered out by the display filters.

❖ **To control the display of the probability of a PTM occurring on a site**

1. In the Protein Identification Details view, click the **Coverage** tab if it is not already selected.

If the site probability scores are available (that is, if you included the *ptmRS* node in the processing workflow), the application color-codes the scores, and you can select the threshold.

2. From the Threshold list in the PTM Site Probabilities area, select the category of probability for a modification to occur on a sequence site.


For example, when you choose 25 from the Threshold list, the Protein Identification Details view displays modifications that have a low (25–100%) probability of occurrence on a site, as shown in [Figure 257](#).

Note The application does not display modifications with less than a 25% chance of occurring on a site.

9 Searching for Post-Translational Modifications


Viewing PTM Information on the Protein Identification Details View

Figure 257. Protein Identification Details view showing modifications with a low probability of occurring on a site



As another example, to display only those modifications with a very high probability (99–100%) of a modification occurring on a site, select 99 from the Threshold list, as shown in [Figure 258](#).

Figure 258. Protein Identification Details view showing modifications with a high probability of occurring on a site



❖ **To display UniProt annotations of the protein**

1. Set up and run a protein annotation workflow that includes the ProteinCenter Annotation node (see “[Creating a Protein Annotation Workflow](#)” on page 296).

This node retrieves UniProt database information, including PTMs, from ProteinCenter and saves it in Proteome Discoverer results files.

2. Open the resulting .pdResult file, select the protein of interest, and open the Protein Identification Details view.
3. Select the **Annotate PTMs Reported in Uniprot** check box.

You can find out if the protein has modification information by clicking the ProteinCard tab and then clicking on the Features tab. All lines containing the MOD_RES or CARBOHYD keyword should be displayed as PTMs in the sequence.

If annotation information is available, the modification appears as a single-letter abbreviation in black font in a small gray box beneath the amino acids of the protein in the sequence table. Information about the PTMs found in the UniProt database appear under the UniProt PTMs category in the area to the left of the sequence table. In addition, the likelihood of a PTM occurring on an amino acid appears under the PTM Site Probabilities heading in this same area. The Help illustrates this information.

❖ To copy the colored bar on the Sequence page

Right-click the bar and choose **Copy**.

❖ To save the colored bar on the Sequence page

1. Right-click the bar and choose **Save As**.
2. In the Save As dialog box, browse to the location where you want to save the file, and type the file name in the File Name box.
3. Click **Save**.

❖ To export the Modification List page to Excel

1. Right-click the Modification List table and choose **Export to Excel Workbook**.
2. In the Save Modification List as Microsoft Excel File dialog box, browse to the location where you want to save the Excel file, or specify the path and file name in the File Name box.
3. Click **Save**.

Filtering Phosphorylation Site Probabilities

You can filter phosphorylation site probability data by using display filters. For phosphorylation probabilities produced by using the *ptmRS* node in the processing workflow, you can use the Exceed operator in the Display Filters dialog box to display the following:

- Any peptides that contain a phosphorylation of an arbitrary amino acid
- All peptides that contain a phosphorylation of an arbitrary amino acid where each probability exceeds a given probability
- Any peptides that contain a phosphorylation of an amino acid that exceeds a given probability value

You can combine these conditions with other filter conditions by using AND and OR.

For general information about filtering data, see [“Filtering with Display Filters” on page 254](#).

For information about the *ptmRS* node, see [“Using the *ptmRS* Node” on page 347](#) and the Help.

❖ To filter phosphorylation site probabilities


1. In a .pdResult containing phosphorylation modification data (see [Figure 259](#)), choose **View > Display Filter**, or click the **Display Filter** icon,  .

Figure 259. .pdResult file containing phosphorylation modification data in the ptmRS Best Site Probabilities column

#	Checked	Confidence	Identifying Node	PSM Ambiguity	Annotated Sequence	Modifications	# Protein Groups	# Proteins	Master Protein Accessions	Protein Accessions
1	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	yPQKNEFLFTkTNIINsr	Y1(Phospho); K4(Dimethyl);	1	1	Q2NJ15	Q2NJ15
2	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	RKAsGPPVSELITK	S4(Phospho)	2	7	P16402; P02253; P15864	P16402; P10412;
3	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	RKAsGPPVSELITK	S4(Phospho)	2	7	P16402; P02253; P15864	P16402; P10412;
4	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	KAsGPPVSELITK	S3(Phospho)	2	8	P16402; P02253; P15864	P16402; P10412;
5	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	KAsGPPVSELITK	S3(Phospho)	2	8	P16402; P02253; P15864	P16402; P10412;
6	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	HLQLAVNDEELNKLGGVr	R7(Methyl)	2	9	Q8IU66; Q64522; Q7ZUY3	P69141; Q8IU66;
7	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	EIQTAVrLVLPGELAK	R7(Methyl)	0	25	P16867	O65819; P54345;
8	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	STISSrEIQTAVRLLPGELAK	R6(Methyl)	0	4	Q96A08	Q96A08; P16887;
9	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	STISSrEIQTAVRLLPGELAK	R6(Methyl)	0	7		Q27894; P02289;
10	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	STISSrEIQTAVRLLPGELAK	R6(Methyl)	0	4	Q96A08	Q96A08; P16887;
11	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	STISSrEIQTAVRLLPGELAK	R6(Methyl)	0	7		Q27894; P02289;
12	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	KVLRDNIQGITKPAIR	K1(Methyl)	1	20	P27996; P84043	Q41811; P08436;
13	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	KVLRDNIQGITKPAIR	R4(Methyl)	1	20	P27996; P84043	Q41811; P08436;
14	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	GIWLNAAASrR	R4(Dimethyl); R10(Methyl)	1	1	P70339	P70339
15	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	KVLRDNIQGITKPAIR	R4(Dimethyl)	1	20	P27996; P84043	Q41811; P08436;
16	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	KVLRDNIQGITKPAIR	R4(Dimethyl)	1	20	P27996; P84043	Q41811; P08436;
17	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	KVLRDNIQGITKPAIR	R4(Dimethyl)	1	20	P27996; P84043	Q41811; P08436;
18	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	KVLRDNIQGITKPAIR	K1(Dimethyl)	1	20	P27996; P84043	Q41811; P08436;
19	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	HArLNAAALAPak	R3(Dimethyl); T10(Phospho)	1	1	Q9NU22	Q9NU22
20	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	LArEIAQDFKTDLR	R3(Dimethyl)	1	1	Q5RCC9	Q5RCC9
21	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	TKTETITGFQVDAIPANGQT	R24(Methyl)	0	2	P11276; P04937	P11276; P04937
22	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	GLSEDITTEETLKESFDGVSr	R20(Methyl)	1	3	Q5RF26	Q4R4J7; P19338;
23	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	SrSSRAGLGFVGR	R2(Methyl)	4	28	Q96QV6; Q8IU66; Q64522; Q93077; P16104; P02262; Q6F113	P35062; Q64523;
24	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	SrSSRAGLGFVGR	R2(Methyl)	4	28	Q96QV6; Q8IU66; Q64522; Q93077; P16104; P02262; Q6F113	P35062; Q64523;
25	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	LrPTTIVALVIMAVFLVGVPr	R2(Dimethyl); T30(Phospho)	1	1	Q9TUE1	Q9TUE1
26	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VrPASSAASVYAGAGSGS	R2(Dimethyl)	1	2	P05784	P05784; Q5BJY9
27	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VTQDELKEVFEDAElr	R17(Dimethyl)	1	3	Q5RF26	Q4R4J7; P19338;
28	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VPLAGAAGPGIGrAAGrGI	R14(Dimethyl); R18(Dimethyl)	1	5	P14678	P27048; P17136;
29	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	STISSrEIQTAVRLLPGELAK	R6(Methyl)	0	7		Q27894; P02289;
30	<input type="checkbox"/>	●	Sequest HT (A2)	Rejected	STISSrEIQTAVRLLPGELAK	R6(Methyl)	0	4	Q96A08	Q96A08; P16887;
31	<input type="checkbox"/>	●	Sequest HT (A2)	Selected	RkASGPPVSELITK	K2(Methyl)	2	7	P16402; P02253; P15864	P16402; P10412;
32	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
33	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
34	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
35	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
36	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
37	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
38	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	VGAGAPVYmAAVLELYTAE	M9(Oxidation)	1	3	Q64522; Q6F113	Q64523; Q64522;
39	<input type="checkbox"/>	●	Sequest HT (A2)	Unambiguous	TEIDKPSQmQVTDVQDNSI	M9(Oxidation)	1	2	P11276; P07589	P11276; P07589



2. In the pane on the left, select **Peptide Groups**.
3. In the Add Property box, select **ptmRS PhosphoSite Probabilities**.
4. In the operator box, select **Exceed** if it is not already selected.
For Amino Acid appears in the values box.
5. Type the appropriate value in the values box to the right of the For Amino Acid box.
6. Click the **Apply** icon, .

Figure 260 shows an empty filter that you can further configure in the three different ways shown in “Example 1,” “Example 2,” and “Example 3.”

9 Searching for Post-Translational Modifications

Filtering Phosphorylation Site Probabilities

Figure 260. Empty Display Filters pane that addresses *ptmRS* probabilities



To reverse an addition or change immediately after you make it, click **Remove** next to it. For more information on using display filters, see “[Filtering with Display Filters](#)” on [page 254](#).

Example 1




To show only those peptides that contain a phosphorylation of an arbitrary amino acid with positive probability, type **0** in the values box to the right of the operator box and **Any** in the values box to the right of the For Amino Acid box, as shown in [Figure 261](#). Click the **Apply** icon,  **Apply**.

Figure 261. Peptides that contain a phosphorylation of an arbitrary amino acid with positive probability



Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info	Result Statistics					
#	Checked	Confidence	Sequence	Modifications	# Protein Groups	# Proteins	# PSMs	Master Protein Accessions	PhosphoRS: Best Site Probabilities		
1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LCDFGSASHVADNDITPYLVS	1xPhospho [Y18(100)]	1	3	10	IPI00013721.3	Y18(Phospho): 100		
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LRSYEDMIGEEVPSDQYYWAPLAQHER	1xPhospho [Y4(99.6)]	1	29	2	IPI00182540.7	Y4(Phospho): 99.56		
3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YDEEGGGEEDQDYDLSQLQPDVPEDAIKPV	1xPhospho [S16(99)]	1	2	2	IPI00290085.2	S16(Phospho): 99.03		
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	KISENSYLDLLEIGPQLSSSTFDSEK	1xPhospho [S3(98.8)]	1	7	4	IPI00883896.1	S3(Phospho): 98.8		
5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	AYSQEEITQGFEETGDTLYAPYSTHFQLNQPP	1xPhospho [S3(99.2)]	1	1	2	IPI00022624.1	S3(Phospho): 99.23		
6	<input type="checkbox"/>	<input checked="" type="checkbox"/>	RYSDLTTLVAFSSSVYPTK	1xPhospho [S3(99.7)]	1	2	4	IPI00013214.2	S3(Phospho): 99.47		
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	REPGYTPPGAGNQNPMMYPTGPK	1xPhospho [T6(99.6)]	1	3	6	IPI00023704.1	T6(Phospho): 99.51		
8	<input type="checkbox"/>	<input checked="" type="checkbox"/>	REPGGWGAGASAPVEDDSDAETYGEENDEQGG	1xPhospho [S/T/Y]	1	1	3	IPI00006904.1	S18(Phospho): 49.92; T22(Phospho): 49.92		
9	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SYDVPPPMPEDHPFYSNISK	1xPhospho [S/Y]	1	3	10	IPI00549725.6	S1(Phospho): 50; Y2(Phospho): 50		
10	<input type="checkbox"/>	<input checked="" type="checkbox"/>	HSSYPAGTEDEDEGMGEEPSFFR	1xPhospho [S19(100)]	1	2	3	IPI00024291.2	S2(Phospho): 98.72		
11	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SSSSGSDDYAYTQALLLHQR	1xPhospho [S/Y/T]	1	2	5	IPI00166840.4	S3(Phospho): 33.07; S6(Phospho): 33.07; Y9(Phospho): 33.07		
12	<input type="checkbox"/>	<input checked="" type="checkbox"/>	TRSYDNLTACDNTVPLASR	1xPhospho [T/S/Y]	1	4	3	IPI00185018.7	T1(Phospho): 33.33; S3(Phospho): 33.33; Y4(Phospho): 33.33		
13	<input type="checkbox"/>	<input checked="" type="checkbox"/>	AFSSRSYTSRPGSR	2xPhospho [S3(100); S6(97.4)]	1	8	1	IPI01021414.1	S3(Phospho): 99.99; S4(Phospho): 33.34; S6(Phospho): 33.33; Y7(Phospho): 33.33		
14	<input type="checkbox"/>	<input checked="" type="checkbox"/>	ADVQLFMDDDSYSHHSGLLEYADPEK	1xPhospho [S11(99.5)]	1	2	8	IPI00019870.1	S11(Phospho): 99.52		
15	<input type="checkbox"/>	<input checked="" type="checkbox"/>	ADREDESSPYAAMLAQDVQR	1xPhospho [S7(99.1)]	1	3	3	IPI00026271.5	S7(Phospho): 99.11		
16	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EGMNPYSYDEADSDQHDAYLER	2xPhospho [S13(100); S/Y]	1	1	2	IPI00005154.1	S6(Phospho): 50; Y7(Phospho): 50; S13(Phospho): 100		
17	<input type="checkbox"/>	<input checked="" type="checkbox"/>	RGLLYDSDEEERPAR	1xPhospho [S7(99.4)]	1	2	1	IPI00184330.5	S7(Phospho): 99.43		
18	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YTLENKEEGLSDTEADAVSGQLPDPTNPSAC	1xPhospho [Y/T/S]	1	1	4	IPI00170935.1	Y1(Phospho): 19.95; T2(Phospho): 19.95; S10(Phospho): 19.95; S12(Phospho): 19.95		
19	<input type="checkbox"/>	<input checked="" type="checkbox"/>	GGPGSAVSPYPTFNFPSSDVAALHK	1xPhospho [S/Y/T]	1	3	3	IPI00218918.5	S8(Phospho): 50; Y10(Phospho): 50		
20	<input type="checkbox"/>	<input checked="" type="checkbox"/>	TIDLGAAAHYTGDKASPDQNASHTHTPQSSVK	1xPhospho [S16(100)]	1	3	3	IPI00975549.1	T11(Phospho): 49.69; S16(Phospho): 49.69		
21	<input type="checkbox"/>	<input checked="" type="checkbox"/>	AYSFCGTVEYMAPEVWNR	1xPhospho [S3(99.7)]	2	11	9	IPI00477982.1; IPI00020898.1	S3(Phospho): 99.67		
22	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YCTETSGVHGDSYPYSGTMDTHSLESK	1xPhospho [Y14(99.2)]	1	2	3	IPI00018370.4	Y14(Phospho): 99.19		
23	<input type="checkbox"/>	<input checked="" type="checkbox"/>	QQAAYYAQTSPQGMQHPAPQGG	1xPhospho [Y6(97.2)]	2	3	14	IPI00641948.4; IPI00983652.2	Y6(Phospho): 97.16		

Example 2

To show only those peptides that contain a phosphorylation of an arbitrary amino acid where each probability exceeds a given probability (90% in this example), type **90** in the box to the right of the operator box and **Any** in the values box to the right of the For Amino Acid box, as shown in [Figure 262](#). Click the **Apply** icon,  **Apply**.

9 Searching for Post-Translational Modifications


Filtering Phosphorylation Site Probabilities

Figure 262. Peptides that contain a phosphorylation of an arbitrary amino acid where each probability exceeds a given probability

Peptide Groups

AND

PhosphoRS: Phospho Site Probabilities



Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info	Result Statistics				
#	Checked	Confidence	Sequence	Modifications	# Protein Groups	# Proteins	# PSMs	Master Protein Accessions	PhosphoRS Best Site Probabilities	
1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LCDFGSASHVADNDITPYLVS	1xPhospho [Y18(100)]	1	3	10	IPI00013721.3	Y18(Phospho): 100	
2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LRSYEDMIGEEVPSDQYVWAPLAQHER	1xPhospho [Y4(99.6)]	1	29	2	IPI000182540.7	Y4(Phospho): 99.56	
3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SHYADVDPENGFLESNLGK	1xPhospho [S17(100)]	1	2	2	IPI000410034.2	Y3(Phospho): 99.41	
4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	IADPEHDHTGFLTEYVATR	1xPhospho [Y15(99.2)]	1	8	1	IPI00018195.3	Y15(Phospho): 99.17	
5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	QISYNSYDLSDSNVTEETPEGEEHHPVADTENK	1xPhospho [Y6(98.2)]	1	2	11	IPI000410693.4	Y6(Phospho): 98.19	
6	<input type="checkbox"/>	<input checked="" type="checkbox"/>	IGEGTYGVVYK	1xPhospho [Y10(100)]	2	11	3	IPI01025117.1; IPI00026689.4	Y6(Phospho): 99.64	
7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YCTETSGVHGDSYSGSGTMDTHSLESK	1xPhospho [Y14(99.2)]	1	2	3	IPI00018370.4	Y14(Phospho): 99.19	
8	<input type="checkbox"/>	<input checked="" type="checkbox"/>	AELSYRGPVSGTEPEPVYSMEAADYR	1xPhospho [Y5(99.4)]	1	6	2	IPI01013939.1	Y5(Phospho): 99.39	
9	<input type="checkbox"/>	<input checked="" type="checkbox"/>	QQAAYYQTSPOGMQHPAPQGG	1xPhospho [Y6(97.2)]	2	3	14	IPI00641948.4; IPI00983652.2	Y6(Phospho): 97.16	
10	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YMIGVITYGGDDIPLSPYR	1xPhospho [Y17(99.6)]	1	9	1	IPI00900293.1	Y17(Phospho): 99.59	
11	<input type="checkbox"/>	<input checked="" type="checkbox"/>	EELTPGAPQHQLPVPVGSPEPYFGQQAPGPEF	1xPhospho [T4(99.9)]	1	1	4	IPI00012199.1	Y22(Phospho): 98.36	
12	<input type="checkbox"/>	<input checked="" type="checkbox"/>	LIEDNEYTAR	1xPhospho [Y7(99.7)]	1	14	1	IPI00013981.4	Y7(Phospho): 99.66	
13	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YTPTQQGNMQVLVITYGGDPIPK	1xPhospho [Y1(99.2)]	1	9	2	IPI00900293.1	Y1(Phospho): 99.2	
14	<input type="checkbox"/>	<input checked="" type="checkbox"/>	YSPSQNSPIHHSRR	2xPhospho [Y1(98.9); S7(100)]	1	10	3	IPI00006079.1	Y1(Phospho): 98.9; S7(Phospho): 100	
15	<input type="checkbox"/>	<input checked="" type="checkbox"/>	SESVYADIR	1xPhospho [Y6(100)]	1	4	1	IPI00022558.2	Y6(Phospho): 100	
16	<input type="checkbox"/>	<input checked="" type="checkbox"/>	VADPDHDHTGFLTEYVATR	1xPhospho [Y15(100)]	1	3	2	IPI00003479.3	Y15(Phospho): 99.98	
17	<input type="checkbox"/>	<input checked="" type="checkbox"/>	VADPDHDHTGFLTEYVATR	2xPhospho [T13(100); Y15(100)]	1	3	3	IPI00003479.3	T13(Phospho): 100; Y15(Phospho): 99.99	
18	<input type="checkbox"/>	<input checked="" type="checkbox"/>	AAYGDLSSSEEEENEPESLGVVYK	2xPhospho [Y3(99.3); S/Y]	1	2	1	IPI00007343.1	Y3(Phospho): 99.01; S7(Phospho): 50.48; S8(Phospho): 50.48	
19	<input type="checkbox"/>	<input checked="" type="checkbox"/>	VIEDNEYTAR	1xPhospho [Y7(99.7)]	1	8	1	IPI00940936.1	Y7(Phospho): 99.68	
20	<input type="checkbox"/>	<input checked="" type="checkbox"/>	IGEGTYGVVYK	2xPhospho [T5(100); Y6(100)]	2	11	2	IPI01025117.1; IPI00026689.4	T5(Phospho): 100; Y6(Phospho): 100	
21	<input type="checkbox"/>	<input checked="" type="checkbox"/>	KPYRIESDEEEDFENVGK	1xPhospho [Y3(99.3)]	1	3	2	IPI00290652.4	Y3(Phospho): 99.33	
22	<input type="checkbox"/>	<input checked="" type="checkbox"/>	IFQGSFGGFTLYENPHYQSPNMHR	1xPhospho [Y12(97.3)]	1	1	1	IPI00181728.1	Y12(Phospho): 97.29	
23	<input type="checkbox"/>	<input checked="" type="checkbox"/>	GFYKGTGPALMAYVAENSVLFMCYGCQQVF	1xPhospho [Y13(96.3)]	1	1	1	IPI00011190.3	Y13(Phospho): 96.32	

Example 3


To show only those peptides that contain a phosphorylation of a specific amino acid (serine) that exceeds a given probability (90% in this example), type **90** in the box to the right of the operator box and **S** in the values box to the right of the For Amino Acid box, as shown in [Figure 263](#). Click the **Apply** icon,  **Apply**.

Figure 263. Peptides that contain a phosphorylation of an amino acid that exceeds a probability of 90%

Peptide Groups

AND

ptmRS [5]: Phospho Site Probabilities for amino acid

Protein Groups		Proteins	Peptide Groups	PSMs	MS/MS Spectrum Info	Result Statistics					
#	Checked	Confidence	Sequence	Modifications	# Protein Groups	# Proteins	# PSMs	Master Protein Accessions	PhosphoRS: Best Site Probabilities		
1	<input type="checkbox"/>	<input checked="" type="radio"/>	SSGNSSSSGSGSGSTSAGSSSPGAR	1xPhospho [S21(100)]	1	7	7	IP100294834.6	S21(Phospho): 100		
2	<input type="checkbox"/>	<input checked="" type="radio"/>	EAYSGCGPVDECCPPPPSPVHK	1xPhospho [S19(99.5)]	1	2	8	IP100910135.1	S19(Phospho): 99.28		
3	<input type="checkbox"/>	<input checked="" type="radio"/>	EGITGPPADSSKPIGPDDAIDALSSDFTCGSPTA	1xPhospho [S31(99.2)]	1	12	3	IP100761160.1	S31(Phospho): 99.21		
4	<input type="checkbox"/>	<input checked="" type="radio"/>	RNGAAGPHSPDPLLEQAFGDLTLPVVPK	1xPhospho [S9(100)]	1	2	2	IP100015894.1	S9(Phospho): 100		
5	<input type="checkbox"/>	<input checked="" type="radio"/>	IDASKNEEDEGHSSNSPR	1xPhospho [S16(100)]	1	4	7	IP100028888.1	S16(Phospho): 99.52		
6	<input type="checkbox"/>	<input checked="" type="radio"/>	HGAPSPSHPIAPQAAAAAALR	1xPhospho [S5(100)]	2	7	1	IP100183294.4; IP100900325.1	S5(Phospho): 100		
7	<input type="checkbox"/>	<input checked="" type="radio"/>	SPSAGILRTPIQVESSPQGLPAGGLEGLK	2xPhospho [S3(98.1); S16(99.1)]	1	4	4	IP100010428.4	S3(Phospho): 98.14; S16(Phospho): 99.06		
8	<input type="checkbox"/>	<input checked="" type="radio"/>	SPTMEQAVQTASAHLPAVAVGR	1xPhospho [S1(100)]	1	3	8	IP100410590.5	S1(Phospho): 100		
9	<input type="checkbox"/>	<input checked="" type="radio"/>	LIRPEYAWIVQVPVSGAVYDRPGASPK	1xPhospho [S24(100)]	1	2	5	IP100748807.2	S24(Phospho): 100		
10	<input type="checkbox"/>	<input checked="" type="radio"/>	SGGSGHVAEAPASPEQLDQNK	1xPhospho [S13(100)]	1	2	3	IP100004233.2	S13(Phospho): 100		
11	<input type="checkbox"/>	<input checked="" type="radio"/>	SVNVSILTQEELDGLDELSVR	1xPhospho [S5(100)]	1	6	2	IP100032064.1	S5(Phospho): 100		
12	<input type="checkbox"/>	<input checked="" type="radio"/>	EDNSEGEILEEVGGLEEEEDHHMEFCR	1xPhospho [S4(100)]	1	3	2	IP100455210.2	S4(Phospho): 100		
13	<input type="checkbox"/>	<input checked="" type="radio"/>	GFSFVATGLMEDDGKPR	1xPhospho [S3(100)]	1	6	2	IP100477982.1	S3(Phospho): 100		
14	<input type="checkbox"/>	<input checked="" type="radio"/>	TDGFAEAIHSPQVAVGR	1xPhospho [S10(100)]	1	1	5	IP100742682.2	S10(Phospho): 100		
15	<input type="checkbox"/>	<input checked="" type="radio"/>	LAFFSPFGHLSEGFSSHR	1xPhospho [S5(100)]	1	5	6	IP100179473.9	S5(Phospho): 100		
16	<input type="checkbox"/>	<input checked="" type="radio"/>	AALGLQDSDEDAAVDIEQIESMFNSK	1xPhospho [S8(100)]	1	2	1	IP100939523.1	S8(Phospho): 100		
17	<input type="checkbox"/>	<input checked="" type="radio"/>	NGAAGPHSPDPLLEQAFGDLTLPVVPK	1xPhospho [S8(100)]	1	2	1	IP100015894.1	S8(Phospho): 100		
18	<input type="checkbox"/>	<input checked="" type="radio"/>	GAGDGSDEEVDGKADGAEAKPAE	1xPhospho [S6(100)]	1	3	10	IP100019502.3	S6(Phospho): 100		
19	<input type="checkbox"/>	<input checked="" type="radio"/>	VPSPPDEHQEAENAVSSGNR	1xPhospho [S3(100)]	1	7	3	IP100000398.2	S3(Phospho): 100		
20	<input type="checkbox"/>	<input checked="" type="radio"/>	GDLSVDEEEEEEMDVDEATGAVK	1xPhospho [S4(100)]	1	1	1	IP100170596.1	S4(Phospho): 100		
21	<input type="checkbox"/>	<input checked="" type="radio"/>	RRSEVESTTESQDK	1xPhospho [S3(100)]	1	2	1	IP100293845.5	S3(Phospho): 100		
22	<input type="checkbox"/>	<input checked="" type="radio"/>	EGEELGEGERPEEDAALSSDEAVEVEEVIE	2xPhospho [S21(100); S22(100)]	1	3	6	IP100176903.2	S21(Phospho): 100; S22(Phospho): 100		
23	<input type="checkbox"/>	<input checked="" type="radio"/>	SFAHLLPSPSPSEGSYVQGHSGQLGGHYADS	2xPhospho [S8(100); S12(100)]	1	1	3	IP100217661.4	S8(Phospho): 100; S12(Phospho): 100		

9 Searching for Post-Translational Modifications

Filtering Phosphorylation Site Probabilities

Performing Quantification

This chapter describes how to perform precursor ion-, reporter ion-, and peak area-based quantification in the Proteome Discoverer application.

Contents

- [Precursor Ion Quantification](#)
- [Reporter Ion Quantification](#)
- [Precursor Ion Area Detection](#)
- [Performing Quantification](#)
- [Using a Quantification Method](#)
- [Searching for Quantification Modifications with Mascot](#)
- [Generating Quantification Ratios](#)
- [Calculating Quantification Results](#)
- [Filtering Quantification Data with Signal-to-Noise Values](#)
- [Identifying Isotope Patterns in Precursor Ion Quantification](#)
- [Viewing Quantification Results](#)
- [Troubleshooting Quantification](#)

Precursor Ion Quantification

In precursor ion quantification, also called isotopically labeled quantification, protein abundance is determined from the relative MS signal intensities of an isotopically labeled sample and an unlabeled control sample. (However, this type of quantification does not necessarily have to use labeled and unlabeled samples. Some laboratories use heavy and medium labels to eliminate the bias of unlabeled peptides.)

Stable-isotope labeling by amino acid in cell culture (SILAC™) is a proteomics identification and quantification technique that uses in-vivo metabolic labeling to detect differences in the abundance of proteins in multiple samples. It is a type of isotopically labeled quantification, which uses stable (non-radioactive) heavy isotopes as labels. You can also introduce the stable isotopes by chemical labeling at the protein or peptide level with the isotopomeric tags, for example, dimethyl labeling.

The following default quantification methods are available for precursor ion (isotopically labeled) quantification:

- SILAC 2plex (Arg10, Lys6): Uses label 13C(6)15N(4) for arginine and label 13C(6) for lysine.
- SILAC 2plex (Arg10, Lys8): Uses label 13C(6)15N(4) for arginine and 13C(6)15N(2) for lysine.
- SILAC 2plex (Ile6): Uses label 13C(6) on isoleucine.
- SILAC 3plex (Arg6, Lys4|Arg10, Lys8): Uses heavy label 13C(6)15N(4) for arginine and 13C(6)15N(2) for lysine. Uses medium label 13C(6) for arginine and 2H(4) for lysine.
- SILAC 3plex (Arg6, Lys6|Arg10, Lys8): Uses heavy label 13C(6)15N(4) for arginine and 13C(6)15N(2) for lysine. Uses medium label 13C(6) for arginine and 13C(6) for lysine.
- Dimethylation 3plex: Chemically adds isotopically labeled dimethyl groups to the N-terminus and to the ϵ -amino group of lysine.
- Full ^{18}O labeling: Introduces 4 Da mass tags through the enzyme-catalyzed exchange reaction of C-terminal oxygen atoms with ^{18}O . This method is for fully labeled samples.
- Incomplete ^{18}O labeling: Introduces 2 or 4 Da mass tags through the enzyme-catalyzed exchange reaction of C-terminal oxygen atoms with ^{18}O . This method is for incompletely labeled samples.

The discussion of precursor ion quantification includes the following topics:

- [Chemical Elements Supported in Precursor Ion Quantification](#)
- [SILAC 2plex Methods](#)
- [SILAC 3plex Methods](#)
- [Dimethylation 3plex Method](#)
- [18O Labeling Method](#)

Chemical Elements Supported in Precursor Ion Quantification


You can use any chemical element and each of its available isotopes for precursor ion quantification. However, you must adhere to certain syntactical restrictions when you enter chemical compositions. The notation from unimod.org, which sets the nominal mass prefix in parentheses—for example, (13)C4 (15)N2—only applies to the elemental isotopes ^2H , ^{13}C , ^{15}N , and ^{18}O .

SILAC 2plex Methods

In a typical SILAC quantification experiment, two cell populations grow in media that are deficient in lysine and arginine. One population grows in a medium containing normal (“light”) amino acids, such as lysine ($^{12}\text{C}_6^{14}\text{N}_2$). The other population grows in a medium containing amino acids where stable heavy isotopes, such as lysine 6 ($^{13}\text{C}_6^{14}\text{N}_2$) or lysine 8 ($^{13}\text{C}_6^{15}\text{N}_2$), have been substituted for normal atoms. SILAC quantification usually uses “heavy” arginine and lysine, because these are the cleavage sites for the generally used trypsin protease. Both populations incorporate these amino acids into proteins through natural cellular protein synthesis. The cells growing in the medium with the heavy isotopes incorporate these isotopes into all of their proteins.

After altering the proteome in one sample through chemical treatment or genetic manipulation, you combine equal amounts of protein from both cell populations and then digest with trypsin before MS analysis. Because peptides labeled with “heavy” and “light” amino acids are chemically identical, they co-elute during reverse-phase chromatographic separation. This means they are detected simultaneously during MS analysis. To determine the average change in protein abundance in the treated sample, you use the relative peak intensities of multiple isotopically distinct peptides from each protein, as shown in [Figure 264](#).

Figure 264. Schematic workflow for SILAC-based peptide and protein quantification



SILAC can differentiate peptides in single MS mode without requiring you to perform tandem mass spectrometry. However, SILAC cannot identify peptides, so you must use tandem mass spectrometry for that purpose.

You can choose from several SILAC 2plex methods, for example, (Arg10, Lys6) and (Arg10, Lys8), to compare two samples.

SILAC 3plex Methods

SILAC 3plex methods are similar to SILAC 2plex methods except, in addition to a “heavy” sample (containing, for example, Arg10 and Lys8), they also use a “medium” sample (containing, for example, Arg6 and Lys4). Protein abundance is determined from the relative MS signal intensities of the heavy sample, medium sample, and a control sample containing “light” (^{12}C and ^{14}N) arginine and lysine.

Dimethylation 3plex Method

The application also includes the dimethylation 3plex method, which is not metabolomic labeling in cell culture but a form of peptide chemical labeling. This method uses formaldehyde and sodium cyanoborohydride to add dimethyl groups $(\text{CH}_3)_2$ to the N-terminus and to the ϵ -amino group of lysine. By choosing the isotopomers of formaldehyde and sodium cyanoborohydride, you can create light, medium, and heavy labels:

- For the light label, the (natural-isotope) dimethyl group is $^{12}\text{C}_2^1\text{H}_6$.
- For the medium label, the dimethyl group is $^{12}\text{C}_2^2\text{H}_4^1\text{H}_2$, which is 4 Da more massive.
- For the heavy label, the dimethyl group is $^{13}\text{C}_2^2\text{H}_6$, which is an additional 4 Da more massive.

You can use the dimethylation 3plex method to compare up to three samples.

^{18}O Labeling Method

The ^{18}O labeling method introduces 2 or 4 Da mass labels through the enzyme-catalyzed exchange reaction of C-terminal oxygen atoms with ^{18}O .

- Full ^{18}O labeling: Introduces 4 Da mass tags through the enzyme-catalyzed exchange reaction of C-terminal oxygen atoms with ^{18}O . This method is for fully labeled samples.
- Incomplete ^{18}O labeling: Introduces 2 or 4 Da mass tags through the enzyme-catalyzed exchange reaction of C-terminal oxygen atoms with ^{18}O . This method is for incompletely labeled samples.

Reporter Ion Quantification

In contrast to the metabolic labeling used by isotopically labeled precursor ion quantification methods such as SILAC, isobarically labeled reporter ion quantification methods use reagents, or tags, to enzymatically or chemically label proteins and peptides. Reporter ion quantification uses tags that have the same mass. (A reporter ion is a fragment ion produced from the MS/MS of a peptide with an attached tag.)

The application supports reporter ion quantification for Tandem Mass Tag (TMT), Isobaric Tag for Relative and Absolute Quantification (iTRAQ), and any user-defined tags. The MS/MS scan performs identification and quantification with both TMT and iTRAQ.

Note User-defined tags require an additional license key provided by Thermo Fisher Scientific.

For iTRAQ, 4plex and 8plex default methods are available. For TMT, 2plex, 6plex, and 10plex default methods are available. You can also add new methods by either deriving them from an existing method or creating new ones from the beginning.

The discussion of reporter ion quantification includes the following topics:

- [TMT Quantification](#)
- [iTRAQ Quantification](#)

TMT Quantification

TMT quantification is a reproducible, highly accurate quantification method that provides both comparative and absolute MS/MS-based quantification of proteins and peptides in biological samples. TMT tagging produces data to calculate the relative abundances of proteins. You can evaluate differential protein expression in one to six samples in a single experiment.

Each sample is labeled with chemically identical tags before mixing the samples, and a single MS run generates a single peak for each peptide, irrespective of which tag it has been given. Between the normalizer and reporter is a cleavable linker, which breaks during MS/MS. The mass reporter ion is split off and measured by the mass spectrometer.

Only MS/MS fragmentation can differentiate the tagged peptides. The reporter ion, measured by the mass spectrometer, generates a different low-mass peak for each sample. As a result, the peak height/peak integral for each reporter denotes the relative amount of protein originating from each of the labeled samples.

With the quantification functions in the application, you can set filters to see only unique peptides so that every protein associated with the same peptide is not counted, producing a best-results list of peptides. Filtering the number of proteins can give you a more robust final analysis of your experimental set.

Quantification with TMT tags is no different from quantification with iTRAQ (described in [“iTRAQ Quantification”](#) on [page 388](#)), except that it uses the following default methods by Proteome Sciences PLC:

- TMT 2plex
- iodo TMT 6plex
- Low-resolution iodo TMTe 6plex

- TMTe 6plex, which is the standard method when you use a high-enough resolution. In this method, different channel impurities do not affect each other.
- Low-resolution TMTe 6plex, which is the method that you must use when you acquire MS/MS data in an ion trap or in an instrument with low resolution. In this method, different channel impurities can affect each other (mass differences between ¹⁵N and ¹³C are not resolved).
- TMT 10plex

Note The TMT 6plex quantification method, created for an earlier version of 6-plex labels, is no longer listed. The TMTe 6plex methods replace it. If you require an analysis of the earlier 6-plex method, you must derive a new method with the appropriate reporter ions from the 10plex method, deactivating all the “_N” reporter ion tags in the method.

You can use these default methods to create your own quantification templates. For information on adding quantification methods, see [“Setting Up the Quantification Method”](#) on page 408.

[Table 26](#) lists the masses of the reporter ions of the tags available in the different TMT kits. The masses for the original TMT reagents, which are no longer available, are included for reference.

Table 26. Monoisotopic masses of the reporter ions after CID or HCD fragmentation of the tags in the different TMT kits


TMT 2plex		TMT 6plex (Original)		TMTe 6plex (Current)		TMT 10plex		iodo TMT 6plex	
Tag	Mass	Tag	Mass	Tag	Mass	Tag	Mass	Tag	Mass
126	126.12773	126	126.12773	126	126.12773	126	126.127726	126	126.12773
127	127.13108	127	127.13108	127	127.124760	127_N	127.124761	127	127.124760
						127_C	127.131081		
		128	128.13444	128	128.13444	128_N	128.128116	128	128.13444
						128_C	128.134436		
		129	129.13779	129	129.131468	129_N	129.131471	129	129.1311468
						129_C	129.13779		
		130	130.14115	130	130.14115	130_N	130.134825	130	130.14115
						130_C	130.141145		
		131	131.13818	131	131.13818	131	131.13818	131	131.13818

The iodo TMT 6plex includes cysteine reactive TMT reagents.

The TMT 10plex leverages the high resolution of recent mass spectrometers to routinely differentiate the ¹³C isotopes from the ¹⁵N isotopes^{1, 2}. For the 127, 128, 129, and 130 tags, the TMT 10plex contains two reagents: the ¹³C reagent and the ¹⁵N reagent. For the monoisotopic masses of the different reporter ions after CID or HCD fragmentation, see [Table 26](#).

Figure 265 shows the position of the ^{13}C and ^{15}N atoms in the different reagents. In this illustration, the stars indicate the positions of the ^{13}C and the ^{15}N substitutions, the red lines indicate the position of the ETD fragmentation sites, and the blue lines indicate the position of the CID fragmentation sites.

Figure 265. Structures of the TMT reagents contained in the TMT 10plex quantification method



Recent research concludes that, when you avoid applying any correction for isotopic impurities, the quantification results improve for the TMTe 6plex, TMT 10plex, and iodo TMT 6plex kits, so the default methods for these kits turn off the purity correction.

iTRAQ Quantification

iTRAQ is a protein quantification technique that uses isobaric amine-specific, stable isotope reagents to label all peptides in up to eight different samples simultaneously. The labeled peptides from each sample are combined and analyzed with liquid chromatography tandem mass spectrometry (LC/MS/MS). The same peptide from each sample appears as a single peak in the MS spectrum. In single MS mode, the differentially labeled versions of a peptide are

¹ McAlister, G. C.; Huttlin, E. L.; Haas, W.; Ting, L.; Jedrychowski, M. P.; Rogers, J. C.; Kuhn, K.; Pike, I.; Grothe, R. A.; Blethrow, J. D.; Blethrow, G. S. P. Increasing the Multiplexing Capacity of TMTs Using Reporter Ion Isotopologues with Isobaric Masses. *Anal. Chem.* **2012**, *84*, 7469–7478.

² Werner, T.; Becher, I.; Sweetman, G.; Doce, C.; Savitski, M. M.; Savitski, B. M. High-Resolution Enabled TMT 8-plexing. *Anal. Chem.* **2012**, *84*, 7188–7194.

indistinguishable. In tandem MS mode, which isolates and fragments peptides, each tag generates a unique reporter ion. Protein quantification compares the peak intensity of the reporter ions in the MS/MS spectra to assess the relative abundance of the peptides and, therefore, the proteins that they are derived from.

iTRAQ includes two default methods available from AB Sciex that you can use to label all peptides:

- iTRAQ 4plex
- iTRAQ 8plex

The application includes default quantification methods for processing data from iTRAQ 4plex- and iTRAQ 8plex-labeled samples. You can use these methods to create your own workflow templates. For information on adding quantification methods, see [“Changing a Quantification Method”](#) on page 424.

iTRAQ quantification works exactly the same as TMT quantification, except that TMT quantification offers 2plex, 6plex, and 10plex quantification methods, and iTRAQ offers 4plex and 8plex quantification methods.

Precursor Ion Area Detection

If you want to determine the area for any peptides, you can use precursor ion area detection. You might want to use this quantification method to obtain an idea of the relative quantities of all peptides in a sample for datasets where isotopically labeled methods such as SILAC or TMT were not employed. This type of quantification not only produces areas for individual peptides in a dataset but also produces the average of the peptide peak areas identified by the Top *N* Peaks Filter node for proteins in the final .pdResult report. This enables protein quantification across samples. The default quantification methods calculate the protein peak area on the basis of the top three unique peptides for the given protein.

Performing Quantification


Use the instructions in the following topics to perform quantification:

- [Performing Precursor Ion Quantification](#)
- [Performing Reporter Ion Quantification](#)
- [Performing Precursor Ion Area Detection](#)

Performing Precursor Ion Quantification

This topic describes the general procedure to use to perform precursor ion quantification and then describes the processing and consensus workflows in detail.

❖ **To perform precursor ion quantification**

1. Create or open a study and an analysis:
 - To create a study, see “[Creating a Study](#)” on [page 38](#).
 - To open an existing study, see “[Opening an Existing Study](#)” on [page 40](#).
 - To create an analysis, see “[Creating an Analysis](#)” on [page 71](#).
 - To open an existing analysis, see “[Opening an Existing Analysis](#)” on [page 72](#).
2. Select the quantification method. For information, see “[Setting Up the Quantification Method](#)” on [page 408](#). If the appropriate quantification method is not available, create one by following the instructions in “” on [page 419](#).
3. Create a precursor ion quantification workflow. For information, see “[Creating a Processing Workflow for Precursor Ion Quantification](#),” and “[Creating a Consensus Workflow for Precursor Ion Quantification](#)” on [page 394](#).
4. Specify the quantification channels. For information, see “[Specifying the Quantification Channels](#)” on [page 410](#).
5. Save the analysis. See “[Saving an Analysis](#)” on [page 80](#).
6. Save the study. See “[Saving a Study](#)” on [page 62](#).
7. Click the **Run** icon,  **Run**, in the upper right corner of the Analysis window.

The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

The application creates one or more MSF and .pdResult files, depending on whether the data files to be analyzed are batch-processed. It might reuse MSF files to produce a new .pdResult file if a different report for the given data is required.

The .pdResult file has the same name as the input file but with a .pdResult extension. The Analysis Results page displays the name of the .pdResult file.

Creating a Processing Workflow for Precursor Ion Quantification

To create a precursor ion quantification method, you must set up a processing workflow that includes the Precursor Ions Quantification node.

Default processing workflows are available for several instruments. You can use these default workflows as is or modify them to suit your needs. To access them, see “[Using Common Workflow Templates](#)” on [page 119](#). When you use a default workflow, you do not need to perform [step 1](#) through [step 7](#) in the following procedure because these steps are already included in the workflow. However, you must check the node parameters in the workflow to make sure they are appropriate for your data and add modifications to the search for the quantification to be successful.

❖ **To create a processing workflow for precursor ion quantification**

Note This procedure uses a SILAC 2plex example.

1. In the Workflow Editor, follow the general instructions for creating a processing workflow. See “Creating a Processing Workflow” on page 104.

The processing workflow must include the following nodes as a minimum:

- Spectrum Files node
 - Event Detector node
 - Spectrum Selector node
 - A search engine node (for example, Mascot or Sequest HT)
 - Precursor Ions Quantifier node
 - Fixed Value PSM Validator node, Target Decoy PSM Validator node, or Percolator node
2. In the Processing Workflow window, drag the **Spectrum Files** node to the workspace.
 3. Drag the **Spectrum Selector** node and the **Event Detector** node to the workspace.
 4. Connect the Spectrum Files node to the Event Detector node and to the Spectrum Selector node.
 5. Drag the **Precursor Ions Quantifier** node to the Workflow Tree pane and attach it directly to the Event Detector node.


The Precursor Ions Quantifier node performs quantification for isotopically labeled amino acids.

Note You cannot use the Precursor Ions Quantifier node and the Precursor Ions Area Detector node in the same workflow. You cannot use the Reporter Ions Quantifier node in a workflow that includes the Precursor Ions Quantifier node.

6. Drag the appropriate search engine node—for example, **Sequest HT**—to the Workflow Tree pane and attach it to the Spectrum Selector node.
7. Drag the **Fixed Value PSM Validator**, **Target Decoy PSM Validator**, or the **Percolator** node to the Workflow Tree pane and attach it to the search engine node.

Figure 266 illustrates the basic processing workflow for precursor ion quantification.

Figure 266. Basic processing workflow for precursor ion quantification



8. (Optional) Add any other appropriate nodes and connect them to the workflow.
9. In the Parameters pane of the Workflow Editor, click **Show Advanced Parameters**.
10. Click the **Event Detector** node and set the parameters for it in the Parameters pane as follows:
 - a. In the Mass Precision box, specify the expected standard deviation of the mass precision.

To create extracted ion chromatograms, use three times the standard deviation. The minimum value is 0.5 ppm. The maximum value is 4 ppm. The default is 2 ppm.
 - b. In the S/N Threshold box, specify a threshold signal-to-noise value that determines whether the application removes peaks from the spectrum.

The application removes peaks with a signal-to-noise value below this threshold.
The minimum value is 0.0, and there is no maximum value. The default is 1.
11. Click the **Spectrum Selector** node, and set the parameters for it in the Parameters pane:
 - a. Change the setting in the Max. Precursor Mass box to an appropriate setting. For example, for SILAC 2plex (Arg10, Lys6) quantification, set this option to **6500**.
 - b. Change the setting in the S/N Threshold box to an appropriate setting. For example, for SILAC 2plex (Arg10, Lys6) quantification, use the default of **1.5**.

For optional parameters that you can set for the Spectrum Selector node, refer to the Help.

12. Click the search engine node—for example, **Sequest HT**—and set its parameters in the Parameters pane:
 - a. In the Protein Database box, select the **FASTA** database.
 - b. In the Dynamic Modifications area, select the dynamic modifications.

Note If you are using the Mascot search engine node, you can use the **Quan Modifications** parameter rather than the Dynamic Modifications parameters to specify the modifications to search for. For instructions on using this method, see “[Checking the Quantification Method](#)” on [page 424](#).

For example, for SILAC 2plex (Arg10, Lys6) quantification, you might select the following two dynamic modifications:

- **13C(6)/ +6.020 Da (K)**
- **13C(6)/15N(4)/+10.008 Da (R)**

If you do not find these labels, you can enable them by following the instructions in “[Defining Chemical Modifications](#)” on [page 227](#).

- c. In the Static Modifications area, select the static modifications. For example, for SILAC 2plex (Arg10, Lys6) quantification, select **Carbamidomethyl/ +57.021 Da (C)** in the Static Modification box.
 - d. Set any other parameters that you prefer. For information on the search engine parameters available, refer to the Help.
13. Click the **Precursor Ions Quantifier** node and set the parameters for it in the Parameters pane.

For information on the parameters that you can set for the Precursor Ions Quantifier node, refer to the Help.


14. Click the **Fixed Value PSM Validator** node, the **Percolator** node, or the **Target Decoy PSM Validator** node, and set its Maximum Delta Cn parameter.

For information on this parameter, refer to the Help.

15. Set the parameters for all other nodes in the Parameters pane.

For information about all the parameters that you can set for each node, refer to the Help.

16. (Optional) Save the workflow as follows:

- a. In the Name box above the Workflow Tree pane, type a name for the processing workflow.
- b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the processing workflow.
- c. In the Workflow Editor, click the **Save** icon,  **Save** .
- d. In the Save Workflow dialog box, do the following:

- i. Browse to the file to save the template in, or type a file name in the File Name box.
- ii. In the Save As Type box, select **Processing Workflow File (*.pdProcessingWF)**.
- iii. Click **Save**.

The application saves the workflow with a .pdProcessingWF file name extension.

Creating a Consensus Workflow for Precursor Ion Quantification

To use a precursor ion quantification method, you must set up a consensus workflow that includes the Peptide and Protein Quantifier node.

Default consensus workflows are available. You can use these default workflows as is or modify them to suit your needs. To access them, see [“Using Common Workflow Templates” on page 119](#). When you use a default workflow, you do not need to perform [step 1](#) through [step 8](#) in the following procedure because these steps are already included in the workflow. However, you must check the node parameters in the workflow to make sure they are appropriate for your data and add modifications to the search for the quantification to be successful.

❖ To create a consensus workflow for precursor ion quantification

1. Follow the general instructions for creating a consensus workflow with the Workflow Editor. See [“Creating a Consensus Workflow” on page 112](#).


The consensus workflow must include the following nodes as a minimum:

- MSF Files node
 - PSM Grouper node
 - Peptide Validator node
 - Peptide And Protein Filter node
 - Protein Scorer node
 - Protein Grouping node
 - Peptide and Protein Quantifier node
2. In the Consensus Workflow window, drag the **MSF Files** node to the workspace.
 3. Drag the **PSM Grouper** node to the workspace, and connect it to the MSF Files node.
 4. Drag the **Peptide Validator** node to the workspace, and connect it to the PSM Grouper node.
 5. Drag the **Peptide And Protein Filter** node to the workspace, and connect it to the Peptide Validator node.

6. Drag the **Protein Scorer** node to the workspace, and connect it to the Peptide And Protein Filter node.
7. Drag the **Protein Grouping** node to the workspace, and connect it to the Protein Scorer node.
8. Drag the **Peptide and Protein Quantifier** node to the Workflow Tree pane, and connect it to the Protein Grouping node.
9. (Optional) Drag the **Data Distributions** node to the Post-Processing Nodes pane.
This node adds heat maps, which aid in data review.


Figure 267 illustrates the basic consensus workflow for precursor ion quantification.

Figure 267. Basic consensus workflow for precursor ion quantification




10. (Optional) Add any other appropriate nodes and connect all the nodes together.
11. Set the parameters for each node.

For information about all the parameters that you set for each node, refer to the Help.

12. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the consensus workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the consensus workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save** .
 - d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **ConsensusWorkflow File (*.pdConsensusWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdConsensusWF file.

13. Save the analysis. See “[Saving an Analysis](#)” on [page 80](#).
14. Save the study. See “[Saving a Study](#)” on [page 62](#).
15. Click the **Run** icon,  **Run** , in the upper right corner of the Analysis window.


The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

Performing Reporter Ion Quantification

This topic describes the general procedure to use to perform reporter ion quantification and then describes the processing and consensus workflows in detail.

❖ To perform reporter ion quantification

1. Create or open a study and an analysis:
 - To create a study, see “[Creating a Study](#)” on [page 38](#).
 - To open an existing study, see “[Opening an Existing Study](#)” on [page 40](#).
 - To create an analysis, see “[Creating an Analysis](#)” on [page 71](#).
 - To open an existing analysis, see “[Opening an Existing Analysis](#)” on [page 72](#).
2. Select the quantification method. For information, see “[Setting Up the Quantification Method](#)” on [page 408](#). If the appropriate quantification method is not available, create one by following the instructions in “” on [page 419](#).
3. Create a reporter ion quantification workflow. For information, see “[Creating a Processing Workflow for Reporter Ion Quantification](#),” and “[Creating a Consensus Workflow for Reporter Ion Quantification](#)” on [page 400](#).

4. Specify the quantification channels. For information, see “[Specifying the Quantification Channels](#)” on [page 410](#).
5. Save the analysis. See “[Saving an Analysis](#)” on [page 80](#).
6. Save the study. See “[Saving a Study](#)” on [page 62](#).
7. Click the **Run** icon,  **Run**, in the upper right corner of the Analysis window.

The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

The application creates a .pdResult file, which has the same name as the input file but with a .pdResult extension, and displays it on the Analysis Results page.

Creating a Processing Workflow for Reporter Ion Quantification

To create a reporter ion quantification method, you must set up a processing workflow that includes the Reporter Ions Quantification node.

Default processing workflows are available for several instruments. You can use these default workflows as is or modify them to suit your needs. To access them, see “[Using Common Workflow Templates](#)” on [page 119](#). When you use a default workflow, you do not need to perform [step 1](#) through [step 7](#) in the following procedure because these steps are already included in the workflow. However, you must check the node parameters in the workflow to make sure they are appropriate for your data and add modifications to the search for the quantification to be successful.

Setting up the workflow for TMT and iTRAQ quantification is basically the same.

❖ To create a processing workflow for reporter ion quantification

Note This procedure uses a TMTe 6plex example.

1. In the Workflow Editor, follow the general instructions for creating a processing workflow. See “[Creating a Processing Workflow](#)” on [page 104](#).

The processing workflow must include the following nodes as a minimum:

- Spectrum files node
 - Spectrum Selector node
 - A search engine node (Mascot or Sequest HT)
 - Reporter Ions Quantifier node
 - PSM Validation node
2. In the Processing Workflow window, drag the **Spectrum Files** node from the Data Input area of the Workflow Nodes pane to the Workflow Tree pane.

3. Drag the **Spectrum Selector** node to the workspace.
4. Connect the Spectrum Selector node to the Spectrum Files node.
5. Drag the **Reporter Ions Quantifier** node to the Workflow Tree pane, and attach it directly to the Spectrum Files node.


The Reporter Ions Quantifier node performs quantification for isobarically labeled amino acids.

Note You cannot use the Reporter Ions Quantifier node in a workflow that includes the Precursor Ions Quantifier node.

6. Drag the appropriate search engine node—for example, **Sequest HT**—to the Workflow Tree pane, and attach it to the Spectrum Selector node.
7. Drag the **Fixed Value PSM Validator** node, **Target Decoy PSM Validator** node, or the **Percolator** node to the Workflow Tree pane and attach it to the search engine node.

Figure 268 illustrates the basic processing workflow for reporter ion quantification.

Figure 268. Basic processing workflow for reporter ion quantification



8. (Optional) Add any other nodes that you want and connect them to the workflow.
9. In the Parameters pane of the Workflow Editor, click **Show Advanced Parameters**.
10. Click the search engine node—for example, **Sequest HT**—and set the parameters for it in the Parameters pane:
 - a. In the Protein Database box, select the appropriate FASTA database.



- b. In the Dynamic Modifications area, select the dynamic modifications from these options:

Note If you are using the Mascot search engine node, you can use the Quan Modifications parameter rather than the Dynamic Modifications parameters to specify the modifications to search for. For instructions on using this method, see [“Checking the Quantification Method”](#) on [page 424](#).

- TMT 2plex:
 - TMT 2plex for lysine and N-terminal (Use these as static or dynamic modifications.)
 - Static or dynamic TMT 2plex for threonine
- TMTe 6plex or TMT 6plex:
 - TMT 6plex for lysine and N-terminal (Use these as static or dynamic modifications.)
 - Static or dynamic TMT 6plex for threonine
- TMT 10plex: the same modifications as for TMT 6plex
- iodo TMT 6plex: iodo TMT 6plex for cysteine (Use these as static or dynamic modifications.)

For example, for TMTe 6plex quantification, you would select a static or dynamic modification of **TMT6plex / +229.163 Da (K)**. If you do not find this label, you can enable it by following the instructions in [“Defining Chemical Modifications”](#) on [page 227](#).

- c. In the Static Modifications area, select the static modifications. For example, for TMTe 6plex quantification, you would select **TMT6plex / +229.163 Da** in the Peptide N-Terminus box.
- d. Set any other parameters that you prefer. For information on the available search engine parameters, refer to the Help.
11. Click the **Reporter Ions Quantifier** node and set the parameters for it in the Parameters pane.
- For information on the parameters that you can set for the Reporter Ions Quantifier node, refer to the Help.
12. Set the parameters for all other nodes in the Parameters pane.
- For information about the parameters that you can set for each node, refer to the Help.
13. (Optional) Save the workflow:
- a. In the Name box above the Workflow Tree pane, type a name for the processing workflow.

- b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the processing workflow.
- c. In the Workflow Editor, click the **Save** icon,  **Save** , or the **Save Common** icon,  **Save Common** .
- d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Processing Workflow File (*.pdProcessingWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdProcessingWF file.

Creating a Consensus Workflow for Reporter Ion Quantification

To use a reporter ion quantification method, you must use a consensus workflow that includes the Peptide and Protein Quantifier node. This workflow is the same as that for precursor ion quantification but uses different parameters of the Peptide and Protein Quantifier node.

Default consensus workflows are available. You can use these default workflows as is or modify them to suit your needs. To access them, see [“Using Common Workflow Templates”](#) on [page 119](#). When you use a default workflow, you do not need to perform [step 1](#) through [step 8](#) in the following procedure because these steps are already included in the workflow. However, you must check the node parameters in the workflow to make sure they are appropriate for your data and add modifications to the search for the quantification to be successful.

Setting up the workflow for TMT and iTRAQ quantification is basically the same.

❖ To create a consensus workflow for reporter ion quantification

1. Follow the general instructions for creating a consensus workflow. See [“Creating a Consensus Workflow”](#) on [page 112](#).

The consensus workflow must include the following nodes as a minimum:


- MSF Files node
- PSM Grouper node
- Peptide Validator node
- Peptide and Protein Filter node
- Protein Scorer node
- Protein Grouping node
- Peptide and Protein Quantifier node

2. In the Processing Workflow window, drag the **MSF Files** node to the workspace.
3. Drag the **PSM Grouper** node to the workspace, and connect it to the MSF Files node.
4. Drag the **Peptide Validator** node to the workspace, and connect it to the PSM Grouper node.
5. Drag the **Peptide and Protein Filter** node to the workspace, and connect it to the Peptide Validator node.
6. Drag the **Protein Scorer** node to the workspace, and connect it to the Peptide and Protein Filter node.
7. Drag the **Protein Grouping** node to the workspace, and connect it to the Protein Scorer node.
8. Drag the **Peptide and Protein Quantifier** node to the Workflow Tree pane, and connect it to the Protein Grouping node.
9. (Optional) Drag the **Data Distributions** node to the Post-Processing Nodes pane.

This node adds heat maps, which aid in data review.

Figure 266 illustrates the basic consensus workflow for reporter ion quantification.

Figure 269. Basic consensus workflow for reporter ion quantification




10. (Optional) Add any other appropriate nodes and connect all the nodes together.

11. Set the parameters for each node.

For information about all the parameters that you set for each node, refer to the Help.

12. (Optional) Save the workflow:


- a. In the Name box above the Workflow Tree pane, type a name for the consensus workflow.
- b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the consensus workflow.
- c. In the Workflow Editor, click the **Save** icon,  **Save** .
- d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **ConsensusWorkflow File (*.pdConsensusWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdConsensusWF file.

Performing TMT Quantification on HCD and CID Scans

If a raw data file contains both CID scans for identification and HCD scans for quantification, you can use the following workflow to both quantify the HCD scans and identify peptides in the CID scans, the HCD scans, or both.


❖ To perform TMT Quantification on HCD and CID scans

1. Drag the **Reporter Ions Quantifier** node to the Workflow Tree pane and connect it to the workflow.
2. Set the Activation Type parameter for the Reporter Ions Quantifier node to **HCD**.
3. Set the Activation Type parameter for the Spectrum Selector node to **Any, Is CID, HCD**, or **Is CID**, depending on your method setup and identification strategy.
4. Set all other parameters—modifications, tolerances, FASTA files, and so forth.
5. Click the **Run** icon,  **Run** .

Performing Precursor Ion Area Detection

This topic describes the general procedure to use to perform precursor ion area detection and then describes the processing and consensus workflows in detail.

❖ To perform precursor ion area detection

1. Create or open a study and an analysis:
 - To create a study, see “[Creating a Study](#)” on [page 38](#).
 - To open an existing study, see “[Opening an Existing Study](#)” on [page 40](#).
 - To create an analysis, see “[Creating an Analysis](#)” on [page 71](#).
 - To open an existing analysis, see “[Opening an Existing Analysis](#)” on [page 72](#).
2. Select the quantification method. For information, see “[Setting Up the Quantification Method](#)” on [page 408](#).
3. Create a precursor ion area detection workflow. For information, see “[Creating a Processing Workflow for Precursor Ion Area Detection](#)” on [page 403](#) and “[Creating a Consensus Workflow for Precursor Ion Area Detection](#)” on [page 406](#).
4. Specify the quantification channels. For information, see “[Specifying the Quantification Channels](#)” on [page 410](#).
5. Save the analysis. See “[Saving an Analysis](#)” on [page 80](#).
6. Save the study. See “[Saving a Study](#)” on [page 62](#).
7. Click the **Run** icon,  , in the upper right corner of the Analysis window.

The job queue appears, as shown in [Figure 79](#) on [page 112](#), displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

The application creates a .pdResult file, which has the same name as the input file but with a .pdResult extension, and displays it on the Analysis Results page.

Creating a Processing Workflow for Precursor Ion Area Detection

To create a precursor ion area detection method, you must set up a processing workflow that includes the Precursor Ions Area Detector node.

❖ To create a processing workflow for precursor ion area detection

1. Follow the general instructions for creating a processing workflow in “[Creating a Processing Workflow](#)” on [page 104](#).

The processing workflow must include the following nodes as a minimum:

- Spectrum Files node
- Spectrum Selector node
- Event Detector node
- Precursor Ions Area Detector node


- Search engine node (Sequest HT or Mascot)
 - Fixed Value PSM Validator or Percolator node
2. From the Data Input area of the Workflow Nodes pane, drag the **Spectrum Files** node to the Workflow Tree pane.
 3. Drag the **Spectrum Selector** node and the **Event Detector** node to the workspace.
 4. Connect the Spectrum Selector node and the Event Detector node to the Spectrum Files node.
 5. Drag the **Precursor Ions Area Detector** node to the Workflow Tree pane and attach it directly to the Event Detector node.

The Precursor Ions Area Detector node calculates the area of each precursor ion.

6. Drag the search engine node that you prefer—for example, **SequestHT**—to the Workflow Tree pane and attach it to the Spectrum Selector node.
7. Drag the **Fixed Value PSM Validator**, **Target Decoy PSM Validator**, or the **Percolator** node to the Workflow Tree pane and attach it to the search engine node.

Figure 270 illustrates the basic processing workflow for precursor ion area detection.

Figure 270. Basic processing workflow for precursor ion area detection



8. Add any other appropriate nodes and connect all the nodes together.
9. Click the **Spectrum Files** node and specify the raw data file in the Parameters pane.
10. Click the **Event Detector** node and set the parameters for it in the Parameters pane:

- a. In the Mass Precision box, specify the expected standard deviation of the mass precision.

To create extracted ion chromatograms, use three times the standard deviation. The minimum value is 0.5 ppm. The maximum value is 4 ppm. The default is 2 ppm.

- b. In the S/N Threshold box, specify a threshold signal-to-noise value that determines whether the application removes peaks from the spectrum.

The application removes peaks with a signal-to-noise value below this threshold. The minimum value is 0.0, and there is no maximum value. The default is 1.

11. Click the **Spectrum Selector** node, and set the parameters for it in the Parameters pane:

- a. Change the setting in the Max. Precursor Mass box to an appropriate setting. For example, you can set this option to **6500**.
- b. Change the setting in the S/N Threshold box to an appropriate setting. For example, you can use the default of **1.5**.

For other optional parameters that you can set for the Spectrum Selector node, refer to the Help.

12. Click the search engine node—for example, **Sequest HT**—and set the parameters for it in the Parameters pane:

- a. In the Protein Database box, select an appropriate FASTA database.
- b. In the Dynamic Modifications area, select the dynamic modifications.

Note If you are using the Mascot search engine node, you can use the Quan Modifications parameter rather than the Dynamic Modifications parameters to specify the modifications to search for. For instructions on using this method, see [“Checking the Quantification Method”](#) on page 424.

For example, you might select the **Oxidation+15.995 Da (M)** dynamic modification. If you do not find this label, you can enable it by following the instructions in [“Defining Chemical Modifications”](#) on page 227.

- c. In the Static Modifications area, select the static modifications.



For example, you might select **Carbamidomethyl/ +57.021 Da (C)** in the Static Modification box.

- d. Set any other parameters that you prefer. For information on the available search engine parameters, refer to the Help.

13. Set the parameters for all other nodes in the Parameters pane.

For information about all the parameters that you can set for each node, refer to the Help.

The Precursor Ions Area Detector node has no parameters.

14. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the processing workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the processing workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save** , or the **Save Common** icon,  **Save Common** .
 - d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Processing Workflow File (*.pdProcessingWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdProcessingWF file.

Creating a Consensus Workflow for Precursor Ion Area Detection

To create a precursor ion area detection method, you must set up a consensus workflow that includes the Peptide and Protein Quantifier node.

❖ To create a consensus workflow for precursor ion area detection

1. Follow the general instructions for creating a consensus workflow with the Workflow Editor. See “[Creating a Consensus Workflow](#)” on [page 112](#).

The consensus workflow must include the following nodes as a minimum:


- MSF Files node
 - PSM Grouper node
 - Peptide Validator node
 - Peptide And Protein Filter node
 - Protein Scorer node
 - Protein Grouping node
 - Peptide and Protein Quantifier node
2. In the Processing Workflow window, drag the **MSF Files** node to the workspace.
 3. Drag the **PSM Grouper** node to the workspace, and connect it to the MSF Files node.
 4. Drag the **Peptide Validator** node to the workspace, and connect it to the PSM Grouper node.

5. Drag the **Peptide and Protein Filter** node to the workspace, and connect it to the Peptide Validator node.
6. Drag the **Protein Scorer** node to the workspace, and connect it to the Peptide and Protein Filter node.
7. Drag the **Protein Grouping** node to the workspace, and connect it to the Protein Scorer node.
8. Drag the **Peptide and Protein Quantifier** node to the Workflow Tree pane, and connect it to the Protein Grouping node.
9. (Optional) Drag the **Data Distributions** node to the Post-Processing Nodes pane.

This node adds heat maps, which aid in data review.



Figure 271 illustrates the basic consensus workflow for precursor ion area detection.

Figure 271. Basic consensus workflow for precursor ion area detection



10. (Optional) Add any other nodes that you want and connect all the nodes together.
11. Set the parameters for each node.

For information about all the parameters that you set for each node, refer to the Help.

12. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the consensus workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the consensus workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save** , or the **Save Common** icon,  **Save Common** .
 - d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **ConsensusWorkflow File (*.pdConsensusWF)**.
 - iii. Click **Save**.

The application saves the workflow in the *file_name*.pdConsensusWF file.

Using a Quantification Method

This topic explains how to set up, add, check, change, remove, import, and export a quantification method.

- [Setting Up the Quantification Method](#)
- [Specifying the Quantification Channels](#)
- [Checking the Quantification Method](#)
- [Changing a Quantification Method](#)
- [Removing a Quantification Method](#)
- [Importing a Quantification Method](#)
- [Exporting a Quantification Method](#)
- [Restoring Quantification Method Template Defaults](#)

Setting Up the Quantification Method

Setting up the quantification method is similar for both precursor ion quantification and reporter ion quantification. Both methods use values called quantification (quan) channels as the basis for the ratio reporting. You do not need to set up a quantification method for precursor ion area detection quantification.

For reporter ion quantification, a quantification channel is one of several masses, states, or tags (depending on which quantification method you use) for which you measure a quantification value. The application calculates the reported quantification ratios from the quantification values of the different quantification channels. For example, for iTRAQ 4plex, the different reporter tags (114, 115, 116, 117) are the four quantification channels of the iTRAQ 4plex method. The application calculates the ratios from the detected quantification values of the four quantification channels.

For precursor ion quantification, a quantification channel is one of the different possible labeling states of a peptide corresponding to the different heavy amino acids used in the cell cultures. For example, the SILAC 2plex methods are normally used with two quantification channels named “light” and “heavy.” The light quantification channel uses the natural isotopes of lysine ($^{12}\text{C}_6^{14}\text{N}_2$) and arginine ($^{12}\text{C}_6^{14}\text{N}_4$). In the heavy quantification channel, arginine 10 ($^{13}\text{C}_6^{15}\text{N}_4$) replaces all arginines, and either lysine 6 ($^{13}\text{C}_6^{14}\text{N}_2$) or lysine 8 ($^{13}\text{C}_6^{15}\text{N}_2$) replaces all lysines.


❖ **To set up the quantification method**

On the Study Definition page of the study, select the check box of the appropriate quantification method or methods in the Quantification Methods area, shown in [Figure 272](#).

The example shows the iodo TMT 6plex method selected.

If the quantification method that you want to use is not in the list, you must create it. For instructions, see [“Specifying the Quantification Channels”](#) on [page 410](#).

Figure 272. Quantification Methods area of the Study Definition page



If you use an input file in one processing step and another input file in another processing step, you can select a different quantification method for each file.

Note Once you add a quantification method to a study, it is completely separated from the original method on the server. Edits that you make to the quantification method in the study do not propagate to the method on the server and vice versa, although they seem to be the same at first glance. You must be able to copy a study in the study folder to another computer that might contain different quantification methods or none at all. For the study to work on the new computer, the application saves the quantification method with the study rather than on the server.

Specifying the Quantification Channels

The first step in setting up the quantification is to specify the quantification channels to use. The process of specifying label modifications is similar for precursor ion quantification and reporter ion quantification, but it also has some differences:

- For precursor ion quantification, you specify the label modifications for each quantification channel. For reporter ion quantification, you set the label modification that produces all reporter ion quantification channels.
- For precursor ion quantification, specifying the label modifications for quantification channels other than the unlabeled channel is mandatory. For reporter ion quantification, specifying the label modifications is optional because the information about the modification of the peptides is not necessary for processing the data. It is only used to verify the peptides when the application loads the reports.

When you specify at least one of the label modifications in the quantification method, the application verifies that each identified peptide has at least one of the specified modifications during processing.

- When the application identifies a peptide with none of the specified label modifications, this peptide cannot be the source of reporter peaks in the MS/MS spectra. As a result, the application marks the peptide “No Quan Labels” in the .pdResult file. It does not use these peptides when it calculates the protein quantification values.
- When the application finds a peptide that does not have an iTRAQ or TMT label as a modification, even though reporter ions were present, it leaves the *Ratio* columns in the .pdResult file blank.


When you install the application, the default methods for TMT and iTRAQ do not include the correct label modification. When you set up a workflow, ensure that these modifications are selected as either fixed or variable. The application does not automatically update already existing reporter methods; you must manually specify the label modifications.

When you open old .pdResult files that contain reporter quantification data, the label modifications of the quantification method of the .pdResult file appear as None on the Quan Channels page of the Quantification Method Editor dialog box. You can manually specify the label modification, which then triggers the validation of the peptides, and save the change in the quantification method in the .pdResult file.

When you do not set the label modifications on the Quan Channels page, the application does not perform the validation.


The process of specifying quantification channels for precursor ion quantification is slightly different from the process of specifying label modifications for reporter ion quantification.

❖ To specify quantification channels for precursor ion quantification

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon, .

The Quantification Methods view opens, as shown in [Figure 273](#). It lists all of the available methods for both precursor ion and reporter ion quantification.

Figure 273. Quantification Methods view





The Status column indicates whether the quantification method is valid for use in quantification:

- A green check mark means that the quantification method is valid and can be used for quantification.
- An exclamation point in a yellow triangle means that the quantification method is not valid. Double-click this mark to view a message that describes the error and provides information on how to fix it.

Figure 274 gives examples of these symbols in the Status column.


Figure 274. Method validity symbols in the Quantification Methods view



2. Double-click the row for the appropriate method in the Method Name or Description column, or click the column to the left of the Status column for the method and click the **Add** icon,  **Add**.

The Create New Quantification Method dialog box now appears, as shown in Figure 275.


Figure 275. Create New Quantification Method dialog box



3. Select one of the following methods of creating a quantification method:
 - (Default) From Factory Defaults: Creates a new method using the same settings from one of the default settings that come with a new application installation.
 - From Existing Method: Uses the same settings as those of the existing quantification method that you select from the list. The list of methods is the same as that given at the beginning of “[Setting Up the Quantification Method](#)” on [page 408](#).
 - From Scratch: Uses one of the following templates so that you can build a new processing method from the beginning:
 - Reporter Ion Quan Method: Provides a template for reporter ion quantification. This method requires an additional license.
 - Precursor Ion Quan Method: Provides a template for precursor ion quantification.
4. Click **Create**.

The Quantification Method Editor dialog box appears, as shown in [Figure 276](#).

Figure 276. Quantification Method Editor dialog box for precursor ion quantification




5. If you want the application to ignore modifications on the N terminus, clear the **Quan Labels Themselves Are Modifications on a Side Chain or the N-Terminus** check box.

By default, this check box is selected. Select this check box for dimethylation and clear it for SILAC modifications.

6. In the box at the left, select the name of the labeling method to use:
 - Heavy: Refers to amino acid labels that use heavy isotopes, for example, Arg10 and Lys8.
 - Medium (3plex methods only): Refers to amino acid labels that use less massive isotopes, for example, Arg6 and Lys4.
 - Light: Refers to amino acid labels that use normal isotopes.
7. To add a quantification channel, click + beneath the list of quantification channels in the box on the left.

The default name of New *number* now appears in the list of quantification channels and in the Channel Name box, as shown in [Figure 277](#).

Figure 277. New quantification channel on the Quan Channels page for precursor ion quantification



To remove a quantification channel, select the quantification channel in the list of quantification channels and click – beneath the list.

8. To specify a name for the new quantification channel, backspace over the default name in the Channel Name box and type the new name.

The new name now appears in the quantification channel (left) box.

9. To specify a quantification label to assign to a quantification channel, click + beneath the Quantification Labels box.

A default quantification label of New *number* now appears in the Quantification Labels box and the Label Name box.


To remove an existing quantification label, select the label in the Quantification Labels box and click – beneath the box.

10. To change the default quantification channel name, backspace over the name in the Label Name box and type the new name. The example in [Figure 278](#) uses Arg6, Lys6.
11. In the Modification Target area, select the location of the label on the peptide:
 - Side Chain Modification: Indicates that the label occurs on a side chain.
 - N-Terminal Modification: Indicates that the label occurs on the N terminus.
 - C-Terminal Modification: Indicates that the label occurs on the C terminus.

- From the Modification list, select the modification for the amino acid label. The example in [Figure 278](#) shows Label:13C(6) / +6.020 Da.
- From the list adjacent to the Modification list, select the abbreviation of the amino acid selected in the Quantification Labels box that the modification should occur on. In this example, K is selected.


The completed Quan Channels page resembles [Figure 278](#).

Figure 278. Completed Quan Channels page



- When you finish setting the parameters in the Quantification Method Editor, click **OK**. The Save Quantification Method dialog box opens, as shown in [Figure 279](#).


Figure 279. Save Quantification Method dialog box




- In the Save as New Method box, type the name of the new method.
- Click **Save**.

The new quantification method now appears in the Quantification Methods view.

❖ **To specify quantification channels for reporter ion quantification**

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon, .

The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification. For information about the Status column in this view, see [“To specify quantification channels for precursor ion quantification”](#) on [page 411](#).


2. Double-click the row for the appropriate method in the Method Name or Description column, or click the column to the left of the Status column for the method and click the **Add** icon, .

The Create New Quantification Method dialog box now appears, as shown in [Figure 275](#) on [page 413](#).

3. Select one of the following methods of creating a quantification method:
 - (Default) From Factory Defaults: Creates a new method using the same settings from one of the default settings that come with a new application installation.
 - From Existing Method: Uses the same settings as those of the existing quantification method that you select from the list. The list of methods is the same as that given at the beginning of [“Setting Up the Quantification Method”](#) on [page 408](#).
 - From Scratch: Uses one of the following templates so that you can build a new processing method from the beginning:
 - Reporter Ion Quan Method: Provides a template for reporter ion quantification. This method requires an additional license.
 - Precursor Ion Quan Method: Provides a template for precursor ion quantification.
4. Click **Create**.


The Quantification Method Editor dialog box appears, as shown in [Figure 280](#).

Figure 280. Quantification Method Editor dialog box for reporter ion quantification



5. From the Residue Modification list, select the label modification to be found on the target amino acid residue. From the adjacent list, select the appropriate letter to indicate that the modification should occur on the indicated residue and will have an increased mass.
6. From the N-Terminal Modification list, select the label modification to be found on the N terminus of each peptide.
7. In the TMT Reporter Ion Isotope Distributions or iTRAQ Software Correction Factors area, place the isotope intensities for specific isotope shifts in the following columns:
 - a. In the -2 column, type the value that is 2 Da lower than that of the reporter ion.
 - b. In the -1 column, type the value that is 1 Da lower than that of the reporter ion.
 - c. In the +1 column, type the value that is 1 Da higher than that of the reporter ion.
 - d. In the +2 column, type the value that is 2 Da higher than that of the reporter ion.
 - e. Select the check box in the Active column if you want to use that reporter ion in the method.
 - f. If only some, but not all, of the tags are used—for example, all but 129—clear the active check box for that reporter ion mass.
 - g. Click **OK**.

The Save Quantification Method dialog box opens.

Figure 281. Save Quantification Method dialog box

8. In the Save as New Method box, type the name of the new method.
9. Click **Save**.

The new quantification method now appears in the Quantification Methods view.

To create a new quantification method to correct for isotopic impurities in TMT 10plex kits, see [“Excluding PSMs with High Levels of Coisolation”](#) on page 435.

Correcting Reporter Ion Quantification Results for Isotopic Impurities

This topic explains how to set up a quantification method for isotopic impurities.

When you buy a TMT kit from Pierce, you receive a product data sheet called the Certificate of Analysis (COA), which summarizes the essential information about the kit and the lot that it was taken from. The COA includes a table with the exact masses of the different labels and their isotopic impurities. For reporter ion methods, you can add quantification channel data like monoisotopic mass, isotopic impurities, and associated isotopic variant relations as they are stated in the COA to the corresponding quantification method. With existing reporter ion quantification methods, you can only change the values of the isotopic impurities. The impurity values of the default methods are zero, so you must first enter the correct values from the COA to make the purity corrections work.

You can create a new method for correcting reporter ion quantification results through any of the following means:


- Add a new method based on a factory default method as a template and fill in the reporter ion correction factor table.
- Use an existing method as a template.
- Create a new method from the beginning. This method requires an additional license.

In this new method, you can use all of the available labels in a kit or only a subset.

When you create a new quantification method from the factory default template or an existing method, you cannot change the reporter ion masses, change the mass tags, or edit the relations among the impurities. For instructions on creating a new method from the first two templates, see [“Specifying the Quantification Channels”](#) on page 410. To create a new method from the beginning, see the next topic.

❖ **To create a quantification method for correcting reporter ion quantification results**

Note The following procedure uses a TMT 10plex kit as an example.


1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon, .

The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.

2. Click the **Add** icon,  **Add**.

The Create New Quantification Method dialog box now opens, as shown in [Figure 282](#).

Figure 282. Create New Quantification Method dialog box



3. Select the **From Scratch** option.


Note The From Scratch option requires a separate license.

4. Select **Reporter Ion Method**.

A Channels box opens to the right of the method selection list.

5. In the Channels box, specify the number of channels, in this case, **10**, as shown in [Figure 283](#).

Figure 283. Create New Quantification Method dialog box for creating TMT 10plex method




6. Click **Create**.

The Quantification Method Editor dialog box opens, as shown in [Figure 284](#). The Quantification Method Editor dialog box for reporter ion quantification is different from the Quantification Method Editor dialog box for precursor ion quantification (see [Figure 276](#) on page 414).

The Quantification Method Editor dialog box for reporter ion quantification displays channels with no masses in red until you input the masses (see [Figure 284](#)).

Figure 284. Quantification Method Editor dialog box showing channels with no masses



7. From the Residue Modification list, select **TMT6plex/+229.163 Da**.
8. From the list next to the Residue Modification box, select **K**.
9. From the N-Terminal Modification list, select **TMT6plex/+229.163 Da**.
10. Type the data from the TMT 10plex kit into the appropriate columns of the Quantification Method Editor dialog box, as shown in [Figure 285](#).

Note You must enter data in the Mass Tag and Reporter Ion Mass columns before you can select items from the lists in the Affects columns.

10 Performing Quantification

Using a Quantification Method

Figure 285. Quantification Method Editor dialog box with values entered from a TMT 10plex kit

Quantification Method Editor: newMethod1

Quan Channels

Residue Modification: TMT6plex / +229.163 Da K

N-Terminal Modification: TMT6plex / +229.163 Da

Reporter Ion Isotopic Distribution Correction Factors

Mass Tag	Reporter Ion Mass	- 2	Affects	- 1	Affects	Main	+ 1	Affects	+ 2	Affects	Active
126	126.127726	0	N/A	0	N/A	100	4.7	127C	0	128C	<input checked="" type="checkbox"/>
127N	127.124761	0	N/A	0.4	126	100	6.5	128N	0	129N	<input checked="" type="checkbox"/>
127C	127.131081	0	N/A	0.2	126	100	4.6	128C	0.3	129C	<input checked="" type="checkbox"/>
128N	128.128116	0	126	0.9	127N	100	4.7	129N	0.2	130N	<input checked="" type="checkbox"/>
128C	128.134436	0.1	126	0.5	127C	100	2.6	129C	0	130C	<input checked="" type="checkbox"/>
129N	129.131471	0	127N	0.7	128N	100	2.5	130N	0	131	<input checked="" type="checkbox"/>
129C	129.13779	0	127C	1.3	128C	100	2.5	130C	0	N/A	<input checked="" type="checkbox"/>
130N	130.134825	0	128N	1.2	129N	100	2.8	131	2.7	N/A	<input checked="" type="checkbox"/>
130C	130.141145	0.1	128C	2.9	129C	100	2.9	N/A	0	N/A	<input checked="" type="checkbox"/>
131	131.13818	0	129N	2.4	130N	100	1.4	N/A	0	N/A	<input checked="" type="checkbox"/>

Reporter quan method

OK Cancel Help

The application displays channels that are not active with a gray background, as shown in Figure 286.

Figure 286. Quantification Method Editor dialog box showing inactive channels

Quantification Method Editor: newMethod1

Quan Channels

Residue Modification: TMT6plex / +229.163 Da K

N-Terminal Modification: TMT6plex / +229.163 Da

Reporter Ion Isotopic Distribution Correction Factors

Mass Tag	Reporter Ion Mass	- 2	Affects	- 1	Affects	Main	+ 1	Affects	+ 2	Affects	Active
126	126.127726	0	N/A	0	N/A	100	4.7	127C	0	128C	<input checked="" type="checkbox"/>
127N	127.124761	0	N/A	0.4	126	100	6.5	128N	0	129N	<input checked="" type="checkbox"/>
127C	127.131081	0	N/A	0.2	126	100	4.6	128C	0.3	129C	<input checked="" type="checkbox"/>
128N	128.128116	0	126	0.9	127N	100	4.7	129N	0.2	130N	<input checked="" type="checkbox"/>
128C	128.134436	0.1	126	0.5	127C	100	2.6	129C	0	130C	<input checked="" type="checkbox"/>
129N	129.131471	0	127N	0.7	128N	100	2.5	130N	0	131	<input type="checkbox"/>
129C	129.13779	0	127C	1.3	128C	100	2.5	130C	0	N/A	<input type="checkbox"/>
130N	130.134825	0	128N	1.2	129N	100	2.8	131	2.7	N/A	<input checked="" type="checkbox"/>
130C	130.141145	0.1	128C	2.9	129C	100	2.9	N/A	0	N/A	<input checked="" type="checkbox"/>
131	131.13818	0	129N	2.4	130N	100	1.4	N/A	0	N/A	<input checked="" type="checkbox"/>


Reporter quan method

OK Cancel Help

11. Click **OK**.

The Save Quantification Method dialog box appears, as shown in [Figure 287](#).

Figure 287. Save Quantification Method dialog box




12. Type the name of the new method in the Save as New Method box, and click **Save**.

The name of the new method now appears in the Method Name column in the Quantification Methods view.

When you edit factory-provided methods, the application automatically sets the mass tag relations and stores them in the method definition. You cannot edit them when you create a new method from the beginning. If you want to edit a mass tag, you must create a new method.

When you double-click the name of the new method in the Quantification Methods view, the Quantification Method Editor opens, but now a usage status column replaces the Active column. It displays text indicating whether each channel is used or not used, as shown in [Figure 288](#).

Figure 288. Quantification Method Editor showing usage status column



Usage status column

Checking the Quantification Method

The application checks the parameters that you have set for the quantification method. For reporter ion quantification, it verifies that the method has at least two channels. For precursor ion quantification, it checks for the following:

- At least one quantification channel.
- At least one label for each quantification channel. The label can be None, which the default SILAC methods use.
- Unique label names in a channel.
- The modification of each label applied to at least one amino acid, unless you chose None for a modification.
- Each amino acid labeled only once in a channel (labels must have a defined elemental composition).
- Each label mass used only once (label masses vary by at least 1 Da).

You cannot apply changes to a quantification method unless the method meets all these criteria.

Changing a Quantification Method


You can change a quantification method before you generate the results report. You can also change a quantification method after you generate a results report, but you can make only limited changes.

Changing a Quantification Method Before Generating a Results Report


Follow these steps to change a quantification method before generating a results report.

❖ To change an existing quantification method before generating a results report


1. Open the Quantification Method Editor dialog box by doing one of the following:

In the Quantification Methods area of the Study Definition page of a study, select the check box for the appropriate quantification method, and click the **Edit Quantification Method** icon, , for the selected quantification method.

—or—

- a. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon, .

The Quantification Methods view appears, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.

- b. Double-click the row for the appropriate method in the Method Name or Description column, or click the column to the left of the Status column for the method and click the **Edit** icon,  **Edit**.

The Quantification Method Editor dialog box appears, as shown in [Figure 276](#) on [page 414](#) for precursor ion quantification or [Figure 280](#) on [page 418](#) for reporter ion quantification.

2. In the Quantification Method Editor dialog box, make any appropriate modifications.

For information on setting the parameters in this dialog box, see “[Specifying the Quantification Channels](#)” on [page 410](#).

3. Click **OK**.

The application checks the parameters that you have changed to be sure that they conform to the guidelines given in “[Checking the Quantification Method](#)” on [page 424](#). It does not apply the changes to a quantification method unless the method meets all these criteria.

Changing a Quantification Method After Generating a Results Report


After you perform quantification, you can change the quantification method and reprocess the data.

❖ To change a quantification method after generating a results report


1. On the Input Files page of the study, choose a new quantification method in the Quan Method box.

–or–

Open the Quantification Method Editor by doing the following.


- a. In the Quantification Methods area of the Study Definition page, click the **Edit Quantification Method** icon, , for the selected quantification method.


–or–

Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon, , and then double-click one of the methods shown in the Method Name column of the Quantification Methods view.

- b. In the Quantification Method Editor dialog box, add, or change the modifications for the preexisting channels.

However, you cannot define mass tags or labels as you can when setting up the initial quantification method, because they have already been measured.


2. Click the **Analysis Results** tab.
3. Click the **Reprocess** icon,  **Reprocess**.

4. Choose **All Analysis Steps**.
5. In the Analysis window, click the **Run** icon,  .

Removing a Quantification Method


You can delete a quantification method if it is no longer useful, or make a quantification method temporarily unavailable to new workflows.


❖ To remove a quantification method

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon,  .


The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.

2. Select the row of the method that you want to remove.

The Remove icon,  **Remove** , now becomes available.

3. Click the **Remove** icon,  **Remove** .
4. In the Delete Methods dialog box, click **OK**.

❖ To deactivate a quantification method

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon,  .

The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.


2. Clear the check box in the Is Active column on the line containing the quantification method that you want to make inactive.

To make the quantification method active again, select the same check box.

Importing a Quantification Method

You can import a new quantification method from another computer.

❖ To import a quantification method

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon,  , either in the toolbar or on the Administration page.

The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.


2. Click the **Import** icon,  **Import** .

3. In the Import Quan Method dialog box, select the METHOD file containing the method that you want to import, and click **Open**.
 - If the new method is valid, the Quantification Method Editor dialog box opens, showing the new method.
 - If the new method is not valid, a message box appears that describes the error.
4. For a valid method, click **OK** in the Quantification Method Editor dialog box.
5. Change the name of the imported quantification method by changing it in the Method Name column of the Quantification Methods view.


Exporting a Quantification Method

You can save a quantification method to use on another computer.

❖ To export a quantification method

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon,  , either in the toolbar or on the Administration page.

The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.


2. Select the method that you want to export in the Quantification Methods view by clicking in the column to the left of the Status column.
3. Click the **Export** icon,  **Export** .
4. In the Export Quan Method dialog box, browse to the folder where you want to store the METHOD file containing the quantification method, and click **Save**.

The application saves the method in a file with the name of the method and a .method suffix, for example, TMT 2plex.method.


Restoring Quantification Method Template Defaults

If you have altered one of the quantification method templates listed at the beginning of “[Setting Up the Quantification Method](#)” on [page 408](#), you can restore the original template in effect when the application was newly installed.

❖ To restore the original template

1. Choose **Administration > Maintain Quantification Methods**, or click the **Maintain Quantification Methods** icon,  .

The Quantification Methods view opens, as shown in [Figure 273](#) on [page 411](#). It lists all of the available methods for both precursor ion and reporter ion quantification.

2. To open the Quantification Method Editor dialog box, click the **Add** icon,  **Add**, in the Quantification Methods view.

The Create Quantification Method dialog box opens, as shown in [Figure 275](#) on [page 413](#).

3. Select the appropriate template from the Create from Factory Defaults list.
4. Set up the quantification method according to the instructions in “[Setting Up the Quantification Method](#)” on [page 408](#).

Searching for Quantification Modifications with Mascot

When you use the Mascot node as the search engine in a quantification workflow, you can set the dynamic and static modifications as parameters. For samples with isotopic labels and several PTMs, you might need to specify several dynamic modifications usable within a single search, but the current number that you can specify is limited to nine.

To avoid this limitation, you can configure quantification methods on the Mascot server. In a quantification method, modifications are organized into groups classified as fixed, variable, or exclusive. You can define modification groups as variable or exclusive at the component level, where they usually characterize the component. You can also define them at the method level, but only as fixed or variable. Defining modifications at the method level is convenient for modifications that are important to the method and saves having to choose them in the Workflow Editor. Exclusive groups are effectively a choice of fixed modifications, so the restrictions that apply to fixed modifications also apply to them.

With the Mascot node, you can use the modification groups specified as part of a quantification method on the Mascot server. You can use the node’s From Quan Method parameter in the Parameters pane to select the dynamic modifications to search for, rather than manually specifying each modification with a Dynamic Modifications parameter.

In the editor in the Mascot server window, you can specify that these groups be variable, fixed, or exclusive. You can also define them directly for the method in report ion quantification or for each component in precursor ion quantification.

❖ To specify the quantification modifications to search for

1. Choose **Administration > Configuration > Mascot**, and configure Mascot by following the instructions in “[Configuring the Mascot Search Engine](#)” on [page 25](#). In the Mascot Server URL box, be sure that you enter the URL of the Mascot server to be used for Mascot searches.
2. Set up a processing workflow that includes, at a minimum, the nodes shown in these figures:
 - For precursor ion quantification, see [Figure 266](#) on [page 392](#).
 - For reporter ion quantification, see [Figure 268](#) on [page 398](#).

- For precursor ion area detection, see [Figure 270](#) on [page 404](#).

Follow the instructions in these topics to set up the workflow.

3. For the search engine node, drag the **Mascot** node to the workflow, and connect it to the other nodes.
4. Set the parameters of the Mascot node.
5. Select the dynamic modifications to search for from the following:
 - Select a dynamic modification from the list in each Dynamic Modification parameter.

You can select up to nine modifications.

–or–

- Click the **From Quan Method** parameter in the Parameters pane under Modification Groups, and from the list (see [Figure 289](#) for an example), select the modifications that you want to search for.


You can select more than nine modifications.

Note Do not use the modifications that you specify as part of the modification groups in the selected quantification method as additional dynamic or static modifications.

10 Performing Quantification

Searching for Quantification Modifications with Mascot

Figure 289. From Quan Method list



- (Optional) If you want to group these modifications, go to the editor in the Mascot server window and choose **Configuration Editor > Quantitation**.

Once you group the modifications, you can define them as fixed, variable, or exclusive. You can also define them directly for the method in reporter ion quantification or for each component in precursor ion quantification. Refer to the Mascot documentation for information on grouping modifications and defining the groups.

For the final search results, it does not matter whether you explicitly specified a modification as either a dynamic or a static modification, or indirectly specified a modification from the chosen quantification method. As an exception, when you select an exclusive modification group, Mascot modifies all or none of the affected residues of a peptide sequence. Therefore, peptide matches with inconsistent labeling no longer occur.

Note Using a Mascot quantification method to retrieve the modification groups to use does not affect how the application performs the quantification. The application itself exclusively performs the quantification. You must specify in the application's methods any quantification labels used for the quantification.

Generating Quantification Ratios

To generate quantification ratios, see [“Specifying Quantification Ratios from Selected Sample Groups”](#) on page 82.

Calculating Quantification Results

This topic gives an overview of the methodology that the Proteome Discoverer application uses to perform quantification and then examines each step in detail.

Quantification Methodology

The application uses the following general methodology to perform quantification:

1. It calculates the raw quantification values and associates them with the PSMs identified by Sequest HT, Mascot, or other search engines. The quantification values are peak areas for precursor ion quantification, and intensities or S/N values of the mass or reporter tags in reporter ion quantification.
2. It sums the quantification channel values of the PSMs to the peptide groups and proteins. PSMs that do not meet the criteria set by the method, such as minimum S/N threshold or isolation purity, are not included in this sum. Otherwise, this methodology aggregates all contributing signals of the protein, irrespective of the charge and modification state of the peptides. The result is quantification channel abundance values for all channels from all files for the peptides and proteins.
3. The application optionally normalizes the abundance values so that all channels have the same total abundance if you selected a normalization method in the consensus workflow.
4. It optionally scales the normalized abundance values for each protein and peptide so that the average abundance is 100. For instance, if you select a TMT 6plex quantification method, the application scales the calculated abundance values of each protein so that their sum is 600 (100×6). If the abundances are equal, each one is 100.
5. It groups the optionally normalized and scaled abundance values according to the sample grouping that you specified when you set up the analysis. Grouping in this case means effectively averaging the abundance values of all replicates.
6. As the last step, the application calculates the quantification ratios from the grouped (averaged) abundance values according to the ratios selected when you set up the analysis.

The Peptide and Protein Quantifier node in the consensus workflow calculates quantification abundances and ratios from the raw quantification values used in the processing workflows. It calculates quantification abundances from the data processed by the Precursor Ions Quantifier node and the Reporter Ions Quantifier node.

The Peptide and Protein Quantifier node performs the following steps to calculate the quantification ratios:

- [Calculating PSM Abundances](#)
- [Classifying PSMs](#)
- [Calculating Peptide Group Abundances](#)
- [Classifying Peptide Groups](#)
- [Calculating Protein Abundances](#)
- [Normalizing Peptide Groups and Protein Abundances](#)
- [Calculating Group Abundances and Ratios](#)
- [Calculating Standard Errors](#)
- [Calculating Protein Group Ratios](#)

The Peptide and Protein Quantifier node also calculates areas from the data of the Precursor Ions Quantifier node and the Precursor Ions Area Detector node (see [“Calculating Areas”](#) on page 447).

Finally, the node calculates the emPAI values for all processing workflows (see [“Calculating emPAI Values”](#) on page 448).

Calculating PSM Abundances

As a first step in calculating the quantification ratios, the Peptide and Protein Quantifier node determines the minimum detected quantification value for each spectrum file. Then it assigns the raw quantification values to the PSMs.

A quantification value is the intensity, S/N, or area detected for a given quantification channel. For reporter ion quantification, a quantification channel is one of the mass or reporter tags, and for precursor ion quantification, it is one of the different possible labeling states of a PSM corresponding to the different heavy amino acids used in the cell cultures. “Intensity” refers to both the intensity of the reporter peaks in reporter ion quantification and to the areas detected in precursor ion quantification. The S/N corresponds to the intensity of the peak divided by the noise values calculated for each peak during acquisition.

You can set the following parameters to control this step:

- **Replace Missing Values with Minimum Value:** Determines whether the application replaces missing values with the minimum detected quantification value in the spectrum file.
- **Reject Quan Results with Missing Channels:** Determines whether the application rejects all quantification values of the PSM when a value is missing but not replaced with the minimum value.
- **Use Single-Peak Quan Channels:** Determines whether the application uses quantification channels identified with only a single peak for the precursor ion pattern for precursor quantification.
- **Average Reporter S/N Threshold:** Specifies the minimum average reporter signal-to-noise threshold value to use. The application ignores PSMs with quantification results below this threshold.

[Table 27](#) lists some of the different circumstances that can arise in calculating quantification ratios for PSMs from the selected quantification values. The parameters in the table belong to the Peptide and Protein Quantifier node.

[Table 27](#) does not include cases resulting from PSM uniqueness and protein grouping. It focuses on cases where some or all of the quantification channels are zero. In these cases, the application detects nothing on a channel because the spectrum does not contain one of the reporter peaks, the heavy or light isotope pattern is missing, a quantification value falls below a specified minimum threshold, or the calculated ratios are very high or very low.

[Table 27](#) lists the different possible cases exemplified by arbitrary values. The values in the tables have [counts] as units if the cases are presented for reporter ion quantification. For precursor ion quantification, 114 and 115 are replaced by Light and Heavy, and the quantification values have [counts × seconds] as units.

Table 27. Ratios calculated from grouped abundances

Case	Minimum detected quan. value	Replace Missing Quan Values with Minimum Value parameter	Reject Quan Results with Missing Quan Channels parameter	Detected quantification values			Displayed or used quantification values		
				114	115	116	114	115	116
All quan. values detected	33	Irrelevant	Irrelevant	100	50	300	100	50	300
Quan. value missing for a quan. channel	33	No	No	100	0	300	100	0	300
	33	No	Yes	100	0	300	0	0	0
	33	Yes	Irrelevant	100	0	300	100	33	300
Quan. value missing for all quan. channels	33	Irrelevant	Irrelevant	0	0	0	0	0	0

Using Reporter Ion Isotopic Distribution Values to Correct for Impurities

iTRAQ and TMT kits consist of labels that contain different numbers of ^{13}C atoms, ^{15}N atoms, or both. For simplicity, assume that a 4plex kit yields peaks at 114, 115, 116, and 117 m/z , which correspond to $^{13}\text{C}1$, $^{13}\text{C}2$, $^{13}\text{C}3$, and $^{13}\text{C}4$, respectively. Because the label substances are not 100 percent isotopically pure, each label contains a certain number of other atoms. For example, the 116 label would not consist only of label molecules having three ^{13}C atoms but might also contain label molecules with only one or two ^{13}C atoms or even four or five ^{13}C atoms. As a result, these impurities lead to an observed peak at 116 m/z , which is smaller than might be expected if the tag were 100 percent isotopically pure, and additional peaks at positions -2 , -1 , $+1$, $+2$ Da apart from 116 m/z . The intensities of the latter peaks are proportional to the amount of the described isotopic impurities. When the 116 label and the 114, 115, and 117 labels are used, these latter three labels contribute to the peak at 116 m/z because of their isotopic impurities.

The intensity of the peak at 116 m/z effectively includes the following contributions:

$$(\text{observed intensity } 116) = (\text{true intensity } 116) - (\text{intensity loss because of } 116 \text{ impurities}) + (\text{intensity gain because of other label impurities})$$

To obtain the true intensity value of the 116 label—that is, the amount of the substance initially labeled with the 116 tag—you must correct the experimentally observed peak for the impurity of the labels.

For a 4plex sample, four formulas use the above equation for each of the labels, and the proper correction would consider both contributions in the formula by solving the system of coupled linear equations:

*(intensity_of_loss_because_of_116_impurity and
intensity_of_gain_because_of_other_label_impurities)*

For this correction, you must enter the isotopic distribution of each of the labels used in the quantification method, as described in “[Excluding PSMs with High Levels of Coisolation.](#)” The values are part of each of the iTRAQ or TMT label kits that were used.

You can also deconvolve the overlapping labels using other methods. Compatible with the Mascot search engine, the application uses a first-order approximation to the solution. The error made is small when the intensities of all possible contributing labels are of similar height, and it becomes larger if the intensity differences become larger. You can find more information on this topic, at www.matrixscience.com/help/quant_config_help.html. In the Sections menu, choose **Component > Correction > Type AB Certificate**.

In TMT kits, the values of impurities in reporter ion labels are normalized to the highest peak. In the application, they are normalized to 100 percent, which the application considers to be the sum of all impurities.

For TMTe 6plex kits, using correction factors is optional, because the resolution in Q-Exactive or Orbitrap Elite instruments is enough to separate ¹³C and ¹⁵N isotopes.

TMT 10plex kits use no correction factors by default. You can add correction factors to these kits by following the instructions in “[Excluding PSMs with High Levels of Coisolation.](#)”

In the kits, correction factors are the same for the same batch of isotopes and are usually valid for three to five years. See “[Correcting Reporter Ion Quantification Results for Isotopic Impurities](#)” on [page 419](#).

Excluding PSMs with High Levels of Coisolation

To create a fragment spectrum, you select a precursor mass for isolation, isolate and fragment the ions within a mass window that you define, and record the product ion masses created.

Ideally, you would isolate and fragment only the precursor ions of a single selected component. However, in practice you isolate the precursor ions within a user-specified window—typically 1 or 2 daltons around the isolation mass. Coeluting components with a mass falling into this isolation window are also isolated and fragmented. This process is called coisolation. The coisolating components are likely to be peptides whose fragments are observed in the created fragment spectra. The coisolation can prevent the identification of the selected peptide or lower the identification confidence.

Coisolation is also an issue for reporter ion quantification. In this type of quantification, the peptides from different samples—for example, different treatment states—are modified with isobaric labels. The isobaric labels fragment during precursor ion fragmentation and create reporter tags that appear in the low-mass region of the fragment spectra. The intensity ratio of the observed fragment tags is used for relative quantification of the peptides from the different sample charges.

The coisolating peptides also create reporter tags with the same masses as those from the selected peptide. The intensities of the reporter ions are therefore the sum of the intensities of the reporter ions for all coisolated peptides rather than the target peptide. As a result, the intensities are perturbed and are not accurate representations of the true abundance of the selected peptide. Furthermore, the perturbed ratios of the selected peptides that are greatly affected by coisolation can also adversely affect the ratios that the application calculates for the proteins that include these peptides.

Determining the extent to which the real reporter tag ratios of the selected peptides are perturbed depends on the level of coisolation and the isolation characteristics of the instrument. The application calculates and displays the percentage of interference within the precursor isolation window. This percentage is the relative amount of ion current within the isolation window that is not attributed to the precursor itself:

$$\%_isolation_interference = 100 \times \left[1 - \left(\frac{precursor_intensity_in_isolation_window}{total_intensity_in_isolation_window} \right) \right]$$

The application displays the calculated interference value in the % Isolation Interference column on the PSM and MS/MS Spectrum Info pages. For reporter ion quantification, a high isolation interference value could indicate that a calculated peptide ratio is skewed by the presence of coisolated peptide species.

Note The application only calculates the % Isolation Interference value if the precursor scans are high-resolution, high-mass-accuracy scans.

You can use the Co-Isolation Threshold parameter of the Peptide and Protein Quantification parameter to specify a threshold of between 0 and 100 percent for the allowed coisolation interference. The default value is 100 percent, which means that no PSM is excluded. This parameter is only used for reporter ion quantification.

Classifying PSMs


Seven different scenarios can occur when you derive protein quantification ratios from peptide quantification ratios. These cases show how the validity of using a given quantification result for the quantification of a certain protein depends on whether this particular quantification result is unique or shared among other peptides.

The peptide quantification ratios are taken from the associated quantification results. The term *quantification result* in this topic refers to MS/MS reporter intensities taken from the same scan as the identification (for example, ID-CID) or from a separate quantification scan (for example, Quan-HCD). The term also refers to intensities derived from the precursor scans in precursor ion quantification. A quantification result here is a general quantity associated with one or more peptides that are, in turn, associated with one or more proteins.

Case 1: Quantification Result Associated with One Spectrum, One Peptide, and One Protein

Case 1, shown in [Figure 290](#), is the simplest case. The quantification result is associated with one identification spectrum—whether the quantification results come from the same identification spectrum, from a different quantification spectrum, or from the precursor ion—and one peptide that is contained in one protein. The quantification result is unique for this protein. The application can mark peptide A “Unique” in the Quan Info column of the PSMs page if the quantification result meets other criteria.


Figure 290. Case 1: Quantification result associated with one identification spectrum, one peptide, and one protein



Case 2: Two Quantification Results Associated with Two Spectra, One Peptide, and One Protein

Case 2, shown in [Figure 291](#), is a variant of case 1. Each of two different quantification results is associated with a different identification spectrum. Both identification spectra identify peptide A, which is a peptide with the same sequence. Peptide A is only contained in one protein. Each of the two different quantification results is unique for just one protein. The peptides are redundantly identified and quantified, and you can use both for the quantification of protein A.

Figure 291. Case 2: Two different quantification results associated with two identification spectra, one peptide, and one protein




Case 3: Quantification Result Associated with Two Spectra, Two Peptides, and One Protein

Case 3, shown in [Figure 292](#), is similar to case 2 but varies from it in a slight but important way. In case 3, the two identification spectra are associated with the same quantification result rather than with two different quantification results. For example, you might obtain these results if you trigger the same precursor two times for MS/MS. It does not matter whether peptide A and peptide B are the same peptides (redundantly identified) or different peptides that are accidentally contained in the same protein. It also does not matter whether they are identified by the same search engine or by two different search engines, for example, a CID spectrum and an ETD spectrum. The quantification result is still unique for just one protein.

However, you cannot use the quantification ratio of both peptides A and B to calculate the quantification ratio of protein A, because it is the same quantification result, and you do not want to use the same quantification result multiple times for the same protein. In this case, the application marks peptide A—the peptide with the better identification—as “Unique” and the other peptide as “Redundant” for quantification (rather than redundant for identification).


Figure 292. Case 3: Quantification result associated with two identification spectra, two peptides, and one protein



Case 4: Quantification Result Associated with One Spectrum, Two Peptides, and One Protein

In case 4, shown in [Figure 293](#), the two peptides can be identified by the same search engine and have different ranks, or they can be identified by different search engines and both have rank 1. It does not matter whether peptide A and B have the same sequence with different PTM states or different sequences. The quantification result is unique for protein A. You can use it to calculate the protein ratio, but you must only use it once. The application marks the “better” peptide as “Unique” and the other as “Redundant” for quantification.


Figure 293. Case 4: Quantification result associated with one identification spectrum, two peptides, and one protein



Case 5: Quantification Result Associated with One Spectrum, One Peptide, Two Proteins

In case 5, shown in [Figure 294](#), the quantification result is associated with one identification spectrum and one peptide, but this peptide is contained in more than one protein. The quantification result is potentially shared between these proteins, and you do not know how to share it. If the quantification method specifies using only unique peptides for protein quantification, you would not use peptide A in this case. If the quantification method specifies using all peptides for protein quantification, the quantification result of peptide A would be divided equally between both proteins.


Figure 294. Case 5: Quantification result associated with one identification spectrum, one peptide, and two proteins



Case 6: Quantification Result Associated with One Spectrum, Two Peptides, and Two Proteins

In case 6, shown in [Figure 295](#), the quantification result is associated with one identification spectrum from which two different peptides are identified either by the same search engine as different ranks or by different search engines. The two different peptides are contained in two different proteins. The two different peptides are both unique to just one protein. Nevertheless, the associated quantification result is the same, and you do not want to use it for the calculation of the protein ratios if you specified in the quantification method to use only unique peptides. Only if you specify using all peptides can you use them for protein quantification. This case illustrates the discrepancy between the uniqueness of peptides and the uniqueness of the quantification results.


Figure 295. Case 6: Quantification result associated with one identification spectrum and two peptides unique to one protein



Case 7: Quantification Result Associated with Two Spectra, Two Peptides, and Two Proteins

Case 7, shown in [Figure 296](#), is a variant of case 6. Either the same search engine or different search engines identify different identification spectra, for example, CID and ETD spectra. As in case 6, the peptides are unique, but the quantification result is not. The result depends on whether you specified in the quantification method to use only unique peptides or all peptides.

Figure 296. Case 7: Quantification result associated with two identification spectra and two different peptides unique to one protein



Classification Flow Chart

The application calculates the quantification values for the protein as the sum of the intensity or S/N values for all PSMs belonging to the proteins that are marked as being usable. Whether the application considers a PSM usable is determined by the parameters of the Peptide and Protein Quantifier node, including those labeled Peptides to Use and Consider Protein Groups for Peptide Uniqueness. The Use Only Unique Peptides parameter includes the quantification peptides that do not occur in other proteins. The Proteins Groups for Peptide Uniqueness parameter defines peptide uniqueness on the basis of protein groups rather than individual proteins.

10 Performing Quantification

Calculating Quantification Results

When it determines peptide uniqueness for classification in the PSM Ambiguity column on the PSMs page, the application only considers the PSMs that it considered when creating the protein groups, if you select the Use Only Unique Peptides parameter. For example, it does not use for quantification a PSM of low confidence that it did not use to create the protein groups.

The flowchart in [Figure 297](#) shows how the application classifies PSMs for protein quantification. It displays this classification in the Quan Info column of the results report.

Figure 297. Classifying peptides for protein quantification


Depending on the settings of the parameters of the Peptide and Protein quantifier node, the application excludes the following peptides:

- No Quan Labels: No reporter label.
- Indistinguishable Channels: Not all defined channels can be distinguished.
- Inconsistently Labeled: Labels are from different channels.
- Excluded by Method: Quan channels are missing, single-peak channels are missing, ratios exceed limits, and so forth.

The application does not consider high-scoring peptides if they have no protein links.

The application marks peptides that are not unique as Redundant.

This check depends on the setting of the Consider Protein Groups for Peptide Uniqueness parameter of the Peptide and Protein Quantifier node. It classifies the peptides as Unique if they match the proteins within the same protein group.



Calculating Peptide Group Abundances

After classifying PSMs, the Peptide and Protein Quantifier node calculates peptide group abundances for the different samples from their associated PSMs. It calculates the abundance for a peptide group as a simple summation of its associated and used PSM abundances.

When the node calculates peptide group abundances, it creates the following columns:

- Peptide Quan Usage column on the PSMs page, which shows the PSMs that are used for calculating peptide group abundances. This column displays all PSMs, whether or not they are unique.
- Abundances column on the Peptide Groups page, which displays the peptide group abundances.
- Abundance Counts column on the Peptide Groups page, which displays the number of PSM abundances used for calculating the peptide group abundances.

Classifying Peptide Groups

The Peptide and Protein Quantifier node classifies which peptide groups are used for protein quantification.

Use the following parameters to control the classification of peptide groups:

- The Consider Protein Groups for Peptide Uniqueness parameter determines which peptide groups are unique. If you set it to True, the Peptide and Protein Quantifier node considers a peptide unique if it is included in only one protein group. If you set it to False, the node considers a peptide unique if it is included in only one protein.
- The Peptides to Use parameter determines which peptides to use.
 - Unique: Uses only unique peptide groups.
 - Unique + Razor: Uses both unique peptide groups and peptide groups containing razor peptides for the best associated master protein. Razor peptides are shared among multiple protein groups or proteins. When you set the Consider Protein Groups for Peptide Uniqueness parameter to True, the node uses peptide groups that are not unique for all proteins belonging to the protein group of the best master protein. The best master protein is the protein with the largest value in the # Protein Unique Peptides column on the Proteins page and with the smallest value in the Coverage column (the longest protein).
 - All: Uses all peptide groups, whether they are unique or not.

When it classifies peptide groups, the Peptide and Protein Quantifier node creates the following columns in the .pdResult file:

- Quan Info column on the Peptide Groups page, which displays No Quan Values for peptide groups that have no abundance value. It displays Not Unique or Unique for all other peptide groups.

- Protein Quan Usage column on the Peptide Groups page, which shows which peptide group abundances are used for protein quantification. This column is only visible when you look at the peptide groups in the associated tables.
- # Razor Peptides column on the Proteins page, which displays the number of razor peptides for proteins when you use razor peptides for quantification.

Calculating Protein Abundances

After the Peptide and Protein Quantifier node classifies peptide groups, it calculates protein abundances for the samples from the associated peptide group abundances. It calculates a protein abundance as the simple summation of its associated and used peptide group abundances.

Precursor ion area detection uses a different approach. It presents the average peak areas of the top *N* unique peptides for each protein, usually set to 3 by default.

The node creates the Abundance Counts column on the Proteins page, which shows the number of peptide group abundances used for calculating abundances.

Normalizing Peptide Groups and Protein Abundances

After the Peptide and Protein Quantifier node calculates the protein abundances for the samples, it normalizes the peptide groups and protein abundances, scales them, or both. In general, the node applies the normalization of the total abundance values for each channel across all files, equalizing the total abundance between different runs.

You can use the following parameters of the Peptide and Protein Quantifier node to normalize peptide groups and protein abundances.

Normalization Mode Parameter

Use the settings of the node's Normalization Mode parameter to specify how the node should perform normalization:

- (Default) None: Does not perform normalization.
- Total Peptide Amount: Sums the peptide group abundances for each sample and determines the maximum sum for all files. The normalization factor is the factor of the sum of the sample and the maximum sum in all files.
- Specific Protein Amount: Calculates the normalization factor from the abundances of selected proteins in the specified FASTA file. You can specify a FASTA file (use the Proteins for Normalization parameter to specify the name of the FASTA file) that can contain one or multiple proteins. The Peptide and Protein Quantifier node uses all proteins in the FASTA file that are contained in the result file and that have any protein abundance. It calculates the maximum sum for all files. The normalization factor is the factor of the sum of a sample and the maximum sum in all files.

If you set the Normalization parameter to Specific Protein Amount but do not select a FASTA file, or if the .pdResult file contains no proteins that appear in the FASTA file, the application cannot perform normalization.

When normalization is based on specific proteins, the node sums the abundance values of all proteins in the .pdResult file that are contained in the selected FASTA file. It compares the protein sequences in the FASTA file to the sequences in the .pdResult file rather than to the accession strings. You can therefore use the same FASTA file when the title line is slightly different.

After calculating the normalization factors, the Peptide and Protein Quantifier node normalizes peptide group abundances and protein abundances by dividing abundances with the normalization factor over all samples.

Scaling Mode Parameter

After normalization, the Peptide and Protein Quantifier node scales the normalized abundances if you set the Scaling Mode parameter to On Channels Average (Per File) or On Control Channels Avg. (Per File). By default, the node does not scale abundances.

- **On Channels Average (Per File):** After it aggregates all the abundance or normalized abundance values per quantification channel, the node scales the abundance values of each channel so that the average of all channels is 100. It scales per file.
- **On Control Channels Avg. (Per File):** After it aggregates all the abundance or normalized abundance values per quantification channel, the node scales the abundance values of each channel so that the average of all control channels in a file is 100 and then scales all other channels up or down relative to 100. It scales per file.

The application displays normalized abundances in the Abundances (Normalized) column and scaled abundances in the Abundances (Scaled) column. The Abundances (Normalized) column appears on the Proteins and Peptide Groups pages of the .pdResult report when you set the Normalization Mode parameter of the Peptide and Protein Quantifier node to Total Peptide Amount or Specific Protein Amount. The Abundances (Scaled) column appears on the same pages when you set the Scaling Mode parameter of the Peptide and Protein Quantifier node to On Channels Average (Per File) or On Control Channels Average (Per File). These columns do not appear if you do not apply normalization and scaling.

Sample Information Used to Calculate and Display Quantification Results

The application displays the Abundance (Scaled) column and the Abundance (Grouped) column in blue, white, and red tones, as shown in [Figure 298](#). For information about how the application calculates the values in these columns, see [“Calculating Group Abundances and Ratios”](#) on [page 446](#).

Figure 298. Colors used to display the Abundance (Scaled) and Abundance (Grouped) columns

Fold-Change Color Threshold	Scaled Abundance Color Threshold	Color
0.10	18.2	Dark Blue
0.13	22.2	Blue
0.17	28.6	Light Blue
0.25	40.0	Very Light Blue
0.50	66.7	White
0.00	100.0	White
2.00	133.3	Light Red
4.00	160.0	Red
6.00	171.4	Dark Red
8.00	177.8	Very Dark Red
10.00	181.8	Dark Red

The application uses the settings of the 1st Fold Change Threshold, 2nd Fold Change Threshold, 3rd Fold Change Threshold, 4th Fold Change Threshold, and 5th Fold Change Threshold parameters of the Peptide and Protein Quantifier node to color the scaled abundance values and the ratios calculated from the abundance values shown in the Abundance (Scaled) column and the Abundance (Grouped) columns. For the ratios, the application directly uses the fold changes as thresholds for applying colors from dark blue (downregulation) to white (no regulation) to dark red (upregulation). For the scale abundance values, the application colors the values according to their deviation from 100. It calculates the thresholds of the scaled and scaled-and-grouped abundances as follows:


$$abundance\ coloring\ threshold = 100 \frac{fold\ change\ threshold - 1}{fold\ change\ threshold + 1}$$

T is the given fold change, and 100.0 is used to scale the average of all channels to 100.0 (per file or across files).

Figure 298 shows the thresholds of the scaled abundance values for the default fold changes.

Figure 299 shows the fold-change coloring scheme for the grouped abundance and scaled abundance values.

Figure 299. Abundances (Grouped) and Abundances (Scaled) columns on the Proteins page



Calculating Group Abundances and Ratios

After normalization, the Peptide and Protein Quantifier node calculates the grouped abundances for peptide groups and proteins. It uses the following preprocessed abundances as a basis, if they are available:

- Ungrouped raw abundances, where no normalization and no scaling are applied
- Ungrouped normalized abundances
- Ungrouped scaled abundances (optionally also normalized)

The node calculates these as the average of the preprocessed sample abundances that belong to the same sample group. It calculates group abundances only when there is at least one sample group consisting of multiple samples.

Finally, the node calculates ratios by dividing the abundance values of the associated sample groups. It can calculate ratios only when there are group abundances.

The Maximum Allowed Fold Change parameter defines the maximum fold change that the node calculates. All ratios are between 1/fold change and fold change. When the numerator or the denominator is 0, the ratio is the fold change when the denominator is 0 or 1/fold change when the numerator is 0.

Table 27 shows how the application calculates ratios from grouped abundances.

Table 28. PSM ratios calculated when quantification values are missing

Case	Maximum Allowed Fold Change parameter	Grouped abundances		Calculated ratio	Displayed ratio
		114	115	115/114	115/114
Ratio is within the limits.	100	1000	500	0.5	0.5
	100	1000	2000	2	2
Ratio exceeds the limits.	100	1000	1	0.001	0.01
	100	1	1000	1000	100
Ratio is outside limit because one quan. channel is missing.	100	1000	0	0	0.01
	100	0	1000	∞	100
Ratio is outside the limit because both abundance values are missing.	100	0	0	0	-

Calculating Standard Errors

The Peptide and Protein Quantifier node calculates the standard error of the grouped abundances as follows:

$$standard_error = 100 \times \frac{1.483 \times MAD(a, \dots, a)}{b \times \sqrt{b}}$$

Calculating Protein Group Ratios

The Peptide and Protein Quantifier node does not display quantification information for Protein Groups.

Calculating Areas

The Peptide and Protein Quantifier node calculates areas for peptide groups and proteins when the application calculates PSM areas in the processing workflow. It calculates areas for each sample group that you defined on the Grouping & Quantification page of the analysis.

The node calculates areas for peptide groups and proteins only when the processing workflow includes the Precursor Ions Area Detector node. It does not calculate areas in this case, because the abundances of peptide groups and proteins are equivalent to areas. The Peptide and Protein Quantifier node calculates peptide group areas for each sample group. The area of a peptide group is the maximum PSM area that is associated with this peptide group and that belongs to the sample group.

The node also calculates protein areas for each sample group by calculating the area for each sample of the sample group as the average of a specified number of the most abundant distinct peptides identified for the protein and the sample. The area for the sample group is the average of the areas of each sample.

You can specify the number of peptides used in the protein area calculation with the Top *N* Peptides Used for Area Calculation parameter of the Peptide and Protein Quantifier node in the consensus workflow. The peptides must have different sequences to be considered distinct. The application considers peptides with different charge states and modification variants of the same sequence to be the same peptide.

You can direct the application to use only unique peptides in protein area calculation. Unique peptides are specific to a single protein or a protein group, so they are not affected by the abundance of other proteins in the sample. Without the effects of other proteins, the quantification is more precise. However, the application might quantify fewer proteins, because the list of peptide candidates for quantification is shorter when restricted to unique peptides.

You can specify the usage of only unique peptides for area calculation with the Peptides to Use and the Consider Protein Groups for Peptide Uniqueness parameters of the Peptide and Protein Quantifier node in the consensus workflow. If you set the Peptides to Use parameter to Unique or Unique + Razor, the application uses only unique peptides.

The Consider Protein Groups for Peptide Uniqueness parameter of the Peptide and Protein Quantifier node determines if the peptide must only be unique for a protein group. If you set this parameter to True, the application considers a peptide unique if it is included in only one protein group. This setting therefore is less restrictive. If you set this parameter to False, the application considers a peptide unique if it is included in only one protein. This setting is therefore more restrictive.

Calculating emPAI Values

The exponentially modified protein abundance index (emPAI) is a simple measurement of protein abundance that is based on the number of found peptides³. It correlates to the absolute amount of protein in a sample.

The PAI is calculated as follows:

³ Ishihama, Y. et al. *Molecular & Cellular Proteomics*. **2005**, *4*, 1265.

$$\text{PAI} = \frac{N_{obsd}}{N_{obsbl}}$$

where:

- N_{obsd} is the number of observed, or found, peptides.
- N_{obsbl} is the number of observable peptides.

The emPAI itself is an exponential transformation of the PAI:

$$\text{emPAI} = 10^{\text{PAI}} - 1$$

The application calculates the emPAI only for tryptic digests.

See the emPAI column on the Proteins page for the emPAI values.

Filtering Quantification Data with Signal-to-Noise Values

For Orbitrap data, you can filter reporter ion quantification spectra that have too much variability in the intensities of the reporter ions by using signal-to-noise values instead of intensities.

Using Signal-to-Noise Values as Quantification Channel Values

Because of ion statistical effects, measurements of intensities are more variable and therefore less accurate when the mass peaks arise from a smaller population of ions. When you perform precision quantification, you might want to filter out reporter peaks or spectra that have a variability that is too high because of ion statistical effects. However, with an Orbitrap analyzer, the intensity of a signal is not proportional to the number of ions corresponding to the signal. Instead, the number of ions is proportional to the S/N value of a peak ($\text{number_of_ions} = 6 \times \text{S/N}$ for D30 Orbitraps, and $\text{number_of_ions} = 4 \times \text{S/N}$ for D 20 Orbitraps).

To filter reporter ion quantification input data that is detected from a number of ions that is too small, you use the signal-to-noise (S/N) values of the reporter ion peaks instead of their intensities to provide such a measurement when ions are acquired in the Orbitrap. You can use the settings of the Reporter Abundance Based On parameter of the Peptide and Protein Quantifier node to select signal-to-noise intensity values. For information on this parameter, see the Help.

Filtering Quantification Data with Average Reporter Ion Signal-to-Noise Values

For Orbitrap data, you can obtain more reliable TMT quantification results by filtering out quantification spectra with too few ions for the reporter ion signals. The Peptide and Protein Quantifier node includes a parameter, Average Reporter S/N Threshold, that you can use to specify an average reporter S/N threshold value that determines which PSMs the application excludes from quantification. The application filters quantification values by excluding PSMs that have an average reporter ion S/N value smaller than this threshold value.

In addition, the Average Reporter S/N threshold parameter adds the Average Reporter S/N column, which displays average reporter ion S/N values on the PSMs and Quan Spectra pages. The application calculates the values as the sum of S/N values found, divided by the number of defined tags.

In the following example, the application calculates the average S/N values for a .pdResult file created from an MSF file containing TMT 6plex data. The quantification method has six channels, and the application calculates the average S/N value by $\text{sum}(S/N)$ divided by 6. [Table 29](#) shows the results.

Table 29. Average S/N values for a single MSF file containing TMT 6plex data

Peptide	126	127	128	129	130	131	Average S/N
1	10	10	12	8	7	3	$50/6=8.333$
2	10		8		2		$20/6=3.333$
3							$0/6=0$

In the next example, the application calculates the average S/N values for a .pdResult file created from two MSF files containing TMT 6plex data. The quantification method for both files has six channels, and the application calculates the average S/N by $\text{sum}(S/N)$ divided by 6. [Table 30](#) shows the results.

Table 30. Average S/N values for two combined MSF files containing TMT 6plex data


Peptide	126	127	128	129	130	131	Average S/N
A1	10	10	12	8	7	3	$50/6=8.333$
A2	10		8		2		$20/6=3.333$
B1	20	21	17	22		10	$80/6=13.333$
B2							$0/6=0$

When you combine MSF files containing reporter ion quantification data with MSF files lacking reporter ion quantification data, the PSMs without reporter ion quantification data display no values in the Average Reporter S/N column on the PSMs and Quan Spectra pages.

Identifying Isotope Patterns in Precursor Ion Quantification

The quantification spectra on the pages of the .pdResult file show the isotope pattern used for quantifying the peptides. The algorithm used in precursor ion quantification finds isotope patterns by identifying target components—that is, known elemental compositions from event lists. It identifies the peptides and searches in the event lists for the isotope patterns of these identified peptides. After peptide identification, the algorithm follows the steps shown in [Figure 300](#) to identify the isotope patterns.

Figure 300. Identifying isotope patterns



Viewing Quantification Results

You can view the results of quantification by using the Quantification page of the Result Summaries, data distribution quantification maps, the Quantification Channel Values chart, the Quantification Spectrum chart, the Quan Spectra chart, and the Quantification Ratio chart. For details, see the following topics:

- [Using the Quantification Page of the Result Summaries](#)
- [Displaying Data Distribution Maps for Quantification](#)
- [Displaying the Quantification Channel Values Chart](#)
- [Displaying the Quantification Spectrum Chart](#)
- [Displaying the Quan Spectra Page](#)
- [Displaying Quantification Ratio Charts](#)
- [Treatment of Missing Reporter Ions for Quantification](#)

Using the Quantification Page of the Result Summaries

For information on the Quantification page of the Result Summaries pane, refer to the Help.

Displaying Data Distribution Maps for Quantification

In workflows that include quantification nodes, you can use the Peptide and Protein Quantifier node in the consensus workflow to display ratio, standard error, abundances, and count columns as data distribution maps in .pdResult reports. These maps show the distribution of values across sample groups and ratios. They can help you visually validate the results of the quantification.

You can also use the Data Distributions node in the consensus workflow to display the Found in Samples, Found in Sample Groups, Found in Files, and Found in Fractions columns in .pdResult reports as data distribution maps. These maps show the distribution of values across the available files, samples, and sample groups.

The application combines all values for each generated group of quantification values (ratios, standard errors, abundances, and counts) into one data distribution map. You can see all grouped ratios in relation to the ratios of the individual replicates and the associated variabilities.

For detailed information about the Peptide and Protein Quantifier node, refer to the Help.


Sample Information Used to Calculate and Display Quantification Results

When the application processes the data, it first extracts the raw quantification values during the processing step. For reporter ion quantification, these values are the intensities of the reporter peaks for each MSⁿ spectrum. For precursor ion quantification, these values are the integrated peak areas of the light, medium, or heavy forms of the peptides. In the consensus workflow, the Peptide and Protein Quantifier node refers to the information generated during the analysis setup to determine which ratios to calculate for each file and which of these ratios to group by calculating an average ratio from them.

The ratio changes are currently color-coded. The area values are also color-coded, with the color ranging from yellow to red with increasing area.

Figure 301 shows the data distribution maps of the found ratio changes for the specified sample groups.


Figure 301. Distribution maps with color-coded ratio change for the specified sample groups



You can use the *n*th Fold Change Threshold parameters of the Peptide and Protein Quantifier node to adjust the thresholds of the color highlighting in the Abundance and Abundance Ratios (log₂) column of the Proteins page. For more information on these parameters and the color highlighting, refer to the Help.

The display filters provide a means to filter data on the basis of the data distribution maps (for more information on display filters, see “Filtering with Display Filters” on page 254). The example in Figure 302 filters for the proteins that have a fold change of at least four in at least five of the ratio groups present.


Figure 302. Display filters applied to the distribution map columns



Sample Information Used to Display Identifications and Quantifications Among Files and Samples

The application uses the information about the grouping of samples in the data distribution maps to show which protein group, protein, or peptide was identified in which file, sample, and sample group, as shown in Figure 303. The color-coding refers to the confidence level of the identification—that is, the confidence of the best PSM found for the particular file, sample, or sample group.

Figure 303. Data distribution maps with color-coded identification confidence levels for proteins found in sample groups, samples, and files





In precursor ion identification, there is frequently no confidence information from an identified MS n spectrum, but the application found a quantification value for a particular sample. For example, it identified only the light peptide in a SILAC experiment, but it also found an integrated extracted ion current (XIC) chromatogram for the heavy peptide. The application displays such cases in blue on the distribution map.

Generating Data Distribution Maps


❖ To generate a distribution map

1. Create or open a study and an analysis:
 - To create a study, see “[Creating a Study](#)” on [page 38](#).
 - To open an existing study, see “[Opening an Existing Study](#)” on [page 40](#).
 - To create an analysis, see “[Creating an Analysis](#)” on [page 71](#).
 - To open an existing analysis, see “[Opening an Existing Analysis](#)” on [page 72](#).
2. Create a quantification processing workflow. For information on this procedure, see “[Creating a Processing Workflow for Precursor Ion Quantification](#)” on [page 390](#) or “[Creating a Processing Workflow for Reporter Ion Quantification](#)” on [page 397](#).
3. Create a quantification consensus workflow. For information on this procedure, see “[Creating a Consensus Workflow for Precursor Ion Quantification](#)” on [page 394](#) or “[Creating a Consensus Workflow for Reporter Ion Quantification](#)” on [page 400](#).
4. Set the appropriate parameters of the Peptide and Protein Quantifier node to display the appropriate columns in the output .pdResult file.

For information about these parameters, refer to the Help.
5. (Optional) If you want the .pdResult report to display the Found in Samples, Found in Sample Groups, Found in Files, and Found in Samples columns, drag the **Data Distributions** node to the to the Post-Processing area of the Workflow Editor.
6. (Optional) Save the workflow:
 - a. In the Name box above the Workflow Tree pane, type a name for the consensus workflow.
 - b. (Optional) In the Description box above the Workflow Tree pane, type a brief description of the consensus workflow.
 - c. In the Workflow Editor, click the **Save** icon,  **Save** , or the **Save Common** icon,  **Save Common** .
 - d. In the Save Workflow dialog box, do the following:
 - i. Browse to the file to save the template in, or type a file name in the File Name box.
 - ii. In the Save As Type box, select **Consensus Workflow File (*.pdConsensusWF)**.

iii. Click **Save**.

The application saves the workflow in the *file_name*.pdConsensusWF file.

7. Save the analysis. See “Saving an Analysis” on page 80.
8. Save the study. See “Saving a Study” on page 62.
9. Click the **Run** icon, , in the upper right corner of the Analysis window.

The job queue appears, as shown in Figure 79 on page 112, displaying the status of your search. Use the job queue to check the status of your search as the search progresses. For information about the job queue, refer to the Help.

Displaying Data Distribution Maps

❖ To display a distribution map

1. Open the .pdResult file and use the Field Chooser to make the appropriate columns visible in the .pdResult file.
2. Click the **Protein Groups**, **Proteins**, or **Peptide Groups** tab, as appropriate.
3. In the Field Chooser, select the check boxes of the appropriate data distribution columns. Examples of these columns include the following:
 - Found in Samples
 - Found in Files
 - Found in Sample Groups
 - Found in Fractions

These columns display the best confidence of the PSMs of the protein that the application identified in the files.

The colors of the boxes in the Found in Samples, Found in Sample Groups, Found in Files, and Found in Fractions columns represent the following:

- Green: High confidence
- Yellow: Medium confidence
- Red: Low confidence
- Blue: Found but unidentified PSM. Only the results from precursor ion quantification searches contain blue boxes.

Each column pertaining to data distribution features an Expand icon, , to the right of the column name, as shown in Figure 304.

10 Performing Quantification

Viewing Quantification Results


Figure 304. Data Distribution columns marked by Expand icons

Checked	Master	Accession	Description	Found in Samples	Found in Sample Groups	Areas	Biological Process	Cellular Component	Molecular Function	Ratios	Ratio E	Ratio V	Ratio C
		P20336	Ras-related protein Ra			6.1e6							
		O75223-2	Isoform 2 of Gamma-g			7.7e6 (3.2e7)				4.54	11.29	83.2	2
		P19367-2	Isoform 2 of Hexokinase			1.3e6 (2.5e6)				2.14	2.04	7.6	3
		P09972	Fructose-bisphosphate			6.2e6 (3.3e7)				4.98	1.97	19.8	4
		O00429-6	Isoform 6 of Dynamini-			4.2e7 (4.6e7)				4.03	0.75	5.2	3
		P07741-1	Adenine phosphoribos			6.1e6 (1.9e6)				2.68	0.69	4.1	5
		Q14240-2	Isoform 2 of Eukaryoti			1.9e6 (4.2e6)				2.25	11.13	36.5	2
		P62491-1	Ras-related protein Ra			1.9e6 (4.9e6)				2.44	0.91	15.5	48
		P10809	60 kDa heat shock pro			1.1e7 (4.9e7)				4.43			1
		O75914-3	Isoform 3 of Serine thr			1.0e7							
		Q08257-3	Isoform 3 of Quinine c			4.0e7 (1.0e8)				2.58			1
		P15154-1	Ras-related C3 botulin			3.9e6 (6.1e6)				1.83	2.13	13.0	11
		O43175	D-3-phosphoglycerate			4.9e6 (9.9e6)				2.20	0.66	7.2	24
		P06733-1	alpha-enolase [OS-Hc			2.9e6 (5.9e6)				2.38	1.89	13.5	8
		P22234	multifunctional protein			1.2e6 (1.3e7)				4.13			1
		Q9NQW7-1	xaα-Pro aminopeptid										
		O50331	Phosphatidylinositol 4-							1.02			1
		Q14232	Translation initiation fa			7.2e7 (7.3e7)				2.23	4.16	18.7	4
		P37837	Transaldolase [OS-Hc			2.4e6 (5.0e6)							
		Q15418-1	Ribosomal protein S6										
		P36873-1	serine/threonine-prote			3.7e7 (7.7e7)				3.86	18.75	136.1	2
		Q02218-2	Isoform 2 of 2-oxogluc			3.4e7 (6.7e7)				1.94	18.26	63.5	2
		P62491-2	Isoform 2 of Ras-relat			1.8e6 (4.2e6)				2.25	11.13	36.5	2
		P48507	Glutamate-cysteine li			2.3e7 (1.6e7)				0.76			1
		P35998-2	Isoform 2 of 26S prote			2.9e7 (3.6e7)				5.12	12.28	151.0	3
		Q964N2-2	Isoform 2 of Putative a							2.18	1.26	11.0	16
		P00558	phosphoglycerate kina			1.9e6 (2.3e6)							
		P08237-3	Isoform 3 of ATP-depe							2.58	1.89	18.4	14
		Q15046	Lysine-4RNA ligase [C			1.6e7 (5.1e7)				4.10	5.37	51.5	2
		Q15008-2	Isoform 2 of 26S prote			6.2e7 (2.9e8)				3.18			1
		P63208	α-phase kinase-associ			1.8e6 (4.3e6)				2.60	4.29	27.8	6
		P54819	Adenylate kinase 2, m										
		Q43865-2	Isoform 2 of Putative a			1.9e6 (2.3e6)				2.18	1.26	11.0	16
		P00558-2	Isoform 2 of Phosphog			9.9e7 (2.4e8)				3.07	1.93	14.6	6
		O95831-1	Apoptosis-inducing fac			6.9e7 (1.6e8)				1.94			1
		P67837	ribulose-bisphosphate										

- (Optional) In the appropriate columns pertaining to data distribution, click the Expand icon, , to open a legend that displays the identification number of the input file (shown in the ID column on the Input Files page) and the quantification label.

Figure 305 shows the legends for the Found in Samples, Found in Sample Groups, Areas, Biological Process, and Ratios columns of the .pdResult report.

Figure 305. Legends in data distribution columns



❖ **To sort the data in distribution map columns**

See the Help.

Displaying the Quantification Channel Values Chart

When you select a PSM or a quantification spectrum for reporter ion quantification, the Quantification Channel Values chart displays the absolute intensity of the quantification values detected for the available quantification channels. When you select a protein or a peptide group, the chart displays the abundance of the item.

When you select a PSM or a quantification spectrum for precursor ion quantification, the Quantification Channel Values chart displays the area of the quantification values detected for the available quantification channels. When you select a protein or a peptide group, the chart displays the abundance of the item.

The chart shows the normalized abundances, if available, or the raw abundances. It can also display grouped abundances, if available. Grouped abundances, which are the abundances of the sample groups defined on the Grouping and Quantification page of the study, appear as blue bars. Ungrouped abundances, which are the abundances of the samples defined on the Grouping and Quantification page of the study, are shown as pastel-colored bars. Ungrouped


abundances are ordered. When the results contain grouped and ungrouped abundances, you can use the shortcut menu in the chart to choose the data to display. The Show Legend command on the shortcut menu labels the ungrouped bars with the sample identifier. It also displays a legend for the colored bars.

The processing workflow must include a quantification node for this chart to appear in the results report.

❖ To display the Quantification Channel Values chart

1. Click the appropriate tab (**Proteins**, **Peptide Groups**, **PSMs**, or **Quan Spectra**) in the .pdResult file.
2. Select the row of the protein, peptide group, PSM, or quantification spectrum that interests you.

To obtain meaningful results when you display the chart on the PSMs page, “Used” must appear in the Quan Info column of the report.

3. Choose **View > Quan Channel Values**, or click the **Quan Channel Values** icon, .

To view the results, see the following topics:

- [Displaying Quantification Channel Values for Reporter Ion Quantification](#)
- [Displaying Quantification Channel Values for Precursor Ion Quantification](#)

Displaying Quantification Channel Values for Reporter Ion Quantification

The Quantification Channel Values chart for reporter ion quantification displays the absolute intensity of the reporter ions detected for the available quantification channels when you select a PSM or a quantification spectrum. When you select a protein or a peptide group, the chart displays the abundance of the item.


Displaying the Quantification Channel Values Chart for PSMs and Quantification Spectra

Reporter ions, also called reporters, are produced from MS/MS of the labels affixed to peptide samples in reporter ion quantification. You can use the quantification value intensity to calculate the relative ratio of a peptide. You might also want to view the absolute quantification value intensity to verify that the peptide ratio calculation is correct.

The *x* axis of the chart shows the names of the quantification channels, and the *y* axis shows the intensity of the reporter ions, in counts.

The TMTe 6plex quantification method has six reporter ions. Suppose that they are used to label six biological samples: 126, 127, 128, 129, 130, and 131. Figure 306 shows the Quantification Channel Values chart created by the Show Quan Channel Values command for these samples. It shows the relative intensities of the samples labeled with the 126, 127, 128, 129, 130, and 131 reporter ions. Clearly, the sample labeled 130 is the sample with the greatest reporter ion intensity.

Figure 306. Quantification Channel Values chart showing the absolute intensity of a PSM in reporter ion quantification




A legend is not available in the Quantification Channel Values chart when you display PSMs and quantification spectra.

Displaying the Quantification Channel Values Chart for Proteins and Peptide Groups

Figure 307 shows the grouped and ungrouped abundances of a selected protein. The x axis of the chart shows the quantification channels, and the y axis shows the abundance of the selected protein for the quantification channel, defined by arbitrary units. In this example, the quantification channels are labeled as study factors rather than as reporter ion masses. Error bars denoting the standard error are displayed above the blue bar that shows the average abundance of the ungrouped abundances. An optional legend at the bottom displays the grouped and ungrouped abundances.

Figure 307. Quantification Channel Values chart showing the abundance of a protein in reporter ion quantification



❖ **To display the legend**

Right-click the Quan Channel Values pane and choose **Show Legend**.

❖ **To display grouped abundances**

Right-click and choose **Show Grouped Abundances**.

❖ **To display ungrouped abundances**

Right-click and choose **Show Ungrouped Abundances**.

Displaying Quantification Channel Values for Precursor Ion Quantification


The Quantification Channel Values chart for precursor ion quantification displays the area of the isotopes detected for the available quantification channels when you select a PSM or a quantification spectrum. When you select a protein or a peptide group, the chart displays the abundance of that item.

Displaying the Quantification Channel Values Chart for PSMs and Quantification Spectra

Heavy isotopes are incorporated into proteins in precursor ion quantification. You can use the quantification value area to calculate the relative ratio of a peptide. You might also want to view the quantification value area to verify that the peptide ratio calculation is correct.

The x axis of the chart shows the quantification channels, and the y axis shows the detected area for the given quantification channel, defined by counts \times seconds.

The dimethylation 3plex quantification method in SILAC has a mixture of samples labeled with Light, Medium, and Heavy isotopic labels. [Figure 308](#) shows the chart created by the Show Quan Channel Values command for these samples. It shows the relative peak area for the samples with the Light, Medium, and Heavy isotope labels. The sample labeled Medium is the sample with the greatest area.


Figure 308. Quantification Channel Values chart for precursor ion quantification

A legend is not available in the Quantification Channel Values chart when you display PSMs and quantification spectra.

Displaying the Quantification Channel Values Chart for Proteins and Peptide Groups

Figure 309 shows the grouped and ungrouped abundances of a selected protein. The x axis of the chart shows the quantification channels, and the y axis shows the abundance of the selected protein for the quantification channel, defined by arbitrary units. In this example, the quantification channel names are replaced by the names of the study factors. An optional legend at the bottom displays the grouped and ungrouped abundances.

Figure 309. Quantification Channel Values chart showing the abundance of a protein in precursor ion quantification



❖ **To display the legend**

Right-click the Quan Channel Values pane and choose **Show Legend**.

❖ **To display grouped abundances**

Right-click and choose **Show Grouped Abundances**.

❖ **To display ungrouped abundances**

Right-click and choose **Show Ungrouped Abundances**.

Displaying the Quantification Spectrum Chart

You can generate a chart showing the spectrum used for quantification. This chart is available for every peptide with an associated quantification result. You must conduct a search with a workflow that includes a quantification node for this page to appear in the results report.

❖ To display the Quantification Spectrum chart

1. Select the PSM of interest.

The peptide must be labeled “Used” in the Quan Usage column of the PSMs page.

2. Choose **View > Show Quantification Spectrum**, or click the **Show Quantification Spectrum** icon, .


To see the results, see the following topics:

- [Displaying the Quantification Spectrum Chart for Precursor Ion Quantification](#)
- [“Displaying the Quantification Spectrum Chart for Reporter Ion Quantification” on page 471](#)

Displaying the Quantification Spectrum Chart for Precursor Ion Quantification

For precursor ion quantification, the Quantification Spectrum chart displays a quantification spectrum for each peptide. It also displays the different abundances of the identified Light, Medium, and Heavy isotopic peak patterns used to quantify a peptide. The abundances are measured by calculating the area of the extracted ion chromatogram of each isotope of a pattern. The chart highlights the corresponding isotope pattern peaks and labels them with the quantification channel names, as shown in [Figure 310](#). It also includes any peaks that are not part of an isotope pattern.

Figure 310. Quantification Spectrum chart for precursor ion quantification




The *x* axis of the chart displays the mass-to-charge ratio of the isotopes, and the *y* axis displays the area of the extracted ion chromatogram for the isotopes. Filled blue circles mark the isotope pattern peaks that were used for calculating the quantification values for the different quantification channels. Unfilled blue circles mark the isotope pattern peaks that were identified but not used. The Quantification Spectrum chart always compares the exact same isotopic pattern peaks for each label. For example, the chart in Figure 311 compares the first two isotopic pattern peaks among all three types: Light, Medium, and Heavy. But the chart also contains an additional Light isotopic pattern peak and an additional Medium isotopic pattern peak that are not used, so these two peaks are represented by unfilled circles.

10 Performing Quantification


Viewing Quantification Results

Figure 311. Extra isotopic pattern peaks represented by unfilled circles in the Quantification Spectrum chart



The Quantification Spectrum chart can also indicate whether an expected quantification pattern peak is absent. Regions in pink indicate where a quantification pattern peak was expected but is absent (see [Figure 312](#)). This ion pattern peak is not used in calculating the quantification values for the different quantification channels.

Figure 312. Expected but absent peak in the Quantification Spectrum chart



Regions in blue, shown in Figure 313, indicate where a quantification pattern peak was expected but is unsuitable. Pattern peaks might be unsuitable because of the wrong centroid retention time, out-of-delta mass range, wrong intensity, or a peak that has been used by another isotopic pattern. This ion pattern peak is not used in calculating the quantification values for the different quantification channels.

10 Performing Quantification

Viewing Quantification Results

Figure 313. Expected but unsuitable peaks in the Quantification Spectrum chart




Table 31 shows what the various colors mean on the Quantification Spectrum charts in Figure 312 and Figure 313.

Table 31. The meaning of colors in the Quantification Spectrum chart (Sheet 1 of 2)

Color	Meaning
Filled blue circle	Indicates the isotope pattern peaks that are used in calculating the quantification values for the different quantification channels.
Unfilled blue circle	Indicates the isotope pattern peaks that are not used in calculating the quantification values for the different quantification channels.

Table 31. The meaning of colors in the Quantification Spectrum chart (Sheet 2 of 2)

Color	Meaning
Yellow box	Indicates that the pattern includes peaks from indistinguishable channels. This ion pattern peak is not used in calculating the quantification values for the different quantification channels.
Pink bar	Indicates that a quantification pattern peak is expected but is missing. This ion pattern peak is not used in calculating the quantification values for the different quantification channels.
Blue bar	Indicates that a quantification pattern peak is present but is unsuitable because of errors in peptide labeling or because of the wrong centroid retention time, the out-of-delta mass range, the wrong intensity, or a peak that has been used by another isotopic pattern. This ion pattern peak is not used in calculating the quantification values for the different quantification channels.

Displaying the Quantification Spectrum Chart for Reporter Ion Quantification

For reporter ion quantification, the Quantification Spectrum chart displays the intensity of the reporter ions, in counts. It shows a spectrum for each peptide, except for those peptides labeled with “No Quan Values” in the Quan Info column of the PSMs page.

The width of the bars shown in the Quantification Spectrum is determined by the setting of the Integration Tolerance parameter in the Reporter Ions Quantifier node. [Figure 314](#) shows an example of a quantification spectrum from a TMTe 6plex sample quantified with an Integration Tolerance setting of 20 ppm for extracting the reporter peaks from the quantification spectrum.

10 Performing Quantification

Viewing Quantification Results

Figure 314. Quantification Spectrum chart for a TMTe 6plex showing the intensity of reporter ions





Figure 315 shows an example of a quantification spectrum from a TMT 10plex sample quantified with an Integration Tolerance setting of 10 ppm.

Figure 315. TMT 10plex sample quantified with an Integration Tolerance setting of 10 ppm

The Quantification Spectrum chart includes the following features:

- The light blue boxes represent the integration windows for the reporter tags. The boxes are centered on the masses of the reporter tags, as specified in the quantification method. The width of boxes is the integration window used for extracting the reporter tags. It is ± 20 ppm, as specified by the settings of the parameters in the Reporter Ions Quantifier node (you can look up all these values on the Quantification Summary page). The height of the line in the box represents the actual tag intensity used for calculating the peptide ratios. The height of the box represents the corrected tag intensity. The height depends on the setting of the Integration Method parameter specified in the Reporter Ions Quantifier node. It is always the value that results from correction for isotopic impurities, as specified in the TMT Reporter Ion Isotope Distributions area of the Quan Channels page of the Quantification Method Editor dialog box, shown in [Figure 280](#) on [page 418](#).
- To calculate the actual intensity of a particular tag, the application chooses the blue fragment peaks from the spectrum and considers only peaks in the integration window.

- The black fragment peaks represent peaks that are present in the spectrum but that are not chosen for calculating the tag intensities. They might not be chosen because the peaks lie outside of any integration window, or because the setting of the Integration Method parameter specified in the Reporter Ions Quantifier node determined that only one peak per integration window should be chosen from any integration window. A different peak was picked for this integration window according to the criterion specified by the Integration Method setting.

Displaying the Quan Spectra Page

The Quan Spectra page of the .pdResult report lists all the spectra that were used in reporter ion quantification. It displays quantification values for identified spectra and spectra that were quantified but not identified. The Proteome Discoverer application generates the Quan Spectra page only if you include a Reporter Ion Quantification node in your processing workflow. When the consensus workflow includes the Peptide and Protein Quantifier node, the Quan Spectra page contains columns listing the raw quantification values and the average reporter signal-to-noise values, if available.

Figure 316. Quan Spectra page

	Checked	Best PSM Ambiguity	Mass Analyzer	Activation Type	MS Order	# PSMs	Isolation Interference [%]	Ion Inject Time [ms]	Precursor m/z [Da]	Precursor MH+ [Da]	Precursor Charge	RT [min]	First Scan	Spectrum File	(127_N) / ((127_C) / (127_S))
1	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	39	50.000	908.51385	1816.02043	2	172.8601	70520	141007KSAM00279TMT20-R3.raw	1.647 2.248
2	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	37.589	804.97198	1608.93669	2	125.9169	50438	141007KSAM00279TMT20-R3.raw	1.256 0.571
3	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	40.958	940.50275	1879.98222	2	125.9146	50437	141007KSAM00279TMT20-R3.raw	1.049 0.802
4	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	8	50.000	890.48187	2659.43106	3	191.8664	78591	141007KSAM00279TMT20-R3.raw	0.745 0.833
5	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	10	2.732	655.70703	1965.10654	3	171.7154	70035	141007KSAM00279TMT20-R3.raw	1.053 1.023
6	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	33.829	712.37689	1423.74651	2	88.4541	34317	141007KSAM00279TMT20-R3.raw	2.973 1.355
7	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	0	13.328	636.09607	2541.36245	4	125.9054	50433	141007KSAM00279TMT20-R3.raw	0.812 1.535
8	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	3	50.000	1005.60071	2010.19414	2	212.4780	87338	141007KSAM00279TMT20-R3.raw	0.466 0.591
9	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	0	50.000	639.02148	1915.04990	3	156.6958	63629	141007KSAM00279TMT20-R3.raw	1.208 1.325
10	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	10	50.000	700.76596	2100.28214	3	156.6981	63630	141007KSAM00279TMT20-R3.raw	1.624 2.756
11	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	26	50.000	731.39825	1461.78923	2	156.7004	63631	141007KSAM00279TMT20-R3.raw	1.001 0.980
12	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	0	50.000	946.48291	2837.43418	3	156.7050	63633	141007KSAM00279TMT20-R3.raw	1.252 0.465
13	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	23.064	685.71216	1995.12192	3	191.8756	78595	141007KSAM00279TMT20-R3.raw	0.937 0.995
14	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	50.000	881.02405	1761.04082	2	212.4224	87314	141007KSAM00279TMT20-R3.raw	1.498 2.313
15	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	42.646	746.72522	2238.16111	3	125.8773	50421	141007KSAM00279TMT20-R3.raw	1.518 1.969
16	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	7	1.540	657.06702	1969.18650	3	125.9008	50431	141007KSAM00279TMT20-R3.raw	1.512 0.950
17	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	69	9.121	575.65497	1724.95035	3	109.5093	43382	141007KSAM00279TMT20-R3.raw	1.071 1.199
18	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	4	24	3.900	611.84570	2444.36098	4	156.6752	63620	141007KSAM00279TMT20-R3.raw	1.719 1.860
19	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	4.952	740.45117	1479.89507	2	109.5024	43379	141007KSAM00279TMT20-R3.raw	1.351 1.061
20	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	70	4.073	497.29605	1489.87360	3	68.3020	25677	141007KSAM00279TMT20-R3.raw	1.637 0.972
21	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	50.000	946.48309	2837.43473	3	156.6264	63699	141007KSAM00279TMT20-R3.raw	1.204 0.516
22	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	20	50.000	920.03003	1839.06278	2	212.5592	87373	141007KSAM00279TMT20-R3.raw	0.520 1.163
23	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	0	36.077	768.39968	3070.57290	4	88.5123	34342	141007KSAM00279TMT20-R3.raw	1.922 1.239
24	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	5	50.000	1005.60077	2010.19426	2	212.5524	87370	141007KSAM00279TMT20-R3.raw	0.515 0.841
25	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	44.121	844.98328	1688.95928	2	125.9756	50463	141007KSAM00279TMT20-R3.raw	1.763 8.178
26	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	2	1.522	436.25235	1306.74250	3	88.5026	34338	141007KSAM00279TMT20-R3.raw	0.955 1.065
27	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	29	1.874	497.29553	1489.87204	3	68.3331	25690	141007KSAM00279TMT20-R3.raw	1.759 1.147
28	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	9.538	619.87189	1236.73650	2	68.3354	25691	141007KSAM00279TMT20-R3.raw	1.307 0.640
29	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	12	50.000	629.03241	1885.08268	3	171.7499	70050	141007KSAM00279TMT20-R3.raw	0.953 1.041
30	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	50.000	832.46582	1663.92436	2	171.7408	70046	141007KSAM00279TMT20-R3.raw	1.268 1.239
31	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	15.800	940.50886	1880.00444	2	125.9572	50455	141007KSAM00279TMT20-R3.raw	1.139 0.595
32	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	2	0	10.538	627.33984	1880.00498	3	125.9549	50454	141007KSAM00279TMT20-R3.raw	1.031 1.164
33	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	15	9.698	983.05670	1965.10613	2	171.7384	70045	141007KSAM00279TMT20-R3.raw	0.981 0.865
34	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	50.000	791.96692	1582.92656	2	156.6631	63615	141007KSAM00279TMT20-R3.raw	0.956 0.871
35	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	11	2.679	575.65497	1724.95035	3	109.5352	43393	141007KSAM00279TMT20-R3.raw	1.011 1.167
36	<input type="checkbox"/>	Unambiguous	FTMS	HCD	MS2	1	0	50.000	768.90350	1538.78973	2	156.6218	63697	141007KSAM00279TMT20-R3.raw	0.856 0.864

When you select a row on the Quan Spectra page, you can view data in the Quantification Channels view and the Quantification Spectrum view to see the details of the spectra used for quantification.

❖ To display the Quan Spectra page

1. Perform reporter ion quantification.
2. Open the .pdResult report and click the **Quan Spectra** tab.

Quan Spectra Page Parameters


The columns of the Quan Spectra page are basically the same as those on the MS/MS Spectrum Info page (see the Help). However, they also include reporter ion abundances of the intensity of the fragmented tag in a sample to the intensity of the fragmented tag in the control sample for all spectra, regardless of whether they have been identified.

Displaying Quantification Ratio Charts

For both precursor ion and reporter ion quantification, the application calculates protein and peptide group ratios from the protein and peptide group abundances. The ratios are not calculated from PSM ratios.

You can display a chart of the peptide group ratios. This graph shows the distribution of peptide group ratios for the selected protein, displaying the ratios of the peptide groups associated with the selected protein as a log₂-fold change. You must conduct a search with a workflow that includes a quantification node for this chart to appear in the results report.

❖ To display quantification ratio distribution charts

1. On the Proteins page of the .pdResult file containing quantification results, click the row of the protein that you are interested in.
2. Choose **View > Quan Ratio Distributions** or click the **Quan Ratios Distribution** icon, .

The Quan Ratio Distributions chart appears, as shown in [Figure 317](#) for reporter ion quantification, or [Figure 318](#) on [page 477](#) for precursor ion quantification.

10 Performing Quantification

Viewing Quantification Results

Figure 317. Quan Ratio Distributions chart for reporter ion quantification






Figure 318 shows the Quan Ratio Distributions chart for precursor ion quantification.

Figure 318. Quan Ratio Distributions Chart for precursor ion quantification



The Quantification Ratio Distributions chart shows the distribution and spread of the ratios of all peptides belonging to a particular protein. Figure 319 shows an example for the albumin protein.

Figure 319. Quantification Ratio Distributions chart


The chart shows the distribution of peptide group ratios for each of the ratios reported, as defined in the quantification method for this search. Each of the ratio distribution charts displays the peptide group ratios as the binary logarithm. The logarithmic form is common for such displays, because it provides a reasonable display, even when there is a large spread of the displayed values. In binary logarithmic form, a value of 1 means a two-fold increase, a value of 2 means a four-fold increase, a value of 3 means an eight-fold increase, and so forth.

Each of the separate distribution charts displays the peptide group ratios in three sections. The chart legend explains the meaning of these sections. You can access the chart legend by right-clicking the chart and choosing Show Legend.

The Quan Ratio Distribution charts contain the three sections illustrated in [Figure 320](#):

- The first section displays the distribution of the ratios of all peptides considered for calculating the abundance of this protein as a box-and-whisker plot. A box-and-whisker plot is a convenient way of graphically depicting groups of numerical data through a five-number summary: 5 percent lower bound, lower quartile, median, upper quartile, 95 percent upper bound. The range between the lower and upper quartile (this is the range of the box) is also known as the inter-quartile range (IQR) and, like the standard deviation for normally distributed data, is a measure of the spread of the data.
 - The box represents the peptide group ratios between the 25th and the 75th percentiles.
 - The error bars represent the peptide group ratios below the 5th and the 95th percentiles.
 - The blue lines inside the horizontal bar represent the median of the distribution.
- The second section (blue circles) displays the distribution of the ratios of all peptide groups considered in calculating the protein abundances.
- The third section (red circles) displays the distribution of the ratios of all peptide groups that were not considered in calculating the protein abundance (for example, the peptide group is not unique to this protein or this protein group) according to the rules defined in the quantification method.

Figure 320. Quan Ratio Distributions chart for reporter ion quantification




In addition, each chart displays the median ratio (R) and the inter-quartile range (IQR) in linear and logarithmic format. The header of the chart identifies the protein that the PSM belongs to. Right-click the chart and choose **Show Legend** for the identity of other notations on the chart.

Treatment of Missing Reporter Ions for Quantification

If reporter ions are missing in the quantification spectra, you can use the parameter settings of the Peptide and Protein Quantifier node in the consensus workflow to determine how the application handles this problem. For example, suppose that any one of the eight reporter ions from a TMT 8plex experiment is missing or that the reference ion is missing, as [Figure 321](#) shows for the TMT 8plex method. In this latter case, the 126 ion is missing in the Quan Channel Values view opened by the View > Show Quan Channel Values command. In the default quantification consensus workflows, the application still uses the PSM to calculate the peptide and protein quantification values.

Figure 321. Quantification results with missing reporter ions



You can do one of the following:

- Accept spectra with missing reporter ions, report the minimum or maximum allowed fold change, and use the quantification results (for example, $x/0 = 100.000$, $0/x = 0.010$). This option is used in the default consensus quantification workflows.
- Reject the quantification results.
- Replace the missing reporter ions and use the new quantification results.

The following topics discuss these options in more detail.

Accepting Spectra with Missing Reporter Ions

If one or more of the reporter, or mass, tags are missing in the quantification spectrum, the calculated ratios using these reporter tags are either zero or infinity, depending on which tag intensity is the numerator and which is the denominator. Even if all tags are present, the calculated ratios might be very high or very low. You can use the Maximum Allowed Fold Change parameter of the Peptide and Protein Quantifier node to replace such extremely high or extremely low ratios with the maximum or minimum allowed fold change, respectively.

That is, extremely high ratios (greater than 100) are replaced by 100, and extremely low ratios (less than 0.01) are replaced by 0.01. If the Maximum Allowed Fold Change parameter is set to 1000, the application replaces any ratios higher than 1000 by 1000. It replaces any ratios less than 0.001 by 0.001. By default, the Maximum Allowed Fold Change parameter is set to 100. You can find this parameter in the Quantification - General category in the Parameters area shown in [Figure 322](#). You can expect an inherent dynamic range to be valid or detectable with the given instrumentation and method.


Accept spectra with missing reporter ions, use the default settings of the Peptide and Protein Quantifier node parameters. [Figure 322](#) shows these settings.

Figure 322. Default settings of the Peptide and Protein Quantifier node parameters

Parameters	
Show Advanced Parameters	
1. Quantification - General	
Peptides to Use	Unique + Razor
Consider Protein Groups for Peptide Uniqueness	True
Replace Missing Values with Minimum Value	False
Reject Quan Results with Missing Channels	False
Maximum Allowed Fold Change	100
Top N Peptides Used for Area Calculation	3
2. Reporter Quantification	
Reporter Abundance Based On	Automatic
Apply Quan Value Corrections	True
Co-Isolation Threshold	50
Average Reporter S/N Threshold	0
3. Precursor Quantification	
Use Single-Peak Quan Channels	False
4. Normalization and Scaling	
Normalization Mode	None
Proteins For Normalization	
Scaling Mode	None
5. Display Options	
Show Standard Errors	True
Show Quan Value Counts	False
Show Quan Ratios As	Normal Space Values
6. Quan Ratio Distributions	
1st Fold Change Threshold	2
2nd Fold Change Threshold	4
3rd Fold Change Threshold	6
4th Fold Change Threshold	8
5th Fold Change Threshold	10

In the regenerated results shown in [Figure 323](#), the abundance ratios are all 100.000, because the denominator for each ratio is the reporter ion with mass 126.

Figure 323. Effect of accepting PSMs with missing reporter ions



Rejecting the Quantification Results


When you reject any quantification results from spectra with missing reporter ions, the application does not use the corresponding PSMs to calculate the protein and peptide quantification values. The Peptide Quan Usage column of the .pdResult report displays “Not Used,” and the Quan Info column displays “No Quan Values.”

You can use the Reject Quan Results with Missing Channels parameter of the Peptide and Protein Quantifier node to determine whether the application sets all quantification values from all channels to zero if one or more of the quantification channels have a detected intensity of zero. When you set this parameter to True, as shown in Figure 324, the application does not use a PSM for calculating peptide group and protein quantification values if one or more of the quantification channels has a detected intensity of zero. If you set it to False, which is the default, it calculates and displays the quantification values for those channels that the application found and ignores the missing channels.

Figure 324. Parameters of the Peptide and Protein Quantifier node set to reject quantification results

Parameters	
Show Advanced Parameters	
1. Quantification - General	
Peptides to Use	Unique + Razor
Consider Protein Groups for Peptide Uniqueness	True
Replace Missing Values with Minimum Value	False
Reject Quan Results with Missing Channels	True
Maximum Allowed Fold Change	100
Top N Peptides Used for Area Calculation	3
2. Reporter Quantification	
Reporter Abundance Based On	Automatic
Apply Quan Value Corrections	True
Co-Isolation Threshold	50
Average Reporter S/N Threshold	0
3. Precursor Quantification	
Use Single-Peak Quan Channels	False
4. Normalization and Scaling	
Normalization Mode	None
Proteins For Normalization	
Scaling Mode	None
5. Display Options	
Show Standard Errors	True
Show Quan Value Counts	False
Show Quan Ratios As	Normal Space Values
6. Quan Ratio Distributions	
1st Fold Change Threshold	2
2nd Fold Change Threshold	4
3rd Fold Change Threshold	6
4th Fold Change Threshold	8
5th Fold Change Threshold	10

The regenerated results show no abundance ratios. See [Figure 325](#).

Figure 325. Effect of rejecting quantification results from spectra with missing reporter ions

Replacing the Missing Reporter Ions


The third option is to replace the missing reporter ions. Use the Replace Missing Values with Minimum Value parameter of the Peptide and Protein Quantifier node. This parameter determines whether the application replaces missing values with the minimum detected quantification value in the spectrum file. If you set it to True, as shown in Figure 326, the application replaces any missing quantification values by the minimum intensity value detected. If you set it to False, which is the (default), the application sets the abundance of missing values to 0.

Figure 326. Parameters of the Peptide and Protein Quantifier node set to replace missing reporter ions

Parameters	
Show Advanced Parameters	
1. Quantification - General	
Peptides to Use	Unique + Razor
Consider Protein Groups for Peptide Uniqueness	True
Replace Missing Values with Minimum Value	True
Reject Quan Results with Missing Channels	False
Maximum Allowed Fold Change	100
Top N Peptides Used for Area Calculation	3
2. Reporter Quantification	
Reporter Abundance Based On	Automatic
Apply Quan Value Corrections	True
Co-Isolation Threshold	50
Average Reporter S/N Threshold	0
3. Precursor Quantification	
Use Single-Peak Quan Channels	False
4. Normalization and Scaling	
Normalization Mode	None
Proteins For Normalization	
Scaling Mode	None
5. Display Options	
Show Standard Errors	True
Show Quan Value Counts	False
Show Quan Ratios As	Normal Space Values
6. Quan Ratio Distributions	
1st Fold Change Threshold	2
2nd Fold Change Threshold	4
3rd Fold Change Threshold	6
4th Fold Change Threshold	8
5th Fold Change Threshold	10

In the regenerated results shown in [Figure 327](#), the missing reporter ion is replaced by the value 279, and the abundance ratios vary because this value is used as the denominator.

Figure 327. Effect of replacing missing reporter ions



Troubleshooting Quantification

The following procedures can help you obtain optimal results when performing quantification.


❖ To troubleshoot reporter ion quantification

If you obtain unexpected quantification results, verify that all node settings in your processing workflow are reasonable.

- Make sure that the Integration Tolerance parameter of the Reporter Ions Quantifier node fits the data that you are processing. The default is 20 ppm, which is too low if you are processing PQD data from an ion trap.


- Make sure that the settings of the Mass Analyzer, MS Order, and Activation Type parameters of the Reporter Ions Quantifier node are correct for the data that you are processing. [Figure 328](#) shows the typical settings to use if you want to quantify HCD scans from the Orbitrap mass spectrometer.

Figure 328. Typical settings for quantifying iTRAQ or iodo TMT 6plex tags from HCD scans




To quantify PQD scans from an ion trap, use the typical settings shown in [Figure 329](#).

Figure 329. Typical settings for quantifying ITRAQ or TMT tags from the ion trap PQD scans



- Make sure that you have used the correct set of static and dynamic modifications for the search engine. For example, if you are searching iodo TMT 6plex data with Sequest HT, check that your settings resemble [Figure 330](#).

Figure 330. Modifications required for searching iodo TMT 6plex samples



❖ **To troubleshoot precursor ion quantification**

If you obtain unexpected precursor ion quantification results, verify that all settings of your processing workflow are reasonable.

- Check the dynamic modification parameters in the Mascot or Sequest HT search engines.

These should match your isotope labeling sample.

- Check the node parameters that you set before performing the quantification to see if they are appropriate for your sample.

See “Reporter Ion Quantification” on page 385 for more information.

- Verify that your isotopic labeling is one of the following options in the protein ID/search node (either SequestHT or Mascot):
 - SILAC 2plex (Arg10, Lys6): Uses arginine 10 and lysine 6.
 - SILAC 2plex (Arg10, Lys8): Uses arginine 10 and lysine 8.
 - SILAC 2plex (Ile6): Uses isoleucine 6.
 - SILAC 3plex (Arg6, Lys4|Arg10, Lys8): Uses arginine 10 and lysine 8 for “heavy” labels and arginine 6 and lysine 4 for “medium” labels.
 - SILAC 3plex (Arg6, Lys6|Arg10, Lys8): Uses arginine 10 and lysine 8 for “heavy” labels and arginine 6 and lysine 6 for “medium” labels.
 - Dimethylation 3plex: Chemically adds isotopically labeled dimethyl groups to the N-terminus and to the ϵ -amino group of lysine.
 - ^{18}O labeling: Introduces 2 or 4 Da mass tags through the enzyme-catalyzed exchange reaction of C-terminal oxygen atoms with ^{18}O .

Note Low-mass accuracy MS 1 full-scan data cannot be used for precursor ion quantification or precursor ion area detection.

- Check your tolerance window. If you get too many results, decrease the size of the window. For too few results, increase the size of the window.
- Make sure you chose the right database.
- Check the species listed to make sure the samples came from that species.
- Verify that the activation type used is correct.
- Verify that the instrument type in the Mascot search engine is correct.
- Use only the ETD Spectrum Charger node for low-mass resolution ETD data.

FASTA Reference

This appendix lists the most important FASTA databases and parsing rules that the Proteome Discoverer application uses to obtain protein sequences, accession numbers, and descriptions.

Contents

- [FASTA Databases](#)
- [Custom Database Support](#)

FASTA Databases

These are the most important FASTA databases that the application uses.

- [NCBI](#)
- [UniRef100](#)
- [SwissProt and TrEMBL](#)

Follow the links given for each database if you would like to download the database and save it to your local machine. Some databases are more time-consuming to load than others.

NCBI

NCBI is a nonredundant database compiled by the NCBI (National Center for Biotechnology Information) as a protein database for Blast searches. It contains nonidentical sequences from GenBank CDS translations, Protein Data Bank (PDB), SwissProt, Protein Information Resource (PIR), and Protein Research Foundation (PRF).

<http://blast.ncbi.nlm.nih.gov/Blast.cgi>

<ftp://ftp.ncbi.nih.gov/blast/db/FASTA/nr.gz>

A typical NCBI title line follows:

```
>gi|70561|pir||MYHO myoglobin - horse_i|418678|pir||MYHOZ myoglobin -  
common zebra (tentative sequence) [MASS=16950]
```

FASTA ID:

- Accession#:gi70561
- Description:myoglobin - horse_i

UniRef100

UniRef, also known as UniProt NREF, is a set of comprehensive protein databases curated by the Universal Protein Resource consortium. UniRef100 contains only nonidentical sequences, and UniRef90, and UniRef50 are nonredundant at a sequence similarity level of 90 percent and 50 percent, respectively.

ftp://ftp.uniprot.org/pub/databases/uniprot/current_release/uniref/uniref100/

A typical UniRef100 title line follows:

```
>UniRef100_4U9M9 Cluster: 104 kDa microneme-rhoptry antigen precursor;  
n=1; Theileria annulata|Rep: 104 kDa microneme-rhoptry antigen  
precursor - Theileria annulata
```

FASTA ID:

- Accession#:4U9M9
- Description:Cluster: 104 kDa microneme-rhoptry antigen precursor;
n=1; Theileria annulata|Rep: 104 kDa microneme-rhoptry antigen
precursor - Theileria annulata

SwissProt and TrEMBL

The SwissProt database is developed by the SwissProt groups at the Swiss Institute of Bioinformatics (SIB) and the European Bioinformatics Institute (EBI).

TrEMBL is a computer-annotated supplement of SwissProt that contains all the translations of EMBL nucleotide sequence entries not yet integrated into SwissProt.

<http://www.expasy.org/sprot/>

ftp://ftp.ebi.ac.uk/pub/databases/uniprot/knowledgebase/uniprot_trembl.fasta.gz

A typical SwissProt title line follows:

```
>43495|108_SOLLC Protein 108 precursor - Solanum lycopersicum (Tomato)  
(Lycopersicon esculentum)
```

FASTA ID:108_SOLLC

- Accession#:43495
- Description:Protein 108 precursor - Solanum lycopersicum (Tomato)
(Lycopersicon esculentum)

Custom Database Support

The Proteome Discoverer application has three “general” parsing rules to support custom sequence database formats. The generic parsing rules are applied only if no other parsing rule matches the given FASTA title line.

- [Custom Parsing Rule A](#)
- [Custom Parsing Rule B](#)
- [Custom Parsing Rule C](#)

Custom Parsing Rule A

The application uses custom parsing rule A if the FASTA ID, the accession#, and the description are separated by a pipe (|) symbol. A typical FASTA title line, which matches this parsing rule, would look like this one:

```
>tr|18FC3|18FC3_HALWD IS1341-type transposase - Halouadrum walsbyi (strain DSM 16790).
```

```
FASTA ID:18FC3_HALWD
```

- Accession#:18FC3
- Description:IS1341-type transposase - Halouadrum walsbyi (strain DSM 16790).

Custom Parsing Rule B

The application uses custom parsing rule B if the accession# and the description are separated by a space. A typical FASTA title line, which matches this parsing rule, would look like this one:

```
>HP0001 hypothetical protein {Helicobacter pylori 26695}
```

```
FASTA ID:
```

- Accession#:HP0001
- Description:hypothetical protein {Helicobacter pylori 26695}

Custom Parsing Rule C

The application uses custom parsing rule C if the FASTA title line only contains the accession#. A typical FASTA title line, which matches this parsing rule, would look like this one:

```
>143B_HUMAN
```

```
FASTA ID:
```

- Accession#:143B_HUMAN
- Description:143B_HUMAN

Chemistry References

The tables in this appendix list amino acid symbols and mass values, enzyme cleavage properties, and the fragment ions used in the Proteome Discoverer application.

Contents

- [Amino Acid Mass Values](#)
- [Enzyme Cleavage Properties](#)
- [Fragment Ions](#)

Amino Acid Mass Values

The application uses the amino acid symbols and mass values listed in [Table 32](#) and [Table 33](#).

Table 32. Amino acid mass values (Sheet 1 of 2)

Amino acid	One-letter code	Three-letter code	Monoisotopic mass	Average mass	Sum formula
Glycine	G	Gly	57.02146	57.05177	C ₂ H ₃ NO
Alanine	A	Ala	71.03711	71.07855	C ₃ H ₅ NO
Serine	S	Ser	87.03203	87.07796	C ₃ H ₅ NO ₂
Proline	P	Pro	97.05276	97.11623	C ₅ H ₇ NO
Valine	V	Val	99.06841	99.13211	C ₅ H ₉ NO
Threonine	T	Thr	101.04768	101.10474	C ₄ H ₇ NO ₂
Cysteine	C	Cys	103.00918	103.14464	C ₃ H ₅ NOS
Isoleucine	I	Ile	113.08406	113.15890	C ₆ H ₁₁ NO
Generic Leu/Ile	J	Jli	113.08406	113.15890	C ₆ H ₁₁ NO
Leucine	L	Leu	113.08406	113.15890	C ₆ H ₁₁ NO
Asparagine	N	Asn	114.04293	114.10354	C ₄ H ₆ N ₂ O ₂
Aspartic Acid	D	Asp	115.02694	115.08826	C ₄ H ₅ NO ₃

Table 32. Amino acid mass values (Sheet 2 of 2)

Amino acid	One-letter code	Three-letter code	Monoisotopic mass	Average mass	Sum formula
Glutamine	Q	Gln	128.05858	128.13032	C ₅ H ₈ N ₂ O ₂
Lysine	K	Lys	128.09496	128.17358	C ₆ H ₁₂ N ₂ O
Glutamic Acid	E	Glu	129.04259	129.11504	C ₅ H ₇ NO ₃
Methionine	M	Met	131.04048	131.19820	C ₅ H ₉ NOS
Histidine	H	His	137.05891	137.14062	C ₆ H ₇ N ₃ O
Phenylalanine	F	Phe	147.06841	147.17571	C ₉ H ₉ NO
Arginine	R	Arg	156.10111	156.18707	C ₆ H ₁₂ N ₄ O
Tyrosine	Y	Tyr	163.06333	163.17512	C ₉ H ₉ NO ₂
Tryptophan	W	Trp	186.07931	186.21220	C ₁₁ H ₁₀ N ₂ O

Table 33. Special amino acids

Amino acid	One-letter code	Three-letter code	Monoisotopic mass	Average mass	Sum formula
Avg. N/D	B	Bnd	114.53494	114.59590	
Avg. Q/E	Z	Zqe	128.55059	128.62326	
Unknown acid (X)	X	Xxx	110	110	N/A
Pyrrolysine	O	Pyl	237.14773	237.29874	C ₁₂ H ₁₉ N ₃ O ₂
Seleno Cysteine	U	Sec	150.95364	150.0369	C ₃ H ₅ NOSe

Enzyme Cleavage Properties

Table 34 lists the enzymes and reagents with cleavage properties.

Table 34. Cleavage properties of enzymes and reagents (Sheet 1 of 2)

Enzymes/Reagents	Cleaves after	Cleaves before	Except when
Enzymes for digestion			
AspN		D	
Chymotrypsin	F, W, Y, or L		
Chymotrypsin (FWY)	F, W, or Y		P is after F, W, or Y
Clostripain	R		
Elastase	A, L, I, or V		P is after A, L, I, or V
Elastase/Tryp/Chymo	A, L, I, V, K, R, W, F, or Y		P is after AL, L, I, V, K, R, W, F, or Y

Table 34. Cleavage properties of enzymes and reagents (Sheet 2 of 2)

Enzymes/Reagents	Cleaves after	Cleaves before	Except when
GluC	E or D		
LysC	K		
No-Cleavage			
No-Enzyme			
Proline_Endopept	P		
Staph_protease	E		
Trypsin	K or R		P is after K or R
Trypsin (KRLNH)	K, R, L, N, or H		
Trypsin_K	K		P is after K
Trypsin_R	R		P is after R
Chemicals for degradation			
Cyanogen bromide	M		
Iodobenzoate	W		

Fragment Ions

Several different fragmentation techniques produce fragment ions of peptides, such as ECD, ETD, CID, EThcD, higher-energy C-trap dissociation (HCD), and infrared multi-photon dissociation (IRMPD).

As an example, MS/MS and ESI generate low-energy CID spectra, which are sequence-specific. The fragment ion spectra contain peaks of the fragment ions formed by the cleavage of the peptide bond and are used to determine amino acid sequences. A fragment must have at least one charge for it to be detected.

The fragment ions produced are identified according to where they are fragmented in the peptide. A, b, and c fragment ions have a charge on the N-terminal side, and x, y, and z fragment ions have a charge on the C-terminal side. Fragment ions a*, b*, and y* are ions that have lost ammonia (–17 Da), and fragment ions a^o, b^o, and c^o are ions that have lost water (–18 Da). The subscript next to the letter indicates the number of residues in the fragment ion.¹

Table 35 summarizes the fragment ions used in the Proteome Discoverer application.

¹ For more information about fragment ions and nomenclature, see Roepstorff, P.; Fohlman, J. Proposal for a Common Nomenclature for Sequence Ions in Mass Spectra of Peptides. *Biomed. Mass Spectrom.* **1984**, *11* (11) 601.

Table 35. Fragment ions

Ions	Description
a	A ion with charge on the N-terminal side
b	B ion with charge on the N-terminal side
c	C ion with charge on the N-terminal side
y	Y ion with charge on the C-terminal side
z	Z ion with charge on the C-terminal side
a*	A ion that has lost ammonia (-17 Da)
b*	B ion that has lost ammonia (-17 Da)
y*	Y ion that has lost ammonia (-17 Da)
a ^o	A ion that has lost water (-18 Da)
b ^o	B ion that has lost water (-18 Da)
c ^o	C ion that has lost water (-18 Da)

Technical and Biological Replicates

This appendix explains what replicates are and how to specify them in your experiments.

Contents

- [Definitions](#)
- [Specifying Replicates](#)

Definitions

Replicates are measurements of samples with the same conditions for the study variables that differentiate them within your analysis.

- Biological replicates are samples that should be identical (as much as you can or want to control) but are biologically separated (different cells, different individuals, different organisms, different population, and so forth).
- Technical replicates are samples taken from the same biological source and should be identical but technically you repeated the measurement, the chromatography, or the sample preparation.

Specifying Replicates

In the Proteome Discoverer application, you usually do not explicitly specify that two samples are replicates. Instead, you specify the relevant factors for each sample and the factors you want to use to group your samples.

For example, consider a TMT 10plex file with different time points after drug application and sample preparation methods used. [Table 36](#) summarizes this data.

Table 36. Data from a TMT 10plex raw data file

20140511_TMT10_drugA_0-5-10_sampleprep_1-2.raw			
Sample type	Quantification Channel	Time point	Sample prep method
Control	126	0 min	SP_1
Sample	127_C	5 min	SP_1
Sample	127_N	5 min	SP_1
Sample	128_C	10 min	SP_1
Sample	128_N	10 min	SP_1
Control	129_C	0 min	SP_2
Sample	129_N	5 min	SP_2
Sample	130_C	5 min	SP_2
Sample	130_N	10 min	SP_2
Sample	131	10 min	SP_2

The study variables in this example are quantification channels, time points, and sample preparation methods.

If you do not select any study variables to ignore for sample separation, you would create the following quantification ratios:

- (127_C, 5 min, SP_1)/(126, 0 min, SP_1)
- (127_N, 5 min, SP_1)/(126, 0 min, SP_1)
- (128_C, 10 min, SP_1)/(126, 0 min, SP_1)
- (128_N, 10 min, SP_1)/(126, 0 min, SP_1)
- (129_N 5 min, SP_2)/(129_C, 0 min, SP_2)
- (130_C, 5 min, SP_2)/(129_C, 0 min, SP_2)
- (130_N, 10 min, SP_2)/(129_C, 0 min, SP_2)
- (131, 10 min, SP_2)/(129_C, 0 min, SP_2)

If you group only by time point, you effectively ignore the potential effect of introducing two different sample preparation methods:

- (5 min)/(0 min)
 - (127_C, 5 min, SP_1)/(126, 0 min, SP_1)
 - (127_N, 5 min, SP_1)/(126, 0 min, SP_1)
 - (129_N 5 min, SP_2)/(129_C, 0 min, SP_2)

- (130_C, 5 min, SP_2)/(129_C, 0 min, SP_2)
- (10 min)/(0 min)
 - (128_C, 10 min, SP_1)/(126, 0 min, SP_1)
 - (128_N, 10 min, SP_1)/(126, 0 min, SP_1)
 - (130_N, 10 min, SP_2)/(129_C, 0 min, SP_2)
 - (131, 10 min, SP_2)/(129_C, 0 min, SP_2)

If you group by time point and sample preparation method, you effectively consider both the potential effect of introducing two different sample preparation methods and the effect of applying the drug at two different time points:

- (5 min, SP_1)/(0 min, SP_1)
 - (127_C, 5 min, SP_1)/(126, 0 min, SP_1)
 - (127_N, 5 min, SP_1)/(126, 0 min, SP_1)
- (5 min, SP_2) / (0 min, SP_2)
 - (129_N 5 min, SP_2)/(129_C, 0 min, SP_2)
 - (130_C, 5 min, SP_2)/(129_C, 0 min, SP_2)
- (10 min, SP_1)/(0 min, SP_1)
 - (128_C, 10 min, SP_1)/(126, 0 min, SP_1)
 - (128_N, 10 min, SP_1)/(126, 0 min, SP_1)
- (10 min, SP_2)/(0 min, SP_2)
 - (130_N, 10 min, SP_2)/(129_C, 0 min, SP_2)
 - (131, 10 min, SP_2)/(129_C, 0 min, SP_2)



Index

A

- a* fragment ions 497–498
- a^o fragment ions 497–498
- a fragment ions
 - definition 7, 498
- Abort icon 177
- Abundance Counts column 442–443
- Abundance Ratios (log2) column 454
- Abundance Ratios column 454
- abundances
 - grouped 446, 459
 - normalized 459
 - peptide group 442–443
 - protein 443, 448
 - PSM 432
 - raw 459
 - ungrouped 459
- Abundances column 442
- accession keys 330
- Accession Rules 210
- Accession Rules area 211
- acetylation 357
- activation (fragmentation) types 6
- Activation Type parameter
 - Reporter Ions Quantifier node 402, 487
 - Scan Event Filter node 131–132
- Add Files dialog box 49, 100
- Add Files icon 65, 100
- Add Fraction dialog box 50
- Add Fraction icon 65
- Add icon 177, 201, 211
- Add Protein References page 181, 192
- adding and removing amino acids 232
- Affects All Levels command 257
- Affects First Level Only command 257
- amino acids
 - adding to chemical modifications 232
 - deleting from chemical modifications 233
 - in Chemical Modifications view 229
 - mass values 495
 - symbols 495
- analyses
 - adding input files to 72
 - adding multiple processing steps 73
 - Analysis window
 - features 71
 - parameters 80
 - creating 71
 - creating consensus workflows 78
 - creating processing workflows 74
 - description 33
 - duplicating existing processing step 73
 - factory default templates 119
 - opening existing 72
 - opening from results 97
 - re-adding deleted consensus step 74
 - reprocessing existing 101
 - reprocessing existing MSF results 102
 - saving 63, 80
 - searching analysis results for specified string 98
 - using only the processing step 74
- Analysis Results page 68, 98
- Analysis Sequence Details window 97
- Analysis Validation Issues dialog box 92
- Analysis window
 - features 71
 - parameters 80
- Annotate Flanking Residues of the Peptide parameter 360
- Annotated Sequence column 360
- annotated spectra 6
- Annotation Aspect Editor 305, 315
- annotation aspects
 - Annotation Aspect Editor
 - description 305
 - parameters 315
 - deactivating user-defined 311
 - defining 304–305

- deleting user-defined 314
- displaying 304
- displaying user-defined GO protein 300
- editing user-defined 313
- exporting user-defined 314
- filtering by GO Slim categories 301
- importing user-defined from XML file 314
- parameters on Annotation Aspects view 317
- saving user-defined 310
- updating database before defining new 315
- Annotation Aspects parameter 312
- Annotation Aspects view 304, 317
- antioxidant activity GO Slim category 339
- Any filter 111
- Apply Strict Parsimony Principle parameter 284
- arginine 383
- As Batch check box 80, 95, 100–101
- As Species Map parameter 219
- As Species Names parameter 219
- Auto Layout icon 106, 114
- Automatic parameter 23
- autosamplers 158–159
- Average Reporter S/N column, 450
- Average Reporter S/N Threshold parameter 433, 450
- Avoid Expression Rules area 210–211
- Avoid Expressions parameter 214–215

B

- b* fragment ions 497–498
- b^o fragment ions 497–498
- b fragment ions
 - definition 7, 498
 - generated by CID 3
- batch mode 94
- batch processing
 - creating workflow 145
 - definition 146
 - in Proteome Discoverer Daemon utility 143, 146, 148, 160
- biological process codes 336
- Biological Process column 291, 296
- biological processes 335, 343
- Biological Processes page 290, 327, 335
- biological replicates 499
- blank sample types 54
- Blast searches 491
- Browse for Program dialog box 156
- Bulk Ratio Generation area 84

C

- c^o fragment ions 497–498
- c fragment ions
 - abstracting proton from precursor 7
 - definition 7, 498
 - generated by ETD and ECD 3
- catalytic activity GO Slim category 339
- categorical factor dialog box 44
- categorical study factors 44
- cell communication GO Slim category 343
- cell death GO Slim category 343
- cell differentiation GO Slim category 343
- cell division GO Slim category 344
- cell growth GO Slim category 344
- cell homeostasis GO Slim category 344
- cell motility GO Slim category 344
- cell organization and biogenesis GO Slim category 344
- cell proliferation GO Slim category 344
- cell surface GO Slim category 340
- cellular component codes 335
- Cellular Component column 291, 296
- cellular components 334, 340
- Cellular Components page 290, 327, 334
- Change Instruments In Use dialog box 158
- Check for Updates icon 174, 178
- chemical modifications
 - adding 231
 - adding amino acids 232
 - categories of 229
 - deleting 234
 - deleting amino acids 233
 - displaying 228
 - dynamic 227
 - importing from local file 234
 - importing from UNIMOD 234
 - static 227
 - updating existing 232
- Chemical Modifications view
 - adding modifications 231
 - description 227–228
 - displaying 228
 - Position column 232
- Child Steps bar 71
- Chromosome column 296, 322–323
- chromosome GO Slim category 340
- CID activation type
 - analyzed by Sequest HT 3
 - description 6
 - fragmenting ions 7, 497
- Clear All Filter Criteria icon 187
- Clear All Filters icon 196

- Clear All icon 87, 89, 258
 - Clear icon 258, 267
 - Clear Recent Results File List dialog box 20
 - Clear Reference Filter Criteria icon 186
 - Clear Reference Filter icon 196
 - cleavage enzyme specificities 239, 243
 - cleavage properties 496
 - cleavage reagents
 - adding 243
 - cleavage properties 496
 - deleting 243
 - displaying 238
 - enzyme specificities 239
 - filtering data 244
 - modifying 244
 - Cleavage Reagents view
 - description 237
 - opening 238
 - parameters on 240
 - coagulation GO Slim category 344
 - coisolation 435
 - Co-IsolationThreshold parameter 436
 - Common Templates folder
 - opening
 - directly in Windows Explorer 120
 - through Documents library in Windows Explorer 120
 - workflow templates in 119
 - common workflow templates 122, 126
 - Compile FASTA Database page 187, 191, 193
 - Concatenated FDR Calculation parameter 278
 - confidence indicators
 - distribution in target false discovery rates 278
 - on Peptides page 278
 - Configuration page
 - specifying server to connect to 144
 - Configuration view
 - configuring search engines 21
 - Proteome Discoverer Daemon utility area 150
 - configuring Mascot search engine 25, 27
 - configuring Sequest HT search engine 22
 - conjugation GO Slim category 344
 - Consensus Step box 71
 - consensus workflow
 - creating basic 78, 112, 114
 - decoy database searches 280
 - description 34
 - factory default templates 119
 - marked contaminants searches 217
 - precursor ion area detection 406
 - precursor ion quantification 394
 - protein annotation 297
 - PTM analysis 354
 - reporter ion quantification 400
 - reprocessing existing MSF results by using new 102
 - saving 79, 115
 - specifying number executed in parallel 15, 32
 - spectrum library 225
 - uploading search results to ProteinCenter 324
 - Consider Protein Groups for Peptide Uniqueness parameter 439, 442, 448
 - contacting us xiv
 - Contains icon 61
 - Contaminant column 218
 - contaminants 217
 - control sample types 54, 68
 - Coverage column 442
 - Create New Quantification Method dialog box 41, 412, 417, 420
 - Create New Study/Analysis icon 38
 - Custom Filter dialog box 185, 244
 - Custom option 185
 - custom parsing rule A 493
 - custom parsing rule B 493
 - custom parsing rule C 493
 - custom workflow templates 123, 126
 - cytoplasm GO Slim category 340
 - cytoskeleton GO Slim category 340
 - cytosol GO Slim category 340
- ## D
- data distribution maps
 - displaying 457
 - displaying samples 455
 - displaying species names 219
 - filtering data in 454
 - for files 455
 - for sample groups 454–455
 - for samples 455
 - Found in Files column 453, 457
 - Found in Fractions column 453, 457
 - Found in Sample Groups column 453, 457
 - Found in Samples column 453, 457
 - generating in workflows with quantification nodes 453, 456
 - purpose 5
 - sorting data in 459
 - Data Distributions node
 - description 453
 - generating data distribution maps in workflows with quantification nodes 453
 - precursor ion quantification workflow 395
 - reporter ion quantification workflow 401

deconvolution 435

decoy database searches
 calculating false discovery rates 277
 consensus workflow 280
 processing workflow 279
 setting up FDRs 279

defense response GO Slim category 344

Delete Jobs dialog box 180, 224

Delete Methods dialog box 426

Deleted FASTA Indexes table 206–207

Delta Cn column 246–247

Delta Cn value 246–247

development GO Slim category 344

device drivers xiii

dimethyl modifications 358

dimethylation 3plex quantification method 382, 385, 490

Diseases page 327, 337

Display Filter command/icon 254

Display Filter node
 applying a filter set to another .pdResult file 268
 filtering with 267

Display Filter pane 255

display filters
 applying 257
 applying filters specific to different search engines 276
 clearing 258
 definition 254
 examples 261
 filter sets 254, 258–259, 266–267
 filtering annotation aspects by GO Slim categories 301
 filtering data distribution maps 454
 filtering data in .pdResult files 255
 filtering numerical values by specified precision value 270
 filtering out PSMs on Protein Identification Details view 371
 filtering phosphorylation site probabilities 374
 finding common or unique proteins in multiple searches 272
 finding proteins common to two or more searches 272
 finding proteins unique to each search 274
 loading filter sets 268
 purpose 245
 setting additional conditions 256
 setting additional filters 257
 setting the first condition 256
 temporarily disabling 257
 viewing settings in the Result Summaries view 260

Display Filters command/icon 255

DNA 339–340

DNA binding GO Slim category 339

documentation
 accessing xi
 additional xi

Documents library 120–121

Documents Library Locations dialog box 121

dot score 224

Download from ProteinCenter dialog box 174, 176

Download icon 174

DTA files 12

Dynamic Modification parameter
 Mascot 429
 Sequest HT 75, 393, 399, 405

dynamic modifications
 definition 227

E

EC number 334–336

ECD activation type
 analyzed by Sequest HT 3
 description 7
 fragmenting ions 497

Edit icon 211, 425

Edit Parameter Text dialog box 253

emPAI algorithm 448

endoplasmic reticulum GO Slim category 341

endosome GO Slim category 341

Ensembl Gene ID column 322–323

Ensembl genome database
 description 5, 292, 296
 retrieving annotation results from ProteinCenter 322
 retrieving information from ProteinCenter 5, 290, 292, 296, 299

Ensemble Gene ID column 296

Entrez Gene database
 description 5, 292, 296
 retrieving annotation results from ProteinCenter 321
 retrieving information from ProteinCenter 5, 290, 292, 296, 299
 website 330

Entrez Gene ID column 292, 296, 321

Enzymatic Unspecific cleavage enzyme specificity 239, 243

Enzyme Category (EC) information 333

enzyme regulator activity GO Slim category 339

enzymes 496

ESI
 fragment ions 497
 peptides and fragment ions 7
 with PQD 7

ETciD activation type 131

ETD activation type
 analyzed by Sequest HT 3
 description 7
 fragmenting ions 497

ETD Spectrum Charger node 490

ETHeD activation type
 analyzed by Sequest HT 3
 description 7
 fragmenting ions 497
 searching ETHcD spectra 131
 searching HCD spectra 131

European Bioinformatics Institute (EBI) 5, 291–292, 333

Event Detector node
 precursor ion area detection 404
 precursor ion quantification workflow 391–392

evidence codes 333, 335–336

Excel
 exporting Modification List page to 374

Expand icon 457

exponentially modified protein abundance index. *See* emPAI

Export icon 258, 427

Export Parameter File page 153, 163

Export to ProteinCenter dialog box 325

External Links page 290, 327, 337

extracellular GO Slim category 341

F

false discovery rate
 calculating for
 peptide groups 281
 protein groups 282
 proteins 281
 PSMs 281
 calculating for PSMs 281
 definition 277
 relaxed
 peptide confidence indicators 279
 specifying for decoy database search 278
 setting up 279
 strict
 peptide confidence indicators 279
 specifying for decoy database search 278
 target 278
 validating MS/MS searches 277

FASTA Database Utilities command 181–182, 187

FASTA Database Utilities dialog box
 Add Protein References page 181, 192
 Compile FASTA Database page 187, 191, 193
 filtering protein references 184
 Find Protein References page 182, 195
 opening 181
 parameters on 192

FASTA files
 adding protein references 181
 adding protein sequences 181
 adding to Proteome Discoverer application 174
 compiling 187
 deleting 180
 downloading from ProteinCenter 174
 downloading from sources other than ProteinCenter 177
 excluding protein sequences or references 191
 exporting 180
 filtering protein reference searches 184
 finding protein sequences or references 182
 indexes. *See* FASTA indexes
 listing available 172
 most important databases 491
 parsing rules. *See* FASTA parsing rules
 proteins marked as contaminants 217
 refining search for protein sequences or references 185
 selecting with
 Mascot search engine 108
 Sequest HT search engine 108
 title lines 209–210, 213
 updating from ProteinCenter 178
 use in normalizing protein abundances 443

FASTA Files view
 Abort icon 177
 Add icon 177
 Check for Updates icon 174, 178
 Download icon 174
 opening 172
 Update icon 179

FASTA Index Creator dialog box 201

FASTA indexes
 automatically creating 201
 changing maximum number stored 207
 changing storage location 207
 controlling automatic removal 205
 definition 196
 deleting 206
 deleting from deleted FASTA files 209
 discarding changes from previous session 209
 displaying 199
 manually creating 201
 no-enzyme searches 197
 non-specific digestion 201
 restoring deleted 207
 semi-enzymatic searches 197, 201
 specifying location and number stored 197

FASTA Indexes Options dialog box 208

FASTA Indexes view
 choosing columns to display 201

- displaying after adding FASTA index 204
- opening 199
- FASTA parsing rules
 - A 493
 - adding 209
 - adding regular expressions 210–211
 - applying changes 212
 - B 493
 - C 493
 - categories of 210
 - changing or renaming regular expressions 211
 - displaying available 210
 - modifying 209
 - storage locations 215
 - testing 211
 - XML file 215
- FASTA Parsing Rules view
 - features of 210
 - opening 209
- FASTA title lines 212–213, 215
- FastaTitleParsingRules.xml file 215
- FDRs. *See* false discovery rate
- Features page 290, 327, 331
- FidoCT Protein Validator node
 - calculating FDRs for protein groups 282
 - calculating FDRs for proteins 281
- Field Chooser dialog box
 - FASTA Indexes view 201
- Field Chooser icon 201
- Filter Set parameter 268
- filter sets
 - applying to another .pdResult file 268
 - factory defaults 119
 - loading 259, 266, 268
 - saving 254, 258, 266–267
 - specifying default for Proteins page of .pdResult file 269
- filtering
 - data in .pdResult files 245, 255
 - false discovery rates 277
 - filter sets. *See* filter sets
 - protein reference searches 184
 - proteins 254
 - PSMs 246
 - specifying score names to filter by 253
 - with
 - Display Filter node 267
 - display filters. *See* display filters
 - Fixed Value PSM Validator node 248
 - Maximum Delta Cn parameter of the MSF Files node 245, 249
 - Maximum Delta Cn parameter of the Percolator node 245, 247
 - Maximum Delta Cn parameter of the Target Decoy PSM Validator node 245
 - Maximum Delta Mass parameter of the MSF Files node 251
 - Peptide and Protein Filter node 254
 - Score and Threshold parameters of the MSF Files node 251
 - Target Decoy PSM Validator node 248
- .filterSet file 267
- Find Protein References page 182, 195
- Fixed Value PSM Validator node
 - creating basic processing workflow 105
 - filtering PSMs 248
 - Maximum Delta Cn parameter 393
 - parallel workflow 131
 - precursor ion area detection workflow 404
 - precursor ion quantification workflow 391, 393
 - PTM analysis workflow 352
 - reporter ion quantification 398
 - spectrum library processing workflow 225
- flanking residues 359–360
- Fold Change Threshold parameters 454
- For All Proteins parameter 359
- Found in Files column 453, 457
- Found in Fractions column 453, 457
- Found in Sample Groups column 453, 457
- Found in Samples column 453, 457
- fractions
 - adding to study 50
 - creating subsets of 50
 - definition 49
 - MudPIT files 49
- fragment ions
 - activation types producing 497
 - ammonia loss 497
 - charged on C-terminal side 497
 - charged on N-terminal side 497
 - factors dependent on 7
 - post-translational modifications (PTMs) 347
 - resulting from
 - CID activation 3
 - ECD activation 3
 - ETD activation 3
 - HCD activation 3
 - PQD activation 3
 - types 6–7, 497
 - water loss 497
- fragmentation methods. *See* activation types
- From Quan Method parameter 428–429
- Full ¹⁸O labeling quantification method 382, 385
- Full cleavage enzyme specificity 239, 243

G

GenBank database 491
 Gene Ontology (GO) database. *See* GO database
 General page 290, 327–328
 Generate R-Friendly Headers parameter 16
 Generated Sample Groups area 84
 glycosylation 359
 GO accessions
 description 317
 displaying annotation results from ProteinCenter 318
 GO Accessions column 296, 318
 GO codes 333, 335–336
 GO consortium 291
 GO database
 description 5, 291, 296
 features contrasted with Pfam features 292
 GO Slim annotations 291
 GO terms 291
 retrieving information from ProteinCenter 5, 290–291,
 296, 299
 website 291
 GO protein annotation aspects 300
 GO Slim annotations 291, 296
 GO Slim categories
 biological processes 343
 cellular components 340
 colors 300
 definition 334–336
 filtering annotation aspects by 301
 molecular functions 339
 GO terms 291, 334
 golgi GO Slim category 341
 grouped abundances 446, 459, 462, 465
 Grouping & Quantification page
 areas on 83
 opening 82
 parameters on 70
 purpose 82
 saving settings on 91
 when displayed 72
 grouping samples 84

H

HCD activation type
 analyzed by Sequest HT 3
 description 7
 fragmenting ions 497
 specifying with
 Reporter Ions Quantifier node 487
 Heavy isotopic peak patterns 467
 Hidden Markov Model 292

HTML files
 contents of exported 6
 Human Proteome Organization (HUPO) 11–12

I

Import icon 234, 426
 Import Modifications dialog box 235
 Incomplete ¹⁸O labeling quantification method 382, 385
 infrared multi-photon dissociation. *See* IRMPD activation
 type
 input files
 adding to analysis 72
 adding to study 48
 Input Files page 67
 input files to Proteome Discoverer application 11
 Integration Method parameter 473–474
 Integration Tolerance parameter 471, 486
 intensity 432
 International Union of Biochemistry and Molecular Biology
 334–336
 InterPro database 332
 iodo TMT 6plex quantification method
 as default 386
 reporter ion masses 387
 ion trap mass spectrometer
 processing PQD data from 486–487
 Ions Score value
 PSMs page 135
 IRMPD activation type 7, 497
 Is filter 111
 Is Not filter 111
 isobaric tags for relative and absolute quantification. *See*
 iTRAQ quantification
 isobarically labeled quantification. *See* reporter ion
 quantification
 isoforms
 calculating confidence sequence probability 350
 displaying confidence estimate of 356
 modification 359
 phosphorylation 349
 probability of being correct 356
 pRS score of 350
 putative 348, 351
 isotope patterns 451
 isotopically labeled quantification. *See* precursor ion
 quantification
 isotopomers 382, 385
 iTRAQ 4plex quantification method 4, 389
 iTRAQ 8plex quantification method 4, 389

Index: J

iTRAQ quantification
default methods available in 389
description 4, 388
isotopic distribution values 434
Reporter Ions Quantifier node 4
specifying label modifications 410
See also reporter ion quantification
IUBMB Enzyme Nomenclature 334–336

J

job queue
in creating FASTA indexes 203
Search ID command 135
Job Queue page 148

K

Keys page 290, 327, 330

L

LC/MS 8
Light isotopic peak patterns 467
Load Failures pane 124
Load icon 259
lock symbol 102
low-resolution iodo TMT 6plex quantification method 15, 386
low-resolution TMTe 6plex quantification method 387
LTQ Orbitrap mass spectrometers
troubleshooting 487
LTQ Orbitrap output template 119
lysine 383

M

Magellan server 2, 144, 149
Maintain Chemical Modifications command/icon 228, 231
Maintain Cleavage Reagents command 238
Maintain FASTA Files command/icon 172, 174, 177–178, 180
Maintain FASTA Indexes command 201, 204
Maintain FASTA Parsing Rules command 209
Maintain Quantification Methods command/icon
adding quantification method 41, 420
changing quantification method 424–425
deactivating quantification method 426
exporting quantification method 427
importing quantification method 426
removing quantification method 426

restoring quantification method template defaults 427
specifying channels for precursor ion quantification 411, 417
Maintain Spectrum Libraries command 222, 224
MALDI 7
Manual Ratio Generation area 84
marked contaminants searches
creating consensus workflow 133, 217
creating processing workflow 133, 217
description 217
Mascot Communication Tester 26
Mascot Generic Format files. *See* MGF files
Mascot node
configuring 29
From Quan Method parameter 428–429
parallel processing workflow 131
PTM analysis workflow 352
See also Mascot search engine
Mascot search engine
adding FASTA files before searching 174
applying filters in multiconsensus reports 276
configuring 27
configuring parameters for 21, 25
configuring protein-scoring calculation 30
creating a processing workflow 105
description 4
Ions Score value
PSMs page 135
Mascot server 26
merging results 138
precursor ion quantification workflow 393
reporter ion quantification workflow 398
searching for quantification modifications 428
searching multiple sequence databases 134
selecting FASTA files 108
testing connection between Proteome Discoverer application and Mascot server 26
troubleshooting failed searches 31
Mascot server 26
Mass Analyzer parameter 487
Mass Precision parameter 392, 405
master proteins
displaying modification positions for 362
specifying default filter set for Proteins page of .pdResult file 269
Max. MGF File Size (MB) parameter 28
Max. Number of Consensus Workflows in Parallel Execution parameter 32
Max. Number of Processing Workflows in Parallel Execution parameter 32
Max. Precursor Mass parameter 392, 405
Maximum Allowed Fold Change parameter 446, 480

- Maximum Delta Cn parameter
 - calculating maximum Delta Cn [246](#)
 - MSF Files node [245, 249](#)
 - Percolator node [245, 247](#)
 - Target Decoy PSM Validator node [245](#)
- Maximum Delta Cn value [246](#)
- Maximum Delta Mass parameter [251](#)
- Medium isotopic peak patterns [467](#)
- membrane GO Slim category [342](#)
- Merge Mode parameter
 - determining how the application merges data in multiconsensus reports [137](#)
 - Do Not Merge setting [141](#)
 - Globally by Search Engine Type setting [139](#)
 - Per File and Search Engine Type setting [140](#)
- metabolic process GO Slim category [344](#)
- metal ion binding GO Slim category [339](#)
- methyl modifications [358](#)
- methylation [358](#)
- MGF files
 - contents of [28](#)
 - input file type [11, 36](#)
 - output file type [12](#)
 - splitting [28](#)
- Minimum PSM Confidence parameter [363](#)
- mirror plots [225–226](#)
- missing reporter ions
 - accepting spectra with [480](#)
 - rejecting quantification results [482](#)
 - replacing [484](#)
 - treatment of reporter ion quantification
 - missing reporter ions [479](#)
- mitochondrion GO Slim category [342](#)
- Modification in Master Proteins column [362](#)
- Modification in Proteins column [362](#)
- Modification parameter [363](#)
- Modification Sites Reported parameter [363](#)
- Modification Sites Shown parameter [363–364](#)
- Modifications column
 - Peptide Groups page [364](#)
 - Proteins page [357, 363](#)
- molecular function codes [333](#)
- Molecular Function column [291, 296](#)
- molecular functions [333, 339](#)
- Molecular Functions page [290, 327, 333](#)
- motifs [359](#)
- motor activity GO Slim category [339](#)
- MS Order parameter [487](#)
- MS/MS spectra
 - fragmenting reporter ions in [460](#)
 - generating CID [7, 497](#)
 - types of fragment ions observed in [7](#)
 - validating searches with FDRs [277](#)
- MS/MS Spectrum Info page
 - finding those unique to search or search node [272](#)
 - PSM Ambiguity column [285](#)
 - % Isolation Interference column [436](#)
- MSF files
 - adding to existing study [100](#)
 - adding to studies [53](#)
 - adding to study [100](#)
 - converting to .pdResult files [99](#)
 - reprocessing by using new consensus workflow [102](#)
 - submitting to Workflow Editor [4](#)
- MSF Files node
 - Avoid Expressions parameter [214](#)
 - creating a consensus workflow [113](#)
 - determining how the Proteome Discoverer application parses FASTA title lines [213](#)
 - determining protein title line parsing [212](#)
 - filtering with
 - Maximum Delta Cn parameter [245, 249](#)
 - Maximum Delta Mass parameter [251](#)
 - Score and Threshold parameters [251](#)
 - Merge Mode parameter [137, 139](#)
 - precursor ion area detection [406](#)
 - precursor ion quantification [394](#)
 - Preferred Accession parameter [214](#)
 - Preferred Taxonomy parameter [214](#)
 - PTM analysis workflow [354](#)
 - reporter ion quantification [401](#)
 - Title Line Rule parameter [213](#)
- MSPepSearch node
 - adding spectrum library for searching with [221](#)
 - description [220](#)
 - dot score [224](#)
 - MSPepSearch score [224](#)
 - reverse dot score [224](#)
 - scores reported [224](#)
 - spectrum library search processing workflow [224](#)
 - using as search node [224](#)
- MSPepSearch score [224](#)
- MudPIT
 - creating a search workflow [49](#)
 - creating parameter file to call Proteome Discoverer
 - Daemon from the Xcalibur data system [154](#)
 - creating workflow for processing [145](#)
 - description [8, 146](#)
 - fractions [49](#)
 - processing samples in Proteome Discoverer Daemon
 - utility [143, 146, 166](#)
 - scoring [30](#)
 - using processing method [163](#)
- multicell editing mode [59](#)

Index: N

multiconsensus reports
 applying filters specific to different search engines 276
 creating 135
 creating by using multiple search engines 136
 creating from searches using multiple input files 136
 finding proteins common to two or more searches 272
 finding proteins unique to each search 274
 merging results by search engine node 138
 Proteins page 135
 PSMs page 135
Multidimensional Protein Identification Technology. *See*
 MudPIT
multiple searches 143
MZDATA files
 input file type 11, 36
 output file type 12
MZML files
 input file type 11, 36
 output file type 12
MZXML files
 input file type 36
 output file type 11

N

National Center for Biotechnology Information (NCBI) 5,
 292, 329, 491
National Institute of Standards and Technology (NIST) 220
National Institutes of Health (NIH) 333, 335–336
New Analysis icon 66
New Sequence Template dialog box 160
New Study and Analysis dialog box 38
No Sorting icon 85
No-Cleavage reagent 240
No-Enzyme reagent 240
no-enzyme searches 197
Normalization Mode parameter 443
Normalization parameter 444
normalized abundances 459
normalizing peptide groups and proteins 443
N-Terminal Modification parameter 75
nucleotide binding GO Slim category 339
nucleus GO Slim category 342
Number Flanking Residues in Connection Tables parameter
 361
Number of Attempts to Submit the Search parameter 28
Number of Parallel Tasks parameter 23
Protein Unique Peptides column 442
Razor Peptides column 443
Number of Spectra Processed at Once parameter 23

Unique Peptides column
 Proteins page 287
numeric factor dialog box 48
numeric study factors 47

O

¹⁸O labeling quantification method 382, 490
Only for Master Proteins parameter 360
Open Analysis Template icon 66, 72
Open Common icon 103, 117, 122
Open Containing Folder icon 65
Open Existing Study icon 40
Open icon 127
Open QualBrowser command 236
Open Study icon 127
Options dialog box
 ProteinCenter page 295
Options icon 208–209
Orbitrap Fusion output template 119
Orbitrap mass spectrometers 14, 449
organelle lumen GO Slim category 343
outputs of Proteome Discoverer application 12
oxidation 357

P

Parallel Job Execution parameters 15, 32
parallel workflows 129
parameter file
 creating 153
 MudPIT processing 163
 purpose 153
Parameters pane 35
.pdAnalysis files 80
PDB database 491
.pdConsensusWF files 79, 115
.pdProcessingWF files 77, 111
.pdResult files
 adding to existing study 100
 adding to studies 53
 adding to study 100
 annotating proteins in 5
 deleting from Start Page 20
 deleting from study 96
 example 96
 filtering. *See* filtering
 loading stored filters 268
 opening analysis through 72
 opening associated workflows from 97, 127
 Quan Spectra page 474
 reprocessing existing analysis 101

- saving filter set 267
- searching analysis results for specified string 98
- uploading search results to ProteinCenter 325
- Peptide and Protein Filter node
 - basic consensus workflow 113
 - filtering PSMs 254
 - in marked contaminant searches 217
 - precursor ion area detection workflow 407
 - precursor ion quantification workflow 394
 - protein grouping 285
 - PTM analysis workflow 354
 - reporter ion quantification workflow 401
- Peptide and Protein Quantified node
 - Maximum Allowed Fold Change parameter 446
- Peptide and Protein Quantifier node
 - Average Reporter S/N Threshold parameter 433
 - calculating protein group ratios 447
 - calculating PSM ratios 432
 - calculating quantification ratios 432
 - Co-Isolation Threshold parameter 436
 - Consider Protein Groups for Peptide Uniqueness parameter 439, 448
 - description 432
 - Fold Change Threshold parameters 454
 - Peptides to Use parameter 439, 448
 - precursor ion area detection 113, 128, 406–407
 - precursor ion area detection workflow 407
 - precursor ion quantification 113, 128, 395
 - Reject Quan Results with Missing Channels parameter 433
 - Replace Missing Values with Minimum Value parameter 433, 484
 - reporter ion quantification 113, 128, 400–401
 - Top *N* Peptides Used for Area Calculation parameter 448
 - Use Single-Peak Quan Channels parameter 433
- Peptide and Protein Quantifier node
 - Maximum Allowed Fold Change parameter 480
 - Reject Quan Results with Missing Channels parameter 482
- peptide group ratios 475, 478
- peptide groups
 - abundances 442
 - calculating area in quantification 447
 - calculating FDRs for 281
 - calculating grouped abundances 446
 - classifying 442
 - filtering 245
 - normalizing 443
- Peptide Groups page
 - Annotated Sequence column 360
 - columns added by the PSM Grouper node 357
 - determining how modifications are displayed 365
 - determining positions of PTMs 359
 - displaying peptide modifications 364
 - Modifications column 364
 - Position in Proteins column 359
 - Positions in Master Proteins column 360
- Peptide in Protein Annotation node
 - adding Modification column to Proteins page 357
 - Annotate Flanking Residues of the Peptide parameter 360
 - Minimum PSM Confidence parameter 363
 - Modification parameter 363
 - Modification Sites Reported parameter 363
 - modifications display workflow 359
 - Number Flanking Residues in Connection Tables parameter 361
 - Protein Modifications Reported parameter 362
 - PTM analysis consensus workflow 351, 354
 - Report Only PTMs parameter 363
- Peptide Quan Usage column 442, 482
- Peptide Relevance Factor parameter 24
- Peptide Spectrum Match Identification Details view
 - spectrum library searches 225
- Peptide Validator node
 - basic consensus workflow 113
 - peak area calculation detection 406
 - precursor ion quantification workflow 394
 - PTM analysis workflow 354
 - reporter ion quantification workflow 401
- PeptideAtlas home page 220
- peptides
 - confidence indicators 278
 - ratios 436
- Peptides to Use parameter 439, 442, 448
- pepXML files 12
- % Isolation Interference column 436
- Percolator node
 - calculating FDRs for peptide groups 281
 - calculating FDRs for PSMs 281
 - creating processing workflow 105
 - filtering with Maximum Delta Cn parameter 245, 247, 393
 - peak area calculation detection workflow 404
 - precursor ion quantification workflow 391, 393
 - reporter ion quantification workflow 398
 - validating identifications from different search engines 279
- Pfam database
 - accession identifier 332
 - computational enrichment 332
 - description 5, 291, 296
 - displaying annotation results from ProteinCenter 320
 - features contrasted with GO database features 292

- Hidden Markov Model 292
- retrieving information from ProteinCenter 5, 290–291, 296, 299
- Pfam IDs column 292, 296, 320
- phosphopeptide analysis 247
- phosphorylation
 - calculating site probabilities 351
 - filtering site probabilities 374
 - isoforms 349
- pinning names on Start Page 19
- PIR database 491
- Please Select a FASTA Database dialog box 183
- Position column 232
- Position in Proteins column 359
- Positions in Master Proteins column 360
- Post-Processing Nodes pane 35, 115
- post-translational modifications (PTMs)
 - acetylation 357
 - calculating site probabilities with the ptmRS node 347, 355
 - creating PTM analysis workflow 129, 351
 - determined by dynamic modifications 227
 - displaying in Protein Identification Details view 370–371
 - displaying protein sequences 369
 - filtering phosphorylation site probabilities 374
 - glycosylation 359
 - lysine 358
 - methylation 358
 - oxidation 357
 - Peptide Groups page 357
 - Proteins page 357
 - PSMs page 355
- PQD activation type
 - analyzed by Sequest HT 3
 - description 7
 - specifying with
 - Reporter Ions Quantifier node 486
- precursor ion area detection
 - creating consensus workflow 406
 - creating processing workflow 403
 - creating workflow 128
 - description 4, 389
 - performing 403
 - Precursor Ions Area Detector node 404
- precursor ion quantification
 - calculating PSM ratios 475
 - calculating quantification ratios 477
 - changing a quantification method
 - after generating a results report 425
 - before generating a results report 424
 - checking the quantification method 424
 - chemical elements supported 383
 - creating a consensus workflow for 394
 - creating a processing workflow for 390
 - creating workflow 128
 - data distribution maps 453, 456
 - deactivating a quantification method 426
 - default methods available in 382
 - description 381
 - dimethylation 3plex method 385
 - displaying quantification channel values 459, 463
 - displaying Quantification Spectrum chart 466
 - distribution maps 454
 - exporting a quantification method 427
 - Full ¹⁸O labeling quantification method 385
 - identifying isotope patterns 451
 - importing a quantification method 426
 - Incomplete ¹⁸O labeling quantification method 385
 - intensity 432
 - performing 390
 - removing a quantification method 426
 - restoring quantification method defaults 427
 - setting up quantification method 409
 - SILAC 2plex methods 383
 - SILAC 3plex methods 385
 - specifying quantification channels 410–411
 - troubleshooting 490
- Precursor Ions Area Detector node
 - performing precursor ion area detection 4, 128, 403–404, 406
- Precursor Ions Quantifier node
 - performing precursor ion quantification 128, 390
 - precursor ion quantification workflow 391, 393–394
 - using with Precursor Ions Area Detector node 391
 - using with Reporter Ions Quantifier node 391
- PrediSi database 332
- Preferred Accession parameter 214
- Preferred Taxonomy parameter 214
- PRF database 491
- processing methods 155, 157, 163
- Processing Setup window 155
- Processing Step box 71
- processing workflow
 - creating basic 74, 104, 106
 - decoy database searches 279
 - description 34
 - factory default templates 119
 - marked contaminants searches 217
 - parallel 130
 - precursor ion area detection 403
 - precursor ion quantification 390
 - protein annotation 296
 - PTM analysis 352
 - reporter ion quantification 397
 - saving 77, 111

- specifying number executed in parallel 15, 32
- spectrum library 224
- Programs icon 155
- Programs view 155
- protein abundances 448
- protein annotation
 - configuring URLs for 293
 - creating consensus workflow 129, 297
 - creating processing workflow 129, 296
 - creating workflow for 296
 - Ensembl genome database 5, 290, 292
 - Entrez Gene database 5, 290, 292, 296
 - GO database 5, 290–291, 296
 - Pfam database 5, 290–291, 296
 - retrieving from ProteinCenter 296
 - UniProt database 5, 290, 293, 296, 373
- protein binding GO Slim category 339
- protein complex GO Slim category 342
- Protein Database parameter
 - Mascot node 134
 - processing workflows 75
 - Sequest HT node 108
 - setting in precursor ion area detection 405
 - setting in precursor ion quantification 393
 - setting in reporter ion quantification 398
- Protein Details command/icon 326
- Protein FDR Validator node
 - calculating FDRs for proteins 281
- Protein Grouping node
 - Apply Strict Parsimony Principle parameter 284
 - in marked contaminant searches 217
 - peak area calculation detection workflow 407
 - precursor ion quantification workflow 395
 - PTM analysis workflow 354
 - reporter ion quantification workflow 401
 - workflow used to group proteins 286
- protein groups
 - calculating FDR 282
 - filtering 245
- Protein Identification Details dialog box 327
- Protein Identification Details view
 - controlling probability of site PTM occurrence 371
 - copying colored bar on Sequence page 374
 - displaying
 - filtered-out PSMs 371
 - only PTMs 371
 - ProteinCard 290
 - PTMs 369–371
 - UniProt annotations of protein 373
 - displaying filtered-out PSMs 371
 - exporting Modification List page to Excel 374
 - finding protein PTMs 373
 - modification site probabilities calculated by ptmRS
 - node 358
 - ProteinCard page 290, 326
 - PTM site localization scores 352
 - saving colored bar on Sequence page 374
- Protein Marker node
 - creating consensus workflow with 217
 - in marked contaminant searches 217
- Protein Modifications Reported parameter 362
- Protein Positions for Peptides parameter 359
- Protein Quan Usage column 443
- protein quantification 440
- protein ratios
 - calculating 475
 - calculating from peptide ratios 436
- protein references 181–182, 185, 191
- Protein Scorer node
 - basic consensus workflow 113
 - marked contaminant consensus workflow 217
 - precursor ion area detection workflow 407
 - precursor ion quantification workflow 395
 - PTM analysis workflow 354
 - reporter ion quantification workflow 401
- Protein Scoring Options parameter 24
- protein sequences 181–182, 185, 191
- ProteinCard page
 - Biological Processes page 290, 327, 335
 - Cellular Components page 290, 327, 334
 - Diseases page 327, 337
 - External Links page 290, 327, 337
 - Features page 290, 327, 331
 - General page 290, 327–328
 - information displayed in 290
 - Keys page 290, 327, 330
 - Molecular Functions page 290, 327, 333
 - opening 290, 326
 - pages of 290, 328
- ProteinCenter
 - accessing annotation data in 299
 - accessing ProteinCards for each protein 326
 - configuring the Proteome Discoverer application for
 - uploading results 293, 295
 - configuring URLs for access to 293
 - description 5, 290
 - downloading FASTA files from 174
 - retrieving annotations with ProteinCenter Annotation
 - node 296
 - retrieving information from
 - Ensembl genome 296
 - Ensembl genome database 5, 290, 292
 - Entrez Gene database 5, 290, 292, 296
 - GO database 5, 290–291, 296

Index: P

- Pfam database 5, 290–291, 296
- UniProt database 5, 290, 293, 296, 373
- updating FASTA files from 178
- uploading search results to 323
- ProteinCenter Annotation node
 - retrieving protein annotations from ProteinCenter 129, 290, 296, 373
 - selecting user-defined annotation aspects 312
- ProteinCenter page 295
- proteins
 - accession keys 330
 - annotation. *See* protein annotation
 - biological processes 335
 - calculating abundance 448
 - calculating area in quantification 447
 - calculating FDR 281
 - calculating grouped abundances 446
 - cellular components 334
 - classifying abundances 443
 - diseases associated with 337
 - displaying
 - general information about 328
 - filtering 245, 254
 - finding PTMs 373
 - finding those common to two or more searches 272
 - finding those unique to each search 272
 - finding those unique to each search node 274
 - function of 333
 - group ratios 447
 - grouping
 - algorithm used 283
 - Peptide and Protein Filter node 285
 - Protein Grouping node 284, 286
 - PSMs identified by multiple workflow nodes 287
 - marked as contaminants in FASTA files 217
 - normalizing 443
 - retrieving information from ProteinCenter 290
 - selecting PSMs for protein group 285
 - sequence features 331
 - web links to information 337
- Proteins page
 - # Unique Peptides column 287
 - accessions on 209
 - Chromosome column 322–323
 - Contaminant column 218
 - descriptions on 209
 - determining positions of PTMs 359
 - displaying
 - accessions and descriptions 213
 - annotation aspects 300
 - GO database information from ProteinCenter 290
 - Pfam annotations 292
 - Ensembl Gene ID column 322–323
 - Entrez Gene ID column 321
 - finding those unique to search or search node 272
 - GO Accessions column 318
 - Modifications column 357, 363
 - multiconsensus reports 135
 - Pfam IDs column 320
 - Species column 219
 - Species Map column 219
 - specifying default filter set for 269
- Proteome Discoverer application
 - activation types supported 6
 - closing 18
 - components 2
 - configuring for protein annotation 293
 - configuring URL for uploading results to ProteinCenter 293, 295
 - data export formats 6
 - FASTA databases used 6, 11
 - features 1–2
 - filtering data 6
 - graphical views in 5
 - input files 11
 - Magellan server 2
 - MudPIT experiments 8
 - new features in this release 12
 - outputs 12
 - ProteinCenter access 5
 - Proteome Discoverer client 2
 - Proteome Discoverer Daemon utility 2
 - Qual Browser application 6
 - quantification 4
 - search engines available 3
 - Start Page 17
 - starting 17
 - steps in using 8
 - study management system 3
 - types of mass spectrometer data supported 11
 - Workflow Editor. *See* Workflow Editor
- Proteome Discoverer application icon 144
- Proteome Discoverer client 2
- Proteome Discoverer Daemon utility
 - batch processing 145
 - batch processing with
 - multiple processing methods 160
 - single processing method 157
 - command-line syntax 166
 - command-line syntax examples 168
 - connecting to Magellan server 144
 - connecting to remote server 149
 - creating parameter file for calling utility from Xcalibur data system 153
 - creating processing method 155
 - description 2, 143

- Job Queue page 148
 - monitoring job execution 148
 - MudPIT processing 145
 - processing samples by using a processing method 163
 - running
 - from command line 165
 - from window 145
 - from Xcalibur data system 150
 - on local server 167, 170
 - on remote server 167, 169
 - running on remote server 167
 - specifying sample types to be sent to 157
 - specifying server to connect to 144
 - Start Jobs page 147
 - starting from Xcalibur data system 151
 - starting in window 143
 - Proteome Discoverer Daemon utility window
 - Configuration page 144
 - Start Jobs page 145
 - ProtXML files 12
 - PSM Ambiguity column 285, 440
 - PSM Grouper node
 - adding Modifications column to Peptide Groups page 357
 - basic consensus workflow 113
 - grouping PSMs 288
 - Modification columns added to results report 363–364
 - Modification Sites Shown parameter 363–364
 - peak area calculation detection workflow 406
 - precursor ion quantification workflow 394
 - reporter ion quantification workflow 401
 - Site Probability Threshold parameter 363, 367
 - PSMs
 - calculating abundances in quantification 432
 - calculating FDRs 281
 - calculating peptide group abundances in quantification 442
 - calculating ratios 432, 475, 478
 - classifying 436, 439
 - displaying filtered-out 371
 - excluding from quantification 450
 - excluding those with high levels of coisolation 435
 - filtering
 - with the Fixed Value PSM Validator node 248
 - with the Maximum Delta Cn parameter 246
 - with the MSF Files node 245, 249, 251
 - with the Peptide and Protein Filter node 254
 - with the Percolator node 247
 - with the Target Decoy PSM Validator node 248
 - grouping 288
 - PSMs page
 - columns added by ptmRS node 355
 - Contaminant column 218
 - finding those unique to search or search node 272
 - multiconsensus reports 135
 - PSM Ambiguity column 285
 - ptmRS Best Site Probabilities column 355, 357
 - ptmRS Binomial Peptide Score column 356–357
 - ptmRS Isoform Confidence Probability column 356–357
 - ptmRS Isoform Group Report column 356–357
 - ptmRS *Modification* Site Probabilities column 355, 357
 - Search ID column 135
 - XCorr column 135
 - % Isolation Interference column 436
 - PTM analysis
 - creating consensus workflow 129, 351, 354
 - creating processing workflow 129, 351–352
 - results of 355
 - ptmRS Best Site Probabilities column 355, 357
 - ptmRS Binomial Peptide Score column 356–357
 - ptmRS Isoform Confidence Probability column 356–357
 - ptmRS Isoform Group Report column 356–357
 - ptmRS *Modification* Site Probabilities column 355, 357
 - ptmRS node
 - adding modification site probabilities to Protein Identification Details view 358
 - calculating binomial peptide score 350
 - calculating ptmRS sequence probabilities 351
 - calculating sequence probabilities 348, 350–351
 - calculating site probabilities 348, 351
 - creating PTM analysis processing workflow 351
 - description 347
 - filtering phosphorylation site probabilities 374
 - input to 347
 - internal workflow for phosphorylation 348
 - methylation 358
 - processing workflow 352
 - ptmRS Best Site Probabilities column 355, 357
 - ptmRS Binomial Peptide Score column 356–357
 - ptmRS Isoform Confidence Probability column 356–357
 - ptmRS Isoform Group Report column 356–357
 - ptmRS *Modification* Probabilities column 355, 357
 - See also* PTM analysis 351
 - PubMed database 333
 - purity correction factors 434
 - putative isoforms 348, 351
- ## Q
- Qual Browser application 6
 - Quan Channel Values command/icon 460
 - Quan Info column 482
 - Peptide Groups page 442
 - PSMs page 440

- Quan Ratio Distributions chart
 - displaying 475
 - features 478
- Quan Ratio Distributions command/icon 475
- Quan Spectra page
 - description 474
 - displaying 475
 - parameters 475
- Quan Usage column 466
- quantification
 - calculating areas 447–448
 - calculating grouped abundances and ratios 446
 - calculating peptide group abundances 442
 - calculating protein abundances 443
 - calculating PSM abundances 432
 - calculating PSM ratios 432, 475
 - calculating quantification ratios 432
 - calculating results 432
 - classifying peptide groups 442
 - classifying PSMs 436
 - creating workflow 128
 - custom ratios
 - description 82
 - generating according to channels 90
 - displaying areas and ratios as data distribution maps 453, 456
 - displaying channel values 460
 - displaying spectrum chart 466
 - excluding PSMs with high levels of coisolation 435
 - handling missing reporter ions 479
 - ITRAQ 388
 - normalizing peptide groups and protein abundances 443
 - precursor ion 381
 - precursor ion area detection 389
 - reporter ion 385
 - setting up the quantification method 408
 - specifying the quantification channels 410
 - steps performed in 431
 - TMT 386
 - validation checks 93
- Quantification Channel Values chart
 - description 459
 - displaying 460
 - displaying grouped abundances 460, 462, 465
 - displaying legend 462, 465
 - displaying ungrouped abundances 460, 462, 465
 - grouped abundances 459
 - normalized abundances 459
 - precursor ion quantification
 - description 459
 - displaying for
 - peptide groups 459, 463–464
 - proteins 459, 463–464
 - PSMs 459, 463
 - quantification spectra 459, 463
 - reporter ion quantification
 - description 459
 - displaying for
 - peptide groups 459–460, 462
 - proteins 459–460, 462
 - PSMs 459–460
 - quantification spectra 459–460
 - shortcut menu 460
 - ungrouped abundances 459
- quantification channels
 - definition 432
 - definition for precursor ion quantification 409
 - definition for reporter ion quantification 409
 - displaying values 460
 - precursor ion quantification 411
 - quantification values 432
 - reporter ion quantification 417
- quantification method
 - changing
 - after generating a result report 425
 - before generating a result report 424
 - checking the parameters set 424
 - configuring quantification method on Mascot server 428
 - copying studies to another computer 64
 - creating to correct TMT 10plex kits 420
 - deactivating 426
 - exporting 427
 - importing 426
 - removing 426
 - restoring original template 427
 - selecting 53
 - selecting in study 40
 - setting up
 - in study 409
 - precursor ion quantification 409
 - reporter ion quantification 409
 - specifying quantification channels 410
 - validating 93
 - ways of creating 413, 417
- Quantification Method Editor dialog box
 - changing a method after generating a results report 425
 - changing a method before generating a results report 424
 - changing a quantification method before generating a results report 425

- precursor ion quantification 413
- reporter ion quantification 417, 421
- Quantification Methods view
 - adding quantification method 41, 411, 417, 420
 - changing quantification method 424
 - deactivating quantification method 426
 - exporting quantification method 427
 - importing quantification method 426
 - removing quantification method 426
 - restoring quantification method template defaults 427
- Quantification Ratio Distributions chart
 - purpose 477
- quantification ratios 432
 - custom
 - generating manually 88
 - generating semiautomatically 84
- Quantification Spectrum chart
 - displaying
 - precursor ion quantification 466
 - reporter ion quantification 471
 - opening 466
- quantification value, definition 432
- QuickGO browser 333
- quality algorithm 281

R

- raw data files
 - contents of Xcalibur data system 11
 - in parallel workflows 129–130
 - input file type 36
 - performing multiple searches on multiple 143
 - processing multiple from multiple samples 143
 - processing multiple from one sample 49, 143
 - submitting multiple files to Workflow Editor 4
 - when Qual Browser does not open 236
- razor peptides 442–443
- RDB equivalents 236
- receptor activity GO Slim category 339
- Refresh icon 204
- regular expressions
 - adding 211
 - applying changes 212
 - changing or renaming 211
 - description 209
- regulation of biological process GO Slim category 345
- Reject Quan Results with Missing Channels parameter 433, 482
- Remove FASTA Databases dialog box 180
- Remove FASTA indexes confirmation box 206
- Remove Files icon 65, 96
- Remove icon
 - Chemical Modifications view 233
 - FASTA Files view 180
 - FASTA Indexes view 206
 - Quantification Methods view 426
- Remove Spectrum Libraries Databases dialog box 224
- Replace Missing Values with Minimum Value parameter 433, 484
- replicates
 - biological 499
 - specifying 499
 - technical 499
- Report Only PTMs parameter 363
- Reporter Abundance Based On parameter 449
- reporter ion quantification
 - calculating PSM ratios 475, 479
 - calculating quantification ratios 475
 - changing a quantification method
 - after generating a results report 425
 - before generating a results report 424
 - checking the quantification method 424
 - coisolation 436
 - creating a consensus workflow 400
 - creating consensus workflow 128
 - creating processing workflow 128, 397
 - data distribution maps 453, 456
 - deactivating a quantification method 426
 - description 4, 385
 - displaying quantification channel values 459–460
 - displaying Quantification Spectrum chart 471
 - distribution maps 454
 - exporting a quantification method 427
 - filtering data
 - with average signal-to-noise values 450
 - with signal-to-noise values instead of intensities 449
 - filtering Orbitrap data 449
 - importing a quantification method 426
 - intensity 432
 - isotopic distribution values 434
 - iTRAQ quantification 388
 - Orbitrap data 14
 - performing 396
 - performing on HCD and CID scans 402
 - Quan Spectra page 474
 - removing a quantification method 426
 - restoring quantification method default 427
 - setting up quantification method 409
 - specifying quantification channels 417
 - TMT 386
 - troubleshooting 486
 - iTRAQ quantification *See* iTRAQ quantification
 - TMT quantification *See* TMT quantification

Reporter Ions Quantifier node
 Activation Type parameter 487
 Integration Method parameter 473–474
 Integration Tolerance parameter 471, 486
 Mass Analyzer parameter 487
 MS Order parameter 487
 nodes not used with 398
 performing reporter ion quantification 128
 performing TMT quantification on HCD and CID scans 402
 reporter ion quantification 398–399
 Reprocess icon 72, 101, 425
 reprocessing existing analysis 101
 reprocessing MSF results 102
 reproduction 345
 Restore FASTA indexes confirmation box 207
 Restore icon 207
 Result Summaries view
 display filter settings 260
 results report. *See* .pdResult files
 Reuse icon 101
 reverse dot score 224
 R-friendly headers 16
 ribosome GO Slim category 342
 RNA 328–329, 339, 342
 RNA binding GO Slim category 339
 Run Sequence dialog box 158
 Run Sequence icon 158

S

S/N Threshold parameter 392, 405
 sample groups
 calculating areas for 447
 calculating protein areas for 448
 creating in Discoverer Daemon 146
 displaying in Quantification Channel Values chart 459
 displaying on Grouping & Quantification page 70
 specifying quantification ratios from selected 82
 sample types 68, 157
 samples
 grouping 84
 grouping in data distribution maps 455
 MudPIT 143
 setting factor values for 54–55
 setting values for multiple samples at the same time on
 Input Files page 58
 Samples page 61
 specifying
 blank sample types 54
 control sample types 54, 68
 sample types 68

specifying standard sample types 54, 68
 viewing 55
 Samples page 67
 Save a Parameter File dialog box 154
 Save Common icon 126
 Save icon 126
 Save Modification List as Microsoft Excel File dialog box 374
 Save Quantification Method dialog box 416, 418, 423
 Scan Event Filter node
 parallel workflow 131
 used for Mascot or Sequest analysis 111
 Search ID column 135
 searches
 performing in batch mode 94
 performing in individual mode 92
 Seattle Proteome Center 11–12
 Select a Spectrum Library dialog box 222
 Select Processing Method dialog box 158
 Semi (C-Term) cleavage enzyme specificity 239, 243
 Semi (N-Term) cleavage enzyme specificity 239, 243
 Semi cleavage enzyme specificity 239–240, 243
 Semi-C cleavage enzyme specificity 240
 semi-enzymatic searches 197, 201
 Semi-N cleavage enzyme specificity 240
 Sequence Setup icon 158
 Sequence View icon 158
 Sequest HT node
 parallel processing workflow 131
 PTM analysis workflow 352
See also Sequest HT search engine
 Sequest HT search engine
 adding FASTA files before searching 174
 applying filters in multiconsensus reports 276
 configuration parameters for 21–22
 creating a processing workflow 105
 creating FASTA index 201
 data types analyzed 3
 merging results 138
 precursor ion quantification workflow 393
 protein score configuration parameters 24
 reporter ion quantification workflow 398
 selecting FASTA files 108
 SpScore value 3
 XCorr values
 function 3
 in multiconsensus reports 135
 setting XCorr confidence thresholds 23–24
 Settings Applied dialog box 32
 SGD database 209
 Show Details icon 97
 Show Legend command 460

- Show Protein/Peptide ID Details command/icon 369
- Show Quantification Spectrum command/icon 466
- Show Workflow icon 74, 78, 103
- Show/Hide Advanced Parameters command 118
- signal transducer activity GO Slim category 339
- Significance Middle parameter 29
- SILAC 2plex (Arg10, Lys6) quantification method 382–383, 490
- SILAC 2plex (Arg10, Lys8) quantification method 382–383, 490
- SILAC 2plex (Ile6) quantification method 382–383, 490
- SILAC 3plex (Arg6, Lys4|Arg10, Lys8) quantification method 382, 385, 490
- SILAC 3plex (Arg6, Lys6|Arg10, Lys8) quantification method 382, 385, 490
- SILAC. *See* precursor ion quantification
- site probabilities 374
- Site Probability Threshold parameter 363, 367
- software requirements xiii
- Sort Ascending icon 85
- Species column 219
- Species Map column 219
- species names 219
- spectra exportation
 - creating a consensus workflow 133
 - creating a processing workflow 133
- Spectral Libraries command 220
- Spectrum Files node
 - basic processing workflow 104
 - parallel workflow 130
 - peak area calculation detection workflow 404
 - precursor ion quantification workflow 391
 - PTM analysis workflow 352
 - reporter ion quantification workflow 397
- spectrum libraries
 - generating mirror plots 225
 - searching with MSPepSearch node 220
- Spectrum Libraries view
 - displaying 220
- spectrum library searches
 - adding spectrum library 221
 - creating consensus workflow 225
 - creating processing workflow 224
 - deleting spectrum libraries 224
 - description 220
 - displaying downloaded libraries 220
 - generating mirror plots 226
 - listing available libraries 220
 - visually verifying spectrum library matches 225
- Spectrum Selector node
 - creating a processing workflow 104
 - parallel workflow 130
 - peak area calculation detection 404–405
 - precursor ion quantification workflow 391–392
 - PTM analysis workflow 352
 - reporter ion quantification workflow 398
 - selecting precursor mass to use 105
- spliceosome GO Slim category 342
- SpScore value 3
- standard errors 447
- standard sample types 54, 68
- Start Jobs page 145, 147
- Start Page
 - accessing 17, 19
 - clearing all names from 20
 - closing 21
 - creating a study 38
 - deleting studies or result files from 20
 - description 18
 - opening existing studies 40
 - retaining study and result file names on 19
- Start Page command/icon 19
- Static Modification parameter
 - Sequest HT 75, 393, 399, 405
- static modifications
 - definition 227
 - setting for FASTA indexes 203
- structural molecule activity GO Slim category 340
- studies
 - adding
 - .pdResult files to 53
 - fractions to 50
 - input files to 48
 - MSF files to 53
 - MSF or .pdResult files to existing 100
 - MSF or .pdResult files to new 100
 - multiple unrelated input files to 49
 - single input file to 49
 - study factors to 44
 - copying to another computer 63
 - creating 38
 - creating subsets of fractions 50
 - deleting .pdResult files 96
 - deleting from Start Page 20
 - description 3, 33
 - example 36
 - Grouping and Quantification page 82
 - opening existing 40
 - opening workflows and analyses from results 97
 - retaining names on Start Page 19
 - saving 62, 80
 - searching analysis results for specified string 98
 - selecting quantification method 40
 - setting factor values for samples 54–55
 - setting up the quantification method 409

Index: T

- setting values for multiple samples at the same time on
 - Input Files page 58
 - Samples page 61
- specifying quantification method for input files 53
- Study window parameters 65
- Study Definition page 66
- study management system 3
- Study Variables area 83
- SwissProt database 491–492
- system requirements xii

T

- tab-delimited TXT files
 - export format 12
- TAIR Web pages 209
- tandem mass tag quantification. *See* TMT quantification
- Target Decoy PSM Validator node
 - creating basic processing workflow 105
 - filtering PSMs 248
 - filtering with Maximum Delta Cn parameter 245, 393
 - peak area calculation detection workflow 404
 - precursor ion quantification workflow 391, 393
 - PTM analysis workflow 352
 - reporter ion quantification workflow 398
 - validating identifications from different search engines 279
- target FDRs 278
- Target Name parameter 324
- taxonomy identification number 175
- Taxonomy Rules area 210–211
- technical replicates 499
- third-party nodes 105, 113
- Time Interval Between Attempts to Submit a Search (sec) parameter 28
- Time Interval Between Attempts to Submit the Annotation Request (sec) parameter 324
- Title Line Rule parameter 213
- Title Line Rules area 210–211
- title lines 210, 213
- Tmap database 332
- TMT 10plex quantification method
 - as default 4, 387
 - correcting for isotopic impurities 419
 - correction factors 435
 - reporter ion tag masses 387
- TMT 2plex quantification method
 - as default 4, 386
 - reporter ion tag masses 387
- TMT 6plex quantification method
 - as default 4
 - reporter ion tag masses 387

- TMT quantification
 - default methods available in 386
 - description 4, 386
 - isotopic distribution values 434
 - performing on HCD and CID scans 402
 - reporter ion tag masses 387
 - Reporter Ions Quantifier node 4
 - specifying label modifications 410
 - See also* reporter ion quantification
- TMTe 6plex quantification method
 - as default 387
 - correction factors 435
 - reporter ion tag masses 387
- To ProteinCenter node 324
- Top *N* Peptides Used for Area Calculation parameter 448
- transcription regulator activity GO Slim category 340
- translation regulator activity GO Slim category 340
- transport GO Slim category 345
- transporter activity GO Slim category 340
- TrEMBL database 492
- trimethyl modifications 358
- TXT files. *See* tab-delimited TXT files

U

- U.S. National Library of Medicine (NLM) 333, 335–336
- ungrouped abundances 459, 462, 465
- UNIMOD
 - importing chemical modifications 234
 - updating chemical modifications 227
- UniProt database
 - description 5, 293, 296
 - displaying annotations 323
 - retrieving information from ProteinCenter 5, 290, 293, 296, 299, 373
 - specifying in ProteinCard 332
 - taxonomy identification numbers 175
- UniProt NREF database 492
- UniRef database 492
- UniRef50 database 492
- UniRef90 database 492
- Universal Protein Resource consortium 492
- Update icon 179
- uploading search results to ProteinCenter 323
- Use MudPIT Scoring parameter 30
- Use Single-Peak Quan Channels parameter 433

V

- vacuole GO Slim category 343
- validation errors 92

W

warnings 92

Wellcome Trust Sanger Institute (WTSI) 5, 291–292

Workflow Editor

- aligning nodes 106, 114
- closing analyses 104
- common workflow templates 119
- consensus workflow. *See* consensus workflow
- correcting workflow template errors 123
- creating a quantification workflow 128
- custom workflow templates 123
- deleting workflow templates 127
- description 4, 34
- features 35
- incorporating stand-alone workflow into an analysis 117
- job queue. *See* job queue
- joining nodes 105, 114
- Load Failures pane 124
- opening analyses 103
- opening common workflow templates 122
- opening custom workflow templates 123
- opening workflow associated with .pdResult file 127
- organizing nodes 105, 114
- outdated annotation node symbol 124
- outdated node version symbol 123
- outdated parameter setting symbol 124
- parameter filters 111
- parameters 118
- Parameters pane 107
- Post-Processing Nodes pane 115
- processing workflow. *See* processing workflow
- saving consensus workflow as template 115
- saving processing workflow as template 111
- saving workflow as common template 126
- saving workflow as custom template 126
- setting node parameters 107, 114
- stand-alone
 - closing 104
 - opening 103
 - window 34
- steps to take before using 36
- using third-party nodes 105, 113
- Workflow Nodes pane 104
- workflow templates. *See* workflow templates

Workflow Editor icon 104, 112

Workflow Nodes pane 35, 113

Workflow page 69

workflow templates

- common 119
- correcting errors 123
- custom 123
- definition 119

- deleting 127
- opening common 122
- opening custom 123
- saving common 126
- saving custom 126

Workflow Tree pane 35

Workflows page 72

X

x fragment ions

- definition 7

Xcalibur data system

- creating a parameter file 153
- Home Page window 155
- Qual Browser 236
- running injections sequence that starts the Proteome Discoverer Daemon utility 158
- running Proteome Discoverer Daemon utility from 150
- starting Discoverer Daemon in 2.1.0 158

Xcalibur Sequence Setup dialog box 154

XCorr column

- PSMs page 135

XCorr Confidence Thresholds parameter 23–24

XCorr values

- function 3
- in multiconsensus reports 135

Y

y* fragment ions 497–498

y fragment ions

- abstracting proton from precursor 7
- definition 7, 498
- generated by CID 3

Z

z fragment ions

- definition 7, 498
- generated by ETD and ECD 3

